
Clear-Retriever



(1) Soft Mesh Net

• Less possibility of damaging tissues.

(2) Wide Net Loop

• Enables to hold piecemeal polyps.

Retrieval of polyp and multiple polyps







SPECIFICATION

Model	Loop Size	Type of Retation	Loop Shape	Sheath 0.D	Working Length	Number of Use	Validity of Sterilzation
FM-ER0008	30mm	Rotational	Hexagonal	¢2.5mm	2400mm	Single Use Only	3 years

ClearCap Distal Attachment



FEATURE & BENEFITS

(1) Maintaining the Clear View

• This ClearCap Distal Attachment helps maintain the clear view by keeping a moderate distance from the tissues.

(2) Protecting the tissues

• By adopting the soft raw material, it helps users protect the

(3) Two draining holes

• These holes make the discharge of foreign materials and saline irrigation more easily.





SPECIFICATION

Model	Appilcable Lesion	0.D	I.D	Quantity of Hole for Draining	Number of Use	Validity of Sterilzation
FM-EC0001	Lower GI	¢15mm	¢12.6mm			
FM-EC0002	Upper GI	¢11.8mm	¢9.4mm			
FM-EC0003	Upper GI	¢11.35mm	¢8.9mm			
FM-EC0004	Upper GI	¢12.4mm	¢10mm	2	Single Use Only	3years
FM-EC0005	Lower GI	¢13.4mm	¢11mm		,	
FM-EC0006	Lower GI	¢14mm	¢11.6mm			
FM-EC0007	Lower GI	¢15.7mm	¢13.3mm			

Imagine What a Woderful difference You have 9