

# WANCARE

## ULTRASOUND ECG gel

### USAGE

Ensures clearer and uninterrupted transmission of ultrasonographic waves by disconnecting the air between the skin and the prob during any and all Ultrasonography, Doppler ECG or effort test.

### PROPERTIES

pH balanced between 5,00 - 7,00.

No toxic effect. Not irritate the skin.

Water soluble. Easy to clean.

Not contain formaldehyde and salt. Odorless.



### COMPOSITION

Deionized water, Carbormer, TriEthanol Amine, MonoPropylene Glycol, Benzyl Alcohol, Methylchloroisothiazolinone

### STORAGE CONDITIONS

Please store under 35°C keep away from direct sunlight, store in a cool and dry place.



SIZE	REF NO
250ML	KAF G31-2
500ML	KAF G31-1
1LT	KAF G31
5LT	KAF G31-3

