



Omnilink Elite

Peripheral Stent System

Setting a new standard
in iliac treatment.

6F Sheath Compatible
Across All Sizes

Excellent Deliverability

Flexibility Without
Compromising Strength



Endovascular
Commitment



6F
sheath
compatible

.035

 **Abbott**
A Promise for Life


Omnilink Elite

Peripheral Stent System


.035

Setting a new standard in iliac treatment.

6F Sheath Compatible Across All Sizes




- Less need to change introducer sheaths for faster procedures and reduced cost^{1,2}
- Decreased risk of vascular complications¹
- Shorter compression time^{1,3}



Excellent Deliverability





- Low-profile tip
- Dual-lumen catheter provides strong push
- Dual-layer balloon with excellent stent retention



Flexibility Without Compromising Strength



- Thin CoCr stent struts
- Proven flexible stent design
- High RBP
- No compromise in strength or radiopacity



Omnilink Elite Ordering Information

Stent Diameter (mm)	Stent Length (mm)	STOCK NUMBER Catheter Length		Max Post Dilatation Diameter (mm)	Crimped Stent OD (mm)	Min Sheath Size (ID)	Nominal Pressure/RBP (bar)	
		80 (cm)	135 (cm)					
4.0	12	11000-12	11007-12	6.0	1.85	6F (2.12 mm)	11/16	
	16	11000-16	11007-16					
	19	11000-19	11007-19					
5.0	12	11001-12	11008-12		1.90			11/14
	16	11001-16	11008-16					
	19	11001-19	11008-19					
	29	11001-29	11008-29					
	39	11001-39	11008-39					
	59	11001-59	11008-59					
6.0	12	11002-12	11009-12	8.0	1.95			
	16	11002-16	11009-16					
	19	11002-19	11009-19					
	29	11002-29	11009-29					
	39	11002-39	11009-39					
7.0	59	11002-59	11009-59		2.05			
	12	11003-12	11010-12					
	16	11003-16	11010-16					
	19	11003-19	11010-19					
	29	11003-29	11010-29					
	39	11003-39	11010-39					
8.0	59	11003-59	11010-59	2.10				
	19	11004-19	11011-19					
	29	11004-29	11011-29					
	39	11004-39	11011-39					
	59	11004-59	11011-59					
9.0	19	11005-19	11012-19	11.0	2.15	6F (2.15 mm)		
	29	11005-29	11012-29					
	39	11005-39	11012-39					
	59	11005-59	11012-59					
10.0	19	11006-19	11013-19		2.20	6F (2.20 mm)		
	29	11006-29	11013-29					
	39	11006-39	11013-39					

Maximum guide wire: .035 in (0.89 mm) for all sizes

Omnilink Elite Compliance

Stent Diameter (mm)	Inflation Pressure (bar)					
	11	12	13	14	15	16
4.0	4.00	4.08	4.16	4.24	4.32	4.40
5.0	5.00	5.09	5.19	5.28	5.37	5.47
6.0	6.00	6.08	6.16	6.25	-	-
7.0	7.00	7.12	7.25	7.37	-	-
8.0	8.00	8.14	8.27	8.41	-	-
9.0	9.00	9.13	9.26	9.39	-	-
10.0	10.00	10.15	10.31	10.46	-	-

1. Metz D, Meyer F, Touati C *et al.* *Am Heart J* 1997;134(1):131-137.
2. Talley JD, Mauldin PD, Becker ER. *J Interv Cardiol* 1995;8(4):345-353.
3. Büchler JR, Ribeiro EE, Falcão JL *et al.* *J Interv Cardiol* 2008;21(1):50-55

PALMAZ GENESIS Peripheral Stent is a trademark of Cordis Corporation, a Johnson & Johnson company.
Express Vascular LD Peripheral Pre-Mounted Balloon Expandable Stent is a trademark of Boston Scientific Corporation.
Scuba Co-Cr Peripheral Stent System is a trademark of INVATEC.
Chromaxx Balloon Expandable Peripheral Stent is a trademark of C.R. Bard, Inc.

Please contact your local representative for more information.

Abbott Vascular International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium

All illustrations included are artist's renditions.

All tests performed by and data on file at Abbott Vascular. All photographs taken by and on file at Abbott Vascular.

Multi-Link and Omnilink Elite are trademarks of the Abbott Group of Companies.

This product is intended for use by or under the direction of a physician. Prior to use, it is important to read the package insert thoroughly for Instructions for use, warnings and potential complications associated with use of this device.

Information contained herein for distribution outside the U.S. and Japan only. Please check the regulatory status of the device before distribution in areas where CE marking is not the regulation in force. For more information, visit our web site at abbottvascular.com

©2009 Abbott Laboratories. All rights reserved. EHQ 4001-01 09/2009 Rev A

