

INTRODUCING Axium™ Prime Detachable Framing Coils

Order Number	Diameter (mm)	Length (cm)	O.D. (in)	Calculated Coil Volume (mm ³)
FC-3-6-3D	3	6	0.0115	4.02
FC-3-8-3D	3	8	0.0115	5.36
FC-3-10-3D	3	10	0.0115	6.70
FC-3.5-6-3D	3.5	6	0.0115	4.02
FC-3.5-8-3D	3.5	8	0.0115	5.36
FC-3.5-10-3D	3.5	10	0.0115	6.70
FC-4-6-3D	4	6	0.0125	4.75
FC-4-8-3D	4	8	0.0125	6.33
FC-4-10-3D	4	10	0.0125	7.92
FC-4-12-3D	4	12	0.0125	9.50
FC-4-15-3D	4	15	0.0125	11.88
FC-5-8-3D	5	8	0.0125	6.33
FC-5-10-3D	5	10	0.0125	7.92
FC-5-15-3D	5	15	0.0125	11.88
FC-5-20-3D	5	20	0.0125	15.83
FC-6-10-3D	6	10	0.0125	7.92
FC-6-15-3D	6	15	0.0125	11.88
FC-6-20-3D	6	20	0.0125	15.83
FC-6-25-3D	6	25	0.0125	19.79
FC-7-12-3D	7	12	0.0135	11.08
FC-7-15-3D	7	15	0.0135	13.85
FC-7-20-3D	7	20	0.0135	18.47
FC-7-30-3D	7	30	0.0135	27.70

Order Number	Diameter (mm)	Length (cm)	OD (in)	Calculated Coil Volume (mm ³)
FC-8-15-3D	8	15	0.0135	13.85
FC-8-20-3D	8	20	0.0135	18.47
FC-8-30-3D	8	30	0.0135	27.70
FC-9-20-3D	9	20	0.0135	18.47
FC-9-30-3D	9	30	0.0135	27.70
FC-10-20-3D	10	20	0.0135	18.47
FC-10-30-3D	10	30	0.0135	27.70
FC-10-40-3D	10	40	0.0135	36.94
FC-12-30-3D	12	30	0.0145	31.96
FC-12-40-3D	12	40	0.0145	42.61
FC-12-50-3D	12	50	0.0145	53.27
FC-14-30-3D	14	30	0.0145	31.96
FC-14-40-3D	14	40	0.0145	42.61
FC-14-50-3D	14	50	0.0145	53.27
FC-16-40-3D	16	40	0.0145	42.61
FC-16-50-3D	16	50	0.0145	53.27
FC-18-40-3D	18	40	0.0145	42.61
FC-18-50-3D	18	50	0.0145	53.27
FC-20-50-3D	20	50	0.0145	53.27
FC-22-50-3D	22	50	0.0145	53.27
FC-25-50-3D	25	50	0.0145	53.27

Medtronic International Trading Sàrl.

Route du Molliau 31
Case postale
CH-1131 Tolochenaz
Tel. +41 (0)21 802 70 00
Fax +41 (0)21 802 79 00

medtronic.eu

Calculated coil volume determined using formula for cylinder (volume = $\pi r^2 \times \text{length}$).

The Axium™ and Axium™ Prime detachable coils are not intended for all patients and may not be the appropriate treatment for all clinical scenarios. See Instructions for Use for a complete list of warnings, precautions, and contraindications. The Axium™ Prime detachable coil system is indicated for the endovascular embolization of intracranial aneurysms and other neurovascular abnormalities, such as arteriovenous malformations and arteriovenous fistulae. The Axium™ Prime detachable coils are also indicated for arterial and venous embolizations in the peripheral vasculature.

Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

UC201802978EE © 2017 Medtronic.
All Rights Reserved. Printed in Europe.

Medtronic
Further, Together