


4. Optical bench



Mark	P/n	Designation	Qty	Sales unit	Compatibility
1	XDA838B (1209132838)	OPTICAL, OPT. BENCH LAMP LOCK	1	1	ALL
2	FCA026A (1202631026)	OPTICAL, OPTIC BENCH MIRROR	1	1	ALL
3	1300028830 GBM1326 (1203601326)	OPTICAL, LAMP P400 OPTICAL, LAMP PC400	1	1	P400 PC400
4	XDA858AS (1209139858)	OPTICAL, CUV. EMISSION GUN+LENS	1	1	ALL
5	XDA857BS (1209139857)	OPTICAL, LAMP EMISSION GUN+LENS	1	1	ALL
6	XEC106AS (1209179106)	KIT, OPTICAL BENCH+COVERS	1	1	ALL
7	XDA854D (1209134854)	OPTICAL, SPECTRO P400 >V5.06	1	1	ALL*
8	XAA483CS (1209109483)	PCB, ACQU. BOARD P400 >V5.06	1	1	ALL*
9	XAA457BS (1209109457)	PCB, PREAMPLIFIER BOARD P400	1	1	ALL
10	XDB114A (1209141114)	OPTICAL, SPEC.EMISSION GUN+LENS	1	1	ALL



* For the ABX Pentra 400, the acquisition board ref. XAA483CS (1209109483) and the spectrophotometer ref. XDA854D (1209134854) are compatible only if the instrument is equipped with the software version V5.0.6 or higher.

The spectrophotometer ref. XDA854D (1209134854) is delivered with the acquisition board cable ref. DAD134C (1201923134). This cable is not compatible with the ABX Pentra 400 V1, the acquisition board cable ref. DAD134B (1201922134) must be used instead.