



Str. Grenoble, nr. 106 A, ap. 2, mun. Chişinău Tel: (+373) 22 92 02 20; Fax: (+373) 22 63 28 61 Mobil: 079920220; e-mail: office@nitech.md; www.nitech.md

Anexa A la Termenii și condițiile de furnizare

Detalii specifice și standarde tehnice

LOT "Microcentrifugă si Minicentrifugă":

№	Name of product	Minimal technical specifications	Offered specifications MIKRO 200R with accesories HETTICH- GERMANY
1	Microcentrifugă	The microcentrifuge has rotors that are designed to spin microliter and PCR tubes. - RPM: 500 - 15,000 min -1; Adjustable in increments of 10; Max. RCF: 21,382; Max. capacity: 30 x 2.0 ml; - Choice of 3 rotors; Impulse key for short cycle mode - Easy operation with keypad and control knob; - 4 program memories for more individuality; - 9 individual acceleration and deceleration stages; - coolable from -20 to +40 °C with pre-cooling function.	The microcentrifuge has rotors that are designed to spin microliter and PCR tubes. - RPM: 500 - 15,000 min -1; Adjustable in increments of 10; Max. RCF: 21,382; Max. capacity: 30 x 2.0 ml; - Choice of 4 rotors; Impulse key for short cycle mode - Easy operation with keypad and control knob; - 4 program memories for more individuality; - 9 individual acceleration and deceleration stages; - coolable from -20 to +40 °C with pre-cooling function. Delivered with Rotor with lid – 2434 for 24 tubes 2ml
	Name of product	Minimal technical specifications	Offered specifications MICROSPIN 12 with accesories BIOSAN- LATVIA

		Mini centrifuge with tool-free snap spin	Mini centrifuge with tool-free snap spin rotors
		<u> </u>	
		rotors and quick-spin options.	and quick-spin options.
		1 Safety features help ensure it will not	1 Safety features help ensure it will not operate
		operate when the lid is open, a rotor is	when the lid is open, a rotor is jammed, or the
		jammed, or the unit is out of balance. Ideal	unit is out of balance. Ideal for most protocols
		for most protocols requiring fast spins	requiring fast spins (14500 rpm / 12400 x g), and
		(12,500 rpm / 9,800 x g), and appropriate for	appropriate for PCR, microfilter cell separation,
		PCR, microfilter cell separation, and HPLC	and HPLC protocols.
	Minicentrifugă	protocols.	- Select and set the digital speed display in RCF
2		- Select and set the digital speed display in	or RPM;
		RCF or RPM;	- Backlit display showing time, speed, and
		- Backlit display showing time, speed, and	operational information on the same screen;
		operational information on the same screen;	- UV-Resistant plastic, increased useful life
		- UV-Resistant plastic, increased useful life	under harsh light conditions:
		under harsh light conditions:	- double insulated, low voltage power supply-
		- double insulated, low voltage power	push button lid release;
		supply- push button lid release;	- Manual override allows the user to open the lid
		- Manual override allows the user to open	in case of power failure.
		the lid in case of power failure.	Delivered with rotor + lid for 12x 2ml tubes
		- Original manufacturer packing for all items;	- Original manufacturer packing for all items;
	Mandatory requirements	- Quality Certificates (CE and/or DC) or any	- Quality Certificates (CE and/or DC) or any other
	_	other equivalent;	equivalent;
		- ISO 9001-2015	ISO 9001-2015

Toate bunurile vor fi testate de beneficiar. Acceptarea bunurilor și semnarea actelor de recepție va avea loc după testarea pozitivă efectuată de beneficiar.

NUMELE FURNIZORULUI NITECH SRL	Semnătura și ștampila autorizate	
Locul: CHISINALL Rep. Of Moldova	Data:16.02.2023	