

QUALITY CONTROL DATA SHEET – FICHE DE CONTRÔLE QUALITÉ Brucellosis Antigen for Complement Fixation Test

Product code / Code Produit : AG-BRU
Batch / N° de lot: 010 – Prod.: 02/2020 - Exp.: 02/2023

KIT COMPONENTS / COMPOSITION DU KIT

ID.vet

Brucellosis Antigen for CFT

Antigène Brucellique pour RFC **100 mL**
Brucella abortus Antigen für die KBR **(100X)**

In vitro veterinary diagnosis

 +2°C  +8°C **LOT** AG-BRU-010  02/2023
ZuLNr.: FLIC 049 Für Rinder, Schafe, Ziegen und Schweine

ANALYTICAL SENSITIVITY CONTROL / CONTRÔLE DE LA SENSIBILITE ANALYTIQUE

The analytical sensitivity is standardized according to the European directive CEE 64/432 and the OIE requirements. The French substandard (BRU POS SE 2) diluted at 1/200 must give a fixation to 50%.

La sensibilité analytique est standardisée selon la directive Européenne CEE 64/432 et les exigences de l'OIE. Le sous-étalon national (BRU POS SE 2) dilué au 1/200 doit donner une fixation à 50%.

SENSITIVITY CONTROL / CONTRÔLE DE LA SENSIBILITE

6 positive samples were tested including the reference panel from NRL (Anses): all 6 samples were found positive.

6 échantillons positifs testés, dont les sérums du panel de référence du LNR (Anses) : les 6 échantillons ont été trouvés positifs.

Reference panel from NRL (Anses) Panel de référence du LNR (Anses)	Result obtained (UI/ml) Résultats obtenus (en UI/ml)
SR6	320
SR8	320
SR9	266
SR11	640
SR12	532

SPECIFICITY CONTROL / CONTRÔLE DE LA SPECIFICITE

30 negative samples were tested: all samples were found negative.

30 sérums négatifs testés : tous les échantillons ont été trouvés négatifs.



Quality Control Manager : Anaïs Agnel
Responsable Contrôle Qualité
anaïs.agnel@id-vet.com



Director : Philippe Pourquier
Directeur
philippe.pourquier@id-vet.com