

Diagnostic Monitor RadiForce® GX240







2 megapixel high-resolution diagnostic monochrome monitor ideal for a wide variety of tasks from viewing CR, DR, MRI and CT images to use as a PACS/HIS/RIS terminal. High-definition and multifunctional model featuring precise calibration compliant with the DICOM Part 14 standard.

- LED backlight offers long service life and reduction in power consumption*.
- Calibration simplified with Integrated Front Sensor (IFS) that is housed within the front bezel and measures brightness and grayscale tones according to a user-determined schedule.
- Integrated presence sensor unites convenience with energy savings by ensuring that the monitor automatically conserves power when it is not in use.
- Hybrid Gamma function automatically distinguishes the medical images from the non-medical areas such as the tool palettes and displays each area at its optimum brightness.
- 10-bit (1,024 tones) simultaneous grayscale display extends grayscale fidelity to the boundaries of human visual perception abilities.
- * Based on the EIZO's internal measurement with the predecessor model under the same conditions.

Model Variations	GX240-CL: Clear Base GX240-CL-P: Pairing	OSD Languages	English, French, German, Italian, Japanese, Simplified Chinese, Spanish, Swedish, Traditional
Cabinet Color	Black		Chinese
Panel Type	TFT Monochrome LCD Panel (IPS)	Net Weight (With Stand	10.2 kg / 7.5 kg
Panel Size	54 cm / 21.3" (540 mm diagonal)	/ Without Stand)	
Display Size (H x V)	324.0 x 432.0 mm	Hole Spacing	VESA standard 100 x 100 mm
Pixel Pitch	0.270 x 0.270 mm	Certifications and Standards ¹	CE (Medical Device Directive), EN60601-1, UL60601-1, CSA C22.2 No. 601-1, IEC60601-1, VCCI-B, FCC-B, Canadian ICES-003-B, C-tick,
Native Resolution	1200 x 1600		
Grayscale Tones	10-bit (DisplayPort): 1,024 from a palette of 16,369 tones		RoHS, CCC, GOST-R
	8-bit: 256 from a palette of 16,369 tones	FDA 510(k) Clearance	Pending (for General Radiography ²)
Viewing Angles (H,V)	176°, 176°	Supplied Accessories	AC power cord, signal cable (DVI-D ~ DVI-D, DisplayPort ~ DisplayPort), USB cable, Utility Disk
Brightness (typical)	1,200 cd/m ²		(RadiCS LE, ScreenManager Pro for Medical, user's
Recommended	500 cd/m ²		manual), warranty card
Brightness for Calibration		Warranty	Five Years
Contrast Ratio (typical)	1400:1	,	
Response Time (typical)	40 ms (On/Off)	Dimensions (Unit: mm)	THEOD
Scanning Frequency (H, V)	31 - 100 kHz, 59 - 61 Hz (VGA Text: 69 - 71 Hz) Frame synchronous mode: 59 - 61 Hz	 376	
Dot Clock	164.5 MH z	B B B C C C C C C C C C C C C C C C C C	
Input Terminals	DVI-D x 1, DisplayPort x 1		
USB Ports / Standard	1 upstream, 2 downstream / Rev. 2.0		
Power Requirements	AC 100 - 120 V, 200 - 240 V: 50 / 60 Hz		
Maximum Power Consumption / Save Mode	76 W / Less than 1.6 W		
Power Management	DVI DMPM, DisplayPort 1.1a		
Sensor	Backlight Sensor, Integrated Front Sensor, Presence Sensor, Ambient Light Sensor		

¹ Please contact the EIZO subsidiary or distributor in your country for the latest information.

² Display of mammography images for diagnosis is not supported.



EIZO NANAO CORPORATION

153 Shimokashiwano, Hakusan, Ishikawa 924-8566 Japan Phone +81-76-277-6792 Fax +81-76-277-6793 WWW.eiZO.COM All product names are trademarks or registered trademarks of their respective companies. EIZO, RadiForce, FlexScan, ScreenManager, RadiCS and RadiNET are registered trademarks of Eizo Nanao Corporation. Specifications are subject to change without notice.

Copyright © 2012 EIZO NANAO CORPORATION. All rights reserved. (110810E)