





USER MANUAL [PRODUCT DESCRIPTION SHEET]

Issue No. : 1 Date : 01-09-2022



Rev. No. : 0 Date : 01-09-2022

Product Name	:	Angeltouch Pressure Monitoring Line		
Product Code	:	IT-35A		
CE Mark with No.	:			
Manufacturer's Name	:	Angioplast Private Limited		
Manufacturer's Address	:	Plot No. 4803, Phase IV, Vatva GIDC, Ahmedabad 382 445. [India]		
EC REP Address	:	mdi Europa GmbH, Langenhagener Str. 71, 30855, Langenhagen, Germany		
Intended Use	:	The sterile Set is used for extending the infusion set. It can withstands force of 800psi.		
Product Description	:	Set is made up of Protective Cap, Male Luer Lock, Tubing, Female Luer Lock, Protective Cap, Length 150cm.		
Product Specification	:			
Material List& Make	:	Sr. No.	Name of Components	No.
		1	Male Luer Lock	1
		2	Female Luer Lock	1
		3	Tubing	1
		4	Protective Cap	2
		5	Pouch	1
Instruction For Use	:	<ol style="list-style-type: none"> 1. Follow physician's instructions in opening package and use of contents. 2. Use aseptic technique. 3. Verify gastrostomy tube compatibility prior to use. 4. Attach inlet to outlet of administration set. 		

Prepared by		Approved by
	Page 1 of 2	
S. A Christian		M. M Patel

USER MANUAL [PRODUCT DESCRIPTION SHEET]**Issue No. : 1 Date : 01-09-2022****Rev. No. : 0 Date : 01-09-2022**

		5. Attach extension set outlet to feeding tube. 6. Close clamp of IV set and plug until ready to use or when not in use.
Product Use	:	For Single Use Only
Is product Sterilized	:	Yes
Sterilization Method	:	Ethylene Oxide Gas Sterilization
Storage Conditions	:	Keep away from getting wet. Store in cool & dry place.
Packing Detail		
Unit Pack [Pouch]	:	238mm x 120mm x Paper : 300G and Plastic : 200G
Gross Pack [Size & Qty.]	:	260mm x 520mm x 150G - 25 Nos.
Shipper Pack [Size & Qty.]	:	595mm x 375mm x 325mm - 600 Nos.

Prepared by		Approved by
	Page 2 of 2	
S. A Christian		M. M Patel