

Sales only to commercial customers. No sales to consumers within the meaning of §13 BGB.



Zoom in

Bacterial/Viral Filter, with HME

Bacterial/Viral Filter with HME and CO₂-port. The purpose of the filter is to protect the patient and to humidify the inhaled air. The filter is placed between the endotracheal tube/tracheal cannula and the Y-piece of the breathing system.

Shape (Ventilation Filter): Round

Inner diameter (patient end): 20/15 mm

ISO style (patient end): Male/Female

Sterile: No

Country availability: Please note: Not all products are available in all countries. Feel free to contact our international sales team about any questions relating to our sales and services.

Jutta Dorenberg

Director International Sales

Tel.: +49 2236 3913 - 405

eMail: jutta.dorenberg@dahlhausen.de

shorten text

Article no.	Application area	Inner diameter (machine end)	Particle retention	ISO style (machine end)	Max. tidal volume (Breathing Filter)	minimal quantity
4300906120	Child	13 mm	99,99 %	Male	900 ml	50
4300906020	Adult	22/13 mm	99,999 %	Male/Female	1.500 ml	25
4300905708	Adult	22/13 mm	99,999 %	Male/Female	1.500 ml	50