


INTRODUCERS

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.



INTRODUCER

Kit to introduce and exchange catheters and endocavitary electrodes



Intended purpose / Indications

The Introducer kit is intended to provide access and facilitate the percutaneous introduction of various devices into veins and/or arteries for a variety of diagnostic and therapeutic procedures.

Set components:

- sheath with a valve
- dilator
- guidewire
- angiographic needle


Reference:	Sheath size [F]	Sheath length [mm]	Colour identification	Dilator size [F]	Angiographic needle size [G]	Guidewire size [inch]	Guidewire length [cm]	Guidewire type
INT3F	3	60	white	3	21	0.018"	35	STRAIGHT
INT4F	4	70	red	4	19	0.022"	40	"J"
INT5F	5	110	grey	5	18	0.035"	40	
INT5F23		230					60	
INT6F	6	110	green	6			40	
INT6F23		230					60	
INT7F	7	110	orange	7			40	
INT7F23		230					60	
INT8F	8	110	blue	8			40	
INT8F23		230					60	
INT9F	9	110	black	9			40	
INT9F23		230					60	
INT10F	10	110	white	10			40	
INT10F23		230					60	
INT11F	11	110	blue	11			40	
INT11F23		230					60	

When ordering use Ref. Number, e.g.: INT6F.



INTRODUCER (FOR RADIAL ARTERY)

Kit to introduce and exchange catheters and endocavitary electrodes



Intended purpose / Indications

The Introducer kit is intended to provide access and facilitate the percutaneous introduction of various devices into veins and/or arteries for a variety of diagnostic and therapeutic procedures.

Set components:

- sheath with a valve
- dilator
- guidewire
- angiographic needle

The set contains a guidewire of special construction. The soft distal part of the guidewire enables its introduction to the vessel, and the stiff proximal part provides an scaffold for the dilator.

Reference:	Sheath size [F]	Sheath length [mm]	Colour identification	Dilator size [F]	Angiographic needle size [G]	Guidewire size [inch]	Guidewire length [cm]	Guidewire type
INT4FKR	4	70	red	4	21	0.018"	45	"J"
INT4FKRP						0.018"		STRAIGHT
INT5FKR	5	110	grey	5		0.018"	45	"J"
INT5FKRP						0.018"		STRAIGHT
INT5F21KRP						0.021"		STRAIGHT
INT5FKR23		230				0.018"	60	"J"
INT5FKR23P						0.018"		STRAIGHT
INT5F21KR23P						0.021"		STRAIGHT
INT6FKR	6	110	green	6		0.018"	45	"J"
INT6FKRP						0.018"		STRAIGHT
INT6F21KRP						0.021"		STRAIGHT
INT6FKR23		230				0.018"	60	"J"
INT6FKR23P						0.018"		STRAIGHT
INT6F21KR23P						0.021"		STRAIGHT
INT7FKR	7	110	orange	7		0.018"	45	"J"
INT7FKRP						0.018"		STRAIGHT
INT7F21KRP						0.021"		STRAIGHT
INT7FKR23		230				0.018"	60	"J"
INT7FKR23P						0.018"		STRAIGHT
INT7F21KR23P						0.021"		STRAIGHT

When ordering use Ref. Number, e.g.: INT5FKR.

Special order can contain introducer with radiopaque marker

Model reference number: INT7FKRM – introducer with a marker.


Hydrophilic coating is available on request.

Model reference number: INT7FKRH – introducer with a hydrophilic coating.



INTRODUCER BRAIDED (COIL TYPE)

Kit to introduce and exchange catheters and endocavitary electrodes



Intended purpose / Indications

The Introducer kit is intended to provide access and facilitate the percutaneous introduction of various devices into veins and/or arteries for a variety of diagnostic and therapeutic procedures.

Set components:

- sheath with a valve
- dilator
- angiographic needle

Hydrophilic coating
High flexibility & trackability
Distal platinum marker for high visibility under fluoroscopy

Reference:	Sheath size [F]	Sheath length [mm]	Dilator size [F]	Angiographic needle size [G]	Use with guidewire [inch]
INT4F45T	4	450	4	19	0.022"
INT4F60T		600			
INT4F90T		900			
INT4F100T		1000			
INT5F45T	5	450	5	18	0.035"
INT5F60T		600			
INT5F90T		900			
INT5F100T		1000			
INT6F45T	6	450	6		
INT6F60T		600			
INT6F90T		900			
INT6F100T		1000			
INT7F45T	7	450	7		
INT7F60T		600			
INT7F90T		900			
INT7F100T		1000			
INT8F45T	8	450	8		
INT8F60T		600			
INT8F90T		900			
INT8F100T		1000			
INT9F45T	9	450	9		
INT9F60T		600			
INT9F90T		900			
INT9F100T		1000			

When ordering use Ref. Number, e.g.: INT6F90T.


Other lengths in the range of 23 cm to 100 cm, available on request.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.



INTRODUCER (LARGE SIZES)

Kit to introduce and exchange catheters and endocavitary electrodes



Intended purpose / Indications

The Introducer kit is intended to provide access and facilitate the percutaneous introduction of various devices into veins and/or arteries for a variety of diagnostic and therapeutic procedures.

Set components:

- sheath with a valve
- dilator

Platinum marker for high visibility under fluoroscopy.

Reference:	Sheath size [F]	Sheath length [mm]	Dilator size [F]
INT12F	12	120	12
INT 12F23		230	
INT12F35		350	
INT14F	14	120	14
INT14F23		230	
INT14F35		350	
INT16F	16	120	16
INT16F23		230	
INT16F35		350	


When ordering use Ref. Number, e.g.: INT16F23.

NOTE: Use with guidewire 0.035" and angiographic needle size 18 G.



INTRODUCER (LONG)

Kit to introduce and exchange catheters and endocavitary electrodes



Intended purpose / Indications

The Introducer kit is intended to provide access and facilitate the percutaneous introduction of various devices into veins and/or arteries for a variety of diagnostic and therapeutic procedures.

Set components:

- sheath with a valve
- dilator
- angiographic needle

MOST COMMONLY USED SIZES

Reference:	Sheath size [F]	Sheath length [mm]	Colour identification	Dilator size [F]	Angiographic needle size [G]
INT5F55	5	550	grey	5	18
INT5F90		900			
INT6F90	6	900	green	6	
INT7F55	7	550	orange	7	

ON REQUEST (please contact us for more information)

Reference:	Sheath size [F]	Sheath length [mm]	Colour identification	Dilator size [F]	Angiographic needle size [G]
INT5F50	5	500	grey	5	18
INT5F60		600			
INT5F65		650			
INT5F70		700			
INT5F80		800			
INT6F50	6	500	green	6	
INT6F55		550			
INT6F60		600			
INT6F65		650			
INT6F70		700			
INT6F80		800			
INT7F50	7	500	orange	7	
INT7F60		600			
INT7F65		650			
INT7F70		700			
INT7F80		800			
INT7F90		900			
INT8F50	8	500	blue	8	
INT8F55		550			
INT8F60		600			
INT8F65		650			
INT8F70		700			
INT8F80		800			
INT8F90		900			

When ordering use Ref. Number, e.g.: INT5F55. Special order can contain introducer with radiopaque marker.

Model reference number: INT8F90M- introducer with a marker.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

