


5.4 Changing the Radiography Conditions

Press buttons or keys on the touch panel to change the radiography conditions.

- **Technique**..... ⑧
- **Focus**..... ②
- **kV**..... ①
- **Switch mA/mAs** ④
- **mA/mAs**..... ③
- **sec** ⑤
- **AEC (optional)**
 - AEC ON/OFF ⑥
 - Density ⑦
 - Photo pickup field..... ⑨
 - **Speed** ⑪



Reverts changed radiography conditions to their original settings.

Changed radiography conditions are displayed underlined.