

Product catalog

This document reflects your selection of products from the Hamilton Medical online product catalog. The content is based on the following parameters:

Creation Date: **12 April 2021**
Country settings: **International**

Your selection: **Accessories & consumables**
• Flow sensors

Ventilators
• HAMILTON-C6
• HAMILTON-C1

Search keywords

The online product catalog is subject to change without notice. The latest and most comprehensive product information is available online.

Legend

 Single use

 Reusable

 Autoclavable

 Adult

 Pediatric

 Neonatal



The packaging unit of each single product is marked within this box icon.

Adapters & filters

Product	Description
	281803 Adapter OD15 / ID15 For the connection of the neonatal airway adapter to the neonatal flow sensor. For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1 GTIN: 07630002803762

CO2 measurement

Product	Description
	281803 Adapter OD15 / ID15 For the connection of the neonatal airway adapter to the neonatal flow sensor. For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1 GTIN: 07630002803762

Flow sensors

Product

Description



281803

Adapter



OD15 / ID15



For the connection of the neonatal airway adapter to the neonatal flow sensor.



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1

GTIN: 07630002803762



281637

Flow sensor



OD15 x ID15 / OD22 connectors



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002800860

Dimensions: L: 188 cm / 74 in

Includes: [PN 279937](#) Calibration adapter for flow sensor



282049

Flow sensor



OD15 x ID15 / OD22 connectors



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002800877

Dimensions: L: 260 cm / 102 in

Includes: [PN 279937](#) Calibration adapter for flow sensor



282051

Flow sensor



OD15 x ID15 / OD22 connectors



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002800891

Dimensions: L: 330 cm / 130 in

Includes: [PN 279937](#) Calibration adapter for flow sensor



279937

Calibration adapter for flow sensor



For use with flow sensor PN281637, PN282049, PN282051



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO



GTIN: 07630002800914

Product

Description



260177

Flow sensor



OD15 x ID15 (patient side)



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002802840

Dimensions: L: 160 cm / 63 in

Notes: If using Hamilton Medical Capnostat 5 mainstream airway adapters (PN 281720, PN 281722) or a third-party CO2 monitoring adapter, you must also order the adapter OD15 / ID15, PN 281803.

Includes: [PN 279964](#) Calibration adapter for flow sensor



155500

Flow sensor



OD15 x ID15 (patient side)



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002801362

Dimensions: L: 188 cm / 74 in

Notes: If using Hamilton Medical Capnostat 5 mainstream airway adapters (PN 281720, PN 281722) or a third-party CO2 monitoring adapter, you must also order the adapter OD15 / ID15, PN 281803.

Includes: [PN 279964](#) Calibration adapter for flow sensor



260179

Flow sensor



OD15 x ID15 (patient side)



For ventilator: HAMILTON-S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002802857

Dimensions: L: 310 cm / 122 in

Notes: If using Hamilton Medical Capnostat 5 mainstream airway adapters (PN 281720, PN 281722) or a third-party CO2 monitoring adapter, you must also order the adapter OD15 / ID15, PN 281803.

Includes: [PN 279964](#) Calibration adapter for flow sensor



279964

Calibration adapter for flow sensor



For use with neo flow sensor: PN 155500, PN 260177, PN260179



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002804417

Product

Description



950185

Flow sensor



OD15 x ID15 x OD22 connectors

Includes 3 tubing attachment straps



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1

GALILEO, RAPHAEL



GTIN: 07630002803489

Dimensions: L: 188 cm / 74 in

Notes: The accuracy of the sensor is reduced to +/- 15% or 15 ml whichever is greater, even after calibration when used with the RAPHAEL ventilator.

Includes: [PN 282323](#) Calibration adapter for flow sensor



Intelligent Ventilation since 1983

Manufacturer:

Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

☎ +41 (0)58 610 10 20

info@hamilton-medical.com

www.hamilton-medical.com

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets.
INTELLiVENT-ASV is not available in the US. For all proprietary trademarks and third-party trademarks used by Hamilton Medical AG see
www.hamilton-medical.com/trademarks. © 2017 Hamilton Medical AG. All rights reserved.