

Cat. No.	Category	Test Item	Diseases	Specimen Type	Specimen Volume	Reaction Time	Measuring Range	Qualification	
W290	COVID-19	2019-nCoV RBD Antibody Test	COVID-19	Fingerstick or Venous WB/ S/P	20µL	15 mins	Quantitative	CE	
W295		2019-nCoV RBD IgG Test		Fingerstick or Venous WB/ S/P	20µL	15 mins	Quantitative	CE	
W286		2019-nCoV Antigen Test		NP/ OP	/	15 mins	Qualitative	CE	
W276		2019-nCoV Antibody Test				10 µL	10 mins	Qualitative	CE
W278		2019-nCoV IgM/IgG Test				S/P/WB	10 µL	10 mins	Qualitative
W241P0002	Bone Metabolism	Vitamin D	Bone Metabolism	S/P	75 µL	15 mins	5~100 ng/mL	CE NMPA	
W207	Diabetes	HbA1c	Diabetes Mellitus	WB	10 µL	5 mins	4.0~14.5%	CE IFCC NGSP NMPA	
W206		One Step MAU	Renal Injury	Urine	75 µL	3 mins	5~300 mg/L	CE NMPA	
W219	Kidney Function Markers	CysC	Renal Injury	S/P/WB	10 µL	5 mins	0.2~10 mg/L	CE NMPA	
W228		NGAL	Renal Injury	Urine	75 µL	15 mins	10~1500 ng/mL	CE NMPA	
W229		β2-MG	Proximal Renal Tubule Dysfunction	S/P/WB	10 µL	10 mins	0.3~20 mg/L	CE NMPA	
W231	Thyroid Function	T3	Thyroid Malfunction	S/P/WB	75 µL	15 mins	0.61~9.22 nmol/L	CE NMPA	
W232		T4	Thyroid Malfunction	S/P/WB	75 µL	15 mins	12.87~300 nmol/L	CE NMPA	
W220		TSH	Thyroid Malfunction	S/P/WB	75 µL	15 mins	0.1~100 mIU/L	CE NMPA	
W200		ft3	Thyroid Malfunction	S/P/WB	75 µL	15 mins	0.40~50.00 pmol/L or 0.26~32.55 pg/mL	CE	
W259		ft4	Thyroid Malfunction	S/P/WB	75 µL	15 mins	1.00~100 pmol/L or 0.078~7.77 ng/dL	CE	
W225	Sex Hormone Markers	β-HCG	Fertility	S/P/WB	20 µL	15 mins	2~200,000 mIU/mL	CE NMPA	
W244		LH	Fertility	S/P/WB	75 µL	15 mins	1~100 mIU/mL	CE NMPA	
W245		FSH	Fertility	S/P/WB	75 µL	15 mins	1~100 mIU/mL	CE NMPA	
W233		Prog	Fertility	S/P/WB	75 µL	15 mins	1.5~60 ng/mL	CE NMPA	
W246P0001		PRL	Fertility	S/P/WB	75 µL	15 mins	1~200 ng/mL	CE NMPA	
W247-C4P1-E		E2	Fertility	S/P/WB	75 µL	15 mins	9~3000 pg/mL	CE NMPA	
W248-C4P1-E		Testosterone	Fertility	S/P/WB	75 µL	15 mins	0.2~15 ng/mL	CE NMPA	
W257P0001		AMH	Fertility	S/P/WB	75 µL	15 mins	0.1~16 ng/mL	CE NMPA	
W202	Cardiac Markers	NT-proBNP	Heart Failure	S/P/WB	75 µL	15 mins	18~35000 pg/mL	CE NMPA	
W203		cTn I	Myocardial Infarction	S/P/WB	75 µL	15 mins	0.1~50 ng/mL	CE NMPA	
W204		Myo	Myocardial Infarction	S/P/WB	75 µL	15 mins	2.0~400 ng/mL	CE NMPA	
W205		CK-MB	Myocardial Infarction	S/P/WB	75 µL	15 mins	0.3~100 ng/mL	CE NMPA	
W217		H-FABP	Myocardial Infarction	S/P/WB	75 µL	15 mins	1~120 ng/mL	CE NMPA	
W216		3 in 1 (cTn I/Myo/CK-MB)	Myocardial Infarction	S/P/WB	75 µL	15 mins	Same with single items	CE NMPA	
W222		BNP	Heart Failure	P/WB	75 µL	15 mins	5~5000 pg/mL	CE NMPA	
W223		2 in 1 (cTn I/ NT-proBNP)	Myocardial Infarction	S/P/WB	75 µL	15 mins	Same with single items	CE NMPA	
W240		hsCRP	CVD	S/P/WB	10 µL	3 mins	0.1~10mg/L	CE NMPA	
W254P0001		One Step cTnT	Myocardial Infarction	S/P/WB	100 µL for whole blood 75 µL for plasma/serum	15 mins	0.03~10 ng/mL	CE NMPA	
W242	Lp-PLA2	Atherosclerosis	S/P/WB	75 µL	15 mins	10~1000 ng/mL	NMPA		
W239	Brain Injury Markers	S100β	Cerebral Ischemia Cerebral Hemorrhage	S/P/WB	75 µL	15 mins	0.05~10 ng/mL	NMPA	
W211	Thrombus Markers	D-Dimer	Thrombosis	P/WB	15 µL for whole blood 10 µL for plasma	5 mins	0.1~10 mg/L	CE NMPA	
W227		One step D-Dimer	Thrombosis	P/WB	100 µL for whole blood 75 µL for plasma/serum	15 mins	0.1~10 mg/L	CE NMPA	
W210/ W210P0004		PCT	Sepsis	S/P/WB	75 µL for whole blood 50 µL for plasma/serum	15 mins	0.1~100 ng/mL	CE NMPA	
W201	Inflammation Markers	CRP (hsCRP+CRP)	Infection	S/P/WB	8.5 µL for whole blood 5 µL for plasma/serum	3 mins	0.5~200 mg/L	CE NMPA	
W218		2 in 1 (CRP+PCT)	Infection	S/P/WB	75 µL	15 mins	CRP:0.5~150 mg/L PCT:0.1~100 ng/mL	CE NMPA	
W221		SAA	Infection	S/P/WB	8.5 µL for whole blood 5 µL for plasma/serum	75 µL	3 mins	1 ~ 500 mg/L	CE NMPA
W251P0001		IL-6	Infection	S/P/WB	75 µL	15 mins	3~4000 pg/mL	CE NMPA	
W208	Tumor Markers	AFP	Liver Cancer	S/P/WB	75 µL	15 mins	5~400 ng/mL	CE NMPA	
W209		PSA	Prostate Cancer	S/P/WB	75 µL	15 mins	2~100 ng/mL	CE <sub>0123</sub> NMPA	
W226		CEA	Cancer	S/P/WB	75 µL	15 mins	1~500 ng/mL	CE NMPA	
W230		fPSA	Prostate Cancer	S/P/WB	75 µL	15 mins	0.2~30 ng/mL	CE <sub>0123</sub> NMPA	
W279P0001		iFOB	Gastrointestinal bleeding	Feces	stick to sample at 6 different parts of the feces	5 mins	25~1000 ng/mL	CE	
W282P0001	Infections Disease	Dengue NS1 Ag	Infectious disease	S/P/WB	75 µL	12 mins	Qualitative	CE	
W270P0001		H.pylori Ag	Infectious disease	Feces	6 different parts of the feces	10 mins	Qualitative		
W298P0001	Allergy	Total IgE	Allergic diseases	S/P/WB	75 µL	10 mins	1~1000 IU/mL	CE	
W250	Hypertension	Cortisol	Cortisol disorders	S/P/WB	75 µL	15 mins	50~1000 nmol/L	CE NMPA	
W269P0001	Anemia Markers	Ferritin	Iron metabolism	S/P/WB	75 µL	10 mins	5 ~ 1000 ng/mL	CE	
W297P0001	Rheumatoid Arthritis	ASO	Rheumatic fever and glomerulonephritis	S/P/WB	10 µL	5 mins	20~1000 IU/mL	CE	
W262P0001		RF	Rheumatoid Arthritis	S/P/WB	10 µL	5 mins	10 ~ 650 IU/mL		
/		Anti-CCP	Rheumatoid Arthritis	S/P/WB	10 µL	10 mins	3.5-300 U/mL		



# Finecare™ FIA Meter Plus (FS-113)

REF W915

Fluorescence Immunoassay Analyzer

CE NGSP IFCC

Racing For Life



# Fineware™ FIA METER PLUS

- Convenient**
- Reliable Accuracy**
- Temperature control**
- Rapid Test**

- Portable  
Easy to operate
- Good correlation  
with CLIA method
- Provide better  
reaction condition
- Get results in  
3-15 minutes

## FS-113

- Dimension: 270\*238\*146mm
- Weight: 2.8 Kg

## Connection

- Direct connection to LIS/HIS

## ID chip inlet

- Automatically reading of batch number and calibration information in the ID chip

## Designed Sampler

- Specially-designed sampler for small volume of specimen

Inner Temperature Control



## Built-in Thermal Printer

- Real-time auto-printing
- Support for external printers

## Large Display

- 8" LCD Touch Screen
- Intuitive user interface


## Test Channel

- Multiple testing modes available
- Unique quick mode for large number of samples

## Test Reagent

- Room temperature storage
- Suitable for whole blood, serum, plasma and urine

## TEST PROCEDURE



## APPLICATION

