



SCI - MFI: 101

Mini L Plate 2.0 MM

Shaft with two holes & head with two holes.

Profile

Thickness 1.0 mm
Width 5.0 mm
Hole Spacing 6.0 mm

Cat No.	Side
SCI-101-01	Left
SCI-101-10	Right



SCI - MFI: 102

Mini T Plate 2.0 MM

Head with two holes.

Profile

Thickness 1.0 mm
Width 5.0 mm
Hole Spacing 6.0 mm

Cat No.	No. of Holes in Shaft
SCI-102-02	2
SCI-102-03	3



SCI - MFI: 103

DCP 2.0 MM

Profile

Thickness (4 to 6 holes) 1.0 mm
Thickness (7 & 8 holes) 1.5 mm
Width 5.0 mm
Hole Spacing 5.0 mm

Cat No.	Holes
SCI-103-03	3
SCI-103-04	4
SCI-103-05	5
SCI-103-06	6
SCI-103-07	7
SCI-103-08	8



SCI - MFI: 104

DCP 2.7 MM

Profile

Thickness (2 to 6 holes) 2.0 mm
Thickness (7 to 12 holes) 2.5 mm
Width 8.0 mm
Hole Spacing 8.0 mm & 12.0 mm

Cat No.	Holes
SCI-104-02	2
SCI-104-03	3
SCI-104-04	4
SCI-104-05	5
SCI-104-06	6
SCI-104-07	7
SCI-104-08	8
SCI-104-09	9
SCI-104-10	10
SCI-104-12	12



SCI - MFI: 105

Mini L Plate 2.7 MM

Shaft with three holes & head with two holes.

Profile

Thickness 1.2 mm
Width 7.0 mm
Hole Spacing 8.0 mm

Cat No.	Side
SCI-105-01	Left
SCI-105-10	Right



SCI - MFI:106

Mini T Plate 2.7 MM

Shaft with three holes & head with two holes.

Profile

Thickness 1.2 mm
Width 7.0 mm
Hole Spacing 8.0 mm

Cat No.	No. of Holes in Shaft
SCI-106-02	2
SCI-106-03	3



SCI - MFI: 107

Straight Reconstruction Plate 2.7 MM

Profile

Thickness 3.1 mm
Width 8.0 mm
Hole Spacing 8.0 mm

Cat No.	Holes	Cat No.	Holes
SCI-107-02	2	SCI-107-11	11
SCI-107-03	3	SCI-107-12	12
SCI-107-04	4	SCI-107-14	14
SCI-107-05	5	SCI-107-16	16
SCI-107-06	6	SCI-107-18	18
SCI-107-07	7	SCI-107-20	20
SCI-107-08	8	SCI-107-22	22
SCI-107-09	9	SCI-107-24	24
SCI-107-10	10		



SCI - MFI: 108

Quarter Tubular Plate 2.7 MM

Profile

Thickness 1.0 mm
Width 7.0 mm
Hole Spacing 8.0 mm

Cat No.	Holes	Cat No.	Holes
SCI-108-02	2	SCI-108-07	7
SCI-108-03	3	SCI-108-08	8
SCI-108-04	4	SCI-108-09	9
SCI-108-05	5	SCI-108-10	10
SCI-108-06	6		



SCI - SFI: 300

Cortical Screw 3.5 mm

Profile

Thread Diameter 3.5 mm
 Drill Bit for Threaded Hole 2.5 mm
 Drill Bit for Gliding Hole 3.5 mm
 Core Diameter 2.4 mm
 Hexagonal Socket 2.5 mm
 Head Diameter 6.0 mm

Cat No.	Length	Cat No.	Length
SCI-300-10	10mm	SCI-300-32	32mm
SCI-300-12	12mm	SCI-300-34	34mm
SCI-300-14	14mm	SCI-300-36	36mm
SCI-300-16	16mm	SCI-300-38	38mm
SCI-300-18	18mm	SCI-300-40	40mm
SCI-300-20	20mm	SCI-300-42	42mm
SCI-300-22	22mm	SCI-300-44	44mm
SCI-300-24	24mm	SCI-300-46	46mm
SCI-300-26	26mm	SCI-300-48	48mm
SCI-300-28	28mm	SCI-300-50	50mm
SCI-300-30	30mm		



SCI - SFI: 301

Cancellous Screw 3.5 mm

Profile

Thread Diameter 3.5 mm
 Drill Bit for Threaded Hole 2.5 mm
 Drill Bit for Gliding Hole 2.5 mm
 Core Diameter 1.9 mm
 Hexagonal Socket 2.5 mm
 Head Diameter 6.0 mm

Cat No.	Length	Cat No.	Length
SCI-301-10	10mm	SCI-301-28	28mm
SCI-301-12	12mm	SCI-301-30	30mm
SCI-301-14	14mm	SCI-301-35	35mm
SCI-301-16	16mm	SCI-301-40	40mm
SCI-301-18	18mm	SCI-301-45	45mm
SCI-301-20	20mm	SCI-301-50	50mm
SCI-301-22	22mm	SCI-301-55	55mm
SCI-301-24	24mm	SCI-301-60	60mm
SCI-301-26	26mm		



SCI - SFI: 302

Cancellous Bone Screw 4.0 mm Short Thread

Profile

Thread Diameter	4.0 mm
Drill Bit for Threaded Hole	2.5 mm
Shaft Diameter	2.4 mm
Core Diameter	1.9 mm
Hexagonal Socket	2.5 mm
Head Diameter	6.0 mm
4.0mm Cancellous Bone Screw (Short Thrd.)	

Cat No.	Length	Cat No.	Length
SCI-302-10	10mm	SCI-302-30	30mm
SCI-302-12	12mm	SCI-302-35	35mm
SCI-302-14	14mm	SCI-302-40	40mm
SCI-302-16	16mm	SCI-302-45	45mm
SCI-302-18	18mm	SCI-302-50	50mm
SCI-302-20	20mm	SCI-302-55	55mm
SCI-302-22	22mm	SCI-302-60	60mm
SCI-302-24	24mm	SCI-302-65	65mm
SCI-302-26	26mm	SCI-302-70	70mm
SCI-302-28	28mm		



SCI - SFI: 303

Cancellous Bone Screw 4.0 mm Fully Threaded

Profile

Thread Diameter	4.0 mm
Drill Bit for Threaded Hole	2.5 mm
Core Diameter	1.9 mm
Hexagonal Socket	2.5 mm
Head Diameter	6.0 mm

4.0mm Cancellous Bone Screw (Fully Thrd.)

Cat No.	Length	Cat No.	Length
SCI-303-10	10mm	SCI-303-30	30mm
SCI-303-12	12mm	SCI-303-35	35mm
SCI-303-14	14mm	SCI-303-40	40mm
SCI-303-16	16mm	SCI-303-45	45mm
SCI-303-18	18mm	SCI-303-50	50mm
SCI-303-20	20mm	SCI-303-55	55mm
SCI-303-22	22mm	SCI-303-60	60mm
SCI-303-24	24mm	SCI-303-65	65mm
SCI-303-26	26mm	SCI-303-70	70mm
SCI-303-28	28mm		



SCI - SFI: 304

Cannulated Bone Screw 4.0 mm Short Thread

Profile

Thread Diameter 4.0 mm
 Drill Bit for Threaded Hole 2.5 mm
 Shaft Diameter 2.4 mm
 Core Diameter 1.9 mm
 Hexagonal Socket 2.5 mm
 Head Diameter 6.0 mm

Cat No.	Length	Cat No.	Length
SCI-304-20	20mm	SCI-304-40	40mm
SCI-304-22	22mm	SCI-304-42	42mm
SCI-304-24	24mm	SCI-304-44	44mm
SCI-304-26	26mm	SCI-304-46	46mm
SCI-304-28	28mm	SCI-304-48	48mm
SCI-304-30	30mm	SCI-304-50	50mm
SCI-304-32	32mm	SCI-304-52	52mm
SCI-304-34	34mm	SCI-304-54	54mm
SCI-304-36	36mm	SCI-304-56	56mm
SCI-304-38	38mm	SCI-304-58	58mm



SCI - SFI: 305

Cannulated Bone Screw 4.0 mm Fully Threaded

Profile

Thread Diameter 4.0 mm
 Drill Bit for Threaded Hole 2.5 mm
 Core Diameter 1.9 mm
 Hexagonal Socket 2.5 mm
 Head Diameter 6.0 mm

Cat No.	Length	Cat No.	Length
SCI-305-20	20mm	SCI-305-40	40mm
SCI-305-22	22mm	SCI-305-42	42mm
SCI-305-24	24mm	SCI-305-44	44mm
SCI-305-26	26mm	SCI-305-46	46mm
SCI-305-28	28mm	SCI-305-48	48mm
SCI-305-30	30mm	SCI-305-50	50mm
SCI-305-32	32mm	SCI-305-52	52mm
SCI-305-34	34mm	SCI-305-54	54mm
SCI-305-36	36mm	SCI-305-56	56mm
SCI-305-38	38mm	SCI-305-58	58mm



SCI - SFI: 10

Washer for 4.0 mm Screw


Cat No.	Dia
SCI-10-0	17.0mm



SCI - SFI: 11

Washer for 7.0 mm Screw

Cat No.	Length
SCI-11-10	13.0mm



Small Fragment Implants



SCI - SFI : 9

One Third Tubular Plate

Used as buttress plate for fibula, metatarsals and metacarpals. Fixation with 3.5 & 4.0 mm screws. Compression effect can be achieved by eccentric insertion of remote screws.

Profile

Thickness 1.0 mm
 Width 9.0 mm
 Hole Spacing 12.0 and 16.0 mm
 (One-third tube 12 12.0 and 16.0 mm)

Cat No	Holes	Cat No	Holes
SCI-109-02	2	SCI-109-07	7
SCI-109-03	3	SCI-109-08	8
SCI-109-04	4	SCI-109-09	9
SCI-109-05	5	SCI-109-10	10
SCI-109-06	6	SCI-109-12	12



SCI - SFI : 10

One Third Tubular Plate with Collar

Used as buttress plate for fibula, metatarsals and metacarpals. Fixation with 3.5 & 4.0 mm screws. Compression effect can be achieved by eccentric insertion of remote screws.

Profile

Thickness 1.0 mm
 Width 9.0 mm
 Hole Spacing 12.0 and 16.0 mm
 (One-third tube 12 mm dia)

Cat No	Holes	Cat No	Holes
SCI-110-02	2	SCI-110-07	7
SCI-110-03	3	SCI-110-08	8
SCI-110-04	4	SCI-110-09	9
SCI-110-05	5	SCI-110-10	10
SCI-110-06	6	SCI-110-12	12



SCI - SFI : 11

Small DCP 3.5 mm

For Ulna and Radius

Profile

Thickness 3.5 mm
 Width 10.0 mm
 Hole Spacing 12.0 and 16.0 mm

Cat No	Holes	Cat No	Holes
SCI-111-02	2	SCI-111-07	7
SCI-111-03	3	SCI-111-08	8
SCI-111-04	4	SCI-111-09	9
SCI-111-05	5	SCI-111-10	10
SCI-111-06	6	SCI-111-12	12



SCI - SFI : 12

Straight Reconstruction Plate 3.5 mm

Note: Do not Bend More than 15

Profile

Thickness 2.8 mm
 Width 10.0 mm
 Hole Spacing 12.0 mm

Cat No	Holes	Cat No	Holes
SCI-112-05	5	SCI-112-12	12
SCI-112-06	6	SCI-112-14	14
SCI-112-07	7	SCI-112-16	16
SCI-112-08	8	SCI-112-18	18
SCI-112-09	9	SCI-112-20	20
SCI-112-10	10	SCI-112-22	22



SCI - SFI : 113

T Plate Small, Right Angled

For the Left and Right Radius. Fixation with 3.5 & 4.0mm screws. Head with three holes.

Profile

Thickness 1.5 mm
Width 10.0 mm

Cat No	Holes	Cat No	Holes
SCI-113-03	3	SCI-113-07	7
SCI-113-04	4	SCI-113-08	8
SCI-113-05	5	SCI-113-09	9
SCI-113-06	6	SCI-113-10	10



SCI - SFI : 114

T Plate Small, Oblique Angle

For the Left and Right Radius. Fixation with 3.5 & 4.0mm screws. Head with three holes.

Profile

Thickness 1.5 mm
Width 10.0 mm

Cat No	Holes	Cat No	Holes
SCI-114-03	3	SCI-114-07	7
SCI-114-04	4	SCI-114-08	8
SCI-114-05	5	SCI-114-09	9
SCI-114-06	6	SCI-114-10	10



SCI - SFI : 115

Y Construction Plate 3.5 mm

For Elbow (Compound Fracture) . Fixation with 3.5 screws. Head with three holes.

Profile

Thickness 3.0 mm
Width 10.0 mm
Leaf Thickness 2.0 mm

Cat No	Holes	Cat No	Holes
SCI-115-03	3	SCI-115-07	7
SCI-115-04	4	SCI-115-08	8
SCI-115-05	5	SCI-115-09	9
SCI-115-06	6	SCI-115-10	10



SCI - SFI : 116

Cloverleaf Plate 3.5 mm

For distal Intra-Articular tibial Fractures. Fixation with 3.5 & 4.0 mm screws.

Profile

Thickness 2.0 mm
Width 15.2 mm
Cloverleaf Thickness 1.2 mm

Cat No	Holes	Cat No	Holes
SCI-116-03	3	SCI-116-07	7
SCI-116-04	4	SCI-116-08	8
SCI-116-05	5	SCI-116-09	9
SCI-116-06	6	SCI-116-10	10



SCI - SFI : 17

Semi Tubular Plate 3.5 mm

Fixation with 3.5 screws. Head with three holes.

Profile

Thickness 1.0 mm
Width 11.0 mm
Hole Spacing 16.0 & 26.0 mm

Cat No	Holes	Cat No	Holes
SCI-117-03	3	SCI-117-09	9
SCI-117-04	4	SCI-117-10	10
SCI-117-05	5	SCI-117-12	12
SCI-117-06	6	SCI-117-14	14
SCI-117-07	7	SCI-117-16	16
SCI-117-08	8	SCI-117-18	18



SCI - SFI : 18

Calcaneal Plate 3.5 mm

For Calcaneus Fractures , Fixation with 3.5 and 4.0 mm screws

Profile

Thickness 1.2 mm

Cat No	Length
SCI-118-60	60 mm
SCI-118-70	70 mm
SCI-118-80	80 mm
SCI-118-90	90 mm



SCI - SFI : 19

LC - DCP 3.5 MM

For Ulna & Radius. Fixation with 3.5 & 4.0 mm screws.

Profile

Thickness 3.0 mm
Width 10.0 mm
Leaf Thickness 2.0 mm

Cat No	Holes	Cat No	Holes
SCI-119-03	3	SCI-119-09	9
SCI-119-04	4	SCI-119-10	10
SCI-119-05	5	SCI-119-11	11
SCI-119-06	6	SCI-119-12	12
SCI-119-07	7	SCI-119-12	12
SCI-119-08	8	SCI-119-12	12



SCI - SFI : 20

Straight Reconstruction Plate 4.5 mm

Profile

Thickness 2.8 mm
Width 12.0 mm
Hole Spacing 16.0 mm

Cat No	Holes	Cat No	Holes
SCI-120-05	5	SCI-120-12	12
SCI-120-06	6	SCI-120-14	14
SCI-120-07	7	SCI-120-16	16
SCI-120-08	8	SCI-120-18	18
SCI-120-09	9	SCI-120-20	20
SCI-120-10	10	SCI-120-22	22
SCI-120-11	11		



SCI - LFI : 306

Cortical Screw 4.5 mm



SCI - LFI : 307

Malleolar Screw 4.5 mm

Self tapping Screw is available

Profile

Thread Diameter	4.5 mm
Drill Bit for Threaded Hole	3.2 mm
Drill Bit for Gliding Hole	4.5 mm
Core Diameter	3.0 mm
Hexagonal Socket	3.5 mm
Head Diameter	8.0 mm

Cat No	L	Cat No	L	Cat No	L
SCI-306-12	12 mm	SCI-306-34	34 mm	SCI-306-56	56 mm
SCI-306-14	14 mm	SCI-306-36	36 mm	SCI-306-58	58 mm
SCI-306-16	16 mm	SCI-306-38	38 mm	SCI-306-60	60 mm
SCI-306-18	18 mm	SCI-306-40	40 mm	SCI-306-62	62 mm
SCI-306-20	20 mm	SCI-306-42	42 mm	SCI-306-64	64 mm
SCI-306-22	22 mm	SCI-306-44	44 mm	SCI-306-66	66 mm
SCI-306-24	24 mm	SCI-306-46	46 mm	SCI-306-68	68 mm
SCI-306-26	26 mm	SCI-306-48	48 mm	SCI-306-70	70 mm
SCI-306-28	28 mm	SCI-306-50	50 mm	SCI-306-72	72 mm
SCI-306-30	30 mm	SCI-306-52	52 mm	SCI-306-74	74 mm
SCI-306-32	32 mm	SCI-306-54	54 mm		

Profile

Thread Diameter	4.5 mm
Drill Bit for Threaded Hole	3.2 mm
Shaft Diameter	3.0 mm
Core Diameter	3.0 mm
Hexagonal Socket	3.5 mm
Head Diameter	8.0 mm

Cat No	L	Cat No	L
SCI-307-12	12 mm	SCI-307-34	34 mm
SCI-307-14	14 mm	SCI-307-36	36 mm
SCI-307-16	16 mm	SCI-307-38	38 mm
SCI-307-18	18 mm	SCI-307-40	40 mm
SCI-307-20	20 mm	SCI-307-42	42 mm
SCI-307-22	22 mm	SCI-307-44	44 mm
SCI-307-24	24 mm	SCI-307-46	46 mm
SCI-307-26	26 mm	SCI-307-48	48 mm
SCI-307-28	28 mm	SCI-307-50	50 mm



**INTERNAL FIXATION™
SYSTEMS**



SCI - LFI : 308

**Cancellous Screw 6.5 mm
Thread Length 16 mm**

Profile

Thread Diameter	6.5 mm
Drill Bit for Threaded Hole	3.2 mm
Shaft Diameter	4.5 mm
Drill Bit for Shaft	4.5 mm
Core Diameter	3.0 mm
Hexagonal Socket	3.5 mm
Head Diameter	8.0 mm

Cat No	L	Cat No	L
SCI-308-25	25 mm	SCI-308-70	70 mm
SCI-308-30	30 mm	SCI-308-75	75 mm
SCI-308-35	35 mm	SCI-308-80	80 mm
SCI-308-40	40 mm	SCI-308-85	85 mm
SCI-308-45	45 mm	SCI-308-90	90 mm
SCI-308-50	50 mm	SCI-308-95	95 mm
SCI-308-55	55 mm	SCI-308-100	100 mm
SCI-308-60	60 mm	SCI-308-105	105 mm
SCI-308-65	65 mm	SCI-308-110	110 mm

Large Fragment Implants




SCI - LFI : 309

**Cancellous Screw 6.5 mm
Thread Length 32 mm**

Profile

Thread Diameter	6.5 mm
Drill Bit for Threaded Hole	3.2 mm
Shaft Diameter	4.5 mm
Drill Bit for Shaft	4.5 mm
Core Diameter	3.0 mm
Hexagonal Socket	3.5 mm
Head Diameter	8.0 mm

Cat No	L	Cat No	L
SCI-309-35	35 mm	SCI-309-75	75 mm
SCI-309-40	40 mm	SCI-309-80	80 mm
SCI-309-45	45 mm	SCI-309-85	85 mm
SCI-309-50	50 mm	SCI-309-90	90 mm
SCI-309-55	55 mm	SCI-309-95	95 mm
SCI-309-60	60 mm	SCI-309-100	100 mm
SCI-309-65	65 mm	SCI-309-105	105 mm
SCI-309-70	70 mm	SCI-309-110	110 mm



SCI - LFI : 310

**Cancellous Screw 6.5 mm
(Fully Threaded)**

Profile

Thread Diameter	6.5 mm
Drill Bit for Threaded Hole	3.2 mm
Core Diameter	3.0 mm
Hexagonal Socket	3.5 mm
Head Diameter	8.0 mm

Cat No	Length	Cat No	Length
SCI-310-25	25 mm	SCI-310-70	70 mm
SCI-310-30	30 mm	SCI-310-75	75 mm
SCI-310-35	35 mm	SCI-310-80	80 mm
SCI-310-40	40 mm	SCI-310-85	85 mm
SCI-310-45	45 mm	SCI-310-90	90 mm
SCI-310-50	50 mm	SCI-310-95	95 mm
SCI-310-55	55 mm	SCI-310-100	100 mm
SCI-310-60	60 mm	SCI-310-105	105 mm
SCI-310-65	65 mm	SCI-310-110	110 mm



SCI - LFI : 311

**Cancellous Screw 7.0 mm
(Thread Length 16 mm)**

Profile

Thread Diameter	7.0 mm
Drill Bit for Threaded Hole	4.5 mm
Shaft Diameter	4.5 mm
Core Diameter	4.5 mm
Hexagonal Socket	3.5 mm
Head Diameter	8.0 mm
Cannulation	2.1 mm
Guide Wire	2.0 mm

Cat No	Length	Cat No	Length
SCI-311-50	50 mm	SCI-311-90	90 mm
SCI-311-55	55 mm	SCI-311-95	95 mm
SCI-311-60	60 mm	SCI-311-100	100 mm
SCI-311-65	65 mm	SCI-311-105	105 mm
SCI-311-70	70 mm	SCI-311-110	110 mm
SCI-311-75	75 mm	SCI-311-115	115 mm
SCI-311-80	80 mm	SCI-311-120	120 mm
SCI-311-85	85 mm		




SCI - LFI : 312

**Large Cannulated Screw 7.0 mm
(Thread Length 32 mm)**

Profile

Thread Diameter	7.0 mm
Drill Bit for Threaded Hole	4.5 mm
Shaft Diameter	4.5 mm
Core Diameter	4.5 mm
Hexagonal Socket	3.5 mm
Head Diameter	8.0 mm
Cannulation	2.1 mm
Guide Wire	2.0 mm

Cat No	Length	Cat No	Length
SCI-312-35	35 mm	SCI-312-80	80 mm
SCI-312-40	40 mm	SCI-312-85	85 mm
SCI-312-45	45 mm	SCI-312-90	90 mm
SCI-312-50	50 mm	SCI-312-95	95 mm
SCI-312-55	55 mm	SCI-312-100	100 mm
SCI-312-60	60 mm	SCI-312-105	105 mm
SCI-312-65	65 mm	SCI-312-110	110 mm
SCI-312-70	70 mm	SCI-312-115	115 mm
SCI-312-75	75 mm	SCI-312-120	120 mm




SCI - LFI : 313

**Cannulated Screw 7.0 mm
(Thread Fully)**

Profile

Thread Diameter	7.0 mm
Drill Bit for Threaded Hole	4.5 mm
Shaft Diameter	4.5 mm
Core Diameter	4.5 mm
Hexagonal Socket	3.5 mm
Head Diameter	8.0 mm
Cannulation	2.1 mm
Guide Wire	2.0 mm

Cat No	Length	Cat No	Length
SCI-313-25	25 mm	SCI-313-70	70 mm
SCI-313-30	30 mm	SCI-313-75	75 mm
SCI-313-35	35 mm	SCI-313-80	80 mm
SCI-313-40	40 mm	SCI-313-85	85 mm
SCI-313-45	45 mm	SCI-313-90	90 mm
SCI-313-50	50 mm	SCI-313-95	95 mm
SCI-313-55	55 mm	SCI-313-100	100 mm
SCI-313-60	60 mm	SCI-313-105	105 mm
SCI-313-65	65 mm	SCI-313-110	110 mm



Large Fragment Implants



SCI - LFI : 314

**Cortical Cannulated Screw 4.5 mm
(Fully Threaded)**

Profile

Threaded Diameter 4.5 mm
Drill Bit for Threaded Hole 3.2 mm

Core Diameter 2.7 mm
Hexagonal Socket 3.5 mm
Head Diameter 6.5 mm
Cannulation 1.75 mm
Guide Wire 1.6 mm



SCI - LFI : 315

**Cortical Cannulated Screw 4.5 mm
(Short Threaded)**

Profile

Threaded Diameter 4.5 mm
Drill Bit for Threaded Hole 3.2 mm

Core Diameter 2.7 mm
Hexagonal Socket 3.5 mm
Head Diameter 6.5 mm
Cannulation 1.75 mm
Guide Wire 1.6 mm
Shaft Diameter 3.1 mm

Cat No	Length	Cat No	Length
SCI-314-20	20 mm	SCI-314-38	38 mm
SCI-314-22	22 mm	SCI-314-40	40 mm
SCI-314-24	24 mm	SCI-314-45	45 mm
SCI-314-26	26 mm	SCI-314-50	50 mm
SCI-314-28	28 mm	SCI-314-55	55 mm
SCI-314-30	30 mm	SCI-314-60	60 mm
SCI-314-32	32 mm	SCI-314-65	65 mm
SCI-314-34	34 mm	SCI-314-70	70 mm
SCI-314-36	36 mm		

Cat No	Length	Cat No	Length
SCI-315-20	20 mm	SCI-315-38	38 mm
SCI-315-22	22 mm	SCI-315-40	40 mm
SCI-315-24	24 mm	SCI-315-45	45 mm
SCI-315-26	26 mm	SCI-315-50	50 mm
SCI-315-28	28 mm	SCI-315-55	55 mm
SCI-315-30	30 mm	SCI-315-60	60 mm
SCI-315-32	32 mm	SCI-315-65	65 mm
SCI-315-34	34 mm	SCI-315-70	70 mm
SCI-315-36	36 mm		



SCI - LFI : 121

Narrow LC-DCP 4.5 mm



SCI - LFI : 122

Broad LC- DCP 4.5 mm

Fixation with 4.5 and 6.5 mm screws.
The LC-DCP holes allow compression in both directions of the longitudinal axis.
The end holes can be used for cancellous bone screws
Undercuts : For reduced vascular damage

Fixation with 4.5 and 6.5 mm screws.
The LC-DCP holes allow compression in both directions of the longitudinal axis. The end holes used for cancellous bone screws.

Undercuts : For reduced vascular damage.

Narrow LC-DCP 4.5 Profile

Thickness 4.6 mm
Width 13.5 mm Hole Spacing 18.0 mm

Cat No	Holes	Cat No	Holes
SCI-121-04	4	SCI-121-12	12
SCI-121-05	5	SCI-121-13	13
SCI-121-06	6	SCI-121-14	14
SCI-121-07	7	SCI-121-15	15
SCI-121-08	8	SCI-121-16	16
SCI-121-09	9	SCI-121-17	17
SCI-121-10	10	SCI-121-18	18
SCI-121-11	11		

Broad LC-DCP 4.5 Profile

Thickness 6.0 mm
Width 17.5 mm Hole Spacing 18.0 mm

Cat No	Holes	Cat No	Holes
SCI-122-04	4	SCI-122-12	12
SCI-122-05	5	SCI-122-13	13
SCI-122-06	6	SCI-122-14	14
SCI-122-07	7	SCI-122-15	15
SCI-122-08	8	SCI-122-16	16
SCI-122-09	9	SCI-122-17	17
SCI-122-10	10	SCI-122-18	18
SCI-122-11	11		



SCI - LFI : 123

Spoon Plate

Shaft with V-profile for the fixation on the frontal edge of the tibia. Fixation with 4.5 mm cortical screws

Profile

Thickness 2.0 mm
Head Width 34.0 mm

Cat No	Holes	Length
SCI-123-04	4	100 mm
SCI-123-05	5	120 mm
SCI-123-06	6	140 mm
SCI-123-07	7	140 mm

Fixation with 4.5 mm cortical screws.
A compression effect




SCI - LFI : 124

Semi Tubular Plate 4.5 mm

can be achieved by eccentric positioning of the screws remote from the fracture line. **Profile**

Thickness 1.0 mm
Width 11.0 mm
Hole Spacing 16.0 mm and 26.0 mm

Cat No	Holes	Cat No	Holes
SCI-124-04	4	SCI-124-12	12
SCI-124-05	5	SCI-124-13	13
SCI-124-06	6	SCI-124-14	14
SCI-124-07	7	SCI-124-15	15
SCI-124-08	8	SCI-124-16	16
SCI-124-09	9	SCI-124-17	17
SCI-124-10	10	SCI-124-18	18
SCI-124-11	11		



Large Fragment Implants



SCI - LFI :125

Narrow DCP 4.5 mm

Fixation with 4.5 and 6.5 mm screws. All holes are designed as self-compressing holes. All the end holes are designed for cancellous bone screws.
 Indication: As a neutralisation and tension band plate on the tibia (not intended for the Femur).

Profile

Thickness 3.6 mm
 Width 12.0 mm
 16.0 mm and 25.0 mm

Cat No	Holes	Cat No	Holes
SCI-125-04	4	SCI-125-10	10
SCI-125-05	5	SCI-125-11	11
SCI-125-06	6	SCI-125-12	12
SCI-125-07	7	SCI-125-13	13
SCI-125-08	8	SCI-125-14	14
SCI-125-09	9	SCI-125-16	16Hole Spacing



SCI - LFI : 126

Broad DCP 4.5 mm


Fixation with 4.5 and 6.5 mm screws. All holes are designed as self-compressing holes. All the end holes are designed for cancellous bone screws. Indication: As a tension band plate on the femur and for humeral pseudoarthrosis (not intended for the Tibia).

Use broad plates with six or seven holes in exceptional cases only, e.g. on the Humerus.

Profile

Thickness 4.5 mm
 Width 16.0 mm
 16.0 mm and 25.0 mm

Cat No	Holes	Cat No	Holes
SCI-126-06	6	SCI-126-12	12
SCI-126-07	7	SCI-126-13	13
SCI-126-08	8	SCI-126-14	14
SCI-126-09	9	SCI-126-15	15
SCI-126-10	10	SCI-126-16	16
SCI-126-11	11	SCI-126-18	18Hole Spacing



SCI - LFI :127

T Plate 4.5 mm

Shaft with round profile, for humeral and Tibial head. Fixation with 4.5 & 6.5 mm screws.

Profile

Thickness 3.6 mm
 Width 12.0 mm
 16.0 mm and 25.0 mm

Hole Spacing With double bend, as buttress for Tibial and Humeral head



SCI - LFI : 126

T Buttress Plate 4.5 mm

fractures. Left angled for right leg & right angled for left leg.

Profile

Thickness 2.0 mm
 Width 16.0 mm

T Plate 4.5 mm

Cat No (R)	Holes	Length
SCI-127-03	3	68 mm
SCI-127-04	4	84 mm
SCI-127-05	5	100 mm
SCI-127-06	9	116 mm
SCI-127-07	7	132 mm
SCI-127-08	8	148 mm
SCI-127-09	9	164 mm
SCI-127-10	10	176 mm

T Buttress Plate 4.5 mm Cat No (R) Holes Length

Cat No (R)	Holes	Length
SCI-126-03	3	68 mm
SCI-126-04	4	84 mm
SCI-126-05	5	100 mm
SCI-126-06	9	116 mm
SCI-126-07	7	132 mm
SCI-126-08	8	148 mm
SCI-126-09	9	164 mm
SCI-126-10	10	176 mm



SCI - LFI :129

Y-Reconstruction Plate 4.5 mm

For Elbow (Compound Fracture) . Fixation with 4.5 screws. Head with three holes.

Profile

Thickness 3.0 mm
Width 10.0 mm
2.0 mm

Cat No	Holes in shaft	Cat No	Holes in shaft
SCI-129-03	3	SCI-129-07	7
SCI-129-04	4	SCI-129-08	8
SCI-129-05	5	SCI-129-09	9
SCI-129-06	6	SCI-129-10	10

Leaf Thickness



SCI - LFI : 130


L- Buttress Plate 4.5 mm

For Elbow (Compound Fracture) . Fixation with 4.5 screws. Head with three holes.


Profile

Thickness 3.0 mm
Width 16.0 mm

Cat No (R)	Holes	Cat No (L)	Length
SCI-130-03	3	SCI-131-03	68 mm
SCI-130-04	4	SCI-131-04	84 mm
SCI-130-05	5	SCI-131-05	100 mm
SCI-130-06	6	SCI-131-06	116 mm
SCI-130-07	7	SCI-131-07	132 mm
SCI-130-08	8	SCI-131-08	148 mm
SCI-130-09	9	SCI-131-09	164 mm
SCI-130-10	10	SCI-131-10	176 mm



Large Fragment Implants



SCI - LFI :132

Condylar Buttress Plate 4.5 mm

With DCP dynamic compression holes and slot for tension device. As buttress for the femoral condyles with multiple fragment fractures.

Shaft Profile

Thickness 5.0 mm
Width 16.0 mm
16.0 mm

for Left Leg Cat No	Holes	for Right Leg Cat No
SCI-132-07	7	SCI-133-07
SCI-132-09	9	SCI-133-09
SCI-132-11	11	SCI-133-11
SCI-132-13	13	SCI-133-13
SCI-132-15	15	SCI-133-15

Hole Spacing



SCI - LFI : 134

Lateral Tibial Head Plate 4.5 mm

Fixation with 4.5 & 6.5mm screws shaft with DCP dynamic compression holes and slot for tension device. Left angled for right leg and right angled for left leg. **Shaft Profile**

Thickness 3.8 mm
Width 14.0 mm
16.0 mm

for Left Leg Cat No	Holes	for Right Leg Cat No
SCI-134-07	7	SCI-135-07
SCI-134-09	9	SCI-135-09
SCI-134-11	11	SCI-135-11
SCI-134-13	13	SCI-135-13
SCI-134-15	15	SCI-135-15

Hole Spacing



SCI-LFI: 136

Cobra Head Plate 4.5 mm

Shaft with DCP dynamic compression holes and slot for tension device. For arthrodesis of the hip. In the two shaft holes next to the head, two 4.5 mm cortical screws, length 68 mm, are normally used.

Shaft Profile

Thickness 6.0 mm
Width 16.0 mm
16.0 mm

Cat No	Holes	Cat No	Holes
SCI-136-08	8	SCI-136-12	12
SCI-136-09	9	SCI-136-14	14
SCI-136-10	10	SCI-136-16	16
SCI-136-11	11	Hole Spacing	



SCI-LFI: 137

Hook Plate 4.5 mm

Fixation with 4.5 mm Cortical Screws. Shaft with DCP dynamic compression holes and slot for tension device. **Shaft Profile**

Thickness 6.0 mm
Width 16.0 mm

Cat No	Holes
SCI-137-08	6
SCI-137-09	7
SCI-137-10	8

Profile

Thickness	5.6 mm
Width	16.0 mm
Hole Spacing	16.0 mm
U-profile blade	6.5 mm x 16.0 mm

With DCP dynamic compression holes and slot for tension device. All condylar plates have an angle of 95° between blade and shaft. The two holes next to the blade are suitable for 6.5mm Cancellous bone screws, all other holes are, self-compressing holes. Fixation with 4.5 mm cortical screws. For fractures in the distal and proximal femoral region as well as for intertrochanteric valgus osteotomy.



SCI-LFI: 138

Condylar Plate 95° for Adults

5 Holes		
Shaft length 92 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"138-50	50 mm	SCI"143-50
SCI"138-60	60 mm	SCI"143-60
SCI"138-70	70 mm	SCI"143-70
SCI"138-80	80 mm	SCI"143-80

7 Holes		
Shaft length 124 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"139-50	50 mm	SCI"144-50
SCI"139-60	60 mm	SCI"144-60
SCI"139-70	70 mm	SCI"144-70
SCI"139-80	80 mm	SCI"144-80

9 Holes		
Shaft length 156 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"140-50	50 mm	SCI"145-50
SCI"140-60	60 mm	SCI"145-60
SCI"140-70	70 mm	SCI"145-70
SCI"140-80	80 mm	SCI"145-80

12 Holes		
Shaft length 204 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"141-50	50 mm	SCI"146-50
SCI"141-60	60 mm	SCI"146-60
SCI"141-70	70 mm	SCI"146-70
SCI"141-80	80 mm	SCI"146-80

14 Holes		
Shaft length 236 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"142-50	50 mm	SCI"147-050
SCI"142-60	60 mm	SCI"147-060
SCI"142-70	70 mm	SCI"147-070
SCI"142-80	80 mm	SCI"147-080



SCI-LFI: 148

Angled Blade Plates 130° for Adults

Profile

Thickness	5.6 mm
Width	16.0 mm
Hole Spacing	16.0 mm
U-profile blade	6.5 mm x 16.0 mm

With DCP dynamic compression holes and slot for tension device (from 6 holes up). These plates have an angle of 130° between blade and shaft. Fixation with 4.5 mm cortical screws. For femoral neck and pertrochanteric fractures.

4 Holes		
Shaft length 60 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"148-50	50 mm	SCI"153-50
SCI"148-60	60 mm	SCI"153-60
SCI"148-70	70 mm	SCI"153-70
SCI"148-80	80 mm	SCI"153-80

6 Holes		
Shaft length 104 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"149-50	50 mm	SCI"154-50
SCI"149-60	60 mm	SCI"154-60
SCI"149-70	70 mm	SCI"154-70
SCI"149-80	80 mm	SCI"154-80

9 Holes		
Shaft length 152 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"150-50	50 mm	SCI"155-50
SCI"150-60	60 mm	SCI"155-60
SCI"150-70	70 mm	SCI"155-70
SCI"150-80	80 mm	SCI"155-80

12 Holes		
Shaft length 200 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"151-50	50 mm	SCI"156-50
SCI"151-60	60 mm	SCI"156-60
SCI"151-70	70 mm	SCI"156-70
SCI"151-80	80 mm	SCI"156-80

14 Holes		
Shaft length 232 mm		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"152-50	50 mm	SCI"156-050
SCI"152-60	60 mm	SCI"156-060
SCI"152-70	70 mm	SCI"156-070
SCI"152-80	80 mm	SCI"156-080

Blade Length : 55 mm, 65 mm, 75 mm, 85 mm and 90 mm will be manufactured against special order. # Holes : 5H, 7H, 8H, 10H and 11H will be manufactured against special order.



SCI-LFI: 157

Angled Blade Plate for Intertrochanteric Femoral Osteotomies in Adults 90°

With DCP dynamic compression holes and slot for tension device. Fixation with 4.5 mm cortical screws. The hole in the bend enables fixation of the calcar femoris with a long cortical screw.

Profile

Thickness	5.6 mm
Width	16.0 mm
Hole Spacing	16.0 mm
U-profile blade	6.5 mm x 16.0 mm



SCI-LFI: 160

Angled Blade Plate for Intertrochanteric Femoral Osteotomies in Adults 120°

With DCP dynamic compression holes and slot for tension device. Fixation with 4.5 mm cortical screws. The hole in the bend enables fixation of the calcar femoris with a long cortical screw.

Profile

Thickness	5.6 mm
Width	16.0 mm
Hole Spacing	16.0 mm
U-profile blade	6.5 mm x 16.0 mm

4 Holes		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"157-60	60 mm	SCI"159-60
SCI"157-65	65 mm	SCI"159-65
SCI"157-70	70 mm	SCI"159-70
SCI"157-75	75 mm	SCI"159-75
SCI"157-80	80 mm	SCI"159-80
SCI"157-85	85 mm	SCI"159-85

5 Holes		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"158-60	60 mm	SCI"160-60
SCI"158-65	65 mm	SCI"160-65
SCI"158-70	70 mm	SCI"160-70
SCI"158-75	75 mm	SCI"160-75
SCI"158-80	80 mm	SCI"160-80
SCI"158-85	85 mm	SCI"160-85

4 Holes		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"160-60	60 mm	SCI"162-60
SCI"160-65	65 mm	SCI"162-65
SCI"160-70	70 mm	SCI"162-70
SCI"160-75	75 mm	SCI"162-75
SCI"160-80	80 mm	SCI"162-80
SCI"160-85	85 mm	SCI"162-85

5 Holes		
Cat No. (DCP hole)	Blade Length	Cat No. (Round Hole)
SCI"161-60	60 mm	SCI"163-60
SCI"161-65	65 mm	SCI"163-65
SCI"161-70	70 mm	SCI"163-70
SCI"161-75	75 mm	SCI"163-75
SCI"161-80	80 mm	SCI"163-80
SCI"161-85	85 mm	SCI"163-85



SCI-LFI: 173

DCS Plate 95°

With DCP dynamic compression holes and slot for tension device. Fixation with 4.5 mm cortical screws. For condylar fractures, avoids further damage to the bone.

Profile

Thickness	5.4 mm
Width	16.0 mm
Hole Spacing	16.0 mm
Outer Barrel Diameter	12.5 mm
Barrel Length	25.0 mm

Cat. No.	Holes	Length
SCI-173-06	6	100 mm
SCI-173-07	7	115 mm
SCI-173-08	8	130 mm
SCI-173-09	9	146 mm
SCI-173-10	10	163 mm
SCI-173-12	12	198 mm
SCI-173-14	14	225 mm
SCI-173-16	16	260 mm
SCI-173-18	18	306 mm
SCI-173-20	20	338 mm
SCI-173-22	22	370 mm



SCI-LFI: 316

**Hip Screw with Compression Screw
(DHS/DCS Lag Screw)**

Use with DHS & DCS Plate

Profile

Thread Diameter	12.5 mm
Thread Length	22.0 mm
Shaft Diameter	8.0 mm

Cat. No.	Length	Cat No.	Length
SCI-316-50	50 mm	SCI-316-100	100 mm
SCI-316-55	55 mm	SCI-316-105	105 mm
SCI-316-60	60 mm	SCI-316-110	110 mm
SCI-316-65	65 mm	SCI-316-115	115 mm
SCI-316-70	70 mm	SCI-316-120	120 mm
SCI-316-75	75 mm	SCI-316-125	125 mm
SCI-316-80	80 mm	SCI-316-130	130 mm
SCI-316-85	85 mm	SCI-316-135	135 mm
SCI-316-90	90 mm	SCI-316-140	140 mm
SCI-316-95	95 mm	SCI-316-145	145 mm



SCI - LFI : 1

DHS Barrel Doubled Angled Plate 120°

With DCP dynamic compression holes and slot for tension device..
Fixation with 4.5 mm Cortical Screws.

Profile

Thickness	5.4 mm
Width	16.0 mm
Hole Spacing	16.0 mm
Outer Barrel Diameter	12.5 mm
Barrel Length	25.0 mm

Cat. No.	Holes	Cat. No.	Holes
SCI-174-04	4	SCI-174-07	7
SCI-174-05	5	SCI-174-08	8
SCI-174-06	6	SCI-174-10	10



SCI - LFI : 1

Compression Screw

For use with DHS and DCS plates, after mounting in leg screw its create compression.

Profile

Hexagonal Socket	3.5 mm
------------------	--------

Cat. No.	Length
SCI-317-36	36 mm



SCI - CI : 1

Wainwright Osteotomy Plate

Cat. No.	Holes	Cat. No.	Holes
SCI-175-04	4	SCI-175-07	7
SCI-175-05	5	SCI-175-08	8
SCI-175-06	6	SCI-175-10	10



SCI - CI : 2

Jewett Nail Plate 135°

Cat. No.	Holes	Cat. No.	Holes
SCI-176-04	4	SCI-176-07	7
SCI-176-05	5	SCI-176-08	8
SCI-176-06	6	SCI-176-10	10



SCI - SI : 1

Hart Shell


Cat. No.	Length
SCI- 400-75	75 mm
SCI- 400-100	100 mm
SCI- 400-125	125 mm
SCI- 400-150	150 mm
SCI- 400-175	175 mm
SCI- 400-200	200 mm
SCI- 400-225	225 mm
SCI- 400-250	250 mm



SCI - SI : 2

H Plate for Cervical Spine

Cat. No.	Holes	Hole Spacing	Length
SCI- 401-05	5	16 mm	23 mm
SCI- 402-05	5	18 mm	25 mm
SCI- 403-05	5	21 mm	28 mm
SCI- 404-08	8	16 + 18 mm	41 mm
SCI- 405-08	8	18 + 21 mm	46 mm
SCI- 406-11	11	16+18+18 mm	58 mm



SCI - SI : 3

Steffee Plate

Catalogue No.	Holes
SCI- 410-10	1
SCI- 410-15	1.5
SCI- 410-20	2
SCI- 410-25	2.5
SCI- 410-30	3
SCI- 410-35	3.5
SCI- 410-40	4
SCI- 410-45	4.5
SCI- 410-50	5



SCI - SI : 4

Steffee Screw

Length	Diameter			
	4.5 mm	5.5 mm	6.25 mm	7.0 mm
	Catalogue No.			
24 mm	SCI- 450-24	SCI- 451-24	SCI- 452-24	SCI- 453-24
26 mm	SCI- 450-26	SCI- 451-26	SCI- 452-26	SCI- 453-26
28 mm	SCI- 450-28	SCI- 451-28	SCI- 452-28	SCI- 453-28
30 mm	SCI- 450-30	SCI- 451-30	SCI- 452-30	SCI- 453-30
32 mm	SCI- 450-32	SCI- 451-32	SCI- 452-32	SCI- 453-32
34 mm	SCI- 450-34	SCI- 451-34	SCI- 452-34	SCI- 453-34
36 mm	SCI- 450-36	SCI- 451-36	SCI- 452-36	SCI- 453-36
38 mm	SCI- 450-38	SCI- 451-38	SCI- 452-38	SCI- 453-38
40 mm	SCI- 450-40	SCI- 451-40	SCI- 452-40	SCI- 453-40
42 mm	SCI- 450-42	SCI- 451-42	SCI- 452-42	SCI- 453-42
44 mm	SCI- 450-44	SCI- 451-44	SCI- 452-44	SCI- 453-44
46 mm	SCI- 450-46	SCI- 451-46	SCI- 452-46	SCI- 453-46
48 mm	SCI- 450-48	SCI- 451-48	SCI- 452-48	SCI- 453-48
50 mm	SCI- 450-50	SCI- 451-50	SCI- 452-50	SCI- 453-50
52 mm	SCI- 450-52	SCI- 451-52	SCI- 452-52	SCI- 453-52
54 mm	SCI- 450-54	SCI- 451-54	SCI- 452-54	SCI- 453-54
56 mm	SCI- 450-56	SCI- 451-56	SCI- 452-56	SCI- 453-56
58 mm	SCI- 450-58	SCI- 451-58	SCI- 452-58	SCI- 453-58
60 mm	SCI- 450-60	SCI- 451-60	SCI- 452-60	SCI- 453-60



SCI - SI : 5

Drummers Wire

Cat No.	Dia
SCI- 454-18	18 SWG (1.25 mm)
SCI- 454-20	20 SWG (0.90 mm)



SCI - SI : 6

Harrington Rod & Hooks


Cat No.	Description & Size
SCI-411-25	Harrington Rod : 2
SCI-411-51	Harrington Rod : 51 mm
SCI-411-76	Harrington Rod : 76 mm
SCI-411-102	Harrington Rod : 102 mm
SCI-411-127	Harrington Rod : 127 mm
SCI-411-152	Harrington Rod : 152 mm
SCI-411-178	Harrington Rod : 178 mm
SCI-411-203	Harrington Rod : 203 mm
SCI-411-224	Harrington Rod : 224 mm
SCI-411-254	Harrington Rod : 254 mm
SCI-412-01	Harrington Hook Upper : Sharp
SCI-412-02	Harrington Hook Upper : Blunt
SCI-412-03	Harrington Hook Lower : Sharp
SCI-412-04	Harrington Hook Lower : Blunt
SCI-413-01	Clip Washer
SCI-414-125	Sacral Rod With Six Nuts : 125 mm
SCI-414-175	Sacral Rod With Six Nuts : 175 mm
SCI-414-200	Sacral Rod With Six Nuts : 200 mm



SCI - WP : 1

Stainless Steel Suture Wire


Standard Lng. Roll	Size	10 Meter Roll
Catalogue No.		Catalogue No.
SCI-500-16	16 SWG (1.60 mm)	SCI-501-16
SCI-500-17	17 SWG (1.50 mm)	SCI-501-17
SCI-500-18	18 SWG (1.25 mm)	SCI-501-18
SCI-500-19	19 SWG (1.00 mm)	SCI-501-19
SCI-500-20	20 SWG (0.90 mm)	SCI-501-20
SCI-500-21	21 SWG (0.80 mm)	SCI-501-21
SCI-500-22	22 SWG (0.70 mm)	SCI-501-22
SCI-500-23	23 SWG (0.60 mm)	SCI-501-23
SCI-500-24	24 SWG (0.55 mm)	SCI-501-24
SCI-500-26	26 SWG (0.45 mm)	SCI-501-26
SCI-500-28	28 SWG (0.35 mm)	SCI-501-28
SCI-500-30	30 SWG (0.30 mm)	SCI-501-30



SCI - WP : 2

K-Wires

Diameter	Length				
	70 mm	100 mm	150 mm	225 mm	300 mm
	Cat No.	Cat No.	Cat No.	Cat No.	Cat No.
0.8 mm	SCI-502-08	SCI-503-08	SCI-504-08	-	-
0.9 mm	SCI-502-09	SCI-503-09	SCI-504-09	-	-
1.0 mm	-	SCI-503-10	SCI-504-10	SCI-505-10	-
1.1 mm	-	SCI-503-11	SCI-504-11	SCI-505-11	-
1.2 mm	-	SCI-503-12	SCI-504-12	SCI-505-12	-
1.4 mm	-		SCI-504-14	SCI-505-14	-
1.5 mm	-	-	SCI-504-15	SCI-505-15	SCI-506-15
1.6 mm	-	-	SCI-504-16	SCI-505-16	SCI-506-16
1.8 mm	-	-	SCI-504-18	SCI-505-18	SCI-506-18
2.0 mm	-	-	SCI-504-20	SCI-505-20	SCI-506-20
2.2 mm	-	-	SCI-504-22	SCI-505-22	SCI-506-22
2.3 mm	-	-	SCI-504-23	SCI-505-23	SCI-506-23
2.5 mm	-	-	SCI-504-25	SCI-505-25	SCI-506-25
3.0 mm	-	-	SCI-504-30	SCI-505-30	SCI-506-30



**Kirschner Wire with Threaded TIP
(Guide Wire)**



Length : 225 mm	Diameter	Length : 150 mm
Catalogue No.		Catalogue No.
SCI-507-15	1.5 mm	SCI-508-15
SCI-507-20	2.0 mm	SCI-508-20
SCI-507-25	2.5 mm	SCI-508-25




Steinmann Pin

Length	Dia				
	3.0mm	3.5mm	4.0mm	4.5mm	5.0mm
	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
150mm	SCI-509-15	SCI-510-15	SCI-511-15	SCI-512-15	SCI-513-15
175mm	SCI-509-17	SCI-510-17	SCI-511-17	SCI-512-17	SCI-513-17
200mm	SCI-509-20	SCI-510-20	SCI-511-20	SCI-512-20	SCI-513-20
225mm	SCI-509-22	SCI-510-22	SCI-511-22	SCI-512-22	SCI-513-22
250mm	SCI-509-25	SCI-510-25	SCI-511-25	SCI-512-25	SCI-513-25
275mm	SCI-509-27	SCI-510-27	SCI-511-27	SCI-512-27	SCI-513-27
300mm	SCI-509-30	SCI-510-30	SCI-511-30	SCI-512-30	SCI-513-30




Denham Pin (Steinmann Pin Centrally Threaded)

Cat No.	Size	Length
SCI-514-40	4.0mm	225mm
SCI-515-45	4.5mm	225mm
SCI-516-50	5.0mm	225mm
SCI-517-60	6.0mm	225mm




SCI - WP : 6

Cerclage Wire with Loop



Cat No. (Length 280mm)	Dia	Cat No. (Length 600mm)
SCI-518-17	17 SWG (1.5 mm)	SCI-519-17
SCI-518-18	18 SWG (1.25 mm)	SCI-519-18
SCI-518-19	19 SWG (1.0 mm)	SCI-519-19
SCI-518-20	20 SWG (0.9 mm)	SCI-519-20
SCI-518-21	21 SWG (0.8 mm)	SCI-519-21



SCI - WP : 7

Austin Moore Pin with 2 Nuts




SCI - WP : 7

Knowels Pin

Cat No.	Length	Cat No.	Length
SCI-520-50	50 mm	SCI-520-105	105 mm
SCI-520-55	55 mm	SCI-520-110	110 mm
SCI-520-60	60 mm	SCI-520-115	115 mm
SCI-520-65	65 mm	SCI-520-120	120 mm
SCI-520-70	70 mm	SCI-520-125	125 mm
SCI-520-75	75 mm	SCI-520-130	130 mm
SCI-520-80	80 mm	SCI-520-135	135 mm
SCI-520-85	85 mm	SCI-520-140	140 mm
SCI-520-90	90 mm	SCI-520-145	145 mm
SCI-520-95	95 mm	SCI-520-150	150 mm
SCI-520-100	100 mm		

3.0mm Dia Cat No.	Length	4.0 mm Dia Cat No.
SCI-521-65	65 mm	SCI-521-065
SCI-521-70	70 mm	SCI-522-070
SCI-521-75	75 mm	SCI-522-075
SCI-521-80	80 mm	SCI-522-080
SCI-521-85	85 mm	SCI-522-085
SCI-521-90	90 mm	SCI-522-090
SCI-521-95	95 mm	SCI-522-095
SCI-521-100	100 mm	SCI-522-0100



SCI - S : 1

Staples - U Type

Dia	Width		
	15mm	20mm	25mm
Cat No.			
1.5 mm	SCI-600-15	SCI-601-15	SCI-602-15
2.0 mm	SCI-600-20	SCI-601-20	SCI-602-20
2.5 mm	SCI-600-25	SCI-601-25	SCI-602-25
3.0 mm	SCI-600-30	SCI-601-30	SCI-602-30



SCI - S : 2

Staples - Convex

Cat No.	Dia
SCI-603-05	5 mm Step
SCI-603-10	10 mm Step
SCI-603-15	15 mm Step



SCI - S : 3

**6.5 Dia, Cancellous Traction Screw
Thread Length 32 mm**

Cat No.	Length	Cat No.	Length
SCI-604-90	90 mm	SCI-604-115	115 mm
SCI-604-95	95 mm	SCI-604-120	120 mm
SCI-604-100	100 mm	SCI-604-130	130 mm
SCI-604-105	105 mm	SCI-604-135	135 mm
SCI-604-110	110 mm		



SCI - N : 1



Kuntscher Cloverleaf Nail


Length	Diameter									
	6 mm	7 mm	8 mm	9 mm	10 mm	11 mm	12 mm	13 mm	14 mm	15 mm
	Cat No.	Cat No.	Cat No.	Cat No.	Cat No.	Cat No.	Cat No.	Cat No.	Cat No.	Cat No.
20 cm	SCI-706-20	SCI-707-20	-	-	-	-	-	-	-	-
21 cm	SCI-706-21	SCI-707-21	-	-	-	-	-	-	-	-
22 cm	SCI-706-22	SCI-707-22	-	-	-	-	-	-	-	-
23 cm	SCI-706-23	SCI-707-23	-	-	-	-	-	-	-	-
24 cm	SCI-706-24	SCI-707-24	SCI-708-24	SCI-709-24	-	-	-	-	-	-
25 cm	SCI-706-25	SCI-707-25	SCI-708-25	SCI-709-25	-	-	-	-	-	-
26 cm	SCI-706-26	SCI-707-26	SCI-708-26	SCI-709-26	SCI-710-26	SCI-711-26	SCI-712-26	SCI-713-26	SCI-714-26	SCI-715-26
27 cm	SCI-706-27	SCI-707-27	SCI-708-27	SCI-709-27	SCI-710-27	SCI-711-27	SCI-712-27	SCI-713-27	SCI-714-27	SCI-715-27
28 cm	SCI-706-28	SCI-707-28	SCI-708-28	SCI-709-28	SCI-710-28	SCI-711-28	SCI-712-28	SCI-713-28	SCI-714-28	SCI-715-28
29 cm	SCI-706-29	SCI-707-29	SCI-708-29	SCI-709-29	SCI-710-29	SCI-711-29	SCI-712-29	SCI-713-29	SCI-714-29	SCI-715-29
30 cm	SCI-706-30	SCI-707-30	SCI-708-30	SCI-709-30	SCI-710-30	SCI-711-30	SCI-712-30	SCI-713-30	SCI-714-30	SCI-715-30
31 cm	SCI-706-31	SCI-707-31	SCI-708-31	SCI-709-31	SCI-710-31	SCI-711-31	SCI-712-31	SCI-713-31	SCI-714-31	SCI-715-31
32 cm	SCI-706-32	SCI-707-32	SCI-708-32	SCI-709-32	SCI-710-32	SCI-711-32	SCI-712-32	SCI-713-32	SCI-714-32	SCI-715-32
33 cm	SCI-706-33	SCI-707-33	SCI-708-33	SCI-709-33	SCI-710-33	SCI-711-33	SCI-712-33	SCI-713-33	SCI-714-33	SCI-715-33
34 cm	SCI-706-34	SCI-707-34	SCI-708-34	SCI-709-34	SCI-710-34	SCI-711-34	SCI-712-34	SCI-713-34	SCI-714-34	SCI-715-34
35 cm	SCI-706-35	SCI-707-35	SCI-708-35	SCI-709-35	SCI-710-35	SCI-711-35	SCI-712-35	SCI-713-35	SCI-714-35	SCI-715-35
36 cm	SCI-706-36	SCI-707-36	SCI-708-36	SCI-709-36	SCI-710-36	SCI-711-36	SCI-712-36	SCI-713-36	SCI-714-36	SCI-715-36
37 cm	SCI-706-37	SCI-707-37	SCI-708-37	SCI-709-37	SCI-710-37	SCI-711-37	SCI-712-37	SCI-713-37	SCI-714-37	SCI-715-37
38 cm	SCI-706-38	SCI-707-38	SCI-708-38	SCI-709-38	SCI-710-38	SCI-711-38	SCI-712-38	SCI-713-38	SCI-714-38	SCI-715-38
39 cm	SCI-706-39	SCI-707-39	SCI-708-39	SCI-709-39	SCI-710-39	SCI-711-39	SCI-712-39	SCI-713-39	SCI-714-39	SCI-715-39
40 cm	SCI-706-40	SCI-707-40	SCI-708-40	SCI-709-40	SCI-710-40	SCI-711-40	SCI-712-40	SCI-713-40	SCI-714-40	SCI-715-40
41 cm	SCI-706-41	SCI-707-41	SCI-708-41	SCI-709-41	SCI-710-41	SCI-711-41	SCI-712-41	SCI-713-41	SCI-714-41	SCI-715-41
42 cm	SCI-706-42	SCI-707-42	SCI-708-42	SCI-709-42	SCI-710-42	SCI-711-42	SCI-712-42	SCI-713-42	SCI-714-42	SCI-715-42
43 cm	SCI-706-43	SCI-707-43	SCI-708-43	SCI-709-43	SCI-710-43	SCI-711-43	SCI-712-43	SCI-713-43	SCI-714-43	SCI-715-43
44 cm	SCI-706-44	SCI-707-44	SCI-708-44	SCI-709-44	SCI-710-44	SCI-711-44	SCI-712-44	SCI-713-44	SCI-714-44	SCI-715-44
45 cm	SCI-706-45	SCI-707-45	SCI-708-45	SCI-709-45	SCI-710-45	SCI-711-45	SCI-712-45	SCI-713-45	SCI-714-45	SCI-715-45
46 cm	SCI-706-46	SCI-707-46	SCI-708-46	SCI-709-46	SCI-710-46	SCI-711-46	SCI-712-46	SCI-713-46	SCI-714-46	SCI-715-46
47 cm	-	-	-	-	SCI-710-47	SCI-711-47	SCI-712-47	SCI-713-47	SCI-714-47	SCI-715-47
48 cm	-	-	-	-	SCI-710-48	SCI-711-48	SCI-712-48	SCI-713-48	SCI-714-48	SCI-715-48



SCI - N : 2

Bend Kuntscher Cloverleaf Nail

Length	Diameter			
	7 mm Cat. No.	8 mm Cat. No.	9 mm Cat. No.	10 mm Cat. No.
26 cm	SCI-727-26	SCI-728-26	SCI-729-26	SCI-730-26
27 cm	SCI-727-27	SCI-728-27	SCI-729-27	SCI-730-27
28 cm	SCI-727-28	SCI-728-28	SCI-729-28	SCI-730-28
29 cm	SCI-727-29	SCI-728-29	SCI-729-29	SCI-730-29
30 cm	SCI-727-30	SCI-728-30	SCI-729-30	SCI-730-30
31 cm	SCI-727-31	SCI-728-31	SCI-729-31	SCI-730-31
32 cm	SCI-727-32	SCI-728-32	SCI-729-32	SCI-730-32
33 cm	SCI-727-33	SCI-728-33	SCI-729-33	SCI-730-33
34 cm	SCI-727-34	SCI-728-34	SCI-729-34	SCI-730-34
35 cm	SCI-727-35	SCI-728-35	SCI-729-35	SCI-730-35
36 cm	SCI-727-36	SCI-728-36	SCI-729-36	SCI-730-36
37 cm	SCI-727-37	SCI-728-37	SCI-729-37	SCI-730-37
38 cm	SCI-727-38	SCI-728-38	SCI-729-38	SCI-730-38
39 cm	SCI-727-39	SCI-728-39	SCI-729-39	SCI-730-39
40 cm	SCI-727-40	SCI-728-40	SCI-729-40	SCI-730-40
41 cm	SCI-727-41	SCI-728-41	SCI-729-41	SCI-730-41
42 cm	SCI-727-42	SCI-728-42	SCI-729-42	SCI-730-42
43 cm	SCI-727-43	SCI-728-43	SCI-729-43	SCI-730-43
44 cm	SCI-727-44	SCI-728-44	SCI-729-44	SCI-730-44



SCI - N : 3

Square Nail for Radius

Dia.	Dia.	Dia.	Dia.	Dia.	Dia.	Length
1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0mm	
Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	
SCI-735-17	SCI-736-17	SCI-737-17	-	-	-	17 cm
SCI-735-18	SCI-736-18	SCI-737-18	-	-	-	18 cm
SCI-735-19	SCI-736-19	SCI-737-19	-	-	-	19 cm
SCI-735-20	SCI-736-20	SCI-737-20	SCI-738-20	SCI-739-20	SCI-740-20	20 cm
SCI-735-21	SCI-736-21	SCI-737-21	SCI-738-21	SCI-739-21	SCI-740-21	21 cm
SCI-735-22	SCI-736-22	SCI-737-22	SCI-738-22	SCI-739-22	SCI-740-22	22 cm
SCI-735-23	SCI-736-23	SCI-737-23	SCI-738-23	SCI-739-23	SCI-740-23	23 cm
SCI-735-24	SCI-736-24	SCI-737-24	SCI-738-24	SCI-739-24	SCI-740-24	24 cm
SCI-735-25	SCI-736-25	SCI-737-25	SCI-738-25	SCI-739-25	SCI-740-25	25 cm
SCI-735-26	SCI-736-26	SCI-737-26	SCI-738-26	SCI-739-26	SCI-740-26	26 cm
SCI-735-27	SCI-736-27	SCI-737-27	SCI-738-27	SCI-739-27	SCI-740-27	27 cm
SCI-735-28	SCI-736-28	SCI-737-28	SCI-738-28	SCI-739-28	SCI-740-28	28 cm
SCI-735-29	SCI-736-29	SCI-737-29	SCI-738-29	SCI-739-29	SCI-740-29	29 cm
SCI-735-30	SCI-736-30	SCI-737-30	SCI-738-30	SCI-739-30	SCI-740-30	30 cm



SCI - N : 4

Square Nail for Ulna

Dia.	Dia.	Dia.	Dia.	Dia.	Dia.	Length
1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0mm	
Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	
SCI-745-17	SCI-746-17	SCI-747-17	-	-	-	17 cm
SCI-745-18	SCI-746-18	SCI-747-18	-	-	-	18 cm
SCI-745-19	SCI-746-19	SCI-747-19	-	-	-	19 cm
SCI-745-20	SCI-746-20	SCI-747-20	SCI-748-20	SCI-749-20	SCI-750-20	20 cm
SCI-745-21	SCI-746-21	SCI-747-21	SCI-748-21	SCI-749-21	SCI-750-21	21 cm
SCI-745-22	SCI-746-22	SCI-747-22	SCI-748-22	SCI-749-22	SCI-750-22	22 cm
SCI-745-23	SCI-746-23	SCI-747-23	SCI-748-23	SCI-749-23	SCI-750-23	23 cm
SCI-745-24	SCI-746-24	SCI-747-24	SCI-748-24	SCI-749-24	SCI-750-24	24 cm
SCI-745-25	SCI-746-25	SCI-747-25	SCI-748-25	SCI-749-25	SCI-750-25	25 cm
SCI-745-26	SCI-746-26	SCI-747-26	SCI-748-26	SCI-749-26	SCI-750-26	26 cm
SCI-745-27	SCI-746-27	SCI-747-27	SCI-748-27	SCI-749-27	SCI-750-27	27 cm
SCI-745-28	SCI-746-28	SCI-747-28	SCI-748-28	SCI-749-28	SCI-750-28	28 cm
SCI-745-29	SCI-746-29	SCI-747-29	SCI-748-29	SCI-749-29	SCI-750-29	29 cm
SCI-745-30	SCI-746-30	SCI-747-30	SCI-748-30	SCI-749-30	SCI-750-30	30 cm



SCI - N : 5

Rush Nail

Length	Diameter						
	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm	5.0 mm	6.0 mm
	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
5 cm	SCI-755-05	-	-	-	-	-	-
6 cm	SCI-755-06	-	-	-	-	-	-
7 cm	SCI-755-07	-	-	-	-	-	-
8 cm	SCI-755-08	-	-	-	-	-	-
9 cm	SCI-755-09	-	-	-	-	-	-
10 cm	SCI-755-10	SCI-756-10	SCI-757-10	SCI-758-10	-	-	-
11 cm	SCI-755-11	SCI-756-11	SCI-757-11	SCI-758-11	-	-	-
12 cm	SCI-755-12	SCI-756-12	SCI-757-12	SCI-758-12	-	-	-
13 cm	SCI-755-13	SCI-756-13	SCI-757-13	SCI-758-13	-	-	-
14 cm	SCI-755-14	SCI-756-14	SCI-757-14	SCI-758-14	-	-	-
15 cm	SCI-755-15	SCI-756-15	SCI-757-15	SCI-758-15	SCI-759-15	-	-
16 cm	SCI-755-16	SCI-756-16	SCI-757-16	SCI-758-16	SCI-759-16	-	-
17 cm	SCI-755-17	SCI-756-17	SCI-757-17	SCI-758-17	SCI-759-17	-	-
18 cm	SCI-755-18	SCI-756-18	SCI-757-18	SCI-758-18	SCI-759-18	-	-
19 cm	SCI-755-19	SCI-756-19	SCI-757-19	SCI-758-19	SCI-759-19	-	-
20 cm	SCI-755-20	SCI-756-20	SCI-757-20	SCI-758-20	SCI-759-20	SCI-760-20	SCI-761-20
21 cm	SCI-755-21	SCI-756-21	SCI-757-21	SCI-758-21	SCI-759-21	SCI-760-21	SCI-761-21
22 cm	SCI-755-22	SCI-756-22	SCI-757-22	SCI-758-22	SCI-759-22	SCI-760-22	SCI-761-22
23 cm	SCI-755-23	SCI-756-23	SCI-757-23	SCI-758-23	SCI-759-23	SCI-760-23	SCI-761-23
24 cm	SCI-755-24	SCI-756-24	SCI-757-24	SCI-758-24	SCI-759-24	SCI-760-24	SCI-761-24
25 cm	SCI-755-25	SCI-756-25	SCI-757-25	SCI-758-25	SCI-759-25	SCI-760-25	SCI-761-25
26 cm	-	SCI-756-26	SCI-757-26	SCI-758-26	SCI-759-26	SCI-760-26	SCI-761-26
27 cm	-	SCI-756-27	SCI-757-27	SCI-758-27	SCI-759-27	SCI-760-27	SCI-761-27
28 cm	-	SCI-756-28	SCI-757-28	SCI-758-28	SCI-759-28	SCI-760-28	SCI-761-28
29 cm	-	SCI-756-29	SCI-757-29	SCI-758-29	SCI-759-29	SCI-760-29	SCI-761-29
30 cm	-	SCI-756-30	SCI-757-30	SCI-758-30	SCI-759-30	SCI-760-30	SCI-761-30
31 cm	-	-	-	-	SCI-759-31	SCI-760-31	SCI-761-31
32 cm	-	-	-	-	SCI-759-32	SCI-760-32	SCI-761-32
33 cm	-	-	-	-	SCI-759-33	SCI-760-33	SCI-761-33
34 cm	-	-	-	-	SCI-759-34	SCI-760-34	SCI-761-34
35 cm	-	-	-	-	SCI-759-35	SCI-760-35	SCI-761-35
36 cm	-	-	-	-	SCI-759-36	SCI-760-36	SCI-761-36
37 cm	-	-	-	-	SCI-759-37	SCI-760-37	SCI-761-37
38 cm	-	-	-	-	SCI-759-38	SCI-760-38	SCI-761-38
39 cm	-	-	-	-	SCI-759-39	SCI-760-39	SCI-761-39
40 cm	-	-	-	-	SCI-759-40	SCI-760-40	SCI-761-40



SCI - N : 6

Profile

Proximal Shaft
Distal Shaft
Hole for 4.9 mm
Locking Bolt

Two round holes
Three round holes
5.0 mm

Interlocking Nail (Universal Nail Femur)

Length	Diameter			
	9 mm	10 mm	11 mm	12 mm
	Cat. No.	Cat. No.	Cat. No.	Cat. No.
26 cm	-	-	-	-
28 cm	-	-	-	-
30 cm	SCI-765-30	SCI-766-30	SCI-767-30	SCI-768-30
32 cm	SCI-765-32	SCI-766-32	SCI-767-32	SCI-768-32
34 cm	SCI-765-34	SCI-766-34	SCI-767-34	SCI-768-34
36 cm	SCI-765-36	SCI-766-36	SCI-767-36	SCI-768-36
38 cm	SCI-765-38	SCI-766-38	SCI-767-38	SCI-768-38
40 cm	SCI-765-40	SCI-766-40	SCI-767-40	SCI-768-40
42 cm	SCI-765-42	SCI-766-42	SCI-767-42	SCI-768-42
44 cm	SCI-765-44	SCI-766-44	SCI-767-44	SCI-768-44



SCI - N : 7

Profile

Proximal Shaft

Distal Shaft
Hole for 4.9 mm
Locking Bolt

One round and one
compression hole
Three round holes
5.0 mm

Interlocking Nail (Universal Nail Tibia)

Length	Diameter			
	8 mm	09 mm	10 mm	11 mm
	Cat. No.	Cat. No.	Cat. No.	Cat. No.
26 cm	SCI-775-26	SCI-776-26	SCI-777-26	SCI-778-26
28 cm	SCI-775-28	SCI-776-28	SCI-777-28	SCI-778-28
30 cm	SCI-775-30	SCI-776-30	SCI-777-30	SCI-778-30
32 cm	SCI-775-32	SCI-776-32	SCI-777-32	SCI-778-32
34 cm	SCI-775-34	SCI-776-34	SCI-777-34	SCI-778-34
36 cm	SCI-775-36	SCI-776-36	SCI-777-36	SCI-778-36
38 cm	SCI-775-38	SCI-776-38	SCI-777-38	SCI-778-38
40 cm	SCI-775-40	SCI-776-40	SCI-777-40	SCI-778-40



SCI - N : 8

Supra Condylar Nail (IMSC) - Five Holes

Profile

Proximal Shaft Dia 13.0 mm
Proximal Holes for 6.5 mm Screws
Distal Hole for 4.9 mm Locking Bolts 5.0 mm

Length	Diameter			
	10 mm	11 mm	12 mm	13 mm
	Cat. No.	Cat. No.	Cat. No.	Cat. No.
15 cm	SCI-780-15	SCI-781-15	SCI-782-15	SCI-783-15
20 cm	SCI-780-20	SCI-781-20	SCI-782-20	SCI-783-20
25 cm	SCI-780-25	SCI-781-25	SCI-782-25	SCI-783-25
30 cm	SCI-780-30	SCI-781-30	SCI-782-30	SCI-783-30




SCI - N : 9

Supra Condylar Nail (IMSC) - Multiple Holes

Profile

Proximal Shaft Dia 13.0 mm
Proximal Holes for 6.5 mm Screws
Distal Hole for 4.9 mm Locking Bolts 5.0 mm

Length	Diameter			
	10 mm	11 mm	12 mm	13 mm
	Cat. No.	Cat. No.	Cat. No.	Cat. No.
15 cm	SCI-785-15	SCI-786-15	SCI-787-15	SCI-788-15
20 cm	SCI-785-20	SCI-786-20	SCI-787-20	SCI-788-20
25 cm	SCI-785-25	SCI-786-25	SCI-787-25	SCI-788-25
30 cm	SCI-785-30	SCI-786-30	SCI-787-30	SCI-788-30




SCI - N : 10

Profile
Proximal Shaft Dia 8.0 mm
Hole Dia for 3.4 mm Locking Bolts 3.5 mm

Un-Reamed Humeral Nail (UHN)

Length	Diameter	
	6.7 mm	7.5 mm
	Catalogue No.	
19 cm	SCI-790-19	SCI-791-19
20 cm	SCI-790-20	SCI-791-20
21 cm	SCI-790-21	SCI-791-21
22 cm	SCI-790-22	SCI-791-22
23 cm	SCI-790-23	SCI-791-23
24 cm	SCI-790-24	SCI-791-24
25 cm	SCI-790-25	SCI-791-25
26 cm	SCI-790-26	SCI-791-26
27 cm	SCI-790-27	SCI-791-27
28 cm	SCI-790-28	SCI-791-28
29 cm	SCI-790-29	SCI-791-29
30 cm	SCI-790-30	SCI-791-30
31 cm	SCI-790-31	SCI-791-31
32 cm	SCI-790-32	SCI-791-32



SCI - N : 11

Reconstruction Nail (Recon Nail)

Length	Left Leg					Right Leg				
	Diameter					Diameter				
	9 mm	10 mm	11 mm	12 mm	13 mm	9 mm	10 mm	11 mm	12 mm	13 mm
	Catalogue No.					Catalogue No.				
34 cm	SCI-795-34	SCI-796-34	SCI-797-34	SCI-798-34	SCI-799-34	SCI-800-34	SCI-801-34	SCI-802-34	SCI-803-34	SCI-804-34
36 cm	SCI-795-36	SCI-796-36	SCI-797-36	SCI-798-36	SCI-799-36	SCI-800-36	SCI-801-36	SCI-802-36	SCI-803-36	SCI-804-36
38 cm	SCI-795-38	SCI-796-38	SCI-797-38	SCI-798-38	SCI-799-38	SCI-800-38	SCI-801-38	SCI-802-38	SCI-803-38	SCI-804-38
40 cm	SCI-795-40	SCI-796-40	SCI-797-40	SCI-798-40	SCI-799-40	SCI-800-40	SCI-801-40	SCI-802-40	SCI-803-40	SCI-804-40
42 cm	SCI 795 42	SCI 796 42	SCI 797 42	SCI 798 42	SCI 799 42	SCI 800 42	SCI 801 42	SCI-802-42	SCI-803-42	SCI-804-42
44 cm	SCI-795-44	SCI-796-44	SCI-797-44	SCI-798-44	SCI-799-44	SCI-800-44	SCI-801-44	SCI-802-44	SCI-803-44	SCI-804-44
46 cm	SCI-795-46	SCI-796-46	SCI-797-46	SCI-798-46	SCI-799-46	SCI-800-46	SCI-801-46	SCI-802-46	SCI-803-46	SCI-804-46



SCI - N : 12

Proximal Femoral Nail (PFN)

Length	Diameter			
	9 mm	10 mm	11 mm	12 mm
	Catalogue No.			
22 cm	SCI- 809-22	SCI- 810-22	SCI- 811-22	SCI- 812-22
24 cm	SCI- 809-24	SCI- 810-24	SCI- 811-24	SCI- 812-24
26 cm	SCI- 809-26	SCI- 810-26	SCI- 811-26	SCI- 812-26



INTERNAL FIXATION
SYSTEMS™

Locking Screws



SCI - LS : 1

Locking Bolt 6.5 mm Self Tapping

Solid Catalogue No.	Length	Cannulated Catalogue No.
SCI-318-50	50 mm	SCI-319-50
SCI-318-55	55 mm	SCI-319-55
SCI-318-60	60 mm	SCI-319-60
SCI-318-65	65 mm	SCI-319-65
SCI-318-70	70 mm	SCI-319-70
SCI-318-75	75 mm	SCI-319-75
SCI-318-80	80 mm	SCI-319-80
SCI-318-85	85 mm	SCI-319-85
SCI-318-90	90 mm	SCI-319-90
SCI-318-95	95 mm	SCI-319-95
SCI-318-100	100 mm	SCI-319-100
SCI-318-105	105 mm	SCI-319-105
SCI-318-110	110 mm	SCI-319-110
SCI-318-115	115 mm	SCI-319-115
SCI-318-120	120 mm	SCI-319-120



SCI - LS : 2

Locking Bolt 7.5 mm Self Tapping

Solid Catalogue No.	Length	Cannulated Catalogue No.
SCI-320-50	50 mm	SCI-321-050
SCI-320-55	55 mm	SCI-321-055
SCI-320-60	60 mm	SCI-321-060
SCI-320-65	65 mm	SCI-321-065
SCI-320-70	70 mm	SCI-321-070
SCI-320-75	75 mm	SCI-321-075
SCI-320-80	80 mm	SCI-321-080
SCI-320-85	85 mm	SCI-321-085
SCI-320-90	90 mm	SCI-321-090
SCI-320-95	95 mm	SCI-321-095
SCI-320-100	100 mm	SCI-321-0100
SCI-320-105	105 mm	SCI-321-0105
SCI-320-110	110 mm	SCI-321-0110
SCI-320-115	115 mm	SCI-321-0115
SCI-320-120	120 mm	SCI-321-0120



SCI - LS : 3

Locking Bolt 4.9 mm Self Tapping

Profile	
Thread Diameter	4.9 mm
Fully Threaded	
Drill Bit for Threaded Hole	4.0 mm
Core Diameter	4.3 mm
Hexagonal Socket	3.5 mm
Head Diameter	8.0 mm



SCI - LS : 4

Locking Bolt 3.9 mm Self Tapping

Profile	
Thread Diameter	3.9 mm
Fully Threaded	
Drill Bit for Threaded Hole	3.2 mm
Core Diameter	3.2 mm
Hexagonal Socket	2.5 mm
Head Diameter	8.0 mm

Locking Bolt 4.9 mm Self Tapping

Cat No.	Length	Cat No.	Length	Cat No.	Length
SCI-321-24	24 mm	SCI-321-42	42 mm	SCI-321-60	60 mm
SCI-321-26	26 mm	SCI-321-44	44 mm	SCI-321-65	65 mm
SCI-321-28	28 mm	SCI-321-46	46 mm	SCI-321-70	70 mm
SCI-321-30	30 mm	SCI-321-48	48 mm	SCI-321-75	75 mm
SCI-321-32	32 mm	SCI-321-50	50 mm	SCI-321-80	80 mm
SCI-321-34	34 mm	SCI-321-52	52 mm	SCI-321-85	85 mm
SCI-321-36	36 mm	SCI-321-54	54 mm	SCI-321-90	90 mm
SCI-321-38	38 mm	SCI-321-56	56 mm	SCI-321-95	95 mm
SCI-321-40	40 mm	SCI-321-58	58 mm	SCI-321-100	100mm

Locking Bolt 3.9 mm Self Tapping

Cat No.	Length	Cat No.	Length
SCI-322-20	20 mm	SCI-322-42	42 mm
SCI-322-22	22 mm	SCI-322-44	44 mm
SCI-322-24	24 mm	SCI-322-46	46 mm
SCI-322-26	26 mm	SCI-322-48	48 mm
SCI-322-28	28 mm	SCI-322-50	50 mm
SCI-322-30	30 mm	SCI-322-52	52 mm
SCI-322-32	32 mm	SCI-322-54	54 mm
SCI-322-34	34 mm	SCI-322-56	56 mm
SCI-322-36	36 mm	SCI-322-58	58 mm
SCI-322-38	38 mm	SCI-322-60	60 mm
SCI-322-40	40 mm		



End Caps for Un-reamed Humeral Nail (UHN)

Extension	Outer Dia.	
	5.0mm	8.0mm
0 mm	SCI-324-00	SCI-325-00
5 mm	SCI-324-05	SCI-325-05
10 mm	SCI-324-10	SCI-325-10
15 mm	SCI-324-15	SCI-325-15

Locking Bolt 3.4 mm Self Tapping

Cat No.	Length	Cat No.	Length
SCI-323-14	14 mm	SCI-323-30	30 mm
SCI-323-16	16 mm	SCI-323-32	32 mm
SCI-323-18	18 mm	SCI-323-34	34 mm
SCI-323-20	20 mm	SCI-323-36	36 mm
SCI-323-22	22 mm	SCI-323-38	38 mm
SCI-323-24	24 mm	SCI-323-40	40 mm
SCI-323-26	26 mm	SCI-323-45	45 mm
SCI-323-28	28 mm	SCI-323-50	50 mm



SCI - LS : 6

**Austin Moore Hip Prosthesis:-
Excel, Standard & Narrow Stem**

Diameter	Standard Stem (Non-Sterile)	Standard Stem (Sterile)	Narrow Stem (Non-Sterile)	Narrow Stem (Sterile)
	Cat No.	Cat No.	Cat No.	Cat No.
38 mm	SCI-900-38	SCI-901-38	SCI-902-38	SCI-903-38
39 mm	SCI-900-39	SCI-901-39	SCI-902-39	SCI-903-39
40 mm	SCI-900-40	SCI-901-40	SCI-902-40	SCI-903-40
41 mm	SCI-900-41	SCI-901-41	SCI-902-41	SCI-903-41
42 mm	SCI-900-42	SCI-901-42	SCI-902-42	SCI-903-42
43 mm	SCI-900-43	SCI-901-43	SCI-902-43	SCI-903-43
44 mm	SCI-900-44	SCI-901-44	SCI-902-44	SCI-903-44
45 mm	SCI-900-45	SCI-901-45	SCI-902-45	SCI-903-45
46 mm	SCI-900-46	SCI-901-46	SCI-902-46	SCI-903-46
47 mm	SCI-900-47	SCI-901-47	SCI-902-47	SCI-903-47
48 mm	SCI-900-48	SCI-901-48	SCI-902-48	SCI-903-48
49 mm	SCI-900-49	SCI-901-49	SCI-902-49	SCI-903-49
50 mm	SCI-900-50	SCI-901-50	SCI-902-50	SCI-903-50
51 mm	SCI-900-51	SCI-901-51	SCI-902-51	SCI-903-51
52 mm	SCI-900-52	SCI-901-52	SCI-902-52	SCI-903-52
53 mm	SCI-900-53	SCI-901-53	SCI-902-53	SCI-903-53
54 mm	SCI-900-54	SCI-901-54	SCI-902-54	SCI-903-54
55 mm	SCI-900-55	SCI-901-55	SCI-902-55	SCI-903-55
56 mm	SCI-900-56	SCI-901-56	SCI-902-56	SCI-903-56
57 mm	SCI-900-57	SCI-901-57	SCI-902-57	SCI-903-57
58 mm	SCI-900-58	SCI-901-58	SCI-902-58	SCI-903-58



SCI - LS : 7

**Austin Moore Hip Prosthesis:-
Excel - Broad Stem**



SCI - LS : 8

**Austin Moore Hip Prosthesis:-
Excel - Extra Long Stem
Stem Length - 300 mm**


Non Sterile	Dia.	Sterile
Cat No.		Cat. No.
SCI-905-37	37 mm	SCI-906-37
SCI-905-38	38 mm	SCI-906-38
SCI-905-39	39 mm	SCI-906-39
SCI-905-40	40 mm	SCI-906-40
SCI-905-41	41 mm	SCI-906-41
SCI-905-42	42 mm	SCI-906-42
SCI-905-43	43 mm	SCI-906-43
SCI-905-44	44 mm	SCI-906-44
SCI-905-45	45 mm	SCI-906-45
SCI-905-46	46 mm	SCI-906-46
SCI-905-47	47 mm	SCI-906-47
SCI-905-48	48 mm	SCI-906-48
SCI-905-49	49 mm	SCI-906-49
SCI-905-50	50 mm	SCI-906-50
SCI-905-51	51 mm	SCI-906-51
SCI-905-52	52 mm	SCI-906-52
SCI-905-53	53 mm	SCI-906-53
SCI-905-54	54 mm	SCI-906-54
SCI-905-55	55 mm	SCI-906-55
SCI-905-56	56 mm	SCI-906-56
SCI-905-57	57 mm	SCI-906-57
SCI-905-58	58 mm	SCI-906-58

Cat No.	Dia.
SCI-907-37	37 mm
SCI-907-38	38 mm
SCI-907-39	39 mm
SCI-907-40	40 mm
SCI-907-41	41 mm
SCI-907-42	42 mm
SCI-907-43	43 mm
SCI-907-44	44 mm
SCI-907-45	45 mm
SCI-907-46	46 mm
SCI-907-47	47 mm
SCI-907-48	48 mm
SCI-907-49	49 mm
SCI-907-50	50 mm
SCI-907-51	51 mm
SCI-907-52	52 mm
SCI-907-53	53 mm
SCI-907-54	54 mm
SCI-907-55	55 mm
SCI-907-56	56 mm
SCI-907-57	57 mm
SCI-907-58	58 mm



**Thompson Hip Prosthesis - Excel
Standard & Narrow Stem**

Diameter	Standard Stem (Non-Sterile)	Standard Stem (Sterile)	Narrow Stem (Non-Sterile)	Narrow Stem (Sterile)
	Cat No.	Cat No.	Cat No.	Cat No.
37 mm	SCI-908-37	SCI-909-37	SCI-910-37	SCI-911-37
38 mm	SCI-908-38	SCI-909-38	SCI-910-38	SCI-911-38
39 mm	SCI-908-39	SCI-909-39	SCI-910-39	SCI-911-39
40 mm	SCI-908-40	SCI-909-40	SCI-910-40	SCI-911-40
41 mm	SCI-908-41	SCI-909-41	SCI-910-41	SCI-911-41
42 mm	SCI-908-42	SCI-909-42	SCI-910-42	SCI-911-42
43 mm	SCI-908-43	SCI-909-43	SCI-910-43	SCI-911-43
44 mm	SCI-908-44	SCI-909-44	SCI-910-44	SCI-911-44
45 mm	SCI-908-45	SCI-909-45	SCI-910-45	SCI-911-45
46 mm	SCI-908-46	SCI-909-46	SCI-910-46	SCI-911-46
47 mm	SCI-908-47	SCI-909-47	SCI-910-47	SCI-911-47
48 mm	SCI-908-48	SCI-909-48	SCI-910-48	SCI-911-48
49 mm	SCI-908-49	SCI-909-49	SCI-910-49	SCI-911-49
50 mm	SCI-908-50	SCI-909-50	SCI-910-50	SCI-911-50
51 mm	SCI-908-51	SCI-909-51	SCI-910-51	SCI-911-51
52 mm	SCI-908-52	SCI-909-52	SCI-910-52	SCI-911-52
53 mm	SCI-908-53	SCI-909-53	SCI-910-53	SCI-911-53
54 mm	SCI-908-54	SCI-909-54	SCI-910-54	SCI-911-54
55 mm	SCI-908-55	SCI-909-55	SCI-910-55	SCI-911-55
56 mm	SCI-908-56	SCI-909-56	SCI-910-56	SCI-911-56
57 mm	SCI-908-57	SCI-909-57	SCI-910-57	SCI-911-57
58 mm	SCI-908-58	SCI-909-58	SCI-910-58	SCI-911-58



SCI - LS : 10

**Bi-Polar Prosthesis - Talwakar Type -
Hard Top Non Fenestrated**

Profile

Stem Length 150 mm
Angle 135°
Inner Ball Dia 26mm for 39 to 43mm
 28mm for 45 to 53mm

Inner Material of Cup UHMWPE
Stainless Steel Grade 316L

Cat No.	Dia.
SCI-912-37	37 mm
SCI-912-38	38 mm
SCI-912-39	39 mm
SCI-912-40	40 mm
SCI-912-41	41 mm
SCI-912-42	42 mm
SCI-912-43	43 mm
SCI-912-44	44 mm
SCI-912-45	45 mm
SCI-912-46	46 mm
SCI-912-47	47 mm
SCI-912-48	48 mm
SCI-912-49	49 mm
SCI-912-50	50 mm
SCI-912-51	51 mm
SCI-912-52	52 mm
SCI-912-53	53 mm
SCI-912-54	54 mm
SCI-912-55	55 mm
SCI-912-56	56 mm
SCI-912-57	57 mm
SCI-912-58	58 mm



SCI - LS : 11

**Bi-Polar Prosthesis - Talwakar Type -
Hard Top Fenestrated**

Profile

Stem Length 150 mm
Angle 135°
Inner Ball Dia 26mm for 39 to 43mm
 28mm for 45 to 53mm

Inner Material of Cup UHMWPE
Stainless Steel Grade 316L

Cat No.	Dia.
SCI-913-37	37 mm
SCI-913-38	38 mm
SCI-913-39	39 mm
SCI-913-40	40 mm
SCI-913-41	41 mm
SCI-913-42	42 mm
SCI-913-43	43 mm
SCI-913-44	44 mm
SCI-913-45	45 mm
SCI-913-46	46 mm
SCI-913-47	47 mm
SCI-913-48	48 mm
SCI-913-49	49 mm
SCI-913-50	50 mm
SCI-913-51	51 mm
SCI-913-52	52 mm
SCI-913-53	53 mm
SCI-913-54	54 mm
SCI-913-55	55 mm
SCI-913-56	56 mm
SCI-913-57	57 mm
SCI-913-58	58 mm