



WUHU RUIJIN MEDICAL INSTRUMENT
& DEVICE CO.,LTD

Address/Tel/Fax: No.33, EastWanchun Road, Wuhu Economic
& Technology Development Zone, Wuhu Anhui China;
Tel: 86-553-2672511 Fax: 86-553-2672510
Website: www.whruijin.com



RJ-PD Micro high speed drill II ***Product performance and application***

Motor handpiece: Apply for ENT, Spine, Microsurgery.
Adapter handpiece: Apply for Neurosurgery, hand and
foot surgery, maxillofacial, joint surgery, etc

>> Feature

- 01 *Refined design*
- 02 *Imported moror, autoclavable*
- 03 *Low noise, high speed*
- 04 *Long service life.*




>> Main performance parameters

Sequence	Model No.:	Performance parameters
1	RJ-MZ-II	Power Console Input Power: AC 100-240V 600VA 50/60Hz Speed and frequency of electric grinding, drilling, milling, sawing, planing: i.Grinding and milling 0-40000r/min ii. Drilling 0-1000r/min iii. Saw 0-12000 times/min iv. Planing 0-5000r/min reciprocating rotation or 0-12000r/min unidirectional rotation v. Temperature rise: temperature rise for accessible surface is less than amibient temperature 20°C after non-loading operating for 3 minutes. vi. No-load noise: rotating ≤75db, non-rotating ≤90db vii. The blades of electric grinding, drilling, sawing and planing should be heat treated, and the hardness more than HRC30

Standard configuration >>

Power Console



Pedal switch



Straight attachment


Angled attachment

Motor handpiece A

Motor handpiece B

0-40000r/min

0-40000r/min



Cutting drill bit

Diamond drill bit

Cranial Mill Adapter

Oscillating Saw Adapter

Reciprocating Saw Adapter


Cranial Drill Adapter

0-40000r/min

0-12000r/min

0-12000r/min

0-1000r/min



Sterilization case for power tool

*Power tools excluded

