

Comprehensive solutions for healthcare providers

Product catalogue



tradition since

ZARYS is a family-owned company with Polish capital, operating on the market since 1989

Evolution from a small enterprise to a leader in the sourcing industry is a process that led to the creation of more comprehensive offers and solutions that facilitate the work of medics and save the lives of patients. ZARYS is currently hiring over 250 people, who are working to serve thousands of clients operating in the field of health care worldwide.

comprehensive medical supplies for hospitals

The company's customers have many product groups at their disposal, starting from medical gloves, through dressing materials, to operating kits and treatment packages. Products from different categories have individual characteristics and specific purposes. In order to facilitate the identification of materials and improve the ordering process, 13 different categories were created in addition to the latest category "ZARYS kids".



| | | W |
|---------------------------------|------------------|------------------------|
| surgical drapes and gowns | surgical sets | perso prote equi |

0



basic medical devices

medical

gloves

dressings

medical clothing

disposable devices

medical

onal ective pment patient hygiene

O omegapack

🛆 deltaset



Griolity

safe

BETAtex ALPHAtex

stila **• care**











sports & physiotherapy products



vitahealth





being a manufacturer of disposable medical equipment, ZARYS meets the highest required standards

Part of ZARYS production plants are located in Silesia and combine sterile production - a cleanroom with an area of over 600 m2 - a place for manufacturing dedicated operating sets, a sewing room and an automated production line.

Zarys logistic center

The ZARYS logistic center is the largest warehouse of disposable medical equipment in Central Europe. Located in Zabrze area as part of Katowice Special Economic Zone. Modern, high-bay warehouse, occupies an area of over 22,850 m2 that contains over 22,500 pallet spaces.

oro cucu

Topologicalwarehouse
vice Special
varehouse
at containsAugust
Descension</



ZARYS

TABLE OF

CONTENTS

WWW.ZARYS.COM

MEDICALGLOVES pp. 10-23

NITRILE GLOVES

15

27

36

- 13 | easyCARE nitrile 13 easyCARE nitrile long NEW
- easyCARE nitrile pink
- easvCARE nitrile black
- easvCARE nitrile Hvdra+
- mediCARE nitrile SafeDon
- 17 mediCARE nitrile SafeDon Compact

BASIC MEDICAL DEVICES pp. 24-85

BASIC INJECTION EQUIPMENT

- dispoFINE injection needle
- 27 dispoGUARD safety injection needle dispoSULIN insulin pen needle
- 28 28
- side hole needle
- dicoSPIKE, dicoSPIKE Chemo withdrawal cannula 29 31 duoNEX single-use syringe, 2-part, Luer
 - dicoNEX SN single-use syringe, 3-part, Luer,
- with needle alongside the syringe
- 32 dicoNEX single-use syringe, 3-part, Luer
- 32 dicoNEX single-use syringe, 3-part, Luer-Lock 33
- dicoNEX single-use syringe, 3-part, Luer-Lock + amber 33
- dicoNEX single-use catheter syringe, 3-part
- dicoTUBER tuberculin syringe 34 dicoSULIN insulin syringe

INFUSION THERAPY

- 36 | SANVIIIon I.V. cannula
- SANVIFION S I.V. cannula without port
- 37 easyFLOW IS PREMIUM infusion set NEW
- easvFLOW IS infusion set 37
- easyFLOW IS infusion set (dual chamber) NEW 38
- easyFLOW IS ECO infusion set 38
- 39 easyFLOW IS SAFE infusion set NEW
- easyFLOW IS AMBER infusion set, amber NEW 39 easvFLOW IS REG infusion set with precision flow
- 40 rate regulator NEW
- 40 precision flow regulator NEW
- easyFLOW TS PREMIUM transfusion set NEW 41
- easyFLOW TS transfusion set 41
- 42
- easyFLOW LINE extension tube for infusion pump
- 42 needle-free valve
- 43 needle free valve with extension line NEW
- 43 3-way stopcock
- 44 3-way stopcock with extension tube
- 44 | cannula stopper

OXYGEN THERAPY

- 46 | nasal oxygen cannula
- 46 oxygen tubing
- 46 oxygen mask
- 47 non-rebreather mask
- 47 multi-vent mask
- 47 nebulizer mask
- 48 nebulizer

gauze and non-woven dressings

- 89 GAZA lux dressing gauze, zig-zag folded
- 89 GAZA lux cutting gauze
- GAZA lux S cutting gauze 89
- 90 KOMPRI lux gauze swabs BESTSELLER
- 91 KOMPRI lux S gauze swabs
- 92 SERVI lux gauze lap sponges
- SERVI lux S gauze lap sponges 92
 - TUPFER lux gauze balls
- TUPFER lux S gauze balls

93

- 93 NONVI lux non-woven swabs
- 94 94 NONVI lux S non-woven swabs
- 94 NONVI lux S non-woven swabs, with Y-incision
- 95 MULTIabsorb ABD pads, multi-layer
- 95 96
- 96
- 96

- SENSIplast adhesive fabric tape, white
- 99 SOFTplast adhesive non-woven tape
- SILKplast adhesive silk tape 99
- 100 FILMplast adhesive PE tape
 - 100 POREplast adhesive perforated non-woven tape

plasters with dressing

- 101 | NONVIplast non-woven plaster with dressing

- - 101 **FABRIplast** fabric plaster with dressing

- MULTIabsorb S ABD pads, multi-layer
- elastoLUMENAL eye pads
- elastoLUMENALS eye pads
- cellPAD cellulose swabs NEW

adhesive tapes

- PLASTIplast adhesive fabric tape, skin colour 98



LATEX GLOVES AND VINYL GLOVES

19 | HIGH RISK NEW

20easyCARE latex PH21easyCARE vinyl PF easyCARE latex PF

20

50

50

50

51

51

52

52

52

53

53

54

54

55

55

56

64

64

65

65

66

66 67

67

68

URINE BAGS

elbow connector

STOMACH TUBES

Nelaton catheter

Tiemann catheter

urine bag hanger

urine bag

easyCARE latex

SURGICAL GLOVES 23 | SafeCARE powdered NEW

23 SafeCARE pf powder-free NEW



ANAESTHESIA PRODUCTS endotracheal tube uncuffed 70 endotracheal tube uncuffed NEW 70 Yankauer surgical suction cannula, bulb tip CH23 endotracheal tube cuffed 71 surgical suction tube endotracheal tube cuffed NEW 71 surgical suction setwith Yankauer reinforced endotracheal tube cuffed cannula, standard intubation stylet 72 surgical suction set with Yankauer tracheal tube introducer cannula, bulb tip endotracheal tube holder oropharyngeal airway/guedel airway **BASIC PRODUCTS** laryngeal mask PVC FOR GINECOLOGY laryngeal mask silicone 74 | vaginal speculum, cusco type tracheostomy tube cuffed 74 cervical brush, standard air cushion anaesthetic mask 75 cervical brush, special anaesthetic mask with open seal 75 umbilical cord clamp catheter mount with double swivel

56 dual-limb breathing circuit

SUCTION CATHETERS,

| 58 | suction catheter with vacuum control |
|----|--|
| 58 | suction catheter |
| 59 | suction catheter with frozen surface |
| 59 | feeding tube graduated, with X-ray contrast stomach tube/duodenal tube |
| 60 | stomach tube/duodenal tube |

UROLOGICAL CATHETERS,

| 62 | Foley catheter 2-way, silicone-coated latex |
|----|---|
| | Foley catheter 3-way, silicone-coated latex |
| 63 | Foley catheter 2-way with tiemann tip, |

silicone-coated latex Foley catheter 2-way, 100% silicone

- Foley catheter 3-way, 100% silicone
- Pezzer catheter silicone-coated latex samples taking urine bag
- urine bag with sample port
- 68 | urine bag with sample port 2W

GENERAL SURGERY

Yankauer surgical suction cannula, standard CH21

OTHER PRODUCTS

- 77 | ID wristband
- entrance mat 77
- 77 enema bag
- 78 enema bag
- 78 enema set
- 78 vomit baq
- 79 safeLANCE, pressure-activated safety lancet
- 79 dental brush with suction
- 79 oral swab
- 80 tactical tourniquet CAT
- 80 automatic tourniquet
- 81 latex-free tourniquet
- 81 tongue depressor
- 82 electronic thermometer
- 82 wash glove
- 83 PE gloves
- 83 PE oversleeves **NEW**
- 84 medical waste container
- 85 | hygienic pet underpad NEW



| SP | ECI | AL | PU | OSE |
|----|-----|-----|----|-----|
| DR | ESS | SIN | GS | |

self-adhesive non-woven dressings

| elastopor, elastopor E adhesive non-woven |
|---|
| dressing tape |
| elastopor STERIL non-woven dressing |
| with absorbent pad BESTSELLER |
| elastoKIDS STERIL non-woven dressing |
| with absorbent pad |
| elastopor IV non-woven IV cannula |
| dressing BESTSELLER |
| elastoKIDS IV non-woven IV cannula dressing |
| elastopor STERILD non-woven dressing |
| with absorbent pad, incision and O-hole |
| elastoKIDS STERILD non-woven dressing |
| with absorbent pad NEW |
| elastopor EYE, elastoKIDS EYE non-woven |
| eye dressing with absorbent pad |
| elastoKIDS colorful kids plasters |
| set of injection plasters NEW |
| |

stabilizing tapes

- 108 I elastoNASAL non-woven nasal cannula fixation plaste
- elastoNASAL KIDS non-woven nasal cannula 108 fixation plaster NEW
- 109 elastoFIXAL adhesive tape for catheter tubing fixation
- elastoFIXAL KIDS adhesive tape for catheter 109 tubing fixation NEW
- 109 elastoFIXAL Set adhesive tape for endotracheal tube
- elastoFIXAL Set KIDS adhesive tape 110
- for endotracheal tube NEW 110 elastoFIXALTUBE band for tracheostomy
- tube fixation
- 110 adhesive dressing pouch for dialysis catheters NEW

- 112 | elastoBAND BASIC knitted supporting bandage 112 elastoBAND FLEX elastic supporting bandage wove
- elastoBAND FIX elastic supporting 113
- bandage, woven **NEW** elastoBAND CREPP cohesive supporting 113
- bandage, non-woven NEW elastoBAND KNITT elastic supporting 114
- bandage knitted NEW 114
- 115 elastoNET elastic tubular net bandage

immobilizing dressings

- orthoSYNTHETIC orthopedic padding, synthetic
- 117 elastoSLING triangular sling
 - orthoBAND FAST P.O.P. bandages
- 118
- 118

133 | urinary bladder catheterization kit 2

dialysis on kit 1a

dialysis off kit 1b

- yellowBAND cohesive bandage

- 116 | orthoNATURAL orthopedic padding, natural
- 116
- 118

- orthoBAND NORMAL P.O.P. bandages

- orthoBAND LONG P.O.P. bandages
- ALKOPAD DUO gauze pad moistened with alcohol
 - 125 thermCARE rescue blanket, foil, termic

121 | paraffiNET paraffin gauze dressing

elastoDERM F foil dressing with frame

elastoDERM F-IV IV cannula dressing

elastoDERM IV IV cannula dressing

elastoDERM PAD foil dressing

124 | thermPAD cold/warm compress NEW

125 ALKOPAD gauze pad moistened

gynPAD gynaecological-maternity pad NEW

Alcohol Prep Pad gauze pad moistened

121 | elastoDERM foil dressing

with absorbent pad

with frame

with alcoho

with alcoho

134 | dialysis on/off kit 2

anasthesia kit 1

135 suture application kit 1

sterile circular banded covers

sterile rectangular banded covers

sterile adhesive fluid collection pouch

for medical devices

for medical devices

163 sterile armboard cover set

164 | sterile surgical tapes

SURGICAL SETS

universal set 1

universal set 2

universal set 5-G

universal set 7

universal set 8

universal set 5-GP

universal set 4 BESTSELLER

universal set 5 **BESTSELLER**

universal set 6 **BESTSELLER**

universal set 5-M NEW

universal set 6-M NEW

pediatric universal set 9 . universal set 10

universal sets

167 | basic set 1

sterile c-arm cover set

sterile adhesive pouches

sterile absorbent drape

sterile absorbent drape

165 sterile incise film elastoFILM

foil dressings

122

122

122

122

124

125

135

163

163

163

16/

164

164

165

167

167

167

168

168

169

169

170

170

171

171

172

172

accessories

orthopaedics and traumatology

general surgery

174

174

175

175

173 | abdominal set 1-IF

173 abdo-perineal set 1

laparoscopy set 1-IF

laparoscopy set 3

varicose vein set 1

175 varicose vein set 2

laparoscopy set 2 BESTSELLER

- 176 | extremity set 1
- 176 177 extremity set 2 arthroscopy set 1
- 177 arthroscopy set 2
- 178 arthroscopy set 2-IP
- 178 arthroscopy set 3
- 179 arthroscopy set 4-IP 179 arthroscopy set 5
- 180 hip set 1
- 180 hip set 2
- 181 hip set 3
- 181 shoulder set 1
- 182 shoulder set 2

MEDICAL GOWNS

non-woven

204

204

205

207

207

207

208

208

208

SURGICAL CLOTHING

CAPS, PROTECTIVE MASKS

DONA medical bouffant cap

LILI medical cap with rubber

SOL medical cap with ties

ABI medical cap with ties

POLA medical cap with rubber

- 182 dynamic hip screw set 1-IF NEW
- vertical isolation set 1-IP NEW 182

MEDICAL CLOTHING

201 | medical gown with elastic cuffs, non-woven

201 medical gown with knitted cuffs, non-woven

203 | surgical suit, short-sleeved shirt+trousers,

surgical shirt, short-sleeved, non-woven

surgical shirt, long-sleeved, non-woven

surgical trousers, with ties, non-woven

VITAHEALTH pp. 222-223

medical bouffant cap, non-woven

delivery set 4 189 delivery set 5

delivery set 1

delivery set 2

189 delivery set 3

gynaecology and urology

gynaecology set 3

uro/gynaecology set 2

uro/gynaecology set 3

185 tur set 2-IP **BESTSELLER**

c-section set 1-IF

c-section set 2-IF

c-section set 3-IF

cystoscopy set 1

183 | gynaecology set 1

tur set 1-IP

186 | c-section set 1

183

184

184

184

185

185

187

187

188

188

188

189

obstetrics

ophthalmology

190 | ophthalmic set 2-IF 190 ophthalmic set 3-IF

pp. 198-221

209 | NINO medical cap

BEN medical cap with ties TOM medical cap with ties

ROS medical cap with ties

medical mask with ties

PATIENT CLOTHING

209

209

210

210

210

210

211

211

211

211

213

213

214

21/

214

215

6 Π∄

> 80 ШÓ

pp. 126-135

STANDARD DELTASET SURGICAL KITS

131 | puncture kit 1

DELTASET

- 132 suture removal kit 1
- 132 urinary bladder catheterization kit 1

OMEGAPACK pp. 136-139

PROCEDURE SETS

138 advantages of using omegapack operating sets139 components of our set

PERSONAL PROTECTIVE EQUIPMENTS pp. 140-141

133

134

BETASAFE

141 | protective coverall SAFE+

SURGICAL SETS, DRAPES, GOWNS pp. 142-197

SURGICAL DRAPES

2-layer surgical drapes

- 151 | sterile surgical drapes with cellulose layer
- sterile surgical drapes (economic version)
- 151 sterile surgical drapes 151 sterile surgical drapes with adhesive edge
- (economic version) sterile surgical drapes with adhesive edge 152
- sterile surgical drapes with central fenestration 152 (economic version)
- 152 sterile surgical drapes with central fenestration sterile surgical drapes with adhesive central 152
- fenestration (economic version) 152 sterile surgical drapes with adhesive central fenestration
- sterile surgical drapes with adhesive central 153 adjustable fenestration (economic version)
- sterile surgical drapes with adhesive central 153 adiustable fenestration sterile adhesive epidural drapes 153
- 153 | sterile instrument table covers

3-layer surgical drapes

- 154 | sterile surgical drapes
- sterile surgical drapes with adhesive edge 154
- sterile surgical drapes with central fenestration 154
- sterile surgical drapes with adhesive central 155 fenestration
- 155 sterile surgical drapes with adhesive fenestration, reinforced

specialistic surgical drapes 156 | abdominal drape 1-IF

orthopaedic drape 2

orthopaedic drape 3

orthopaedic drape 4

extremity drape 1

extremity drape 2

extremity drape 3

extremity drape 4

extremity drape 5

delivery drape 1

SURGICAL COVERS.

162 | sterile mayo stand cover,

sterile limb cover

sterile set of limb cover

with an adhesive tape

reinforced - red BESTSELLER

sterile camera cabls cover - perforated,

sterile camera cable cover - unperforated,

with 2 adhesive tapes **BESTSELLER**

ACCESSORIES

surgical covers

gynaecology drape 1

gynaecology drape 2

ophthalmic drape 1-IF

ophthalmic drape 2-IF

ophthalmic drape 3-IF

ophthalmic drape 4-IF

- 156 abdo-perineal drape 1
- 156 laparoscopy drape 3 157 orthopaedic drape 1

157

157

157

158

158

158

158

159

159

159

159

160

160

160

160

162

162

162

16.3

neurosurgery and otolaryngology

gynaecology set 2 BESTSELLER

- 191 | otolaryngology set 1 191 otolaryngology set 2
- 191 | craniotomy set 1-IF

186 | percutaneous lithotripsy set 1-IP NEW

- 192 | ablation set 1-IF 192 angiography set 2
- 192 angiography set 3 NEW

SURGICAL GOWNS

| 194 194 194 195 195 195 196 | sterile gown NORMAL - 22 surgical gown STANDARD - 35 BESTSELLER surgical gown STANDARD - 35 (economic version) surgical gown STANDARD - 40 NEW surgical gown STANDARD - 45 surgical gown STANDARD PLUS - 35 BESTSELLER surgical gown STANDARD PLUS - 35 (economic version) surgical gown STANDARD PLUS - 40 NEW surgical gown STANDARD PLUS - 40 NEW |
|---|--|
| 196 | |
| 196 | surgical gown STANDARD PLUS - 45 |
| 196 | surgical gown EXTRA SAFE - 35 BESTSELLER |
| 196 | surgical gown COMFORT - 68 |
| 197 | surgical gown COMFORT PLUS - 68 |

non-sterile surgical gowns

- 197 | surgical gown STANDARD 35
- 197 | surgical gown PREMIUM 35



SHOE COVERS

| h | sr | 10 | oc | ł | |
|---|----|----|----|---|--|

KOKO long medical cap with DORA long medical cap with snood

OLAF medical cap with ties NEW

MAYA medical cap with rubber NEW

medical mask with earloops medical mask with earloops NEW

213 | patient's gown, tied-up, non-woven patient's gown, Y-neck, non-woven patient's gown, V-neck, non-woven patient's gown, for delivery, non-woven patient's gown, for children , non-woven NEW colonoscopy boxers, SMS non-woven colonoscopy boxers, PP non-woven 215 gynecological skirt, non-woven

- 217 | shoe covers
- 217 | shoe covers, calf-high

CAPS, PROTECTIVE MASKS

- 219 | disposable medical bed sheet disposable medical pillow case 219
- 219 disposable medical quilt cover
- 219 disposable medical bedding set
- 220 disposable absorbent underpad
- disposable absorbent underpad for operating table 220
- 220 absorbent underpad for children **NEW**
- 221 disposable medical bed sheet, reinforced
- 221 thermal patient cover soft
- 221 thermal patient cover



EXAMINATION GLOVES

The device is intended for patient diagnostics, providing a barrier to the transmission of microorganisms between the operator and the patient, minimizing the risk of cross-contamination.

NITRILE GLOVES are made of synthetic acrylonitrile butadiene rubber (NBR), they are an alternative to latex gloves. They do not contain latex proteins, and the content of chemical accelerators is minimized, so they are friendly to the skin.

They are characterized by:

- high mechanical strength (also for puncture),
- chemical and biological resistance (constitute a barrier for viruses and bacteria),
- flexibility, thanks to which they adjust to the hand, providing the work comfort to the user.

VINYL GLOVES made of PVC (polyvinyl chloride) are intended for low-risk procedures due to limited mechanical, chemical and biological resistance. The glove material does not cause Type I allergic reactions, ZARYS vinyl gloves are free of harmful phthalates.

LATEX GLOVES made of natural rubber are recognized as the best barrier against blood borne pathogens and part of hazardous chemical substances. They are also characterized by high flexibility, providing the work comfort to the user by fitting to the hand like a "second skin". Latex gloves, due to the natural raw material, are biodegradable.

STERILE, SURGICAL AND PROTECTIVE GLOVES

Made from natural rubber latex, intended to be worn on hands in surgical settings to provide barrier against potentially infectious fluids or other contaminants, and to protect the patient by ensuring sterility of the wound area. These gloves are intended for single use only, and are classified as Medical Devices Class IIa as well as a Personal Protective Equipment Category III. The indicated Gloves are designed to protect against substances and mixtures which are hazardous and against harmful biological agents.





medical **gloves**

nitrile gloves I latex gloves I vinyl gloves I surgical gloves

ZARYS

WWW.ZARYS.COM

easyCARE nitrile

| type | examination glove, nitrile |
|-------------------------------|--|
| external surface | micro-textured with additional texturing at the fir |
| internal surface | chlorinated |
| cuff | beaded |
| AQL | 1.0 |
| classification | class I medical device personal protective equipment category III |
| length | 245 mm |
| features | increased tactile sensitivity |
| packing unit/ bulk packing | 100 pcs / 10 x 100 pcs 200 pcs / 10 x 200 pcs |
| REF | RNBXX10001 (XX:XS-XL) RNBXX20001 (XX:XS-XL) |
| easyCARE nitrile | Researchered water wat |

easyCARE Ditrile

| | EW | |
|---|----|--|
| N | | |
| | | |

| type | examination glove, nitrile , with extended cuff |
|-------------------------------|--|
| external surface | micro-textured with additional texturing at the f |
| internal surface | chlorinated |
| cuff | beaded |
| AQL | 1.0 |
| classification | class I medical device personal protective equipment category III |
| length | 300 mm |
| features | increased tactile sensitivity, extended cuff |
| packing unit/ bulk packing | 100 pcs / 10 x 100 pcs |
| REF | RNBPXX10001 (XX:S-XL) |
| 6666 | |
| | LING AND |
| | sizes |
| | |

NITRILE GLOVES



easyCARE nitrile **pink**

| type | examination glove, nitrile | |
|-------------------------------|--|--------------|
| external surface | micro-textured with additional texturing at the fingertips | |
| internal surface | examination glove, nitrile micro-textured with additional texturing at the fingertips chlorinated beaded | |
| cuff | beaded Billion | |
| AQL | 1.0 | |
| classification | class I medical device personal protective equipment category III | |
| length | 240 mm | |
| features | intense colour | |
| packing unit/ bulk packing | 100 pcs / 10 x 100pcs | |
| REF | RNBRXX10001 (XX:XS-L) | 100 201 |
| easyCARE nitrile | | |
| | ISON ASTM | 5604 F167 |

1ARYS AQL 1.0 XS S M L

ceptable quali level EN 455-1 EN ISO 374-2

mediCARE nitrile Hydra+

Dermatology solution for skin problems

gloves that protect and nourish the hands, while providing long-lasting hydration



Almond oil

egenerates, moisturizes, slows the aging process of the skin, smoothens, emoving dead cells, brighten the skin, soothens irritations and accelerates trengthens the hydro-lipid layer of the skin



Vitamin E



Anti-wrinkle effects proven in clinical test* *mediCARE Hydra+ glove presents a better anti-wrinkle effect compared to a traditional nitrile glove with visible reduction of wrinkles from the first day of use with progressive increase in improvment until the end of the tests.

Nitrile gloves with moisturizing and nourishing coating

| type | examination glove, nitrile |
|-------------------------------|---|
| external surface | micro-textured with additional texturing at the fi |
| internal surface | chlorinated |
| cuff | beaded |
| AQL | 1.0 |
| classification | class I medical device personal protective equipment category III |
| length | 240 mm |
| features | moisturizing and nourishing coating |
| packing unit/ bulk packing | 100 pcs / 10 x 100 pcs |
| REF | RNBHXX10003 (XX:XS-XL) |



easyCARE nitrile **black**

(11.

external internal

| | | and the state of the |
|-------------------------------|--|--|
| type | examination glove, nitrile | and the second sec |
| external surface | micro-textured with additional texturing at the fingertips | Harding and |
| internal surface | chlorinated | urite Basylillt nitrile |
| cuff | beaded | terile - seiffitt in- |
| AQL | chlorinated beaded 1.0 BESHMIT | 6003- |
| classification | class I medical device personal protective equipment category III | And the form |
| length | 240 mm | |
| features | increased endurance | |
| packing unit/ bulk packing | 100 pcs / 10 x 100pcs | (Mins currile |
| REF | RNBCXX10001 (XX:XS-XL) | 100 10 easyCARE nitrile |
| ISYCARE nitrile | Sizes Si | r k k k k k k k k k k k k k |



tural antioxidant, vitamin of youth – flights off free radicals and prevents inkles, moisturizes for a long time



Glicerin

moisturizes, improves skin appearance, smoothens and soothes irritations











rie Hydrat

mediCARE nitrile SafeDon

| type | examination glove, nitrile |
|-------------------------------|--|
| external surface | micro-textured with additional texturing at the fingertips |
| internal surface | chlorinated |
| cuff | beaded |
| AQL | 1.0 |
| classification | class I medical device personal protective equipment category III |
| length | 240 mm |
| features | single glove removal system always by the cuff |
| packing unit/ bulk packing | 250 pcs (240 for XL) / 8 x 250 pcs (240 pcs) |
| REF | RNBSXX25001 (RNBSXL24001 for XL) (XX:XS-XL) |







SafeDon^m

HOSPITAL INFECTION PREVENTION SYSTEM

SINGLE GLOVE REMOVAL SYSTEM always by the cuff

reduced risk of cross-infection by




Gloves packaged under controlled conditions of **microbiological purity** reduced glove consumption by



Space-saving wall-mounted holder keeps gloves away from potential contamination

SafeDon[™] holder types



SafeDon 250 pcs SafeDon compact 50 pcs

SafeDon 250 pcs SafeDon compact **50 pcs**

mediCARE nitrile **SafeDon**[®] **COMPACT**

examination glove, nitrile external surface micro-textured with additional texturing at the fingertips internal surface chlorinated beaded AQL 1.0 classification class I medical device personal protective equipment category III length 240 mm features single glove removal system always by the cuff 50 pcs / 20 x 50 pcs packing unit/ bulk packing RFF RNBSXX5001 (XX:XS-XL)



MODURA RAIL



SafeDon 250 pcs SafeDon compact 50 pcs

COMPACT





compatible with

SafeDon compact 50 pcs













HGHRISK

INTENDED FOR: DIAGNOSTIC PROCEDURES WHERE THERE IS AN INCREASED RISK OF INFECTION

LENGTH OF GLOVE: 30 CM **ADDITIONAL TEXTURE** ON FINGERTIPS **POWDER - FREE GLOVE, INTERNALLY CHLORINATED**





ZARYS

WWW.ZARYS.COM

LATEX **AND VINYL** GLOVES













easyCARE latex

| type | examination glove, latex |
|-------------------------------|--|
| external surface | micro-textured with additional texturing at the fingertips |
| internal surface | powdered |
| cuff | beaded |
| AQL | 1.0 |
| classification | class I medical device personal protective equipment category III |
| length | 245 mm |
| features | biodegradable |
| packing unit/ bulk packing | 100 pcs / 10 x 100 pcs |
| REF | RLPXX10007 (XX:XS-XL) |
| ols the | |



easyCARE vinyl PF

| type | examination glove, vinyl |
|-------------------------------|--|
| external surface | smooth |
| internal surface | powder-free |
| cuff | beaded |
| AQL | 1.5 |
| classification | class I medical device personal protective equipment category III |
| length | 240 mm |
| packing unit/ bulk packing | 100 pcs / 10 x 100 pcs |
| REF | RWBXX10001 (XX:XS-XL) |



easyCARE latex PF

| | | | All Control of the Co |
|-------------------------------|--|------------------------|---|
| type | examination glove, latex | | |
| external surface | micro-textured with additional texturing at t | hefingertips | - MALL LOICH PF |
| internal surface | powder-free, polymer coated | assignt assignment | - aculture |
| cuff | beaded | agestion. | Suce |
| AQL | 1.0 | | |
| classification | class I medical device personal protective equipment category III | | |
| length | 245 mm | Car | |
| features | easier donning due to inner polymer coating / biodegradable | 1 | Refs and lates pt |
| packing unit/ bulk packing | 100 pcs / 10 x 100 pcs | 100 | easy can or (1) |
| REF | RLBXX10007 (XX:XS-XL) | R | 2000 E 14875 |
| EasyGARE LOCKER PF | | NON E LATEX RHT SYEARS | |
| e Cite | 100 MK | Sizes | protection against blood-borne pathogens Protection 3 against blood-borne pathogens Protection 3 against Blood-borne pathogens Protection 3 against Blood 3 ag |

20 www.zarys.com











ZARYS

WWW.ZARYS.COM



internal surf

| type | natural rubber latex , sterile |
|-------------------------------|--|
| ternal surface | textured |
| ternal surface | powdered |
| cuff | extended cuff, beaded |
| AQL | 0.65 |
| classification | class lla medical device personal protective equipment category III |
| features | standard powdered latex surgical glove |
| packing unit/ bulk packing | 1 pair / 1 x 50 pair |
| REF | RCHLP-XX-50-03 (XX:5.5-9.0) |



5.5 - 9.0



internal

| type | natural rubber latex, sterile |
|-------------------------------|--|
| ternal surface | textured |
| iternal surface | powder-free |
| cuff | extended cuff, beaded |
| AQL | 0.65 |
| classification | class lla medical device personal protective equipment category III |
| features | standard powder-free latex surgical glove |
| packing unit/ bulk packing | 1 pair / 1 x 50 pair |
| REF | RCHLB-XX-50-03 (XX:5.5-9.0) |
| sizes | 5.5 - 9.0 |



SURGICAL GLOVES









www.zarys.com 23

ZARYS International Group prides itself on having a long lasting experience gained in the medical industry by being a leading manufacturer and supplier of single use medical devices.

A wide range of products and complex solutions makes it possible to successfully respond to our customers' needs and meet their expectations.

High standards of design, production, storing and transportation allow to maintain quality, safety and comfort of use to both the medical personnel and patients. A wide range of single use products include basic medical devices and those dedicated to specialized procedures.

ZARYS branded products find a wide application in anesthesiology, urology, gynecology, oxygen and fluid therapy and in many other medical specialities. They make it possible to perform medical procedures including patient diagnostics, prophylaxis, treatment and care.

basic medical devices

basic injection equipment I infusion therapy I oxygen therapy I anaesthesia products I suction catheters, stomach tubes I urological catheters, urine bags I general surgery I basic products for gynecology I other products



ZARYS

dispoFINE

injection needle

- needle tube is made of smooth and durable stainless steel
- sheath and hub made of polypropylene
- thin walls assure adequate stream pace • tip sharpened in three surfaces (back cut)
- bevel type: LB/BL standard long bevel
- surface of needle tube covered with sliding agent (silicone)
- glued needle elements to prevent their separation (glue epoxy resin)
- hubs coloured in accordance with ISO code allowing fast recognition of needle size
- sinale-use
- non-toxic, non-pyrogenic
- packaging: 1 pc./paper-foil

indication

· device designed for intravenous, intramuscular and hypodermic injections

| REF | size | colour | G | intermediate packaging | REF |
|-----------|-----------------|-------------|----|---------------------------|-------------|
| IN-03X13 | 0,3 mm x 13 mm | yellow | 30 | 100 pcs | |
| IN-04X19 | 0,4 mm x 19 mm | medium grey | 27 | 100 pcs | IN-08X40 |
| IN-045X12 | 0,45 mm x 12 mm | brown | 26 | 100 pcs | *IN-08X50 |
| IN-045X16 | 0,45 mm x 16 mm | brown | 26 | 100 pcs | IN-09X25 |
| IN-045X22 | 0,45 mm x 22 mm | brown | 26 | 100 pcs | IN-09X40 |
| IN-05X16 | 0,5 mm x 16 mm | orange | 25 | 100 pcs | IN-11X40 |
| *IN-05X20 | 0,5 mm x 20 mm | orange | 25 | 100 pcs | IN-12X40 |
| IN-05X25 | 0,5 mm x 25 mm | orange | 25 | 100 pcs | IN-12X50 |
| IN-05X40 | 0,5 mm x 40 mm | orange | 25 | 100 pcs | IN-16X40 |
| IN-06X25 | 0,6 mm x 25 mm | deep blue | 23 | 100 pcs | IN-18X40 |
| IN-06X30 | 0,6 mm x 30 mm | deep blue | 23 | 100 pcs | IN-21X40 |
| IN-06X40 | 0,6 mm x 40 mm | deep blue | 23 | 100 pcs | |
| IN-07X30 | 0,7 mm x 30 mm | black | 22 | 100 pcs | |
| IN-07X40 | 0,7 mm x 40 mm | black | 22 | 100 pcs | bulk packag |
| *IN-07X50 | 0,7 mm x 50 mm | black | 22 | 100 pcs | 50 x 100 p |
| | | | | | |



dispoGUARD

safety injection needle

- · safety cover can be activated with one hand and it shields needle immediately after injection
- cover integrated with needle
- tip sharpened in three surfaces (back cut)
- · colour coding system for size, marked on hub
- activation of safety mechanism signalized by audible click
- needle is securely closed inside activated safety cover
- single-use
- non-toxic, non-pyrogenic
- packaging: 1 pc./paper-foil

indication

· device designed for intravenous, intramuscular and hypodermic injections

| REF | size | colour | G | intermediate packaging |
|------------|-----------------|------------|----|---------------------------|
| *INB-03X13 | 0.3 mm x 13 mm | vellow | 30 | 50 pcs |
| INB-045X12 | 0.45 mm x 12 mm | brown | 26 | 50 pcs |
| INB-05X25 | 0.5 mm x 25 mm | orange | 25 | 50 pcs |
| INB-05X20 | 0,5 mm x 40 mm | 0 | 25 | |
| | | orange | | 50 pcs |
| INB-06X25 | 0,6 mm x 25 mm | deep blue | 23 | 50 pcs |
| *INB-06X30 | 0,6 mm x 30 mm | deep blue | 23 | 50 pcs |
| *INB-07X30 | 0,7 mm x 30 mm | black | 22 | 50 pcs |
| INB-07X40 | 0,7 mm x 40 mm | black | 22 | 50 pcs |
| INB-08X40 | 0,8 mm x 40 mm | deep green | 21 | 50 pcs |
| INB-09X40 | 0,9 mm x 40 mm | yellow | 20 | 50 pcs |
| INB-11X40 | 1,1 mm x 40 mm | cream | 19 | 50 pcs |
| INB-12X40 | 1,2 mm x 40 mm | pink | 18 | 50 pcs |



Ø bulk packaging: 50 x 50 pcs

BASIC INJECTION EQUIPMENT

| | | 11 | ۱. | |
|----|----------------|------------|----|--------------|
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | size | colour | G | intermediate |
| | | | | packaging |
| | | | | |
| 10 | 0,8 mm x 40 mm | deep green | 21 | 100 pcs |
| 50 | 0,8 mm x 50 mm | deep green | 21 | 100 pcs |
| 25 | 0,9 mm x 25 mm | yellow | 20 | 100 pcs |
| 10 | 0,9 mm x 40 mm | yellow | 20 | 100 pcs |
| 10 | 1,1 mm x 40 mm | cream | 19 | 100 pcs |
| 10 | 1,2 mm x 40 mm | pink | 18 | 100 pcs |
| 50 | 1,2 mm x 50 mm | pink | 18 | 100 pcs |
| 10 | 1,6 mm x 40 mm | white | 16 | 100 pcs |
| 10 | 1,8 mm x 40 mm | blue-grey | 15 | 100 pcs |
| 10 | 2,1 mm x 40 mm | pale green | 14 | 100 pcs |
| | | | | |













dispoSULIN

insulin pen needle

features

- made of polypropylene, polyethylene and stainless steel
 triple-point needle blade allows easy and atraumatic skin penetration
 covering the needle with silicone lubricant reduces resistance during
- injection • thin walls of needle ensure easy flow of medicine and the right flow rate
- coloured cap and paper cover of packaging make it easy to quickly recognize size of needle
- available in various sizes, allowing you to adapt to individual needs • universal, compatible with all pen models
- single-use
- non-pyrogenic, non-toxic packaging: 100 pcs/box

indication

- for subcutaneous injection
- dedicated for use with pre-filled pens (insulin pens)

| REF | size c | olour | G | intermediate packaging |
|------------|-----------------|--------|----|---------------------------|
| INP-033X12 | 0,33 mm x 12 mm | pink | 29 | 100 pcs |
| INP-030X08 | 0,30 mm x 8 mm | yellow | 30 | 100 pcs |
| INP-025X06 | 0,25 mm x 6 mm | purple | 31 | 100 pcs |
| INP-025X08 | 0,25 mm x 8 mm | blue | 31 | 100 pcs |



\$ bulk packaging: 50 x 100 pcs



features

- easy to use, ergonomic shape
 body with knurled surface for a secure grip
 sharp receiving spike enabling tight and safe connection with containers of different capacity

dicoSPIKE, dicoSPIKE Chemo

- or online enclosed
 receiving spike groove guarantees complete emptying of the container
 hermetic self-Locking cap securing the connector against contamination
 available in version with air filter and in version with air and particle filter
 air filter protects against microbiological contamination and against harmful aerosols
- particle filter is an effective barrier to protect drugs and liquids from , contamination
- female Luer-Lock connector compatible with products with either male or female Luer-Lock
- non-pyrogenic
- single-usepackaging: 1 pc./ paper-foil

indication

- · device for needle-free dissolving and / or withdrawing drugs from multi-
- dose packages (bottles, vials, etc.) in both directions dicoSPIKE Chemo model is designed for medicines releasing dangerous aerosols such as cytostatics

| REF | cap colour | air filter | particle filter | purpose | interme packag |
|------|---------------|---------------|--------------------|--------------------|-------------------|
| PPLZ | green | 0,1 µm | - | standard medicines | 100 pc |
| PPLN | blue | 0,1 µm | 5 µm | standard medicines | 100 pcs |

| dicoSPIKE | Chemo | | | | |
|-----------|---------------|---------------|--------------------|-------------|--------------------|
| REF | cap colour | air filter | particle filter | purpose | interme packagi |
| PPLC | red | 0,1 µm | 5 µm | cytostatics | 100 pcs |



side hole needle

features

- made of stainless steel
- the needle hub ideally suited for Luer and Luer-Lock (to ensure full tightnessof the connection with syringes)
- CE sign comply with PN-EN-ISO 7864
- individual packaging comply with PN-EN 556 and PN-EN 868, type blister-pack coloured cap and paper cover of packaging make it easy to quickly recognize single-use

indication

- made of high quality materials, using the latest technology, which guarantees its high quality and failure-free, intended usage
 special design needle with a pencil cut with a side hole, designed to
- withdraw medication from a vial, including vials with a rubber stopper
 needle construction prevents from accidental, unintended production
- of rubber scrubs, in order to avoid contamination of the drug
- side hole with special cut, placed on the top of the needle enables free dissolution of the drug and limits foaming by directing the fluid directly onto the walls of the vial

| REF | size | colour | G | intermediate packaging | |
|----------|-------------|--------|----|---------------------------|----------------------------------|
| HPP-1830 | 1,2 x 30 mm | pink | 18 | 100 pcs | bulk packaging: 100 x 100 pcs |








bulk packaging: 10 x 100 pcs

nediate iging CS



construction of the syringe



GET TO KNOW OUR SYRINGES

catheter tip 90°

catheter tip 45°

syringes

the drug/liquid residues

type of connection used for rinsing (cleaning)


insulin and tuberculin

special design of the plunger seal reduces

catheters, stomach tubes and feeding purposes

nozzle cut at 45 degrees increases the lumen surface

nozzle types:



Luer slip, centric

tip of the syringe is pushed into the hub of the needle or other device (e.g. catheters) and twisted. It allows to remove the connection immediately if necessary



Luer slip, eccentric

eccentric location of the tip is used in procedures involving the administration of drugs and fluids where closer proximity to skin is necessary



Luer-Lock

designed to secure connection between the syringe and the device (needle, tube, catheter, stopcock, etc.). The tip is equipped with thread to screw in the syringe to prevent disconnection

duoNEX

syringe type:



2-part - without plunger sealing



3-part - with plunger sealing

graduation type:

.



extended graduation



dicoNEX

regular graduation



package colour and shape code:





duoNEX

single-use syringe, 2-part, Luer

features • white line

- transparent barrel allows full visual control of content
- white contrasting plunger
 Luer nozzle centrally located in 2 ml syringe, eccentric in 5, 10, 20 ml syringes smooth plunger sliding
- limitation flange securing against accidental removal of plunger
 clear, legible and permanent black graduation facilitating dosage,
- adequately:
- every 0,1 ml for capacity of 2 ml
- every 0,2 ml for capacity of 5 ml
- every 0,5 ml for capacity of 10 ml
- every 1 ml for capacity of 20 ml
- extended graduation that includes a 10% extension of the nominal capacity • name of syringe and manufacturer's logo placed on barrel
- non-pyrogenic, non-toxic
- packaging: 1 pc./paper-foil

indication

· 2-piece syringe for subcutaneous, intramuscular and intravenous infusions intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient

| REF | type | capacity | position of nozzle | intermediate packaging | bulk packag |
|-------------|------|----------|-----------------------|---------------------------|----------------|
| | | | | | |
| 002ML-2CZ-B | Luer | 2 ml | centric | 100 pcs | 30 x 100 pcs |
| 005ML-2CZ-B | Luer | 5 ml | eccentric | 100 pcs | 24 x 100 pcs |
| 010ML-2CZ-B | Luer | 10 ml | eccentric | 100 pcs | 16 x 100 pcs |
| 020ML-2CZ-B | Luer | 20 ml | eccentric | 50 pcs | 16 x 50 pcs |

STERILEEO

dicoNEX SN

single-use syringe, 3-part, Luer, with needle alongside the syringe

features


- in set with needle alongside the syringe
- transparent barrel allows full visual control of content, green contrasting plunger
- · Luer nozzle placed centrally for 2 ml capacity, eccentric for 5, 10, 20, 50 ml capacity • sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- · limitation flange securing against accidental removal of plunge • name of syringe and manufacturer's logo placed on barrel
- · clear, legible and permanent black graduation facilitating dosage, adequately - every 0,1 ml for capacity of 2 ml
- every 0,2 ml for capacity of 5 ml (extended to 6 ml)
- every 0,2 ml for capacity of 10 ml
- every 0,5 ml for capacity of 20 ml
- every 1 ml for capacity of 50 ml • needle tube is made of smooth and durable stainless steel
- needle tip sharpened in three surfaces
- bevel type: LB/BL standard long bevel
- hubs coloured in accordance with ISO code
- non-toxic, non-pyrogenic

indicatio

· 3-piece syringe for subcutaneous, intramuscular and intravenous infusions

| REF | type | capacity | position of nozzle | needle size | needle colour | intermediate packaging | bulk packaging |
|--------------|------|----------|-----------------------|--------------|------------------|---------------------------|-------------------|
| 002ML-3CZ-SN | Luer | 2 ml | centric | 0,6mm x 30mm | deep blue | 100 pcs | 30 x 100 pcs |
| 005ML-3CZ-SN | Luer | 5 ml | eccentric | 0,7mm x 30mm | black | 100 pcs | 24 x 100 pcs |
| 010ML-3CZ-SN | Luer | 10 ml | eccentric | 0,8mm x 40mm | deep green | 100 pcs | 16 x 100 pcs |
| 020ML-3CZ-SN | Luer | 20 ml | eccentric | 0,8mm x 40mm | deep green | 50 pcs | 16 x 50 pcs |
| 050ML-3CZ-SN | Luer | 50 ml | eccentric | 0,8mm x 40mm | deep green | 25 pcs | 16 x 25 pcs |





extended graduation

.....

Luer slip and Luer-Lock syringes

...........



aging

DCS pcs pcs



dicoNEX

single-use syringe, 3-part, Luer

features

- transparent barrel allows full visual control of content
- green contrasting plunger • Luer nozzle
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- · limitation flange securing against accidental removal of plunger
- special design of the 1 ml syringe plunger stopper reduces dead space
- name of syringe and manufacturer's logo placed on barrel
- clear, legible and permanent black graduation facilitating
- dosage, adequately:
- every 0,01ml for capacity of 1ml
- every 0,1 ml for capacity of 2 ml
- every 0,2 ml for capacity of 5 ml
- every 0,2 ml for capacity of 10 ml - every 0,5 ml for capacity of 20 ml
- every 1 ml for capacity of 50/60 ml
- non-pyrogenic, non-toxic
- packaging: 1 pc./paper-foil

indication

• 3-piece syringe intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient

| REF | type | capacity | position of nozzle | intermediate packaging | bulk packaging |
|--------------|------|----------|-----------------------|---------------------------|-------------------|
| 001ML-3CZ-BL | Luer | 1 ml | centric | 100 pcs | 32 x 100 pcs |
| 002ML-3CZ-BL | Luer | 2 ml | centric | 100 pcs | 30 x 100 pcs |
| 005ML-3CZ-BL | Luer | 5 ml | centric | 100 pcs | 24 x 100 pcs |
| 010ML-3CZ-BL | Luer | 10 ml | centric | 100 pcs | 16 x 100 pcs |
| 020ML-3CZ-BL | Luer | 20 ml | eccentric | 50 pcs | 16 x 50 pcs |
| 050ML-3CZ-BL | Luer | 50/60 ml | eccentric | 25 pcs | 16 x 25 pcs |



dicoNEX single-use syringe, 3-part, Luer-Lock

features

- transparent barrel allows full visual control of content
- green contrasting plunger
- Luer-Lock nozzle
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- · limitation flange securing against accidental removal of plunger
- name of syringe and manufacturer's logo placed on barrel • clear, legible and permanent black graduation facilitating dosage, adequately:
- every 0,1 ml for capacity of 3 ml
- every 0,2 ml for capacity of 5 ml and of 10 ml
- every 0,5 ml for capacity of 20 ml
- every 1 ml for capacity of 30 ml
- non-pyrogenic, non-toxic
- packaging: 1 pc./paper-foil

indication

· 3-piece syringe intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient

| REF | type | capacity | position of nozzle | colour of plunger | intermediate packaging | bulk packaging |
|-----------------|-----------|----------|-----------------------|----------------------|---------------------------|-------------------|
| 003ML-3CZ-LL-BL | Luer-Lock | 3 ml | centric | green | 100 pcs | 30 x 100 pcs |
| 005ML-3CZ-LL-BL | Luer-Lock | 5 ml | centric | green | 100 pcs | 24 x 100 pcs |
| 010ML-3CZ-LL-BL | Luer-Lock | 10 ml | centric | green | 100 pcs | 16 x 100 pcs |
| 020ML-3CZ-LL-BL | Luer-Lock | 20 ml | centric | green | 50 pcs | 16 x 50 pcs |
| 030ML-3CZ-LL-BL | Luer-Lock | 30 ml | centric | green | 50 pcs | 16 x 50 pcs |





dicoNEX

single-use syringe, 3-part, Luer-Lock + amber

features

- transparent cylinder with a transparent or amber color
- green contrasting plunger (amber-coloured syringe white contrasting plunger) Luer-Lock type termination
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
 - limitation flange securing against accidental removal of plunger
 name of syringe and manufacturer's logo placed on barrel

 - clear, legible and permanent black graduation facilitating dosage scaling extension to 60 ml
 - non-pyrogenic, non-toxic
 - packaging: 1 pc./paper-foil

indication

- 3-piece syringe intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient
- amber-coloured syringe is intended for drawing up and administering of photosensitive medicines

| | REF | type | capacity | position of nozzle | colour of plunger | intermediate packaging |
|---|-------------------|-----------|----------|-----------------------|----------------------|---------------------------|
| | 050ML-3CZ-LL-PI | Luer-Lock | 50/60 ml | centric | green | 25 pcs |
| a | mber type | | | | | |
| | 050ML-3CZ-LL-B-PI | Luer-Lock | 50/60 ml | centric | white | 25 pcs |
| | | | SYEARS | | | |

dicoNEX

single-use catheter syringe, 3-part

features

- 3-piece syringe with centrally located nozzle for catheter
- made of polypropylene
- transparent barrel allows full visual control of content
- green contrasting plunger (capacity 50/60 ml)
- white contrasting plunger (capacity 100 ml)
- two Luer connectors attached
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- limitation flange securing against accidental removal of plunger
 nozzle cut at 90° (capacity 50/60 ml)
- nozzle cut at 45° (capacity 100 ml)
- catheter nozzle secured with cap (capacity 50/60 ml)
- name of syringe and manufacturer's logo placed on barrel
- extended graduation (capacity 50/60 ml) clear, legible and permanent black graduation facilitating dosage
- non-pyrogenic, non-toxic • packaging: 1 pc./paper-foil

indication

- 3-piece syringe with catheter nozzle
- intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient

| REF | type | capacity | position of nozzle | intermediate packaging | bulk packa |
|------------------|--------------|----------|-----------------------|---------------------------|---------------|
| 050ML-3CZ-CEW-BL | catheter tip | 50/60 ml | centric | 25 pcs | 16 x 25 pcs |
| 100ML-3CZ-CEW-BL | catheter tip | 100 ml | centric | 25 pcs | 4 x 25 pcs |
| | | | | | |











kaging

dicoTUBER

tuberculin syringe

features

- with bypacked needle 0,45 x 13 mm
- made entirely of polypropylene
- white contrasting plunger • transparent barrel allows full visual control of content
- Luer nozzle
- smooth plunger sliding thanks to rubber sealing
- · limitation flange securing against accidental removal of plunger • clear, legible and permanent black graduation every 0,01 ml that allows precise dosing of medicine/liquid
- · special design of plunger stopper reduces the residual volume of medicine/liquid
- name of syringe and manufacturer's logo placed on barrel single-use
- non-pyrogenic, non-toxic • packaging: 1 pc./paper-foil


indication

• 3-piece syringe intended for precise drawing up and administering of medicines/liquids

| REF | no. of units | capacity | intermediate packaging | |
|-------------|--------------|----------|---------------------------|---------------------------------|
| TU-U20-1-BL | 20 | 1 ml | 100 pcs | bulk packaging: 32 x 100 pcs |

Ø bulk packaging: 32 x 100 pcs





0.2 0.3 0.4 0.4 0.5 0.6 0.7 0.8 0.9 1.0 m

dicoSULIN insulin syringe

features

- with bypacked needle 0,4 x 13 mm
- transparent barrel made of polypropylene allows full visual control of content
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- · limitation flange securing against accidental removal of plunger special design of plunger stopper reduces the residual volume
- of medicine/liquid clear, legible and permanent black graduation facilitating dosage
- name of syringe and manufacturer's logo placed on barrel
- sinale-use
- non-pyrogenic, non-toxic
- packaging: 1 pc./paper-foil

indication

· 3-piece syringe for subcutaneous, intramuscular and intravenous infusions for insulin shots

| REF | no. of units | capacity | colour of plunger | intermediate packaging |
|--------------|--------------|----------|----------------------|---------------------------|
| IN-U40-1-BL | 40 | 1 ml | red | 100 pcs |
| IN-U100-1-BL | 100 | 1 ml | orange | 100 pcs |



ZARYS

34 www.zarys.com

INFUSION THERAPY

SANVIflon

I.V. cannula

features

- needle (stylet) made of stainless steel with back-cut blade
 cannula catheter made of Teflon (FEP), equipped with 4 radiopaque lines • cannula equipped with wings and hand-closed cap of the upper port colour coded according to size
- handle that allows you to insert with one hand
 valve that prevents blood leakage to outside after inserting the cannula into the vessel
- flashback chamber secured with a Luer-Lock cap available in various sizes from 14 G to 26 G
- single-use
- packaging: 1 pc. hard blister (PVC + lacquered paper) with a colour print in accordance with the colour identification of the size

indication

- for cannulation of intravenous, peripheral vessels to:
- multiple injections of intravenous drugs · intravenous infusion of solutions and parenteral nutrition suitable for peripheral veins
- maintain the patient's hydration and / or reduce dehydration if the patient is unable to take sufficient fluid orally
- blood transfusions and blood-like preparations
- ensuring vascular access in the event of the need for rapid drug delivery,
- especially during diagnostic and therapeutic procedures

| REF | size [G] | size [mm] | colour | flow rate [ml/min] | intermediate packaging |
|-----------|-------------|--------------|--------|-----------------------|---------------------------|
| KDF-06-19 | 26 x ¾″ | 0,6 x 19 | violet | 17 | 100 pcs |
| KDF-07-19 | 24 x ¾″ | 0,7 x 19 | yellow | 23 | 100 pcs |
| KDF-09-25 | 22 x 1" | 0,9 x 25 | blue | 36 | 100 pcs |
| KDF-11-32 | 20 x 1¼" | 1,1 x 32 | pink | 56 | 100 pcs |
| KDF-13-45 | 18 x 1¾″ | 1,3 x 45 | green | 90 | 100 pcs |
| KDF-15-45 | 17 x 1¾″ | 1,5 x 45 | white | 142 | 100 pcs |
| KDF-17-45 | 16 x 1¾″ | 1,7 x 45 | grey | 200 | 100 pcs |
| KDF-20-45 | 14 x 1¾″ | 2,0 x 45 | orange | 305 | 100 pcs |
| | | | | | |



SANVIfion S

I.V. cannula without port

features

- needle (stylet) made of stainless steel with back-cut blade
- cannula catheter made of Teflon (FEP), equipped with 4 radiopaque lines
- cannula equipped with coloured fixing wings, without upper port
- handle that allows you to insert with one hand
- valve that prevents blood leakage to outside after inserting the cannula into the vessel
- flashback chamber secured with a Luer-Lock cap
- available in two sizes: 24 G and 26 G
- single-use
- packaging: 1 pc. hard blister (PVC + lacquered paper) with a colour print in accordance with the colour identification of the size

indication

for cannulation of intravenous, peripheral vessels to:

- mvultiple injections of intravenous drugs
- intravenous infusion of solutions
- maintain the patient's hydration and / or reduce dehydration if the patient is unable to take sufficient fluid orally
- ensuring vascular access in the event of the need for rapid drug delivery, especially during diagnostic and therapeutic procedures

| REF | size [G] | size [mm] | colour | flow rate [ml/min] | intermediate packaging |
|------------|-------------|--------------|--------|-----------------------|---------------------------|
| KDFS-06-19 | 26 x ¾″ | 0,6 x 19 | violet | 17 | 100 pcs |
| KDFS-07-19 | 24 x ¾″ | 0,7 x 19 | yellow | 23 | 100 pcs |





easyFLOW IS PREMIUM

infusion set

features

- · 2-way, sharp spike of the drip chamber with anti-slip surface that guarantees tight and reliable connection with fluid containers
- air-vent with anti-bacterial filter and blue flap
- flexible droplet chamber, size 6 cm, equipped with additional pressure wings to facilitate insertion into containers with liquids
- dropper: 20 drops = 1 ml ± 0,1 ml
- iquid filter protecting against larger particles with pore size 15 µm
- precise safe roll clamp with integrated spike holder and with a holder for securing the drain
- name of the manufacturer placed on the roll clamp
- soft, flexible drain
- universal Luer-Lock tip secured by an Air Pass plug with a hydrophobic filter, which allows the drain to be filled without accidental contamination and prevents fluid leakagesingle-use single-use
- non-pyrogenic, non-toxic
- individual packaging: paper-foil with a colored code (blue) identifying the device type

indication

device for fluid transfer in infusion therapy

| REF | version | length | bulk packaging |
|---------|---------------------|--------|----------------------------|
| IS-BF-P | without phthalates | 150 cm | bulk packaging: 250 pcs |
| 13-DF-F | without pricialates | 150 cm | 250 pcs |
| STERILE | | 5YEARS | |



infusion set

features

- 2-way sharp closure-piercing device of the drip chamber with anti-slip
- surface that guarantees tight and reliable connection with fluid containers/bags
 air inlet with anti-bacterial filter and coloured (blue) flap
- pliable drip chamber, size 6 cm, equipped with additional pressure wings to facilitate insertion into containers with liquids
- dropper: 20 drops = 1 ml +/- 0,1 ml
- oropper: 20 orops = 1 mi +7-0, mi
 fluid filter protecting against larger particles with pore size 15 µm
 soft, flexible tubing with Luer-Lock conical fitting
 precise safe flow regulator with integrated holder for the closure-piercing device and with hook for securing the tubing
 name of the manufacturer placed on the flow regulator

- · both tips of the device additionally secured with protective caps
- single-use
- non-pyrogenic, non-toxic packaging: 1 pc./paper-foil with a coloured code (blue) identifying the device

indication

device for fluid transfer in infusion therapy

| | REF | version | length | |
|-----|--------------|--------------------|--------|----------------|
| | IS* | with phthalates | 150 cm | |
| | IS-BF ** | without phthalates | 150 cm | bulk packaging |
| NEW | IS-BF-180 ** | without phthalates | 180 cm | 250 pcs |





10 x 100 pcs

bulk packaging:

> bulk packaging: 10 x 100 pcs





easyFLOW IS

infusion set (dual chamber)



features

- made of high quality materials, using the latest technology, which guarantees
 very high quality and trouble-free, intended use
- Very firght quality and trouble-free, Intended USE
 2-way, 2-plane sharp closure-piercing device of the drip chamber with anti-slip surface that guarantees tight and reliable connection with fluid containers
 air-inlet with anti-bacterial filter and coloured (green) flap
 dropler chamber divided into a flexible lower part and a hard upper part with an additional pressure ring for easy insertion into fluid containers
 page of the manufacturer placed on the dropler chamber
- name of the manufacturer placed on the dropler chamber
- device in its entirety including the spike and tip transparent, allowing visual inspection of the flow along its entire length
- dropper: 20 drops = 1 ml ± 0.1 ml
 fluid filter protecting against larger particles with pore size 15 µm
 soft, flexible tubing with Luer-Lock conical fitting
- precise safe roll clamp with integrated spike holder and with a holder for securing the drain
- both tips of the device additionally secured with protective caps
- single-use
- non-pyrogenic, non-toxic
- individual packaging: paper-foil with a coloured code (blue) identifying the device type

indication

• device intended for connection to an intravenous line for the administration of fluids, photosensitive medicinal products solutions and other body fluids during infusion therapy



easyFLOW IS ECO

infusion set

features

- made of high quality materials
- · 2-way sharp closure-piercing device of the drip chamber that guarantees tight and reliable connection with fluid containers
- air inlet with anti-bacterial filter and colored (blue) flap
- pliable drip chamber, size 4cm • dropper: 20 drops = 1 ml +/- 0.1 ml
- fluid filter protecting against larger particles with pore size 15 µm
- soft, flexible tubing
- universal Luer-Lock conical fitting
- precise safe flow regulator
- · both tips of the device additionally secured with protective caps • packaging: 1pc./foil bag
- single-use
- non-pyrogenic, non-toxic

indication

device for fluid transfer in infusion therapy

REF lenath IS-E* with phthalates 150 cm IS-BF-E** without phthalates 150 cm bulk packaging NEW IS-BF-E-180 ** without phthalates 180 cm

Ø

400 pcs





easyFLOW IS SAFE

infusion set

features

- 2-way, sharp spike of the drip chamber with anti-slip surface that guarantees tight and reliable connection with fluid containers
- air-vent with anti-bacterial filter and blue flap
- flexible droplet chamber, size 6 cm, equipped with additional pressure wings to facilitate insertion into containers with liquids
- dropper: 20 drops = 1 ml ± 0.1 ml
- hydrophilic filter Liquid Lock with a mesh size of 8 µm, which automatically stops the flow when the chamber is emptied. Prevents air from entering the tubing to prevent air embolism
- soft, flexible drain with additional port for drug application • precise safe roll clamp with integrated spike holder and with a holder for securing the drain
- universal Luer-Lock tip
- · Air Pass plug with hydrophobic filter, which allows filling of the tubing without accidental contamination and protects against fluid leakage
- single-use
- non-pyrogenic, non-toxic
- individual packaging: paper-foil with a colored code (blue) identifying the device type

indication

• device for fluid transfer in infusion therapy





easyFLOW IS AMBER

infusion set, amber



features

- amber colouring of the tubing and drip chamber
- · double-channel, sharp spike of the drip chamber with a matted surface with non-slip Surface guaranteeing tight and reliable connection with fluid containers/bags
- air inlet with anti-bacterial filter and coloured (blue) flap
- flexible drip chamber with length 6 cm, equipped with additional pressure wings to facilitating puncture into containers with liquids
- dropper: 20 drops = 1 ml +/- 0,1 ml
- fluid filter protecting against larger particles with pore size 15 µm
- soft and flexible drains with additional port for drug application
- precise safe roll clamp with integrated spike holder and with a holder for securing the drain with name of the manufacturer placed on the roll clamp
- single-use
- non-pyrogenic, non-toxic
- universal Luer-Lock tip
- packaging: 1 pc./paper-foil

indication

· device for fluid transfer photosensitive medicinal and other body fluids in infusion therapy









easyFLOW IS REG

infusion set with precision flow rate regulator

features

- 2-way, sharp closure-piercing device of the drip chamber with anti-slip surface that guarantees tight and reliable connection with fluid containers
 adjustable from 5 to 250 ml/h (dark blue digits on white background) for viscosities below 10% and adjustable from 5 to 200 ml/h (white digits on dark blue background)
- for viscosities 10%-40% ability to completely close (OFF) and open (ON) the flow rate
 air-inlet with anti-bacterial filter and coloured (blue) flap
- dropler chamber and drain made of phthalate-free PVC
- dropper: $20 \text{ drops} = 1 \text{ ml} \pm 0.1 \text{ ml}$
- fluid filter protecting against larger particles with pore size 15 µm
- soft, flexible drain with additional port for drug application
 universal Luer-Lock tip secured by an Air Pass plug with a hydrophobic filter, which allows the drain to be filled without accidental contamination and prevents fluid leakage
- precise safe roll clamp with integrated spike holder and with a holder for securing the drain
- name of the manufacturer placed on the roll camp
- both tips of the device additionally secured with protective caps
- single-use
- non-pyrogenic, non-toxic
- individual packaging: paper-foil with a coloured code (blue) identifying the device type

indication

· device intended for connection to an intravenous line for the administration of fluids, photosensitive medicinal products solutions and other body fluids during infusion therapy with adjustable flow control



precision flow regulator

features


- · precision regulator to set the exact desired flow rate within a range of 5-250 ml/h with the possibility to completely close or open the flow rate
- two adjustment scales, differentiated by colour • adjustment from 5 to 250 ml/h (dark blue digits on white background) for viscosities below 10% and adjustment from 5 to 200 ml/h (white digits on dark

NEW

- blue background) for viscosities of 10%-40%.
- additional Y port for drug administration
- luer lock termination • both ends of the tube protected by protective caps
- single-use
- non-pyrogenic, non-toxic
- unit packaging: 1 piece paper/foil

indication

- dual scale regulator for precise administration of drugs and fluids as well as parenteral nutrition
- · compatible with any infusion set model







conallhealth 岱 Your Partner In Healthcar

easyFLOW TS PREMIUM

transfusion set

features

2-way sharp spike of the drip chamber with matted surface that guarantees tight and reliable connection with fluid containers

NEW

- air-vent with anti-bacterial filter and a colored (red) flap

- flexible droplet chamber
 dropper: 20 drops = 1ml +/- 0.1ml
 special blood filter with large area, pore size 200µm
- soft, flexible drain
- universal luer-lock drain tip secured by an Air Pass stopper with a hydrophobic filter that allow the drain to be filled without accidental contamination and protects against fluid leakage
- precise safe roll clamp with integrated spike holder and with a holder for securing the drain
- name of the manufacturer placed on the roll clamp
 both tips of the device additionally secured with protective caps
- single-use
- non-pyrogenic, non-toxic • individual packaging: blister-pack with a colored code (red) identifying the device type

indication

 device intended for connection to an intravenous line for the administration of blood and blood products

| REF | version | length | |
|----------|--------------------|--------|----------------------------|
| TS-BF-P | without phthalates | 150 cm | bulk packaging: 250 pcs |
| STERILEE | | Syears | |

easyFLOW TS

transfusion set

features

- 2-way sharp closure-piercing device of the drip chamber with anti-slip surface that guarantees tight and reliable connection with containers/bags
- air inlet with anti-bacterial filter and a coloured (red) flap
- pliable drip chamber
- dropper: 20 drops = 1 ml +/- 0,1 ml
- \bullet special blood filter with large area, pore size 200 μm
- soft, flexible tubing
- universal Luer-Lock fitting at the end of tubing
- precise safe flow regulator with integrated holder for the closure-piercing device and with hook for securing the tubing
- name of the manufacturer placed on the flow regulator · both tips of the device additionally secured with protective caps
- sinale-use • non-pyrogenic, non-toxic
- packaging: 1 pc./paper-foil with a coloured code (red) identifying the device

indication

device used for the administration of blood and blood products

| ~~~ |
|-------------------------------|
| bulk packaging 250 pcs |
| ates 150 cm halates 150 cm |









easyFLOW LINE

extension tube for infusion pump

features

- tubing made of PVC
- device available in a transparent and amber version for
- photosensitive medicines
- diameter of the tubing outer: 2,4 mm; inner: 1,24 mm
 Luer-Lock fitting male and female
- single-use
- non-pyrogenic, non-toxic
- packaging: 1 pc./paper-foil

indication

• device used to extend infusion line

| REF | colour of tube | length | intermediate packaging |
|-----------|-------------------|--------|---------------------------|
| P450-BF | transparent | 45 cm | 50 pcs |
| P1500-BF | transparent | 150 cm | 50 pcs |
| PB1500-BF | amber | 150 cm | 50 pcs |







needle free valve

with extension line

features

- housing made of transparent material polycarbonate
- high-quality silicone membrane allows multiple port injections
- flat injection surface smooth and evenly ended port membrane that facilitates disinfection extension line made of transparent polyurethane
- metal-free, suitable for MRI examination
- extension lines fitted with sliding clamps
- male connector protected with applicator
- total length 15 cm
- can be used for 7 days or 350 activations
 filling volume of needle-free port 0.09 ml
- dead space 0.001 ml
- resistance to injection fluid pressure 3 bar (44 PSI)
- resistance to back pressure 2 bar (29 PSI)
- single-use packaging: 1 pc. / paper-foil

indication

| REF | diameters of extension lines | flow | number of needle-free ports | |
|-----------|---------------------------------|------------|--------------------------------|---------|
| PB-T/1-27 | 1,0 x 2,7 mm | 115 ml/min | 1 | CGD |
| PB-T/1-40 | 3,0 x 4,0 mm | 185 ml/min | 1 | |
| PB-T/2-27 | 1,0 x 2,7 mm | 100 ml/min | 2 | and the |
| PB-T/2-40 | 3,0 x 4,0 mm | 165 ml/min | 2 | |
| PB-T/3-27 | 1,0 x 2,7 mm | 100 ml/min | 3 | 0000 |
| PB-T/3-40 | 3,0 x 4,0 mm | 160 ml/min | 3 | |
| | | | | |

| PB-T/4-27 | 1,0 x 2,7 mm | 65 ml/min | 4 |
|-----------|--------------|------------|---|
| PB-T/4-40 | 3,0 x 4,0 mm | 135 ml/min | 4 |

needle-free valve

features

- body made of transparent material polycarbonate
- high-quality silicone membrane allows multiple port injections • flat injection surface - smooth and evenly ended port membrane
- that facilitates disinfection
- metal-free, suitable for MRI examination
- male connector of the transparent port is secured with a cap allowing for aseptic removal of the item from the packaging* • may be used for 7 days or up to 350 activations
- priming capacity: 0,09 ml • high flow rate: 350 ml/min
- withstands injection fluid pressure of 3 bars (44 PSI)
- withstands back pressure of 2 bars (29 PSI)
- single-use
- packaging: 1 pc./paper-foil

indication

- intended for closing infusion ports
- device intended for multiple, needle-free injections
- (liquid drug administration, blood sample taking)
- protects infusion lines from contamination and infection
- during infusion therapy
- does not require using safety caps • dedicated for standard Luer and Luer Lock connectors

Ø

20 x 100 pcs

REF intermediate colour packaging PB-N blue 100 pcs bulk packaging: PB-T 100 pcs transparent





3-way stopcock

features

- made of polycarbonate
- transparent housing
- clear optical identifier of open/close position (marked flow directions) operating in 360 degrees
- consistent, straight flow channels and an optimal internal diameter ensure correct pressure measurement
- rotating male Luer-Lock
- all inlets protected with caps
- resistant to lipids
- single-use
- non-pyrogenic, non-toxic • packaging: 1 pc./paper-foil

indication

- device designed for managing the flow of infusion fluids during administration
- enables switching the infusion flow path
- · for direct medicine/fluid supply

| REF | colour | intermediate packaging |
|--------|--------|---------------------------|
| KT-002 | blue | 50 pcs |







- for blood pressure monitoring

- device used to close infusion ports for multiple, needle-free injections and to protect the infusion line from contamination does not require using safety caps
- dedicated for standard Luer and Luer Lock connectors
 - 1-1)
 - 0-3 25-20
- C 13 0 5 5







3-way stopcock with extension tube

feature

- three-way stopcock with an extension tube of various lengths (from 7 cm to 100 cm)
- transparent housing
- clear optical identifier of open/close position (marked flow directions) • three-way stopcock made of polycarbonate, extension tube made of PVC

bulk packaging: 20 x 50 pcs

- all inlets protected with caps
- resistant to lipids
- single-use
- non-pyrogenic, non-toxic packaging: 1 pc./paper-foil

indication

- device designed for managing the flow of infusion fluids during administration
- enables switching the infusion flow path
- for blood pressure monitoring
- for direct medicine/fluid supply

| REF | length of tube | intermediate packaging |
|---------|-------------------|---------------------------|
| KTP-7 | 7 cm | 50 pcs |
| KTP-10 | 10 cm | 50 pcs |
| KTP-25 | 25 cm | 50 pcs |
| KTP-50 | 50 cm | 50 pcs |
| KTP-75 | 75 cm | 50 pcs |
| KTP-100 | 100 cm | 50 pcs |





ZARYS

cannula stopper

features

- precisely and safely closes standard ports
- Luer-Lock stopper has only one male tip
- combi stopper with male and female tip especially dedicated for securing Luer Lock fitting in infusion and transfusion sets
- single-use
- non-pyrogenic • packaging: 1 pc./paper-foil
- indication
- as a stopper of intravenous cannulas, ports and catheters
- intended to be used in combination with other devices during infusion therapy, protecting against the outflow of blood through the cannula and protecting against contamination and infection

| REF | type | intermediate packaging | bulk packaging |
|-------|-----------|---------------------------|-------------------|
| KDKLL | Luer-Lock | 250 pcs | 10 x 250 pcs |
| KDKC | combi | 100 pcs | 50 x 100 pcs |









OXYGEN THERAPY

nasal oxygen cannula

features

- made of medical PVC
- very soft tips with smooth edges

- very soft ups with smooth edges
 fastening under the chin with a soft, sliding adjustment ring
 various spacing of the nasal part of the cannula depending on the size
 universal connector compatible with any source of oxygen
 tubing with star cross-section that is resistant to bending
 available in 3 sizes: for adults, for children and for infants with 140 cm, 200 cm, 300 cm and 500 cm long tubing
- single-use
- packaging: 1 pc./foil

indication

designed to deliver oxygen through the nose

| REF | length of tubing | purpose | intermediate packaging | bulk packaging |
|------------|---------------------|--------------|---------------------------|-------------------|
| CTNN-200 | 200 cm | for infants | 2E pcc | 10 x 25 pcc |
| CTININ-200 | 200 Cm | TOF IMAMUS | 25 pcs | 10 x 25 pcs |
| CTNZ-200 | 200 cm | for children | 25 pcs | 10 x 25 pcs |
| CTND-140 | 140 cm | for adults | 30 pcs | 10 x 30 pcs |
| CTND-200 | 200 cm | for adults | 25 pcs | 10 x 25 pcs |
| CTND-300 | 300 cm | for adults | 15 pcs | 10 x 15 pcs |
| CTND-500 | 500 cm | for adults | 10 pcs | 10 x 10 pcs |





oxygen tubing

features

- made of transparent, non-toxic PVC
- tubing with star cross-section that is resistant to bending,
- with reinforcing straps along the entire length • standardized connectors at the end of the tubing allowing connection to oxygen devices
- single-use
- packaging: 1 pc./paper-foil

indication

 for transporting oxygen from an oxygen source to another medical device for oxygen supply, e.g. oxygen mask

| REF | length of tubing | bulk packaging |
|--------|---------------------|-------------------|
| DT2100 | 210 cm +/- 2% | 100 x 1 pc. |
| DT4260 | 426 cm +/- 2% | 50 x 1 pc. |
| DT7620 | 762 cm +/- 2% | 25 x 1 pc. |

oxygen mask

with tubing

features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
- equipped with a 210 cm (+/- 5%) long tubing with universal connectors
- tubing with star cross-section that is resistant to bending
- revolving connector that allows adjusting to patient's position • single-use
- packaging: 1 pc./paper-foil

indication

· designed to deliver oxygen with precise regulation of concentration, for medium oxygen concentration

| REF | length of tubing | size | |
|---------|---------------------|------|-----------------|
| MTS-S | 210 cm +/- 5% | ç | |
| 10113-3 | 210 011 17- 370 | 2 | \sim |
| MTS-M | 210 cm +/- 5% | Μ | |
| MTS-L | 210 cm +/- 5% | L | bulk packaging: |
| MTS-XL | 210 cm +/- 5% | XL | 50 x 1 pc. |









non-rebreather mask

with tubing

features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
 equipped with a 210 cm (+/- 5%) long tubing with universal connectors
- tubing with star cross-section that is resistant to bending
 1000 ml oxygen reservoir
- revolving connector that allows adjusting to patient's position • silicone valves on the connector and ventilation openings
- single-use
- packaging: 1 pc./paper-foil

indication

• designed to deliver oxygen with precise regulation of concentration, for high oxygen concentration

| REF | length of tubing | size | |
|--------|---------------------|------|-----------------|
| MTW-S | 210 cm +/- 5% | S | ~~ |
| MTW-M | 210 cm +/- 5% | Μ | |
| MTW-L | 210 cm +/- 5% | L | bulk packaging: |
| MTW-XL | 210 cm +/- 5% | XL | 50 x 1 pc. |



multi-vent mask

with tubing

features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
 equipped with a 210 cm (+/- 5%) long tubing with universal connectors
 tubing with star cross-section that is resistant to bending
- equipped with 2 jets white to set mid flow concentration (35%, 40%
- and 50%) and green to set low concentration (24%, 26%, 28% and 30%) single-use
- packaging: 1 pc./foil

indication

· designed to deliver oxygen with precise regulation of concentration, regardless of way of patient's breathing

| REF | length of tubing | size | | |
|---------|---------------------|------|-----------------|-----|
| MTVR-L | 210 cm +/- 5% | L | bulk packaging: | STE |
| MTVR-XL | 210 cm +/- 5% | XL | 100 x 1 pc. | SIL |

nebulizer mask

with tubing

features

- made of transparent, non-toxic PVC

- single-use
 packaging: 1 pc./paper-foil

indication

 for producing aerosol containing a comminuted medicine, that can be delivered directly to the respiratory system as the patient breathes

| REF | length of tubing | size |
|--------|---------------------|------|
| MTN-S | 210 cm +/- 5% | S |
| MTN-M | 210 cm +/- 5% | Μ |
| MTN-L | 210 cm +/- 5% | L |
| MTN-XL | 210 cm +/- 5% | XL |



- - adjustable nose metal plate and rubber fastening
 equipped with a 210 cm (+/- 5%) long tubing with universal connectors
 tubing with star cross-section that is resistant to bending
 - 8 ml nebulizer (graduated every 1 ml)
 average size of the molecule at a flow of 6l/min 3,23 μm



nebulizer

with TF/M adaptor, mouthpiece,

corrugated tube and tubing

features

- made of non-toxic materials
- 210 cm long, bend-resistant tubing with star-shaped cross-section • nebuliser with a 6ml capacity (graduated every 1ml), producing particles with an average size of MMAD = $3.91 \,\mu m$ at a flow rate of 6 l/min.
- nebuliser, tubing, T-adaptor, mouthpiece, one-way valve and corrugated tube included in a kit
- universal T-adaptor with a diameter of 22M/15F 22F enables the kit to be connected to endotracheal, tracheostomy tube or breathing circuit and perform nebulisation on mechanically ventilated patients
- the universal (detachable) mouthpiece allows nebulisation to be given to patients who are breathing on their own
- silicone one-way valve (detachable) reduces medication waste during nebulisation sinale-use

• packaging: 1 pc./foil

indication

- · device is designed to produce an aerosol containing fragmented medication
- which can be delivered deep into the respiratory system with the patient's inhalation • enables nebulisation to be carried out in mechanically ventilated patients by connecting directly to the endotracheal tube, tracheostomy tube or inspiratory arm of the ventilator breathing system
- · allows you to perform nebulisation in patient "on one's own breath"





ZARYS

ANAESTHESIA PRODUCTS

endotracheal tube

| uncuffed | RIB-20 | I.D. (mm) 2.0 | 20 pcs |
|---|--------|---------------|---------|
| | RIB-25 | I.D. (mm) 2.5 | 20 pcs |
| • features | RIB-30 | I.D. (mm) 3.0 | 20 pcs |
| made of flexible PVC | RIB-35 | I.D. (mm) 3.5 | 20 pcs |
| radiopaque line over the entire tube length equipped with lateral Murphy eye | RIB-40 | I.D. (mm) 4.0 | 20 pcs |
| graduated every 1 cm | RIB-45 | I.D. (mm) 4.5 | 20 pcs |
| double depth marker to control the position of the tube | RIB-50 | I.D. (mm) 5.0 | 20 pcs |
| • tube size marked in 3 places (on the connector and in | RIB-55 | I.D. (mm) 5.5 | 20 pcs |
| wo places on the tube body) | RIB-60 | I.D. (mm) 6.0 | 20 pcs |
| detachable connector with an outer diameter of 15 mm single-use | RIB-65 | I.D. (mm) 6.5 | 20 pcs |
| packaging: 1 pc./paper-foil | RIB-70 | I.D. (mm) 7.0 | 20 pcs |
| | RIB-75 | I.D. (mm) 7.5 | 20 pcs |
| • indication | RIB-80 | I.D. (mm) 8.0 | 20 pcs |
| for both oral and nasal intubation | RIB-85 | I.D. (mm) 8.5 | 20 pcs |
| inserted into the trachea to prevent airway obstruction and respiratory failure | RIB-90 | I.D. (mm) 9.0 | 20 pcs |
| obstruction and respiratory failure | | LD (mm) 0 F | 20 0.00 |

REF

RIB-95

RIB-100

REF

ETT20

ETT25

ETT30

ETT35

ETT40

ETT45

ETT50

ETT55

ETT60

ETT65

ETT70

ETT75

ETT80

ETT85

ETT90

ETT95

ETT100

REF

RIM-30

RIM-35

RIM-40

RIM-45

RIM-50

RIM-55

RIM-60

RIM-65

RIM-70

RIM-75

RIM-80

RIM-85

RIM-90

RIM-95

RIM-100

size

I.D. (mm) 3.0

I.D. (mm) 3.5

I.D. (mm) 4.0

I.D. (mm) 4.5

I.D. (mm) 5.0

I.D. (mm) 5.5

I.D. (mm) 6.0

I.D. (mm) 6.5

l.D. (mm) 7.0

I.D. (mm) 7.5

I.D. (mm) 8.0

I.D. (mm) 8.5

I.D. (mm) 9.0

I.D. (mm) 9.5

I.D. (mm) 10.0

NEV

size

I.D. (mm) 9.5

I.D. (mm) 10.0

size

I.D. (mm) 2.0

I.D. (mm) 2.5

I.D. (mm) 3.0

I.D. (mm) 3.5

I.D. (mm) 4.0

I.D. (mm) 4.5

I.D. (mm) 5.0

I.D. (mm) 5.5

packaging

20 pcs

20 pcs



uncuffed

STERILE EO

features

• made of thermoplastic flexible PVC

endotracheal tube

- · radiopaque line over the entire tube length
- equipped with lateral Murphy eye
- graduated every 1 cm
- · double depth marker to control the position of the tube – 2 rings
- tube size marked in three places (on the connector and in two places on the body of the tube)
- name of the manufacturer affixed to the product · detachable connector with an external diameter
- of 15 mm
- banana-type packaging, maintaining the anatomical shape of the tube
- single-use
- packing: 1 pc / paper-foil

indication

- for both oral and nasal intubation
- inserted into the trachea to prevent airway obstruction and respiratory failure



endotracheal tube

cuffed


features

- made of flexible PVC
- equipped with lateral Murphy eye and radiopaque line
- over the entire tube length
- graduated every 1 cm, double depth marker to control
- the position of the tube • low-pressure, high-volume cuff
- tube size marked in 4 places (on the connector, on the
- pilot balloon and in two places on the tube body)
- detachable connector with an outer diameter of 15 mm
- pilot balloon with LOT and size indication
- single-use
- packaging: 1 pc./paper-foil

indication

- for both oral and nasal intubation
- inserted into the trachea to prevent airway obstruction and respiratory failure





endotracheal tube cuffed

features

- made of thermoplastic flexible PVC
- radiopaque line over the entire tube length
- equipped with lateral Murphy eye
- graduated every 1 cm
- double depth marker to control the position of the tube 2 rings
- · equipped with a low-pressure, high-volume sealing cuff
- tube size marked in three places (on the connector and in two places on the body of the tube)

NEV

- name of the manufacturer affixed to the product
- control balloon, marked with batch and size number
- detachable connector with an external diameter of 15 mm
- banana-type packaging, maintaining the anatomical shape of the tube sinale-use
- packing: 1 pc / paper-foil

indication

- for both oral and nasal intubation
- inserted into the trachea to prevent airway obstruction and respiratory failure

| REF | size | intermediate packaging |
|---------|----------------|---------------------------|
| ETTC25 | I.D. (mm) 2.5 | 10 pcs |
| ETTC30 | I.D. (mm) 3.0 | 10 pcs |
| ETTC35 | I.D. (mm) 3.5 | 10 pcs |
| ETTC40 | I.D. (mm) 4.0 | 10 pcs |
| ETTC45 | I.D. (mm) 4.5 | 10 pcs |
| ETTC50 | I.D. (mm) 5.0 | 10 pcs |
| ETTC55 | I.D. (mm) 5.5 | 10 pcs |
| ETTC60 | I.D. (mm) 6.0 | 10 pcs |
| ETTC65 | I.D. (mm) 6.5 | 10 pcs |
| ETTC70 | I.D. (mm) 7.0 | 10 pcs |
| ETTC75 | I.D. (mm) 7.5 | 10 pcs |
| ETTC80 | I.D. (mm) 8.0 | 10 pcs |
| ETTC85 | I.D. (mm) 8.5 | 10 pcs |
| ETTC90 | I.D. (mm) 9.0 | 10 pcs |
| ETTC95 | I.D. (mm) 9.5 | 10 pcs |
| ETTC100 | I.D. (mm) 10.0 | 10 pcs |



reinforced endotracheal tube

cuffed

features

- made of flexible PVC
- reinforcing within the wall of the entire tube
- · equipped with lateral Murphy eye and radiopaque line over the entire tube length
- graduated every 1 cm, double depth marker to control the position of the tube • tube size marked in 4 places (on the connector, on the pilot balloon and in two
- places on the tube body) • indissiociable connector with an outer diameter of 15 mm
- single-use
- packaging: 1 pc./paper-foil

indication

- for both oral and nasal intubation
- inserted into the trachea to prevent airway obstruction and respiratory failure
- for intubation of the patient in an unusual position, with a high risk of kinking the endotracheal tube

| | REF | size | intermediate packaging |
|--|-----------|----------------|---------------------------|
| | RIMZP-30 | I.D. (mm) 3.0 | 10 pcs |
| | RIMZP-35 | I.D. (mm) 3.5 | 10 pcs |
| | RIMZP-40 | I.D. (mm) 4.0 | 10 pcs |
| | RIMZP-45 | I.D. (mm) 4.5 | 10 pcs |
| | RIMZP-50 | I.D. (mm) 5.0 | 10 pcs |
| | RIMZP-55 | I.D. (mm) 5.5 | 10 pcs |
| | RIMZP-60 | I.D. (mm) 6.0 | 10 pcs |
| | RIMZP-65 | I.D. (mm) 6.5 | 10 pcs |
| | RIMZP-70 | I.D. (mm) 7.0 | 10 pcs |
| | RIMZP-75 | I.D. (mm) 7.5 | 10 pcs |
| | RIMZP-80 | I.D. (mm) 8.0 | 10 pcs |
| | RIMZP-85 | I.D. (mm) 8.5 | 10 pcs |
| | RIMZP-90 | I.D. (mm) 9.0 | 10 pcs |
| | RIMZP-95 | I.D. (mm) 9.5 | 10 pcs |
| | RIMZP-100 | I.D. (mm) 10.0 | 10 pcs |



Ø bulk packaging 10 x 20 pcs

24 · 24 · 24 · 24 · 107.2

10 pcs 10 pcs 10 pcs

packaging

20 pcs

I.D. (mm) 6.0 I.D. (mm) 6.5 10 pcs I.D. (mm) 7.0 10 pcs I.D. (mm) 7.5 10 pcs I.D. (mm) 8.0 10 pcs I.D. (mm) 8.5 10 pcs I.D. (mm) 9.0 I.D. (mm) 9.5 I.D. (mm) 10.0 10 x 10 pcs













intubation stylet

features

- made of PVC-coated aluminium
- · special material for both stiffness and pliability
- intended for use during difficult intubation
- size and diameter markings
- soft distal tip
- single-use
- packaging: 1 pc./paper-foil

indication

- allows you to give the endotracheal tube the right curvature • introduced to the endotracheal tube facilitates intubation in routine
- and difficult cases, when the laryngeal inlet is not visible

| REF | size | diameter | length | tube size | intermediate packaging |
|--------|-------|----------|--------------|------------------|---------------------------|
| PRI-06 | 6 Fr | 2,0 mm | 310 mm (±10) | from 2,0 to 4,0 | 20 pcs |
| PRI-10 | 10 Fr | 3,3 mm | 390 mm (±10) | from 4,5 to 6,5 | 20 pcs |
| PRI-14 | 14 Fr | 4,7 mm | 390 mm (±10) | from 7,0 to 10,0 | 20 pcs |







curved

features

- made of medical grade low density polyethylene
- · curved and rounded tip enables safe intubation without
- damaging the soft tissues · optimum flexibility and rigidity facilitate correct insertion
- and safe removal of the introducer after intubation • imprint of the size, diameter and manufacturer's name
- on the device · bright yellow body for increased visibility of the introducer in the tracheal tube canal
- clear, easy to read graduation every 1 cm
- available in three sizes and lengths
- single-use
- PVC free
- packaging: 1 pc./paper-foil

indication

52


www.zarys.com

- for use during difficult intubation procedures
- as an aid to facilitate the insertion of the endotracheal tube during standard intubation

| REF | size | diameter O.D. | total length | tube size | graduation | intermediate packaging |
|-------------|-------|------------------|-----------------|------------------|------------|---------------------------|
| PTIZ-06-535 | 6 Fr | 2,0 mm | 535 mm | from 2,0 to 4,0 | 10 – 45 mm | 10 pcs |
| PTIZ-10-600 | 10 Fr | 3,3 mm | 600 mm | from 4,5 to 6,0 | 10 – 50 mm | 10 pcs |
| PTIZ-10-800 | 10 Fr | 3,3 mm | 800 mm | from 4,5 to 6,0 | 10 – 60 mm | 10 pcs |
| PTIZ-15-600 | 15 Fr | 5,0 mm | 600 mm | from 6,5 to 10,0 | 10 – 50 mm | 10 pcs |
| PTIZ-15-800 | 15 Fr | 5,0 mm | 800 mm | from 6,5 to 10,0 | 10 – 60 mm | 10 pcs |

 \bigotimes bulk packaging: 20 x 10 pcs





oropharyngeal airway/guedel airway

- made of high quality PE (polyethylene), free of PVC
- distinctive anatomical curve of airway tube allows for the opening
- of the upper respiratory tract and the free flow of anaesthetic gases, air colour coded sizes
- numeric marking of the size on the airway
- single-use
- packaging: 1 pc./foil

indication

- used during artificial ventilation procedures to ensure the airway patency as a tool for paramedic to assist in the management of a casualty's
- airway which increases the effectiveness of the rescue operation
- facilitates the stabilization of the endotracheal tube and protects it against biting

| REF | size | number | colour | intermediate packaging |
|--------|--------|--------|-----------|---------------------------|
| RG-40 | 40 mm | 000 | pink | 50 pcs |
| RG-50 | 50 mm | 00 | dark blue | 50 pcs |
| RG-60 | 60 mm | 0 | black | 50 pcs |
| RG-70 | 70 mm | 1 | white | 50 pcs |
| RG-80 | 80 mm | 2 | green | 50 pcs |
| RG-90 | 90 mm | 3 | yellow | 50 pcs |
| RG-100 | 100 mm | 4 | red | 50 pcs |
| RG-110 | 110 mm | 5 | blue | 50 pcs |
| RG-120 | 120 mm | 6 | orange | 50 pcs |
| | | | | |



laryngeal mask

PVC

features

• made of medical PVC

- elastic mask body with description (size, cuff volume, suggested patient's weight)
- elliptic shape, flexible sealing cuff and universal, transparent connector: 15 mm colour-coded valve protection makes it easy to quickly recognize size
- single-use
- compatible with magnetic resonance (MR)
 no metal elements
- packaging: 1 pc./paper-foil

indication

- device inserted through the mouth, up to the back of the throat, forms a seal at the top of the glottis in order to:
- direct protecting the patient's airway patency
- artificial ventilation also during cardiopulmonary resuscitation maintaining airway patency in: short-term surgical procedures,
- without the need of muscles relaxation, when endotracheal intubation is not necessary and under anesthesia for diagnostic procedures

| REF | size | colour of valve protection | weight of patient | max. inflation | intermediate packaging |
|-------------|------|----------------------------|----------------------|-------------------|---------------------------|
| MKPJ-10 | 1 | pink | ≤ 5 kg | 4 ml | 5 pcs |
| MKPJ-15 | 1.5 | orange | 5-10 kg | 7 ml | 5 pcs |
| MKPJ-20 2 | 2 | dark green | 10-20 kg | 10 ml | 5 pcs |
| MKPJ-25 2.5 | 2.5 | brown | 20-30 kg | 14 ml | 5 pcs |
| MKPJ-30 | 3 | green | 30-50 kg | 20 ml | 5 pcs |
| MKPJ-40 | 4 | purple | 50-70 kg | 30 ml | 5 pcs |
| MKPJ-50 | 5 | dark blue | 70-100 kg | 40 ml | 5 pcs |
| | | | | | |











laryngeal mask

silicone

features

- made of biocompatible silicone, which ensures excellent flexibility of the device • elastic mask body with description (size, cuff volume, suggested patient's weight) • elliptic shape, flexible sealing cuff and universal, transparent connector: 15 mm
- single-use
- compatible with magnetic resonance (MR) no metal elements
- packaging: 1 pc./paper-foil

indication

- device inserted through the mouth, up to the back of the throat, forms a seal at the top of the glottis in order to:
- direct protecting the patient's airway patency
- artificial ventilation also during cardiopulmonary resuscitation - maintaining airway patency in: short-term surgical procedures, without the need of muscles relaxation, when endotracheal intubation is not necessary and under anesthesia for diagnostic procedures

| REF | size | weight of patient | max. inflation | intermediate packaging |
|---------|------|----------------------|-------------------|---------------------------|
| MKSJ-10 | 1 | ≤ 5 kg | 4 ml | 5 pcs |
| MKSJ-15 | 1.5 | 5-10 kg | 7 ml | 5 pcs |
| MKSJ-20 | 2 | 10-20 kg | 10 ml | 5 pcs |
| MKSJ-25 | 2.5 | 20-30 kg | 14 ml | 5 pcs |
| MKSJ-30 | 3 | 30-50 kg | 20 ml | 5 pcs |
| MKSJ-40 | 4 | 50-70 kg | 30 ml | 5 pcs |
| MKSJ-50 | 5 | 70-100 kg | 40 ml | 5 pcs |
| | | | | |

flexible neckplates with slits for strap insertion to secure tube position

• pilot balloon in contrasting colour clearly indicates the level of the cuff's



tracheostomy tube

• made of medical, flexible PVC • swivel 15 mm connector

• soft, low-pressure, thin-walled cuff

• 2 fixation straps protected with foil

clear size marking on the neckplates

• atraumatic, soft rounded tip, curved cannula

• to control respiratory tracts in people after tracheostomy • introduced to the airways by a surgical incision of the tracheal's

• provides free breathing bypassing the nose, throat and larynx

intermediate

packaging

10 pcs

10 pcs

10 pcs

10 pcs

10 pcs

cuffed

features

inflation

radiopaque line

single-use

indication

REF

RTM-40

RTM-45

RTM-50

RTM-60

RTM-65

RTM-70

RTM-75

RTM-80

RTM-85

RTM-90

RTM-95

RTM-100

• equipped with obturator

• packaging: 1 pc./rigid blister

front wall to allow ventilation

size

4.0

4.5

5.0

6.0

6.5

7.0

7.5

8.0

8.5

9.0

9.5

10.0

bulk packaging 10 x 5 pcs



features

- due to the high quality of materials used, ventilation is conducted in an atraumatic and safe way for the patient
- anatomically shaped inflatable cuff, ensures perfect tightness of the mask with minimal pressure
- device, when applied to the face, tightly covers the patient's mouth and nose without compressing the eyes
- clear shell provides continuous visual observation of the patient
- · colour-coded retaining rings, easy for size recognition
- retaining ring can be removed when it is not needed • full size range allows to properly choose the size of the mask depending on the age of the patient
- · all sizes are equipped with inflatable check valve to adjust air amount for best cushion fitting auto patient's face
- single-use
- packaging: 1 pc./foil

indication

- intended for the administration of a gas mixture coming through devices for artificial ventilation
- enables performing medical procedures including artificial ventilation, administration of anesthetic gas and carrying out emergency treatment of an unconscious patient

| REF | size | age of patient/ size | colour | bulk packaging |
|-------|------|-------------------------|--------|-------------------|
| MAP-0 | 0 | neonate | white | 200 x 1 pc. |
| MAP-1 | 1 | infant | pink | 200 x 1 pc. |
| MAP-2 | 2 | pediatric | yellow | 150 x 1 pc. |
| MAP-3 | 3 | adult S | green | 150 x 1 pc. |
| MAP-4 | 4 | adult M | red | 100 x 1 pc. |
| MAP-5 | 5 | adult L | blue | 100 x 1 pc. |
| MAP-6 | 6 | adult XL | orange | 80 x 1 pc. |

NON DATEX STERLE

anaesthetic mask with open seal

- due to the high quality of materials used, ventilation is conducted in an atraumatic an safe way for the patient
- anatomically shaped open seal made of TPE ensures perfect tightness of the mask with minimal pressure
- polypropylene, transparent mask dome allows continuous visual observation of the patient
- · device, when applied to the face, tightly covers the patient's mouth and nose without compressing the eyes
- size embossed on the body and colour-coded seal
- full size range allows to properly choose the size of the mask depending on the age of the patient
- sinale-use
- packaging: 1 pc./medical foil

indication

- device intended for the administration of a gas mixture coming through devices for artificial ventilation
- enables performing medical procedures including artificial ventilation, administration of anesthetic gas and carrying out emergency treatment of an unconscious patient
- · device used in hospitals (operating theaters, intensive care units, emergency departments) as well as during pre-hospital care by rescue teams carrying out rescue operations in an unconscious patient

| catalogue number | size | age of patient | seal colour | connector | bulk pac |
|---------------------|------|------------------|----------------|-----------|-------------|
| MAO-1 | 1 | infant | blue | 15 mm | 200 x 1 pc |
| MAO-2 | 2 | pediatric | orange | 15 mm | 150 x 1 pc |
| MAO-3 | 3 | small adult | pink | 22 mm | 150 x 1 pc |
| MAO-4 | 4 | medium adult | white | 22 mm | 100 x 1 pc |
| MAO-5 | 5 | large adult | yellow | 22 mm | 100 x 1 pc |
| MAO-6 | 6 | very large adult | green | 22 mm | 80 x 1 pc. |













ckaging C. C. C.

catheter mount with double swivel elbow connector

features

- smooth bore type made of medical PVC
- expandable type made of high quality PE, EVA
- universal patient connection: 22M/15F • ventilation device side connector: 15M or 22F
- expandable version characterized by shape memorizing
- and by the ability to control the position of the connector double swivel elbow connector
- · port for suctioning and bronchoscope secured with a stopper
- available in versions: smooth bore and expandable
- single-use • packaging: 1 pc./paper-foil

indication


- used to connect the tracheostomy tube or endotracheal tube with the breathing circuit of the ventilation device provides greater mobility and rotation of the connection.
- increasing patient comfort and safety



ZARYS

| REF | type | connection | length (without connectors) | intermediate packaging |
|----------|------------------------------|-------------|--------------------------------|---------------------------|
| M-15-G | smooth bore (not expandable) | 22M/15F-15M | 15 cm | 10 pcs |
| PM-22-G | smooth bore (not expandable) | 22M/15F-22F | 15 cm | 10 pcs |
| PMR-15-K | expandable | 22M/15F-15M | 8-15 cm | 10 pcs |
| PMR-22-K | expandable | 22M/15F-22F | 8-15 cm | 10 pcs |





dual-limb breathing circuit

features


- made of high quality EVA material, without PVC
- 15F/22M and 22F universal connector fits compatible products
- and devices for mechanical ventilation or anesthesia machine
- · expandable arms with maximum length of 160 cm
- · adjustable length of circuit enabling adjustment of length of respiratory system and optimization of dead space
- equipped with latex-free breathing bag: 3I (for adults) or 1I (for children)
- additional expandable tube, length: 100 cm
- elbow with Luer-Lock port
- single-use
- packaging: 1 pc./foil

indication

• for patient's mechanical ventilation • it enables transporting breathing mixture from equipment (ventilator, anesthesia machine) to respiratory protection device (e.g. tracheal tube, tracheostomy mask) and vice versa

| REF | purpose | capacity of bag | bulk packaging |
|------------|--------------|--------------------|-------------------|
| OBWD-1601L | for children | 1 liter | 50x1 pc |
| OBWD-1603L | for adults | 3 liters | 40x1 pc |














SUCTION CATHETERS, **STOMACH TUBES**

suction catheter with vacuum control

features

- made of soft and flexible PVC, that is resistant to bending and kinking
- frozen surface prevents the catheter from sticking to the wall of the
- endotracheal tube
- full control over the suction
- atraumatic, soft rounded, open tip • tip with two alternating lateral eyes with smooth edges
- colour-coded connectors for size recognition, numerical size identification
- on the packaging • single-use
- packaging: 1 pc./paper-foil

indication

• device designed for suction of secretions from the upper respiratory tract

| REF | size | length | colour | intermediate packaging |
|-----------|-------|--------|-------------|---------------------------|
| CO-A06-60 | CH 6 | 600 mm | light green | 50 pcs |
| CO-A08-60 | CH 8 | 600 mm | blue | 50 pcs |
| CO-A10-60 | CH 10 | 600 mm | black | 50 pcs |
| CO-A12-60 | CH 12 | 600 mm | white | 50 pcs |
| CO-A14-60 | CH 14 | 600 mm | green | 50 pcs |
| CO-A16-60 | CH 16 | 600 mm | orange | 50 pcs |
| CO-A18-60 | CH 18 | 600 mm | red | 50 pcs |



\bigotimes bulk packaging: 10 x 50 pcs

Ø bulk packaging: 10 x 50 pcs





suction catheter

features

- made of soft and flexible PVC, that is resistant to bending and kinking • frozen surface prevents the catheter from sticking to the wall of the
- endotracheal tube
- atraumatic, soft rounded, open tip
- tip with two alternating lateral eyes with smooth edges
- colour-coded connectors for size recognition, numerical size identification on the packaging
- single-use
- packaging: 1 pc./paper-foil

indication

· device designed for suction of secretions from the upper respiratory tract ad for agurea of quations

| recommended for | source or s | uction with | its own suctio | n force regulation |
|-------------------------------------|-------------|-------------|----------------|--------------------|
| | | | | |

| REF | size | length | colour | intermediate packaging |
|-----------|-------|--------|-------------|---------------------------|
| CO-B06-40 | CH 6 | 400 mm | light green | 50 pcs |
| CO-B08-40 | CH 8 | 400 mm | blue | 50 pcs |
| CO-B10-40 | CH 10 | 400 mm | black | 50 pcs |
| CO-B12-50 | CH 12 | 500 mm | white | 50 pcs |
| CO-B12-60 | CH 12 | 600 mm | white | 50 pcs |
| CO-B14-50 | CH 14 | 500 mm | green | 50 pcs |
| CO-B14-60 | CH 14 | 600 mm | green | 50 pcs |
| CO-B16-50 | CH 16 | 500 mm | orange | 50 pcs |
| CO-B16-60 | CH 16 | 600 mm | orange | 50 pcs |
| CO-B18-50 | CH 18 | 500 mm | red | 50 pcs |
| CO-B18-60 | CH 18 | 600 mm | red | 50 pcs |
| CO-B20-50 | CH 20 | 500 mm | yellow | 50 pcs |
| CO-B20-60 | CH 20 | 600 mm | yellow | 50 pcs |
| CO-B22-50 | CH 22 | 500 mm | purple | 50 pcs |
| CO-B22-60 | CH 22 | 600 mm | purple | 50 pcs |
| | | | | |



suction catheter with frozen surface

features

- made of soft and flexible PVC, that is resistant to bending and kinking
- frozen surface enabling easy insertion
- atraumatic, soft rounded, open tip
- tip with two opposite lateral eyes with smooth edges • colour-coded connectors for size recognition, numerical size identification on the packaging
- single-use
- packaging: 1 pc./paper-foil

indication

- · device designed for suction of secretions from the upper respiratory tract
- recommended for source of suction with its own suction force regulation

| REF | size | length | colour | intermediate packaging |
|-----------------|-------|--------|-------------|---------------------------|
| CO-BF-ZP-B06-40 | CH 6 | 400 mm | light green | 50 pcs |
| CO-BF-ZP-B06-60 | CH 6 | 600 mm | light green | 50 pcs |
| CO-BF-ZP-B08-40 | CH 8 | 400 mm | blue | 50 pcs |
| CO-BF-ZP-B08-60 | CH 8 | 600 mm | blue | 50 pcs |
| CO-BF-ZP-B10-40 | CH 10 | 400 mm | black | 50 pcs |
| CO-BF-ZP-B10-60 | CH 10 | 600 mm | black | 50 pcs |
| CO-BF-ZP-B12-60 | CH 12 | 600 mm | white | 50 pcs |
| CO-BF-ZP-B14-60 | CH 14 | 600 mm | green | 50 pcs |
| CO-BF-ZP-B16-60 | CH 16 | 600 mm | orange | 50 pcs |
| CO-BF-ZP-B18-60 | CH 18 | 600 mm | red | 50 pcs |
| CO-BF-ZP-B20-60 | CH 20 | 600 mm | yellow | 50 pcs |
| CO-BF-ZP-B22-60 | CH 22 | 600 mm | purple | 50 pcs |
| CO-BF-ZP-B24-60 | CH 24 | 600 mm | dark blue | 50 pcs |
| | | | | |



feeding tube graduated, with X-ray contrast

features

- made of soft PVC with a frosted surface, resistant to bending and kinking due to the appropriate degree of hardness • atraumatic, soft rounded, closed tip facilitates precise insertion of catheter

 - two alternate end holes with gently stretched edges
 - connector compatible with the conical tip of Luer syringe
 - integrated stoper allows the tube to be closed tightly • colour-coded connectors for size recognition, numerical size identification
 - on packaging
 - radiopaque line along the entire length of the catheter
 digital depth graduation every 1 cm (from 5 cm to 40 cm)
- length 50 cm
- name of the manufacturer placed on the catheter
- individually packed straight
- single-use
- packaging: 1 pc. / paper-foil

indication

· device inserted through patient's nose or mouth to perform enteral nutrition

| REF | size | length | colour | intermediate packaging |
|------------|-------|-----------|------------|---------------------------|
| | | | | |
| CKS-04-50 | CH 4 | 500 mm | red | 100 pcs |
| CKS-05-50 | CH 5 | 500 mm | grey | 100 pcs |
| CKS-06-50 | CH 6 | 500 mm | pink | 100 pcs |
| CI(3-00-30 | CITO | 500 11111 | ршк | roo pes |
| CKS-08-50 | CH 8 | 500 mm | pale green | 100 pcs |
| CKS-10-50 | CH 10 | 500 mm | black | 100 pcs |
| | | | | |











stomach tube/duodenal tube

features

- made of PVC
- made of material that is resistant to bending and kinking due to the appropriate degree of hardness
- atraumatic, soft rounded closed tip
- four lateral eyes with smooth edges colour-coded connectors for size recognition, numerical
- size identification on packaging • stomach tube is available in 3 lengths: 80 cm, 100 cm, 125 cm and duodenal tube is available in 150 cm length
- single-use
- packaging: 1 pc./paper-foil

indication

- inserted through nose or mouth into stomach or duodenum
- drainage of stomach/intestinal content decompression of stomach (air or food evacuation)
- gastric lavage
- taking samples of stomach content
- · control of bleeding from the upper section of digestive tract
- evaluation of secretory function • administration of medicines, fluids and/or food content into stomach or duodenum

| REF | size | length | colour | intermediate packaging | bulk packaging |
|---------------|-------|---------|-------------|---------------------------|-------------------|
| * ZZ-80/CH05 | CH 5 | 800 mm | grey | 35 pcs | 20 x 35 pcs |
| * ZZ-80/CH06 | CH 6 | 800 mm | light green | 35 pcs | 20 x 35 pcs |
| ZZ-80/CH08 | CH 8 | 800 mm | blue | 35 pcs | 20 x 35 pcs |
| ZZ-80/CH10 | CH 10 | 800 mm | black | 35 pcs | 20 x 35 pcs |
| ZZ-80/CH12 | CH 12 | 800 mm | white | 35 pcs | 20 x 35 pcs |
| ZZ-80/CH14 | CH 14 | 800 mm | green | 30 pcs | 20 x 30 pcs |
| ZZ-80/CH16 | CH 16 | 800 mm | orange | 30 pcs | 20 x 30 pcs |
| ZZ-80/CH18 | CH 18 | 800 mm | red | 25 pcs | 20 x 25 pcs |
| ZZ-80/CH20 | CH 20 | 800 mm | yellow | 15 pcs | 20 x 15 pcs |
| ZZ-80/CH22 | CH 22 | 800 mm | purple | 15 pcs | 20 x 15 pcs |
| ZZ-80/CH24 | CH 24 | 800 mm | dark blue | 15 pcs | 20 x 15 pcs |
| * ZZ-100/CH05 | CH 5 | 1000 mm | grey | 30 pcs | 20 x 30 pcs |
| * ZZ-100/CH06 | CH 6 | 1000 mm | light green | 30 pcs | 20 x 30 pcs |
| ZZ-100/CH08 | CH 8 | 1000 mm | blue | 30 pcs | 20 x 30 pcs |
| ZZ-100/CH10 | CH 10 | 1000 mm | black | 30 pcs | 20 x 30 pcs |
| ZZ-100/CH12 | CH 12 | 1000 mm | white | 30 pcs | 20 x 30 pcs |
| ZZ-100/CH14 | CH 14 | 1000 mm | green | 25 pcs | 20 x 25 pcs |
| ZZ-100/CH16 | CH 16 | 1000 mm | orange | 20 pcs | 20 x 20 pcs |
| ZZ-100/CH18 | CH 18 | 1000 mm | red | 15 pcs | 20 x 15 pcs |
| ZZ-100/CH20 | CH 20 | 1000 mm | yellow | 15 pcs | 20 x 15 pcs |
| ZZ-100/CH22 | CH 22 | 1000 mm | purple | 10 pcs | 20 x 10 pcs |
| ZZ-100/CH24 | CH 24 | 1000 mm | dark blue | 10 pcs | 20 x 10 pcs |
| * ZZ-125/CH05 | CH 5 | 1250 mm | grey | 30 pcs | 20 x 30 pcs |
| * ZZ-125/CH06 | CH 6 | 1250 mm | light green | 30 pcs | 20 x 30 pcs |
| ZZ-125/CH08 | CH 8 | 1250 mm | blue | 30 pcs | 20 x 30 pcs |
| ZZ-125/CH10 | CH 10 | 1250 mm | black | 30 pcs | 20 x 30 pcs |
| ZZ-125/CH12 | CH 12 | 1250 mm | white | 20 pcs | 20 x 20 pcs |
| ZZ-125/CH14 | CH 14 | 1250 mm | green | 25 pcs | 20 x 25 pcs |
| ZZ-125/CH16 | CH 16 | 1250 mm | orange | 20 pcs | 20 x 20 pcs |
| ZZ-125/CH18 | CH 18 | 1250 mm | red | 20 pcs | 20 x 20 pcs |
| ZZ-125/CH20 | CH 20 | 1250 mm | yellow | 15 pcs | 20 x 15 pcs |
| ZZ-125/CH22 | CH 22 | 1250 mm | purple | 10 pcs | 20 x 10 pcs |
| ZZ-125/CH24 | CH 24 | 1250 mm | dark blue | 10 pcs | 20 x 10 pcs |
| | | | | | |
| * ZD-150/CH5 | CH 5 | 1500 mm | grey | 30 pcs | 20 x 30 pcs |
| * ZD-150/CH6 | CH 6 | 1500 mm | light green | 30 pcs | 20 x 30 pcs |
| * ZD-150/CH8 | CH 8 | 1500 mm | blue | 30 pcs | 20 x 30 pcs |
| * ZD-150/CH10 | CH 10 | 1500 mm | black | 30 pcs | 20 x 30 pcs |
| * ZD-150/CH12 | CH 12 | 1500 mm | white | 30 pcs | 20 x 30 pcs |
| * ZD-150/CH14 | CH 14 | 1500 mm | green | 20 pcs | 20 x 20 pcs |
| * ZD-150/CH16 | CH 16 | 1500 mm | orange | 15 pcs | 20 x 15 pcs |
| * ZD-150/CH18 | CH 18 | 1500 mm | red | 15 pcs | 20 x 15 pcs |
| ZD-150/CH20 | CH 20 | 1500 mm | yellow | 10 pcs | 20 x 10 pcs |
| ZD-150/CH22 | CH 22 | 1500 mm | purple | 10 pcs | 20 x 10 pcs |
| ZD-150/CH24 | CH 24 | 1500 mm | dark blue | 10 pcs | 20 x 10 pcs |
| | | | | | |




2000

> *< available

on request

ZARYS

60 www.zarys.com



UROLOGICAL CATHETERS, **URINE BAGS**

Foley catheter

2-way, silicone-coated latex

features

- made of soft and flexible latex, that is resistant to bending and kinking • atraumatic, soft rounded, special tip facilitating insertion
- of the catheter to the bladder
- tip with two lateral eyes with smooth edges
 universal urine drainage port of the catheter fits into standard connectors in urine bags
- silicone-coated surface that allows atraumatic insertion of the catheter
 easy to fill and to empty, tear-resistant balloon of various capacities
- Plastic or rubber valve to ensure balloon tightness
 CH 6, CH 8, CH 10 size has a guide to facilitate catheterization
- sinale-use
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil

indication

• to drain urine from the bladder

with plastic valve

| REF | length | balloon | size | colour | intermediate packaging |
|-------------------|--------|----------|-------|-------------|---------------------------|
| CFNL-2D-06-03-P | 270 mm | 3 ml | CH 6 | light green | 10 pcs |
| CFNL-2D-08-05-P | 270 mm | 3-5 ml | CH 8 | blue | 10 pcs |
| CFNL-2D-10-05-P | 270 mm | 3-5 ml | CH 10 | black | 10 pcs |
| CFNL-2D-12-10-P | 400 mm | 5-10 ml | CH 12 | white | 10 pcs |
| CFNL-2D-14-10-P | 400 mm | 5-10 ml | CH 14 | green | 10 pcs |
| CFNL-2D-16-10-P | 400 mm | 5-10 ml | CH 16 | orange | 10 pcs |
| CFNL-2D-18-10-P | 400 mm | 5-10 ml | CH 18 | red | 10 pcs |
| CFNL-2D-20-10-P | 400 mm | 5-10 ml | CH 20 | yellow | 10 pcs |
| CFNL-2D-22-10-P | 400 mm | 5-10 ml | CH 22 | purple | 10 pcs |
| CFNL-2D-24-10-P | 400 mm | 5-10 ml | CH 24 | dark blue | 10 pcs |
| CFNL-2D-12-30-P | 400 mm | 30 ml | CH 12 | white | 10 pcs |
| * CFNL-2D-26-30-P | 400 mm | 30 ml | CH 26 | pink | 10 pcs |
| CFNL-2D-28-30-P | 400 mm | 30 ml | CH 28 | brown | 10 pcs |
| * CFNL-2D-30-30-P | 400 mm | 30 ml | CH 30 | grey | 10 pcs |
| CFNL-2D-14-50-P | 400 mm | 30-50 ml | CH 14 | green | 10 pcs |
| CFNL-2D-16-50-P | 400 mm | 30-50 ml | CH 16 | orange | 10 pcs |
| CFNL-2D-18-50-P | 400 mm | 30-50 ml | CH 18 | red | 10 pcs |
| CFNL-2D-20-50-P | 400 mm | 30-50 ml | CH 20 | yellow | 10 szt. |
| CFNL-2D-22-50-P | 400 mm | 30-50 ml | CH 22 | purple | 10 pcs |
| CFNL-2D-24-50-P | 400 mm | 30-50 ml | CH 24 | dark blue | 10 pcs |
| CFNL-2D-26-50-P | 400 mm | 30-50 ml | CH 26 | pink | 10 pcs |
| CFNL-2D-28-50-P | 400 mm | 30-50 ml | CH 28 | brown | 10 pcs |
| CFNL-2D-30-50-P | 400 mm | 30-50 ml | CH 30 | grey | 10 pcs |
| | | | | | |

STERILEEO LATEX STERILE

with rubber valve

| REF | length | balloon | size | colour | intermediate packaging |
|-----------------|--------|---------|-------|-------------|---------------------------|
| CFNL-2D-06-03-L | 270 mm | 3 ml | CH 6 | light green | 10 pcs |
| CFNL-2D-08-05-L | 270 mm | 3-5 ml | CH 8 | blue | 10 pcs |
| CFNL-2D-10-05-L | 270 mm | 3-5 ml | CH 10 | black | 10 pcs |
| CFNL-2D-12-10-L | 400 mm | 5-10 ml | CH 12 | white | 10 pcs |
| CFNL-2D-14-10-L | 400 mm | 5-10 ml | CH 14 | green | 10 pcs |
| CFNL-2D-16-10-L | 400 mm | 5-10 ml | CH 16 | orange | 10 pcs |
| CFNL-2D-18-10-L | 400 mm | 5-10 ml | CH 18 | red | 10 pcs |
| CFNL-2D-20-10-L | 400 mm | 5-10 ml | CH 20 | yellow | 10 pcs |
| CFNL-2D-22-10-L | 400 mm | 5-10 ml | CH 22 | purple | 10 pcs |
| CFNL-2D-24-10-L | 400 mm | 5-10 ml | CH 24 | dark blue | 10 pcs |
| CFNL-2D-26-10-L | 400 mm | 5-10 ml | CH 26 | pink | 10 pcs |
| CFNL-2D-16-30-L | 400 mm | 30 ml | CH 16 | orange | 10 pcs |
| CFNL-2D-18-30-L | 400 mm | 30 ml | CH 18 | red | 10 pcs |
| CFNL-2D-20-30-L | 400 mm | 30 ml | CH 20 | yellow | 10 pcs |
| CFNL-2D-22-30-L | 400 mm | 30 ml | CH 22 | purple | 10 pcs |
| CFNL-2D-24-30-L | 400 mm | 30 ml | CH 24 | dark blue | 10 pcs |
| CFNL-2D-26-30-L | 400 mm | 30 ml | CH 26 | pink | 10 pcs |
| | | | | | |
















Foley catheter

3-way, silicone-coated latex

features

- made of soft and flexible latex, that is resistant to bending and kinking • universal urine drainage port of the catheter fits into standard connectors
- in urine bags
- additional inlet for flushing out the bladder
 atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
 silicone-coated surface that allows atraumatic insertion of the catheter
- easy to fill and to empty, tear-resistant balloon
 plastic valve to ensure balloon tightness
- single-use
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil indication
- to drain urine from the bladder

| REF | length | balloon | size | colour | intern packa |
|-------------------|--------|---------|-------|-----------|-----------------|
| CFNL-3D-16-30-P | 400 mm | 30 ml | CH 16 | orange | 10 pcs |
| CFNL-3D-18-30-P | 400 mm | 30 ml | CH 18 | red | 10 pcs |
| CFNL-3D-20-30-P | 400 mm | 30 ml | CH 20 | yellow | 10 pcs |
| CFNL-3D-22-30-P | 400 mm | 30 ml | CH 22 | purple | 10 pcs |
| CFNL-3D-24-30-P | 400 mm | 30 ml | CH 24 | dark blue | 10 pcs |
| * CFNL-3D-26-30-P | 400 mm | 30 ml | CH 26 | pink | 10 pcs |



Foley catheter

2-way with tiemann tip,

silicone-coated latex

features


- made of soft and flexible latex, that is resistant to bending and kinking universal urine drainage port of the catheter fits into standard connectors
- in urine bags
- angled tip makes it easier to overcome obstacles in the urinary tract
- tip with two lateral eyes with smooth edges
 silicone-coated surface that allows atraumatic insertion of the catheter
- easy to fill and to empty, tear-resistant balloon
- plastic valve to ensure balloon tightness
- single-use
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil indication

• to drain urine from the bladder in the case of a narrowed urethra

| REF | length | balloon | size | colour | intermedi packaging |
|-----------------|--------|---------|-------|-----------|------------------------|
| CFTL-2D-12-10-P | 400 mm | 5-10 ml | CH 12 | white | 10 pcs |
| CFTL-2D-14-10-P | 400 mm | 5-10 ml | CH 14 | green | 10 pcs |
| CFTL-2D-16-10-P | 400 mm | 5-10 ml | CH 16 | orange | 10 pcs |
| CFTL-2D-18-10-P | 400 mm | 5-10 ml | CH 18 | red | 10 pcs |
| CFTL-2D-20-10-P | 400 mm | 5-10 ml | CH 20 | yellow | 10 pcs |
| CFTL-2D-22-10-P | 400 mm | 5-10 ml | CH 22 | purple | 10 pcs |
| CFTL-2D-24-10-P | 400 mm | 5-10 ml | CH 24 | dark blue | 10 pcs |
| CFTL-2D-12-30-P | 400 mm | 30 ml | CH 12 | white | 10 pcs |
| CFTL-2D-14-30-P | 400 mm | 30 ml | CH 14 | green | 10 pcs |
| CFTL-2D-16-30-P | 400 mm | 30 ml | CH 16 | orange | 10 pcs |
| CFTL-2D-18-30-P | 400 mm | 30 ml | CH 18 | red | 10 pcs |
| CFTL-2D-20-30-P | 400 mm | 30 ml | CH 20 | yellow | 10 pcs |







Foley catheter

2-way, 100% silicone

features

- made of soft and flexible silicone, that reduces the allergic
- reaction of the mucous membranes and increases the patient's comfort • atraumatic, soft rounded, special tip facilitating insertion of the
- catheter to the bladder
- tip with two lateral eyes with smooth edges
 universal urine drainage port of the catheter fits into standard connectors in urine bags
- radiopaque line along the entire length of the catheter
- easy to fill and to empty, tear-resistant balloon of various capacities plastic valve to ensure balloon tightness
- CH 6, CH 8, CH 10 size has a nylon guide to facilitate catheterization • available in version with 10 ml sterile water pre-filled syringe with 10% glycerin to seal the balloon
- separately packed Catheter Spigot inside the packaging
- single-use
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil
- intermediate packaging: 10 pcs for REF: CFNS-2D-XX-XX-P, 5 pcs for REF: CFNSG-2D-XX-XX-P
 bulk packaging: 40 x 10 pcs for REF: CFNS-2D-XX-XX-P, 50 x 5 pcs for REF: CFNSG-2D-XX-XX-P

indication

- to drain urine from the bladder
- recommended for patients with blood clots or sediment in their urine

| REF | intermediate packaging | REF for pre-filled syringe | intermediate packaging | length | balloon | size | colour |
|------------------|---------------------------|-------------------------------|---------------------------|--------|---------|-------|-----------|
| CFNS-2D-06-03-P* | 10 pcs | CFNSG-2D-06-03-P | 5 pcs | 310 mm | 3 ml | CH 6 | pink |
| CFNS-2D-08-05-P* | 10 pcs | CFNSG-2D-08-05-P | 5 pcs | 310 mm | 3-5 ml | CH 8 | black |
| CFNS-2D-10-05-P* | 10 pcs | CFNSG-2D-10-05-P | 5 pcs | 310 mm | 3-5 ml | CH 10 | grey |
| CFNS-2D-12-10-P* | 10 pcs | CFNSG-2D-12-10-P | 5 pcs | 400 mm | 5-10 ml | CH 12 | white |
| CFNS-2D-14-10-P* | 10 pcs | CFNSG-2D-14-10-P | 5 pcs | 400 mm | 5-10 ml | CH 14 | green |
| CFNS-2D-16-10-P* | 10 pcs | CFNSG-2D-16-10-P | 5 pcs | 400 mm | 5-10 ml | CH 16 | orange |
| CFNS-2D-18-10-P* | 10 pcs | CFNSG-2D-18-10-P | 5 pcs | 400 mm | 5-10 ml | CH 18 | red |
| CFNS-2D-20-10-P* | 10 pcs | CFNSG-2D-20-10-P | 5 pcs | 400 mm | 5-10 ml | CH 20 | yellow |
| CFNS-2D-22-10-P* | 10 pcs | CFNSG-2D-22-10-P | 5 pcs | 400 mm | 5-10 ml | CH 22 | purple |
| CFNS-2D-24-10-P* | 10 pcs | CFNSG-2D-24-10-P | 5 pcs | 400 mm | 5-10 ml | CH 24 | dark blue |



anda Maria

14%-Up

Foley catheter

3-way, 100% silicone

features

- made of soft and flexible silicone, that reduces the allergic
- reaction of the mucous membranes and increases the patient's comfort • universal urine drainage port of the catheter fits into standard
- connectors in urine bags · additional inlet for flushing out the bladder
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- radiopaque line along the entire length of the catheter
- easy to fill and to empty, tear-resistant balloon
- plastic valve to ensure balloon tightness
- separately packed Catheter Spigot inside the packaging
- single-use • packaging: 1 pc./double packaging: internal- foil, external- paper-foil

indication

- to drain urine from the bladder
- recommended for patients with blood clots or sediment in their urine

| REF | length | balloon | size | colour | intermediate packaging |
|------------------|-----------|---------|-------|------------|---------------------------|
| | | • | | | |
| CFNS-3D-16-30-P | 400 mm | 30 ml | CH 16 | orange | 10 pcs |
| CENS-3D-18-30-P | 400 mm | 30 ml | CH 18 | red | 10 pcs |
| CFIN3-3D-10-30-F | 400 11111 | 50 111 | CH IO | reu | io pes |
| CFNS-3D-20-30-P | 400 mm | 30 ml | CH 20 | vellow | 10 pcs |
| | | | | J | - 1 |
| CFNS-3D-22-30-P | 400 mm | 30 ml | CH 22 | purple | 10 pcs |
| CFNS-3D-24-30-P | 400 mm | 30 ml | CH 24 | dark blue | 10 pcs |
| CFN3-3D-24-30-F | 400 11111 | 30 111 | CI124 | uai k biue | TO pes |
| | | | | | |



TATEX STERILE EO

18Fe/Ch

Nelaton catheter

features

- made of soft and flexible PVC, that is resistant to bending and kinking
- frozen surface for easier insertion
- semi-transparent tube (catheters available in a transparent version for better visual inspection)
- atraumatic, soft rounded, closed tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- · colour-coded connectors for size recognition, numerical size identification on the packaging
- single-use
- packaging: 1 pc./paper-foil

indication

• for urinary bladder catheterization

| REF for frozen tube | intermediate packaging | REF for NEW! transparent tube | intermediate packaging | length |
|------------------------|---------------------------|---|---------------------------|--------|
| CN-06-40** | 50 pcs | CNP-06-40* | 50 pcs | 400 mm |
| | | | | |
| CN-08-40** | 50 pcs | CNP-08-40* | 50 pcs | 400 mm |
| CN-10-40** | 50 pcs | CNP-10-40* | 50 pcs | 400 mm |
| CN-12-40** | 50 pcs | CNP-12-40* | 50 pcs | 400 mm |
| CN-14-40** | 50 pcs | CNP-14-40* | 50 pcs | 400 mm |
| CN-16-40** | 50 pcs | CNP-16-40* | 50 pcs | 400 mm |
| CN-18-40** | 50 pcs | CNP-18-40* | 50 pcs | 400 mm |
| CN-20-40** | 50 pcs | CNP-20-40* | 50 pcs | 400 mm |
| CN-22-40** | 50 pcs | CNP-22-40* | 50 pcs | 400 mm |
| CN-24-40** | 50 pcs | CNP-24-40* | 50 pcs | 400 mm |



Tiemann catheter

features

- made of soft and flexible PVC, that is resistant to bending and kinking
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- cone-shaped angled tip, bent at an angle of 45°
- transparent tube allowing visual inspection of excreted urine
- colour-coded connectors for size recognition, numerical size identification on the packaging
- single-use
- packaging: 1 pc./paper-foil

indication

• to drain urine from the bladder in the case of a narrowed urethra

| REF | length | size | colour | intermediate packaging | bulk packaging |
|----------|--------|-------|-----------|---------------------------|-------------------|
| CT-08-40 | 400 mm | CH 8 | blue | 60 pcs | 10 x 60 pcs |
| CT-10-40 | 400 mm | CH 10 | black | 60 pcs | 10 x 60 pcs |
| CT-12-40 | 400 mm | CH 12 | white | 50 pcs | 10 x 50 pcs |
| CT-14-40 | 400 mm | CH 14 | green | 50 pcs | 10 x 50 pcs |
| CT-16-40 | 400 mm | CH 16 | orange | 50 pcs | 10 x 50 pcs |
| CT-18-40 | 400 mm | CH 18 | red | 40 pcs | 10 x 40 pcs |
| CT-20-40 | 400 mm | CH 20 | yellow | 40 pcs | 10 x 40 pcs |
| CT-22-40 | 400 mm | CH 22 | purple | 40 pcs | 10 x 40 pcs |
| CT-24-40 | 400 mm | CH 24 | dark blue | 40 pcs | 10 x 40 pcs |
| | | | | | |







Ø bulk packaging: 10 x 50 pcs

| SIZE | COIOUI |
|-------|-------------|
| CH 6 | light green |
| CH 8 | blue |
| CH 10 | black |
| CH 12 | white |
| CH 14 | green |
| CH 16 | orange |
| CH 18 | red |
| CH 20 | yellow |
| CH 22 | purple |
| CH 24 | dark blue |
| | |
| | |

| size | colour |
|----------|-------------|
| CH 6 | light green |
| CH 8 | blue |
| CH 10 | black |
| CH 12 | white |
| CH 14 | green |
| CH 16 | orange |
| CH 18 | red |
| CL 1 2 0 | |

Pezzer catheter

silicone-coated latex

features

- made of soft and flexible silicone-coated latex, resistant to bending and kinking
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with three lateral eyes with smooth edges single-use
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil

indication

- to secure proper urinary drainage directly from the bladder or kidney in patients after surgical procedures
- introduced by the surgical wound, it is used for permanent urinary drainage, facilitating wound healing

| REF | length | size | intermediate packaging |
|----------|--------|-------|---------------------------|
| CP-12-40 | 400 mm | CH 12 | 10 pcs |
| CP-14-40 | 400 mm | CH 14 | 10 pcs |
| CP-16-40 | 400 mm | CH 16 | 10 pcs |
| CP-18-40 | 400 mm | CH 18 | 10 pcs |
| CP-20-40 | 400 mm | CH 20 | 10 pcs |
| CP-22-40 | 400 mm | CH 22 | 10 pcs |
| CP-24-40 | 400 mm | CH 24 | 10 pcs |
| CP-26-40 | 400 mm | CH 26 | 10 pcs |
| CP-28-40 | 400 mm | CH 28 | 10 pcs |
| CP-30-40 | 400 mm | CH 30 | 10 pcs |
| CP-32-40 | 400 mm | CH 32 | 10 pcs |
| CP-34-40 | 400 mm | CH 34 | 10 pcs |
| CP-36-40 | 400 mm | CH 36 | 10 pcs |
| CP-38-40 | 400 mm | CH 38 | 10 pcs |
| CP-40-40 | 400 mm | CH 40 | 10 pcs |
| | | | |



Ø

bulk packaging:

25 x 100 pcs

samples taking urine bag

features

- latex-free bag with a capacity of 100 ml, size 19 cm x 7 cm
- made of delicate, durable PE foil
- easy-to-read scale, adequate for diuresis measurement, highly accurate every 10 ml
 adhesive part of the bag coated with anti-allergic, high-quality glue
 additional sealing made of anti-allergic foam that protects against the wall sticking of the bag, thanks to which free flow of urine is possible
- single-use
- packaging: 1 pc./foil

indication

• for urine sample collection in infants and small children

| REF | capacity | purpose | version | intermediate packaging |
|----------|----------|-----------|----------------|---------------------------|
| WMCH-001 | 100 ml | for boys | with sponge | 100 pcs |
| WMCH-002 | 100 ml | for boys | without sponge | 100 pcs |
| WMDZ-001 | 100 ml | for girls | with sponge | 100 pcs |
| WMDZ-002 | 100 ml | for girls | without sponge | 100 pcs |
| | | | | |





urine bag

features

- latex-free 2000 ml bag made of 0,22 mm thick PVC (double wall)
- double sealed PVC for maximum protection against tearing the bag
- effective non-return valve
- tight T-valve, that can be easily operated with one hand
- · connecting tubing with universal stepped connector, standard tubing length 90 cm and 150 cm
- tubing made of material that is resistant to bending and kinking and which ensures a free and effective urine outflow
- precise, easy-to-read scale (every 100 ml), adequate for diuresis measurement • strengthened holes for hanging, matching standard hangers
- white back wall making urine easy to notice
- single-use
- packaging: 1 pc./foil

indication

- for 24-hour extracorporeal collection of urine and other body fluids
- maximum time of use 24 hours

| REF | capacity | length of tubing | intermediate packaging | |
|-------------|----------|---------------------|---------------------------|-----------------|
| WMD2000-90 | 2000 ml | 90 cm | 10 pcs | bulk packaging: |
| WMD2000-150 | 2000 ml | 150 cm | 10 pcs | 25 x 10 pcs |



urine bag

with sample port

features

- latex-free 2000 ml bag made of 0,40 mm thick PVC (double wall)
- effective non-return valve
- tight T-valve, that can be easily operated with one hand
- connecting tubing with universal stepped connector, standard tubing length 120 cm (other lengths available: 90 cm, 150 cm)
- tubing made of material that is resistant to bending and kinking and which ensures a free and effective urine outflow
- self-tightening port for sample taking
- precise, easy-to-read scale, adequate for diuresis measurement, with high measurement accuracy from 25 ml, every 50 ml from 50 ml to 200 ml and every 100 ml to 2000 ml
- strengthened holes for hanging, matching standard hangers
- · white back wall making urine easy to notice
- single-use • packaging: 1 pc./paper-foil

indication

• for extracorporeal collection of urine and other body fluids • maximum time of use - 7 days

| REF | capacity | length of tubing | intermediate packaging |
|-------------|----------|---------------------|---------------------------|
| WMT2000-90 | 2000 ml | 90 cm | 10 pcs |
| WMT2000-120 | 2000 ml | 120 cm | 10 pcs |
| WMT2000-150 | 2000 ml | 150 cm | 10 pcs |







24 Fr/Ch





urine bag hanger

features

- made of strong and durable plastic
- special clamp that prevents tube bending
- suitable for round and square bed rails
- two handles on each side for hanging urine bags of various sizes · manufacturer's name placed on the device
- single-use

indication

• for hanging standard urine collection bags







ZARYS

urine bag


with sample port 2W

features

- latex-free 2000 ml bag made of medical grade PVC
- · white back wall making urine easy to notice
- effective non-return valve • tight T-valve, that can be easily operated with one hand and which can be folded and tucked in a hard sheath
- Pasteur type droplet chamber providing a "dry" air barrier
- two hydrophobic air filters (in a bag and chamber), which prevent sucking,
- balancing the pressure inside the system
- 120 cm long, large diameter tubing, ensuring a free and effective urinary outflow with a universal stepped connector
- self-tightening port for sample taking at the end of tubing • tubing equipped with a sliding clamp and a fastening pin to stabilize
- the tubing • at the place where the tubing connects to the chamber, an additional
- reinforcing element prevents bending and kinking the tubing • precise, easy-to-read scale, adequate for diuresis measurement, with high measurement accuracy, every 25 ml from 25 ml to 100 ml and every 100 ml
- from 100 ml to 2000 ml
- tag to describe the patient's details and date
- integrated universal hanger system matching the round and square bed rails • additional fixing with a string
- single-use
- packaging: 1 pc./paper-foil

indication

• device is intended for 14-days urine collection in a closed system







GENERAL SURGERY

Yankauer surgical suction cannula

standard CH21

features

- completely transparent for full possibility of observing the sucked secretion • light and perfectly balanced handle ensures high comfort of work and reduces hand fatigue
- perfectly smooth inner surface prevents bLocking or retaining the sucked content, ensuring correct and optimal flow
- · cannula with double bend of curvature • four lateral eyes, which prevent the tip from sticking to the tissues and one
- central opening • universal stepped connector ensures tight connection with a wide range of suction connecting tubes
- sinale-use
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil

indication

• evacuating secretions and body fluids from the operating field in general surgery





Yankauer surgical suction cannula

bulb tip CH23

features

- cannula with bulb tip with four lateral eyes, which prevent the tip from sticking to the tissues and one central opening
- completely transparent for full possibility of observing the sucked secretion
- · light and perfectly balanced handle ensures high comfort of work and reduces hand fatigue
- · perfectly smooth inner surface prevents bLocking or retaining the sucked content, ensuring correct and optimal flow
- · cannula with double bend of curvature
- universal stepped connector ensures tight connection with a wide range of suction connecting tubes
- single-use • packaging: 1 pc./double packaging: internal- foil, external- paper-foil

indication

• evacuating secretions and body fluids from the operating field in general surgery





surgical suction tube

features

- made of soft and flexible PVC
- along the tube a special longitudinal reinforcement to prevent bending and sucking of the tube
- perfectly smooth inner surface of tube prevents the sucked content from settling on the walls, ensuring correct and optimal flow
- connector type funnel-funnel cut-to-fit with the possibility of cutting in the marked places of the connector and matching to any type of medical suction device (for the diameter of 8 mm to 18 mm)
- funnel connector equipped with a special bending system, enabling convenient manipulation of the tube
- tube connectors compatible with the standard suction cannulas, enabling their reliable, effective and tight connection
- single-use
- packaging: 1 pc./double packaging: internal- foil, external- paper-foil indication

- · device designed for connecting the suction cannula to medical suction device • Kapkon connector - designed to connect the catheter for aspiration of upper respiratory tract to medical suction device

| REF | size | lenght of tube | connector |
|----------|------|-------------------|--------------------------|
| DO-24201 | CH24 | 210 cm | funnel-funnel |
| DO-24301 | CH24 | 300 cm | funnel-funnel |
| DO-30201 | CH30 | 210 cm | funnel-funnel |
| DO-30301 | CH30 | 300 cm | funnel-funnel |
| DO-24204 | CH24 | 210 cm | funnel-funnel cut-to-fit |
| DO-24202 | CH24 | 210 cm | funnel-Kapkon |

STERILEEO

surgical suction set

with Yankauer cannula, standard

CANNULA

- completely transparent for full possibility of observing the sucked secretion
- · light and perfectly balanced handle ensures high comfort of work and reduces hand fatigue
- · perfectly smooth inner surface prevents bLocking or retaining the sucked content, ensuring correct and optimal flow
- cannula with double bend of curvature
- four lateral eyes, which prevent the tip from sticking to the tissues and one central opening • universal stepped connector ensures tight connection with a wide range
- of suction connecting tubes

TUBE

- made of soft and flexible PVC
- along the tube a special longitudinal reinforcement to prevent bending and sucking of the tube
- perfectly smooth inner surface of tube prevents the sucked content
- from settling on the walls, ensuring correct and optimal flow connector type funnel-funnel cut-to-fit with the possibility of cutting in the marked places of the connector and matching to any type of medical suction device (for the diameter of 8 mm to 18 mm)
- funnel connector equipped with a special bending system, enabling convenient manipulation of the tube
- tube connectors compatible with the standard suction cannulas, enabling their reliable, effective and tight connection

features

- available with or without suction control
- single-use
- packaging: 1 set/double packaging: internal- foil, external- paper-foil

indication

• evacuating secretions and body fluids from the operating field in general surgery

| REF | size | lenght of tube | connector |
|---------------|-----------|-------------------|---------------------------------------|
| ZO-212400104 | CH21/CH24 | 210 cm | without suction control, funnel-funn |
| ZO-212400114 | CH21/CH24 | 210 cm | with suction control, funnel-funnel |
| ZO-402423104 | CH21/CH24 | 400 cm | without suction control, funnel-funn |
| ZO-402423114 | CH21/CH24 | 400 cm | with suction control, funnel-funnel |
| ZO-2124001041 | CH21/CH24 | 210 cm | without suction control, funnel-funn |
| ZO-2124001141 | CH21/CH24 | 210 cm | with suction control, funnel-funnel c |







nel cut-to-fit bulk packaging: cut-to-fit 50 x 1 set

surgical suction set

with Yankauer cannula, bulb tip

CANNULA

- · cannula with bulb tip with four lateral eyes, which prevent the tip from sticking to the tissues and one central opening
- completely transparent for full possibility of observing the sucked secretion • light and perfectly balanced handle ensures high comfort of work and reduces hand fatigue
- perfectly smooth inner surface prevents bLocking or retaining the sucked content, ensuring correct and optimal flow
- · cannula with double bend of curvature
- universal stepped connector ensures tight connection with a wide range of suction connecting tubes

TUBE

- made of soft and flexible PVC
- along the tube a special longitudinal reinforcement to prevent bending and sucking of the tube
- perfectly smooth inner surface of tube prevents the sucked content from settling on the walls, ensuring correct and optimal flow
- connector type funnel-funnel cut-to-fit with the possibility of cutting in the marked places of the connector and matching to any type of medical suction device (for the diameter of 8 mm to 18 mm)
- funnel connector equipped with a special bending system, enabling convenient manipulation of the tube
- tube connectors compatible with the standard suction cannulas, enabling their reliable, effective and tight connection

features

- available with or without suction control
- single-use • packaging: 1 set/double packaging: internal- foil, external- paper-foil


indication

· evacuating secretions and body fluids from the operating field in general surgery

| | REF | size | lenght of tube | connector |
|-----|---------------|-----------|-------------------|---|
| - 1 | | | | |
| | ZO-2124002041 | CH23/CH24 | 210 cm | without suction control, funnel-funnel cut-to-fit |
| | ZO-2124002141 | CH23/CH24 | 210 cm | with suction control, funnel-funnel cut-to-fit |
| | ZO-3024002041 | CH23/CH24 | 300 cm | without suction control, funnel-funnel cut-to-fit |
| | 70 2024002141 | CH32/CH34 | 200 cm | with suction control funnal funnal cut to fit |



STERILEEO



BASIC PRODUCTS FOR GYNECOLOGY





vaginal speculum

cusco type

features

- both blades made of transparent polystyrene
- both blades made of transparent polystyrene
 locking bolt made of polyethylene
 blades position blocked by a quarter turn of the Locking bolt (without screwing in the bolt)
 perfectly smooth and streamlined edges
 colour coded bolts depending on the size
 cipalouse
- single-use
- packaging: 1 pc./foil

indication

• device designed for gynecological examination of the cervix and vagina

| REF | width | colour | size | |
|--------|-------|--------|------|-----------------|
| WG-XXS | 16 mm | yellow | XXS | |
| WG-XS | 20 mm | white | XS | |
| WG-S | 24 mm | blue | S | |
| WG-M | 26 mm | red | Μ | |
| WG-L | 30 mm | green | L | bulk packaging: |
| WG-XL | 32 mm | black | XL | 100 pcs |





cervical brush

special

features

- broom type enables obtaining samples from the cervix, cervical canal and transformation zone
- soft bristle tips enables precise collection of material for examination
- made entirely of polyethylene
 optimal outreach thanks to the long handle

- single-use
 packaging: 1 pc./paper-foil

indication

device designed for obtaining cell samples from the cervix, cervical canal and transformation zone in smear tests

| | packaging | |
|----------------|-----------|---------------------------------|
| SC-002 special | 100 pcs | bulk packaging: 10 x 100 pcs |



umbilical cord clamp

features

- made of polypropylene
 safe clamping teeth with rounded edges do not cut
 the umbilical cord
- the umbilical cord
 serrated clamping surface of the arms stabilize the clamp on the umbilical cord, after closing it is not possible to open the clamp again
 double Lock catch prevents accidental opening
 universal size
 single-use

- non-pyrogenicpackaging: 1 pc./paper-foil

indication

 device is designed for clamping the umbilical cord of a newborn in a labour ward in order to close the lumen of the umbilical blood vessels after cutting it off of the placenta

| REF | size | colour | intermediate packaging | bulk assisted |
|------|-----------|--------|---------------------------|--------------------------------|
| ZDPB | universal | white | 50 pcs | bulk packaging: 50 x 50 pcs |



cervical brush standard

features

- straight shape enables obtaining samples from the cervical canal
 soft bristle tips to prevent irritation
 optimal outreach thanks to the long handle

- single-use
- packaging: 1 pc./paper-foil


indication

 device designed for obtaining cell samples from the cervix, cervical canal and transformation zone in smear tests









ZARYS

ID wristband

features

- durable wristband coating characterized by high humidity resistance and resistance to the activity of soaps, hand disinfectants and other liquids used in hospitals
- made of soft and delicate vinyl
- flexible, delicate, very light and easy to put on
- rounded edges prevent from injuring
- high tensile strength
- · colour marking depending on the purpose
 - · single-use fastening resistant to accidental break
 - hypoallergenic, latex-free, phthalate-free • available in a version with a description field directly on the band
 - and with an insert card
 - packaging: 100 pcs/box

indication

• for identification of patients in hospital wards

| REF | purpose | colour | version | intermediate packaging |
|----------|-------------------------|--------|-------------|---------------------------|
| OPD-B-01 | for children and adults | white | insert card | 100 pcs |
| OPN-R-01 | for infants | pink | insert card | 100 pcs |
| OPN-N-01 | for infants | blue | insert card | 100 pcs |
| 0 | | r. | | |

entrance mat

features

- blue mat consisting of 30 layers thin sticky sheets
- each of the sticky sheets made of polyethylene resin
- · does not contain heavy metals
- excellent adhesive properties
- the first and last sheets secured with a non-adhesive foil
- low mat profile: app. 2 mm thick
- the bottom surface adheres very well to the floor over the entire surface, preventing adverse movement of the mat, allowing trouble-free passage of trolleys with equipment or beds, without the risk of curling up the mat or tearing • no need to clean the mat - as soon as the surface of the mat is contaminated,
- it should be torn off to uncover a new, clean sheet
- there is numbering of sheets on a white non-adhesive material in the mat corner facilitate quick identification and easy removal
- of subsequent sheets
- light, easy to spread

indication

preventing spreading of contamination brought on shoe soles and wheelchairs to operating rooms





features

- bag with a capacity of 1750 ml with two holes for hanging made of medical PVC
- scaled from 50 ml to 1750 ml
- linear scale every 100 ml, digital scale every 200 ml
- tube with a length of 150 cm with a sliding clamp, ended
- with an atraumatic opening and one lateral opening • tube end is secured with a tube cap
- upper part of the bag is reinforced with an additional flange
- stabilizing the inlet during filling with liquid • flange is equipped with a non-return valve
- sinale-use
- packaging: 1 pc./foil

indication


• bag that allows introduction of fluid into the large bowel through the anus to remove fecal material and stimulate peristaltic bowel movement

 \bigotimes





| | 0 | ΓH | Ε | R |
|----|----|----|---|---|
| PR | Dl | JC | | S |









enema bag

features

- bag with a capacity of 1500 ml with a hole for hanging
- made of medical PVC
- digitally scaled every 250 ml
 tube with a length of 115 cm with a sliding clamp, ended with an atraumatic opening and one lateral opening
- greased end of the tube secured with a cap

single-use

indication

• bag that allows introduction of fluid into the large bowel through the anus to remove fecal material and stimulate peristaltic bowel movement



enema set

features

- bag with a capacity of 1500 ml with a hole for hanging
- made of medical PVC • digitally scaled every 250 ml
- tube with a length of 115 cm with a sliding clamp, ended with
- an atraumatic opening and one lateral opening
- · greased end of the tube secured with a cap
- drape 42 cm x 45 cm
- liquid soap (9 ml)
- PE gloves 1 pair
- single-use
- · packaging: 1 set/foil

indication

• set that allows introduction of fluid into the large bowel through the anus to remove fecal material and stimulate peristaltic bowel movement

Ø REF NON STERRIE CATERY PHT OF STERRIE bulk packaging: ZLN 50 x 1 set



vomit bag

features

- bag is made of durable, transparent blue foil, enabling the observation of secretions • total bag capacity: 2000 ml
- exact measuring graduation (from 10 ml to 100 ml every 10 ml and from 100 to 2000 ml every 100 ml) placed on the bag allows accurate estimation of fluid volume
- direction of graduation reading for small volumes is marked with the tip
- wide bag inlet protected by a plastic collar
- circle-shaped collar ensures a firm grip, makes it easier to maneuver the bag, reducing the risk of content contamination
- collar equipped with a special incision to close the bag and hygienic disposal of the vomit - "twist and hook"
- · device name, manufacturer's name and pictorial instructions for use printed directly on the bag
- single-use
- · class I non-sterile medical device
- packaging: 20 pcs/foil bag

indication

- recommended in uncomfortable situations, in cases of gastroesophageal reflux requiring immediate action
- used in hospitals, outpatient departments, ambulances, etc.

| REF | colour | intermediate packaging | bulk assissing | (LALER) (NHT) |
|-------|--------|---------------------------|--------------------------------|----------------|
| WW-01 | blue | 20 pcs | bulk packaging: 30 x 20 pcs | NON STERILE |



safeLANCE

pressure-activated safety lancet

features

- easy to use, minimizes unpleasant feelings
- special construction allows lancet tip to be automatically inserted into finger pad and then to be retracted so that needle is not visible
- before or after puncturing · correctly matched size of lancet allows you to get optimal amount
- of blood with minimal patient's pain design ensures safety for medical personnel and patient by eliminating
- accidental puncture and infection · perfect guidance of lancet needle inside housing reduces vibration,
- thus decreases patient's pain and prevents damage to skin tissue • available in two sizes:
- 23 G, puncture depth 1,8 mm green, dedicated for women and children - 21 G, puncture depth 2,4 mm - pink, dedicated for men
- single-use • packaging: 100 pcs/box

indication

· capillary blood sampling from fingertip or ear flap for testing, for example, glucose level

| REF | size / puncture depth | colour | intermediate packaging | |
|-------|--------------------------|--------|---------------------------|---------------|
| NA-18 | 23 G / 1,8 mm | green | 100 pcs | bulk packagir |
| NA-24 | 21 G / 2,4 mm | pink | 100 pcs | 20 x 100 pc |

dental brush

with suction

features

- · made of polypropylene
- soft bristles on one side, foam on the other side
- total length 18 cm
- length of the cleaning part 2,5 cm
- suction opening from both the bristle side and the space between the foam and the bristles
- controlled suction connector at 45° angle for convenient manipulation single-use
- packaging: 1 pc./foil

indication

- hygiene of the patient's oral cavity
- function for simultaneous brushing, rinsing and suction

| REF | intermediate packaging | | \wedge | (ATEX) | (Sul) |
|--------|---------------------------|--------------------------------|----------|--------|-------|
| HJM-SO | 50 pcs | bulk packaging: 10 x 50 pcs | NON | XX | X |

oral swab

features

- overall length 15 cm (+/- 2 mm)
- length of the foam part 2.5 cm handle made of polyester
- swab made of polypropylene
- single-use
- packaging: 1 pc./foil

indication

- intended for oral hygiene in patients partially or completely
- dependent on medical personnel
- allows the oral mucosa to be kept in proper humidity allows you to perform a procedure to prevent
- pathological changes in the mouth
- \bigotimes REF intermediate packaging bulk packaging: HJM-GB 50 pcs 10 x 50 pcs



















tactical tourniquet CAT

features

- made of soft and breathable material to prevent chafing and skin burns
- · additional element to secure the tourniquet and control blood flow
- locking system to prevent accidental opening of the tourniquet
- one-handed mechanism
- durable loop system with internal, free-moving strap that provides circumferential compression to the limb
- the unique double fastening system allows to avoid the use of screws and clips
- high quality of workmanship
- place to write down the time of establishment
- weight: 82g
- folded dimensions: 16 cm x 6 cm x 4cm
- unfolded dimensions: 96 cm x 4 cm
- · individually wrapped in foil packaging single-use

indication

· tactical tourniquet used in case of limb injuries with massive arterial or venous bleeding that cannot be controlled by traditional techniques

| REF | dimensions | intermediate packaging | |
|--------|--------------|---------------------------|--------------------------------|
| ST-001 | 96 cm x 4 cm | 1 рс. | bulk packaging: 100 x 1 pcs |



automatic tourniquet

features

- components: blue stretchy strap and white ABS plastic buckle
- length 480 mm +/- 10 mm, width 25 mm +/- 10 mm • wide, stretchable band put above the venipuncture site reduces
- the pressure intensity • protects against the congestion of the arm, minimizes unpleasant sensations, improves the quality of the blood sample
- reusable
- packaging: 1 pc./double packaging: internal- foil, external- box

indication

• used to induce venous stasis and to make a vein selected for puncture more prominent







for children

features

- components: polyester stretch strap and ABS plastic items
- length 450 mm +/- 5 mm, width 25 mm +/- 1 mm
- wide stretch strap to minimise pressure build-up during compression
- significantly protects against arm congestion, improves the quality of the blood sample for
- examination • reduces unpleasant sensations during compression
- a simple one-handed mechanism for easy fastening and unfastening with gradual release of pressure
- reusable

indication

• device used to induce venous stasis and to make a vein selected for puncture


more prominent

RFF intermediate intermediate packaging packaging



Ø bulk packaging: yellow + imprint 20 x 21 pcs







latex-free tourniquet

features

- made of a wide stretch TPE strap (thermoplastic elastomer)
- · latex-free material prevents allergic reactions and skin irritation
- high tensile strength
- single-use bands reduce the risk of infection • total roller length 11,25 m
- single-use
- individual packaging contains: 25 tourniquets
- · packaging: 1 roll/paper dispenser that allows convenient sharing of perforated tourniquets

indication

• used to induce venous stasis and to make a vein selected for puncture more prominent



tongue depressor

wooden

features

- made of birch wood
- rounded edges preventing injuries
- natural colour
- · long outreach thanks to appropriate length
- durable
- single-use

packaging: 1 pc./paper

indication

diagnostic examination of the oral cavity

| SLS-001 sterile 150 mm x 17 mm 20 x 100 pcs | REF | version | dimensions | bulk packaging |
|---|---------|---------|----------------|-------------------|
| | SLS-001 | sterile | 150 mm x 17 mm | 20 x 100 pcs |



tongue depressor


wooden

features

- made of birch wood
- rounded edges preventing injuries • natural colour
- · long outreach thanks to appropriate length
- durable
- single-use
 - packaging: 100 pcs/box

indication


· diagnostic examination of the oral cavity











electronic thermometer

features

- digital thermometer
- fast, safe and reliable
- temperature readings
- waterproof
- measurement in 10 seconds
 range 32.0°C-42.9°C (90.0°F-109.9°F)
 accuracy ±0.1°C (±0.2°F)
- has a soft, flexible measuring tip
- sound signal
- automatic turning off
- memory of the last measurement
- low battery indicator
- battery life: approximately 200 hours
- package: 1 piece / plastic cover and paper box

indication

• for measuring body temperature in the mouth, armpit and rectal area

Ø

bulk packaging:

20 x 50 pcs





wash glove


molton

features

- wash glove in the form of a glove to put on a hand
- made of MOLTON material
- ultrasonic sewing
- soft, smooth, gentle and absorbent
- dry, without cleaning agents • tear resistant
- single-use



vitahealth





wash glove

Spunlace

features

- wash glove in the form of a glove to put on a hand
 made of SPUNLACE material
- ultrasonic sewing • very soft, gentle and absorbent
- dry, without cleaning agents
- tear resistant
- single-use







wash glove

with film lining

features

- wash glove with film lining in the form of a glove to put on a hand
 made of Airlaid material + PE
- ultrasonic sewing
- gentle and absorbent • dry, without cleaning agents
- tear resistant
- sinale-use

| 5 | | | |
|---|--|--|--|
| | | | |
| | | | |
| | | | |
| | | | |

| REF | grammage | dimensions | unit packaging | |
|--------|----------|-----------------|----------------|--------------------------------|
| MB13PE | 80 g/m² | 160 mm x 220 mm | 50 pcs | bulk packaging: 20 x 50 pcs |
| | | _ | | |



PE gloves

features

- made of polyethylene
- coarse surface and tear-resistant weld
- ambidextrous
- elongation min. 130%, tensile strength min. 10 MPa
 approved for contact with food
- single-use • packaging: 100 pcs/foil bag enabling individual removal of gloves

indication

• intended for use as a contribution to other gloves (double-gloving) and for minimally invasive procedures such as drug preparation, feeding the patient

| REF | size | length | |
|---------|-----------|--------|-----------------|
| BS-RF-M | M (women) | 280 mm | bulk packaging: |
| BS-RF-L | L (men) | 300 mm | 100 x 100 pcs |



| PE ov | versleev | /es | NEW |
|--|------------------|--------|----------------|
| ended w colour: b | d for contact wi | | |
| REF | dimensions | weight | unit packaging |



BS-ZF01 22 x 40 cm 2,6 g 100 pcs

82 www.zarys.com









medical waste container

features

- blow-moulded container made of high density polyethylene (HDPE), lid made of polypropylene (PP)
- lid equipped with a rosette insertion hole and two holes to separate blades from handles, needles from syringes, etc.
- dual-function lid for temporary or permanent closure of the container
- high resistance to mechanical deformation
- manufactured in red warning colour, other colours (yellow, blue) available on request • supplied with a BML-type label, injected during blowing, with an international warning sign and place for easily writing the waste data, even with a pen • resistant to UV radiation, formalin and ethyl alcohol
- it is approved by the National Institute of Public Health
- according to the requirements of the standards: ASTM F2132-01 (puncture resistance) and EN ISO 175 (resistance to immersion in liquid chemicals) sinale-use
- MOQ 5000 pcs

indication

- intended for the collection of solid, hazardous medical waste such as needles, intravenous cannulas, syringes, disposable gloves, etc.
- the container enables the safe storage, transport and disposal of medical waste, ensuring safety in terms of hygiene and sanitary standards, as well as mechanical safety to prevent unwanted events such as punctures and cuts when medical staff dispose of sharp waste

\$ bulk packaging: 20 pcs

| REF | capacity | colour | bottom diameter | upper diameter | height |
|-------------|----------|--------|--------------------|-------------------|--------|
| *POM-00700R | 0,7L | red | 93 mm | 109 mm | 97 mm |
| *POM-00700Y | 0,7L | yellow | 93 mm | 109 mm | 97 mm |
| *POM-00700B | 0,7L | blue | 93 mm | 109 mm | 97 mm |
| *POM-02000R | 2L | red | 95 mm | 109 mm | 246 mm |
| *POM-02000Y | 2L | yellow | 95 mm | 109 mm | 246 mm |
| *POM-02000B | 2L | blue | 95 mm | 109 mm | 246 mm |





medical waste container

features

- ontainer entirely made of polypropylene (PP)
- · lid equipped with a rosette insertion hole and two holes to separate blades from handles, needles from syringes, etc.
- dual-function lid for temporary or permanent closure of the container
- high resistance to mechanical deformation
- manufactured in red warning colour, other colours (yellow, blue) available on request
- supplied with label with an international warning sign and place for easily writing the waste data
- single-use
- MOQ 5000 pcs

indication

- intended for the collection of solid, hazardous medical waste such as needles, intravenous cannulas, syringes, disposable gloves, etc.
- the container enables the safe storage, transport and disposal of medical waste, ensuring safety in terms of hygiene and sanitary standards, as well as mechanical safety to prevent unwanted events such as punctures and cuts when medical staff dispose of sharp waste

| REF | capacity | colour | bottom diameter | upper diameter | height |
|-------------|----------|--------|--------------------|-------------------|--------|
| *POM-01000R | 1L | red | 92,5 mm | 109 mm | 140 mm |
| *POM-01000Y | 1L | yellow | 92,5 mm | 109 mm | 140 mm |
| *POM-01000B | 1L | blue | 92,5 mm | 109 mm | 140 mm |
| *POM-05000R | 5L | red | 190 mm | 230 mm | 210 mm |
| *POM-05000Y | 5L | yellow | 190 mm | 230 mm | 210 mm |
| *POM-05000B | 5L | blue | 190 mm | 230 mm | 210 mm |







features

- underpad with guilted absorbent pad containing superabsorbent
- white colour
- equipped with 4 adhesives for sticking the underpad to the substrate
- · absorbs fluid and holds it inside
- foil outer layer prevents leakage
- protects against contamination and soiling of the substrate
- soft
- disposable

| REF | dimensions | packaging unit | bulk packaging |
|------|---------------|-------------------|-------------------|
| PHZ1 | 60 cm x 60 cm | 30 pcs | 10 x 30 pcs |
| PHZ2 | 60 cm x 90 cm | 15 pcs | 10 x 15 pcs |

NEW









dressings

GAUZE AND NON-WOVEN DRESSINGS

DRESSING GAUZE

features

- 100% hydrophilic cotton
- bleached with the use of hydrogen peroxide
- soft and very absorbent
- increase in absorption capacity with increasing number of gauze threads
- homogeneous structure through a strong twist of fibers eliminates the risk of loose fibers on the surface

indication

• general wound dressing in out-patient clinics, surgical sites, procedure department

dressing gauze, zig-zag folded

GAZA lux

| in accordar | on open wounds, the de ace with a validated steri ith steam, ethylene oxid |
|-------------|--|
| REF | dimensions |

NEW GR1709100-XR 0.9 m x 100 m GR1309100 0,9 m x 100 m GR1709100 0,9 m x 100 m

of can be cut

cutting gauze

GAZA lux

and formed

| Before use on open wounds, the de in accordance with a validated steri sterilized with steam, ethylene oxid | | | | | |
|---|--------------------|--|--|--|--|
| REF | size | | | | |
| * CG1325SQM | 1/4 m ² | | | | |
| CG1305SQM | 1/2 m ² | | | | |
| CG131SQM | 1 m ² | | | | |
| * CG1725SQM | 1/4 m ² | | | | |
| *CG1705SQM | 1/2 m ² | | | | |
| CG171SQM | 1 m² | | | | |
| | | | | | |

NON DATES STEALS

GAZA lux S

| | REF | size | no. of threads | packaging unit paper-foil | bulk packaging carton |
|-----|--------------|--------------------|-------------------|------------------------------|--------------------------|
| NEW | CG1748255-SS | 25 cm x 25 cm | 17 | 5 pcs | 25 x 5 pcs |
| | *CG1325SQM-S | 1/4 m ² | 13 | 1 pc. | 100 x 1 pc. |
| | CG1305SQM-S | 1/2 m ² | 13 | 1 pc. | 50 x 1 pc. |
| | CG131SQM-S | 1 m ² | 13 | 1 pc. | 25 x 1 pc. |
| | *CG1725SQM-S | 1/4 m ² | 17 | 1 pc. | 100 x 1 pc. |
| | CG1705SQM-S | 1/2 m ² | 17 | 1 pc. | 50 x 1 pc. |
| | CG171SQM-S | 1 m ² | 17 | 1 pc. | 25 x 1 pc. |
| | | | | | |



TRADITIONAL DRESSING MATERIALS











available on request

vice must be sterilized ilization method. It can be le or radiation

| no. of threads | packaging unit foil bag | bulk packaging carton |
|-------------------|----------------------------|-----------------------|
| 17 | 1 pc. | 5 x 1 pc. |
| 13 | 1 pc. | 5 x 1 pc. |
| 17 | 1 pc. | 5 x 1 pc. |
| | | |



evice must be sterilized ilization method. It can be de or radiation

| no. of threads | packaging unit paper-foil | bulk packaging carton |
|-------------------|------------------------------|-----------------------|
| 13 | 1 pc. | 100 x 1 pc. |
| 13 | 1 pc. | 50 x 1 pc. |
| 13 | 1 pc. | 25 x 1 pc. |
| 17 | 1 pc. | 100 x 1 pc. |
| 17 | 1 pc. | 50 x 1 pc. |
| 17 | 1 pc. | 25 x 1 pc. |





GAUZE SWABS

features

- 100% hydrophilic cotton

- bleached with the use of hydrogen peroxide
 edges ES-type folded, which prevents fraying of loose threads
 available with a sewed in X-ray detectable thread, which enables fast and non-invasive identification
- of the swab left in the operating field

indication


• general wound dressing in out-patient clinics, surgical sites, procedure departments



Before use on open wounds, the device must be sterilized in accordance with a validated sterilization method. It can be sterilized with steam or ethylene oxide.









without X-ray detectable thread

| REF | dimensions | no. of threads | no. of plies | packaging unit paper bag | bulk packaging carton |
|-------------|-----------------|-------------------|-----------------|-----------------------------|-----------------------|
| GS1305081 | 5 cm x 5 cm | 13 | 8 | 100 pcs | 250 x 100 pcs |
| GS1375081 | 7,5 cm x 7,5 cm | 13 | 8 | 100 pcs | 60 x 100 pcs |
| GS1310081 | 10 cm x 10 cm | 13 | 8 | 100 pcs | 40 x 100 pcs |
| GS1305121 | 5 cm x 5 cm | 13 | 12 | 100 pcs | 80 x 100 pcs |
| GS1305125 | 5 cm x 5 cm | 13 | 12 | 500 pcs | 4 x 500 pcs |
| GS1375121 | 7,5 cm x 7,5 cm | 13 | 12 | 100 pcs | 40 x 100 pcs |
| GS1375125 | 7,5 cm x 7,5 cm | 13 | 12 | 500 pcs | 4 x 500 pcs |
| GS1310121 | 10 cm x 10 cm | 13 | 12 | 100 pcs | 20 x 100 pcs |
| GS1310125 | 10 cm x 10 cm | 13 | 12 | 500 pcs | 4 x 500 pcs |
| GS1705081 | 5 cm x 5 cm | 17 | 8 | 100 pcs | 80 x 100 pcs |
| GS1775081 | 7,5 cm x 7,5 cm | 17 | 8 | 100 pcs | 90 x 100 pcs |
| GS1710081 | 10 cm x 10 cm | 17 | 8 | 100 pcs | 40 x 100 pcs |
| GS171020081 | 10 cm x 20 cm | 17 | 8 | 100 pcs | 40 x 100 pcs |
| GS1705121 | 5 cm x 5 cm | 17 | 12 | 100 pcs | 96 x 100 pcs |
| GS1705125 | 5 cm x 5 cm | 17 | 12 | 500 pcs | 4 x 500 pcs |
| GS1775121 | 7,5 cm x 7,5 cm | 17 | 12 | 100 pcs | 60 x 100 pcs |
| GS1775125 | 7,5 cm x 7,5 cm | 17 | 12 | 500 pcs | 4 x 500 pcs |
| GS1710121 | 10 cm x 10 cm | 17 | 12 | 100 pcs | 20 x 100 pcs |
| GS1710125 | 10 cm x 10 cm | 17 | 12 | 500 pcs | 4 x 500 pcs |
| GS171020121 | 10 cm x 20 cm | 17 | 12 | 100 pcs | 20 x 100 pcs |
| GS1705161 | 5 cm x 5 cm | 17 | 16 | 100 pcs | 96 x 100 pcs |
| GS1775161 | 7,5 cm x 7,5 cm | 17 | 16 | 100 pcs | 60 x 100 pcs |
| GS1710161 | 10 cm x 10 cm | 17 | 16 | 100 pcs | 40 x 100 pcs |
| GS171020161 | 10 cm x 20 cm | 17 | 16 | 100 pcs | 20 x 100 pcs |

with X-ray detectable thread

| REF | dimensions | no. | no. | packaging unit | bulk packaging |
|--------------|-----------------|------------|---------|----------------|----------------|
| | | of threads | ofplies | paper bag | carton |
| *GS1705081X | 5 cm x 5 cm | 17 | 8 | 100 pcs | 80 x 100 pcs |
| GS1775081X | 7,5 cm x 7,5 cm | 17 | 8 | 100 pcs | 90 x 100 pcs |
| GS1710081X | 10 cm x 10 cm | 17 | 8 | 100 pcs | 40 x 100 pcs |
| *GS1705121X | 5 cm x 5 cm | 17 | 12 | 100 pcs | 96 x 100 pcs |
| GS1775121X | 7,5 cm x 7,5 cm | 17 | 12 | 100 pcs | 60 x 100 pcs |
| GS1710121X | 10 cm x 10 cm | 17 | 12 | 100 pcs | 20 x 100 pcs |
| GS171020121X | 10 cm x 20 cm | 17 | 12 | 100 pcs | 20 x 100 pcs |
| GS1705161X | 5 cm x 5 cm | 17 | 16 | 100 pcs | 96 x 100 pcs |
| GS1775161X | 7,5 cm x 7,5 cm | 17 | 16 | 100 pcs | 60 x 100 pcs |
| GS1710161X | 10 cm x 10 cm | 17 | 16 | 100 pcs | 40 x 100 pcs |
| GS171020161X | 10 cm x 20 cm | 17 | 16 | 100 pcs | 20 x 100 pcs |





vitahealth



*<

available

on request

 \checkmark self-adhesive labels

for medical

documentation

| REF | dimensions | no. of threads | no. of plies | packaging unit paper-foil / paper-paper | intermediate packaging box | bulk packaging carton |
|-------------------------------|-----------------|-------------------|-----------------|--|-------------------------------|--------------------------|
| GS1705081-S* | 5 cm x 5 cm | 17 | 8 | 1 pc. | 50 x 1 pc. | 30 x 50 pcs |
| GS1705082-S | 5 cm x 5 cm | 17 | 8 | 2 pcs | 50 x 2 pcs | 16 x 50 pcs |
| GS1705083-S | 5 cm x 5 cm | 17 | 8 | 3 pcs | 50 x 3 pcs | 16 x 50 pcs |
| GS1705085-S | 5 cm x 5 cm | 17 | 8 | 5 pcs | 20 x 5 pcs | 16 x 20 pcs |
| GS17050810-S | 5 cm x 5 cm | 17 | 8 | 10 pcs | 40 x 10 pcs | 10 x 40 pcs |
| GS17050820-SS | 5 cm x 5 cm | 17 | 8 | 20 pcs | 20 x 20 pcs | 10 x 20 pcs |
| GS1775081-S* | 7,5 cm x 7,5 cm | 17 | 8 | 1 pc. | 50 x 1 pc. | 30 x 50 pcs |
| GS1775082-S | 7,5 cm x 7,5 cm | 17 | 8 | 2 pcs | 50 x 2 pcs | 12 x 50 pcs |
| GS1775083-S | 7,5 cm x 7,5 cm | 17 | 8 | 3 pcs | 50 x 3 pcs | 12 x 50 pcs |
| GS1775085-S | 7,5 cm x 7,5 cm | 17 | 8 | 5 pcs | 20 x 5 pcs | 12 x 20 pcs |
| GS17750810-S | 7,5 cm x 7,5 cm | 17 | 8 | 10 pcs | 40 x 10 pcs | 10 x 40 pcs |
| GS17750820-SS | 7,5 cm x 7,5 cm | 17 | 8 | 20 pcs | 20 x 20 pcs | 10 x 20 pcs |
| GS1710081-S* | 10 cm x 10 cm | 17 | 8 | 1 pc. | 50 x 1 pc. | 30 x 50 pcs |
| GS1710082-SS | 10 cm x 10 cm | 17 | 8 | 2 pcs | 50 x 2 pcs | 8 x 50 pcs |
| GS1710083-S | 10 cm x 10 cm | 17 | 8 | 3 pcs | 50 x 3 pcs | 8 x 50 pcs |
| GS1710085-S | 10 cm x 10 cm | 17 | 8 | 5 pcs | 50 x 5 pcs | 8 x 50 pcs |
| GS17100810-S | 10 cm x 10 cm | 17 | 8 | 10 pcs | 40 x 10 pcs | 10 x 40 pcs |
| GS17100820-SS | 10 cm x 10 cm | 17 | 8 | 20 pcs | 10 x 20 pcs | 6 x 20 pcs |
| GS1705122-SS | 5 cm x 5 cm | 17 | 12 | 2 pcs | 50 x 2 pcs | 16 x 50 pcs |
| GS1705122-55 | 5 cm x 5 cm | 17 | 12 | 3 pcs | 50 x 2 pcs | 16 x 50 pcs |
| GS1705123-S-IP* | 5 cm x 5 cm | 17 | 12 | 3x1 pcs | 20 x 3 pcs | 20 x 20 pcs |
| GS1705125-SS | 5 cm x 5 cm | 17 | 12 | 1 | | |
| GS1705125-SS GS17051210-SS | 5 cm x 5 cm | | | 5 pcs | 20 x 5 pcs | 16 x 20 pcs |
| | | 17 | 12 | 10 pcs | 20 x 10 pcs | 12 x 20 pcs |
| GS1775122-SS | 7,5 cm x 7,5 cm | 17 | 12 | 2 pcs | 50 x 2 pcs | 18 x 50 pcs |
| GS1775123-SS | 7,5 cm x 7,5 cm | 17 | 12 | 3 pcs | 50 x 3 pcs | 12 x 50 pcs |
| GS1775123-S-IP* | 7,5 cm x 7,5 cm | 17 | 12 | 3x1 pcs | 20 x 3 pcs | 20 x 20 pcs |
| GS1775125-SS | 7,5 cm x 7,5 cm | 17 | 12 | 5 pcs | 20 x 5 pcs | 12 x 20 pcs |
| GS17751210-SS | 7,5 cm x 7,5 cm | 17 | 12 | 10 pcs | 40 x 10 pcs | 10 x 40 pcs |
| GS17751220-SS | 7,5 cm x 7,5 cm | 17 | 12 | 20 pcs | 20 x 20 pcs | 10 x 20 pcs |
| GS1710123-SS | 10 cm x 10 cm | 17 | 12 | 3 pcs | 50 x 3 pcs | 8 x 50 pcs |
| GS1710123-S-IP* | 10 cm x 10 cm | 17 | 12 | 3x1 pcs | 20 x 3 pcs | 20 x 20 pcs |
| GS1710125-S | 10 cm x 10 cm | 17 | 12 | 5 pcs | 20 x 5 pcs | 8 x 20 pcs |
| GS17101210-SS | 10 cm x 10 cm | 17 | 12 | 10 pcs | 20 x 10 pcs | 6 x 20 pcs |
| GS17101220-SS | 10 cm x 10 cm | 17 | 12 | 20 pcs | 20 x 20 pcs | 4 x 20 pcs |
| GS171020123-SS | 10 cm x 20 cm | 17 | 12 | 3 pcs | 50 x 3 pcs | 8 x 50 pcs |
| GS1710201210-SS | 10 cm x 20 cm | 17 | 12 | 10 pcs | 10 x 10 pcs | 8 x 10 pcs |
| GS171020125-SS | 10 cm x 20 cm | 17 | 12 | 5 pcs | 20 x 5 pcs | 8 x 20 pcs |
| GS1705163-SS | 5 cm x 5 cm | 17 | 16 | 3 pcs | 50 x 3 pcs | 16 x 50 pcs |
| GS1705165-SS | 5 cm x 5 cm | 17 | 16 | 5 pcs | 30 x 5 pcs | 12 x 30 pcs |
| GS1775163-SS | 7,5 cm x 7,5 cm | 17 | 16 | 3 pcs | 50 x 3 pcs | 12 x 50 pcs |
| GS1775165-SS | 7,5 cm x 7,5 cm | 17 | 16 | 5 pcs | 20 x 5 pcs | 12 x 20 pcs |
| GS17751610-SS | 7,5 cm x 7,5 cm | 17 | 16 | 10 pcs | 16 x 10 pcs | 12 x 16 pcs |
| GS1710163-SS | 10 cm x 10 cm | 17 | 16 | 3 pcs | 50 x 3 pcs | 8 x 50 pcs |
| GS1710165-SS | 10 cm x 10 cm | 17 | 16 | 5 pcs | 30 x 5 pcs | 6 x 30 pcs |
| GS17101610-SS | 10 cm x 10 cm | 17 | 16 | 10 pcs | 20 x 10 pcs | 8 x 20 pcs |
| GS171020163-SS | 10 cm x 20 cm | 17 | 16 | 3 pcs | 50 x 3 pcs | 4 x 50 pcs |
| GS1710201610-SS | 10 cm x 20 cm | 17 | 16 | 10 pcs | 10 x 10 pcs | 8 x 10 pcs |
| GS179020062-SS | 90 cm x 20 cm | 17 | 6 | 2 pcs | 40 x 2 pc. | 4 x 40 pcs |
| | 90 cm x 20 cm | 17 | 6 | 5 pcs | 16 x 5 pc. | 4 x 16 pcs |
| GS179020065-SS | 90 CHI X ZU CHI | | | | | |

with X-ray detectable thread

| REF | dimensions | no. of threads | no. of plies | packaging unit paper-foil | bulk packaging box |
|---------------------|-----------------|-------------------|-----------------|------------------------------|-----------------------|
| GS17750810-SSXR | 7,5 cm x 7,5 cm | 17 | 8 | 10 pcs | 40 x 10 pcs |
| GS17750820-SSXR | 7,5 cm x 7,5 cm | 17 | 8 | 20 pcs | 20 x 20 pcs |
| ✔ GS17100810-SSXR | 10 cm x 10 cm | 17 | 8 | 10 pcs | 40 x 10 pcs |
| * GS1705125-SSXR | 5 cm x 5 cm | 17 | 12 | 5 pcs | 20 x 5 pcs |
| * GS17051210-SSXR | 5 cm x 5 cm | 17 | 12 | 10 pcs | 40 x 10 pcs |
| ✔ GS1775125-SSXR | 7,5cm x 7,5cm | 17 | 12 | 5 pcs | 20 x 5 pcs |
| ✔ GS17751210-SSXR | 7,5cm x 7,5cm | 17 | 12 | 10 pcs | 10 x 10 pcs |
| ✔ GS17101210-SSXR | 10cm x 10cm | 17 | 12 | 10 pcs | 20 x 10 pcs |
| ✔ GS17101220-SSXR | 10cm x 10cm | 17 | 12 | 20 pcs | 20 x 20 pcs |
| ✔ GS17101240-SSXR | 10cm x 10cm | 17 | 12 | 40 pcs | 10 x 40 pcs |
| ✔ GS1710201210-SSXR | 10cm x 20cm | 17 | 12 | 10 pcs | 10 x 10 pcs |
| ✔ GS1710201220-SSXR | 10cm x 20cm | 17 | 12 | 20 pcs | 6 x 20 pcs |
| * GS17051610-SSXR | 5 cm x 5 cm | 17 | 16 | 10 pcs | 30 x 10 pcs |
| ✔ GS1775165-SSXR | 7,5 cm x 7,5 cm | 17 | 16 | 5 pcs | 20 x 5 pcs |
| ✔ GS17751610-SSXR | 7,5 cm x 7,5 cm | 17 | 16 | 10 pcs | 16 x 10 pcs |
| * GS1710165-SSXR | 10 cm x 10 cm | 17 | 16 | 5 pcs | 20 x 5 pcs |
| ✔ GS17101610-SSXR | 10 cm x 10 cm | 17 | 16 | 10 pcs | 20 x 10 pcs |
| | | | | | |

KOMPRI lux S

without X-ray detectable thread

GAUZE LAP SPONGES

features

- 100% hydrophilic cotton
- bleached with the use of hydrogen peroxide
- with a sewed in X-ray detectable thread or chip, which enables fast and non-invasive identification
- of the lap sponge left in the operating field with a belt to fix the lap sponge beyond the operating field

indication

Before used on open wounds, the device must be sterilized in accordance with a validated

sterilization method. It can be sterilized with steam or ethylene oxide

NONE TATER SHAFT STRATE

• as a surgical invasive dressing material for intraoperative absorption of large amounts of blood and body fluids from surgical wounds

gauze lap sponges

SERVI lux

with X-ray thread

| | REF | dimensions | no. of threads | no. of plies | packaging unit foil bag | bulk packaging carton |
|-----|-------------|---------------|-------------------|-----------------|----------------------------|-----------------------|
| NEW | LS1745458-X | 45 cm x 45 cm | 17 | 8 | 20 pcs | 10 x 20 pcs |
| NEW | LS1745704-X | 45 cm x 70 cm | 17 | 4 | 20 pcs | 10 x 20 pcs |
| NEW | LS1730304-X | 30 cm x 30 cm | 17 | 4 | 50 pcs | 10 x 50 pcs |
| NEW | LS1745454-X | 45 cm x 45 cm | 17 | 4 | 20 pcs | 10 x 20 pcs |
| NEW | LS1775904-X | 75 cm x 90 cm | 17 | 4 | 20 pcs | 6 x 20 pcs |

with X-ray thread and belt

| | REF | dimensions | no. of threads | no. of plies | packaging unit foil bag | bulk packaging carton |
|-----|--------------|---------------|-------------------|-----------------|----------------------------|-----------------------|
| | LS1745454-XB | 45 cm x 45 cm | 17 | 4 | 20 pcs | 10 x 20 pcs |
| | LS1745704-XB | 45 cm x 70 cm | 17 | 4 | 20 pcs | 10 x 20 pcs |
| | LS1745456-XB | 45 cm x 45 cm | 17 | 6 | 20 pcs | 10 x 20 pcs |
| | LS1745706-XB | 45 cm x 70 cm | 17 | 6 | 20 pcs | 10 x 20 pcs |
| | LS1745458-XB | 45 cm x 45 cm | 17 | 8 | 20 pcs | 10 x 20 pcs |
| | LS2045454-XB | 45 cm x 45 cm | 20 | 4 | 20 pcs | 10 x 20 pcs |
| NEW | LS1775904-XB | 75 cm x 90 cm | 17 | 4 | 20 pcs | 6 x 20 pcs |
| NEW | LS1730454-XB | 30 cm x 45 cm | 17 | 4 | 50 pcs | 10 x 50 pcs |
| NEW | LS1730306-XB | 30 cm x 30 cm | 17 | 6 | 50 pcs | 10 x 50 pcs |
| NEW | LS1775752-XB | 75 cm x 75 cm | 17 | 2 | 20 pcs | 10 x 20 pcs |
| NEW | LS1730304-XB | 30 cm x 30 cm | 17 | 4 | 50 pcs | 10 x 50 pcs |

pre-washed with X-ray chip and belt

| REF | dimensions | no. of threads | no. of plies | packaging unit foil bag | bulk packaging carton |
|-----------------|---------------|-------------------|-----------------|----------------------------|-----------------------|
| LS2045454-XB-PW | 45 cm x 45 cm | 20 | 4 | 20 pcs | 10 x 20 pcs |



The advantages of self-adhesive labels for medical documentation: - fast data transport without the need to rewrite, - safe identification of the device -elimination of the possibility of errors, - improvement of the preparation of medical documentation.



92 www.zarys.com

| with X-ray threa | d and belt | for medical documentation | | |
|-------------------|---------------|---------------------------|----------|--|
| REF | dimensions | no. of threads | r s c | |
| ✓ LS17454542-SSXB | 45 cm x 45 cm | 17 | 4 | |
| ✔ LS17454545-SSXB | 45 cm x 45 cm | 17 | 4 | |
| LS17457042-SSXB | 45 cm x 70 cm | 17 | 4 | |

SERVI lux S

| ith X-ray thr | self-adhesive labels for medical documentation | | |
|---------------|--|-------------------|-------------|
| REF | dimensions | no. of threads | no. of p |

| REF | dimensions | no. of threads | no. of plies | packaging unit paper-foil | bulk packaging box |
|-------------------|---------------|-------------------|-----------------|------------------------------|-----------------------|
| ✓ LS17454542-SSXB | 45 cm x 45 cm | 17 | 4 | 2 pcs | 20 x 2 pcs |
| ✓ LS17454545-SSXB | 45 cm x 45 cm | 17 | 4 | 5 pcs | 8 x 5 pcs |
| ✓ LS17457042-SSXB | 45 cm x 70 cm | 17 | 4 | 2 pcs | 20 x 2 pcs |
| ✓ LS17454563-SSXB | 45 cm x 45 cm | 17 | 6 | 3 pcs | 10 x 3 pcs |
| ✓ LS17454582-SSXB | 45 cm x 45 cm | 17 | 8 | 2 pcs | 10 x 2 pcs |
| ✔ LS17454585-SSXB | 45 cm x 45 cm | 17 | 8 | 5 pcs | 5 x 5 pcs |
| ✔ LS17457043-SSXB | 45 cm x 70 cm | 17 | 4 | 3 pcs | 10 x 3 pcs |
| ✓ LS17759041-SSXB | 75 cm x 90 cm | 17 | 4 | 1 pc. | 15 x 1 pcs |

pre-washed with X-ray chip and belt

| | REF | dimensions | no. of threads | no. of plies | packaging unit foil bag | bulk packaging carton |
|-----|--------------------|---------------|-------------------|-----------------|----------------------------|-----------------------|
| NEW | LS17454542-SSXB-PW | 45 cm x 45 cm | 17 | 4 | 2 pcs | 20 x 2 pcs |
| NEW | LS17454543-SSXB-PW | 45 cm x 45 cm | 17 | 4 | 3 pcs | 15 x 3 pcs |
| NEW | LS17454545-SSXB-PW | 45 cm x 45 cm | 17 | 4 | 5 pcs | 10 x 5 pcs |



gauze balls

TUPFER lux

features

indication

NEW

NEW NEW

without X-ray thread

with X-ray thread

| REF | size of gauze | diameter approx. | shape | packaging unit foil bag | bulk packaging carton |
|------------|---------------|---------------------|-------|----------------------------|-----------------------|
| TK1515250 | 15 cm x 15 cm | 2,5 cm | ball | 250 pcs | 24 x 250 pcs |
| TK2020250 | 20 cm x 20 cm | 3 cm | ball | 250 pcs | 24 x 250 pcs |
| TK3030100 | 30 cm x 30 cm | 4 cm | ball | 100 pcs | 10 x 100 pcs |
| TFT1515100 | 15 cm x 15 cm | 2,5 cm | bean | 1000 pcs | 60 x 1000 pcs |
| TF1515100 | 15 cm x 15 cm | 2,5 cm | bean | 100 pcs | 60 x 100 pcs |

| Before use on open wounds, the device must |
|--|
| be sterilized in accordance with a validated |
| sterilization method. It can be sterilized |
| with steam or ethylene oxide. |



| REF | size of gauze | diameter approx. | shape | packaging unit foil bag | bulk packaging carton |
|-------------|---------------|---------------------|-------|----------------------------|-----------------------|
| TK1515250X | 15 cm x 15 cm | 2,5 cm | ball | 250 pcs | 24 x 250 pcs |
| TK2020250X | 20 cm x 20 cm | 3 cm | ball | 250 pcs | 24 x 250 pcs |
| TK2525250X | 25 cm x 25 cm | 3,5 cm | ball | 250 pcs | 8 x 250 pcs |
| TK3030100X | 30 cm x 30 cm | 4 cm | ball | 100 pcs | 10 x 100 pcs |
| TK3040100X | 30 cm x 40 cm | 4,5 cm | ball | 100 pcs | 10 x 100 pcs |
| TK5050100X | 50 cm x 50 cm | 5 cm | ball | 100 pcs | 10 x 100 pcs |
| TFT1515100X | 15 cm x 15 cm | 2,5 cm | bean | 1000 pcs | 10 x 1000 pcs |
| TF1515100X | 15 cm x 15 cm | 2,5 cm | bean | 100 pcs | 60 x 100 pcs |





100% hydrophilic cotton
bleached with the use of hydrogen peroxide
available in a version with or without X-ray detectable thread, which enables fast and non-invasive identification of the gauze ball left in the operating field

• absorption of blood and body fluids from surgical wounds and hard to reach places



| size of gauze | diameter approx. | shape | packaging unit paper-foil | bulk packaging carton |
|-----------------|---------------------|------------|------------------------------|--------------------------|
| 9,5 cm x 9,5 cm | | bean | 5 pcs | 50 x 5 pcs |
| 9,5 cm x 9,5 cm | | bean | 10 pcs | 25 x 10 pcs |
| 15 cm x 15 cm | | bean | 5 pcs | 25 x 5 pcs |
| 15 cm x 15 cm | | bean | 10 pcs | 20 x 10 pcs |
| 12 cm x 12 cm | | triangular | 5 pcs | 15 x 5 pcs |
| 18 cm x 18 cm | | triangular | 10 pcs | 15 x 10 pcs |
| 15 cm x 15 cm | 2,5 cm | ball | 10 pcs | 20 x 10 pcs |
| 20 cm x 20 cm | 3 cm | ball | 10 pcs | 20 x 10 pcs |
| 30 cm x 30 cm | 4 cm | ball | 10 pcs | 10 x 10 pcs |
| | | | | |



non-woven swabs

NONVI lux

Before use on open wounds, the device must be sterilized in accordance with a validated sterilization method. It can be sterilized with steam, ethylene oxide or radiation.

| REF | dimensions | grammage | packaging unit paper bag | bulk packaging carton |
|--------------|-----------------|---------------------|-----------------------------|--------------------------|
| NL50-100 | 5 cm x 5 cm | 30 g/m ² | 100 pcs | 50 x 100 pcs |
| NL50-100-40 | 5 cm x 5 cm | 40 g/m ² | 100 pcs | 50 x 100 pcs |
| NL75-100 | 7,5 cm x 7,5 cm | 30 g/m ² | 100 pcs | 40 x 100 pcs |
| NL75-100-40 | 7,5 cm x 7,5 cm | 40 g/m ² | 100 pcs | 40 x 100 pcs |
| NL100-100 | 10 cm x 10 cm | 30 g/m ² | 100 pcs | 40 x 100 pcs |
| NL100-100-40 | 10 cm x 10 cm | 40 g/m ² | 100 pcs | 20 x 100 pcs |
| NL1020-100 | 10 cm x 20 cm | 30 g/m ² | 100 pcs | 20 x 100 pcs |

NONVI lux S



NONE TERME

NON-WOVEN SWABS

features

- hydrophilic medical non-woven, 30 or 40 g/m²
- 4-plies
- soft, non-fraying, durable

indication

 general wound dressing in out-patient clinics, surgical sites, procedure department





| Υ | *< |
|--------------|------------|
| swab with | available |
| "Y" incision | on request |

| REF | dimensions | grammage | packaging unit paper-foil | intermediate packaging box | bulk packaging box |
|-------------|-----------------|---------------------|------------------------------|-------------------------------|-----------------------|
| NL052-S | 5 cm x 5 cm | 30 g/m ² | 2 pcs | 25 x 2 pcs | 60 x 25 pcs |
| NL053-S | 5 cm x 5 cm | 30 g/m ² | 3 pcs | 25 x 3 pcs | 45 x 25 pcs |
| NL055-S | 5 cm x 5 cm | 30 g/m ² | 5 pcs | 25 x 5 pcs | 36 x 25 pcs |
| NL0510-S | 5 cm x 5 cm | 30 g/m ² | 10 pcs | 10 x 10 pcs | 32 x 10 pcs |
| NL052-S-40 | 5 cm x 5 cm | 40 g/m ² | 2 pcs | 25 x 2 pcs | 60 x 25 pcs |
| NL053-S-40 | 5 cm x 5 cm | 40 g/m ² | 3 pcs | 25 x 3 pcs | 45 x 25 pcs |
| NL055-S-40 | 5 cm x 5 cm | 40 g/m ² | 5 pcs | 25 x 5 pcs | 24 x 25 pcs |
| NL0510-S-40 | 5 cm x 5 cm | 40 g/m ² | 10 pcs | 10 x 10 pcs | 32 x 10 pcs |
| NL752-S | 7,5 cm x 7,5 cm | 30 g/m ² | 2 pcs | 25 x 2 pcs | 48 x 25 pcs |
| NL753-S | 7,5 cm x 7,5 cm | 30 g/m ² | 3 pcs | 25 x 3 pcs | 36 x 25 pcs |
| NL755-S | 7,5 cm x 7,5 cm | 30 g/m ² | 5 pcs | 25 x 5 pcs | 24 x 25 pcs |
| NL7510-S | 7,5 cm x 7,5 cm | 30 g/m ² | 10 pcs | 10 x 10 pcs | 12 x 10 pcs |
| NL752-S-40 | 7,5 cm x 7,5 cm | 40 g/m ² | 2 pcs | 25 x 2 pcs | 48 x 25 pcs |
| NL753-S-40 | 7,5 cm x 7,5 cm | 40 g/m ² | 3 pcs | 25 x 3 pcs | 24 x 25 pcs |
| NL755-S-40 | 7,5 cm x 7,5 cm | 40 g/m ² | 5 pcs | 25 x 5 pcs | 18 x 25 pcs |
| NL7510-S-40 | 7,5 cm x 7,5 cm | 40 g/m ² | 10 pcs | 10 x 10 pcs | 12 x 10 pcs |
| NL7520-S-40 | 7,5 cm x 7,5 cm | 40 g/m ² | 20 pcs | 5 x 20 pcs | 12 x 5 pcs |
| NL102-S | 10 cm x 10 cm | 30 g/m ² | 2 pcs | 25 x 2 pcs | 30 x 25 pcs |
| NL103-S | 10 cm x 10 cm | 30 g/m ² | 3 pcs | 25 x 3 pcs | 24 x 25 pcs |
| NL105-S | 10 cm x 10 cm | 30 g/m ² | 5 pcs | 20 x 5 pcs | 12 x 20 pcs |
| NL1010-S | 10 cm x 10 cm | 30 g/m ² | 10 pcs | 10 x 10 pcs | 12 x 10 pcs |
| NL102-S-40 | 10 cm x 10 cm | 40 g/m ² | 2 pcs | 25 x 2 pcs | 30 x 25 pcs |
| NL103-S-40 | 10 cm x 10 cm | 40 g/m ² | 3 pcs | 25 x 3 pcs | 18 x 25 pcs |
| NL105-S-40 | 10 cm x 10 cm | 40 g/m ² | 5 pcs | 20 x 5 pcs | 12 x 25 pcs |
| NL1010-S-40 | 10 cm x 10 cm | 40 g/m ² | 10 pcs | 10 x 10 pcs | 12 x 25 pcs |
| NL1020-S-40 | 10 cm x 10 cm | 40 g/m ² | 20 pcs | 5 x 20 pcs | 12 x 25 pcs |

ABSORBENT WOUND DRESSINGS

features

- layered structure: non-woven cellulose cellulose pulp
- increased absorption capacity due to the high absorbent cellulose layer
- bottom layer of the dressing covered with waterproof blue non-woven

indication

removing the exudate from strongly oozing wounds
covering wounds after surgical procedures

ABD pads, multi-layer

MULTIabsorb

| Before use on open wounds, the device must be sterilized in accordance with a validated sterilization method. It can be sterilized with ethylene oxide or radiation. | | | | | | |
|--|---------------|----------------------------|--------------------------|--|--|--|
| REF | dimensions | packaging unit foil bag | bulk packaging carton | | | |
| ABD1010 | 10 cm x 10 cm | 50 pcs | 40 x 50 pcs | | | |
| ABD1020 | 10 cm x 20 cm | 50 pcs | 20 x 50 pcs | | | |
| ABD2025 | 20 cm x 25 cm | 50 pcs | 12 x 50 pcs | | | |



MULTIabsorb S

STERILEEO

REF dimensions

| ABD1010-S | 10 cm x 10 cm |
|-----------|---------------|
| ABD1020-S | 10 cm x 20 cm |
| ABD2025-S | 20 cm x 25 cm |
| ABD2040-S | 20 cm x 40 cm |
| | |

NONVI lux S with Y-incision

features

- hydrophilic medical non-woven 40 g/m²
- 4-plies • Y-incision for precise, aseptic protection of the outlet of various types of medical drains

indication


NEW

(-s)

• securing the outlet of drains, catheters, tracheostomy tubes etc.

| REF sterilized by EO | REF sterilized by steam | dimensions | packaging unit paper-foil | bulk packaging box |
|----------------------|-------------------------|---------------|------------------------------|-----------------------|
| NL105-S-40Y | * NL105-SS-40Y | 10 cm x 10 cm | 5 pcs | 20 x 5 pcs |
| | NL104-SS-Y | 10 cm x 10 cm | 4 pcs | 30 x 4 pcs |





| packaging unit paper-foil | intermediate packaging box | bulk packaging carton |
|------------------------------|-------------------------------|-----------------------|
| 1 pc. | 25 x 1 pc. | 24 x 25 pcs |
| 1 pc. | 25 x 1 pc. | 18 x 25 pcs |
| 1 pc. | 15 x 1 pc. | 18 x 15 pcs |
| 1 pc. | 8 x 1 pc. | 12 x 8 pcs |

EYE DRESSINGS

features

- layered structure: non-woven + cellulose pulp
- needs to be fixed with a secondary dressing e.g. with adhesive tape or tubular net bandage
- hypoallergenic

indication

· covering an eye and its area after ophthalmic procedures and surgical procedures sight correction in strabismus therapy



eye pads

elastoLUMENAL

| REF | dimensions | packaging unit paper-paper | intermediate packaging box | bulk packaging carton |
|--------|-----------------|-------------------------------|-------------------------------|--------------------------|
| 803311 | 5,5 cm x 7,5 cm | 1 pc. | 50 x 1 pc. | 24 x 50 pcs |


elastoLUMENALS



| Conforms well to facial contours thanks to the oval shape. | | | | |
|--|-----------------|-------------------------------|-------------------------------|-----------------------|
| REF | dimensions | packaging unit paper-paper | intermediate packaging box | bulk packaging carton |
| 803310 | 5,5 cm x 7,5 cm | 1 pc. | 50 x 1 pc. | 24 x 50 pcs |









| REF | dimensions | packaging unit foil bag | bulk packaging carton |
|--------|-------------|----------------------------|-----------------------|
| KC4050 | 4 cm x 5 cm | 2 rolls (2 x 500 pcs) | 16 x 2 rolls |



features

• regular perforations to facilitate tearing off individual pieces as small swabs • each roll contains 500 pieces of individual swabs

indication

 device used for absorbing minor exudates, dressing minor cuts, abrasions on the skin and injection sites

ADHESIVE TAPES AND PLASTERS WITH DRESSING




ADHESIVE TAPES

features

- flexible, conforms to the body's contours, guarantees freedom of movement
- characterized by adequate air and water vapuor permeability
- good adhesion
- $\cdot\,\tilde{\mathsf{can}}$ be divided lengthwise and widthwise without the use of scissors
- hypoallergenic

indication

- stabilization, protection and fixation of primary, non-adhesive dressings, applied directly to wounds and skin damages
 stabilization, protection and fixation of drains, catheters, medical tubes that do not require specialized fixing






adhesive fabric tape, white

| 1,25 cm x 5 m | 24 pcs | 30 x 24 pcs |
|------------------|--|---|
| 2,5 cm x 5 m | 12 pcs | 30 x 12 pcs |
| 5 cm x 5 m | 6 pcs | 30 x 6 pcs |
| 1,25 cm x 9,14 m | 24 pcs | 30 x 24 pcs |
| 2,5 cm x 9,14 m | 12 pcs | 30 x 12 pcs |
| 5 cm x 9,14 m | 6 pcs | 30 x 6 pcs |
| 2,5 cm x 5 m | 12 pcs | 30 x 12 pcs |
| | 2,5 cm x 5 m 5 cm x 5 m 1,25 cm x 9,14 m 2,5 cm x 9,14 m 5 cm x 9,14 m | 2,5 cm x 5 m 12 pcs 5 cm x 5 m 6 pcs 1,25 cm x 9,14 m 24 pcs 2,5 cm x 9,14 m 12 pcs 5 cm x 9,14 m 6 pcs |





| | REF | dimensions | packaging unit box | bulk packaging carton |
|-----|--------------|---------------|-----------------------|-----------------------|
| | SZO12505S | 1,25 cm x 5 m | 24 pcs | 30 x 24 pcs |
| | SZO25005S | 2,5 cm x 5 m | 12 pcs | 30 x 12 pcs |
| | SZO50005S | 5 cm x 5 m | 6 pcs | 30 x 6 pcs |
| NEW | SZO25005S-L* | 2,5 cm x 5 m | 12 pcs | 30 x 12 pcs |

Ĺ.J

on a roll/

spooled'

 \bigcirc

beige colour











fabric

adhesive silk tape

SILKplast

| REF | dimensions | packaging unit box | bulk packaging carton |
|----------|------------------|-----------------------|--------------------------|
| SS12505W | 1,25 cm x 5 m | 24 pcs | 30 x 24 pcs |
| SS25005W | 2,5 cm x 5 m | 12 pcs | 30 x 12 pcs |
| SS50005W | 5 cm x 5 m | 6 pcs | 30 x 6 pcs |
| SS12509W | 1,25 cm x 9,14 m | 24 pcs | 30 x 24 pcs |
| SS25009W | 2,5 cm x 9,14 m | 12 pcs | 30 x 12 pcs |
| SS50009W | 5 cm x 9,14 m | 6 pcs | 30 x 6 pcs |

on a roll



 \sim \sim serrated



adhesive fabric tape, beige colour

PLASTIplast

adhesive non-woven tape

| dimensions | packaging unit box | bulk packaging carton |
|------------------|-----------------------|-----------------------|
| 1,25 cm x 5 m | 24 pcs | 30 x 24 pcs |
| 2,5 cm x 5 m | 12 pcs | 30 x 12 pcs |
| 5 cm x 5 m | 6 pcs | 30 x 6 pcs |
| 1,25 cm x 9,14 m | 24 pcs | 30 x 24 pcs |
| 2,5 cm x 9,14 m | 12 pcs | 30 x 12 pcs |
| 5 cm x 9,14 m | 6 pcs | 30 x 6 pcs |










edge





adhesive PE tape

FILMplast





PLASTERS WITH DRESSING

features

- self-adhesive plasters with absorbent pads
- centrally situated absorbent pad absorbs a small and medium amount of exudate and protects from external factors
- absorbent pad coated with a polyethylene mesh that prevents sticking to the wound
- characterized by adequate air and water vapuor permeability
- flexible, conforms to the body's contours, guarantees freedom of movement
- silicone backing paper facilitates precise, painless and efficient application

• protection of small skin injuries, cuts, abrasions and lacerations

hypoallergenic

indication

dimensions

6 cm x 1 m

6 cm x 5 m

8 cm x 1 m

8 cm x 5 m

 $\uparrow \uparrow \uparrow$

structure

NONVIplast

REF



NP6100 NP6500 NP8100 NP8500 Ð

NON DATES STATES

made of hydrophobic non-woven

fabric plaster with dressing

FABRIplast

| REF | dimensions | packaging unit box | bulk packaging carton |
|--------|------------|-----------------------|--------------------------|
| FP6100 | 6 cm x 1 m | 1 pc. | 200 x 1 pc. |
| FP6500 | 6 cm x 5 m | 1 pc. | 100 x 1 pc. |
| FP8100 | 8 cm x 1 m | 1 pc. | 200 x 1 pc. |
| FP8500 | 8 cm x 5 m | 1 pc. | 100 x 1 pc. |



ተተተ air permeable structure




adhesive perforated non-woven tape

POREplast

| - | REF | dimensions | colour | packaging unit box / cover* | bulk packaging carton |
|-----|-----------|-----------------|--------|--------------------------------|-----------------------|
| | SPW25009W | 2,5 cm x 9,14 m | white | 12 pcs | 30 x 12 pcs |
| NEW | VH-POPP * | 2,5 cm x 5 m | pink | 1 рс | 24 x 24 pcs |
| NEW | VH-POPG* | 2,5 cm x 5 m | green | 1 pc | 24 x 24 pcs |



NONE TATEX STRATES

















vitahealth

non-woven plaster with dressing

| packaging unit box | bulk packaging carton |
|-----------------------|-----------------------|
| 1 pc. | 200 x 1 pc. |
| 1 pc. | 100 x 1 pc. |
| 1 pc. | 200 x 1 pc. |
| 1 рс. | 100 x 1 pc. |

 \Diamond œ air permeable possibility white of cutting colour





SELF-ADHESIVE NON-WOVEN DRESSINGS

DRESSING TAPE

features

- made of hydrophobic non-woven
- covered with hypoallergenic acrylic glue
- micropores in the structure of the non-woven provide adequate vapor permeability, thanks to which
- a proper gas exchange between the dressing and the skin takes place • flexible, conforms to the body's contours, guarantees freedom of movement
- backing layer made of silicone paper

indication

• stabilization, protection and fixation of primary, non-adhesive dressings, applied directly to wounds and skin damages • stabilization, protection and fixation of drains, catheters, medical tubes

elastopor

| REF | dimensions |
|--------|---------------|
| 806009 | 2,5 cm x 10 m |
| 806010 | 5 cm x 10 m |
| 806011 | 10 cm x 10 m |
| 806012 | 15 cm x 10 m |
| 806013 | 20 cm x 10 m |
| 806014 | 25 cm x 10 m |
| 806015 | 30 cm x 10 m |
| | |







of backing paper

elastopor E

| REF | dimensions | packaging unit box | bulk packaging carton |
|-----------|---------------|-----------------------|-----------------------|
| * 806009E | 2,5 cm x 10 m | 1 pc. | 200 x 1 pc. |
| 806010E | 5 cm x 10 m | 1 pc. | 100 x 1 pc. |
| 806011E | 10 cm x 10 m | 1 pc. | 50 x 1 pc. |
| 806012E | 15 cm x 10 m | 1 pc. | 50 x 1 pc. |
| 806013E | 20 cm x 10 m | 1 pc. | 25 x 1 pc. |
| 806014E | 25 cm x 10 m | 1 pc. | 25 x 1 pc. |
| 806015E | 30 cm x 10 m | 1 pc. | 25 x 1 pc. |

| wavy | |
|------------------|--|
| intersection | |
| of backing paper | |
| | |



SPECIA DRESSINGS











cutting





• backing paper with metric scale facilitates the division of the tape and precise, painless and efficient application • packed in a cardboard box - a dispenser that allows dosing without removing the tape from the packaging

adhesive non-woven dressing tape

| packaging unit box | bulk packaging carton |
|-----------------------|-----------------------|
| 1 pc. | 200 x 1 pc. |
| 1 pc. | 100 x 1 pc. |
| 1 pc. | 50 x 1 pc. |
| 1 pc. | 50 x 1 pc. |
| 1 pc. | 25 x 1 pc. |
| 1 pc. | 25 x 1 pc. |
| 1 pc. | 25 x 1 pc. |





SELF-ADHESIVE DRESSINGS WITH ABSORBENT PAD

features

Ŷ

made of

hydrophobic

non-woven

 \Diamond

white

colour

 $\uparrow\uparrow\uparrow$

air permeable

structure

- self-adhesive dressings with centrally situated absorbent pad, made of hydrophobic non-woven, covered with acrylic glue • absorbent pad coated with a polyethylene mesh that prevents sticking to the wound, absorbs a small and medium amount of exudate and protects against the influence of external factors
- rounded edges prevent accidental removal of the dressing
- silicone backing paper facilitates the precise, painless and efficient application
- micropores in the structure of the non-woven provide adequate vapor permeability, thanks to which a proper gas exchange
- between the dressing and the skin takes place
- hypoallergenic

non-woven dressing with absorbent pad

elastopor STERIL 📀 BESTSELLER

indication

• especially recommended for covering surgical wounds • for treating wounds with a small and medium amount of exudate

| REF | | dimensions | packaging unit paper-paper | intermediate packaging box | bulk packaging carton |
|--------|---------|---------------|-------------------------------|-------------------------------|-----------------------|
| 801001 | 801001E | 5 cm x 7,2 cm | 1 pc. | 100 x 1 pc. | 24 x 100 pcs |
| 801012 | 801012E | 6 cm x 10 cm | 1 pc. | 50 x 1 pc. | 45 x 50 pcs |
| 801002 | 801002E | 8 cm x 10 cm | 1 pc. | 30 x 1 pc. | 36 x 30 pcs |
| 801013 | 801013E | 10 cm x 10 cm | 1 pc. | 30 x 1 pc. | 24 x 30 pcs |
| 801021 | 801021E | 10 cm x 12 cm | 1 pc. | 30 x 1 pc. | 30x 30 pcs |
| 801014 | 801014E | 10 cm x 15 cm | 1 pc. | 30 x 1 pc. | 24 x 30 pcs |
| 801003 | 801003E | 6 cm x 15 cm | 1 pc. | 50 x 1 pc. | 24 x 50 pcs |
| 801004 | 801004E | 8 cm x 15 cm | 1 pc. | 30 x 1 pc. | 24 x 30 pcs |
| 801005 | 801005E | 8 cm x 20 cm | 1 pc. | 30 x 1 pc. | 18 x 30 pcs |
| 801015 | 801015E | 10 cm x 20 cm | 1 pc. | 25 x 1 pc. | 18 x 25 pcs |
| 801016 | 801016E | 10 cm x 25 cm | 1 pc. | 25 x 1 pc. | 18 x 25 pcs |
| 801017 | 801017E | 10 cm x 30 cm | 1 pc. | 25 x 1 pc. | 12 x 25 pcs |
| 801018 | 801018E | 10 cm x 35 cm | 1 pc. | 25 x 1 pc. | 12 x 25 pcs |
| | | | | | |

indexes without the ending E - cutting the securing paper along the long edge (exception 801001; 801012) indexes with the ending E - the intersection of the securing paper for all indexes along the short edge



non-woven dressing with absorbent pad

elastoKIDS STERIL

indication • especially recommended for covering surgical wounds • treating wounds waith a small and medium amount of exudate

| packaging | |
|-----------|--|
| | |
| 50 pcs | |
| 30 pcs | |
| 25 pcs | |
| | 50 pcs 30 pcs 25 pcs |

SELF-ADHESIVE DRESSINGS FOR KIDS

features

 \Diamond

white

colour

- self-adhesive dressing with absorbent pad, made of hydrophobic non-woven material with color overprint, covered by hypoallergenic glue
- absorbent pad coated with polyethylene mesh to prevent sticking to the wound, absorbing small and medium amounts of exudate and protecting against external influences
- rounded edges to prevent accidental removal of the dressing • protective layer of silicone paper for precise, painless and effective application
- micropores in the non-woven material structure provide an adequate rate of moisture vapour permeability so correct gas exchange between the dressing and skin is guaranteed



non-woven IV cannula dressing

elastoKIDS IV

features

sticking to the injection site from causing damage to the skin

indication

| REF | dimensions | packaging unit paper-paper | intermediate packaging box |
|-----|------------|-------------------------------|-------------------------------|
| | | | |

• centrally situated absorbent pad protects the non-woven material from • non-woven lining should be placed under cannula's wings, which prevents

• stabilization, protection and fixation of intravenous cannulas

| | dimensions | packaging unit paper-paper | intermediate packaging box | bulk packaging carton |
|---|-----------------|-------------------------------|-------------------------------|-----------------------|
| | | | | |
| E | 6 cm x 8 cm | 1 pc. | 100 x 1 pc. | 20 x 100 pcs |
| Е | 5,1 cm x 7,6 cm | 1 pc. | 100 x 1 pc. | 20 x 100 pcs |
| E | 6 cm x 8 cm | 1 pc. | 50 x 1 pc. | 40 x 50 pcs |
| Е | 5,1 cm x 7,6 cm | 1 pc. | 50 x 1 pc. | 40 x 50 pcs |
| | | | | |

• centrally situated absorbent pad protects the non-woven material from

• non-woven lining should be placed under cannula's wings, which prevents it

• application of the dressing around the cannula

• stabilization, protection and fixation of intravenous cannulas





non-woven dressing with absorbent pad, incision and O-hole

elastopor STERILD

indication

- dressing application around the inserted drain
- securing the outlet of drains and catheters · for treating wounds with a small and medium amount of exudate



| REF dimensions | | dimensions | packaging unit paper-paper | intermediate packaging box | bulk packaging carton |
|-----------------------------------|--------------------|---------------|-------------------------------|-------------------------------|-----------------------|
| 801031 | 801031E | 9 cm x 10 cm | 1 pc. | 30 x 1 pc. | 24 x 30 pcs |
| 801030 | 801030E | 12 cm x 14 cm | 1 pc. | 25 x 1 pc. | 16 x 25 pcs |
| to all a constants and the second | t the second sec E | | | | |

NEW

xes without the ending E - cutting the securing paper along the long edge (exception 801001; 801012) indexes with the ending E - the intersection of the securing paper for all indexes along the short edge



vitahealth

NON DATEX STEALE

wound

features

elastoKIDS

- and proper application
- hypoallergenic
- for wound healing
- each plaster individually protected by paper-paper packaging

indication

| REF | dimensions | packaging unit paper-paper | intermediate packaging box |
|------------|-----------------|-------------------------------|-------------------------------|
| VH-PDP7219 | 7,2 cm x 1,9 cm | 1 pc. | 100 x 1 pc. |

elastoKIDS STERIL D

indication



| protoction | | | | | | |
|------------|--------------|-------------------------------|-------------------------------|--|--|--|
| | | | | | | |
| REF | dimensions | packaging unit paper-paper | intermediate packaging box | | | |
| 801031A | 9 cm x 10 cm | 1 pc. | 30 x 1 pc. | | | |




| 801031A 9 cm x 10 cm 1 pc. | | | 1.1.1.1 |
|----------------------------|---------|--------------|---------|
| 801031A 9 cm x 10 cm 1 pc. | | | |
| | 801031A | 9 cm x 10 cm | 1 pc. |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

protection of the site of drain or catheter application









set of injection plasters

features

- made of polyurethane film transparent
- and protecting from external influences
- removal of the dressing
- and proper application
- hypoallergenic
- covered with hypoallergenic acrylic glue
 good air and steam permeability provides the right conditions to promote wound healing processes

indication

of exudate

| REF | dimensions | packaging unit paper-paper | intermediate packaging box |
|--------|-----------------|-------------------------------|-------------------------------|
| PI2572 | 2,5 cm x 7,2 cm | 1 pc. | 100 x 1 pc. |

HOW TO CHOOSE A DRESSING PROPERLY?

When choosing the proper type of dressing, pay special attention to:

- absorption rate and absorption capacity, - vapor permeability,
- possibility of an atraumatic dressing change,
- safety of use.

 $\Diamond \Diamond \Diamond$

white, beige

or colourful

Colorful kids plasters

made of polyethylene film with colourful print

- centrally situated absorbent pad that absorbs small and medium amounts of exudate and protects against external factors
- absorbent pad coated with polyethylene mesh to prevent sticking to the
- flexible, conforms to the body contours, guarantees freedom of movement
- good air and water vapour permeability provides the right conditions
- for protection against contamination of small wounds, skin injuries, cuts, abrasions, lacerations with small to medium amount of exudate







- centrally placed absorbent pad absorbs small and medium wound exudate
- absorbent pad coated with polyethylene mesh protects against external factors' influence and does not stick to the wound as well as facilitates
- fexible, adapting to the shape of the body, provides freedom of movement
- silcone paper backing facilitates precise, painless and efficient application
- for protection against contamination of minor wounds, injuries, cuts, abrasions, cuts, lacerations on the skin with a small and medium amount



elastoFIXAL

features a part fastening the drain stabilization • covered with hypoallergenic glue indication

REF dimensions 816101 9 cm x 4 cm

elastoFIXAL KIDS

| ea | tu | re | S | |
|-----|----|-----|---|--|
| ~~~ | | ••• | ~ | |

- two-part fixing tape consisting of a part attached to the skin and a part fastening the drain
- stabilization • possibility of reposition and correction of the drain position thanks to Velcro micropores in the non-woven fabric structure provide appropriate vapour-permeability rate, thanks to which proper gas exchange between the dressing and the skin takes place
- elastic adheres well, conforms to the body contours

indication

catheters and tubes

| REF | dimensions | packaging unit paper-foil | intermediate packaging box | |
|---------|-------------|------------------------------|-------------------------------|--------------------------------|
| 816101A | 9 cm x 4 cm | 1 pc. | 50 x 1 pc. | bulk packaging: 20 x 50 pcs |

adhesive tape for endotracheal tube elastoFIXAL Set

- part with Velcro on the other side
- the possibility of repositioning and correcting the position of the tube thanks to
- pads and strap covered with hypoallergenic glue • has excellent adhesive properties and high mechanical strength which prevents premature and accidental detachment from the skin

indication

· long-term stabilization of the endotracheal or oropharyngeal tube

REF kit composition 816102 2 self-adhesive non 1 fastening strap

108 www.zarys.com

non-woven nasal cannula fixation plaster

elastoNASAL

features

- flexible fixation plaster made of non-woven
- covered with acrylic glue • adequate vapor permeability, thanks to which a proper gas exchange between the dressing and the skin takes place
- thin and flexible plaster is comfortable for the patient, conforms to the body's contours
- "two-stage" application type facilitates precise and efficient application self-adhesive
- hypoallergenic

indication

• stabilization of catheters and stomach tubes fixed in the nasal area

| REF | dimensions | purpose | packaging unit box |
|--------|-----------------|--------------|-----------------------|
| 811014 | 7 cm x 7,1 cm | for adults | 50 x 1 pc. |
| 11015 | 5,9 cm x 6,2 cm | for children | 50 x 1 pc. |
| 311016 | 3,6 cm x 4 cm | for infants | 50 x 1 pc. |









and the skin takes place • elastic - adheres well, conforms to the body contours • does not absorb X-rays

- acrylic adhesive applied continuously over the entire surface of the plaster
- backing layer of siliconised paper
- two-step application system facilitating precise and effective application hypoallergenic

indication

• stabilization of catheters and stomach tubes fixed in the nasal area

| REF | dimensions | purpose | packaging unit box | |
|---------|-----------------|--------------|-----------------------|--------------------|
| 811015A | 5,9 cm x 6,2 cm | for children | 50 x 1 pc. | bulk pac 40 x 5 |













features





elastoNASALKIDS







<u>^</u> air pe str









| ermeable | blue-white | | |
|----------|------------|--|--|
| ructure | colour | | |
| | SYEARS | | |

| <u>↑</u> | | |
|-------------|------------|--|
| : neable | blue-white | |
| ure | colour | |
| | | |























↑STABILIZING↓TAPES





air permeable structure





- pads are covered
- the Velcro

adhesive tape for catheter tubing fixation

- two-part fixing tape consisting of a part attached to the skin and
- part attached to the patient's skin made of non-woven fabric
- drain fixing part is integrated with a part adhered to the patient's skin and has an additional adhesive tape - Velcro and adhesive field for better drain
- possibility of reposition and correction of the drain position thanks to Velcro
- stabilization of various types of drains and catheters

NEW

| packaging unit paper-foil | intermediate packaging box | |
|------------------------------|-------------------------------|--------------------------------|
| 1 рс. | 50 x 1 pc. | bulk packaging: 20 x 50 pcs |
| | | |

- part attached to the patient's skin made of adhesive colorful non-woven drain fixing part is integrated with a part adhered to the patient's skin and has an additional adhesive tape - Velcro and adhesive field for better drain
- acrylic adhesive applied continuously over the entire surface of the tape
- device intended for stabilisation and protection of all types of drains,

- set consisting of 2 self-adhesive non-woven pads and a fastening strap • self-adhesive non-woven pads attached to the patient's cheeks with an adhesive
- a strap for stabilization of the endotracheal or oropharyngeal tube made of non-woven fabric which should be attached to the velcro layer with which the
- adhesive area in the middle of the strip for better stabilization of the medical tube

| | packaging unit paper-foil | intermediate packaging box | |
|-------------|------------------------------|-------------------------------|---------------------------------|
| -woven pads | 1 set | 30x1 set | bulk packaging: 24 x 30 sets |

elastoFIXAL Set KIDS



features

- micropores in the non-woven fabric structure provide appropriate vapourpermeability rate, thanks to which proper gas exchange between the dressing and the skin takes place
- elastic adheres well, conforms to the body contours
- does not absorb X-rays
- acrylic adhesive applied continuously over the entire surface of the plaster backing layer of siliconised paper
- two-step application system facilitating precise and effective application
- hypoallergenic

indication

• device intended for stabilisation of catheters and stomach tubes fixed in the nasal area

| REF | kit composition | packaging unit paper-foil | intermediate packaging box | |
|---------|---|------------------------------|-------------------------------|--------------------|
| 816102A | 2 self-adhesive non-woven pads 1 fastening strap | 1 set | 30 x 1 set | bulk pac 24 x 3 |



NEW



NON CALEX AND GYEARS



band for tracheostomy tube fixation elastoFIXAL TUBE

features

- stabilising band consisting of two parts made of gentle hydrophobic nonwoven fabric in flesh-coloured
- on each part there are fastening strips with Velcro on one side enabling to attac it to the fastening wings of the tracheostomy tube
- shorter part of the stabilising band is additionally equipped with a second Velcro strap attached to an elastic white rubber band which makes it
- possible to adjust and determine the total length of the band stabilising the tracheostomy tube hypoallergenic
- single use

indication

device intended for stabilisation of tracheostomy tube

| REF | dimensions | packaging unit paper-foil | intermediate packaging box | bulk packaging box |
|---------|---------------|------------------------------|-------------------------------|-----------------------|
| EFT2520 | 25cm x 2,0 cm | 1 set | 30 x 1 set | 20 x 30 sets |
| EFT4325 | 43cm x 2,5cm | 1 set | 20 x 1 set | 20 x 20 sets |



conallhealth 欲 Your Partner In Healthcare



adhesive dressing pouch for dialysis catheters

features

- self-adhesive dressing with a pocket to protect dialysis catheters • the upper fixing part is made of self-adhesive, breathable, waterproof polyurethane film which is a barrier against contact with water and contamination
- breathable, non-sticky pocket with an outer layer resistant to contact with water and microorganisms provides adequate protection for the tips of dialysis catheters
- lower part stabilizing the pocket made of hydrophobic non-woven fabric

indication

• protective dressing is used to fix dialysis catheters quickly and safely dressing protect dialysis catheters against mechanical stress and contamination

| REF | dimensions | packaging unit paper-foil | intermediate packaging box | |
|--------|---------------|------------------------------|-------------------------------|-------------------------------|
| 811501 | 24 cm x 10 cm | 1 pc. | 25 x 1 pc. | bulk packaging: 4 x 25 pcs |



STABILIZING DRESSINGS



SUPPORTING BANDAGES

features

• delicate, light, air and X-ray permeable • conforms good to the body's contours hypoallergenic

indication

NON THERE AND ALL AND

fixation of all types of non-adhesive dressings



knitted supporting bandage elastoBAND BASIC

| features • contains viscose • single-use | | | | | | |
|--|-------------|-------------------------|-----------------------|--|--|--|
| REF | dimensions | packaging unit paper | bulk packaging carton | | | |
| OD05400 | 5 cm x 4 m | 1 pc. | 250 x 1 pc. | | | |
| OD10400 | 10 cm x 4 m | 1 pc. | 150 x 1 pc. | | | |
| OD15400 | 15 cm x 4 m | 1 pc. | 100 x 1 pc. | | | |



NONE TATES STATES

elastic supporting bandage, woven \bigcirc elastoBAND FLEX reusable

features

- supporting bandage, elastic, woven
 high cotton content makes the bandage safe for the skin and well-tolerated by people with allergic tendency
- equipped with a clip for fixing

indication

- fixation of all types of non-adhesive dressings, especially where a non-elastic bandages can slide down
- in ession therapy to facilitate correct blood circulation
- stabilization at bone-articular system damage, e.g. dislocation, joint twisting, weakness of muscles, tendons and ligaments
- as a sports bandage
 - securing and strengthening integument after surgery

| REF - cardboard box | REF - paper pouch | dimensions | packaging unit box/paper | bulk packaging carton |
|------------------------|----------------------|-------------|-----------------------------|-----------------------|
| | | | | |
| OE08400-E | OE08400P-E | 8 cm x 4 m | 1 pc. | 124 x 1 pc. |
| OE10400-E | OE10400P-E | 10 cm x 4 m | 1 pc. | 108 x 1 pc. |
| OE12400-E | OE12400P-E | 12 cm x 4 m | 1 pc. | 96 x 1 pc. |
| OE15400-E | OE15400P-E | 15 cm x 4 m | 1 pc. | 72 x 1 pc. |
| ОЕ08500-Е | OE08500P-E | 8 cm x 5 m | 1 pc. | 124 x 1 pc. |
| OE10500-E | OE10500P-E | 10 cm x 5 m | 1 pc. | 108 x 1 pc. |
| OE12500-E | OE12500P-E | 12 cm x 5 m | 1 pc. | 96 x 1 pc. |
| OE15500-E | OE15500P-E | 15 cm x 5 m | 1 pc. | 72 x 1 pc. |



NON DETERMINE STERILE STERILE STERILE

NON UNITED AUTO

elastic supporting bandage, woven elastoBAND FIX

features

- one-layer, woven elastic bandage breathable
- healing process

indication

after surgery

| REF | dimensions | packaging unit paper | intermediate packaging box | bulk packaging carton |
|----------|-------------|-------------------------|-------------------------------|--------------------------|
| BET06400 | 6 cm x 4 m | 1 pc. | 20 x 1 pc. | 54 x 20 pcs |
| BET08400 | 8 cm x 4 m | 1 pc. | 20 x 1 pc. | 45 x 20 pcs |
| BET10400 | 10 cm x 4 m | 1 pc. | 20 x 1 pc. | 36 x 20 pcs |
| BET12400 | 12 cm x 4 m | 1 pc. | 20 x 1 pc. | 27 x 20 pcs |
| BET15400 | 15 cm x 4 m | 1 pc. | 20 x 1 pc. | 27 x 20 pcs |

cohesive supporting bandage, non-woven elastoBAND CREPP

features

- hair, etc.
- breathable
- healing process

indication

after surgery

| REF | dimensions | packaging unit foil/box* | intermediate packaging box | bulk packaging carton |
|----------|----------------|-----------------------------|-------------------------------|--------------------------|
| OK04400 | 4 cm x 4 m | 1 mg | 10 × 1 mg | 29 × 10 mm |
| UK04400 | 4 CI I X 4 I I | 1 pc. | 10 x 1 pc. | 28 x 10 pcs |
| OK06400 | 6 cm x 4 m | 1 pc. | 10 x 1 pc. | 20 x 10 pcs |
| OK062000 | 6 cm x 20 m | 1 pc.* | - | 48 x 1 pc. |
| OK08400 | 8 cm x 4 m | 1 pc. | 10 x 1 pc. | 16 x 10 pcs |
| OK082000 | 8 cm x 20 m | 1 pc.* | - | 36 x 1 pc. |
| OK10400 | 10 cm x 4 m | 1 pc. | 10 x 1 pc. | 12 x 10 pcs |
| OK102000 | 10 cm x 20 m | 1 pc. * | - | 18 x 1 pc. |
| OK12400 | 12 cm x 4 m | 1 pc. | 10 x 1 pc. | 8 x 10 pcs |
| OK122000 | 12 cm x 20 m | 1 pc.* | - | 12 x 1 pc. |
| | | | | |



- bandage does not have therapeutic properties, it is a product supporting the
- packaged in foil-coated paper packaging from underneath

 device intended for holding dressings and for supporting and compressing body parts to treat sprains and strains, and to strengthen skin layers

one-layer, non-woven bandage
cohesive - adheres to itself (no need for fasteners), does not adhere to skin,

• bandage does not have therapeutic properties, it is a product supporting the

can be divided into smaller sections, no need to use scissors, the bandage
 can be torn in the hands

• packaged in foil package or cardboard box

device intended for holding dressings and for supporting and compressing body parts to treat sprains and strains, and to strengthen skin layers







elastoBAND KNITT

features

- one-layer, woven elastic bandage
- breathable
- bandage does not have therapeutic properties, it is a product supporting the healing process
- packaged in foil-coated paper packaging from underneath

indication

 device intended for holding dressings and for supporting and compressing body parts to treat sprains and strains, and to strengthen skin layers after surgery

| REF | dimensions | packaging unit paper | intermediate packaging box | bulk packaging carton |
|----------|-------------|-------------------------|-------------------------------|-----------------------|
| BED06400 | 6 cm x 4 m | 1 pc. | 20 x 1 pc. | 54 x 20 pcs |
| BED08400 | 8 cm x 4 m | 1 pc. | 20 x 1 pc. | 45 x 20 pcs |
| BED10400 | 10 cm x 4 m | 1 pc. | 20 x 1 pc. | 36 x 20 pcs |
| BED12400 | 12 cm x 4 m | 1 pc. | 20 x 1 pc. | 27 x 20 pcs |
| BED15400 | 15 cm x 4 m | 1 pc. | 20 x 1 pc. | 27 x 20 pcs |

cohesive bandage

yellowBAND

features

• 1-layer, breathable, made of non-woven • cohesive - it coheres to itself, does not adhere to skin, hair or other materials

indication

- fixation of all types of non-adhesive dressings,
- especially where the elastic band slides fixation of splinting rails
- in compression therapy

NON CATEX AND GYEARS

• reducing joints mobility after an injury

AVAILABLE SIZES AND COLOURS:

bulk packaging

| 2,5cm x 4,5m | | | | YCB-2545-X | | YCB-2545-B | | | | | 16 x 36 pcs |
|--------------|------------|------------|------------|------------|------------|------------|------------|-------------|-------------|-------------|-------------|
| 5cm x 4,5m | YCB-545-S | YCB-545-Y | YCB-545-G | YCB-545-X | YCB-545-R | YCB-545-B | YCB-545-W | YCB-545-BC | YCB-545-IO | YCB-545-IG | 24 x 12 pcs |
| 7,5cm x 4,5m | YCB-7545-S | YCB-7545-Y | YCB-7545-G | YCB-7545-X | YCB-7545-R | YCB-7545-B | YCB-7545-W | YCB-7545-BC | YCB-7545-IO | YCB-7545-IG | 16 x 12 pcs |
| 10cm x 4,5m | YCB-1045-S | YCB-1045-Y | YCB-1045-G | YCB-1045-X | YCB-1045-R | YCB-1045-B | YCB-1045-W | YCB-1045-BC | YCB-1045-IO | YCB-1045-IG | 12 x 12 pcs |
| 15cm x 4,5m | | | | YCB-1545-X | YCB-1545-R | YCB-1545-B | | | | | 8 x 12 pcs |

bulk packaging

| 2,5cm x 4,5m | YCB-2545-LB YCB-2545 | -IP YCB-2545-JG | YCB-2545-V | | | | | | 16 x 36 pcs |
|--------------|----------------------|-----------------|--------------|-------------|-------------|-------------|-------------|-------------|-------------|
| 5cm x 4,5m | YCB-545-LB YCB-545- | P YCB-545-JG | YCB-545-V | YCB-545-GC | YCB-545-GP | YCB-545-VP | YCB-545-XP | YCB-545-PH | 24 x 12 pcs |
| 7,5cm x 4,5m | YCB-7545-LB YCB-7545 | -IP YCB-7545-JC | 9 YCB-7545-V | YCB-7545-GC | YCB-7545-GP | YCB-7545-VP | YCB-7545-XP | YCB-7545-PH | 16 x 12 pcs |
| 10cm x 4,5m | YCB-1045-LB YCB-1045 | -IP YCB-1045-JO | G YCB-1045-V | | YCB-1045-GP | YCB-1045-VP | YCB-1045-XP | YCB-1045-PH | 12 x 12 pcs |
| 15cm x 4,5m | | | | | | | | | |



NON DATER AT STERILE

elastic tubular net bandage elastoNET

features

- elastic net bandage in the form of a tube
- contains polyamide and elastomer
- elastic and very soft
- easy and quick to use length before stretching:
- -1m, stretched: 2,5 m
- 10 m, stretched: 25 m

indication

it easier to fix the dressing to the patient's body

| REF | lenght before stretching | size | indication | packaging unit box - dispenser | bulk packaging carton |
|---------|-----------------------------|------|-----------------------------------|-----------------------------------|--------------------------|
| SO-0101 | 1 m | 1 | finger | 1 pc. | 64 x 1 pc. |
| SO-0110 | 10 m | 1 | finger | 1 pc. | 6 x 1 pc. |
| SO-0201 | 1 m | 2 | hand, toe | 1 pc. | 64 x 1 pc. |
| SO-0210 | 10m | 2 | hand, toe | 1 pc. | 6 x 1 pc. |
| SO-0301 | 1 m | 3 | hand, foot | 1 pc. | 64 x 1 pc. |
| SO-0310 | 10 m | 3 | hand, foot | 1 pc. | 6 x 1 pc. |
| SO-0401 | 1 m | 4 | lower leg, knee, arm, foot, elbow | 1 pc. | 64 x 1 pc. |
| SO-0410 | 10 m | 4 | lower leg, knee, arm, foot, elbow | 1 pc. | 6 x 1 pc. |
| SO-0601 | 1 m | 6 | head, shoulder, lower leg, knee | 1 pc. | 64 x 1 pc. |
| SO-0610 | 10 m | 6 | head, shoulder, lower leg, knee | 1 pc. | 6 x 1 pc. |
| SO-0801 | 1 m | 8 | thigh, head, hips | 1 pc. | 64 x 1 pc. |
| SO-0810 | 10 m | 8 | thigh, head, hips | 1 pc. | 6 x 1 pc. |
| SO-1001 | 1 m | 10 | hips, trunk | 1 pc. | 32 x 1 pc. |
| SO-1010 | 10 m | 10 | hips, trunk | 1 pc. | 4 x 1 pc. |
| SO-1401 | 1 m | 14 | chest, trunk | 1 pc. | 32 x 1 pc. |
| SO-1410 | 10 m | 14 | chest, trunk | 1 pc. | 4 x 1 pc. |

Hippocrates' achievements in treating head wounds were so widespread and advanced that they are still used today. The Hippocrates cap also called "mitra is a headgear that allows you to hold the dressing and stop bleeding of the head wound. It is made of two bandage tapes. The main purpose of using the Hippocrates cap is to hold the dressing directly protecting the bleeding head wound, seal and warm the head, hold the compress applied to the skull area, compression, strengthening and immobilization after surgery and additional protection against the external environment. This solution is used in the absence of access to dressings such as elastic dressings and scarves.

RULES FOR BANDAGING

How to make a

HIPOCRATES CAP?

- 1. Prepare a two-head bandage about 5-6 cm wide by sewing together two bands or by twisting the free end to the middle.
- 2. Put the two-headed bandage on the back of the head and drive bandage to the forehead a circular move.
- 3. Cross the bandages on the forehead and drive one through the top of the head towards the occiput by a longitudinal move.
- 4. Drive the other bandage above the ear to the occiput also (half-round). 5. Scroll the longitudinal course through the circular course bandage.
- 6. Drive back the longitudinal course bandage across the top of the head to the forehead.
- 7. Drive the circular course bandage around the other ear to the forehead
- 8. Continue longitudinal and circular moves, the circular course bandage wraps around the scalp
- 9. Little by little, move the longitudinal course bandage towards the temples 10. After covering the whole area with a two-head bandage,
- fix the ends of bandage

114 www.zarys.com

high stretchability and strength of returning to pre-stretch state

• fixation of all dressings on the patient's body, replacing the traditional bandage • the use of tubular net bandage reduces the time of wound dressing and makes

FIRST AID FOR DRESSING WOUNDS USING BANDAGES





ORTHOPEDIC PADDINGS

features

- very soft and fluffy
- adhere well and adjusts to the shape of the body
- coherent structure resistant to stratification makes it easier placement of the immobilizing dressing
- minimize the risk of skin abrasions and irritations under the P.O.P. bandage
- do not absorb moisture and thus do not allow for its contact with the skin
- hypoallergenic

indication

- as a protective layer put under plaster of Paris bandages
- in compression therapy as a padding under the bandages of various elasticity



white colour

triangular sling

elastoSLING

features • triangle shaped sling made of non-woven or cotton of the body, provides freedom of movement hypoallergenic indication or hand bone





orthopedic padding, natural



 \Diamond

white

colour

orthoNATURAL

| REF | dimensions | packaging unit paper cover | intermediate packaging foil bag | bulk packaging box |
|---------|-------------|-------------------------------|------------------------------------|-----------------------|
| *ON-63 | 6 cm x 3 m | 1 pc. | 12 x 1 pc. | 40 x 12 pcs |
| *ON-83 | 8 cm x 3 m | 1 pc. | 12 x 1 pc. | 40 x 12 pcs |
| ON-103 | 10 cm x 3 m | 1 pc. | 12 x 1 pc. | 20 x 12 pcs |
| ON-123 | 12 cm x 3 m | 1 pc. | 6 x 1 pc. | 40 x 6 pcs |
| ON-153 | 15 cm x 3 m | 1 pc. | 6 x 1 pc. | 30 x 6 pcs |
| *ON-203 | 20 cm x 3 m | 1 pc. | 6 x 1 pc. | 20 x 6 pcs |
| *ON-253 | 25 cm x 3 m | 1 pc. | 6 x 1 pc. | 10 x 6 pcs |
| | | | | |



- of the sling is at the level of the elbow of the immobilized hand. 2. Place the second corner of the sling on the opposite arm in relation to the immobilized hand. 3. Carefully place the sick / injured limb on the sling.
- 4. Wrap the sick/injured limb with dangling, third corner of the sling. 5. Make sure that the hand is in the right position, and then tie the corners of the sling around the neck.





orthopedic padding, synthetic

orthoSYNTHETIC

NONE TATEX STERILE

| REF | dimensions | packaging unit paper cover | intermediate packaging foil bag | bulk packaging box |
|--------|-------------|-------------------------------|------------------------------------|-----------------------|
| OS-63 | 6 cm x 3 m | 1 pc. | 12 x 1 pc. | 40 x 12 pcs |
| OS-83 | 8 cm x 3 m | 1 pc. | 12 x 1 pc. | 40 x 12 pcs |
| OS-103 | 10 cm x 3 m | 1 pc. | 12 x 1 pc. | 20 x 12 pcs |
| OS-123 | 12 cm x 3 m | 1 pc. | 6 x 1 pc | 40 x 6 pcs |
| OS-153 | 15 cm x 3 m | 1 pc. | 6 x 1 pc. | 30 x 6 pcs |
| OS-203 | 20 cm x 3 m | 1 pc. | 6 x 1 pc. | 20 x 6 pcs |
| OS-253 | 25 cm x 3 m | 1 pc. | 6 x 1 pc. | 10 x 6 pcs |





• thin, light and elastic - comfortable for the patient, adjusts to the shape

• possibility of connecting several slings to achieve bigger immobilization

• short-term immobilization and support of the upper limb in case of various injuries e.g. dislocation or fracture of collar bone, humerus, forearm

| | description | packaging unit zipper pouch | intermediate packaging foil bag |
|----------|-------------|--------------------------------|------------------------------------|
| | | | |
| x 136 cm | non-woven | 1 pc. | 12 x 1 pc. |
| x 130 cm | cotton | 1 pc. | 12 x 1 pc. |
| | | | |

HOW TO PROPERLY PUT ON A TRIANGULAR SLING?

- 1. Place the triangular sling on the patient's body so that one of the corners

bulk packaging: 20 x 12 pcs

P.O.P. BANDAGES

features

• made of cotton gauze coated with natural plaster of Paris on both sides • on a plastic roll

• enable fast and durable immobilization

 ${\scriptstyle \bullet}$ before fixing need to be soaked in a container with water at 20°C temperature

• do not absorb X-rays

indication

• large and small immobilizing dressings in orthopedic, ambulatory and surgical procedures

000




P.O.P. bandage



| REF | dimensions | type of edge | packaging unit foil bag | bulk packaging carton |
|-----------|-------------|--------------|----------------------------|--------------------------|
| | | | | |
| OB06300-F | 6 cm x 3 m | straight | 1 pc. | 48 x 1 pc. |
| OB08300-F | 8 cm x 3 m | straight | 1 pc. | 48 x 1 pc. |
| OB10300-F | 10 cm x 3 m | serrated | 1 pc. | 48 x 1 pc. |
| OB12300-F | 12 cm x 3 m | serrated | 1 pc. | 48 x 1 pc. |
| OB14300-F | 14 cm x 3 m | serrated | 1 pc. | 36 x 1 pc. |
| OB15300-F | 15 cm x 3 m | serrated | 1 pc. | 36 x 1 pc. |
| OB20300-F | 20 cm x 3 m | serrated | 1 pc. | 24 x 1 pc. |
| | | | | |

orthoBAND NORMAL





NONE TERES STEELE



orthoBAND LONG

| slow-setting slow- | | | | |
|--|-------------|--------------|----------------------------|-----------------------|
| REF | dimensions | type of edge | packaging unit foil bag | bulk packaging carton |
| OB10400-L | 10 cm x 4 m | serrated | 1 pc. | 48 x 1 pc. |
| OB12400-L | 12 cm x 4 m | serrated | 1 pc. | 36 x 1 pc. |
| OB15400-L | 15 cm x 4 m | serrated | 1 pc. | 24 x 1 pc. |





RULES FOR BANDAGING





Select a bandage of proper type and size for the part of the body to be bandaged.

With your right hand, grasp the beginning of the band, and with your left hand hold the rolled part so called "head".



For proper stabilization of dressing, start bandaging with a circular wrap, each subsequent layer covers exactly the previous one in a durable and aesthetic way. Little by little unroll the band during bandaging. Finish the bandaging with a circular wrap, end of the dressing secure with a clip or tie it.



The bandage should adhere correctly to the patient's body and be rolled up with adequate power - too much pressure can lead to impaired circulation, too loose it will slide and fall.







Put the patient in a comfortable position supporting the bandaged lower or upper limb.



Remove the bandage by transferring it from one hand to other to minimize patient movement.







paraffin gauze dressing

paraffiNET

features

- - secondary dressing
- and allows a painless change of dressing
- preparations containing medicinal substances or antiseptics
- can be shaped and cut according to the needs aseptically

- hypoallergenic

indication

- keeping proper moisture in the wound
- recommended for uninfected wounds
- treatment of chronic wounds: various ulcers, bedsores
- protection against radiation skin damage

| REF | dimensions | packaging unit paper-aluminum foil | intermediate packaging box | |
|--------|---------------|---------------------------------------|-------------------------------|--|
| PAR-01 | 5 cm x 5 cm | 1 pc. | 10 x 1 pc. | |
| PAR-02 | 10 cm x 10 cm | 1 pc. | 10 x 1 pc. | |
| PAR-06 | 10 cm x 40 cm | 1 pc. | 10 x 1 pc. | |
| PAR-07 | 15 cm x 20 cm | 1 pc. | 10 x 1 pc. | |

FOIL DRESSINGS

FOIL DRESSINGS

features

- made of polyurethane film
- covered with hypoallergenic acrylic glue
- transparent they allow observation of the dressed place without the need to remove the dressing
- vapour permeability guarantees correct gas exchange between the dressing and the skin, and prevents it from
- excessive moisture accumulation and maceration
- barrier against microorganisms and bacteria, minimizing the risk of infection
- very thin and elastic adjust to the shape of the body, guarantee freedom of movement
- waterproof
- self-adhesive hypoallergenic

elastoDERM



SPECIALIST DRESSINGS





- wide-mesh dressing made of cotton gauze covered with paraffin
- structure of the carrier material allows exudate to be freely drained to the covering
- covering with a layer of paraffin reduces risk of adhesion to the wound
- neutral properties of paraffin allow for using the dressing in combination with other
- needs to be covered with a dressing with absorbent properties, e.g. compress,
- high-absorption dressing, polyurethane foam
- non-adhesive, necessary fixing with a tape, bandage or net tubular dressing
- preventing the wound edges from dying and maceration
- treatment of surgical and burn wounds, skin donor and grafting sites, skin damage, cuts



dressing without absorbent pad, with "three-stage" application type

| ре | indication | packaging unit paper-paper | intermediate packaging box |
|----|---|-------------------------------|-------------------------------|
| | fixing and securing intravenous cannulas secondary dressing bedsore prevention, works as so-called, second skin" for small epidermal abrasions and damages of the skin surface | 1 pc. | 100 x 1 pc. |





features

- "frame" or three-stage" application system makes dressing placement precise and accurate
- available in several sizes and options: - with a central absorbent pad
- with notch enabling the dressing to fit the cannulas
- without notch and absorbent pad

elastoDERM F

dressing without absorbent pad, with "frame" application type

1. Remove the back paper protection.



2. Apply the dressing to the desired place.



3. Remove the paper protection on the outer side of the dressing starting from the cut.



4. Put the notes on the tape to describe the dressing and apply it on the dressing without covering the wound/ the place of needle insertion.

| REF | dimensions | shape | indication | packaging unit paper-paper | intermediate packaging box | bulk packaging box |
|-------------------------------|-----------------|------------|--|-------------------------------|-------------------------------|-----------------------|
| 812003 | 4,4 cm x 4,4 cm | | | 1 pc. | 100 x 1 pc. | 36 x 100 pcs |
| | | | fixing and securing intravenous cannulas secondary dressing bedges argumented uselo | | | |
| 812005 | 6 cm x 7 cm | | bedsore prevention, works as so-called "second skin" | 1 pc. | 100 x 1 pc. | 30 x 100 pcs |
| with special to write note | | | for small epidermal abrasions and damages of the skin surface | | | |
| 812008 | 10 cm x 12 cm | - | - fixing and securing central | 1 pc. | 100 x 1 pc. | 20 x 100 pcs |
| with special to write note | | | venous catheters - fixing epidural catheters - secondary dressing | | | |
| 812011 | 10 cm x 15 cm | | bedsore prevention, works as so-called "second skin" | 1 pc. | 100 x 1 pc. | 20 x 100 pcs |
| with special to write note | | | for small epidermal abrasions and damages of the skin surface | | | |
| STERILE | EO 😿 🤇 | SHELF LIFE | | | | |

elastoDERM F-IV

IV cannula dressing, with "frame" application type packaging unit intermediate RFF dimensions shape indication paper-paper packaging box 812006 6 cm x 7 cm 100 x 1 pc. 1 pc. fixing and securing peripheral cannulas in adults securing cannulas inserted into central veins with special area to write notes Ø bulk packaging: 30 x 100 pcs

elastoDERM IV

IV cannula dressing, with "three-stage" application type

| REF | dimensions | shape | indication | packaging unit paper-paper | intermediate packaging box |
|-------------------------------------|---------------|---|----------------------------------|-------------------------------|---------------------------------|
| 811301 | 5,8 cm x 8 cm | | - fixing and securing peripheral | | |
| with special area to write notes | | 2 | cannulas in adults | 1 pc. | 100 x 1 pc. |
| STERILE | EO 😿 🛇 | SHELF UPE | | | bulk packaging: 20 x 100 pcs |

elastoDERM PAD

dressing with central absorbent pad, with "three-stage" application type

| REF | dimensions | shape | indication | packaging unit paper-paper | intermediate packaging box | bulk packaging box |
|--------|---------------|-------|--|-------------------------------|-------------------------------|-----------------------|
| 802012 | 6 cm x 7 cm | | | 1 pc. | 100 x 1 pc. | 12 x 100 pcs |
| 802014 | 9 cm x 10 cm | | | 1 pc. | 50 x 1 pc. | 24 x 50 pcs |
| 802015 | 9 cm x 15 cm | ** | protection of post-operative wounds and stitches | 1 pc. | 50 x 1 pc. | 24 x 50 pcs |
| 802016 | 10 cm x 20 cm | | wounds and stitches | 1 pc. | 50 x 1 pc. | 18 x 50 pcs |
| 802018 | 10 cm x 25 cm | | | 1 pc. | 50 x 1 pc. | 15 x 50 pcs |





NONE TATES ALL STREET

cold/warm compress



features

- ${\boldsymbol{\cdot}}$ cold/warm gel compress made of non-toxic durable film with gel insert
- presence of a fabric cover to prevent external damage
 ability to store and return accumulated heat or cold for a long time
- remains flexible after freezing, allowing the compress to be easily shaped to the body area

reusable

 indication ${\boldsymbol{\cdot}}$ device used for the local cooling or heating of specific areas of the body to combat pain/neuralgia and reduce swelling and puffiness

| REF | dimensions | packaging unit foil bag | bulk packaging carton | |
|---------|----------------|----------------------------|-----------------------|--|
| KCZ7513 | 7,5 cm x 13 cm | 1 pc. | 100 x 1 pc. | |
| KCZ7552 | 7,5 cm x 52 cm | 1 pc. | 30 x 1 pc. | |
| KCZ1229 | 12 cm x 29 cm | 1 pc. | 25 x 1 pc. | |
| KCZ1314 | 13 cm x 14 cm | 1 pc. | 50 x 1 pc. | |
| KCZ1626 | 16 cm x 26 cm | 1 pc. | 20 x 1 pc. | |
| KCZ2138 | 21 cm x 38 cm | 1 pc. | 10 x 1 pc. | |
| KCZ3040 | 30 cm x 40 cm | 1 pc. | 8 x 1 pc. | |
| | | | | |



NONE TATEX STATES

gynaecological-maternity pad gynPAD

features

an absorbent pad with an adhesive layer

NEW

- it has a layered construction
- pink foil sheet placed in the lower part of the pad as an insulating part to protect against wetness
- provides increased absorption capacity due to the highly absorbent cellulose pad
- structure of the dressing allows free passage of exudate to the absorbent pad inside
- can be used as a primary as well as a secondary dressing • can be sterilised with ethylene oxide

indication

• device used to absorb bodily fluids resulting from surgery, childbirth to protect the patient from contamination and infection

| REF | dimensions | packaging unit foil bag | intermediate packaging box | { |
|---------|--------------|----------------------------|-------------------------------|----------------|
| PG34090 | 34 cm x 9 cm | 1 pc. | 25 x 1 pc. | bulk p 50 x |



gauze pad moistened with alcohol **Alcohol Prep Pad**

| features made of non-woven soaked with isopropyl alcohol (70%) two pouches connected in one-piece packaging can be separated thanks to perforation indication disinfection of small surfaces | | | | | | |
|---|--|-----------------|---------------------------------------|---------------------------|--|--|
| REF | dimensions | no. of plies | packaging unit paper-aluminum foil | bulk packaging box | | |
| P1-IA-2 | 30 mm x 65 mm (unfolded) 30 mm x 32,5 mm (folded) | 2 | 2 x 1 pc. | 100 x 1 pc. 50 x 2 pcs | | |





P2-IA-1 M 65 mm x 55 mm (ur 32,5 mm x 27,5 mm P2-IA-2 XL 110 mm x 90 mm (u

- to perforation (without L-06)

indication

disinfection low surface

REF dimensions

L-06 XXL 190 mm x 140 mm



ALKOPAD DUO

features

- made of non-woven

indication

disinfection low surface

| REF | dimensions | packaging unit paper-aluminum foil | bulk packaging box | |
|---------|-----------------|---------------------------------------|-----------------------|--|
| L-05DUO | 125 mm x 125 mm | 2 x 1 pc. | 50 x 2 pcs | |

• made of metallized PE foil

indication

| REF | dimensions | purpose | packaging unit zipper pouch | intermediate packaging foil bag | |
|----------|-----------------|------------|--------------------------------|------------------------------------|--------------------------------|
| KR160210 | 160 cm x 210 cm | for adults | 1 pc. | 12 х 1 рс. | bulk packaging: 20 x 12 pcs |

features

- made of non-woven

indication

• disinfection low surface



55 mm x 45 mm (fold

ALKOPAD

features

made of non-woven


• soaked with isopropyl alcohol (70%)

conallhealth 🕸





NON CALEX STATES



rescue blanket, foil, termic thermCARE

features

• resistant to damage













• each gauze pad is packed separately, in aluminum foil - paper soaked with isopropyl alcohol (70%) • two pouches connected in one-piece packaging can be separated thanks to perforation (without L-06)

| | packaging unit paper-aluminum foil | bulk packaging box | |
|------------|---------------------------------------|---------------------------|--|
| 6 I I II | | | |
| nfolded) | 2 x 1 pc. | 100 x 1 pc. (50 x 2 pcs.) | |
| n (folded) | | | |
| unfolded) | 2 x 1 pc. | 100 x 1 pc. (50 x 2 pcs.) | |
| lded) | | | |

gauze pad moistened with alcohol

• each gauze pad is packed separately, in aluminum foil - paper

• two pouches connected in one-piece packaging can be separated thanks

| | packaging unit paper-aluminum foil | bulk packaging box | |
|------------|---------------------------------------|--------------------------|--|
| | | | |
| (unfolded) | 2 x 1 pc. | 50 x 1 pc. (25 x 2 pcs.) | |
| | | | |

• each gauze pad is packed separately, in aluminum foil - paper • gauze pad wet soaked with isopropyl alcohol (70%) and gauze pad dry • two pouches connected in one-piece packaging can be separated thanks to perforation

protects against adverse weather conditions – cold, rain, heat, wind, depending on the side, which is applied towards the victim's body, prevents the life-threatening heat disorders and reduces the risk of shock

protection against hypothermia and overheating of the victim's body
in emergency cases such as car accidents, ski accidents, in mountain tourism, in sport


surgical kits

The company's offer also includes small procedure packs under the deltaset brand. They are used in treatment rooms, emergencies, hospital wards or dialysis units. They contain all disposable products necessary to perform the procedure in aseptic manner.

Deltaset packages are dedicated to a specific procedure, e.g. for changing a dressing, catheterization, dialysis or anaesthesia. The set includes: disposable equipment, dressings, surgical drapes, surgical instruments, accessories, medical gloves and many more.










🛆 deltaset standard surgical kits

deltaset label











2-chamber blister



3-chamber blister









ZARYS

WWW.ZARYS.COM

DS-WKL-O-01 PUNCTURE KIT1



features

- single-use
- for short-term use
- · shelf life: 3 years
- label with 2 adhesive tag labels for data archiving;
- the following information is available on each adhesive tag label: REF number, expiry date, batch number, manufacturer's data and barcode
- packaging: double-chamber hard blister indication

- intended for use in treatment rooms, emergency services or hospital wards
- it contains auxiliary disposable medical devices necessary to perform the central puncture procedure

| 10. | quantity | component | specification |
|--------------|----------|---|---|
| L. | 1рс. | 2-layer surgical drape | 90 cm x 75 cm cellulose layer + waterproof polyethylene film; blue colour |
| 2. | 1рс. | 2-layer surgical drape with central adhesive fenestration | 45 cm x 75 cm fenestration dia. 8 cm, cellulose layer + waterproof polyethylene film; blue colour |
| 3. | 10 pcs | gauze ball without X-ray thread | 20 cm x 20 cm 17-thread cotton gauze; white colour |
| ł. | 10 pcs | gauze swab | 7,5 cm x 7,5 cm 17-thread cotton gauze, 8-ply; white colour |
| 5. | 1 pc. | plastic tweezers | length 12,5 cm, polypropylene; blue colour |
| 5. | 1 pc. | metal grasper (Pean) | length 13 cm, stainless steel |
| 7. | 1рс. | metal needle holder MAYO-HEGAR | length 13 cm, stainless steel |
| 3. | 1рс. | two-piece syringe in paper-foil packaging | capacity 10 ml, Luer |
|). | 1рс. | two-piece syringe in paper-foil packaging | capacity 20 ml, Luer |
| 10. | 1рс. | injection needle in paper-foil packaging | 1,2 mm x 40 mm |
| L 1 . | 1рс. | injection needle in paper-foil packaging | 0,8 mm x 40 mm |
| ۱2. | 1рс. | surgical blade | no. 11 |
| L3. | 1рс. | foil dressing with frame and tape for description in paper-foil packaging | 10 cm x 15 cm, polyethylene film |


4

9.

ų

STANDARD DELTASET **SURGICAL KITS**





individual packaging

1 kit double-chamber hard blister



DS-USZ-O-01



features

- single-use
- for short-term use
- shelf life: 3 years
- label with 2 adhesive tag labels for data archiving; the following information is available on each adhesive tag label:
- REF number, expiry date, batch number, manufacturer's data and barcode • packaging: single-chamber hard blister

indication

- intended for use in treatment rooms, emergency services, hospital wards or dialysis stations
- it contains auxiliary disposable medical devices necessary
- to perform the suture removal procedure

surgical blade

plastic tweezers

anatomic metal tweezers

1.

4.4 cm x 1

individual packaging

1 kit single-chamber

hard blister

50x60 cm x 1

5

20x20 cm x 6 7,5x7,5 cm x 8 12,5 cm x 1

4

- 12.5 cm x 1
- 12 cm x 1



50x60 cm x 1

ø5cm

7


 \cap

14 cm x 1

Mx2

120 ml x 1

- 20x20 cm x 3



- the following information is available on each adhesive tag label: REF number, expiry date, batch number, manufacturer's data and barcode
- packaging: single-chamber hard blister

indication

features

single-use

for short-term use

shelf life: 3 years

- intended for use in treatment rooms, emergency services or hospital wards
- it contains auxiliary disposable medical devices necessary to perform the urinary bladder catheterisation procedure

| no. | quantity | component | spe |
|-----|----------|--|--------------------|
| 1. | 1рс. | 2-layer surgical drape with central, fenestration and incision | 75c nor poly |
| 2. | 1рс. | 2-layer surgical drape | 45c poly |
| 3. | 2 pcs | examination gloves, nitrile, powderfree, in paper packaging | size |
| 4. | 4 pcs | gauze ball without X-ray thread | 20c |
| 5. | 5 pcs | gauze swab | 7,50 |
| 6. | 1 pc. | two-piece syringe in paper-foil packaging | сар |
| 7. | 1 pc. | injection needle in paper-foil packaging | 1,2n |
| 8. | 1 pc. | plastic tweezers | lenç |
| 9. | 1 pc. | plastic grasper (Pean) | len |

DS-DIA-O-01A DIALYSIS ON KIT 1A



features

- single-use
- for short-term use
- shelf life: 3 years
- label with 2 adhesive tag labels for data archiving;
- the following information is available on each adhesive tag label
- REF number, expiry date, batch number, manufacturer's data and barcode
- packaging: soft blister paper-foil

indication

- intended for use in treatment rooms, emergency services,
- hospital wards or dialysis stations
- it contains auxiliary disposable medical devices necessary to perform the dialysis procedure
- no. quantity component specification 2-layer surgical drape 50cm x 35cm pc. cellulose layer + waterproof polyethylene film; blue colour size. M, 100% nitrile; light blue colour examination gloves, nitrile, 2 pcs powderfree in paper packaging IV cannula dressing, self-adhesive 6cm x 8cm, hydrophobic non-woven 2 pcs 7,5cm x 7,5cm, non-woven 30G, 4 ply; white colour pcs non-woven swab 2 pcs dressing tape, self-adhesive 2,5cm x 15cm

DS-CEW-O-01

no. quantity component

1 pc.

3 pcs

1 pc.

l nc

URINARY BLADDER CATHETERIZATION KIT1

specification

length 12,5cm, polypropylene, blue colour

length 12cm, stainless steel

no. 11

non-woven ball without X-ray thread 20 cm x 20 cm, non-woven 30 G, white colour



features

- single-use
- for short-term use
- shelf life: 3 years
- label with 2 adhesive tag labels for data archiving; the following information is available on each adhesive tag label:
- REF number, expiry date, batch number, manufacturer's data and barcode • packaging: single-chamber hard blister

indication


- intended for use in treatment rooms, emergency services or hospital wards
- it contains auxiliary disposable medical devices necessary to perform the
- urinary bladder catheterisation procedure


| no. | quantity | component | specification | individual packaging |
|-----|----------|---|---|---|
| 1. | 1рс. | 2-layer surgical drape | 50cm x 60cm non-woven polypropylene + waterproof, polyethylene film; blue colour | |
| 2. | 1рс. | 2-layer surgical drape with central, fenestration and incision | 50cm x 60cm fenestration: diameter 5cm, non-woven polypropylene + waterproof, polyethylene film; blue colour | |
| 3. | 2 pcs | examination gloves, nitrile, powderfree, in paper packaging | size. M, 100% nitrile; light blue colour | 1 kit single-chamber hard blister |
| 4. | 6 pcs | gauze ball without X-ray thread | 20cm x 20cm, 17-thread cotton gauze; white colour | |
| 5. | 8 pcs | gauze swab | 7,5cm x 7,5cm, 17-thread cotton gauze, 8-ply; white colour | |
| 6. | 1рс. | plastic tweezers | length 12,5cm, polypropylene; green colour | |
| 7. | 1рс. | plastic grasper (Pean) | length 14cm, polypropylene; blue colour | |
| 8. | 1рс. | sponge bowl | capacity 120ml, polypropylene; transparent colour | |





· label with 2 adhesive tag labels for data archiving;






1 kit soft blister paper-foil DS-DIA-O-01B DIALYSIS OFF KIT 1B



features

- single-use
- for short-term use
- shelf life: 3 years
- · label with 2 adhesive tag labels for data archiving;
- the following information is available on each adhesive tag label:
- REF number, expiry date, batch number, manufacturer's data and barcode
- packaging: soft blister paper-foil
- indication
- intended for use in treatment rooms, emergency services, hospital wards or dialysis stations
- it contains auxiliary disposable medical devices necessary to perform the dialysis procedure



| no. | quantity | component | specification | individual packaging |
|-----|----------|---|---|----------------------------|
| 1. | 2 pcs | examination gloves, nitrile, powderfree in paper packaging | size. M 100% nitrile; light blue colour | _ |
| 2. | 1рс. | examination gloves, nitrile | size. L 100% nitrile; white colour | 1 kit |
| 3. | 2 pcs | post-injection dressing, self-adhesive | 3cm x 7,2cm polyethylene film; transparent | soft blister paper-foil |
| 4. | 6 pcs | non-woven swab | 7,5cm x 7,5cm non-woven 30G, 4 ply; colour white | - |

1

50x35 cm x 1

2.

6.

nIII/n

Mx2

nM/

Mx2

3



features

- single-use
- for short-term use
- shelf life: 3 years
- label with 2 adhesive tag labels for data archiving; the following information is available on each adhesive tag label: REF number, expiry date, batch number, manufacturer's data and barcode
- packaging: two-piece soft blister with perforation, paper-foil

indication

- intended for use in treatment rooms, emergency services or hospital wards
- it contains auxiliary disposable medical devices necessary to perform the
- puncture procedure e.g. central venepuncture, subarachnoid puncture

| no. | quantity | component | specification | individual packaging |
|-----|----------|---|---|---|
| 1. | 1рс. | 2-layer surgical drape | 50cm x 35cm; cellulose layer + waterproof polyethylene film; blue colour | |
| 2. | 2 pcs | examination gloves, nitrile, powderfree in paper packaging | | |
| 3. | 2 pcs | IV cannula dressing, self-adhesive | 6cm x 8cm, hydrophobic non-woven | |
| 4. | 6 pcs | non-woven swab | 7,5cm x 7,5cm, non-woven 30G, 4 ply; white colour | 1 kit two-piece soft |
| 5 | 4 pcs | dressing tape, self-adhesive | 2,5cm x 15cm | blister with perforation, paper-foil |
| 6. | 2 pcs | examination gloves, nitrile, powderfree in paper packaging | size. M, 100% nitrile; light blue colour | |
| 7. | 1рс. | examination gloves, nitrile | size. L, 100% nitrile; white colour | |
| 8. | 2 pcs | dressing with absorbent pad, selfadhesive | 5cm x 7,2cm | |
| 9. | 4 pcs | non-woven swab | 7,5cm x 7,5cm, non-woven 30G, 4 ply; colour white | |

DS-ZNI-O-01 ANASTHESIA KIT 1



features

- single-use
- for short-term use
- shelf life: 3 years
- label with 2 adhesive tag labels for data archiving;
- the following information is available on each adhesive tag label:
- REF number, expiry date, batch number, manufacturer's data and barcode • packaging: three-chamber hard blister

indication

- intended for use in treatment rooms, emergency
- services or hospital wards
- it contains auxiliary disposable medical devices necessary to perform procedures of conduction anesthesia

| no. | quantity | component | specification |
|-----|----------|---|---|
| 1. | 1 рс. | 2-layer surgical drape with central fenestration and adhesive corners | 75cm x 90cm; fenestra polypropylene + water |
| 2. | 10 pcs | gauze swab | 7,5cm x 7,5cm, 17-three |
| 3. | 1pc. | plastic window forceps, Korcang type, wrapped in a 2-layer drape with a cellulose layer | length 19 cm, 45cm x 4 polyethylene film, blue |

DS-SZY-O-01 SUTURE APPLICATION KIT1


features

- single-use
- for short-term use
- shelf life: 3 years
- label with 2 adhesive tag labels for data archiving;
- the following information is available on each adhesive tag label:
- REF number, expiry date, batch number, manufacturer's data and barcode
- packaging: double-chamber hard blister


indication

- intended for use in treatment rooms, emergency services or hospital wards
- it contains auxiliary disposable medical devices necessary to perform the
- puncture procedure e.g. central venepuncture, subarachnoid puncture

| no. | quantity | component | specification |
|-----|----------|---|---|
| 1. | 1рс. | 2-layer surgical drape with central adhesive fenestration | 50cm x 60cm, fenestration: dia non-woven polypropylene + wat polyethylene film; blue colour |
| 2. | 1рс. | 2-layer surgical drape | 45cm x 75cm, non-woven polyp polyethylene film; blue colour |
| 3. | 3 pcs | gauze ball without X-ray thread | 20cm x 20cm 17-thread cotton gauze; white c |
| 4. | 5 pcs | non-woven swab | 7,5cm x 7,5cm non-woven 30G, 4 ply; white col |
| 5. | 1рс. | anatomic metal tweezers | length 12cm stainless steel |
| 6. | 1рс. | plastic tweezers | length 12,5cm polypropylene; green colour |
| 7. | 1рс. | metal needle holder MAYO-HEGAR | length 13 cm stainless steel |
| 8. | 1рс. | metal scissors sharp/sharp | length 11,4cm stainless steel |



ackaging







O omegapack

surgical sets

Surgical sets are medical devices personalized according to customer needs, dedicated to a selected medical procedure.



PROCEDURE SETS

Operating theatre is a specific environment where time is the most valuable factor influencing the success of the operation. It is a place where every minute - or even a second - counts, because it can affect the patient's life. To meet the expectations of medicine, ZARYS introduces operating sets to its offer that meet the highest standards, helping to improve the preparation of the operating theatre before the procedure

ADVANTAGES of using omegapack operating sets

 \odot



30% less time than selecting individual products and improving work at operating bLock

> REDUCING COST





OUALITY

STANDARDIZATION OF PROCEDURES

properly selected composition of the set according to individual needs for a specific medical procedure in line with the Customer's order

reduced risk of patient infection by minimizing the time of contact of medical staff with devices during the preparation of the theatre and the patient

correct and schematic conduct during the procedure, one package - one treatment

EPIDEMIOLOGICAL SAFETY

EFFECTIVENESS

I

increasing work productivity, so that you can perform more surgical procedures in one day, improvement of work of medical staff



Ş optimization of record keeping - one self-adhesive label with the identification of the used set, instead of each used product separately

REDUCTION OF ERRORS



WIDE RANGE OF PRODUCTS 000

Å

SAVING

THE STOCKS

reduced number of supplies, saving space in the warehouse

a wide range of medical devices selected individually in accordance with the Customer's preferences

The sets are assembled on the basis of individual



Medical devices under the omegapack brand by ZARYS are operating sets containing complete equipment with all disposable products necessary to perform a given medical procedure to which the set is dedicated. The operating set contains carefully selected disposable components, such as:

• surgical drapes, pads

- accessories, e.g. fluid collection
- bags, tapes
- surgical gowns
- surgical gloves
- covers for equipment
- shoe covers, masks, caps

gauze balls

- gauze and non-wover
- bandages, dressings
- catheters, urine colle
- trays, bowls
- taps, three-way stopc
- manifolds, connectors

necessary for the surgical operation, which eliminates the possibility of errors during the selection of individual products

the complete set contains all disposable medical devices

F

requirements and preferences of Customers

| | • infusion sets |
|-----------|---|
| swabs | suction sets |
| | disposable surgical instruments |
| tion bags | injection needles and syringes |
| | disposable electrosurpicel electrodes |
| ocks | and accessories |
| | and many more |



safe

personal protective

equipment



PROTECTIVE COVERALL SAFE+ TYPE 4b - 5b - 6b

Personal protective equipments (PPE) are products that protect against one or more risks to the user's health or safety. These products must comply with Regulation (EU) 2016/425 of the European Parliament and of the Council. PPE is classified into three categories according to the degree of risk.

The coverall is a PPE category III in accordance with the Regulation (EU) 2016/425 of the European Parliament and of the Council



The category III covers products of complex construction, intended to protect health or life against risks that may cause death or irreversible damage to health.





Coverall is personal protective equipment of category III in conformity with Regulation (EU) 2016/425 of the European Parliament and of the Council

EN ISO 13688:2013 Protective clothing. General requirements

EN 14605:2005+A1:2009 Protective clothing against liquid chemicals. Performance requirements for clothing with liquid-tight (Type 3) or spray-tight (Type 4) connections, including items providing protection to parts of the body only (Types PB [3] and PB [4])

EN ISO 13982-1:2004+A1:2010 Protective clothing for use against solid particulates. Part 1: Performance requirements for chemical protective clothing providing protection to the full body against airborne solid particulates (type 5 clothing)



EN ISO 13982-1:2004 +A1:2010

EN 13034:2005+A1:2009 Protective clothing against liquid chemicals. Performance requirements for chemical protective clothing offering limited protective performance against liquid chemicals (Type 6 and Type PB [6] equipment)



EN 14126:2003 + AC:2004

EN 1073-2:2002 Protective clothing against radioactive contamination. Part 2: Requirements and test methods for non-ventilated protective clothing against particulate radioactive contamination

EN 14126:2003 + AC:2004 Protective clothing. Performance requirements and tests methods for protective clothing against infective agents

class 6 Resistance to penetration of contaminated liquids under the influence of hydrostatic pressure class 6 Resistance to penetration by infective agents due to mechanical contact with substances class 3 Resistance to penetration of contaminated liquid aerosols class 3 Resistance to penetration of contaminated solid particles

<u>4</u> EN 1149-5:2018

EN 1149-5:2018 Protective clothing. Electrostatic properties. Part 5: Material performance and design requirements

EN 14325:2004 Protective clothing against chemicals. Test methods and performance classification of chemical protective clothing materials, seams, joins and assemblages







AGAINST CHEMICAL AND INFECTIVE AGENTS

features:

- made of a laminate (non-woven polypropylene and polyethylene) with a weight of 63 g/m²
- polyester seams additionally covered with tape
- equipped with a hood with elastic trim
- equipped with an elastic waist band, wrists and ankles
- zipper covered with a self-adhesive flap
- available in version with boots



| | SIZE ROZMIAR | | | | | | |
|---|-----------------|-------|--------|---------|---------|---------|---------|
| | cm | 86-94 | 94-102 | 102-110 | 110-118 | 118-129 | 129-141 |
| - | 158-166 | S | | | | | |
| T | 166-174 | | | | | | |
| | 174-182 | | | . L | | | |
| | 182-190 | | | | XL | | |
| | 190-198 | | | | | XXL | |
| ै | 198-206 | | | | | | XXXL |

| REF | size | version | packaging unit |
|-----------------|------|---------------|-------------------|
| BS-KCB02-S | S | without boots | 1 pc. |
| BS-KCB02-M | Μ | without boots | 1 pc. |
| BS-KCB02-L | L | without boots | 1 pc. |
| BS-KCB02-XL | XL | without boots | 1 pc. |
| BS-KCB02-XXL | XXL | without boots | 1 pc. |
| BS-KCB02-XXXL | XXXL | without boots | 1 pc. |
| BS-KCB02-B-S | S | with boots | 1 pc. |
| BS-KCB02-B-M | Μ | with boots | 1 pc. |
| BS-KCB02-B-L | L | with boots | 1 pc. |
| BS-KCB02-B-XL | XL | with boots | 1 pc. |
| BS-KCB02-B-XXL | XXL | with boots | 1 pc. |
| BS-KCB02-B-XXXL | XXXL | with boots | 1 pc. |
| | | | |

Preventing infections in the operating field is a priority for all medicals. Thanks to our complex offer of singleuse protective solutions of ALPHAtex STANDARD, STANDARD+, COMFORT, the post-procedural complications and risk of infection in the operating field are minimized while increasing the protection barrier against the potentially pathogenic germs residing in the human skin layers. We offer:

- patient and equipment drapes
- sets of surgical drapes
- incise films
- surgical covers
- adhesive tapes
- accessories such as fluid collection pouches
- surgical gowns

Properly selected features of surgical drapes provide safety and comfort of use to medical staff and patients. A variety of multi-layer non-woven and laminates reinforced in critical areas in lines STANDARD+ and COMFORT allow to adjust our solutions to particular surgical procedures with maximal cost efficiency.

ALPHAtex products are characterized by parameters frequently exceeding the requirements of EN 13795 for 'Surgical drapes, gowns and clean air suits, used as medical devices for patients, clinical staff and equipment'.

Moreover, they provide a successful protection against the risk of infection in the operating field. Practical solutions and composition of our sets enable their quick and simple application while performing medical procedures.

C ALPHAtex surgical sets, drapes, gowns

surgical drapes I surgical covers and accessories I surgical sets I surgical gowns



Prevention of surgical site infection



The wide range of ALPHAtex products is a result of over 30 years of ZARYS International Group existence on medical market. Sharing our experience and knowledge we are convinced that using ALPHAtex products will provide comfort, safety and economic effectiveness during daily work in the operating room.

Surgical site infection (SSI) -

is defined as infection occurring up to 30 days after surgery, when there is a possibility of microorganisms penetration into the open surgical wound. Microorganisms could have various sources EXOGENOUS (carried by an external factor as medical staff or medical devices) or ENDOGENOUS (patient's bacterial flora from the closest area surrounding the surgical postoperative wound). The most important routes of infections are airway or contact path. Various type of barriers are used to minimalize the surgical site infections: surgical drapes, surgical gowns and medical equipment covers.



Surgical drapes (the names are based on the EN 13795 standard)

Surgical drapes are used to ensure clean microbiological surgical wound surrounding, minimalize the microorganisms migration from patient skin to surgical wound and to supervise the spreading of body fluids from the wound area. Surgical drapes meet the requirements of medical device in accordance with the Directive 93/42/EEC.

Surgical gowns (the names are based on the EN 13795 standard)

Surgical gowns are used to prevent direct microorganisms transmission from medical staff to patient. Moreover, surgical gowns prevent spreading of exfoliated epidermis in the air of the operating room.






ZONES OF SURGICAL DRAPES ON PATIENT

on the operating room

ZONE 4 NOT HAVING ANY CONTACT WITH PATIENT'S BODY performs the function of anaesthetic cover and operating table cover

ZONE 3

PARTIAL CONTACT WITH NOT OPERATED PARTS OF PATIENT'S BODY occasionally exposed to body fluids



disinfected patient's skin in direct surroundings of surgical wound, mostly protected by incise film



Materials legend:

detailed information and material structures available in **material template ALPHAtex**



STAN CHARACTERISTICS STANDARD UNIT CRIT PRODU microbial penetration EN ISO 22612 CFU not requ - dry microbial penetration EN ISO 22610 ≥2,8 I_ - wet cleanliness microbial EN ISO 11737-1 CFU/100 cm² ≤ 300 bioburden log₁₀ (lint count) particulate release EN ISO 9073-10 ≤ 4,0 liquid penetration EN ISO 811 $cm H_2O$ ≥30 bursting strength EN ISO 13938-1 kPa ≥40 - dry bursting streangth EN ISO 13938-1 kPa ≥40 - wet tensile strength ≥15 EN 29073-3 Ν - dry tensile strength EN 29073-3 ≥15 Ν - wet

 $\mathbf{>}$

SURGICAL DRAPES

STANDARDS GUIDE

> SURGICAL GOWNS

| | | | REQUIREMENT | | | |
|---|----------------|--------------------------|--------------------------|-------------------------------|--------------------------|-------------------------------|
| CHARACTERISTICS | STANDARD | UNIT | STANDARD R | EQUIREMENT | HIGH REQUIREMENT | |
| | | | CRITICAL PRODUCT AREA | LESS CRITICAL PRODUCT AREA | CRITICAL PRODUCT AREA | LESS CRITICAL PRODUCT AREA |
| microbial penetration - dry | EN ISO 22612 | CFU | not required | ≤ 300 ª | not required | ≤ 300 ª |
| microbial penetration - wet | EN ISO 22610 | I _B | ≥2,8 ^b | not required | 6,0 ^{bc} | not required |
| cleanliness - microbial bioburden | EN ISO 11737-1 | CFU/ 100 cm ² | ≤ 300 | ≤ 300 | ≤ 300 | ≤ 300 |
| particulate release | EN ISO 9073-10 | log₁₀ (lint count) | ≤ 4,0 | ≤ 4,0 | ≤ 4,0 | ≤ 4,0 |
| liquid penetration | EN ISO 811 | $\rm cm~H_2O$ | ≥20 | ≥10 | ≥100 | ≥10 |
| bursting strength - dry | EN ISO 13938-1 | kPa | ≥40 | ≥40 | ≥40 | ≥40 |
| bursting streangth - wet | EN ISO 13938-1 | kPa | ≥40 | not required | ≥40 | not required |
| tensile strength - dry | EN 29073-3 | Ν | ≥20 | ≥20 | ≥20 | ≥20 |
| tensile strength - wet | EN 29073-3 | Ν | ≥20 | not required | ≥20 | not required |

Standard EN 13795-1:2019 -

Essential requirements for surgical drapes and surgical gowns results from directive 93/42/EEC and have been included in EN 13795 standard "Surgical drapes, gowns and clean air suits, used as medical devices for patients, clinical staff and equipment". The purpose of EN 13795 standard is to maintain stable level of surgical drapes and surgical gowns safety for the whole using period as well as minimalization of direct microorganisms transmission from medical staff to patient during surgical procedures.



| REQUIREMENT | | | | | |
|---------------|-------------------------------|--------------------------|-------------------------------|--|--|
| IDARD RI | EQUIREMENT | HIGH REQUIREMENT | | | |
| CAL T AREA | LESS CRITICAL PRODUCT AREA | CRITICAL PRODUCT AREA | LESS CRITICAL PRODUCT AREA | | |
| quired | ≤ 300 ª | not required | ≤ 300 ª | | |
| 3 ь | not required | 6,0 ^{bc} | not required | | |
| 00 | ≤ 300 | ≤ 300 | ≤ 300 | | |
| 0 | ≤ 4,0 | ≤ 4,0 | ≤ 4,0 | | |
| 0 | ≥10 | ≥100 | ≥10 | | |
| 0 | ≥40 | ≥40 | ≥40 | | |
| C | not required | ≥40 | not required | | |
| 5 | ≥15 | ≥20 | ≥20 | | |
| 5 | not required | ≥20 | not required | | |

ZARYS

BASIC TYPES OF MATERIALS USED IN SURGICAL DRAPES:

PF DRAPE - two layer material of the drape containing absorbency layer of celulose on top and PE foil layer of the drape on back

NF, NFA, NFF, NFFA, NFFAR DRAPE - two layer material of the drape containing absorbency layer of non-woven material on top

NFN, NFNA, NFNF, NFNFA DRAPE - three layer material on the drape containing absorbency non-woven layer on top, PE foil layer

STERILE SURGICAL DRAPES - 2-LAYERED

| REF | dimensions | packaging unit | bulk packaging |
|--------------|---------------|----------------|----------------|
| AT-NF-S 1 E | 45cm x 45cm | 100 pcs | 4 x 100 pcs |
| AT-NF-S 3 E | 45cm x 60cm | 100 pcs | 4 x 100 pcs |
| AT-NF-S 4 E | 45cm x 75cm | 100 pcs | 4 x 100 pcs |
| AT-NF-S 5 E | 75cm x 75cm | 50 pcs | 4 x 50 pcs |
| AT-NF-S 9 E | 75cm x 90cm | 50 pcs | 4 x 50 pcs |
| AT-NF-S 13 E | 100cm x 150cm | 50 pcs | 4 x 50 pcs |
| AT-NF-S 15 E | 150cm x 180cm | 20 pcs | 4 x 20 pcs |
| AT-NF-S 21 E | 37cm x 45cm | 100 pcs | 4 x 100 pcs |

STERILE SURGICAL DRAPES - 2-LAYERED

| REF | dimensions | packaging unit | bulk packaging |
|------------|---------------|----------------|----------------|
| AT-NF-S 1 | 45cm x 45cm | 1 pc. | 2 x 240 pcs |
| AT-NF-S 2 | 50cm x 50cm | 1 pc. | 2 x 200 pcs |
| AT-NF-S 3 | 50cm x 60cm | 1 pc. | 2 x 280 pcs |
| AT-NF-S 4 | 50cm x 75cm | 1 pc. | 2 x 160 pcs |
| AT-NF-S 5 | 75cm x 75cm | 1 рс. | 2 x 150 pcs |
| AT-NF-S 8 | 80cm x 90cm | 1 рс. | 2 x 150 pcs |
| AT-NF-S 9 | 75cm x 90cm | 1 pc. | 2 x 150 pcs |
| AT-NF-S 10 | 90cm x 120cm | 1 рс. | 2 x 80 pcs |
| AT-NF-S 12 | 75cm x 140cm | 1 рс. | 2 x 50 pcs |
| AT-NF-S 13 | 100cm x 150cm | 1 pc. | 2 x 50 pcs |
| AT-NF-S 14 | 120cm x 150cm | 1 pc. | 2 x 50 pcs |
| AT-NF-S 15 | 150cm x 180cm | 1 pc. | 2 x 20 pcs |
| AT-NF-S 16 | 150cm x 200cm | 1 pc. | 2 x 20 pcs |
| AT-NF-S 19 | 160cm x 220cm | 1 pc. | 2 x 20 pcs |
| AT-NF-S 20 | 150cm x 240cm | 1 pc. | 2 x 20 pcs |
| AT-NF-S 22 | 170cm x 300cm | 1 pc. | 2 x 20 pcs |


STERILE SURGICAL DRAPES WITH ADHESIVE EDGE - 2-LAYERED


| REF | dimensions | packaging unit | bulk packaging |
|---------------|---------------|----------------|----------------|
| AT-NFA-S 4 E | 45cm x 75cm | 100 pcs | 4 x 100 pcs |
| AT-NFA-S 9 E | 75cm x 90cm | 50 pcs | 4 x 50 pcs |
| AT-NFA-S 10 E | 90cm x 120cm | 50 pcs | 4 x 50 pcs |
| AT-NFA-S 13 E | 100cm x 150cm | 50 pcs | 4 x 50 pcs |
| AT-NFA-S 15 E | 150cm x 180cm | 20 pcs | 4 x 20 pcs |
| AT-NFA-S 16 E | 150cm x 200cm | 20 pcs | 4 x 20 pcs |
| AT-NFA-S 20 E | 150cm x 240cm | 20 pcs | 4 x 20 pcs |


SURGICAL DRAPES

STERILE DRAPES WITH CELLULOSE LAYER - 2-LAYERED

| REF | dimensions | packaging unit | bulk packaging |
|-----------|--------------|----------------|----------------|
| AT-PF-S 1 | 50cm x 50cm | 1 pc. | 2 x 250 pcs |
| AT-PF-S 3 | 50cm x 70cm | 1 pc. | 2 x 200 pcs |
| AT-PF-S 5 | 80cm x 80cm | 1 pc. | 2 x 100 pcs |
| AT-PF-S 6 | 75cm x 90cm | 1 pc. | 2 x 100 pcs |
| AT-PF-S 8 | 90cm x 150cm | 1 pc. | 2 x 50 pcs |
| | | | |












STERILE SURGICAL DRAPES WITH ADHESIVE EDGE - 2-LAYERED

| REF | dimensions | packaging unit | bulk packaging |
|-------------|---------------|----------------|----------------|
| AT-NFA-S 1 | 45cm x 45cm | 1 pc. | 2 x 230 pcs |
| AT-NFA-S 3 | 50cm x 60cm | 1 pc. | 2 x 100 pcs |
| AT-NFA-S 4 | 50cm x 75cm | 1 pc. | 2 x 210 pcs |
| AT-NFA-S 5 | 75cm x 75cm | 1 pc. | 2 x 135 pcs |
| AT-NFA-S 9 | 75cm x 90cm | 1 pc. | 2 x 125 pcs |
| AT-NFA-S 10 | 90cm x 120cm | 1 pc. | 2 x 80 pcs |
| AT-NFA-S 13 | 100cm x 150cm | 1 pc. | 2 x 50 pcs |
| AT-NFA-S 15 | 150cm x 180cm | 1 pc. | 2 x 20 pcs |
| AT-NFA-S 16 | 150cm x 200cm | 1 pc. | 2 x 20 pcs |
| AT-NFA-S 20 | 150cm x 240cm | 1 pc. | 2 x 20 pcs |
| AT-NFA-S 22 | 170cm x 300cm | 1 pc. | 2 x 20 pcs |



0

 \bigcirc

economic version

41ª 40 g/m²


BHAT SYEARS

STERILE SURGICAL DRAPES WITH CENTRAL FENESTRATION - 2-LAYERED

| BEF | dimensions | fenestration | packaging unit | bulk packaging |
|--------------|-------------|--------------|----------------|----------------|
| AT-NFF-S 2 E | 45cm x 75cm | 6cm x 8cm | 100 pcs | 4 x 100 pcs |
| AT-NFF-S 3 E | 45cm x 60cm | ø 6cm | 100 pcs | 4 x 100 pcs |
| AT-NFF-S 4 E | 75cm x 90cm | ø 7cm | 50 pcs | 4 x 50 pcs |
| AT-NFF-S 5 E | 75cm x 90cm | 6cm x 8cm | 50 pcs | 4 x 50 pcs |




| REF | dimensions | fenestration | packaging unit | bulk packaging |
|-------------|-------------|--------------|----------------|----------------|
| AT-NFF-S 1 | 45cm x 45cm | ø 6cm | 1 pc. | 2 x 200 pcs |
| AT-NFF-S 2 | 50cm x 75cm | 6cm x 8cm | 1 pc. | 2 x 120 pcs |
| AT-NFF-S 3 | 50cm x 75cm | ø 7cm | 1 pc. | 2 x 120 pcs |
| AT-NFF-S 4 | 75cm x 90cm | ø 7cm | 1 pc. | 2 x 120 pcs |
| AT-NFF-S 5 | 75cm x 90cm | 6cm x 8cm | 1 pc. | 2 x 120 pcs |
| AT-NFF-S 21 | 50cm x 60cm | ø 7cm | 1 pc. | 2 x 300 pcs |




STERILE SURGICAL DRAPES WITH ADHESIVE CENTRAL FENESTRATION - 2-LAYERED

| REF | dimensions | fenestration | packaging unit | bulk packaging |
|---------------|--------------|--------------|----------------|----------------|
| AT-NFFA-S 1 E | 45cm x 45cm | ø 5cm | 100 pcs | 4 x 100 pcs |
| AT-NFFA-S 2 E | 45cm x 75cm | 6cm x 8cm | 100 pcs | 4 x 100 pcs |
| AT-NFFA-S 3 E | 50cm x 60cm | ø 6cm | 100 pcs | 4 x 100 pcs |
| AT-NFFA-S 4 E | 75cm x 90cm | ø 7cm | 50 pcs | 4 x 50 pcs |
| AT-NFFA-S 5 E | 75cm x 90cm | 6cm x 8cm | 50 pcs | 4 x 50 pcs |
| AT-NFFA-S 6 E | 90cm x 120cm | ø 8cm | 50 pcs | 4 x 50 pcs |



STERILE SURGICAL DRAPES WITH ADHESIVE CENTRAL FENESTRATION - 2-LAYERED

| REF | dimensions | fenestration | packaging unit | bulk packaging |
|--------------|---------------|--------------|----------------|----------------|
| AT-NFFA-S 1 | 45cm x 45cm | ø 6cm | 1 pc. | 2 x 200 pcs |
| AT-NFFA-S 2 | 50cm x 75cm | 6cm x 8cm | 1 pc. | 2 x 120 pcs |
| AT-NFFA-S 3 | 50cm x 75cm | ø 7cm | 1 pc. | 2 x 120 pcs |
| AT-NFFA-S 4 | 75cm x 90cm | ø 7cm | 1 pc. | 2 x 100 pcs |
| AT-NFFA-S 5 | 75cm x 90cm | 6cm x 8cm | 1 pc. | 2 x 100 pcs |
| AT-NFFA-S 6 | 90cm x 120cm | ø 10cm | 1 pc. | 2 x 60 pcs |
| AT-NFFA-S 7 | 120cm x 150cm | ø 7cm | 1 pc. | 2 x 50 pcs |
| AT-NFFA-S 8 | 120cm x 150cm | ø 10cm | 1 pc. | 2 x 50 pcs |
| AT-NFFA-S 11 | 50cm x 60cm | 6cm x 8cm | 1 pc. | 2 x 200 pcs |




STERILE SURGICAL DRAPES WITH ADHESIVE CENTRAL **ADJUSTABLE FENESTRATION - 2-LAYERED**

| REF | dimensions | fenestration | packaging unit | bulk pack |
|----------------|-------------|--------------|----------------|-------------|
| AT-NFFAR-S 2 E | 45cm x 75cm | 6cm x 8cm | 100 pcs | 4 x 100 pcs |
| AT-NFFAR-S 4 E | 75cm x 90cm | ø 7cm | 50 pcs | 4 x 50 pcs |
| | | | | |



| REF | dimensions | fenestration | packaging unit | bulk pack |
|--------------|-------------|--------------|----------------|-------------|
| AT-NFFAR-S 2 | 45cm x 75cm | 6cm x 8cm | 1 pc. | 2 x 150 pcs |
| AT-NFFAR-S 4 | 75cm x 90cm | ø 7cm | 1 pc. | 2 x 90 pcs |
| | | 、 、 | | |

| | REF | dimensions | fenestration | packaging unit | bulk packa |
|-----|---------------|-------------|--------------|----------------|------------|
| | AT-NFZRS-S 1 | 75cm x 80cm | 6cm x 15cm | 1 pc. | 2x60 pcs |
| NEW | AT-NFZRS-S 2 | 60cm x 75cm | 6cm x 9,5cm | 1 pc. | 2x120 pcs |
| | AT-NFZRS-S 11 | 75cm x 90cm | 10cm x 15cm | 1 pc. | 2x60 pcs |
| | | | | | |



| REF | dimensions | packaging unit | bulk packaging |
|-----------|---------------|----------------|----------------|
| AT-SD-S 3 | 100cm x 150cm | 1 pc. | 100 pcs |
| AT-SD-S 4 | 150cm x 190cm | 1 pc. | 60 pcs |
| AT-SD-S 5 | 150cm x 260cm | 1 pc. | 50 pcs |



414 70 g/m²

STERILE SURGICAL DRAPES WITH ADHESIVE CENTRAL FENESTRATION - 3-LAYERED

| REF | dimensions | fenestration | packaging unit | bulk package |
|--------------|--------------|--------------|----------------|--------------|
| AT-NFNFA-S 2 | 50cm x 70cm | 6cm x 8cm | 1 pc. | 200 pcs |
| AT-NFNFA-S 3 | 50cm x 75cm | ø 7cm | 1 pc. | 100 pcs |
| AT-NFNFA-S 4 | 75cm x 90cm | ø 7cm | 1 pc. | 100 pcs |
| AT-NFNFA-S 5 | 75cm x 90cm | 6cm x 8cm | 1 pc. | 100 pcs |
| AT-NFNFA-S 6 | 90cm x 120cm | ø 10cm | 1 pc. | 63 pcs |
| | | | | |



| ERILE R | ATEX | 5YEARS |
|---------|------|------------|
| | | SHELF LIFE |

STERILE SURGICAL DRAPES WITH ADHESIVE FENESTRATION, REINFORCED

| REF | dimensions | fenestration | packaging unit | bulk packa |
|--------------|---------------|--------------|----------------|------------|
| AT-NFZRF-S 1 | 120cm x 170cm | ø 12cm | 1 pc. | 50 pcs |
| AT-NFZRF-S 2 | 200cm x 280cm | ø 10cm | 1 pc. | 20 pcs |
| AT-NFZRF-S 3 | 200cm x 320cm | ø 15cm | 1 pc. | 20 pcs |
| | | | | |





STERILE SURGICAL DRAPES - 3-LAYERED

| | REF | dimensions | packaging unit | bulk packaging |
|-----|-------------|---------------|----------------|----------------|
| | AT-NFN-S 1 | 45cm x 45cm | 1 pc. | 300 pcs |
| | AT-NFN-S 2 | 50cm x 50cm | 1 pc. | 300 pcs |
| | AT-NFN-S 3 | 50cm x 60cm | 1 pc. | 150 pcs |
| | AT-NFN-S 4 | 50cm x 70cm | 1 pc. | 160 pcs |
| | AT-NFN-S 5 | 75cm x 75cm | 1 pc. | 100 pcs |
| NEW | AT-NFN-S 8 | 80cm x 90cm | 1 pc. | 80 pcs |
| | AT-NFN-S 9 | 75cm x 90cm | 1 pc. | 80 pcs |
| | AT-NFN-S 10 | 90cm x 120cm | 1 pc. | 75 pcs |
| | AT-NFN-S 13 | 100cm x 150cm | 1 pc. | 50 pcs |
| NEW | AT-NFN-S 14 | 120cm x 150cm | 1 pc. | 50 pcs |
| | AT-NFN-S 15 | 150cm x 180cm | 1 pc. | 50 pcs |
| | AT-NFN-S 16 | 150cm x 210cm | 1 pc. | 25 pcs |
| NEW | AT-NFN-S 19 | 160cm x 220cm | 1 pc. | 25 pcs |
| | AT-NFN-S 20 | 150cm x 240cm | 1 pc. | 25 pcs |

STERILE R 4<u>1</u>4 75 g/m²

STERILE SURGICAL DRAPES WITH ADHESIVE EDGE - 3-LAYERED

| | REF | dimensions | packaging unit | bulk packaging |
|-----|--------------|---------------|----------------|----------------|
| | AT-NFNA-S 2 | 50cm x 50cm | 1 pc. | 256 pcs |
| | AT-NFNA-S 4 | 50cm x 70cm | 1 pc. | 160 pcs |
| | AT-NFNA-S 5 | 75cm x 75cm | 1 pc. | 100 pcs |
| IEW | AT-NFNA-S 8 | 80cm x 90cm | 1 pc. | 80 pcs |
| | AT-NFNA-S 9 | 75cm x 90cm | 1 pc. | 90 pcs |
| | AT-NFNA-S 13 | 100cm x 150cm | 1 pc. | 80 pcs |
| | AT-NFNA-S 15 | 150cm x 180cm | 1 pc. | 25 pcs |
| | AT-NFNA-S 16 | 150cm x 210cm | 1 pc. | 50 pcs |
| IEW | AT-NFNA-S 19 | 160cm x 220cm | 1 pc. | 25 pcs |
| | AT-NFNA-S 20 | 150cm x 240cm | 1 pc. | 25 pcs |





STERILE SURGICAL DRAPES WITH CENTRAL FENESTRATION - 3-LAYERED


| REF | dimensions | fenestration | packaging uni | t 🖗 bulk packaging |
|-------------|--------------|--------------|---------------|--------------------|
| AT-NFNF-S 1 | 45cm x 45cm | ø 6cm | 1 pc. | 300 pcs |
| AT-NFNF-S 2 | 50cm x 75cm | 6cm x 8cm | 1 pc. | 160 pcs |
| AT-NFNF-S 3 | 50cm x 75cm | ø 7cm | 1 pc. | 200 pcs |
| AT-NFNF-S 4 | 75cm x 90cm | ø 7cm | 1 pc. | 80 pcs |
| AT-NFNF-S 5 | 75cm x 90cm | 6cm x 8cm | 1 pc. | 80 pcs |
| AT-NFNF-S 6 | 90cm x 120cm | ø 10cm | 1 pc. | 100 pcs |



154 www.zarys.com

aging 0

kaging



| S | URGICAL DRAPES |
|---|---------------------|
| S | SPECIALISTIC DRAPES |

| AT-SD-ABD1-IF-C | ABDOMINAL | ORAPE 1- | IF |
|-----------------|---------------|----------|---|
| REF | dimensions | quantity | component |
| AT-SD-ABD1-IF-C | 200cm x 320cm | 1 pc. | drape with adhesive fenestration 25cm x 30cm filled with an incise film, reinforced |
| AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, reinforced |

V bulk packaging: 20 pcs

AT-SD-S 4


STERILE R

| AT-SD-ABPER1-S | ABDO-PERINEAL DRAPE 1 | | | |
|----------------|-----------------------|----------|---|--------|
| REF | dimensions | quantity | component | |
| AT-SD-ABPER1-S | 180/280cm x 220cm | 1 pc. | drape integrated with limb covers, adhesive trapezium-shaped fenestration 30cm x 20cm x and perineal fenestration 10cm x 15cm with covering flap | < 20cm |

1 pc.

| bulk packaging: 16 pcs | | SHELF LIFE |
|---------------------------|--|------------|
|---------------------------|--|------------|


150cm x 190cm



150x190cm

AT-SD-LAP3-SP LAPAROSCOPY DRAPE 3 REF dimensions quantity component drape integrated with limb covers, adhesive fenestration 25cm x 30cm, reinforced, AT-SD-LAP3-SP 200/270cm x 300cm 1 pc. with 8 cable holders Velcro type AT-SD-S 4 150cm x 190cm instrument table cover, reinforced 1 pc.

| | SYEARS |
|---------------------------|----------------|
| bulk packaging: 12 pcs | SHELF LIFE |



150x190cm

200/270x300cm

111

200x320cm



| AT-SD-ORT1-SP | ORTHOPAEDI | C DRAPE 1 | |
|---------------------------|---------------|-----------|--|
| REF | dimensions | quantity | component |
| AT-SD-ORT1-SP | 150cm x 250cm | 1 pc. | drape with adhesive U-sp 7cm x 70cm, reinforced |
| AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, r |
| bulk packaging: 20 pcs | | | ure - |

| AT-SD-ORT2-SP | ORTHOPAEDIC DRAPE 2 | | | |
|---------------|---------------------|----------|--|--|
| REF | dimensions | quantity | component | |
| AT-SD-ORT2-SP | 230cm x 260cm | 1 pc. | drape with adhesive L 20cm x 100cm, reinfor | |
| AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cove | |
| | | | rs , | |

bulk packaging: 16 pcs

| AT-SD-ORT3-C | ORTHOPAEDIC DRAPE 3 | | | | |
|--------------|---------------------|----------|--|--|--|
| REF | dimensions | quantity | component | | |
| AT-SD-ORT3-C | 150cm x 200cm | 1 pc. | drape with adhesive U-sp 6cm x 65cm, reinforced | | |
| AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, r | | |
| | | | | | |



| AT-SD-ORT4-C | ORTHOPAEDIC | C DRAPE 4 | |
|---------------------------|---------------|-----------|--|
| REF | dimensions | quantity | component |
| AT-SD-ORT4-C | 200cm x 260cm | 1 pc. | drape with adhesive U-s 8.5cm x 85cm, reinforce |
| AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, |
| bulk packaging: 15 pcs | | | - |






150x190cm 200x260cm

| AT-SD-EXT1-SP | EXTREMITY DR | APE 1 | | | | | AT-SD-EXT5-C | EXTREMITY DR | APE 5 | |
|----------------------------|---------------|-------------|--|----------|----------------------------|----------------------------|-------------------------------|---------------------|-----------|--|
| REF | dimensions | quantity | component | | | | REF | dimensions | quantity | component |
| AT-SD-EXT1-SP | 150cm x 180cm | 1 pc. | drape with elastic fenestration ø 3,5cm, reinforced | | | • • | AT-SD-EXT5-C | 200cm x 300cm | 1 pc. | drape with elastic fenestration ø 3,5cm, reinforced |
| AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, reinforced | | | | AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, reinforced |
| | | | | | 150 x 190 cm | 150 x 180 cm | | | | |
| bulk packaging: 24 pcs | STERILE R | SYEA | 55 197 | | | | bulk packaging: 20 pcs | STERILE R | SHEF UP | τ. |
| AT-SD-EXT2-SP | EXTREMITY DR | APE 2 | | * | | | | | | |
| REF | dimensions | quantity | component | | | | | | | |
| AT-SD-EXT2-SP AT-SD-S 4 | 230cm x 300cm | 1 pc. | drape with elastic fenestration ø 3,5cm, reinforced instrument table cover, reinforced | | | | G | YNAECOLOGY A | ND UROLO | DGY |
| AT-50-54 | 150cm x 190cm | 1 pc. | instrument table cover, reinforced | | | _ <u></u> | | | | |
| bulk packaging: 18 pcs | STERILE R | | | | | | AT-SD-GYN1-S | GYNAECOLOG | Y DRAPE 1 | |
| | | | - | | | | REF | dimensions | quantity | component |
| | | | | | | | AT-SD-GYN1-S | 160/240cm x 180cm | 1 pc. | drape integrated with limb covers, ac perineal fenestration 10cm x 15cm |
| | | | | | 150x190cm | 230x300cm | AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, reinforced |
| AT-SD-EXT3-SP REF | EXTREMITY DR | APE 3 | component | * | | | bulk packaging: 16 pcs | STERILE R | SHELF U | • |
| AT-SD-EXT3-SP | 230cm x 300cm | 1 pc. | drape with elastic fenestration ø 7cm, rei | inforced | | | | | | |
| AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, reinforced | | | | AT-SD-GYN2-S | GYNAECOLOG | Y DRAPE 2 | |
| | | | | | | | REF | dimensions | quantity | component |
| bulk packaging: 18 pcs | STERILE R | SYEAR | are the second sec | | | | AT-SD-GYN2-S | 160/240cm x 180cm | 1 pc. | drape integrated with limb covers, ac perineal fenestration 10cm x 15cm, i fluid collection pouch 35cm x 50cm v formable molder, filter and tubing po |
| | | | | | | | AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, reinforced |
| | | | | | 150x190cm | 230x300cm | bulk packaging: 8 pcs | STERILE R | SYEARS | e e |
| AT-SD-EXT4-C | EXTREMITYDR | APE 4 | | | | | | | | |
| REF | dimensions | quantity | component | | | | AT-SD-DEL1-SP | DELIVERY DRA | DF 1 | |
| AT-SD-EXT4-C | 200cm x 320cm | 1 pc. | drape with elastic fenestration 6cm x 8cm, reinforced | | | | | | | |
| AT-SD-S 4 | 150cm x 190cm | 1 pc. | instrument table cover, reinforced | | | | | dimensions | quantity | component |
| | | | | | | | AT-SD-DEL1-SP | 114cm x 150cm | 1 pc. | delivery drape integrated with fluid c pouch 78cm x 80cm with formable m filter and tubing port, reinforced |
| hully nackaging | | X NHT SYEAT | | | | | | | | |






150x190cm

200x320cm

158 www.zarys.com

bulk packaging: 15 pcs



STERILE R STERILE R

without instrument table cover wrapping

 \bigotimes

bulk packaging: 45 pcs

www.zarys.com 159


114x150cm

| Орн | THALMOLOG | ίΥ | | | |
|--|--|--|---|-----------|----------------------------|
| | | | | | |
| T-SD-OPHT1-IF-S | OPHTHALMIC | DRAPE 1 | -IF | | |
| F | dimensions | quantity | component | | |
| T-SD-OPHT1-IF-S | 100cm x 100cm | 1 pc. | ophthalmic drape with fenestration 10cm x 10cm filled with an incise film | | <i>*</i> |
| r-SD-S 3 | 100cm x 150cm | 1 pc. | instrument table cover, reinforced | | |
| | | | | 100x150cm | 100 x 100 cm |
| packaging: 50 sets | | SHELF | RS V | | |
| | | | | | |
| -SD-OPHT2-IF-C | OPHTHALMIC | DRAPE 2 | -IF | | |
| F | dimensions | quantity | component | | |
| -SD-OPHT2-IF-C | 100cm x 100cm | 1 pc. | ophthalmic drape with fenestration 10cm x 10cm filled with an incise film, integrated with fluid collec- tion pouch and with 2 formable molders, reinforced | | |
| T-SD-S 3 | 100cm x 150cm | 1 pc. | instrument table cover, reinforced | | |
| | <u>ش</u> | | | | |
| packaging: hyd | drophobic S material | RILE R | | 100x150cm | 100 x 100 cm |
| 00 SELS | Sindlendi | | | | |
| T-SD-OPHT3-IF-C | OPHTHALMIC | | -IF | | |
| | | | | | |
| F T-SD-OPHT3-IF-C | dimensions 120cm x 170cm | quantity 1 pc. | component ophthalmic drape with fenestration 10cm x 10cm filled | | |
| | | | with an incise film, integrated with 2 fluid collection | | |
| | | | pouches and with 1 formable molder, reinforced | | |
| r-sd-s 3 | 100cm x 150cm | 1 pc. | pouches and with 1 formable molder, reinforced | | |
| r-sD-s 3 | | 1 pc. | | 100x150cm | 120x170cm |
| packaging: hyd | က drophobic | 1 pc. | | 100x150cm | 120 x 170 cm |
| packaging: hyd | ᠭ | 1 pc. | | 100x150cm | 120x170cm |
| packaging: hyc 40 sets SM | က drophobic | 1 pc. | | 100x150cm | 120 x 170 cm |
| packaging: hyc 40 sets SM | frophobic S material | | instrument table cover, reinforced | 100x150cm | 120 x 170 cm |
| packaging: hyc 40 sets SM | frophobic S material | | | 100x150cm | 120 x 170 cm |
| packaging: 40 sets SM | frophobic S material | formable molde | instrument table cover, reinforced | 100x150cm | 120x170cm |
| packaging: 40 sets SM ERILE R | Trophobic S material | formable molde | instrument table cover, reinforced | 100x150cm | 120x170cm |
| packaging: 40 sets SM | Grophobic S material | formable molde | instrument table cover, reinforced | 100x150cm | 120x170cm |
| | Trophobic S material Trophobic S material OPHTHALMIC dimensions | formable molder DRAPE 4 quantity | instrument table cover, reinforced | 100x150cm | 120x170cm |
| packaging: 40 sets SM | Trophobic S material Trophobic S material OPHTHALMIC dimensions | formable molder DRAPE 4 quantity | instrument table cover, reinforced | 100×150cm | 120x170cm |
| packaging: to sets hyde SM ERILE R SD-OPHT4-IF-C F T-SD-OPHT4-IF-C | Trophobic S material | formable molder DRAPE 4 quantity 1 pc. | instrument table cover, reinforced | 100x150cm | 120x170cm |

ZARYS

SURGICAL COVERS, ACCESSORIES







STERILE CAMERA CABLE COVERS - PERFORATED, WITH 2 ADHESIVE TAPES 🙆 BESTSELLER

| REF | dimensions | packaging unit | bulk packaging |
|-----------|--------------|----------------|----------------|
| AT-SP-S 3 | 14cm x 250cm | 1 pc. | 160 pcs |
| AT-SP-S 4 | 16cm x 200cm | 1 pc. | 200 pcs |
| AT-SP-S 5 | 16cm x 250cm | 1 pc. | 160 pcs |
| AT-SP-S 6 | 16cm x 300cm | 1 pc. | 200 pcs |





sterile camera cover

STERILE CAMERA CABLE COVER - UNPERFORATED, WITH AN ADHESIVE TAPE

| dimensions | packaging unit | bulk packaging |
|-------------|----------------|-------------------|
| 7cm x 120cm | 1 pc. | 200 pcs |
| 7cm x 170cm | 1 pc. | 200 pcs |
| | 7cm x 120cm | 7cm x 120cm 1 pc. |

STERILE CIRCULAR BANDED COVERS FOR MEDICAL DEVICES

| REF | dimensions in stretched state | maximum size of covered device | packaging unit | bulk p |
|------------|-------------------------------|-----------------------------------|-------------------|---------|
| AT-SP-S 9 | 65cm x 85cm | 35cm x 45cm | 1 pc. | 200 pcs |
| AT-SP-S 10 | 85cm x 85cm | 50cm x 50cm | 1 pc. | 100 pcs |
| AT-SP-S 11 | 105cm x 105cm | 60cm x 60cm | 1 pc. | 100 pcs |
| AT-SP-S 12 | 120cm x 120cm | 70cm x 70cm | 1 pc. | 100 pcs |
| AT-SP-S 13 | 140cm x 140cm | 80cm x 80cm | 1 pc. | 50 pcs |

STERILE R

STERILE RECTANGULAR BANDED COVERS FOR MEDICAL DEVICES

| REF | dimensions | packaging unit | bulk packaging |
|------------|--------------|----------------|----------------|
| AT-SP-S 14 | 70cm x 80cm | 1 pc. | 100 pcs |
| AT-SP-S 15 | 90cm x 100cm | 1 pc. | 50 pcs |
| AT-SP-S 16 | 80cm x 150cm | 1 pc. | 50 pcs |

strong PE film 65 µm



STERILE C-ARM COVER SET



| REF | dimensions | packaging unit | bulk packaging | | |
|-------|-------------|----------------|----------------|--|--|
| AT-AC | 35cm x 75cm | 2 pcs | 80 pcs | | |





100x160cm



80x150cm







STERILE SURGICAL TAPES

| REF | dimensions | description | packaging unit | bulk packaging | |
|-----------|-------------|----------------------------|----------------|----------------|--|
| AT-SP-S 1 | 10cm x 50cm | surgical tape, OP tape | 1 pc. | 1100 pcs | |
| AT-SP-S 2 | 2cm x 22cm | surgical tape, Velcro type | 1 pc. | 100 pcs | |



STERILE R STERILE R





Velcro surgical tape

STERILE ADHESIVE POUCHES

| REF | dimensions | description | packaging unit | bulk packaging |
|------------|-------------|--|----------------|----------------|
| AT-SP-S 20 | 30cm x 40cm | adhesive pouch, one-chamber | 1 pc. | 180 pcs |
| AT-SP-S 21 | 50cm x 60cm | adhesive pouch, one-chamber | 1 pc. | 180 pcs |
| AT-SP-S 22 | 40cm x 30cm | adhesive pouch, two-chamber 15cm x 30cm left chamber 25cm x 30cm right chamber | | 250 pcs |



STERILE ADHESIVE FLUID COLLECTION POUCH

| | STERILE ADHESIVE FLUID COLLECTION POUCH | | | | | \bigcirc | |
|-----|---|--------------|----------------|---|----------------|------------|--------------------------|
| | REF | dimensions | packaging unit | description | bulk packaging | | |
| NEW | AT-SP-S 18 | 60cm x 100cm | 1 pc. | arthroscopy cone-shaped pouch with self-sealing fenestration and suction port | 36 pcs | • | |
| | AT-SP-S 19 | 50cm x 50cm | 1 pc. | cone-shaped pouch with filter and tubing port | 65 pcs | 60x100cm | 50 x 50 cm |
| | | | | | | | |

STERILE R



STERILE ABSORBENT DRAPE



STERILE INCISE FILM ELASTOFILM

| | REF | total dimensions | adhesive filed dimensions | packaging unit | bulk packaging |
|-----|--------|------------------|------------------------------|-------------------|-------------------|
| NEW | 808068 | 84cm x 56cm | 84cm x 50cm | 1 pc. | 6 pcs |
| Í | 808066 | 40cm x 42cm | 40cm x 36cm | 1 pc. | 20 pcs |
| | 808065 | 45cm x 55cm | 45cm x 49cm | 1 pc. | 20 pcs |
| Í | 808064 | 45cm x 28cm | 45cm x 22cm | 1 pc. | 20 pcs |
| | 808063 | 30cm x 28cm | 30cm x 22cm | 1 pc. | 20 pcs |
| | 808062 | 15cm x 28cm | 15cm x 22cm | 1 pc. | 40 pcs |
| i | 808061 | 14cm x 25cm | 14cm x 21cm | 1 pc. | 40 pcs |
| Í | 805018 | 55cm x 55cm | 55cm x 51cm | 1 pc. | 20 pcs |
| 1 | 805025 | 30cm x 40cm | 30cm x 36cm | 1 pc. | 20 pcs |
| | 805016 | 20cm x 30cm | 20cm x 26cm | 1 pc. | 30 pcs |







- flexible





• thin incise film made of polyurethane

• perfectly fits the patient's body

• double packed: externally paper-foil, internally paper

ZARYS

SURGICAL

SETS



| AT-S-BASIC1-S | BASIC SET 1 |
|---------------|-------------|
| | BROID DEI I |

| quantity | dimensions | component |
|----------|---------------|--|
| 1 pc. | 75cm x 90cm | drape with central adhesive fenestration |
| 1 pc. | 10cm x 17cm | tray |
| 10 pcs | 10cm x 10cm | gauze swab |
| 5 pcs | 30cm x 30cm | gauze ball |
| 1 рс. | 100cm x 150cm | instrument table cover, reinforced |





| quantity | dimensions | component |
|--------------|---------------|------------------------------------|
| 4 pcs | 75cm x 90cm | drape with adhesive edge |
| 1 рс. | 150cm x 190cm | instrument table cover, reinforced |
| bulk packagi | ng: STERI | |



50 sets

UNIVERSAL SET 2 AT-S-UNI2-S

| quantity | dimensions | component | |
|-------------------------|---------------|------------------------------------|------------|
| 2 pcs | 150cm x 180cm | drape with adhesive edge | |
| 2 pcs | 75cm x 90cm | drape with adhesive edge | |
| 4 pcs | 30cm x 30cm | absorbent towel | |
| 1 рс. | 150cm x 190cm | instrument table cover, reinforced | |
| | | | |
| bulk packagi 20 sets | ng: STERI | LE R DATES AND STEAMS | 2 x |



| quantity | dimensions | component |
|----------|---------------|---|
| 1 рс. | 170cm x 240cm | drape with adhesive edge, reinforced |
| 1 pc. | 170cm x 180cm | drape with adhesive edge, reinforced |
| 2 pcs | 75cm x 90cm | highly absorbent drape with adhesive edge |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |






75x**90**cm

Mayo stand cover, reinforced - red















30x**30**cm











170x180cm



150x190cm


www.zarys.com 167







| quantity | dimensions | component |
|----------|---------------|------------------------------------|
| 1 pc. | 150cm x 240cm | drape with adhesive edge |
| 1 pc. | 180cm x 180cm | drape with adhesive edge |
| 2 pcs | 75cm x 90cm | drape with adhesive edge |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |
| 1 pc. | L | surgical gown STANDARD |
| 2 pcs | XL | surgical gown STANDARD |



AT-S-UNI5-M-S UNIVERSAL SET 5-M




| quantity | dimensions | component |
|----------|---------------|------------------------------------|
| 1 pc. | 150cm x 240cm | drape with adhesive edge |
| 1 pc. | 180cm x 180cm | drape with adhesive edge |
| 2 pcs | 75cm x 90cm | drape with adhesive edge |
| 4 pcs | 30 cmx30 cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |
| - | | |






4<u>1</u>4 56 g/m²













AT-S-UNI5-GP-S UNIVERSALSET 5-GP

| quantity | dimensions | component |
|----------|---------------|---|
| 1 рс. | 150cm x 240cm | drape with adhesive edge |
| 1 pc. | 180cm x 180cm | drape with adhesive edge |
| 2 pcs | 75cm x 90cm | drape with adhesive edge |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |
| 1 pc. | XL | surgical gown STANDARD |
| 1 pc. | XXL | surgical gown, reinforced STANDARD PLUS |
| 1 рс. | 30cm x 40cm | adhesive pouch, one-chamber |






AT-S-UNI6-SP UNIVERSALSET 6 🙆 BESTSELLER


| quantity | dimensions | component |
|----------|---------------|--------------------------------------|
| 1 рс. | 150cm x 240cm | drape with adhesive edge, reinforced |
| 1 рс. | 180cm x 180cm | drape with adhesive edge, reinforced |
| 2 pcs | 75cm x 90cm | drape with adhesive edge, reinforced |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 рс. | 10cm x 50cm | OP tape |
| 1 рс. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 рс. | 150cm x 190cm | instrument table cover, reinforced |


4<u>1</u>4

60 g/m²









| AT-S-U | JNI7-C | UNI | VERSAL SET 7 |
|-------------------------|---------|------------------------|------------------------------------|
| quantity | dimensi | ons | component |
| 1 pc. | 150cm > | k 240cm | drape with adhesive edge |
| 1 pc. | 180cm > | < 180cm | drape with adhesive edge |
| 2 pcs | 75cm x | 90cm | drape with adhesive edge |
| 4 pcs | 30cm x | 30cm | absorbent towel |
| 1 pc. | 10cm x | 50cm | OP tape |
| 1 pc. | 80cm x | 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm > | < 190cm | instrument table cover, reinforced |
| | | | |
| bulk packagin 8 sets | g: | 3-layer dra PREMIUI | |
| | | | |

2

75x90cm

AT-S-UNI6-M-SP UNIVERSALSET 6-M

 \bigotimes

bulk packaging: 12 sets



75x**90**cm

| quantity | dimensions | component |
|----------|---------------|--------------------------------------|
| 1 pc. | 150cm x 240cm | drape with adhesive edge, reinforced |
| 1 pc. | 180cm x 180cm | drape with adhesive edge, reinforced |
| 2 pcs | 75cm x 90cm | drape with adhesive edge, reinforced |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |

4<u>1</u>4


56 g/m²


÷

highly absorbent

material

STERILE R REF REF STEAR






150x240cm

80x145cm

180x180cm






| quantity | aimensions | component |
|----------|---------------|--------------------------------------|
| 1 рс. | 150cm x 240cm | drape with adhesive edge, reinforced |
| 1 рс. | 180cm x 180cm | drape with adhesive edge, reinforced |
| 2 pcs | 75cm x 90cm | drape with adhesive edge, reinforced |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 рс. | 10cm x 50cm | OP tape |
| 1 рс. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |
| | | |



ċĵ highly absorbent material






75x**90**cm

150x190cm

| 170 www.zarys.com |
|-------------------|
|-------------------|






150x240cm






30x30cm **10**x**50**cm

80x145cm




180x180cm




150x190cm






150x240cm






80x145cm



180x180cm



150x190cm

AT-S-UNI9-SP PEDIATRIC UNIVERSAL SET 9




 \bigotimes bulk packaging: 8 sets

STERILE R STERILE R

٢

pediatric set





150x240cm


AT-S-UNI10-S UNIVERSAL SET 10

| quantity | dimensions | component |
|----------|---------------|------------------------------------|
| 1 pc. | 150cm x 230cm | drape with adhesive edge |
| 1 pc. | 150cm x 180cm | drape with adhesive edge |
| 2 pcs | 75cm x 90cm | drape with adhesive edge |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |









150x190cm

180x180cm

80x145cm

10x50cm



| τ pc. | 200CIII x 320CIII | with an incise film, reinforced |
|-------|-------------------|------------------------------------|
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 рс. | 10cm x 50cm | OP tape |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |
| | | |




AT-S-ABPER1-S ABDO-PERINEAL SET 1


| quantity | dimensions | component |
|----------|-------------------|--|
| 1 pc. | 180/270cm x 220cm | drape integrated with limb covers, adhesive trapezium-shaped fenestration 30cm x 20cm x 20c and perineal fenestration 10cm x 15cm with cover |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 рс. | 10cm x 50cm | OP tape |
| 1 рс. | 2cm x 22cm | surgical tape, Velcro type |
| 1 рс. | 14cm x 250cm | camera cable cover |
| 1 рс. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 рс. | 150cm x 190cm | instrument table cover, reinforced |


 \bigotimes bulk packaging: 12 sets





flap to secure the fenestration







AT-S-LAP1-IF-C LAPAROSCOPY SET 1-IF


| quantity | dimensions | component | | |
|--------------------------|---------------|---|-----------|---------------|
| 1 pc. | 200cm x 260cm | drape T-shaped with adhesive fenestration 30cm x 3 filled with an incise film | 0cm | |
| 4 pcs | 30cm x 30cm | absorbent towel | | |
| 1 pc. | 10cm x 50cm | OP tape | | |
| 1 pc. | 14cm x 250cm | camera cable cover | | |
| 2 pcs | 70cm x 120cm | limb cover | 11 | |
| 1 pc. | 80cm x 140cm | Mayo stand cover, reinforced - red | | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | | |
| bulk packagin 12 sets | SMS material | STERS STERFE 10x50cm 14x25 | 200x260cm | 4x 30x30cm |

AT-S-LAP2-SP LAPAROSCOPY SET 2 🙆 BESTSELLER


| quantity | dimensions | component |
|----------|-------------------|--|
| 1 pc. | 200/270cm x 300cm | drape integrated with limb covers, adhesive fenestration 25cm x 30cm, reinforced, with 8 cable holders Velcro type |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 2cm x 22cm | surgical tape, Velcro type |
| 1 pc. | 14cm x 250cm | camera cable cover |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |








cable holders, Velcro type



200/270x300cm





| quantity | dimensions | component |
|----------------|-------------------|--|
| 1 рс. | 200/270cm x 300cm | drape integrated with limb covers, adhesive fenestr 25cm x 30cm, reinforced, with 8 cable holders Velcr |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 2cm x 22cm | surgical tape, Velcro type |
| 1 pc. | 14cm x 250cm | camera cable cover |
| 1 pc. | 30cm x 40cm | adhesive pouch, one-chamber |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |
| bulk packaging | g: hydrophobic | STERILE R |

| 8 sets | SMS material | | SHELF LIFE | _ |
|------------------------------|--------------|--------------------------|------------|----|
| ¢. j | | | | |
| highly absorbent material | | | <u>_</u> | |
| | | 10 x 50 cm | 2x22cm | 14 |

AT-S-VAR1-S VARICOSE VEIN SET 1

| quantity | dimensions | component |
|----------|---------------|--|
| | | |
| 1 pc. | 200cm x 260cm | drape with adhesive U-split 8,5cm x 85cm |
| 1 pc. | 160cm x 180cm | drape with adhesive edge |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 80cm x 140cm | Mayo stand cover, reinforced - red |
| 1 рс. | 150cm x 190cm | instrument table cover, reinforced |
| | | |






| AT-S- | VAR2-C VARIO | COSE VEIN SET 2 |
|-------------------------|---------------------------------|---|
| quantity | dimensions | component |
| 1 рс. | 200cm x 260cm | drape with adhesive U-split 8,5cm x 95cm, integra genital area cover, reinforced |
| 1 pc. | 170cm x 180cm | drape with adhesive edge, reinforced |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 80cm x 140cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |
| | | |
| bulk packagi 10 sets | ng: hydrophobic SMS material | |
| U | A | STERILE R FATER STAR |










AT-S-EXT1-C EXTREMITY SET 1

material

SMS material

surgery



200x300cm 80x145cm 150x190cm

0

AT-S-EXT2-C EXTREMITY SET 2

STERILE R STERILE R

12 sets

| quantity | dimensions | component |
|----------|---------------|---|
| 1 pc. | 200cm x 300cm | drape with elastic fenestration ø 3,5cm, reinforced |
| 1 pc. | 160cm x 180cm | drape with adhesive edge |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |










7

limb

surgery



30x**30**cm



80x145cm


10x50cm

200x300cm

14cm x 250cm camera cable cover 30cm x 60cm elastic limb cover 60cm x 100cm molder and tubing port 80cm x 145cm Mayo stand cover, reinforced - red 150cm x 190cm instrument table cover, reinforced

bulk packaging: 10 sets

÷ highly absorbent material



200cm x 320cm drape with elastic fenestration 6cm x 8cm, reinforced

component

AT-S-ART1-C ARTHROSCOPY SET 1

| 4 pcs | 30cm x 30cm | absorbent towel | |
|-------|---------------|------------------------------------|--|
| 1 рс. | 10cm x 50cm | OP tape | |
| 1 рс. | 30cm x 60cm | elastic limb cover | |
| 1 рс. | 80cm x 145cm | Mayo stand cover, reinforced - red | |
| 1 рс. | 150cm x 190cm | instrument table cover, reinforced | |
| | | | |

 $\widehat{}$ bulk packaging: hydrophobic SMS material 14 sets

quantity dimensions

÷ highly absorbent material

SYEARS STERILE R


> 4x 10x50cm 30x30cm

quantity dimensions component drape with elastic fenestration 6cm x 8cm, reinforced 200cm x 320cm 30cm x 30cm absorbent towel 10cm x 50cm OP tape fluid collection pouch with elastic fenestration, formable

AT-S-ART2-C ARTHROSCOPY SET 2

 \bigotimes $\widehat{\mathcal{T}}$ hydrophobic SMS material

176 www.zarys.com



30x60cm

80x145cm






200x320cm









150x190cm

60x100cm

| quantity | dimensions | component | _ | |
|--------------------------|---------------------------------|--|------------------|-----------|
| 1 pc. | 200cm x 320cm | drape with elastic fenestration 6cm x 8cm, integrated fluid collection pouch, reinforced | | |
| 2 pcs | 30cm x 30cm | absorbent towel | - | |
| 3 pcs | 10cm x 50cm | OP tape | | |
| 1 pc. | 14cm x 250cm | camera cable cover | | |
| 1 pc. | 30cm x 60cm | elastic limb cover | _ | |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red | | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | _ | |
| bulk packagir 10 sets | ng: hydrophobic SMS material | highly absorbent material | | |
| STERILE | R 😿 🕅 (| SYEARS | | 200x320cm |
| | | | - - - | |



AT-S-ART3-C ARTHROSCOPY SET 3

AT-S-ART2-IP-C ARTHROSCOPY SET 2-IP

quantity dimensions component

| 200cm x 320cm | drape with elastic fenestration 6cm x 8cm, reinforced |
|---------------|---|
| 30cm x 30cm | absorbent towel |
| 10cm x 50cm | OP tape |
| 14cm x 250cm | camera cable cover |
| 30cm x 60cm | elastic limb cover |
| 80cm x 145cm | Mayo stand cover, reinforced - red |
| 150cm x 190cm | instrument table cover, reinforced |
| | 30cm x 30cm 10cm x 50cm 14cm x 250cm 30cm x 60cm 80cm x 145cm |

÷

2x


30x**30**cm

10x50cm





highly absorbent material



30x30cm

10x50cm

150x190cm

AT-S-ART5-SP ARTHROSCOPY SET 5

| quantity | aimensions | component |
|----------|---------------|---|
| 1 рс. | 230cm x 300cm | drape with elastic fenestration ø 7cm, reinforced |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 3 pcs | 10cm x 50cm | OP tape |
| 1 pc. | 22cm x 75cm | elastic limb cover |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 2 pcs | 150cm x 190cm | instrument table cover, reinforced |

----+



30x30cm

AT-S-ART4-IP-S ARTHROSCOPY SET 4-IP

| quantity | dimensions | component |
|----------|---------------|--|
| 1 pc. | 230cm x 320cm | drape with elastic fenestration ø 7cm, integrated flui collection pouch |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 3 pcs | 10cm x 50cm | OP tape |
| 1 рс. | 22cm x 75cm | elastic limb cover |
| 1 рс. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 2 pcs | 150cm x 190cm | instrument table cover, reinforced |





ontitu (



150x190cm

200x320cm


80x145cm

80x145cm


30x60cm

14x250cm

178 www.zarys.com



30x60cm



22x75cm

80x145cm

0

150x190cm


230x300cm








22x75cm







www.zarys.com 179

AT-S-HIP1-C HIP SET 1


| quantity | dimensions | component | | |
|-------------------------|-----------------------------|--|--|--|
| 1 pc. | 200cm x 260cm | drape with adhesive U-split 8,5cm x 85cm, reinforced | | |
| 1 pc. | 170cm x 240cm | drape with adhesive edge, reinforced | | |
| 4 pcs | 30cm x 30cm | absorbent towel | | |
| 2 pcs | 10cm x 50cm | OP tape | | |
| 1 pc. | 35cm x 120cm | elastic limb cover | | |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red | | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | | |
| | | | | |
| bulk packagii 8 sets | ng: hydrophob SMS materi | ic highly absorbent drape with | | |



200x260cm



AT-S-HIP2-C HIP SET 2



AT-S-HIP3-SP HIP SET 3

| quantity | dimensions | component |
|----------|---------------|---|
| 1 рс. | 230cm x 260cm | drape with adhesive U-split 20cm x 100cm, reinforce |
| 1 pc. | 150cm x 240cm | drape with adhesive edge, reinforced |
| 1 pc. | 180cm x 180cm | drape with adhesive edge, reinforced |
| 2 pcs | 75cm x 90cm | drape with adhesive edge, reinforced |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 3 pcs | 10cm x 50cm | OP tape |
| 1 pc. | 35cm x 120cm | elastic limb cover |
| 1 pc. | 10cm x 450cm | elastic bandage |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |
| | | |

 \bigotimes ÷ bulk packaging: highly absorbent 5 sets material





AT-S-SHLR1-C SHOULDER SET 1

 $\widehat{}$

hydrophobic


SMS material


| quantity | dimensions | component |
|----------|---------------|---|
| 1 pc. | 200cm x 260cm | drape with adhesive U-split 8,5cm x 85cm, reinfor |
| 1 pc. | 170cm x 180cm | drape with adhesive edge, reinforced |
| 1 pc. | 150cm x 180cm | drape with adhesive edge |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 pc. | 10cm x 50cm | OP tape |
| 1 pc. | 30cm x 60cm | elastic limb cover |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |
| | | |

 \bigotimes bulk packaging: 8 sets


÷ highly absorbent material












200x260cm



170x180cm






30x60cm



80x145cm



150x190cm

AT-S-SHLR2-S

| 5 | SHOULDER SET 2 |
|---|----------------|
|---|----------------|

| quantity | dimensions | component | | | |
|-------------------------|-------------------------------|---|--------------------------|----------------------------|---------------------------------|
| 1 pc. | 263cm x 406cm | drape with elastic fenestration 13cm x fluid collection pouch | 14cm, integrated | | |
| 4 pcs | 30cm x 30cm | absorbent towel | | | |
| 2 pcs | 10cm x 50cm | OP tape | | | |
| 1 pc. | 22cm x 75cm | elastic limb cover | | | |
| 1 pc. | 80cm x 145cm 150cm x 190cm | Mayo stand cover, reinforced - red instrument table cover, reinforced | | | |
| 1 рс. | | instrument table cover, reinforced | | | |
| | | | \frown | | |
| bulk packagir | | | 5years | | |
| 8 sets | | | SHELF UFE | | U |
| | | | | | |
| | | | | | |
| | | | | | 263x406 cm |
| | | | | | |
| | 4x | 2x m 10x50cm 22x75cm | 80x145cm | 150 x 190 cm | |
| | 30x300 | m 10x50cm 22x75cm | 80x145cm | 150x190cm | |
| | | | | | |
| | | | | | |
| AT-S-DH | IS1-IF-SP DYN | AMIC HIP SCREW SET 1-IF | NEW | | |
| quantity | dimensions | component | | | |
| 1 pc. | 250cm x 300cm | drape integrated with limb covers and | | | |
| 1 pc | 50cm x 60cm | fenestration 14cm x 40cm filled with an adhesive pouch, one-chamber | n incise film, reinf | orced | |
| 1 рс. 4 рсs | 30cm x 30cm | absorbent towel | | | |
| l pc. | 10cm x 50cm | adhesive tape | | | |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red | | | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | | | · · · · · · · · · · · · · · · · |
| | | | | | 250x300cm |
| \bigotimes | | | | | |
| ulk packagir 10 sets | ng: | | | | |
| | | | | | |
| STERIL | D (ALER ALER | SVEARS | | | |
| | | SHELF LIFE | | | |
| | | | | | |
| | | | | 4x | |
| | | 10 x 50 cm | 50 x 60 cm | 30 x 30 cm | 80x145cm 150x190cm |
| | | | | | |
| | | | | | |
| AT-S-V | ER1-IP-S VER | TICAL ISOLATION SET 1-IP | NEW | | |
| quantity | dimensions | component | | | |
| | | | | | |
| 1 pc. | 240cm x 340cm | drape integrated with limb covers and fenestration 14cm x 40cm filled with an | | | |
| 4 pcs | 30cm x 30cm | absorbent towel | , | | |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red | | | $\overline{\mathbf{v}}$ |
| 1 рс. | 150cm x 190cm | instrument table cover, reinforced | | | Ť |
| ~~ | | | | | 340x240cm |
| | | | | | |
| ulk packagir 12 sets | ng: foil | | | | |
| 12 Sets | | | | | |
| | D (LATEX) (RHT) | \bigcirc | | | |
| STERILE | | Svears | | | |
| | | | | | |
| | | SHEFERE | | | |
| | | SHED UPE | | | |

80x145cm

30x30cm


150x190cm

GYNAECOLOGY AND UROLOGY

| | AT-S-GYN1-S | | GYNA | ECOLOGY SET 1 |
|---|-------------|------------|-------|--|
| _ | quantity | dimensions | 8 | component |
| | 1 рс. | 160/240cm | 180cm | drape integrated with limb covers, adhesive period fenestration 10cm x 15cm, integrated fluid collecti 35cm x 50cm with formable molder, filter and tubi |
| | 1 рс. | 75cm x 900 | :m | drape without adhesive edge |
| | 4 pcs | 30cm x 300 | :m | absorbent towel |
| | 1 рс. | 10cm x 500 | m | OP tape |
| | 1 рс. | 150cm x 19 | 90cm | instrument table cover, reinforced |








fluid collection pouch with filter and tubing port

| AT-S-GYN2-C GYNA | | GYNA | AECOLOGY SET 2 OBESTSELLER | |
|------------------|----------|------------|----------------------------|---|
| _ | quantity | dimensions | \$ | component |
| | 1 pc. | 160/240cm | n x 180cm | drape integrated with limb covers, adhesive perinea 10cm x 15cm, integrated fluid collection pouch 35cr with formable molder, filter and tubing port, reinfor |
| | 1 pc. | 75cm x 900 | m | drape without adhesive edge |
| | 4 pcs | 30cm x 300 | m | absorbent towel |
| | 1 рс. | 10cm x 500 | m | OP tape |
| | 1 рс. | 80cm x 145 | 5cm | Mayo stand cover, reinforced - red |
| | 1 pc. | 150cm x 19 | 90cm | instrument table cover, reinforced |


 $\widehat{}$ \bigotimes bulk packaging: 6 sets hydrophobic SMS material















160/240x180cm










150x190cm

10x**50**cm

80x145cm


AT-S-GYN3-C GYNAECOLOGY SET 3

| quantity | dimensions | component | | |
|----------|---------------|--------------------------------------|------------------------------|-----------|
| 1 pc. | 170cm x 180cm | drape with adhesive edge, reinforced | | |
| 2 pcs | 100cm x 150cm | drape with adhesive edge | | |
| 1 pc. | 75cm x 90cm | drape with adhesive edge | · | |
| 1 pc. | 100cm x 150cm | drape without adhesive edge | | |
| 4 pcs | 30cm x 30cm | absorbent towel | · | |
| 1 pc. | 10cm x 50cm | OP tape | | |
| 2 pcs | 70cm x 120cm | limb cover | | 2× |
| 1 pc. | 80cm x 140cm | Mayo stand cover, reinforced - red | - 170 x 180 cm | 100x150cr |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | 170X100Cm | 100/1300 |
| - 1 pc. | 150cm x 190cm | | | |






 \bigotimes bulk packaging: 20 sets




AT-S-UROG2-S URO/GYNAECOLOGY SET 2

| quantity | dimensions | component |
|----------|---------------|---|
| 1 pc. | 160cm x 180cm | drape with adhesive edge |
| 1 pc. | 75cm x 100cm | drape with adhesive fenestration 6cm x 14cm |
| 1 pc. | 10cm x 50cm | OP tape |
| 2 pcs | 70cm x 120cm | limb cover |
| 1 pc. | 80cm x 140cm | Mayo stand cover, reinforced - red |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced |







AT-S-UROG3-S URO/GYNAECOLOGY SET 3

| quantity | dimensions | component |
|----------|---------------|--|
| 1 pc. | 75cm x 200cm | drape with adhesive diamond-shaped fenestration 8cm x 12cm |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 2 pcs | 70cm x 120cm | limb cover |
| 1 рс. | 100cm x 150cm | instrument table cover, reinforced |
| | | |



| • |
|---|


Syears


| AT-S-TUR1-IP-C TUR SET 1-IP | | | | |
|-----------------------------|---------------|--|--|--|
| quantity | dimensions | component | | |
| 1 pc. | 190cm x 230cm | drape integrated with limb covers, adhesive abdor fenestration ø 6cm and adhesive perineal fenestra ø 6cm, finger coat, fluid collection pouch with form molder, filter and tubing port | | |
| 4 pcs | 30cm x 30cm | absorbent towel | | |
| 1 pc. | 10cm x 50cm | OP tape | | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | | |




AT-S-TUR2-IP-C BESTSELLER TUR SET 2-IP quantity dimensions component drape with integrated leggings, abdominal fene-175cm x 240cm stration, perineal fenestration with finger coat and fluid collection pouch 30cm x 30cm absorbent towel 2cm x 22cm surgical tape, Velcro type 1 pc. 150cm x 190cm instrument table cover, reinforced










AT-S-PCNL1-IP-S PERCUTANEOUS LITHOTRIPSY SET 1-IP



10x50cm

NEW

AT-S-CSEC1-IF-C C-SECTION SET 1-IF

hydrophobic

SMS material

AT-S-CSEC2-IF-C C-SECTION SET 2-IF

200cm x 320cm

30cm x 30cm

10cm x 50cm

2cm x 22cm

75cm x 80cm

40cm x 40cm

30cm x 30cm

80cm x 140cm 150cm x 190cm

XL

component

absorbent towel

surgical tape, Velcro type

surgical gown STANDARD

surgical gown STANDARD

 \mathbf{O}

fluid collection

pouch 360º

absorbent baby towel

OP tape

gauze ball

quantity dimensions

| quantity | dimensions | component |
|----------|---------------|---|
| 1 pc. | 200cm x 320cm | drape with adhesive fenestration 25cm x 30cm fille an incise film, integrated with fluid collection pouch 75cm x 85cm with tubing ports |
| 4 pcs | 30cm x 30cm | absorbent towel |
| 1 рс. | 10cm x 50cm | OP tape |
| 1 рс. | 75cm x 80cm | absorbent baby towel |
| 1 рс. | 80cm x 145cm | Mayo stand cover, reinforced - red |
| 1 рс. | 150cm x 190cm | instrument table cover, reinforced |
| | <u></u> | Ð |

fluid collection

pouch 360º



bulk packaging:

10 sets





AT-S-CSEC1-C C-SECTION SET 1

| quantity | dimensions | component | | |
|----------|---------------|--|--|--|
| 1 pc. | 200cm x 320cm | drape with adhesive fenestration 25cm x 30cm, integrated with fluid collection pouch 75cm x 85cm with tubing ports | | |
| 4 pcs | 30cm x 30cm | absorbent towel | | |
| 1 pc. | 10cm x 50cm | OP tape | | |
| 1 pc. | 75cm x 80cm | absorbent baby towel | | |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red | | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | | |









 \bigotimes bulk packaging: 6 sets







ATEX Not STERILE R

10x50cm


10x50cm





AT-S-CSEC3-IF-SP C-SECTION SET 3-IF


| quantity | dimensions | component | | |
|----------|-------------------|---|--|--|
| 1 рс. | 180/240cm x 215cm | drape integrated with limb covers, adhesive fene- stration 23,5cm x 23,5cm x 22,5cm x 33cm filled with an incise film, integrated with fluid collection pouch 56cm x 80cm with tubing port, reinforced | | |
| 4 pcs | 30cm x 30cm | absorbent towel | | |
| 1 pc. | 10cm x 50cm | OP tape | | |
| 1 pc. | 75cm x 80cm | absorbent baby towel | | |
| 1 pc. | 80cm x 140cm | Mayo stand cover, reinforced - red | | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | | |



 \bigotimes ÷ bulk packaging: 10 sets highly absorbent material








180/240x215cm

AT-S-DEL1-S DELIVERY SET 1


dimension

| quantity | dimensions | component |
|----------|--------------|-----------------------------|
| 1 рс. | 90cm x 120cm | drape without adhesive edge |
| 2 pcs | 30cm x 30cm | absorbent towel |
| 3 pcs | 75cm x 80cm | absorbent baby towel |
| 1 pc. | 60cm x 90cm | absorbent underpad |
| 1 pc. | 80cm x 130cm | foil apron |
| 2 pcs | - | umbilical cord clamp |






75x80cm 60x90cm



AT-S-DEL2-S DELIVERY SET 2 quantity dimensions component 90cm x 120cm drape without adhesive edge 30cm x 30cm 6 pcs absorbent towel 1 pc. 75cm x 80cm absorbent baby towel 2 pcs 60cm x 90cm absorbent underpad

RHA bulk packaging: 20 sets



60x90cm

75x80cm

| | quantity | dimensions | component |
|---|----------|--------------|-----------------------------|
| | 1 рс. | 90cm x 120cm | drape without adhesive edge |
| | 2 pcs | 30cm x 30cm | absorbent towel |
| | 3 pcs | 75cm x 80cm | absorbent baby towel |
| | 1 pc. | 60cm x 90cm | absorbent underpad |
| | 1 pc. | 80cm x 130cm | foil apron |
| | 2 pcs | - | umbilical cord clamp |
| | 10 pcs | 30cm x 30cm | gauze ball |
| bulk packaging: 20 sets STERILE R STERILE R STERILE R | | | |

AT-S-DEL3-S DELIVERY SET 3

| AT-S-DEL4-S DELIVERY SET 4 | | | |
|----------------------------|--------------|-----------------------------|--|
| quantity | dimensions | component | |
| 2 pcs | 90cm x 120cm | drape without adhesive edge | |
| 6 pcs | 30cm x 30cm | absorbent towel | |
| 1 pc. | 75cm x 80cm | absorbent baby towel | |
| 1 pc. | 60cm x 90cm | absorbent underpad | |
| 1 pc. | 80cm x 130cm | foil apron | |
| 2 pcs | - | umbilical cord clamp | |





| AT-S-DEL5-SP | DELIVERY SET 5 |
|--------------|-----------------------|
|--------------|-----------------------|

| quantity | dimensions | component | | |
|----------|---------------|--|--|--|
| 1 pc. | 114cm x 150cm | delivery drape integrated with fluid collection pouch 78cm x 80cm with formable molder, filter and tubing port, reinforced | | |
| 2 pcs | 30cm x 30cm | absorbent towel | | |
| 1 pc. | 75cm x 80cm | absorbent baby towel | | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | | |








pouch with a very large scale



188 www.zarys.com










| AT-S-O | OPHT2-IF-C O | PHTHALMIC SET 2-IF | | | | | |
|--------------------------|---------------|--|---------|---------|----------|-----------------|---------|
| quantity | dimensions | component | | _ | | | |
| 1 pc. | 120cm x 170cm | ophthalmic drape with fenestration 10cm x 1 an incise film, integrated with 2 fluid collectic and with 1 formable molder, reinforced | | | | | |
| 2 pcs | 30cm x 30cm | absorbent towel | | _ | | | 2x |
| 1 pc. | 40cm x 65cm | Phaco tray cover | | | 100 | 470 | |
| 2 pcs | 35cm x 75cm | armboard cover | | | 120> | (170 cm | 30x30cm |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red | | | | | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | | | | | |
| bulk packagir 12 sets | SMS materia | | Ţ | 2x | | | |
| SIERILE | | SHELFLIFE | 40x65cm | 35x75cm | 80x145cm | 150x190c | m |

30x30cm

40x65cm

AT-S-OPHT3-IF-C OPHTHALMIC SET 3-IF

| quantity | dimensions | component | | |
|----------|---------------|--|-----------|--|
| 1 pc. | 140cm x 240cm | ophthalmic drape with fenestration 10cm x 10cm filled with an incise film, integrated with 2 fluid collection pouches and with 1 formable molder, reinforced | | |
| 4 pcs | 30cm x 30cm | absorbent towel | | |
| 1 pc. | 40cm x 65cm | Phaco tray cover | | |
| 2 pcs | 35cm x 75cm | armboard cover | | |
| 1 pc. | L | surgical gown STANDARD | | |
| 1 pc. | XL | surgical gown STANDARD | | |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced - red | 140x240cm | |
| 1 pc. | 150cm x 190cm | instrument table cover, reinforced | | |



NEUROSURGERY AND OTOLARYNGOLOGY

AT-S-OTOL1-C OTOLARYNGOLOGY SET 1



| AT-S-OTOL2-C | OTOLARYNGOLOGY SET 2 |
|--------------|----------------------|
|--------------|----------------------|


| quantity | ity dimensions | | component | | | |
|-------------------------|----------------|-----------------------------|---|--|--|--|
| 1 pc. | 150cm x 200cm | | drape with adhesive U-split 6cm x 65cm, reinforce | | | |
| 1 pc. | 150 |)cm x 180cm | drape with adhesive edge | | | |
| 1 pc. | 750 | m x 90cm | drape with adhesive edge | | | |
| 4 pcs | 300 | m x 30cm | absorbent towel | | | |
| 1 pc. | 100 | m x 50cm | OP tape | | | |
| 1 рс. | 800 | m x 145cm | Mayo stand cover, reinforced - red | | | |
| 1 рс. | 150 |)cm x 190cm | instrument table cover, reinforced | | | |
| bulk packagi 12 sets | - | hydrophobic SMS material | | | | |

75x**90**cm


AT-S-CRA1-IF-SP CRANIOTOMY SET 1-IF

| quantity | dimensions | component |
|----------|---------------|---|
| 1 рс. | 230cm x 300cm | drape with adhesive fenestration 20cm x 30cm fille with anincise film, with fluid collection pouch, reinfi |
| 4 pcs | 50cm x 50cm | drape with adhesive edge |
| 2 pcs | 30cm x 30cm | absorbent towel |
| 1 рс. | 10cm x 50cm | OP tape |
| 1 рс. | 2cm x 25cm | surgical tape, velcro type |
| 1 pc. | 80cm x 145cm | Mayo stand cover, reinforced |
| 1 рс. | 150cm x 190cm | instrument table cover, reinforced |


| bulk packaging: 16 sets | head surgery | highly absorbent material | | |
|----------------------------|--------------|------------------------------|--------|------------|
| STERILE R | SPEARS | i I | 2x25cm | 2 x |



150x190cm



230x300cm




50x**50**cm

10x50cm

30x30tcm

80x145cm

150x190cm www.zarys.com 191



ZARYS


192 www.zarys.com

150x190cm



SURGICAL GOWNS

STERILE GOWN NORMAL









STERILE SURGICAL GOWN STANDARD

| AT-SGS-M 1E | M | 1 pc. |
|---------------|-----|-------|
| | | . 1 |
| AT-SGS-L 1E | L | 1 pc. |
| AT-SGS-XL 1E | XL | 1 pc. |
| AT-SGS-XXL 1E | (XL | 1 pc. |



STERILE SURGICAL GOWN STANDARD

| REF | size | trimming | packaging unit |
|--------------|------|----------|----------------|
| AT-SGS-M 2 | М | | 1 pc. |
| AT-SGS-L 2 | L | | 1 pc. |
| AT-SGS-XL 2 | XL | | 1 pc. |
| AT-SGS-XXL 2 | XXL | | 1 pc. |

NEW





194 www.zarys.com





STERILE SURGICAL GOWN STANDARD

| REF | size | trimming | packaging unit |
|---------------------------|------------------------------------|----------|-----------------------|
| AT-SGS-M | М | | 1 pc. |
| AT-SGS-L | L | | 1 pc. |
| AT-SGS-XL | XL | | 1 pc. |
| AT-SGS-XXL | XXL | | 1 pc. |
| bulk packaging: 25 pcs | include 2 absorbent 30cm x 3 | towels | لَمْ الْحَ 45 g/m² |
| STERILE EO | | 5YEARS | |

| REF | size | trimming | packaging unit | |
|----------------|------|----------|---|--|
| AT-SGSP-M 1 | М | | 1 pc. | |
| AT-SGSP-L 1 | L | | 1 pc. | |
| AT-SGSP-XL 1 | XL | | 1 pc. | |
| AT-SGSP-XXL 1 | XXL | | 1 pc. | |
| AT-SGSP-XXXL 1 | XXXL | | 1 pc. | |
| | | 7 | ~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~ | |



STERILE SURGICAL GOWN STANDARD PLUS

30cm x 30cm

| REF | size | trimming | packaging uni | t |
|---------------------------|--|----------|---------------------|------------------|
| AT-SGSP-M 1E | М | | 1 pc. | |
| AT-SGSP-L 1E | L L | | 1 pc. | |
| AT-SGSP-XL 1E | XL | | 1 pc. | |
| AT-SGSP-XXL 1E | XXL | | 1 pc. | |
| bulk packaging: 25 pcs | without 2 absorbent to 30cm x 30cr | | economic version | کی کے 35 g/m² |
| STERILEEO | \times \times \setminus | | | |










STERILE SURGICAL GOWN STANDARD PLUS

| REF | size | trimming | packaging unit | | | |
|---------------------------|------------------------------|-----------|------------------|--|--|--|
| AT-SGSP-M 2 | М | | 1 рс. | | | |
| AT-SGSP-L 2 | L | | 1 pc. | | | |
| AT-SGSP-XL 2 | XL | | 1 pc. | | | |
| AT-SGSP-XXL 2 | XXL | | 1 pc. | | | |
| bulk packaging: 25 pcs | inclu 2 absorbe 30cm > | nt towels | کے کے 40 g/m² | | | |
| STERILEEO | STERILE EO | | | | | |



STERILE SURGICAL GOWN STANDARD PLUS

| REF | size | trimming | packaging unit |
|-------------|------|----------|----------------|
| AT-SGSP-M | М | | 1 pc. |
| AT-SGSP-L | L | | 1 pc. |
| AT-SGSP-XL | XL | | 1 pc. |
| AT-SGSP-XXL | XXL | | 1 pc. |
| | _ | | |






30x30cm



STERILE SURGICAL GOWN COMFORT

| REF | size | trimming | packaging unit | | | |
|---------------------------|------------------------------|-----------|----------------|--|--|--|
| AT-SGC-M | М | | 1 рс. | | | |
| AT-SGC-L | L | | 1 pc. | | | |
| AT-SGC-XL | XL | | 1 pc. | | | |
| AT-SGC-XXL | XXL | | 1 pc. | | | |
| bulk packaging: 25 pcs | inclu 2 absorbe 30cm x | nt towels | ۿؚ 68 g/m² | | | |
| STERILE EO | | | | | | |





STERILE SURGICAL GOWN COMFORT PLUS

| REF | size | trimming | packaging unit | | | |
|---------------------------|----------|------------------------------|-----------------------------|--|--|--|
| AT-SGCP-M | М | | 1 pc. | | | |
| AT-SGCP-L | L | | 1 pc. | | | |
| AT-SGCP-XL | XL | | 1 pc. | | | |
| AT-SGCP-XXL | XXL | | 1 pc. | | | |
| bulk packaging: 25 pcs | 2 absorb | udes ent towels x 30cm | کے آ 68 g/m ² | | | |
| STERILEEO | | | | | | |

NON-STERILE SURGICAL GOWNS

NON-STERILE SURGICAL GOWN STANDARD


| size | trimming | packaging unit |
|------|--------------|----------------|
| М | | 1 pc. |
| L | | 1 рс. |
| XL | | 1 рс. |
| XXL | | 1 рс. |
| | M L XL | |



NON-STERILE SURGICAL GOWN PREMIUM

| REF | size | trimming | packaging unit | |
|--------------------------------|--------------------------|----------|---------------------------------|-------------------|
| AT-NS-SGP-L 1E | L | | 10 pcs | |
| AT-NS-SGP-XL 1E | XL | | 10 pcs | |
| bulk packaging: 10 x 10 pcs | fully reinforce PP+PE | | et the standards of EN 14126 | کے آگے 35 g/m² |
| NON CATER A | SYEARS | | | |





It is extremely important for us to offer our clients products that ensure high standards of working conditions and comfort of use.

In order to achieve this goal, we have created a line of high quality products made of different types of non-woven

- polypropylene non-woven

- viscose non-woven
- Spunlace Mesh non-woven

A combination of several factors such as optimal design, perfect fit, breathability and structure makes BETAtex products comfortable and certainly meet our clients' individual comfortand safety for the Medical Personnel and the Patients.



BETAtex medical clothing

medical gowns I surgical clothing I caps and protective masks I patient clothing I shoe covers I bedclothes and absorbent underpads

ZARYS

MEDICAL

GOWNS

MEDICAL GOWN

with elastic cuffs, non-woven

features

made of polypropylene non-woven 20 g/m², 25 g/m², 30 g/m² or 40 g/m²

• ties at waist and neck · long sleeves with elastic cuffs

breathable

available sizes: L. XL

single-use

| REF | grammage | size | colour | packaging unit |
|-----------------|---------------------|------|--------|-------------------|
| BT-019-20B3-L | 20 g/m² | L | blue | 10 pcs |
| BT-019-20B4-XL | 20 g/m² | XL | blue | 10 pcs |
| BT-019-20G3-L | 20 g/m² | L | green | 10 pcs |
| BT-019-20G4-XL | 20 g/m ² | XL | green | 10 pcs |
| *BT-019-20W3-L | 20 g/m² | L | white | 10 pcs |
| BT-019-25B3-L | 25 g/m² | L | blue | 10 pcs |
| BT-019-25B4-XL | 25 g/m² | XL | blue | 10 pcs |
| BT-019-25G3-L | 25 g/m² | L | green | 10 pcs |
| BT-019-25G4-XL | 25 g/m² | XL | green | 10 pcs |
| *BT-019-25W3-L | 25 g/m² | L | white | 10 pcs |
| *BT-019-30W3-L | 30 g/m² | L | white | 10 pcs |
| *BT-019-30G3-L | 30 g/m ² | L | green | 10 pcs |
| *BT-019-30G4-XL | 30 g/m ² | XL | green | 10 pcs |
| *BT-019-30B3-L | 30 g/m² | L | blue | 10 pcs |
| *BT-019-40B3-L | 40 g/m ² | L | blue | 10 pcs |
| *BT-019-40G3-L | 40 g/m ² | L | green | 10 pcs |
| | | | | |

bulk packaging: 10 x 10 pcs



*

available

on request

MEDICAL GOWN

with knitted cuffs, non-woven

features

- made of polypropylene non-woven 20 g/m², 25 g/m², 30 g/m² or 40 g/m²
- ties at waist and neck
- · long sleeves with knitted cuffs
- breathable
- available sizes: L, XL, XXL
- single-use

| REF | grammage | size | colour | packaging unit |
|-----------------|---------------------|------|--------|-------------------|
| BT-015-20G3-L | 20 g/m ² | L | green | 10 pcs |
| BT-015-20G4-XL | 20 g/m ² | XL | green | 10 pcs |
| BT-015-25G3-L | 25 g/m² | L | green | 10 pcs |
| BT-015-25B3-L | 25 g/m² | L | blue | 10 pcs |
| BT-015-25G4-XL | 25 g/m² | XL | green | 10 pcs |
| BT-015-30G3-L | 30 g/m² | L | green | 10 pcs |
| BT-015-30G4-XL | 30 g/m² | XL | green | 10 pcs |
| BT-015-30G5-XXL | 30 g/m² | XXL | green | 10 pcs |
| BT-015-40G3-L | 40 g/m ² | L | green | 10 pcs |
| BT-015-40B3-L | 40 g/m ² | L | blue | 10 pcs |
| BT-015-40G4-XL | 40 g/m ² | XL | green | 10 pcs |

8

available colours:

bulk packaging 10 x 10 pcs









WWW.ZARYS.COM

SURGICAL SUIT

PP short-sleeved shirt+trousers, non-woven

features

• made of SMS non-woven 35 g/m²

- short-sleeved shirt with V-neck and 3 pockets
- trousers with guiding ties
- blue shirt and trousers available also in individual packaging

single-use

| SET (SHIRT + TROUSERS) | | | | |
|------------------------|-------|--------|-------------------|--|
| REF | size | colour | packaging unit | |
| BT-017-S35B1-S-SMS | S | blue | 1 set | |
| BT-017-S35B2-M-SMS | Μ | blue | 1 set | |
| BT-017-S35B3-L-SMS | L | blue | 1 set | |
| BT-017-S35B4-XL-SMS | XL | blue | 1 set | |
| BT-017-S35B5-XXL-SMS | XXL | blue | 1 set | |
| BT-017-S35B6-XXXL-SMS | XXXL | blue | 1 set | |
| BT-017-S35B7-XXXXL-SMS | XXXXL | blue | 1 set | |
| BT-017-S35G1-S-SMS | S | green | 1 set | |
| BT-017-S35G2-M-SMS | М | green | 1 set | |
| BT-017-S35G3-L-SMS | L | green | 1 set | |
| BT-017-S35G4-XL-SMS | XL | green | 1 set | |
| BT-017-S35G5-XXL-SMS | XXL | green | 1 set | |
| BT-017-S35V1-S-SMS | S | purple | 1 set | |
| *BT-017-S35V2-M-SMS | Μ | purple | 1 set | |
| *BT-017-S35V3-L-SMS | L | purple | 1 set | |
| *BT-017-S35V4-XL-SMS | XL | purple | 1 set | |
| *BT-017-S35V5-XXL-SMS | XXL | purple | 1 set | |

| SHIRT | | | |
|-----------------------|------|--------|-------------------|
| REF | size | colour | packaging unit |
| *BT-018-S35B1-S-SMS | S | blue | 1 pc. |
| *BT-018-S35B2-M-SMS | Μ | blue | 1 pc. |
| *BT-018-S35B3-L-SMS | L | blue | 1 pc. |
| *BT-018-S35B4-XL-SMS | XL | blue | 1 pc. |
| *BT-018-S35B5-XXL-SMS | XXL | blue | 1 pc. |

| TROUSERS | | | |
|-----------------------|------|--------|-------------------|
| REF | size | colour | packaging unit |
| *BT-018-T35B1-S-SMS | S | blue | 1 pc. |
| *BT-018-T35B2-M-SMS | Μ | blue | 1 pc. |
| *BT-018-T35B3-L-SMS | L | blue | 1 pc. |
| *BT-018-T35B4-XL-SMS | XL | blue | 1 pc. |
| *BT-018-T35B5-XXL-SMS | XXL | blue | 1 pc. |





SURGICAL CLOTHING



SURGICAL SHIRT SHORT-SLEEVED

non-woven

features

• made of soft polypropylene non-woven 47 g/m²

- short sleeves
- 3 big pockets

 \bigotimes

bulk packaging: 50 x 1 pc.

• snap fastener at the neck single-use

| F | REF | size | packaging unit |
|---|--------------------------|-------|-------------------|
| E | 3T-018-S47B0-XS-SOFT2 | XS | 1 pc. |
| E | BT-018-S47B1-S-SOFT2 | S | 1 pc. |
| E | BT-018-S47B2-M-SOFT2 | Μ | 1 pc. |
| E | BT-018-S47B3-L-SOFT2 | L | 1 pc. |
| E | BT-018-S47B4-XL-SOFT2 | XL | 1 pc. |
| E | BT-018-S47B5-XXL-SOFT2 | XXL | 1 pc. |
| E | BT-018-S47B6-XXXL-SOFT2 | XXXL | 1 pc. |
| E | BT-018-S47B7-XXXXL-SOFT2 | XXXXL | 1 pc. |
| | | | |

available colours:



SURGICAL TROUSERS, WITH TIES

non-woven

features

- made of soft polypropylene non-woven 47 g/m^2
- trousers with guiding ties
- 1 big pocket with snap fastener
- single-use

| REF | size | packaging unit |
|--------------------------|-------|-------------------|
| BT-018-T47B0-XS-SOFT1 | XS | 1 pc. |
| BT-018-T47B1-S-SOFT1 | S | 1 pc. |
| BT-018-T47B2-M-SOFT1 | Μ | 1 pc. |
| BT-018-T47B3-L-SOFT1 | L | 1 pc. |
| BT-018-T47B4-XL-SOFT1 | XL | 1 pc. |
| BT-018-T47B5-XXL-SOFT1 | XXL | 1 pc. |
| BT-018-T47B6-XXXL-SOFT1 | XXXL | 1 pc. |
| BT-018-T47B7-XXXXL-SOFT1 | XXXXL | 1 pc. |
| | | |



SURGICAL SHIRT LONG-SLEEVED

non-woven

features

made of soft polypropylene non-woven 47 g/m²
long sleeves with knitted cuffs
2 big pockets

NON STERILE STERILE SHEP IFE

snap fasteners at the front
single-use

| REF | size | packaging unit |
|----------------------|------|-------------------|
| BT-018-047B1-S-SOFT3 | S | 1 pc. |
| BT-018-047B2-M-SOFT3 | Μ | 1 pc. |

| BT-018-047B2-M-SOFT3 | M | 1 pc. | |
|------------------------|-----|-------|--|
| BT-018-047B3-L-SOFT3 | L | 1 pc. | |
| BT-018-047B4-XL-SOFT3 | XL | 1 pc. | |
| BT-018-O47B5-XXL-SOFT3 | XXL | 1 pc. | |
| | | | |







ZARYS

MEDICAL BOUFFANT CAP

PP non-woven

features

- · cap with an elastic rubber band
- diameter: approx. 53 cm
- made of airy polypropylene non-woven
- excellent breathability and wearing comfort • packaging: ergonomic cardboard box to remove individual pieces
- single-use



DONA medical bouffant cap

features

- · cap with an elastic rubber band
- diameter: approx. 47 cm
- made of perforated viscose non-woven
- excellent breathability and wearing comfort
- · packaging: ergonomic cardboard box to remove individual pieces single-use

| REF | colour | packaging unit |
|--------------------|--------|-------------------|
| BT-008-25G-B1-V | green | 100 pcs |
| BT-008-25B-B1-V | blue | 100 pcs |
| BT-008-25G-B1-V-XL | green | 100 pcs |

÷

bulk packaging: 4 x 100 pcs





LILI

medical cap with rubber


features

• cap with an elastic rubber band in the back part

made of perforated viscose non-woven

• excellent breathability and wearing comfort

- packaging: ergonomic cardboard box to remove individual pieces
- single-use



CAPS, PROTECTIVE MASKS



POLA medical cap with rubber

features

- cap with an elastic rubber band in the back part
- made of durable Spunlace and polypropylene non-woven
- the used non-woven is characterized by good absorption properties and durability and additionally it is pleasant to the touch, which increases comfort while using
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces single-use

REF packaging unit BT-001-45B-SPN 100 pcs NONE TATEX STRATES available colours: bulk packaging: 6 x 100 pcs



medical cap with ties

features

- cap with ties at the back
- made of perforated viscose non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces single-use

| REF | colour | packaging unit | |
|--------------------------------|--------------------|-------------------|--------|
| BT-007-25G-T1-V | green | 100 pcs | |
| BT-007-25B-T1-V | blue | 100 pcs | |
| bulk packaging: 6 x 100 pcs | available colours: | NON STERLE | SYEARS |



ABI

medical cap with ties

features

- deep cap with ties at the back
- being folded out, the front part creates an additional layer protecting against sweat
- made of perforated viscose non-woven
- excellent breathability and wearing comfort • packaging: ergonomic cardboard box to remove individual pieces
- single-use

| REF | colour | packaging unit |
|-----------------|--------|-------------------|
| BT-007-25G-T2-V | green | 100 pcs |
| BT-007-25B-T2-V | blue | 100 pcs |





NINO

medical cap

features

- cap without rubber or ties, designed for short-haired people
- two layers of non-woven fabric in the side part
- made of perforated viscose non-woven excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use

| REF | colour | packaging unit |
|--------------------------------|--------------------|-------------------|
| BT-001-25G-D-V | green | 100 pcs |
| BT-001-25B-D-V | blue | 100 pcs |
| bulk packaging: 6 x 100 pcs | available colours: | STERIES STREET (|

KOKO long

medical cap with snood

features

- cap with snood and elastic rubber band in the back part
- optimal protection for long-haired people
- made of perforated viscose non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use

| REF | colour | packaging unit |
|--------------------------------|---------------|--------------------|
| BT-001-25G-S1-V | green | 100 pcs |
| BT-001-25B-S1-V | blue | 100 pcs |
| bulk packaging: 6 x 100 pcs | for long hair | available colours: |


DORA long

medical cap with snood

features

- cap with anti-sweat insert in the forehead area
- special design allows the cap to adjust to the shape of the head
- and keep one's hair in an excellent manner
- optimal protection for long-haired people
- packaging: ergonomic cardboard box to remove individual pieces
- single-use













- made of Spunlace and polypropylene non-woven
- excellent breathability and wearing comfort



BEN medical cap with ties

features

• forage cap with ties at the back • made of perforated Spunlace Mesh non-woven • excellent breathability and comfort of use single-use REF packaging unit BT-007-35B-SPNM 100 pcs



<u>***</u>

anti-sweat

insert

anti-sweat

available colours:

TOM

medical cap with ties

features

- cap with ties at the back
- anti-sweat insert on the whole circumference • made of polypropylene non-woven 20 g/m²
- packaging: ergonomic cardboard box to remove
- individual pieces
- single-use

| REF | packaging unit |
|-------------|-------------------|
| BT-007-P20G | 100 pcs |



 \bigotimes

bulk packaging:

10 x 100 pcs





features • cap with ties at the back • anti-sweat insert in the forehead area • made of Spunlace and polypropylene non-woven • packaging: ergonomic cardboard box to remove individual pieces • single-use REF packaginę unit

50 pcs



Ø



bulk packaging:

6 x 100 pcs





BT-007-P45B


medical cap with ties


features

- made of polypropylene non-woven 30 g/m²
- forage cap
- cap with ties at the back
- packaging: ergonomic cardboard box to remove
- individual pieces
- single-use

| REF | packaging unit | |
|------------------|-------------------|----------------|
| BT-007-30C-T1-PP | 100 pcs | NON STERILE |







| REF | colour | packaging unit | |
|------------|--------|-------------------|-----------|
| BT-046-32G | green | 50 pcs | A GATER (|
| BT-046-32B | blue | 50 pcs | |
| | | | |

| REF | colour | type in accordance with the standard EN 14683 | packaging unit |
|-------------|--------|---|-------------------|
| BT-046-31G | green | type II | 50 pcs |
| BT-046-31B | blue | type II | 50 pcs |
| BT-046-R31B | blue | type II R | 50 pcs |



| MEDICA with earloo | | NEW | | |
|---|---------|---|-------------------|----|
| features made of three I modeling nose | -piece | | three laye | rs |
| bacteria filtratii packaging: ergo pieces single-use | , | FE ≥ 98% ard box to remove individual | NON | |
| REF | colour | type in accordance with the standard EN 14683 | packaging unit | |
| BT-046-31C | pattern | type II | 50 pcs | |



WWW.ZARYS.COM

PATIENT'S GOWN

tied-up, non-woven

features

- made of SMS non-woven
- short sleeves
- white ties at waist and neck
- easy putting on and taking off in bed-ridden patients
- universal size
- single-use

| REF | material | size | grammage | colour | packaging unit |
|--------------------------------|--------------------|----------------|------------|--------|-------------------|
| BT-0022-33B-SMS | SMS non-woven | universal | 33 g/m² | blue | 10 pcs |
| bulk packaging: 10 x 10 pcs | available colours: | NON STERILE | PER SYLARS | • | |

PATIENT'S GOWN

Y-neck, non-woven

features

- made of polypropylene non-woven
- short sleeves and Y-neck
- slipover
- single-use



PATIENT'S GOWN

V-neck, non-woven

features

• made of Spunlace non-woven

- the used non-woven is characterized by good absorption properties and durability and additionally it is pleasant to the touch, which increases
- , comfort while wearing
- short sleeves and V-neck
- slipover
- single-use

| REF | material | size | grammage | colour | packaging unit | |
|-----------------|--------------------|-----------|----------|--------|-------------------|--------------|
| BT-0022-45W-SPN | Spunlace non-woven | universal | 45 g/m² | white | 10 pcs | pleasant for |
| BT-0022-45B-SPN | Spunlace non-woven | universal | 45 g/m² | blue | 10 pcs | the skin |

bulk packaging: 10 x 10 pcs





PATIENT CLOTHING



| olour | packaging unit |
|----------|-------------------|
| ark blue | 10 pcs |



PATIENT'S GOWN

for delivery, SMS non-woven

features

• made of SMS non-woven

• short sleeves, ties at waist and neck (other models available on request)

• available in 5 sizes: S, M, L, XL, XXL

single-use

 \bigotimes

bulk packaging: 10 x 10 pcs

| REF | material | size | grammage | packaging unit |
|----------------------------|---------------|------|----------|-------------------|
| BT-0022-35B-SMS-DM-001-S | SMS non-woven | S | 35 g/m² | 10 pcs |
| BT-0022-35B-SMS-DM-001-M | SMS non-woven | Μ | 35 g/m² | 10 pcs |
| BT-0022-35B-SMS-DM-001-L | SMS non-woven | L | 35 g/m² | 10 pcs |
| BT-0022-35B-SMS-DM-001-XL | SMS non-woven | XL | 35 g/m² | 10 pcs |
| BT-0022-35B-SMS-DM-001-XXL | SMS non-woven | XXL | 35 g/m² | 10 pcs |

NEW



NON STERLE VATER STERLE

PATIENT'S GOWN, FOR CHILDREN

available colours:

non-woven

features

made of SMS non-woven 35 g/m²
 white ties at waist and neck
 single-use

| REF | size | grammage | packaging unit |
|--------------------|----------------|----------|-------------------|
| BT-0022-35B-SMS-D1 | 72 cm x 100 cm | 35 g/m² | 10 pcs |



COLONOSCOPY BOXERS

SMS non-woven

features

- made of SMS non-woven 35 g/m²
- boxers with a slit in the back part
- completely opaque
- single-use

| REF | material | size | grammage | packaging unit |
|--------------------|---------------|-----------|----------|-------------------|
| BT-0023-35B-SMS | SMS non-woven | universal | 35 g/m² | 10 pcs |
| BT-0023-35B-SMS-XL | SMS non-woven | XL | 35 g/m² | 10 pcs |





COLONOSCOPY BOXERS

PP non-woven

features

- made of polypropylene non-woven 40 g/m²
- boxers with a slit in the back part
- single-use

| REF | material | size | grammage | packa unit |
|--------------------------------|--------------------|-----------|----------|---------------|
| BT-0023-40DB | PP non-woven | universal | 40 g/m² | 10 pc |
| bulk packaging: 10 x 10 pcs | available colours: | | | |
| | SHELFUFE | | | |

GYNECOLOGICAL SKIRT

non-woven

features

made of polypropylene non-woven 40 g/m²

• dimensions: 160 cm (at the waist after stretching)

single-use

| REF | material | size | grammage | ۲ ر |
|--------------------------------|--------------------|-----------|----------|--------|
| BT-0025-40DB | PP non-woven | universal | 40 g/m² | 1 |
| bulk packaging: 10 x 10 pcs | available colours: | | | |



| - | 1/2 |
|------|-----|
| 1 | 11 |
| 1773 | 77) |
| | 24 |
| | |
| | |

kaging

ICS



packaging unit 10 pcs



ZARYS

SHOE

COVERS

SHOE COVERS

features

- made of polyethylene or polypropylene non-woven
- 30 g/m² or 40 g/m² (for shoe covers with anti-skid layer)
- elastic rubber band at the opening
- protection against contaminating rooms where higher cleanliness standards are required
- increased durability
- single-use

| REF | material | grammage | size | colour | packaging unit | bulk packaging |
|------------|-----------------------------------|---------------------|---------------|--------|-------------------|-------------------|
| BT-010-G | non-woven | 30 g/m² | 38 cm x 17 cm | green | 100 pcs | 10 x 100 pcs |
| BT-010-B | non-woven | 30 g/m ² | 38 cm x 17 cm | blue | 100 pcs | 10 x 100 pcs |
| BT-011-B | non-woven with anti-skid layer | 40 g/m² | 38 cm x 17 cm | blue | 100 pcs | 10 x 100 pcs |
| BT-002-PEB | polyethylene | | 41 cm x 15 cm | blue | 100 pcs | 20 x 100 pcs |



NONE TATES STERILE

SHOE COVERS

calf-high

features

made of polypropylene non-woven 30 g/m²

- elastic rubber band at the opening
- protection against contaminating rooms where higher cleanliness
- standards are required
- increased durability
- single-use

| REF | material | grammage | size | packa unit |
|----------|--------------|----------|-----------------|---------------|
| BT-013-B | PP non-woven | 30 g/m² | 47 cm x 38,5 cm | 100 po |
| BT-013-G | PP non-woven | 30 g/m² | 47 cm x 38,5 cm | 100 po |

 \bigotimes bulk packaging: 10 x 100 pcs











aging pcs





ZARYS

DISPOSABLE MEDICAL BED SHEET

features

• made of polypropylene non-woven, laminated

- polypropylene non-woven or SMS non-woven
- · dimensions: according to customer's specification • unit packaging: foil, bulk packaging: cardboard box

| REF | material | grammage |
|---|------------------------|----------|
| BT-030-PP-17-"X" x "Y" - G (X - width , Y - length) | PP non-woven | 17 g/m² |
| BT-030-PP-25-"X" x "Y" - G (X - width, Y - length) | PP non-woven | 25 g/m² |
| BT-030-PP-35-"X" x "Y"- G (X - width, Y - length) | PP non-woven | 35 g/m² |
| BT-030-PF-25-"X" x "Y" - G (X - width, Y - length) | laminated PP non-woven | 25 g/m² |
| BT-030-PF-45-"X" x "Y" - G (X - width, Y - length) | laminated PP non-woven | 45 g/m² |
| BT-030-SMS-35-"X" x "Y" - B (X - width, Y - length) | SMS non-woven | 35 g/m² |

*according to customer's specification

DISPOSABLE MEDICAL PILLOW CASE

| features made of polyprop | ylene non-woven | | bulk packaging: 10 x 50 pcs |
|---|---------------------|--------------|--------------------------------|
| REF | grammage | size | packaging unit |
| BT-PP-25G1 | 25 g/m² | 60cm x 80 cm | 50 pcs |
| BT-PP-25G6 | 25 g/m ² | 80cm x 90 cm | 50 pcs |

DISPOSABLE MEDICAL QUILT COVER

| features • made of polyprop | bulk packaging: 10 x 10 pcs | | |
|-----------------------------|--------------------------------|----------------|-------------------|
| REF | grammage | size | packaging unit |
| BT-PZ-25G1 | 25 g/m ² | 160cm x 210 cm | 10 pcs |

DISPOSABLE MEDICAL BEDDING SET

| components: bed sheet: 150 cm; quilt cover: 160 cm | de of polypropylene non-woven | | bulk packaging: 25 x 1 set |
|--|-------------------------------|--------|-------------------------------|
| REF | grammage | colour | packaging unit |
| BT-ZES-25G3 | 25 g/m² | green | 1 set |
| BT-ZES-40G3 | 40 g/m ² | green | 1 set |

white

1 set

40 g/m²

BT-ZES-40W3

BEDCLOTHES, ABSORBENT UNDERPADS

















www.zarys.com 219

DISPOSABLE ABSORBENT UNDERPAD

features

· composed of five layers:

- laminate + cellulose wadding + cellulose pulp + cellulose wadding + PP nonwoven absorbent insert made of cellulose pulp
- impermeable foiled bottom
- single use

| REF | dimensions | packaging unit | |
|-----------|---------------|-------------------|--|
| BT-PC4060 | 40 cm x 60 cm | 25 pcs | |
| BT-PC6060 | 60 cm x 60 cm | 25 pcs | |
| BT-PC6090 | 60 cm x 90 cm | 25 pcs | |

 \bigotimes bulk packaging: 4 x 25 pcs

NEW



DISPOSABLE ABSORBENT UNDERPAD

for operating table

features

• quilted absorbent layer with addition of SAP (superabsorbent)

- · composed of five layers: polypropylene non-woven, cellulose wadding, cellulose pulp with superabsorbent SAP, cellulose wadding, PE foil
- impermeable foiled bottom
- single use

| REF | dimensions | packaging unit |
|-------------------------------|-----------------|----------------------|
| BT-PC100225 | 100 cm x 225 cm | 10 pcs |
| | | |
| bulk packaging: 2 x 10 pcs | NON CATER ALL | SYEARS SHELF LIFE |

ABSORBENT UNDERPAD FOR CHILDREN

features

- quilted absorbent insert with SAP (superabsorbent)
- white absorbent insert
- film layer with pattern
- single-use

| REF | dimensions | packaging unit | bulk packaging |
|------------|---------------|-------------------|-------------------|
| BT-PC6060C | 60 cm x 60 cm | 30 pcs | 10 x 30 pcs |
| BT-PC6090C | 60 cm x 90 cm | 15 pcs | 10 x 15 pcs |

features • made of Spunlace non-woven with polyester filling • ultrasonic seams over the entire surface

• can be heated up in the heating device

• maintaining patient temperature at a stable level

• single use

| REF dir | | outer material | packaging unit |
|---------------|----------------|--------------------|-------------------|
| BT-TOP-21S 11 | 10 cm x 120 cm | Spunlace non-woven | 1 pc. |
| BT-TOP-22S 11 | 10 cm x 210 cm | Spunlace non-woven | 1 pc. |

THERMAL PATIENT COVER

features

• made of polypropylene non-woven with polyester filling

- ultrasonic seams over the entire surface
- can be heated up in the heating device
- maintaining patient temperature at a stable level
- single use

| REF | dimensions | outer material | packa unit |
|------------|-----------------|-------------------------|---------------|
| BT-TOP-11P | 110 cm x 120 cm | polypropylene non-woven | 1 pc. |
| BT-TOP-12P | 110 cm x 210 cm | polypropylene non-woven | 1 pc. |
| BT-TOP-13P | 150 cm x 210 cm | polypropylene non-woven | 1 pc. |
| | | | |







DISPOSABLE MEDICAL BED SHEET, REINFORCED

features

- made of foil laminated paper
- reinforced with 8 or 48 threads
- grammage 40 g/m²
- soft, waterproof, absorbent, tough
- white with blue threads • product put on mattresses, couches, armchairs
- and other elements that are the equipment
- of health care units, ie hospitals, outpatient clinics, clinics, etc. • the products are designed to protect the patient during
- medical procedures, to prevent the transmission of contaminants and potentially pathogenic microorganisms located in medical rooms treatment rooms and equipment that furnish these rooms and that may lead to
- complications of the patient's treatment process. They are a barrier to contamination non-sterile
- single use

| REF | type | dimensions | packaging unit |
|----------|------------|----------------|-------------------|
| BT-PW101 | 8 threads | 80 cm x 210 cm | 25 pcs |
| BT-PW401 | 48 threads | 80 cm x 210 cm | 25 pcs |
| ~~ | | | |

bulk packaging: 4 x 25 pcs



THERMAL PATIENT COVER SOFT









vitahealth



vitahealth line of adhesive plaster sets are designed to provide optimum protection in every situation and to satisfy the needs of the most demanding users. They are an indispensable element for everyday life.



LIST OF ADHESIVE PLASTERS

| name | ref | number of pieces in a set | individual packaging | bulk packaging |
|---|--------|---------------------------|----------------------|------------------------------|
| ACTIVE Pro set of adhesive plasters for active people | VH-ZAC | 12 | 1 | 60 |
| AQUA DTO set of adhesive waterproof plasters | VH-ZAQ | 18 | 1 | 60 |
| UNIVERSAL set of adhesive universal plasters | VH-ZUN | 24 | 1 | 60 |
| TRANSPARENT set of adhesive transparent plasters | VH-ZTR | 24 | 1 | 60 |
| HANDS pro set of adhesive hands plasters | VH-ZHA | 12 | 1 | 60 |
| FAMILY pro set of adhesive family plasters | VH-ZFA | 28 | 1 | 60 |
| set of adhesive children's plasters - girl | VH-ZKG | 16 | 1 | 60 |
| KIDS pro B set of adhesive children's plasters - boy | VH-ZKB | 16 | 1 | 60 |
| SENSITIVE set of adhesive plasters for sensitive skin | VH-ZSE | 12 | 1 | 30 |
| S.O.S set of adhesive handbag plasters | VH-ZSO | 8 | 1 | 1000 |
| NONVIPIOST non-woven adhesive plaster for cutting | VH-NOP | 1 + scissors | 1 | 120 |
| FABRIDIOST fabric adhesive plaster for cutting | VH-FAP | 1+ scissors | 1 | 120 |
| SENSIPIOST adhesive fabric tape for cutting | VH-SEP | 1 | 1 | 12x1 (box) 30x12 (carton) |
| SOFT plast adhesive non-woven tape for cutting | VH-SOP | 1 | 1 | 12x1 (box) 30x12 (carton) |









www.vita-health.eu













ZARYS International Group sp. z o.o. sp.k.

We create good solutions

| \bigcirc | | $\bigcirc \bigcirc \bigcirc$ |
|--------------|---------------------------------------|------------------------------|
| headquarters | | ZARYS worldwide |
| Poland | | ▷ Czech Republic ▷ Po |
| | · · · · · · · · · · · · · · · · · · · | ▷ Romania |
| | · · · · · · · · · · · · · · · · · · · | ⊳ Italy En |

Diversify United Arab Emirates

Portugal

Rev. XI 30.09.2023

Please contact us if you require further details:

| Area | E-mail | Direct Line |
|----------------|------------------|-----------------|
| Worldwide | export@zarys.pl | +48 327970579 |
| Czech Republic | info@zarys.cz | +420 737954388 |
| Italy | ufficio@zarys.it | +39 334 6591844 |
| Portugal | geral@zarys.pt | +351 914183173 |
| Romania | office@zarys.ro | +40 371485051 |
| UAE | office@zarys.ae | +971529047753 |

Export Department

- Pod Borem 18, 41-808 Zabrze, Poland
- **48** 32 797 05 79

@ e-mail: export@zarys.pl



EC Certificate Directive 93/42/EEC Annex V Production Quality Assurance Medical Devices

Registration No.: DD 60139535 0001

Report No.:

26300232 017

Manufacturer: **ZARYS** International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products:

(see attachments for products and sites included) Replaces EC Certificate, Registration No.: DD 60117020 0001

Expiry Date:

2024-05-27

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

Effective Date:

Date:

2019-06-09

2019-05-27

Notified Body TÜVRheini Tifizierung

Rafal Byczkowski

TUV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.



Doc. 1/6, Rev. 0

TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

- Sterile and non-sterile cutting gauze
- Non-sterile dressing gauze
- Sterile and non-sterile gauze swabs (with or without X-ray thread)
- Sterile and non-sterile gauze lap sponges (with X-ray thread/ with X-ray chip)
- Sterile and non-sterile gauze balls (with or without X-ray thread)
- Sterile and non-sterile gauze rolls (with or without X-ray thread)
- Sterile and non-sterile non-woven swabs (with or without X-ray thread)
- Sterile paraffin gauze dressings
- Sterile three-way stopcocks
- Sterile transfusion sets for single use
- Sterile infusion sets for single use
- Sterile extension tubes for infusion pump

Notified Body



Rafal Byczkowski

Date: 2019-05-27

0/020 h 04.08 @ TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior approva



Doc. 2/6, Rev. 0

TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

- Sterile endotracheal tubes
- Sterile tracheostomy tubes
- Sterile breathing circuits
- Sterile catheter mounts
- Non-sterile anaesthetic masks
- Sterile laryngeal masks
- Sterile oxygen masks
- Sterile Multi-Vent masks
- Sterile non-rebreath masks
- Sterile nebulizer masks
- Sterile nasal oxygen cannulas
- Sterile nebulizer sets
- Sterile oxygen tubing
- Sterile suction catheters
- Sterile abdominal drains
- Sterile feeding tubes
- Sterile stomach and duodenal tubes
- Sterile urology catheters
- Sterile surgical suction sets

TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior app

- Sterile surgical suction cannulas
- Sterile syringes for single use

Notified Body



Rafal Byczkowski

Date: 2019-05-27



Doc. 3/6, Rev. 0

TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

- Sterile insulin syringes
- Sterile tuberculin syringes
- Sterile hypodermic needles
- Sterile insulin pen needles
- Sterile blood lancets
- Sterile IV cannulas
- Sterile needle free valves
- Sterile surgical gloves

TÜVRheinlan **Notified Body** Tinzierungs

Date: 2019-05-27

10/020 h 04.08 (8) TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior approval

Rafal Byczkowski



Doc. 4/6, Rev. 0

TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate **Registration No.: Report No.:**

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

For the following medical devices the scope covers only the aspects of manufacture concerned with securing and maintaining sterile conditions:

- Adhesive cannulla fixation dressings
- Adhesive wound dressings
- Eye pads
- Incise films
- Transparent film dressings
- Foam dressings
- Absorbent wound dressings
- Surgical gowns
- Surgical drapes
- Sets of surgical drapes
- Fluid collection pouches
- Nelaton catheters



Rafal Byczkowski

Date: 2019-05-27

10/020 h 04.08
TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior approval



Doc. 5/6, Rev. 0

TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

For the following medical devices the scope covers only the aspects of manufacture concerned with securing and maintaining sterile conditions:

- Vaginal speculums
- Cervical brushes
- Urine bags
- Tongue depressors
- Guedel airways
- Intubation stylets
- Endotracheal tube holders
- Suction tubes
- Withdrawal cannulas
- Alginate dressings
- Cannula stoppers
- Umbilical cord clamps



Rafal Byczkowski

Date: 2019-05-27

0/020 h 04.08 ®

TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior approva



Doc. 6/6, Rev. 0

TÜV Rheinland LGA Products GmbH Tillystraße 2, 90431 Nürnberg

Attachment to Certificate Registration No.: Report No.:

DD 60139535 0001 26300232 017

Manufacturer:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Pod Borem 18 41-808 Zabrze Polska

Products included:

ZARYS International Group Spolka z ograniczona odpowiedzialnoscia, spolka komandytowa ul. Gustawa Eiffel'a 15 44-109 Gliwice Poland

Activity: Production

Notified Body TÜVRheinla Tifizierung

Date: 2019-05-27

04.08 ®

TUV, TUEV and TUV are registered trademarks. Utilisation and application requires prior app

Rafal Byczkowski