

# Cordis **Avanti**<sup>®</sup>+

Sheath Introducer in Transradial Kit

Designed for Atraumatic Sheath  
Insertion in Transradial Procedures



# Designed for **Atraumatic** Sheath Insertion in **Transradial** Procedures

Building on a reputation for technological leadership, and developed in close collaboration with leading cardiologists, Cordis presents sheath introducers specifically designed for use with the transradial technique.

## The Cordis Advantage:

- ultra-smooth insertion
- outstanding kink resistance
- lower in-patient costs resulting from savings in staff time
- lower costs associated with potential use in the out-patient setting
- improved patient comfort
- a complete product range (sizes ranging from 4F-7F)
- reduced risk of radial artery spasm and its associated discomfort.



## Ordering information

Each Cordis AVANTI+ Transradial Kit includes an AVANTI+ Sheath Introducer with vessel dilator, mini-guidewire and 21G needle.

French size	Cannula Usable Length (cm)	Vessel Dilator	Mini-Guidewire length (cm)	Mini-Guidewire Diameter (inch)	Needle	Catalogue Number
4F	11	18	45	.021	21G (38mm)	504-614Z
4F	23	30	70	.021	21G (38mm)	504-624Z
5F	11	18	45	.021	21G (38mm)	504-615Z
5F	23	30	70	.021	21G (38mm)	504-625Z
6F	11	19.5	45	.021	21G (38mm)	504-616Z
6F	23	31.5	70	.021	21G (38mm)	504-626Z
7F	11	19.5	45	.021	21G (38mm)	504-617Z

### Important Information:

Prior to use, refer to the "Instructions for use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous product development we reserve the right to change product specifications without prior notification.

Europe:  
Tel: +32 2 352 14 11 • Fax: +32 2 352 15 90

© Johnson & Johnson Medical NV/SA • February 2006 • 2E-800-0454-6