## **JJJ DAHLHAUSEN**°



## Nebulizer Kits, funnel/funnel

Nebulizer Kit for inhalation of medicines and liquids by mask and safety tubing with funnel. The set is for single patient use only and should be disposed of after 24 hours.

Aerosol capacity of nebulizer according to EN 13544-1:2007.

Length (oxygen tube):	210 cm		
Set content:	chamber 6 ml with late connection, safety con	Mask with Lips style, nebulizer chamber 6 ml with lateral oxygen connection, safety connecting tube with funnel/funnel connector	
Sterile:	No		
· Art.no.		Application area	
0100008852 0100008853		Child	Committee of
		Adult	- Contraction of the Contraction

P. J. Dahlhausen & Co. GmbH Emil-Hoffmann-Str. 53 50996 Köln

Tel.: (02236) 39 13 - 0 Fax: (02236) 39 13 - 109