


Trocars and Accessories

Size 6 mm




30160MP



30160MP

Trocar, with pyramidal tip,
with insufflation stopcock, size 6 mm,
working length 10.5 cm,
color code: black,
including:
Cannula, without valve
Trocar only
Multifunctional Valve

Size 11 mm



30103MP


30103MP

Trocar, with pyramidal tip,
with insufflation stopcock, size 11 mm,
working length 10.5 cm,
color code: green,
including::
Cannula, without valve
Trocar only
Multifunctional Valve



30140DB

Reduction Sleeve, reusable,
instrument diameter 5 mm,
trocar cannula outer diameter 11 mm,
color code: green



30142HB

Double Reducer,
13/10 mm, 13.5/10 mm, 13/5 mm
and 13.5/5 mm

Accessories

Reducers and Reduction Sleeves

Reduction Sleeves



30140 DB

Instrument Diameter	Cannula Outer Diameter – Color Code:	
	6 mm black	11 mm green
3 mm	30140 KA	30140 DA
5 mm	–	30140 DB

Reducers



62141 AA



62141 DB

62141 AA

Reducer, 6/3 mm

62141 DB

Reducer, 11/5 mm, flip style
(optional for larger ovaries)

Reducers for use with trocars and TERNAMIAN EndoTIP cannulas with silicone leaflet valve



30160 RE

NEW

30160 RE

Reducer, 6/3 mm

NEW

30103 RE

Reducer, 11/5 mm

NEW

30141 HB

Reducer, 13/5 mm and 13.5/5 mm

Thread Sleeves for Trocars



30160 S

30160 S

Thread Sleeve, for trocar size 6 mm,
color code: black

30103 S

Thread Sleeve, for trocar size 11 mm,
color code: green

Distance holders to limit the depth of puncture



30160 D

30160 D

Distance Holder, for trocars size 6 mm,
color code: black

30103 D

Distance Holder, for trocars size 11 mm,
color code: green

Sealing Caps, with color code



30100 XA

30100 XA

Sealing Cap, size 5 mm, autoclavable, package of 5,
for use with trocars size 6 mm and extractors/reduction
sleeves in conjunction with instruments,
color code: black

30100 XB

Sealing Cap, size 10 mm, autoclavable, package of 5,
for use with trocars size 11 mm and extractors/reduction
sleeves in conjunction with instruments,
color code: green

Coagulation and Dissection Electrodes

without suction channel, insulated sheath,
with connector pin for unipolar coagulation

Size 5 mm

Operating instruments, **lengths 30 and 36 cm**,
for use with **trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm



Length	Instrument
30 cm	
36 cm	
43 cm	

Distal Tip	Instrument	
	26772 UF	Coagulation and Dissection Electrode, L-shaped
	26775 UF	
	26778 UF	
	26775 UE	Coagulation and Dissection Electrode, spatula-shaped, blunt
	26778 UE	
	26775 S	Coagulation and Dissection Electrode, knife-shaped
	26775 R	MANGESHIKAR Coagulation Electrode, with ball end

Special Features:

- The hook electrode is suitable for resection, insulating structures as well as dissection and coagulation.
- The distal tip is semicircular: The thick outer curvature enables safe dissection.
- The thin inner curvature enables fast work with minimal coagulation.
- The striated handle allows easier gripping.
- The hook length enables ergonomic work.


Length	Instrument
36 cm	
43 cm	

Distal Tip	Instrument	
	26775 C	CADIERE Coagulation and Dissection Electrode, L-shaped, with cm-marking, distal tip tapered
	26775 CL	

Units for use with Coagulation and Dissection Electrodes see chapter 20, UNITS

Scissors

CLICK'line® – rotating, dismantling,
with and without connector pin for unipolar coagulation



Size 5 mm


Operating instruments, **lengths 30, 36 and 43 cm**,
for use **with trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm

Length	Handle						
	33151	33121	33125	33149	33161	33131	33147
30 cm							
36 cm							
43 cm							

Double action jaws

Working Insert	Complete Instrument						
34210 MS	34251 MS	34221 MS	34225 MS	34249 MS	34261 MS	34231 MS	34247 MS
34310 MS	34351 MS	34321 MS	34325 MS	34349 MS	34361 MS	34331 MS	34347 MS
34410 MS	34451 MS	34421 MS	34425 MS	34449 MS	34461 MS	34431 MS	34447 MS




CLICK'line® METZENBAUM Scissors, curved

34210 MA	34251 MA	34221 MA	34225 MA	34249 MA	34261 MA	34231 MA	34247 MA
34310 MA	34351 MA	34321 MA	34325 MA	34349 MA	34361 MA	34331 MA	34347 MA
34410 MA	34451 MA	34421 MA	34425 MA	34449 MA	34461 MA	34431 MA	34447 MA



CLICK'line® Scissors, serrated blades, spoon-shaped, curved

34210 MW	34251 MW	34221 MW	34225 MW	34249 MW	34261 MW	34231 MW	34247 MW
34310 MW	34351 MW	34321 MW	34325 MW	34349 MW	34361 MW	34331 MW	34347 MW
34410 MW	34451 MW	34421 MW	34425 MW	34449 MW	34461 MW	34431 MW	34447 MW



CLICK'line® Scissors, serrated, curved, conical

34210 MD	34251 MD	34221 MD	34225 MD	34249 MD	34261 MD	34231 MD	34247 MD
34310 MD	34351 MD	34321 MD	34325 MD	34349 MD	34361 MD	34331 MD	34347 MD




CLICK'line® Scissors, straight

Handles for Suction and Irrigation



37112A

37112A **Handle**, straight, with clamping valve, for suction and irrigation, **autoclavable**



37113A

37113A **Handle**, pistol grip, with clamping valve, for suction and irrigation, **autoclavable**



37360LH

Suction and Irrigation Tube, with lateral holes, size 5 mm, length 36 cm, for use with suction and irrigation handles

Dissecting Electrodes

with Exchangeable Electrode Tip

Insulated Sheath with Connector Pin for Unipolar Coagulation



30775UF



30775UF

Coagulation and Dissection Electrode, L-shaped, with connector pin for unipolar coagulation, size 5 mm, length 36 cm including:

Outer Sheath, insulated

Plastic Handle

Electrode, L-shaped

30775UFE

Exchangeable Electrode Tip, L-shaped, autoclavable, package of 6

Dissecting and Grasping Forceps

CLICK'line – rotating, dismantling, insulated,
with connector pin for unipolar coagulation



Size 5 mm


Operating instruments, **lengths 30 and 36 cm**,
for use **with trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm

Length	Handle					
	33151	33152	33153	33156	33121	33149
30 cm						
36 cm						
43 cm						

Double action jaws

Working Insert	Complete Instrument					
33310 MH	33351 MH	33352 MH	33353 MH	33356 MH	33321 MH	33349 MH



CLICK'line Dissecting and Grasping Forceps, atraumatic

33310 DF	33351 DF	33352 DF	33353 DF	33356 DF	33321 DF	33349 DF
33410 DF	33451 DF	33452 DF	33453 DF	33456 DF	33421 DF	33449 DF



CLICK'line Dissecting and Grasping Forceps, atraumatic

33210 DN	33251 DN	33252 DN	33253 DN	33256 DN	33221 DN	33249 DN
33310 DN	33351 DN	33352 DN	33353 DN	33356 DN	33321 DN	33349 DN
33410 DN	33451 DN	33452 DN	33453 DN	33456 DN	33421 DN	33449 DN



CLICK'line Dissecting and Grasping Forceps, "dolphin nose"

33310 KJ	33351 KJ	33352 KJ	33353 KJ	33356 KJ	33321 KJ	33349 KJ
33410 KJ	33451 KJ	33452 KJ	33453 KJ	33456 KJ	33421 KJ	33449 KJ



CLICK'line REDDICK-OLSEN Dissecting and Grasping Forceps


Dissecting and Grasping Forceps

CLICK*line* – rotating, dismantling, insulated,
with connector pin for unipolar coagulation

Size 5 mm

Operating instruments, **lengths 30 and 36 cm**,
for use **with trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm



Length	Handle						
	33151	33152	33153	33156	33121	33125	33149
30 cm							
36 cm							
43 cm							

Double action jaws

Working Insert	Complete Instrument						
33210 AF	33251 AF	33252 AF	33253 AF	33256 AF	33221 AF	33225 AF	33249 AF
33310 AF	33351 AF	33352 AF	33353 AF	33356 AF	33321 AF	33325 AF	33349 AF
33410 AF	33451 AF	33452 AF	33453 AF	33456 AF	33421 AF	33425 AF	33449 AF



CLICK*line* Grasping Forceps, atraumatic, fenestrated

33310 OD	33351 OD	33352 OD	33353 OD	33356 OD	33321 OD	33325 OD	33349 OD
----------	----------	----------	----------	----------	----------	----------	----------



CLICK*line* Grasping Forceps, fenestrated,
with especially fine atraumatic serration

33310 FG	33351 FG	33352 FG	33353 FG	33356 FG	33321 FG	33325 FG	33349 FG
33410 FG	33451 FG	33452 FG	33453 FG	33456 FG	33421 FG	33425 FG	33449 FG



CLICK*line* Grasping Forceps, atraumatic

33210 D	33251 D	33252 D	33253 D	33256 D	33221 D	33225 D	33249 D
33310 D	33351 D	33352 D	33353 D	33356 D	33321 D	33325 D	33349 D



CLICK*line* Grasping Forceps, atraumatic, spoon-shaped, with multiple teeth

33310 K	33351 K	33352 K	33353 K	33356 K	33321 K	33325 K	33349 K
---------	---------	---------	---------	---------	---------	---------	---------



CLICK*line* Grasping Forceps, atraumatic, fenestrated


Dissecting and Grasping Forceps

CLICK'line – rotating, dismantling, insulated,
with connector pin for unipolar coagulation

Size 5 mm

Operating instruments, **lengths 30 and 36 cm**,
for use **with trocars size 6 mm**


Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm



Length	Handle					
	33151	33152	33153	33156	33121	33149
30 cm						
36 cm						
43 cm						


Double action jaws

Working Insert	Complete Instrument					
33210 A	33251 A	33252 A	33253 A	33256 A	33221 A	33249 A
33310 A	33351 A	33352 A	33353 A	33356 A	33321 A	33349 A
33410 A	33451 A	33452 A	33453 A	33456 A	33421 A	33449 A



CLICK'line BABCOCK Grasping Forceps, atraumatic, fenestrated

33210 C	33251 C	33252 C	33253 C	33256 C	33221 C	33249 C
33310 C	33351 C	33352 C	33353 C	33356 C	33321 C	33349 C
33410 C	33451 C	33452 C	33453 C	33456 C	33421 C	33449 C



CLICK'line Bowel Grasper, fenestrated

33310 CK	33351 CK	33352 CK	33353 CK	33356 CK	33321 CK	33349 CK
33410 CK	33451 CK	33452 CK	33453 CK	33456 CK	33421 CK	33449 CK



CLICK'line Bowel Grasper, fenestrated, short


33310 UK	33351 UK	33352 UK	33353 UK	33356 UK	33321 UK	33349 UK
----------	----------	----------	----------	----------	----------	----------



CLICK'line KOH Urethra Grasping Forceps, insulated jaws, concave inwardly curved jaws, continuously toothed

Dissecting and Grasping Forceps

CLICK'linē – rotating, dismantling, insulated,
with connector pin for unipolar coagulation



Size 5 mm

Operating instruments, **lengths 30 and 36 cm**,
for use **with trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm

Length	Handle					
	33151	33152	33153	33156	33121	33149
30 cm						
36 cm						
43 cm						

Single action jaws

Working Insert	Complete Instrument					
33310 M	33351 M	33352 M	33353 M	33356 M	33321 M	33349 M
	CLICK'linē Grasping Forceps , with multiple teeth, width of jaws 3 mm, for atraumatic and accurate grasping					

33310 FM	33351 FM	33352 FM	33353 FM	33356 FM	33321 FM	33349 FM
33410 FM	33451 FM	33452 FM	33453 FM	33456 FM	33421 FM	33449 FM

	CLICK'linē Claw Forceps , 2 x 3 teeth					
--	--	--	--	--	--	--

33210 ME	33251 ME	33252 ME	33253 ME	33256 ME	33221 ME	33249 ME
33310 ME	33351 ME	33352 ME	33353 ME	33356 ME	33321 ME	33349 ME
33410 ME	33451 ME	33452 ME	33453 ME	33456 ME	33421 ME	33449 ME


	CLICK'linē Grasping Forceps , with multiple teeth, width of jaws 4.8 mm, for atraumatic and accurate grasping					
--	--	--	--	--	--	--

33210 MG	33251 MG	33252 MG	33253 MG	33256 MG	33221 MG	33249 MG
33310 MG	33351 MG	33352 MG	33353 MG	33356 MG	33321 MG	33349 MG
33410 MG	33451 MG	33452 MG	33453 MG	33456 MG	33421 MG	33449 MG

	CLICK'linē Grasping Forceps , “tiger jaws”, 2 x 4 teeth, for grasping and removal of solid organs, in particular during adhesiolyses					
--	---	--	--	--	--	--

Dissecting and Grasping Forceps

CLICK'linē – rotating, dismantling, insulated,
with connector pin for unipolar coagulation



Size 5 mm

Operating instruments, **lengths 30 and 36 cm**,
for use **with trocars size 6 mm**

Operating instruments, **length 43 cm**,
for use with telescopes with inbuilt
working channel and trocars size 6 mm

Length	Handle					
	33151	33152	33153	33156	33121	33149
30 cm						
36 cm						
43 cm						


Double action jaws

Working Insert	Complete Instrument					
33210 MD	33251 MD	33252 MD	33253 MD	33256 MD	33221 MD	33249 MD
33310 MD	33351 MD	33352 MD	33353 MD	33356 MD	33321 MD	33349 MD
33410 MD	33451 MD	33452 MD	33453 MD	33456 MD	33421 MD	33449 MD



CLICK'linē KELLY Dissecting and Grasping Forceps

33210 ML	33251 ML	33252 ML	33253 ML	33256 ML	33221 ML	33249 ML
33310 ML	33351 ML	33352 ML	33353 ML	33356 ML	33321 ML	33349 ML
33410 ML	33451 ML	33452 ML	33453 ML	33456 ML	33421 ML	33449 ML



CLICK'linē KELLY Dissecting and Grasping Forceps, long

33310 DK	33351 DK	33352 DK	33353 DK	33356 DK	33321 DK	33349 DK
----------	----------	----------	----------	----------	----------	----------



CLICK'linē DeBAKEY Grasping Forceps,
curved and slender jaws, atraumatic

33310 UY	33351 UY	33352 UY	33353 UY	33356 UY	33321 UY	33349 UY
----------	----------	----------	----------	----------	----------	----------



CLICK'linē UYAMA Dissecting and Grasping Forceps,
long jaws, with reservoir, fine, for dissection

33310 NO	33351 NO	33352 NO	33353 NO	33356 NO	33321 NO	33349 NO
----------	----------	----------	----------	----------	----------	----------




CLICK'linē NOSHIRO Dissecting and Grasping Forceps,
short jaws, fine, for dissection

Applicators for Titanium Clips

dismantling, rotating

Size 10 mm

Operating instruments, **length 36 cm**,
for use **with trocars size 11 mm**



30444 LR

Clip Applicator, dismantling, rotating, size 10 mm, length 36 cm, for PILLING-WECK Titanium Clips 30460 AL (medium-large), with ratchet to lock the jaw holding the clip

including:

Metal Handle, with ratchet

Metal Outer Sheath

Clip Applicator Insert



actual size

30460 AL

PILLING-WECK Titanium Clips, medium-large, box with 16 sterile cartridges, 10 clips each, for use with Clip Applicator 30444 LR



30444 GA

Clip Applicator, dismantling, rotating, size 10 mm, length 36 cm, for PILLING-WECK Titanium Clips 30460 AP (medium), with ratchet to lock the jaw holding the clip

including:

Metal Handle, with ratchet

Metal Outer Sheath

Clip Applicator Insert



actual size

30460 AP

PILLING-WECK Titanium Clip, medium, box with 10 sterile cartridges, 25 clips each, for use with Clip Applicator 30444 GA



30444 ET

Clip Applicator, dismantling, rotating, size 10 mm, length 36 cm, for Ethicon Titanium Clips LT-300 (medium-large), with ratchet to lock the jaw holding the clip

including:

Metal Handle, with ratchet

Metal Outer Sheath

Clip Applicator Insert

Caution:


The use of clips other than those indicated can lead to damage of the jaws.

Components/Spare Parts see chapter 21

Knot Tiers


Size 5 mm

Operating instruments, **lengths 36 and 46 cm**,
for use **with trocars size 6 or 11 mm with reduction sleeves**



26596 CL

CICE Knot Tier, CLERMONT-FERRAND model,
for extracorporeal knotting, size 5 mm, length 36 cm



26596 SK

KÖCKERLING Knot Tier, for extracorporeal knotting,
size 5 mm, length 36 cm

26596 SKL

Same, length 46 cm