



PAGINA 0

**CREATA PENTRU A
PUTEA SEMNA
ELECTRONIC
DOCUMENTUL**

Launcher™ Guide Catheter



MID-RANGE
FLEXIBILITY

SUPPORTIVE
SECONDARY
CURVE

MEDTRONIC FULL-WALL TECHNOLOGY:
Encapsulated flat-wire braid enables thinner,
more robust walls without compromising support,
retention, or torque response.

Encapsulated Wire Braid

BALANCED PERFORMANCE THAT'S AHEAD OF THE CURVE

The Launcher™ guide catheter offers you a blend of flexibility, support and visualisation — the capability you need to respond to your challenging cases.

Advanced Platform Designed for Multiple Interventional Approaches

- Designed for radial and femoral approaches
- Flexible distal segment enables you to engage the guide for backup support
- Full-Wall Technology construction provides kink resistance and stable torque control
- Supportive secondary curve for increased backup support and curve retention

Enhanced Visualisation

- Large 6 F.I.D. (0.071 inches) allows for increased contrast flow for enhanced visualisation
- Vest-Tech Nylon outer jacket provides excellent hub to tip radiopacity

Full Range of Traditional and Specialty Curves

- Extensive collection of traditional and specialty coronary curves available

LAUNCHER™ 6 F GUIDE CATHETERS


LEFT CORONARY CURVES					Backup Support Left
Left Standard (JL)	Short Left (SL)	JCL	Amplatz Left	Short Amplatz Left	EBU
LA6JL25	LA6SL30	LA6JCL30	LA6AL75	LA6SAL75	LA6EBU30
LA6JL30	LA6SL35	LA6JCL35	LA6AL10	LA6SAL10	LA6EBU325
LA6JL35	LA6SL40	LA6JCL40	LA6AL15	LA6SAL15	LA6EBU35
LA6JL40	LA6SL45	LA6JCL45	LA6AL175	LA6SAL20	LA6EBU375
LA6JL40P^	LA6SL50	LA6JCL50	LA6AL20	LA6SAL25	LA6EBU40
LA6JL45	LA6SL60		LA6AL25	LA6SAL30	LA6EBU45
LA6JL50			LA6AL30	LA6SAL40	LA6EBU50
LA6JL60			LA6AL40		
			LA6AL50		
Left Standard (FL)	Short Left (FLST)				
LA6FL35	LA6FL35ST				
LA6FL40	LA6FL40ST				
LA6FL40G^^	LA6FL45ST				
LA6FL40P^	LA6FL50ST				
RIGHT CORONARY CURVES					
Right Standard	Short Right	JCR	Amplatz Right	Shepherd's Crook Right	Backup Support Right
LA6JR30	LA6SR30	LA6JCR35	LA6AR10	LA6SCR35	LA6RBU35
LA6JR35	LA6SR35	LA6JCR40	LA6AR20	LA6SCR40	LA6RBU40
LA6JR40	LA6SR40	LA6JCR45	LA6ALR12	LA6SCR50	LA6RBU45
LA6JR40P^	LA6SR45	LA6JCR50			
LA6JR45	LA6SR50		Short Amplatz Right		ECR
LA6JR50	LA6SR60		LA6SAR10		LA6ECR30
LA6JR60			LA6SAR20		LA6ECR35
					LA6ECR40
					LA6ECR45
Right Standard (FR)					
LA6FR35					
LA6FR40					3-Dimensional
					LA63DRC
					LA63DRIGHT
OTHER CURVES					
Multipurpose	MAC (Multi-Aortic)	Bypass	Champ	Peripheral	Radial Specialty
LA6MB1	LA6MAC30	LA6LCB	LA6CHAMP05	LA6CHAMP15K+	LA6LARA
LA6MB2	LA6MAC35	LA6LCBD*	LA6CHAMP10	LA6CHAMP20K+	LA6MRADIAL
LA6MP1	LA6MAC375	LA6RCB	LA6CHAMP15	LA6CHAMP25K+	LA6MRESS
LA6MP2	LA6MAC40	LA6RCBD*	LA6CHAMP20	LA6IMAK+	LA6RRAD
LA6MPST	LA6MAC45	LA6RCBII	LA6CHAMP25	LA6JR40K+	
LA6HSI	LA6MAC50	LA6IMA	LA6CHAMP30	LA6MP1K+	NOTO
LA6HSII		LA6IMAD*	LA6CHAMP35	LA6MPHK+	LA6NOTO
LA6RMHSII			LA6CHAMP40	LA6PK1W++	
LA6HSIII				LA6RDCK+	EL GAMAL
					LA6ELGAMAL

D* = 90 cm shaft length Shaft Diameter: 6 F K+ = 55 cm shaft length Inner Diameter: 0.071 in. W+++ = 47 cm shaft length P^ = Posterior
 To order curves with sideholes, please add "SH" to end of model number G^^ = Anterior 5 packs available in select curves

Indications for Use: The Medtronic Guiding Catheter is designed to provide a pathway through which therapeutic devices are introduced. The guiding catheter is intended to be used in the coronary or peripheral vascular system.
Important Information: Refer to *Instructions for Use* supplied with the device for indications, warnings, precautions, and adverse effects.
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
 The curve recommendations in this guide are not meant to replace patient-specific clinical judgement.



Medtronic
 3576 Unocal Place
 Santa Rosa, CA 95403
 USA
 Tel: 707.525.0111

LifeLine Customer Support
 Tel: 877.526.7890
 Tel: 763.526.7890

Product Services
 Tel: 888.283.7868
 Fax: 800.838.3103

Europe
Medtronic Intl. Trading S.A.R.L.
 Route du Molliiau 31
 Case Postale
 CH-1131 Tolochenaz
 Switzerland
 Tel: 41.21.802.7000
 Fax: 41.21.802.7900

Canada
Medtronic of Canada Ltd.
 99 Hereford Street
 Brampton, Ontario L6Y 0R3
 Canada
 Tel: 905.460.3800
 Fax: 905.460.3998
 Toll-free: 800.268.5346

Asia Pacific
Medtronic Intl. Ltd.
 49 Changi South Avenue 2
 Nasaco Tech Centre
 Singapore 486056
 Singapore
 Tel: 65.6436.5000
 Fax: 65.6776.6335

Latin America
Medtronic
 9850 Doral Blvd., 4th Floor
 Miami, FL 33178
 USA
 Tel: 786.709.4200
 Fax: 786.709.4244

medtronic.com