

Xiamen Winner Medical Co., Ltd.

Address: 4F, No. 98 Huli Industrial Park, Mei Xi Dao, Tong'an, Xiamen, China 361100

Tel:+86-592-5560185

E-Mail: info@winner-medi.com

CPAP Mask

Usage

CPAP mask is used with CPAP machine. CPAP machines use mild air pressure to keep the airways open, and are typically used by patients who have breathing problems during sleep.

Product Description

Non-vented or vented available, PC mask shell with silicone cushion, forehead support with silicone pad, rotational connector, adjustable head strap,

Specification

1. Material PC for hard plastic components, silicone for cushion and valve.
2. Dimension 22mmOD, according to ISO 5356-1
3. Operation temperature: $-18^{\circ}\text{C} \sim 40^{\circ}\text{C}$
4. Storage relative humidity: $< 85\%$
5. Storage temperature: $-20^{\circ}\text{C} \sim 55^{\circ}\text{C}$
6. Validity: 2 years

Size: M/L

Package: a colorful zipper bag

