



[Home](#) » [Products](#) » [Medical Imaging Products](#) » [Angiographic Syringes](#)



Angiographic Syringes

· Angiographic syringe is intended to deliver contrast media, which is compatible with high-pressure injection equipment during CT, DSA and MRI.

[INQUIRY](#)

FEATURES

- Used for CT, MR, DSA-- It could be used for computed tomographic scanning, angiography, Nuclear magnetic resonance and imaging. The contrast media carried is the gold standard for identifying cancerous tissue, presenting high-definition images for lesions;
- Great variety of Syringes; it could match so many injectors of nearly all manufacturers on the market;



special configurations;

- Good materials to provide more resistant to high temperature and high pressure, good permeability;
- Clear body provides excellent clarification and smooth as glass feeling; it could see air bubble and contrast medium.

USE METHOD

This video tutorial shows you how to use our angiographic syringe step by step -



D A D A M E T E D



Product No.		Volume(ml)	Description	Models	
Single	Dual		Used for Medrad	Single	Dual
42.16.10000	42.16.10001	200	Used for CT: MCT & plus; OP-100; VCT610; ECT710	CT-200-A1	CT-200/200-A1
42.16.10002	42.16.10003	200	Used for CT: STELANT	CT-200-A2	CT-200/200-A2
42.16.10004		150	Used for DSA: Mark V&Mark V Provis	DSA-150-A1	
42.16.10005		200	Used for DSA: Mark V& Mark V Provis	DSA-200-A1	
42.16.10007		130	Used for DSA: Mark III & Mark IV	DSA-130-A1	
	42.16.10008	65	Used for MRI: Spectris		MRI-65/65-A1
	42.16.10009	115	Used for MRI: Spectris		MRI-115/115-A1
	42.16.10010	65/115	Used for MRI: Spectris		MRI-65/115-A1
42.16.10011	42.16.10012	190	Used for CT: SALENT	CT-190-A1	CT-190/190-A1
42.16.10016		150	Used for DSA: Mark 7	DSA-150-A2	



Single	Dual		Used for Libel-Flarsheim	Single	Dual
42.16.20000	42.16.20001	200	Used for CT: CT 9000& CT 9000 ADV	CT-200-B1	CT-200/200-B
42.16.20002		150	Used for DSA: Angiomat 6000	DSA-150-B1	
42.16.20003		150	Used for CT & DSA: Angiomat Illumena	CT/DSA-150-B1	
42.16.20005	42.16.20006	60	Used for MRI: Optistar	MRI-60-B1	MRI-60/60-B:

Product No.		Volume(ml)	Description	Models	
Single	Dual		Used for EZEM	Single	Dual
42.16.50000	42.16.50001	200	Used for CT: Empower	CT-200-EM	CT-200/200-EI



Single	Dual		Used for Nemoto Kyorindo	Single	Dual
42,16,30000	42.16.30001	100	Used for CT: A-25	CT-100-NE	CT-100/100-N
42.16.30002	42,16.30003	200	Used for CT: A-25 A-60	CT-200-NE	CT-200/200-N
42.16.30005		120	Used for DSA: Nemoto 120S	DSA-120-NE	
42.16.30006	42.16.30007	60	Nemoto Dual Shot. Sonic Shot 50	MRI-60-NE	MRI-60/60-N
	42,16,30008	100/200	Used for CT.A-25		CT-100/200-N
	42.16.30012	60/200	Used for CT: A-25.A-60		CT-60/200-N

Product No.		Volume(ml)	Description	Models	
Single	Dual		Used for Metron	Single	Dual



	42.16.40004	65/200	Used for MRI: ACCUTRON		MRI-65/200-M
--	-------------	--------	------------------------	--	--------------



CE CERTIFICATION

This is a CE certified product. Here is the picture for your reference:



Manufacturer:

SCW MEDICATH LTD.

NO.4, Baolong 6th Road, Baolong Industrial
Town, Longgang District, Shenzhen,
518116, Guangdong, P.R. China

We, the manufacturer, herewith declare that the products
Angiographic Syringes

Models: DSA-60-A1, DSA-150-A1, CT-100-X1, CT-200-A1, MRI-115-A1,

DSA-100-B1, DSA-150-B1, CT-200-B1, MRI-65-A1, CT-200-A2,
DSA-200-A1, DSA-130-A1, CT/DSA-150-B1, MRI-60-B1, CT-200-NE,
DSA-90-NE, DSA-120-NE, CT-100-NE, CT-190-A1, CT-190/190-A1,
CT-200/200-A2, CT-200/200-B1, CT-100/100-NE, CT-100/200-NE,
CT-200/200-NE, CT-200-EM, CT-200/200-EM, CT-200-MED,
CT-200/200-MED, MRI-65/115-A1, MRI-65/65-A1, MRI-115/115-A1,
MRI-60/60-B1, MRI-60/60-NE, MRI-65-MED, MRI-65/200-MED

meet the provisions of Directive 93/42/EEC which apply to them.

The medical device has been assigned to class IIa according to Annex IX Rule 2 of the
Directive 93/42/EEC. It bears the mark

CE 0197

The product concerned has been manufactured under a quality management system
according to Annex II of Directive 93/42/EEC.

Compliance of the designated product with the Directive 93/42/EEC has been assured via
assessment of the quality management system by the Notified Body

TÜV Rheinland LGA Products GmbH

Tillystraße 2
90431 Nürnberg
Deutschland

Certificate No.: HD 60144232 0001

Issue date: 26.05.2020

Expiry date: 26.05.2024

following the procedure relating to the EC Declaration of Conformity set out in Annex II of
Directive 93/42/EEC.

This Declaration of conformity is valid in connection with the release document for the
respective batch of produced devices.

The above mentioned declaration of conformity is exclusively under the responsibility of

SCW MEDICATH LTD.

NO.4, Baolong 6th Road, Baolong Industrial Town,
Longgang District, Shenzhen, 518116, Guangdong, P.R. China

Shenzhen, 2020/07/16

Place, date

William Xie, RA

Legally binding signature, Function



PRODUCT RECOMMENDATION



Angiographic Syringes Accessories



*Your Email

Your Phone

*Message

*Verify Code

x c R 5

*We respect your privacy. When you submit your contact information, we agree to only contact you in accordance with our [Privacy Policy](#).

Submit

We'd Love to Hear From You

Curious about product features, price quote, registration, delivery time, etc? We're ready to answer any and all questions within one business day.

Contact Us



Email: info@scwmed.com

Mobile Phone: +8613691980679

Telephone: 86-755-89926266 ext 4

Fax: 86-755-89668199 ext 827

NO.4 Baolong 6th Road, Baolong industrial Town,
Longgang District, Shenzhen, P.R. China



Copyright © SCW Medicath Ltd All rights reserved. [Privacy](#) [Site Map](#) [Terms of use](#)