



Original, pre-prepared (Ready-to-Use) system. The system includes a one-direction valve attached to an air filter. Nurse, while preparing the gas, puts the syringe on to the air filter. Nurse holds the bottle with left hand and stabilizes it while the right hand just holds the syringe without any movement. When nurse pushes the right hand from the bottom hardly, the one-direction valve opens. This action allows the movement of gas from the chamber into the syringe. This system controls the movement of gas (not automatically) in only one direction (prevents escape back into the chamber) into the syringe

ml	SulphurHexafluoride	Octafluoropropane
10ml	VT10SFA®	VT10CFA®
30ml	VT30SFA	VT30CFA
50ml	VT50SFA	VT50CFA

Product Description: Presented in Aerosol  
Store at room temperature