



Perfect Ink Pads

Will not dry out when left unused. Ink pad features a slight residue left on hands that is easily removed by rubbing them together or wiping with a towel. Hard ceramic inking surface to prevent over-inking fingers.

Item #	Description	# of Prints
09610	Perfect Ink Sm Rectangular Pad	600 Rolled prints
09614	Perfect Ink Sm Round Pad	1,200 Flat prints
09616	Perfect Ink Lg Round Pad	4,800 prints
09612	Perfect Ink Lg Rectangular Pad	16,000 prints



Fingerprint (Porelon) Pad

For use in the field or office. Simply place fingertips on the pre-inked touch pad, then print. Special permanent ink is non-fading and quick drying. Requires cleaner to remove residue from fingertips.

Item #	Description	# of Prints
05572	Porelon 2" Round Pad	5,000
05575	Porelon Pad Refill	15,000
05576	Porelon Combo Pocket/Desk Pad	15,000

Latent Print Reference Pad

Latent print reference pads create perfect, consistent test latent fingerprints again and again. Simply press your finger on the desired pad and the perfect amount of latent print residue is left behind. Latent print development techniques can then be accurately tested and gauged based on this reference standard, regardless of who placed the print.

Item #	Description
95389	Sebacous Oil Reference Pad
15090	Amino Acid Reference Pad



Fingerprint Ink

Quick-drying, black as coal, glides like velvet and provides excellent ridge detail.

Item #	Description	Size
05570	Fingerprint Ink Tube	4 oz.
05569	Fingerprint Ink Tube	2 oz.

Pre-Inked Slabs

Our pre-inked slabs have a thin, even coat of fingerprint ink between two pieces of polyester film. Just peel the film surfaces apart and you have two perfectly inked slabs ready for dozens of rolled prints! They're ideal for making elimination prints on-site.

Item #	Description	Size	Quantity
55401	Pre-Inked Slabs	1" x 5"	200 pk.
05571	Pre-Inked Slabs	3" x 10"	100 pk.

