




ARTHROSCOPY
IMPLANTS & INSTRUMENTS



- **T315A** | Cannulated Interference 7mm (ACL Screw)
- **T315B** | Cannulated Interference 8mm (ACL Screw)
- **T315C** | Cannulated Interference 9mm (ACL Screw)
- **T315D** | Cannulated Interference 10mm (ACL Screw)

MATERIAL - Titanium

LENGTH - 20mm to 35mm (5mm Increment)



- **T316** | Herbert Screw 3.5mm Cannulated

MATERIAL - Titanium

LENGTH - 14mm to 40mm (2mm Increment)



- **T316A** | Herbert Screw 2.4mm Cannulated

MATERIAL - Titanium

LENGTH - 14mm to 30mm (2mm Increment)



- **T317.01** | Endo Buttons

MATERIAL - Titanium

LENGTH - 15mm to 30mm (5mm Increment)



- **ST317** | Endo Buttons with Wire

MATERIAL - Titanium

LENGTH - 15mm to 30mm (5mm Increment)



- **T318.01** | Suture Disc

MATERIAL - Titanium



- **T323** | Anchor Screw

MATERIAL - Titanium

LENGTH - 3.0mm, 5.0mm



- **ST323** | Suture Anchor Screw

MATERIAL - Titanium

LENGTH - 3.0mm, 5.0mm



• **T384** | Spiked Washer

MATERIAL - Titanium
SIZE - 4.0mm (5 Spikes)
 4.5mm (7 Spikes)
 6.5mm (9 Spikes)



• **T385** | Spiked Ligament Staple

MATERIAL - Titanium
LENGTH - 6.0mm, 8.0mm, 10.0mm

Instruments



673.001 | Main Jig



673.002 | Acl elbow aimer



673.004 | Pcl Femoral



673.005 | Pcl Tibia



673.006 | Angle Bullet



673.007 | Graft passing 2.4mm drill bit



661.005 | Beth Pin Drill tip



663 | ACL Tibia Straight Cannulated Reamer
6mm, 7mm, 8mm, 9mm, 10mm, 11mm



- 662 | ACL Femoral Reamer Cannulated Flower Tip
6mm, 7mm, 8mm, 9mm, 10mm, 11mm



- 661.002 | Endoscopic Cannulated Reamer



- 673.008 | Depth Gauge for ACL/PCL reconstruction



- 666.050 | Femoral Aimer With Offset 5mm
- 666.070 | Femoral Aimer With Offset 7mm



- 661.078 | Tunnel Plug 7-8mm
- 661.089 | Tunnel Plug 8-9mm
- 661.910 | Tunnel Plug 9-10mm



- 661.070 | Tunnel Dilator 7mm
- 661.080 | Tunnel Dilator 8mm
- 661.090 | Tunnel Dilator 9mm



- 659.035 | Screw Driver For Interference Screw



- 673.009 | Impactor for spiked ligament staple



- 673.010 | Arthroscopic Probe



- 674.030 | Microfracture Awl 30°
- 674.045 | Microfracture Awl 45°
- 674.090 | Microfracture Awl 90°



• **661.003** | Graft Master Board, Suturewise with Tensiometer



• **668.001** | Graft Sizing Block



• **665.060** | Tendon Stripper Close 6mm

• **664.070** | Tendon Stripper Open 7mm



• **674.001** | Straight Rasp





• **674.002** | Convex Rasp



• **674.005** | Femoral PCL Elevator



• **674.006** | T-Handle for Tunnel Dilators



Instruments

673.001	Main Jig	1 Each
673.002	Acl elbow aimer	1 Each
673.004	Pcl Femoral	1 Each
673.005	Pcl Tibia	1 Each
673.006	Angle Bullet	1 Each
673.007	Graft passing 2.4mm drill bit	1 Each
661.005	Beth Pin Drill tip	1 Each
663	ACL Tibia Straight Cannulated Reamer 6mm, 7mm, 8mm, 9mm, 10mm, 11mm	1 Each
662	ACL Femoral Reamer Cannulated Flower Tip 6mm, 7mm, 8mm, 9mm, 10mm, 11mm	1 Each
661.002	Endoscopic Cannulated Reamer	1 Each
673.008	Depth Gauge for ACL/PCL reconstruction	1 Each
666.050	Femoral Aimer With Offset 5mm	1 Each
666.070	Femoral Aimer With Offset 7mm	1 Each
661.078	Tunnel Plug 7-8mm	1 Each
661.089	Tunnel Plug 8-9mm	1 Each
661.910	Tunnel Plug 9-10mm	1 Each
661.070	Tunnel Dilator 7mm	1 Each
661.080	Tunnel Dilator 8mm	1 Each
661.090	Tunnel Dilator 9mm	1 Each

Instruments

659.035	Screw Driver For Interference Screw	1 Each
673.009	Impactor for spiked alignment staple	1 Each
673.010	Arthroscopic Probe	1 Each
674.030	Microfracture Awl 30°	1 Each
674.045	Microfracture Awl 45°	1 Each
674.090	Microfracture Awl 90°	1 Each
661.003	Graft Master Board, Suturewise with Tensiometer	1 Each
668.001	Graft Sizing Block	1 Each
665.060	Tendon Stripper Close 6mm	1 Each
664.070	Tendon Stripper Open 7mm	1 Each
674.001	Straight Rasp	1 Each
674.002	Convex Rasp	1 Each
674.005	Femoral PCL Elevator	1 Each
674.006	T-Handle for Tunnel Dilators	1 Each



- **669.001** | Jig for Herbert Screw



- **658.024** | Herbert Screw Tap Cannulated 2.4mm
- **658.035** | Herbert Screw Tap Cannulated 3.5mm
- **658.045** | Herbert Screw Tap Cannulated 4.5mm



- **659.025** | Screw Driver Cannulated 2.5mm Tip



- **659.020** | Screw Driver Cannulated 2.0mm Tip




- **657.020** | Drill Bit Cannulated 2.0mm
- **657.027** | Drill Bit Cannulated 2.7mm
- **657.032** | Drill Bit Cannulated 3.2mm



- **671.100** | Guide Wire 1.0mm x 225mm
- **671.125** | Guide Wire 1.25mm x 225mm



- **670.125** | Direct Measuring Device for 1.25mm Guide Wire



Implants

T316	Herbert Screw 3.5mm Cannulated 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm	4 Each
T316A	Herbert Screw 2.4mm Cannulated 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm	4 Each

Instruments

669.001	Jig for Herbert Screw	1 Each
658.035	Herbert Screw Tap Cannulated 3.5mm for upper thread of Herbert Screw	1 Each
658.045	Herbert Screw Tap Cannulated 4.5mm for Lower thread of Herbert Screw	1 Each
658.024	Herbert Screw Tap Cannulated 2.4mm for Lower thread of Herbert Screw	1 Each
659.025	Screw Driver Cannulated 2.5mm Tip	1 Each
659.020	Screw Driver Cannulated 2.0mm Tip	1 Each
657.027	Drill Bit Cannulated 2.7mm	1 Each
657.032	Drill Bit Cannulated 3.2mm	1 Each
657.020	Drill Bit Cannulated 2.0mm	1 Each
670.125	Direct Measuring Device for 1.25mm Guide Wire	1 Each
671.125	Guide Wire 1.25mm x 225mm	3 Each
671.100	Guide Wire 1.0mm x 225mm	3 Each



House of Ortho Tools

Hardik International Pvt. Ltd.

Plot No. G-2041/42/A, Kishan Gate, Lodhika G.I.D.C.,
Kalawad Road, Metoda, RAJKOT-360 021, (INDIA).

info@hardikinternational.com

www.hardikinternational.com

☎ : +91 97370 47470/1/2/4 | 📞 : [hardikinternational](https://www.hardikinternational.com)

Follow Us : [!\[\]\(f80254b170d0ecdc443847276e625120_img.jpg\)](#) [!\[\]\(f088b1629d3cb765ad7e1a33317c2efd_img.jpg\)](#) [!\[\]\(3f0caf7929933096753ddd9362ac162d_img.jpg\)](#) [!\[\]\(51dcdb92e2aaf6f9e44bb6dfc26ecefa_img.jpg\)](#) [!\[\]\(9871cf0a3c95e5b4689c992ad0529c32_img.jpg\)](#)