



ELBOW CONNECTOR

QUICK RELEASE

Quick release design allows patients to disconnect and connect to the tubing with minimum effort, increasing patient mobility.

MONITORING

Two standard 4.3mm ports are located on the connector to fulfill the demands of patient monitoring.

DOUBLE ROTATION

360 degrees rotatable connection between mask and connector, and between connector and tubing minimizes interruptions to patients.

EXHALATION VENT

Well positioned exhalation vent allows CO2 to be released easily with less effort for patients.



FRAME & MASK

OPTIMAL FIT

The reinforcing frame incorporates with the silicone mask to prevent deformation. The softtouch cushion helps seal the mask without irritating the face.

The connection between elbow connector and frame can be disassembled easily to allow effective cleansing.

SAFETY

Safe vent design ensures the safety of patients.

High quality medical grade materials are autoclavable to 121 °C.



HEADGEAR

SIMPLE CLIP

The clip-in design allows patients to lock the headgear into the frame and remove it efficiently.

The adjustable strap made with flexible fabric is adaptable to various patients. The permeable mesh on the backside of headgear increases patient comfort.

ORDER INFORMATION

| REF | DESCRIPTION | SIZE | PACKAGE |
|------|----------------|------|-------------------------|
| 1601 | CPAP Face Mask | L | 1set/box, 20sets/carton |
| 1602 | CPAP Face Mask | M | 1set/box, 20sets/carton |
| 1603 | CPAP Face Mask | S | 1set/box, 20sets/carton |