



To whom it may concern

Declaration

We hereby declare that PRO300 Automatic Vacuum Tissue Processor presents at least 3 validated low-temperature processing with standard reagents that allow the processing of biopsies:

- 1-3mm (at room temperature) in a maximum of 2 hours,
- 3 mm (at a maximum temperature of 35oC) in a maximum of 5 hours,
- up to 5 mm (at a maximum temperature of 35oC) in a maximum of 7 hours

Thanking you.

A handwritten signature in blue ink, appearing to read 'Fabrizio Illuminati', is written over a faint blue stamp that reads 'Histo-Line Laboratories S.r.l.'.

Mr. Fabrizio Illuminati
Managing Director

06/09/2023