

ORDIN DE PLATA NR.: 1672

TIP.DOC. 1

DATA EMITERII:18 noiembrie 2022

PLATITI: 76000-00  
ani

LEI: Saptezeci si Sase Mii lei 00 b

PLATITOR: (R) "BIOSISTEM  
MLD" S.R.L.

CONTUL DE PLATI/CODUL IBAN  
MD95ML000000002251429243  
CODUL FISCAL :1010600028048 /

PRESTATORUL PLATITOR

CODUL BANCII:

BC"Moldindconbank"S.A. fil."Invest" Chisinau

:MOLDMD2X329:

BENEFICIAR (R) Centrul pen  
tru achizitii publice central  
izate in sanatate

CONTUL DE PLATI/CODUL IBAN  
MD23TRPCCC518430B01859AA  
CODUL FISCAL :1016601000212 /

PRESTATORUL BENEFICIAR

CODUL BANCII:

Ministerul Finantelor - Trezoreria de Stat

:TREZMD2X

DESTINATIA PLATII:/P102/76000,00 Pentru :  
garantia pentru oferta la procedura de a:  
chizi?ie publica nr. ocds-b3wdp1-MD-1666:  
167332209 din 19.11.2022

TIPUL TRANSFERULUI  
NORMAL/URGENT :N:

L.S.

CODUL TRANZACTIEI:101:

DATA PRIMIRII:18/11/2022

SEMNAURILE

DATA EXECUTARII:

EMITENTULUI

CONDUCATOR:Web Poiata Vitalie

MIIGYwYJKoZiHvcNAQcCoIIGVDCCBLACAQExCzAJBgUrDgMCGGUAMAsGCSqGSIB3:  
DQEHAaCCBGwggRSMIIDUKADAgECAhNHAACjbi1rgFksQ0G4AAAAAKNuMA0GCSqG:  
SIB3DQEBCwUAMCIXIDAeBgNVBAMTF0NFULQxLUNBLU1vbGRpbmRjb25iYW5rMB4X:  
DTIxMDEyODExMzgwNVVoXDTI0MDEyODExNDgwNVowZDZ4CzAJBgNVBAYTAk1EMRAw:  
YDVQKQIEwdNb2xkb3ZHMREwDwYDVQKQHEwhDaGlzaw5hdTEWMBQGA1UEChMNQmLv

(semnatura electronica)

CONTABIL-SEF:Web Nasedchin Alexandr

MIIGZwYJKoZiHvcNAQcCoIIGWDCCBLQCAQExCzAJBgUrDgMCGGUAMAsGCSqGSIB3:  
DQEHAaCCBHAwggRSMIIDVKADAgECAhNHAACjcahRKqbJeg8QAAAAAKNxMA0GCSqG:  
SIB3DQEBCwUAMCIXIDAeBgNVBAMTF0NFULQxLUNBLU1vbGRpbmRjb25iYW5rMB4X:  
DTIxMDEyODExMzkwOFoXDTI0MDEyODExNDkxOFowgaMxCzAJBgNVBAYTAk1EMRAw:  
YDVQKQIEwdNb2xkb3ZHMREwDwYDVQKQHEwhDaGlzaw5hdTEWMBQGA1UEChMNQmLv

L.S.

(semnatura electronica)

CONDUCATOR:

(semnatura manuala)

CONTABIL-SEF:

(semnatura manuala)

SEMNAURA PRESTATORUL

L.S.

MOTIVUL REFUZULUI

L.S.

**CERTIFICAT**  
**privind lipsa sau existența restanțelor față de bugetul public național**

Nr.  
№ A2223060

din  
от 09.11.2022

**1. Destinația / Назначение**

Pentru participarea la proceduri de achiziții publice

**2. Date despre contribuabil / Информация о налогоплательщике**

|                                                                                               |                                                                                    |
|-----------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| Denumirea<br>Наименование                                                                     | Codul fiscal / Numărul de identificare<br>Фискальный код / Идентификационный номер |
| BIOSISTEM MLD S.R.L.                                                                          | 1010600028048                                                                      |
| Adresa sediului de bază (strada, numărul)<br>Адрес основного месторасположения (улица, номер) | Codul - Denumirea localității<br>Код - Наименование населенного пункта             |
| Albisoara nr.16 bl.1 of.7                                                                     | 0150-SEC.RISCANI                                                                   |

**3. Atestarea lipsei sau existenței restanțelor conform datelor Sistemului Informațional Automatizat /**  
Подтверждение отсутствия или наличия недоимки согласно данных Информационной автоматизированной системы

La data emiterii prezentului certificat restanța față de bugetul public național constituie/ На дату выдачи данной справки недоимка перед национальным публичным бюджетом составляет:  
**0,00 lei/лей.**

**4. Valabil până la / Действителен до 24.11.2022**

**5. Autentificarea Serviciului Fiscal de Stat / Подтверждение Государственной налоговой службы**

Digitally signed by Catană Carolina  
Date: 2022.11.09 13:25:20 EET  
Reason: MoldSign Signature  
Location: Moldova



STOICOV Ana

Funcția/Должность

Semnătura/Подпись

Numele și prenumele/Фамилия и имя

L.Ș/M.П.

GOJAN Claudia

Executor:

Numele și prenumele/Фамилия и имя

Este extras din Sistemul Informațional al SFS SIA „Contul curent al contribuabilului”// 09.11.2022 ora 10:46:51  
cu aplicarea prevederilor pct. 82-83 Ordin IFPS nr.400 din 14.03.2014 (Monitorul Oficial 72-77/399, 28.03.2014)

NOTA (0,00)



# BC "MOLDINDCONBANK" S.A.

## Filiala "Invest"

Republica Moldova, MD-2068  
mun. Chișinău, bd. Moscovei, 14/1  
Tel. : (373-22) 43-44-81, 43-46-24  
Fax : (373-22) 43-44-22  
cod: MOLDMD2X329

Data 14. IAN. 2016  
Nr. 03/2 - 19/23

Республика Молдова, MD-2068  
мун. Кишинэу, бул. Московской, 14/1  
Тел. : (373-22) 43-44-81, 43-46-24  
Факс : (373-22) 43-44-22  
код: MOLDMD2X329

Filiala „Invest” BC „Moldindconbank” SA confirmă existența contului curent în moneda națională al **“BIOSISTEM MLD” S.R.L. (c/f 1010600028048)**, cu **IBAN MD95ML000000002251429243**.

Codul băncii MOLDMD2X329.

Director

Nina Turcan

Director financiar



Nina Balmuş

Ex. Diana Brinza  
Tel. 43-45-96



REPUBLICA



MOLDOVA

# CERTIFICAT DE ÎNREGISTRARE

**Societatea cu Răspundere Limitată "BIOSISTEM MLD"**  
— ESTE ÎNREGISTRATĂ LA CAMERA ÎNREGISTRĂRII DE STAT —

*Numărul de identificare de stat - codul fiscal*  
**1010600028048**

*Data înregistrării*

**12.08.2010**

*Data eliberării*

**12.08.2010**

**Svirepova Ludmila, registrator**

*Funcția, numele, prenumele persoanei  
care a eliberat certificatul*

*A. Sizy*  
semnătura

**MD 0101250**







**I.P. "AGENȚIA SERVICII PUBLICE"**

Departamentul înregistrare și licențiere a unităților de drept

**EXTRAS**

**din Registrul de stat al persoanelor juridice**

**nr. 8506 din 28.04.2021**

Denumirea completă: **Societatea cu Răspundere Limitată «BIOSISTEM MLD».**

Denumirea prescurtată: **«BIOSISTEM MLD» S.R.L.**

Forma juridică de organizare: **Societate cu Răspundere Limitată.**

Numărul de identificare de stat și codul fiscal: **1010600028048.**

Data înregistrării de stat: **12.08.2010.**

Sediul: **MD-2001, str. Albișoara, 16/1, ap.(of.) 7, mun. Chișinău, Republica Moldova.**

Obiectul principal de activitate:

- 1 Activitatea farmaceutică;**
- 2 Importul, fabricarea, comercializarea, asistența tehnică și (sau) reparația dispozitivelor medicale și (sau) a opticii;**
- 3 Acordarea asistenței medicale de către instituțiile medico-sanitare private;**
- 4 Comerțul cu ridicata al calculatoarelor, echipamentelor periferice și software-ului;**
- 5 Întreținerea și repararea mașinilor de birou și a tehnicii de calcul;**
- 6 Consultații în domeniul sistemelor de calcul.**

Capitalul social: **5400 lei.**

**Administrator: POIATA VITALIE,**

Asociați:

- 1. POIATA VITALIE 33,40 %**
- 2. NASEDCHIN ALEXANDR 33,30 %**
- 3. KOJEVNIKOV DMITRII 33,30 %.**

Prezentul extras este eliberat în temeiul art. 34 al Legii nr. 220-XVI din 19 octombrie 2007 privind înregistrarea de stat a persoanelor juridice și a întreprinzătorilor individuali și confirmă datele din Registrul de stat la data de: 28.04.2021.

Specialist coordonator  
tel. 022-207-840



**Lazari Aliona**



**EB 0358735**



## **Lista fondatorilor Biosistem-mld SRL**

| <b>Nr.</b> | <b>Nume, Prenume</b>      | <b>IDNP</b>          |
|------------|---------------------------|----------------------|
| <b>1.</b>  | <b>Vitalie Poiata</b>     | <b>0983103892591</b> |
| <b>2.</b>  | <b>Alexandr Nasedchin</b> | <b>2002001070747</b> |
| <b>3.</b>  | <b>Dmitrii Kojevnikov</b> | <b>0972305012362</b> |

SITUAȚIILE FINANCIARE  
pentru perioada 01.01.2021 - 31.12.2021

Entitatea: BIOSISTEM MLD S.R.L.  
Cod CUIO: 40717392  
Cod IDNO: 1010600028048

Sediul:  
MD:  
Raionul(municipiul): 106, DDF RISCANI  
Cod CUATM: 0150, SEC.RISCANI  
Strada: SECTORUL RISCANI STR.Albisoara nr.16 bl.1 of.7

Activitatea principală: G4646, Comert cu ridicata al produselor farmaceutice  
Forma de proprietate: 16, Proprietate colectivă  
Forma organizatorico-juridică: 530, Societăți cu răspundere limitată

Date de contact:  
Telefon: +37322808719  
WEB:  
E-mail: zmii13@mail.ru  
Numele și coordonatele al contabilului-șef: DI (dna) Tel.

Numărul mediu al salariaților în perioada de gestiune: 3 persoane.

Persoanele responsabile de semnarea situațiilor financiare\* Nasedchin Alexandr

Unitatea de măsură: leu

BILANȚUL

Anexa 1

| Nr. cpt. | Indicatori                                                           | Cod rd. | Sold la                         |                                 |
|----------|----------------------------------------------------------------------|---------|---------------------------------|---------------------------------|
|          |                                                                      |         | Începutul perioadei de gestiune | Sfârșitul perioadei de gestiune |
| 1        | 2                                                                    | 3       | 4                               | 5                               |
|          | A C T I V                                                            |         |                                 |                                 |
|          | ACTIVE IMOBILIZATE                                                   |         |                                 |                                 |
|          | I. Imobilizări necorporale                                           |         |                                 |                                 |
|          | 1. Imobilizări necorporale în curs de execuție                       | 010     |                                 |                                 |
|          | 2. Imobilizări necorporale în exploatare, total                      | 020     |                                 |                                 |
|          | din care:                                                            | 021     |                                 |                                 |
|          | 2.1. concesiuni, licențe și mărci                                    |         |                                 |                                 |
|          | 2.2. drepturi de autor și titluri de protecție                       | 022     |                                 |                                 |
|          | 2.3. programe informatice                                            | 023     |                                 |                                 |
|          | 2.4. alte imobilizări necorporale                                    | 024     |                                 |                                 |
|          | 3. Fond comercial                                                    | 030     |                                 |                                 |
|          | 4. Avansuri acordate pentru imobilizări necorporale                  | 040     |                                 |                                 |
|          | Total imobilizări necorporale<br>(rd.010 + rd.020 + rd.030 + rd.040) | 050     |                                 |                                 |
|          | II. Imobilizări corporale                                            |         |                                 |                                 |
|          | 1. Imobilizări corporale în curs de execuție                         | 060     |                                 |                                 |
|          | 2. Terenuri                                                          | 070     |                                 |                                 |
|          | 3. Mijloace fixe, total                                              | 080     | 2793637                         | 3559998                         |
|          | din care:                                                            | 081     |                                 |                                 |
|          | 3.1. clădiri                                                         |         |                                 |                                 |
|          | 3.2. construcții speciale                                            | 082     |                                 |                                 |
|          | 3.3. mașini, utilaje și instalații tehnice                           | 083     | 2791637                         | 3533108                         |
|          | 3.4. mijloace de transport                                           | 084     |                                 |                                 |

|    |                                                                                                                          |     |         |         |
|----|--------------------------------------------------------------------------------------------------------------------------|-----|---------|---------|
| A. | 3.5. inventar și mobilier                                                                                                | 085 |         | 26890   |
|    | 3.6. alte mijloace fixe                                                                                                  | 086 | 2000    |         |
|    | 4. Resurse minerale                                                                                                      | 090 |         |         |
|    | 5. Active biologice imobilizate                                                                                          | 100 |         |         |
|    | 6. Investiții imobiliare                                                                                                 | 110 |         |         |
|    | 7. Avansuri acordate pentru imobilizări corporale                                                                        | 120 |         | 1162136 |
|    | <b>Total imobilizări corporale</b><br>(rd.060 + rd.070 + rd.080 + rd.090 + rd.100 + rd.110 + rd.120)                     | 130 | 2793637 | 4722134 |
|    | <b>III. Investiții financiare pe termen lung</b>                                                                         |     |         |         |
|    | 1. Investiții financiare pe termen lung în părți neafiliate                                                              | 140 |         |         |
|    | 2. Investiții financiare pe termen lung în părți afiliate, total                                                         | 150 |         |         |
|    | din care:                                                                                                                |     |         |         |
|    | 2.1. acțiuni și cote de participație deținute în părțile afiliate                                                        | 151 |         |         |
|    | 2.2 împrumuturi acordate părților afiliate                                                                               | 152 |         |         |
|    | 2.3 împrumuturi acordate aferente intereselor de participare                                                             | 153 |         |         |
|    | 2.4 alte investiții financiare                                                                                           | 154 |         |         |
|    | <b>Total investiții financiare pe termen lung</b><br>(rd.140 + rd.150)                                                   | 160 |         |         |
|    | <b>IV. Creanțe pe termen lung și alte active imobilizate</b>                                                             |     |         |         |
|    | 1. Creanțe comerciale pe termen lung                                                                                     | 170 |         |         |
|    | 2. Creanțe ale părților afiliate pe termen lung                                                                          | 180 |         |         |
|    | inclusiv: creanțe aferente intereselor de participare                                                                    | 181 |         |         |
|    | 3. Alte creanțe pe termen lung                                                                                           | 190 |         |         |
|    | 4. Cheltuieli anticipate pe termen lung                                                                                  | 200 |         |         |
|    | 5. Alte active imobilizate                                                                                               | 210 |         |         |
|    | <b>Total creanțe pe termen lung și alte active imobilizate</b><br>(rd.170 + rd.180 + rd.190 + rd.200 + rd.210)           | 220 |         |         |
|    | <b>TOTAL ACTIVE IMOBILIZATE</b><br>(rd.050 + rd.130 + rd.160 + rd.220)                                                   | 230 | 2793637 | 4722134 |
| B. | <b>ACTIVE CIRCULANTE</b>                                                                                                 |     |         |         |
|    | <b>I. Stocuri</b>                                                                                                        |     |         |         |
|    | 1. Materiale și obiecte de mică valoare și scurtă durată                                                                 | 240 | 51978   | 5346    |
|    | 2. Active biologice circulante                                                                                           | 250 |         |         |
|    | 3. Producția în curs de execuție                                                                                         | 260 |         |         |
|    | 4. Produse și mărfuri                                                                                                    | 270 | 7221203 | 9147976 |
|    | 5. Avansuri acordate pentru stocuri                                                                                      | 280 |         |         |
|    | <b>Total stocuri</b><br>(rd.240 + rd.250 + rd.260 + rd.270 + rd.280)                                                     | 290 | 7273181 | 9153322 |
|    | <b>II. Creanțe curente și alte active circulante</b>                                                                     |     |         |         |
|    | 1. Creanțe comerciale curente                                                                                            | 300 | 3912218 | 2182471 |
|    | 2. Creanțe ale părților afiliate curente                                                                                 | 310 |         |         |
|    | inclusiv: creanțe aferente intereselor de participare                                                                    | 311 |         |         |
|    | 3. Creanțe ale bugetului                                                                                                 | 320 | 74631   | 208171  |
|    | 4. Creanțele ale personalului                                                                                            | 330 |         |         |
|    | 5. Alte creanțe curente                                                                                                  | 340 |         |         |
|    | 6. Cheltuieli anticipate curente                                                                                         | 350 | 2       |         |
|    | 7. Alte active circulante                                                                                                | 360 | 5756117 | 1608597 |
|    | <b>Total creanțe curente și alte active circulante</b><br>(rd.300 + rd.310 + rd.320 + rd.330 + rd.340 + rd.350 + rd.360) | 370 | 9742968 | 3999239 |
|    | <b>III. Investiții financiare curente</b>                                                                                |     |         |         |
|    | 1. Investiții financiare curente în părți neafiliate                                                                     | 380 |         |         |
|    | 2. Investiții financiare curente în părți afiliate, total                                                                | 390 |         |         |
|    | din care:                                                                                                                |     |         |         |
|    | 2.1. acțiuni și cote de participație deținute în părțile afiliate                                                        | 391 |         |         |
|    | 2.2. împrumuturi acordate părților afiliate                                                                              | 392 |         |         |
|    | 2.3. împrumuturi acordate aferente intereselor de participare                                                            | 393 |         |         |



|    |                                                                                                       |     |          |          |
|----|-------------------------------------------------------------------------------------------------------|-----|----------|----------|
|    | 2.4. alte investiții financiare în părți afiliate                                                     | 394 |          |          |
|    | <b>Total investiții financiare curente</b><br>(rd.380 + rd.390)                                       | 400 |          |          |
|    | <b>IV. Numerar și documente bănești</b>                                                               | 410 | 3942779  | 9861933  |
|    | <b>TOTAL ACTIVE CIRCULANTE</b><br>(rd.290 + rd.370 + rd.400 + rd.410)                                 | 420 | 20958928 | 23014494 |
|    | <b>TOTAL ACTIVE</b><br>(rd.230 + rd.420)                                                              | 430 | 23752565 | 27736628 |
|    | <b>P A S I V</b>                                                                                      |     |          |          |
| C. | <b>CAPITAL PROPRIU</b>                                                                                |     |          |          |
|    | <b>I. Capital social și neînregistrat</b>                                                             |     |          |          |
|    | 1. Capital social                                                                                     | 440 | 5400     | 5400     |
|    | 2. Capital nevărsat                                                                                   | 450 | ( )      | ( )      |
|    | 3. Capital neînregistrat                                                                              | 460 |          |          |
|    | 4. Capital retras                                                                                     | 470 | ( )      | ( )      |
|    | 5. Patrimoniul primit de la stat cu drept de proprietate                                              | 480 |          |          |
|    | <b>Total capital social și neînregistrat</b><br>(rd.440 + rd.450 + rd.460 + rd.470 + rd.480)          | 490 | 5400     | 5400     |
|    | <b>II. Prime de capital</b>                                                                           | 500 |          |          |
|    | <b>III. Rezerve</b>                                                                                   |     |          |          |
|    | 1. Capital de rezervă                                                                                 | 510 |          |          |
|    | 2. Rezerve statutare                                                                                  | 520 |          |          |
|    | 3. Alte rezerve                                                                                       | 530 |          |          |
|    | <b>Total rezerve</b><br>(rd.510 + rd.520 + rd.530)                                                    | 540 |          |          |
|    | <b>IV. Profit (pierdere)</b>                                                                          |     |          |          |
|    | 1. Corecții ale rezultatelor anilor precedenți                                                        | 550 | X        |          |
|    | 2. Profit nerepartizat (pierdere neacoperită) al anilor precedenți                                    | 560 | 20060126 | 16230339 |
|    | 3. Profit net (pierdere netă) al perioadei de gestiune                                                | 570 | X        | 10403995 |
|    | 4. Profit utilizat al perioadei de gestiune                                                           | 580 | X        | ( )      |
|    | <b>Total profit (pierdere)</b><br>(rd.550 + rd.560 + rd.570 + rd.580)                                 | 590 | 20060126 | 26634334 |
|    | <b>V. Rezerve din reevaluare</b>                                                                      | 600 |          |          |
|    | <b>VI. Alte elemente de capital propriu</b>                                                           | 610 |          |          |
|    | <b>TOTAL CAPITAL PROPRIU</b><br>(rd.490 + rd.500 + rd.540 + rd.590 + rd.600 + rd.610)                 | 620 | 20065526 | 26639734 |
| D. | <b>DATORII PE TERMEN LUNG</b>                                                                         |     |          |          |
|    | 1. Credite bancare pe termen lung                                                                     | 630 |          |          |
|    | 2. Împrumuturi pe termen lung                                                                         | 640 |          |          |
|    | din care:                                                                                             | 641 |          |          |
|    | 2.1. împrumuturi din emisiunea de obligațiuni                                                         |     |          |          |
|    | inclusiv: împrumuturi din emisiunea de obligațiuni convertibile                                       | 642 |          |          |
|    | 2.2. alte împrumuturi pe termen lung                                                                  | 643 |          |          |
|    | 3. Datorii comerciale pe termen lung                                                                  | 650 |          |          |
|    | 4. Datorii față de părțile afiliate pe termen lung                                                    | 660 |          |          |
|    | inclusiv: datorii aferente intereselor de participare                                                 | 661 |          |          |
|    | 5. Avansuri primite pe termen lung                                                                    | 670 |          |          |
|    | 6. Venituri anticipate pe termen lung                                                                 | 680 |          |          |
|    | 7. Alte datorii pe termen lung                                                                        | 690 |          |          |
|    | <b>TOTAL DATORII PE TERMEN LUNG</b><br>(rd.630 + rd.640 + rd.650 + rd.660 + rd.670 + rd.680 + rd.690) | 700 |          |          |
|    | <b>DATORII CURENTE</b>                                                                                |     |          |          |
|    | 1. Credite bancare pe termen scurt                                                                    | 710 |          |          |
|    | 2. Împrumuturi pe termen scurt, total                                                                 | 720 |          |          |

|    |                                                                                                                                    |     |          |          |
|----|------------------------------------------------------------------------------------------------------------------------------------|-----|----------|----------|
| E. | din care:                                                                                                                          | 721 |          |          |
|    | 2.1. împrumuturi din emisiunea de obligațiuni                                                                                      |     |          |          |
|    | inclusiv: împrumuturi din emisiunea de obligațiuni convertibile                                                                    | 722 |          |          |
|    | 2.2. alte împrumuturi pe termen scurt                                                                                              | 723 |          |          |
|    | 3. Datorii comerciale curente                                                                                                      | 730 | 3252667  | 343711   |
|    | 4. Datorii față de părțile afiliate curente                                                                                        | 740 |          |          |
|    | inclusiv: datorii aferente intereselor de participare                                                                              | 741 |          |          |
|    | 5. Avansuri primite curente                                                                                                        | 750 | 188105   | 355528   |
|    | 6. Datorii față de personal                                                                                                        | 760 | 50       | 350      |
|    | 7. Datorii privind asigurările sociale și medicale                                                                                 | 770 |          |          |
|    | 8. Datorii față de buget                                                                                                           | 780 | 187676   | 150263   |
|    | 9. Datorii față de proprietari                                                                                                     | 790 |          |          |
|    | 10. Venituri anticipate curente                                                                                                    | 800 |          |          |
|    | 11. Alte datorii curente                                                                                                           | 810 | 58541    | 247042   |
|    | <b>TOTAL DATORII CURENTE</b><br>(rd.710 + rd.720 + rd.730 + rd.740 + rd.750 + rd.760 + rd.770 + rd.780 + rd.790 + rd.800 + rd.810) | 820 | 3687039  | 1096894  |
| F. | <b>PROVIZIOANE</b>                                                                                                                 |     |          |          |
|    | 1. Provizioane pentru beneficiile angajaților                                                                                      | 830 |          |          |
|    | 2. Provizioane pentru garanții acordate cumpărătorilor/clientilor                                                                  | 840 |          |          |
|    | 3. Provizioane pentru impozite                                                                                                     | 850 |          |          |
|    | 4. Alte provizioane                                                                                                                | 860 |          |          |
|    | <b>TOTAL PROVIZIOANE</b><br>(rd.830 + rd.840 + rd.850 + rd.860)                                                                    | 870 |          |          |
|    | <b>TOTAL PASIVE</b><br>(rd.620 + rd.700 + rd.820 + rd.870)                                                                         | 880 | 23752565 | 27736628 |

SITUAȚIA DE PROFIT ȘI PIERDERE

de la pînă la

Anexa 2

| Indicatori                                                                                                        | Cod rd. | Perioada de gestiune |          |
|-------------------------------------------------------------------------------------------------------------------|---------|----------------------|----------|
|                                                                                                                   |         | precedenta           | curenta  |
| 1                                                                                                                 | 2       | 3                    | 4        |
| Venituri din vânzări, total                                                                                       | 010     | 25963175             | 38680547 |
| din care:                                                                                                         |         |                      |          |
| venituri din vânzarea produselor și mărfurilor                                                                    | 011     | 25044358             | 37724557 |
| venituri din prestarea serviciilor și executarea lucrărilor                                                       | 012     | 918817               | 951393   |
| venituri din contracte de construcție                                                                             | 013     |                      |          |
| venituri din contracte de leasing                                                                                 | 014     |                      |          |
| venituri din contracte de microfinanțare                                                                          | 015     |                      |          |
| alte venituri din vânzări                                                                                         | 016     |                      | 4597     |
| Costul vânzărilor, total                                                                                          | 020     | 15186814             | 24434231 |
| din care:                                                                                                         |         |                      |          |
| valoarea contabilă a produselor și mărfurilor vândute                                                             | 021     | 15186814             | 24433364 |
| costul serviciilor prestate și lucrărilor executate terților                                                      | 022     |                      |          |
| costuri aferente contractelor de construcție                                                                      | 023     |                      |          |
| costuri aferente contractelor de leasing                                                                          | 024     |                      |          |
| costuri aferente contractelor de microfinanțare                                                                   | 025     |                      |          |
| alte costuri aferente vânzărilor                                                                                  | 026     |                      | 867      |
| <b>Profit brut (pierdere brută)</b> (rd.010 - rd.020)                                                             | 030     | 10776361             | 14246316 |
| Alte venituri din activitatea operațională                                                                        | 040     | 247603               | 5189     |
| Cheltuieli de distribuire                                                                                         | 050     | 19740                | 6076     |
| Cheltuieli administrative                                                                                         | 060     | 1259776              | 1788732  |
| Alte cheltuieli din activitatea operațională                                                                      | 070     | 640169               | 1870642  |
| <b>Rezultatul din activitatea operațională: profit (pierdere)</b><br>(rd.030 + rd.040 - rd.050 - rd.060 - rd.070) | 080     | 9104279              | 10586055 |

|                                                                                                             |     |         |          |
|-------------------------------------------------------------------------------------------------------------|-----|---------|----------|
| Venituri financiare, total                                                                                  | 090 | 519239  | 1517765  |
| din care:                                                                                                   | 091 |         |          |
| venituri din interese de participare                                                                        |     |         |          |
| inclusiv: veniturile obținute de la părțile afiliate                                                        | 092 |         |          |
| venituri din dobânzi                                                                                        | 093 | 25612   | 30619    |
| inclusiv: veniturile obținute de la părțile afiliate                                                        | 094 |         |          |
| venituri din alte investiții financiare pe termen lung                                                      | 095 |         |          |
| inclusiv: veniturile obținute de la părțile afiliate                                                        | 096 |         |          |
| venituri aferente ajustărilor de valoare privind investițiile financiare pe termen lung și curente          | 097 |         |          |
| venituri din ieșirea investițiilor financiare                                                               | 098 |         |          |
| venituri aferente diferențelor de curs valutar și de sumă                                                   | 099 | 493627  | 1487146  |
| Cheltuieli financiare, total                                                                                | 100 | 597528  | 249562   |
| din care:                                                                                                   | 101 |         |          |
| cheltuieli privind dobânzile                                                                                |     |         |          |
| inclusiv: cheltuielile aferente părților afiliate                                                           | 102 |         |          |
| cheltuieli aferente ajustărilor de valoare privind investițiile financiare pe termen lung și curente        | 103 |         |          |
| cheltuieli aferente ieșirii investițiilor financiare                                                        | 104 |         |          |
| cheltuieli aferente diferențelor de curs valutar și de sumă                                                 | 105 | 597528  | 249562   |
| <b>Rezultatul: profit (pierdere) financiar(ă)</b> (rd.090 - rd.100)                                         | 110 | -78289  | 1268203  |
| Venituri cu active imobilizate și excepționale                                                              | 120 |         |          |
| Cheltuieli cu active imobilizate și excepționale                                                            | 130 |         |          |
| <b>Rezultatul din operațiuni cu active imobilizate și excepționale: profit (pierdere)</b> (rd.120 - rd.130) | 140 |         |          |
| <b>Rezultatul din alte activități: profit (pierdere)</b> (rd.110 + rd.140)                                  | 150 | -78289  | 1268203  |
| <b>Profit (pierdere) pînă la impozitare</b> (rd.080 + rd.150)                                               | 160 | 9025990 | 11854258 |
| Cheltuieli privind impozitul pe venit                                                                       | 170 | 1051159 | 1450263  |
| <b>Profit net (pierdere netă) al perioadei de gestiune</b> (rd.160 - rd.170)                                | 180 | 7974831 | 10403995 |

SITUAȚIA MODIFICĂRILOR CAPITALULUI PROPRIU

de la pînă la

Anexa 3

| Nr. d/o | Indicatori                                                                                | Cod rd | Sold la începutul perioadei de gestiune | Majorări | Diminuări | Sold la sfîrșitul perioadei de gestiune |
|---------|-------------------------------------------------------------------------------------------|--------|-----------------------------------------|----------|-----------|-----------------------------------------|
| 1       | 2                                                                                         | 3      | 4                                       | 5        | 6         | 7                                       |
| I.      | <b>Capital social și neînregistrat</b>                                                    |        |                                         |          |           |                                         |
|         | 1. Capital social                                                                         | 010    |                                         |          |           |                                         |
|         | 2. Capital nevărsat                                                                       | 020    | ( )                                     | ( )      | ( )       | ( )                                     |
|         | 3. Capital neînregistrat                                                                  | 030    |                                         |          |           |                                         |
|         | 4. Capital retras                                                                         | 040    | ( )                                     | ( )      | ( )       | ( )                                     |
|         | 5. Patrimoniul primit de la stat cu drept de proprietate                                  | 050    |                                         |          |           |                                         |
|         | <b>Total capital social și neînregistrat</b> (rd.010 + rd.020 + rd.030 + rd.040 + rd.050) | 060    |                                         |          |           |                                         |
| II.     | <b>Prime de capital</b>                                                                   | 070    |                                         |          |           |                                         |
| III.    | <b>Rezerve</b>                                                                            |        |                                         |          |           |                                         |
|         | 1. Capital de rezervă                                                                     | 080    |                                         |          |           |                                         |
|         | 2. Rezerve statutare                                                                      | 090    |                                         |          |           |                                         |
|         | 3. Alte rezerve                                                                           | 100    |                                         |          |           |                                         |
|         | <b>Total rezerve</b> (rd.080 + rd.090 + rd.100)                                           | 110    |                                         |          |           |                                         |
|         | <b>Profit (pierdere)</b>                                                                  |        |                                         |          |           |                                         |
|         | 1. Corecții ale rezultatelor anilor precedenți                                            | 120    | X                                       |          |           |                                         |



|     |                                                                                       |     |   |     |     |
|-----|---------------------------------------------------------------------------------------|-----|---|-----|-----|
| IV. | 2. Profit nerepartizat (pierdere neacoperită) al anilor precedenți                    | 130 |   |     |     |
|     | 3. Profit net (pierdere netă) al perioadei de gestiune                                | 140 | X |     |     |
|     | 4. Profit utilizat al perioadei de gestiune                                           | 150 | X | ( ) | ( ) |
|     | <b>Total profit (pierdere)</b><br>(rd.120 + rd.130 + rd.140 + rd.150)                 | 160 |   |     |     |
| V.  | <b>Rezerve din reevaluare</b>                                                         | 170 |   |     |     |
| VI. | <b>Alte elemente de capital propriu</b>                                               | 180 |   |     |     |
|     | <b>Total capital propriu</b><br>(rd.060 + rd.070 + rd.110 + rd.160 + rd.170 + rd.180) | 190 |   |     |     |

SITUAȚIA FLUXURILOR DE NUMERAR

de la pînă la

Anexa 4

| Indicatori                                                                                                                  | Cod rd | Perioada de gestiune |         |
|-----------------------------------------------------------------------------------------------------------------------------|--------|----------------------|---------|
|                                                                                                                             |        | precedentă           | curentă |
| 1                                                                                                                           | 2      | 3                    | 4       |
| <b>Fluxuri de numerar din activitatea operațională</b>                                                                      |        |                      |         |
| Încasări din vânzări                                                                                                        | 010    |                      |         |
| Plăți pentru stocuri și servicii procurate                                                                                  | 020    |                      |         |
| Plăți către angajați și organe de asigurare socială și medicală                                                             | 030    |                      |         |
| Dobînzî plătite                                                                                                             | 040    |                      |         |
| Plata impozitului pe venit                                                                                                  | 050    |                      |         |
| Alte încasări                                                                                                               | 060    |                      |         |
| Alte plăți                                                                                                                  | 070    |                      |         |
| <b>Fluxul net de numerar din activitatea operațională</b><br>(rd.010 - rd.020 - rd.030 - rd.040 - rd.050 + rd.060 - rd.070) | 080    |                      |         |
| <b>Fluxuri de numerar din activitatea de investiții</b>                                                                     |        |                      |         |
| Încasări din vânzarea activelor imobilizate                                                                                 | 090    |                      |         |
| Plăți aferente intrărilor de active imobilizate                                                                             | 100    |                      |         |
| Dobînzî încasate                                                                                                            | 110    |                      |         |
| Dividende încasate                                                                                                          | 120    |                      |         |
| inclusiv: dividende încasate din străinătate                                                                                | 121    |                      |         |
| Alte încasări (plăți)                                                                                                       | 130    |                      |         |
| <b>Fluxul net de numerar din activitatea de investiții</b><br>(rd.090 - rd.100 + rd.110 + rd.120 ± rd.130)                  | 140    |                      |         |
| <b>Fluxuri de numerar din activitatea financiară</b>                                                                        |        |                      |         |
| Încasări sub formă de credite și împrumuturi                                                                                | 150    |                      |         |
| Plăți aferente rambursării creditelor și împrumuturilor                                                                     | 160    |                      |         |
| Dividende plătite                                                                                                           | 170    |                      |         |
| inclusiv: dividende plătite nerezidenților                                                                                  | 171    |                      |         |
| Încasări din operațiuni de capital                                                                                          | 180    |                      |         |
| Alte încasări (plăți)                                                                                                       | 190    |                      |         |
| <b>Fluxul net de numerar din activitatea financiară</b><br>(rd.150 - rd.160 - rd.170 + rd.180 ± rd.190)                     | 200    |                      |         |
| <b>Fluxul net de numerar total</b><br>(± rd.080 ± rd.140 ± rd.200)                                                          | 210    |                      |         |
| Diferențe de curs valutar favorabile (nefavorabile)                                                                         | 220    |                      |         |
| <b>Sold de numerar la începutul perioadei de gestiune</b>                                                                   | 230    |                      |         |
| <b>Sold de numerar la sfîrșitul perioadei de gestiune</b><br>(± rd.210 ± rd.220 + rd.230)                                   | 240    |                      |         |



[Версия для печати](#)  
[Сохранить](#)

## Расписка 2

Респондент

Фискальный код: 1010600028048, наименование: BIOSISTEM MLD S.R.L.

Предоставил отчёт: RSF1\_21

На фискальный период: A/2021

Дата предоставления: 29.03.2022

Временная метка отчёта зарегистрированного в Информационной Системе НБС : 29.03.2022 17:25:45

National Bureau of Statistics (NBS) received the electronic version of the report, sent by you. The data provided is verified by NBS.

Версия для печати

Сохранить

## Расписка

Респондент

Фискальный код: 1010600028048, наименование: BIOSISTEM MLD S.R.L.

Предоставил отчёт: RSF1\_21

На фискальный период: A/2021

Дата предоставления: 29.03.2022

Временная метка отчёта зарегистрированного в Системе Электронной Отчётности и отправленного в Информационную Систему БНС : 29.03.2022 14:51:06



Staplers | Hemostats | Energy Devices | Hernia Repair Meshes | Ligation Clips | Adhesives | Access Devices  
IUDs | Umbilical Cotton Tape | Bonewax



# Infrastructure

Situated on 180 acres of land and built over a 300,000 sq.ft. area, Meril's ISO 13485 & cGMP certified ultra-modern manufacturing facility is constructed to satisfy the stringent needs of our cardiovascular, orthopedics, diagnostics and endo-surgery businesses, as well as our R&D endeavors.

Fully integrated manufacturing systems ensure backward integration, flawless man-material movement and complete control over processes to eliminate production errors. All manufacturing and sterilization processes are conducted in-house, in addition to the analytical and microbiological QA/QC tests required to meet the world-class production standards.

Meril Park, Vapi



## Talent

Meril's 3200+ personnel are a strong, experienced team comprising designers, engineers, chemists, microbiologists, regulatory affair experts, R&D scientists, clinical affairs experts, legal, finance, sales and marketing professionals with innovative capabilities. The team is continuously striving towards improving and saving lives every day by providing revolutionary solutions for diagnosis, prevention and treatment.

Currently Meril has direct sales & representative offices in India, Germany, Turkey, USA, China and Brazil, and a strong distribution network spread over an additional 102 countries in 6 continents.



# Meril Academy

While innovations and advancements in medical technology play a vital role in contributing to the betterment of a healthy society, a larger role also lies in educating and imparting practical training to doctors, clinicians and other members of the medical fraternity. This initiative will help bridge the gaps between the available technology, its adoptions and its implementation for superior patient care.

Also, to add value to medical professionals and young medical interns, Meril Academy functions as world class institution. Spread over a massive area of 200,000 sq. ft. and surrounded with lush, green natural beauty, Meril Academy is situated at the Meril headquarters in Vapi.





# Meril Endo-Surgery

In line with the objective to design and develop novel, clinically relevant, state-of-art surgical devices, Meril is expanding its endo-surgery franchise. At Meril, we strive to create and develop products through innovative thinking, cutting edge technologies and collaboration with medical professionals and healthcare organizations, with an aim to identify clinical needs and translate them into market-leading solutions.

Meril's vision is to become a 'one-stop' solution for surgeons by providing them with an extensive product portfolio comprising absorbable & non-absorbable surgical sutures, tissue sealants, absorbable hemostats, hernia repair meshes, intra uterine devices, energy devices and mechanical closure devices.

## Contents

|                                  |     |
|----------------------------------|-----|
| Synthetic Absorbable Sutures     |     |
| MITSU                            | 15  |
| MITSU C+                         | 43  |
| MITSU FST                        | 49  |
| MEGASORB                         | 55  |
| MEGASORB PLUS                    | 65  |
| MEGASORB RAPID                   | 67  |
| FILAPRON                         | 71  |
| FILAXYN                          | 79  |
| Natural Absorbable Sutures       |     |
| MERISOFT Plain Catgut            | 89  |
| MERISOFT Chromic Catgut          | 93  |
| Natural Non-Absorbable Sutures   |     |
| FILASILK                         | 99  |
| Synthetic Non-Absorbable Sutures |     |
| MERISTEEL                        | 111 |
| FILAMIDE                         | 115 |
| FILAPROP / P                     | 127 |
| MERICRON XL / P                  | 149 |



# Wound Closure

Wounds which do not heal naturally are treated with wound closure devices. The affected tissues need to be kept in opposition until the healing process provides the wound with proper strength to stand against the stress without any assistance. A wound can be treated with the help of various devices such as sutures, staplers, topical adhesives or clips. Different tissues have their own unique functions and requirements. When one decides on a particular wound closure device, one must take care of the healing needs of the particular tissue which is being treated.

The most common device used for wound closure is suture. Sutures, in general term ‘stitches’, are sterilized surgical threads, which are used to ligate blood vessels or to approximate tissues.

**Sutures are generally categorized as:**

- Absorbable and Non-Absorbable
- Natural and Synthetic
- Monofilament and Braided

Absorbable sutures have a tendency to degrade as the wound heals. Non-Absorbable sutures oppose the body’s efforts to dissolve them. On complete healing of the wound, surgeons remove the non-absorbable sutures. Meril Endo-Surgery has a wide range of synthetic absorbable sutures which are made from Polyglactin 910 (MITSU family), Polydioxanone (FILAXYN), Poliglecaprone 25 (FILAPRON), Polyglycolic Acid (MEGASORB family) as well as non-absorbable sutures which are made from Polypropylene (FILAPROP), Polyamide (FILAMIDE) and Polyester (MERICRON XL). In natural absorbable sutures, Meril portfolio covers Chromic and Plain Catgut (MERISOFT) and in natural non-absorbable we have silk (FILASILK).


For placement of sutures, surgical needles are used, which are made of stainless steel. Different types and sizes of needles are used depending on the area of application. We are dedicated to assist surgeons in tissue repair by giving them a broad gamut of products with our unmatched quality and service.

# Meril Endo-Surgery Needles

In pursuit of offering best quality products, Meril sutures are paired with the superior stainless steel needles for an unparalleled surgical experience.

Meril sutures come with sharper and stronger needles (300 series Stainless Steel) specifically designed to meet the ever growing needs in all the surgical disciplines.

With every needle swaged using Multi Stroke Technology (MST) and tested for a needle pull check, Meril needles are distinguished with a high penetration coating and a square body design for a firm grip.



| Needle Profiles                                                                       | Point type                 | Symbol                                                                                |
|---------------------------------------------------------------------------------------|----------------------------|---------------------------------------------------------------------------------------|
|  | Round bodied/Taper point   |  |
|  | Conventional cutting       |  |
|  | Reverse cutting            |  |
|  | Reverse cutting prime      |  |
|  | Tapercutting               |  |
|  | UltraPoint reverse cutting |  |
|  | UltraPoint spatula curved  |  |
|  | CSU spatula                |  |
|  | SBR spatula                |  |



Needle Shapes



Needle Specifications

| Type       | Taper                                                                                  | Reverse cutting | Diamond | Cutting | Premium cutting | Spatula | Blunt | Taper cut |
|------------|----------------------------------------------------------------------------------------|-----------------|---------|---------|-----------------|---------|-------|-----------|
| Design     |                                                                                        |                 |         |         |                 |         |       |           |
| Point      |                                                                                        |                 |         |         |                 |         |       |           |
| Shape      | Typical Applications                                                                   |                 |         |         |                 |         |       |           |
| Straight   | Skin/subcutaneous layer closure; Purse string                                          |                 |         |         |                 |         |       |           |
| 1/4 Circle | Ophthalmic/Microsurgical procedures                                                    |                 |         |         |                 |         |       |           |
| 3/8 Circle | Skin/subcutaneous layer/abdominal wall closure; Ophthalmic/Vascular/Plastic procedures |                 |         |         |                 |         |       |           |
| 1/2 Circle | Peritoneum closure; Gastrointestinal/Respiratory/Urological procedures                 |                 |         |         |                 |         |       |           |
| 5/8 Circle | Pelvic & Oral cavity; Urological procedures                                            |                 |         |         |                 |         |       |           |
| J-Shape    | Laparotomy closure                                                                     |                 |         |         |                 |         |       |           |
| Ski        | Laparoscopic procedures                                                                |                 |         |         |                 |         |       |           |

Understanding the Suture Codes

Meril code logic represents the necessary information needed to identify the required suture-needle combination. The suture codes are alphanumeric in nature to indicate the suture material, the suture size, and the reference code. The alphabetical suffixes are a key to understand the additional features about the suture and/or the needle.

Material Prefixes

| ABBREVIATION | MEANING                                |
|--------------|----------------------------------------|
| PGN          | Polyglactin 910                        |
| PGP          | Polyglactin 910 (with Chlorhexidine)   |
| PGF          | Polyglactin 910 Fast                   |
| MS           | Polyglycolic Acid                      |
| MSP          | Polyglycolic Acid (with Chlorhexidine) |
| MST          | Polyglycolic Acid (with Triclosan)     |
| PCL          | Poliglecaprone 25                      |
| PDX          | Polydioxanone                          |
| CGL          | Chromic Catgut Loop                    |
| SLK          | Silk                                   |
| SLS          | Silk Non Needled                       |
| NYL          | Polyamide                              |
| NYS          | Polyamide Non Needled                  |
| PPL          | Polypropylene                          |
| PPG          | Polypropylene with Pledgets            |
| STC          | Steel Wire                             |

\*In case of

- Polyester (Mericon XL), codes start with ME6XX, ME6XXX, ME4XXX, MEXXXX; and polyester with pledgets (Mericon XL P), codes start with PMEXXX or MMEXXXX
- Chromic Catgut (Merisoft Chromic), codes start with ME40XX, ME42XX, ME50XX
- Chromic Catgut non needled (Merisoft Chromic), codes start with SS22XX
- Plain Catgut (Merisoft Plain), codes start with ME50XX
- Plain Catgut non needled (Merisoft Plain), codes start with SS20XX

# Suture Sizes

- For sutures having a single-digit size, the size is written as: zero followed by the size digit.  
E.g.: Size 01 (read as size 1 according to USP), Size 03 (read as size 3 according to USP).
- For finer sutures, the size is written as: the size digit followed by a zero.  
E.g.: Size 60 (read as size 6-0 according to USP), Size 70 ( read as size 7-0 according to USP), Size 100 (read as size 10-0 according to USP)


# Suffixes

| ABBREVIATION | MEANING                           |
|--------------|-----------------------------------|
| B/BL         | Black Needle                      |
| BX           | Black Needle; Extra Length Suture |
| DN           | Double Needle                     |
| DNL          | Double Needle; Long Length Suture |
| DS           | Doctor Special Code               |
| G            | Green Suture                      |
| ML           | Multi Strands                     |
| SL           | Short Length Suture               |
| U            | Undyed Suture                     |
| V            | Violet Suture                     |
| W            | White Suture                      |
| X            | Extra Pair of Suture              |
| XF           | Extra Fine Needle                 |
| XL           | Extra Length Suture               |


# Abbreviation for Specification

| ABBREVIATION | MEANING         |
|--------------|-----------------|
| HC           | 1/2 CIRCLE      |
| QC           | 1/4 CIRCLE      |
| CU           | 3/8 CIRCLE      |
| FE           | 5/8 CIRCLE      |
| RC           | REVERSE CUTTING |
| CT           | CUTTING         |
| RB           | ROUND BODIED    |
| TC           | TAPER CUT       |
| ST           | STRAIGHT        |
| SM           | SPATULA         |
| TRO          | TROCAR          |
| BP           | BLUNT POINT     |
| TP           | PLEDGET         |
| VIO          | VIOLET          |
| UND          | UNDYED          |

# Decoding the Codes



# Needle Point Symbol



MITSU<sup>TM</sup>  
Polyglactin 910 Suture

Properties

- 75% of the original is retained until initial 14 Days
- 50% of the original is retained until initial 21 Days
- There is a subsequent loss between four to five weeks post implantation
- Complete absorption usually takes place between 56 to 70 Days

Indications


- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Ophthalmic procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- Smooth knot run-down/Knot slide
- Gentle knot repositioning
- Excellent knot security
- Minimal inflammatory reaction



Product Overview

|                      |                                                                                                         |
|----------------------|---------------------------------------------------------------------------------------------------------|
| Structure            | : Multifilament, Braided                                                                                |
| Colour               | : Dyed, Undyed                                                                                          |
|                      |                    |
| Chemical composition | : Poly (glycolide-co-l-lactide 90/10)                                                                   |
| Coating              | : Poly(glycolide-co-l-lactide 30/70) + Calcium Stearate                                                 |
| Sizes                | : Dyed - USP 8/0 (metric 0.4) - USP 2 (metric 5)<br>: Undyed - USP 6/0 (metric 0.7 ) - USP 2 (metric 5) |
| Type of absorption   | : Hydrolysis                                                                                            |
| Sterilization        | : Ethylene oxide                                                                                        |



Meril Endo-Surgery Product Catalogue | **Wound Closure** 17









Meril Endo-Surgery Product Catalogue | **Wound Closure** 25

Meril Endo-Surgery Product Catalogue | **Wound Closure** 27



Meril Endo-Surgery Product Catalogue | **Wound Closure** 29

Meril Endo-Surgery Product Catalogue | **Wound Closure** 31







Meril Endo-Surgery Product Catalogue | **Wound Closure** 37

| Description                                                 |              |      |                                                                                     |                                                          |           |               |     | Description                                                 |              |      |                                                                                       |                                                           |           |               |     |
|-------------------------------------------------------------|--------------|------|-------------------------------------------------------------------------------------|----------------------------------------------------------|-----------|---------------|-----|-------------------------------------------------------------|--------------|------|---------------------------------------------------------------------------------------|-----------------------------------------------------------|-----------|---------------|-----|
| FG Code                                                     |              | Size | Needle Profile Shape, Type & Length                                                 | Suture Length (cm)                                       | Pack Size | Suture Colour |     | FG Code                                                     |              | Size | Needle Profile Shape, Type & Length                                                   | Suture Length (cm)                                        | Pack Size | Suture Colour |     |
| Description                                                 |              |      |                                                                                     |                                                          |           |               |     | Description                                                 |              |      |                                                                                       |                                                           |           |               |     |
| MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910) |              |      |                                                                                     |                                                          |           |               |     | MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910) |              |      |                                                                                       |                                                           |           |               |     |
| PGLA VIO 6-0 X 70-13MM HC RB                                | PGN600134    | 6-0  |    | 13 mm 1/2 Circle Round Body                              | 70        | 12            | VIO | PGLA VIO 6-0 X 70-17MM HC RB                                | PGN600212V   | 6-0  |    | 17 mm 1/2 Circle Round Body                               | 70        | 12            | VIO |
| PGLA VIO 6-0 X 45-8MM HC SM                                 | PGN600562    | 6-0  |    | 8 mm 1/2 Circle Spatula Micropoint                       | 45        | 12            | VIO | PGLA VIO 7-0 X 30-06MM CU SM                                | PGN70 2561   | 7-0  |    | 6 mm 3/8 Circle Spatula Micropoint                        | 30        | 12            | VIO |
| PGLA UND 6-0 X 45-11MM HC RC                                | PGN600510    | 6-0  |    | 11 mm 1/2 Circle Reverse Cutting                         | 45        | 12            | UND | PGLA VIO 7-0 X 45-7MM HC SM                                 | PGN709565    | 7-0  |    | 7 mm 1/2 Circle Spatula Micropoint                        | 45        | 12            | VIO |
| PGLA UND 6-0 X 45-13MM CU RC                                | PGN600492    | 6-0  |    | 13 mm 3/8 Circle Reverse Cutting                         | 45        | 12            | UND | PGLA VIO 7-0 X 30-6.5MM CU SM DN                            | PGN70 9561DN | 7-0  |    | 6.5 mm 3/8 Circle Spatula Micropoint   Spatula Micropoint | 30        | 12            | VIO |
| PGLA UND 6-0 X 70-17MM HC RB                                | PGN600212    | 6-0  |    | 17 mm 1/2 Circle Round Body                              | 70        | 12            | UND | PGLA VIO 7-0 X 90-6.5MM CU SM DN                            | PGN709561    | 7-0  |    | 6.5 mm 3/8 Circle Spatula Micropoint   Spatula Micropoint | 90        | 12            | VIO |
| PGLA UND 6-0 X 45-13MM CU CT                                | PGN60 2831U  | 6-0  |    | 13 mm 3/8 Circle Conventional Cutting                    | 45        | 12            | UND | PGLA VIO 7-0 X 30-5MM HC SM                                 | PGN700575    | 7-0  |    | 5 mm 1/2 Circle Spatula Micropoint                        | 30        | 12            | VIO |
| PGLA VIO 6-0 X 45-08MM QC SM DN                             | PGN60 9552DN | 6-0  |    | 8 mm 1/4 Circle Spatula Micropoint   Spatula Micropoint  | 45        | 12            | VIO | PGLA VIO 7-0 X 45-7MM HC SM DN                              | PGN709565DN  | 7-0  |    | 7 mm 1/2 Circle Spatula Micropoint   Spatula Micropoint   | 45        | 12            | VIO |
| PGLA VIO 6-0 X 30-10MM FE RB                                | PGN609610    | 6-0  |    | 10 mm 5/8 Circle Round Body                              | 30        | 12            | VIO | PGLA VIO 8-0 X 30-9MM CU RB                                 | PGN80 9577   | 8-0  |    | 9 mm 3/8 Circle Round Body                                | 30        | 12            | VIO |
| PGLA UND 6-0 X 45-8MM QC SM DN                              | PGN609752    | 6-0  |    | 8 mm 1/4 Circle Spatula Micropoint   Spatula Micropoint  | 45        | 12            | UND | PGLA VIO 8-0 X 30-6MM CU SM DN                              | PGN80 2348DN | 8-0  |    | 6 mm 3/8 Circle Spatula Micropoint   Spatula Micropoint   | 30        | 12            | VIO |
| PGLA UND 6-0 X 45-13MM HC RB                                | PGN609981    | 6-0  |    | 13 mm 1/2 Circle Round Body                              | 45        | 12            | UND | PGLA VIO 8-0 X 30-08MM CU RB                                | PGN80 2301   | 8-0  |    | 8 mm 3/8 Circle Round Body                                | 30        | 12            | VIO |
| PGLA VIO 6-0 X 45-08MM QC SM DN C1                          | PGN600563    | 6-0  |   | 8 mm 1/4 Circle Spatula Micropoint   Spatula Micropoint  | 45        | 12            | VIO | PGLA VIO 8-0 X 30-6.5MM CU SM DN                            | PGN80 9560DN | 8-0  |   | 6.5 mm 3/8 Circle Spatula Micropoint   Spatula Micropoint | 30        | 12            | VIO |
| PGLA UND 6-0 X 45-13MM CU CT EX                             | PGN6083365   | 6-0  |  | 13 mm 3/8 Circle Conventional Cutting                    | 45        | 12            | UND | PGLA VIO 8-0 X 90-7.5MM HC SM DN                            | PGN809564    | 8-0  |  | 7.5 mm 1/2 Circle Spatula Micropoint   Spatula Micropoint | 90        | 12            | VIO |
| PGLA VIO 6-0 X 45-12MM CU RC                                | PGN60 2671   | 6-0  |  | 12 mm 3/8 Circle Reverse Cutting                         | 45        | 12            | VIO | PGLA VIO 8-0 X 45-6MM CU SM DN                              | PGN809559    | 8-0  |  | 6 mm 3/8 Circle Spatula Micropoint   Spatula Micropoint   | 45        | 12            | VIO |
| PGLA UND 6-0 X 45-8MM HC SM DN                              | PGN609756DN  | 6-0  |  | 8 mm 1/2 Circle Spatula Micropoint   Spatula Micropoint  | 45        | 12            | UND | PGLA VIO 8-0 X 30-6MM CU SM                                 | PGN809558    | 8-0  |  | 6 mm 3/8 Circle Spatula Micropoint                        | 30        | 12            | VIO |
| PGLA UND 6-0 X 45-11MM CU SM DN                             | PGN609760DN  | 6-0  |  | 11 mm 3/8 Circle Spatula Micropoint   Spatula Micropoint | 45        | 12            | UND | PGLA VIO 8-0 X 45-6MM CU SM                                 | PGN800547    | 8-0  |  | 6 mm 3/8 Circle Spatula Micropoint                        | 45        | 12            | VIO |
| PGLA VIO 6-0 X 45-8MM HC SM DN                              | PGN609566DN  | 6-0  |  | 8 mm 1/2 Circle Spatula Micropoint   Spatula Micropoint  | 45        | 12            | VIO | PGLA VIO 2 X 70-48MM HC RB                                  | PGN028063    | 2    |  | 48 mm 1/2 Circle Round Body Standard Needle               | 70        | 12            | VIO |
| PGLA VIO 6-0 X 45-11MM CU RC DN                             | PGN609537    | 6-0  |  | 11 mm 3/8 Circle Reverse Cutting   Reverse Cutting       | 45        | 12            | VIO | PGLA VIO 0 X 70-37MM HC RB                                  | PGN100324    | 0    |  | 37 mm 1/2 Circle Round Body                               | 70        | 12            | VIO |
| PGLA VIO 6-0 X 45-8MM QC SM DN-3D                           | PGN602670DNZ | 6-0  |  | 8 mm 1/4 Circle Spatula Micropoint   Spatula Micropoint  | 45        | 36            | VIO | PGLA VIO 4-0 X 45 -8MM HC SM DN                             | PGN409467DN  | 4-0  |  | 8 mm 1/2 Circle Spatula Micropoint                        | 45        | 12            | VIO |
| PGLA UND 6-0 X 45-16MM CU RC                                | PGN600499U   | 6-0  |  | 16 mm 3/8 Circle Reverse Cutting                         | 45        | 12            | UND |                                                             |              |      |                                                                                       |                                                           |           |               |     |
| PGLA UND 6-0 X 45-11MM CU RC                                | PGN600489U   | 6-0  |  | 11 mm 3/8 Circle Reverse Cutting                         | 45        | 12            | UND |                                                             |              |      |                                                                                       |                                                           |           |               |     |

38

Wound Closure | Meril Endo-Surgery Product Catalogue

Meril Endo-Surgery Product Catalogue | Wound Closure

39





MITSU C+<sup>TM</sup>

Polyglactin 910 Suture with Chlorhexidine

Properties

- Retention of 75% knot tensile strength after 14 days
- Retention of 50% knot tensile strength after 21 days
- Complete mass absorption within 56 to 70 days

Indications


- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- Smooth knot run-down/Knot slide
- Gentle knot repositioning
- Excellent knot security
- Minimal inflammatory reaction



Product Overview

|                      |                                                                                                |
|----------------------|------------------------------------------------------------------------------------------------|
| Structure            | : Multifilament, Braided                                                                       |
| Colour               | : Dyed, Undyed                                                                                 |
|                      |           |
| Chemical composition | : Poly(glycolide-co-l-lactide 90/10)                                                           |
| Coating              | : Poly(glycolide-co-l-lactide 30/70) + Calcium Stearate<br>+ antibacterial agent Chlorhexidine |
| Sizes                | : USP 5/0 (metric 1) - USP 2 (metric 5)                                                        |
| Type of absorption   | : Hydrolysis                                                                                   |
| Sterilization        | : Ethylene oxide                                                                               |







# MITSU FST™

## Polyglactin 910 Fast Suture

### Properties

- 50% of the original is retained until initial 5 Days
- approximately 10-14 Days post implantation, the original tensile strength is completely lost
- Complete absorption usually takes place by 42 Days

### Indications

- Gynaecology/Obstetric procedures (e.g.: episiotomy repair)
- Urology procedures
- Ophthalmic procedures (e.g.: conjunctiva suturing)
- Odontology procedures (e.g.: oral mucosa)
- Skin closure (intracutaneous, subcutaneous)
- Ligatures






### Advantages

- High initial tensile strength
- Predictable and constant degradation rate
- Excellent knot security
- Excellent handling properties
- Quick mass absorption



### Product Overview

|                      |                                           |
|----------------------|-------------------------------------------|
| Structure            | : Multifilament, Braided                  |
| Colour               | : Undyed                                  |
| Chemical composition | : Poly(glycolide-co-l-lactide 90/10)      |
| Sizes                | : USP 8/0 (metric 0.7) - USP 2 (metric 5) |
| Type of absorption   | : Hydrolysis                              |
| Sterilization        | : Gamma radiation                         |

| Description                                                          | FG Code      | Size |                                                                                     | Needle Profile<br>Shape, Type &<br>Length      | Suture<br>Length<br>(cm) | Pack<br>Size | Suture<br>Colour | Description                                                          | FG Code      | Size |                                                                                       | Needle Profile<br>Shape, Type &<br>Length           | Suture<br>Length<br>(cm) | Pack<br>Size | Suture<br>Colour |
|----------------------------------------------------------------------|--------------|------|-------------------------------------------------------------------------------------|------------------------------------------------|--------------------------|--------------|------------------|----------------------------------------------------------------------|--------------|------|---------------------------------------------------------------------------------------|-----------------------------------------------------|--------------------------|--------------|------------------|
| MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST) |              |      |                                                                                     |                                                |                          |              |                  | MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST) |              |      |                                                                                       |                                                     |                          |              |                  |
| PGLA FAST<br>UND 0 X 110-36MM HC TC                                  | PGF10 2763   | 0    |    | 36 mm 1/2 Circle<br>Taper Cut                  | 110                      | 12           | UND              | PGLA FAST<br>UND 2-0 X 90-36MM HC TC                                 | PGF20 2762   | 2-0  |    | 36 mm 1/2 Circle<br>Taper Cut                       | 90                       | 12           | UND              |
| PGLA FAST<br>UND 0 X 90-40MM HC RB                                   | PGF109977    | 0    |    | 40 mm 1/2 Circle<br>Round Body                 | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 140-36MM<br>HC RC RB DN                       | PGF20 2777DN | 2-0  |    | 36 mm 1/2 Circle<br>Reverse Cutting  <br>Round Body | 140                      | 12           | UND              |
| PGLA FAST<br>UND 0 X 90-48MM HC RB                                   | PGF1016951   | 0    |    | 48 mm 1/2 Circle<br>Round Body                 | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 75-31MM HC RB                                 | PGF209975    | 2-0  |    | 31 mm 1/2 Circle<br>Round Body                      | 75                       | 12           | UND              |
| PGLA FAST<br>UND 0 X 90-36MM HC RB                                   | PGF109036    | 0    |    | 36 mm 1/2 Circle<br>Round Body                 | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 90-48MM HC RB                                 | PGF209978    | 2-0  |    | 48 mm 1/2 Circle<br>Round Body                      | 90                       | 12           | UND              |
| PGLA FAST<br>UND 0 X 60-36MM CU RC                                   | PGF10 2721   | 0    |    | 36 mm 3/8 Circle<br>Reverse Cutting            | 60                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 70-26MM HC RB                                 | PGF204171    | 2-0  |    | 26 mm 1/2 Circle<br>Round Body                      | 70                       | 12           | UND              |
| PGLA FAST<br>UND 0 X 110-40MM HC RB                                  | PGF10 2763RB | 0    |    | 40 mm 1/2 Circle<br>Round Body                 | 110                      | 12           | UND              | PGLA FAST<br>UND 2-0 X 90-36MM HC RB                                 | PGF209451    | 2-0  |    | 36 mm 1/2 Circle<br>Round Body                      | 90                       | 12           | UND              |
| PGLA FAST<br>UND 0 X 90-36MM HC TC                                   | PGF10 9963   | 0    |    | 36 mm 1/2 Circle<br>Taper Cut                  | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 135-36MM<br>HC RB RC DN                       | PGF209321    | 2-0  |   | 36 mm 1/2 Circle<br>Round Body  <br>Reverse Cutting | 135                      | 12           | UND              |
| PGLA FAST<br>UND 0 X 90-30MM CU RC-3D                                | PGF102287Z   | 0    |   | 30 mm 3/8 Circle<br>Reverse Cutting            | 90                       | 36           | UND              | PGLA FAST<br>UND 2-0 X 90-36MM HC RC                                 | PGF209236    | 2-0  |  | 36 mm 1/2 Circle<br>Reverse Cutting                 | 90                       | 12           | UND              |
| PGLA FAST<br>UND 1 X 90-40MM HC TC                                   | PGF016190    | 1    |  | 40 mm 1/2 Circle<br>Taper Cut                  | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 70-19MM CU RC                                 | PGF200428    | 2-0  |  | 19 mm 3/8 Circle<br>Reverse Cutting                 | 70                       | 12           | UND              |
| PGLA FAST<br>UND 1 X 90-40MM HC RB                                   | PGF016191    | 1    |  | 40 mm 1/2 Circle<br>Round Body<br>Heavy Needle | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 60-36MM CU RC                                 | PGF20 2720   | 2-0  |  | 36 mm 3/8 Circle<br>Reverse Cutting                 | 60                       | 12           | UND              |
| PGLA FAST<br>UND 1 X 90-36MM HC RB                                   | PGF016192    | 1    |  | 36 mm 1/2 Circle<br>Round Body                 | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 110-36MM HC TC                                | PGF20 2762XL | 2-0  |  | 36 mm 1/2 Circle<br>Taper Cut                       | 110                      | 12           | UND              |
| PGLA FAST<br>UND 1 X 90-48MM HC RB                                   | PGF016193    | 1    |  | 48 mm 1/2 Circle<br>Round Body                 | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 120-36MM HC TC                                | PGF209948    | 2-0  |  | 36 mm 1/2 Circle<br>Taper Cut                       | 120                      | 12           | UND              |
| PGLA FAST<br>UND 1 X 90-36MM HC TC                                   | PGF019190    | 1    |  | 36 mm 1/2 Circle<br>Taper Cut                  | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 75-26MM CU RC                                 | PGF20 9933   | 2-0  |  | 26 mm 3/8 Circle<br>Reverse Cutting                 | 75                       | 12           | UND              |
| PGLA FAST<br>UND 0 X 90-45MM HC TC                                   | PGF109947    | 0    |  | 45 mm 1/2 Circle<br>Taper Cut                  | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 100-36MM<br>HC TC SP                          | PGF20 2761SP | 2-0  |  | 36 mm 1/2 Circle<br>Taper Cut                       | 100                      | 12           | UND              |
| PGLA FAST<br>UND 2-0 X 90-40MM HC TC                                 | PGF209946    | 2-0  |  | 40 mm 1/2 Circle<br>Taper Cut                  | 90                       | 12           | UND              | PGLA FAST<br>UND 2-0 X 90-36MM<br>HC TC SP                           | PGF20 2762SP | 2-0  |  | 36 mm 1/2 Circle<br>Taper Cut                       | 90                       | 12           | UND              |
| PGLA FAST<br>UND 2-0 X 90-35MM CU RC                                 | PGF209938    | 2-0  |  | 35 mm 3/8 Circle<br>Reverse Cutting            | 90                       | 12           | UND              |                                                                      |              |      |                                                                                       |                                                     |                          |              |                  |
| PGLA FAST<br>UND 2-0 X 100-36MM HC TC                                | PGF20 2761   | 2-0  |  | 36 mm 1/2 Circle<br>Taper Cut                  | 100                      | 12           | UND              |                                                                      |              |      |                                                                                       |                                                     |                          |              |                  |



| Description                                                          | FG Code      | Size |                                                                                     | Needle Profile Shape, Type & Length   | Suture Length (cm) | Pack Size | Suture Colour | Description                                                          | FG Code     | Size |                                                                                       | Needle Profile Shape, Type & Length   | Suture Length (cm) | Pack Size | Suture Colour |
|----------------------------------------------------------------------|--------------|------|-------------------------------------------------------------------------------------|---------------------------------------|--------------------|-----------|---------------|----------------------------------------------------------------------|-------------|------|---------------------------------------------------------------------------------------|---------------------------------------|--------------------|-----------|---------------|
| MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST) |              |      |                                                                                     |                                       |                    |           |               | MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST) |             |      |                                                                                       |                                       |                    |           |               |
| PGLA FAST<br>UND 2-0 X 90-26MM CU RC                                 | PGF209941    | 2-0  |    | 26 mm 3/8 Circle Reverse Cutting      | 90                 | 12        | UND           | PGLA FAST<br>UND 3-0 X 75-36MM CU RC                                 | PGF309961   | 3-0  |    | 36 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | UND           |
| PGLA FAST<br>UND 2-0 X 90-26MM CU RC-3D                              | PGF209941Z   | 2-0  |    | 26 mm 3/8 Circle Reverse Cutting      | 90                 | 36        | UND           | PGLA FAST<br>UND 3-0 X 70-26MM HC RB EX                              | PGF304161   | 3-0  |    | 26 mm 1/2 Circle Round Body           | 70                 | 12        | UND           |
| PGLA FAST<br>UND 2-0 X 75-24MM CU RC-3D                              | PGF209275Z   | 2-0  |    | 24 mm 3/8 Circle Reverse Cutting      | 75                 | 36        | UND           | PGLA FAST<br>UND 3-0 X 75-25MM FE CT                                 | PGF309930   | 3-0  |    | 25 mm 5/8 Circle Conventional Cutting | 75                 | 12        | UND           |
| PGLA FAST<br>UND 2-0 X 75-60MM ST CT                                 | PGF209944    | 2-0  |    | 60 mm Straight Conventional Cutting   | 75                 | 12        | UND           | PGLA FAST<br>UND 3-0 X 70-26MM HC RB-3D                              | PGF304161Z  | 3-0  |    | 26 mm 1/2 Circle Round Body           | 70                 | 36        | UND           |
| PGLA FAST<br>UND 3-0 X 75-16MM CU CT                                 | PGF30 9919   | 3-0  |    | 16 mm 3/8 Circle Conventional Cutting | 75                 | 12        | UND           | PGLA FAST<br>UND 3-0 X 75-19MM CU CT                                 | PGF309319   | 3-0  |    | 19 mm 3/8 Circle Conventional Cutting | 75                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 90-60MM ST CT                                 | PGF309943    | 3-0  |    | 60 mm Straight Conventional Cutting   | 90                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 75-19MM CU RC                                 | PGF409922   | 4-0  |    | 19 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 90-25MM FE RB                                 | PGF309929    | 3-0  |    | 25 mm 5/8 Circle Round Body           | 90                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 75-17MM HC RB                                 | PGF409970   | 4-0  |    | 17 mm 1/2 Circle Round Body           | 75                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 90-35MM CU RC                                 | PGF309937    | 3-0  |    | 35 mm 3/8 Circle Reverse Cutting      | 90                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 75-22MM HC RB                                 | PGF409971   | 4-0  |    | 22 mm 1/2 Circle Round Body           | 75                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 70-26MM CU RC                                 | PGF30 2732   | 3-0  |    | 26 mm 3/8 Circle Reverse Cutting      | 70                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 45-13MM CU CT                                 | PGF408351   | 4-0  |    | 13 mm 3/8 Circle Conventional Cutting | 45                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 45-26MM CU RC                                 | PGF30 2732SL | 3-0  |    | 26 mm 3/8 Circle Reverse Cutting      | 45                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 75-13MM CU RC                                 | PGF409913   | 4-0  |    | 13 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 75-19MM CU RC                                 | PGF309923    | 3-0  |  | 19 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 45-16MM CU CT                                 | PGF40 2718  | 4-0  |  | 16 mm 3/8 Circle Conventional Cutting | 45                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 75-22MM HC RB                                 | PGF309974    | 3-0  |  | 22 mm 1/2 Circle Round Body           | 75                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 75-16MM CU CT                                 | PGF40 9918  | 4-0  |  | 16 mm 3/8 Circle Conventional Cutting | 75                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 90-36MM HC RB                                 | PGF309441    | 3-0  |  | 36 mm 1/2 Circle Round Body           | 90                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 45-13MM ST CT                                 | PGF409909SL | 4-0  |  | 13 mm Straight Conventional Cutting   | 45                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 45-22MM HC CT                                 | PGF30 9935   | 3-0  |  | 22 mm 1/2 Circle Conventional Cutting | 45                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 45-16MM CU CT EX                              | PGF408451   | 4-0  |  | 16 mm 3/8 Circle Conventional Cutting | 45                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 70-18MM HC RB                                 | PGF300318    | 3-0  |  | 18 mm 1/2 Circle Round Body           | 70                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 45-22MM CU CT                                 | PGF409926   | 4-0  |  | 22 mm 3/8 Circle Conventional Cutting | 45                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 45-16MM CU CT                                 | PGF30 2719   | 3-0  |  | 16 mm 3/8 Circle Conventional Cutting | 45                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 45-45MM ST CT                                 | PGF409939   | 4-0  |  | 45 mm Straight Conventional Cutting   | 45                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 70-22MM HC CT                                 | PGF30 2735   | 3-0  |  | 22 mm 1/2 Circle Conventional Cutting | 70                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 45-13MM CU RC                                 | PGF404940   | 4-0  |  | 13 mm 3/8 Circle Reverse Cutting      | 45                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 70-26MM HC RB                                 | PGF30 2764   | 3-0  |  | 26 mm 1/2 Circle Round Body           | 70                 | 12        | UND           | PGLA FAST<br>UND 4-0 X 45-22MM HC CT                                 | PGF409926HC | 4-0  |  | 22 mm 1/2 Circle Conventional Cutting | 45                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 45-60MM ST CT                                 | PGF309943SL  | 3-0  |  | 60 mm Straight Conventional Cutting   | 45                 | 12        | UND           | PGLA FAST<br>UND 5-0 X 75-17MM HC RB                                 | PGF50 9969  | 5-0  |  | 17 mm 1/2 Circle Round Body           | 75                 | 12        | UND           |
| PGLA FAST<br>UND 3-0 X 75-26MM CU RC                                 | PGF309932    | 3-0  |  | 26 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | UND           | PGLA FAST<br>UND 5-0 X 45-13MM CU CT                                 | PGF504931   | 5-0  |  | 13 mm 3/8 Circle Conventional Cutting | 45                 | 12        | UND           |
|                                                                      |              |      |                                                                                     |                                       |                    |           |               | PGLA FAST<br>UND 5-0 X 45-11MM CU CT                                 | PGF504901   | 5-0  |  | 11 mm 3/8 Circle Conventional Cutting | 45                 | 12        | UND           |

| Description                                                          | FG Code    | Size |                                                                                   | Needle Profile<br>Shape, Type &<br>Length | Suture<br>Length<br>(cm) | Pack<br>Size | Suture<br>Colour |
|----------------------------------------------------------------------|------------|------|-----------------------------------------------------------------------------------|-------------------------------------------|--------------------------|--------------|------------------|
| MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST) |            |      |                                                                                   |                                           |                          |              |                  |
| PGLA FAST<br>UND 5-0 X 45-11MM CU RC                                 | PGF50 9915 | 5-0  |  | 11 mm 3/8 Circle<br>Reverse Cutting       | 45                       | 12           | UND              |
| PGLA FAST<br>UND 5-0 X 75-16MM CU RC                                 | PGF502296  | 5-0  |  | 16 mm 3/8 Circle<br>Reverse Cutting       | 75                       | 12           | UND              |
| PGLA FAST<br>UND 5-0 X 70-12MM CU RC                                 | PGF502493  | 5-0  |  | 12 mm 3/8 Circle<br>Reverse Cutting       | 70                       | 12           | UND              |
| PGLA FAST<br>UND 5-0 X 45-16MM CU CT                                 | PGF508441  | 5-0  |  | 16 mm 3/8 Circle<br>Conventional Cutting  | 45                       | 12           | UND              |
| PGLA FAST<br>UND 6-0 X 45-11MM CU RC                                 | PGF60 9913 | 6-0  |  | 11 mm 3/8 Circle<br>Reverse Cutting       | 45                       | 12           | UND              |
| PGLA FAST<br>UND 2-0 X 150 (50X3)                                    | PGF200516  | 2-0  |                                                                                   | Non needled                               | 50                       | 12           | UND              |

# MEGASORB™

## Polyglycolic Acid Suture

### Properties

- 70%-80% of the original is retained until initial 14 Days
- 50% of the original is retained until initial 21 Days
- There is a subsequent loss between four to five weeks post implantation
- Complete absorption usually takes place between 60 and 90 Days

### Indications

- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Ophthalmic procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures


### Advantages

- High tensile strength
- Excellent knotting ability
- Smooth passage through tissue
- Easy handling



### Product Overview

Structure : Multifilament, Braided  
 Colour : Dyed, Undyed



Chemical composition : Polyglycolic acid  
 Coating : Polycaprolactone and Calcium Stearate  
 Sizes : Violet - USP 8/0 (metric 0.4) - USP 2 (metric 5)  
           : Undyed - USP 6/0 (metric 0.7) - USP 2 (metric 5)

Type of absorption : Hydrolysis  
 Sterilisation : Ethylene oxide

Meril Endo-Surgery Product Catalogue | **Wound Closure** 61



Meril Endo-Surgery Product Catalogue | **Wound Closure** 63

Meril Endo-Surgery Product Catalogue | **Wound Closure** 65

| Description                                                      |            |      |                                                                                     |                                                         |  |                    |           | Description                                                      |                                |          |      |                                                                                       |                                     |  |                    |           |               |
|------------------------------------------------------------------|------------|------|-------------------------------------------------------------------------------------|---------------------------------------------------------|--|--------------------|-----------|------------------------------------------------------------------|--------------------------------|----------|------|---------------------------------------------------------------------------------------|-------------------------------------|--|--------------------|-----------|---------------|
| FG Code                                                          |            | Size |                                                                                     | Needle Profile Shape, Type & Length                     |  | Suture Length (cm) | Pack Size | Suture Colour                                                    | FG Code                        |          | Size |                                                                                       | Needle Profile Shape, Type & Length |  | Suture Length (cm) | Pack Size | Suture Colour |
| MEGASORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID) |            |      |                                                                                     |                                                         |  |                    |           | MEGASORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID) |                                |          |      |                                                                                       |                                     |  |                    |           |               |
| PGA VIO 4-0 X 75-22MM HC RB                                      | MS409422   | 4-0  |    | 22 mm 1/2 Circle Round Body                             |  | 75                 | 12        | VIO                                                              | PGA VIO 1 x 70 - 65MM HC RB    | MS010650 | 1-0  |    | 65 mm 1/2 Circle Round Body         |  | 70                 | 12        | VIO           |
| PGA VIO 4-0 X 70-16MM HC RB                                      | MS2305     | 4-0  |    | 16 mm 1/2 Circle Round Body                             |  | 70                 | 12        | VIO                                                              | PGA VIO 1 x 90 - 40MM CU RC    | MS2422   | 1    |    | 40 mm 3/8 Circle Reverse Cutting    |  | 90                 | 12        | VIO           |
| PGA VIO 4-0 X 45-12MM CU RC                                      | MS2494     | 4-0  |    | 12 mm 3/8 Circle Reverse Cutting                        |  | 45                 | 12        | VIO                                                              | PGA VIO 0 x 90 - 36MM HC RB    | MS100346 | 0    |    | 36 mm 1/2 Circle Round Body         |  | 90                 | 12        | UND           |
| PGA UND 4-0 X 75-26MM CU CT                                      | MS9918U    | 4-0  |    | 26 mm 3/8 Circle Conventional Cutting                   |  | 75                 | 12        | UND                                                              | PGA VIO 2 x 70 - 40MM CU RC    | MS209241 | 2-0  |    | 40 mm 3/8 Circle Reverse Cutting    |  | 70                 | 12        | VIO           |
| PGA VIO 5-0 X 45-16MM HC RB                                      | MS2303     | 5-0  |    | 16 mm 1/2 Circle Round Body                             |  | 45                 | 12        | VIO                                                              | PGA VIO 2 x 70 - 30MM CU RC    | MS209242 | 2-0  |    | 30 mm 3/8 Circle Reverse Cutting    |  | 70                 | 12        | VIO           |
| PGA VIO 5-0 X 70-17MM HC RB                                      | MS500303   | 5-0  |    | 17 mm 1/2 Circle Round Body                             |  | 70                 | 12        | VIO                                                              | PGA VIO 2 x 70 - 48MM HC RC    | MS209243 | 2-0  |    | 48 mm 1/2 Circle Reverse Cutting    |  | 70                 | 12        | VIO           |
| PGA VIO 5-0 X 70-13MM HC RB                                      | MS502433   | 5-0  |    | 13 mm 1/2 Circle Round Body                             |  | 70                 | 12        | VIO                                                              | PGA VIO 2 x 75 - 30MM HC RB EX | MS2317SL | 2-0  |    | 30 mm 1/2 Circle Round Body         |  | 75                 | 12        | VIO           |
| PGA VIO 5-0 X 70-36MM HC RB                                      | MS502536   | 5-0  |    | 36 mm 1/2 Circle Round Body                             |  | 70                 | 12        | VIO                                                              | PGA VIO 2 x 70 - 48MM HC RB    | MS208060 | 2-0  |    | 48 mm 1/2 Circle Round Body         |  | 70                 | 12        | VIO           |
| PGA VIO 5-0 X 45-13MM CU RB                                      | MS502513   | 5-0  |    | 13 mm 3/8 Circle Round Body                             |  | 45                 | 12        | VIO                                                              | PGA VIO 4-0 X 70 - 24MM CU RB  | MS402424 | 4-0  |    | 24 mm 3/8 Circle Round Body         |  | 70                 | 12        | VIO           |
| PGA VIO 5-0 X 75-16MM HC RB                                      | MS2303XL   | 5-0  |   | 16 mm 1/2 Circle Round Body                             |  | 75                 | 12        | VIO                                                              | PGA VIO 4-0 X 70-23MM CU RB    | MS402423 | 4-0  |  | 23 mm 3/8 Circle Round Body         |  | 70                 | 12        | VIO           |
| PGA VIO 5-0 X 75-16MM CU RC                                      | MS500500XL | 5-0  |  | 16 mm 3/8 Circle Reverse Cutting                        |  | 75                 | 12        | VIO                                                              | PGA VIO 1 X 70-65MM HC RB      | MS010650 | 1    |  | 65 mm 1/2 Circle Round Body         |  | 70                 | 12        | VIO           |
| PGA VIO 5-0 X 45-12MM CU RC                                      | MS2493     | 5-0  |  | 12 mm 3/8 Circle Reverse Cutting                        |  | 45                 | 12        | VIO                                                              | PGA VIO 0 X 45 ( 45 X3)        | MS100616 | 0    |                                                                                       | Non needed                          |  | 45                 | 12        | VIO           |
| PGA VIO 5-0 X 45-15MM HC CT                                      | MS500557   | 5-0  |  | 15 mm 1/2 Circle Conventional Cutting                   |  | 45                 | 12        | VIO                                                              | PGA VIO 0 X 180                | MS2614NS | 0    |                                                                                       | Non needed                          |  | 180                | 12        | VIO           |
| PGA VIO 6-0 X 70-13MM HC RB                                      | MS602134   | 6-0  |  | 13 mm 1/2 Circle Round Body                             |  | 70                 | 12        | VIO                                                              | PGA VIO 0 X 270 (45X6)         | MS100062 | 0    |                                                                                       | Non needed                          |  | 45                 | 12        | VIO           |
| PGA VIO 6-0 X 45-13MM CU RB                                      | MS602613   | 6-0  |  | 13 mm 3/8 Circle Round Body                             |  | 45                 | 12        | VIO                                                              | PGA VIO 0 X 150                | MS100065 | 0    |                                                                                       | Non needed                          |  | 150                | 12        | VIO           |
| PGA VIO 6-0 X 75-16MM CU RC                                      | MS602612   | 6-0  |  | 16 mm 3/8 Circle Reverse Cutting                        |  | 75                 | 12        | VIO                                                              | PGA VIO 1 X 45 ( 45X3)         | MS016173 | 1    |                                                                                       | Non needed                          |  | 45                 | 12        | VIO           |
| PGA VIO 6-0 X 45-8MM QC SM DN                                    | MS2670DN   | 6-0  |  | 8 mm 1/4 Circle Spatula Micropoint   Spatula Micropoint |  | 45                 | 12        | VIO                                                              | PGA VIO 1 X 45 (45X6)          | MS014506 | 1    |                                                                                       | Non needed                          |  | 45                 | 12        | VIO           |
| PGA VIO 6-0 X 45-12MM CU RC                                      | MS602645   | 6-0  |  | 12 mm 3/8 Circle Reverse Cutting                        |  | 45                 | 12        | VIO                                                              | PGA VIO 2 X 270 (45X6)         | MS020082 | 2    |                                                                                       | Non needed                          |  | 45                 | 12        | VIO           |
| PGA VIO 7-0 X 30-6MM QC SM DN                                    | MS2349DN   | 7-0  |  | 6 mm 1/4 Circle Spatula Micropoint   Spatula Micropoint |  | 30                 | 12        | VIO                                                              | PGA VIO 2-0 X 45 ( 45X3)       | MS200615 | 2-0  |                                                                                       | Non needed                          |  | 45                 | 12        | VIO           |
| PGA VIO 8-0 X 30-6MM CU SM DN                                    | MS2348DN   | 8-0  |  | 6 mm 3/8 Circle Spatula Micropoint   Spatula Micropoint |  | 30                 | 12        | VIO                                                              | PGA VIO 2-0 X 180              | MS2615NS | 2-0  |                                                                                       | Non needed                          |  | 180                | 12        | VIO           |
| PGA VIO 2 x 90 - 65MM HC RB                                      | MS020880   | 2    |  | 65 mm 1/2 Circle Round Body                             |  | 90                 | 12        | VIO                                                              | PGA VIO 2-0 X 150              | MS300055 | 2-0  |                                                                                       | Non needed                          |  | 150                | 12        | VIO           |
|                                                                  |            |      |                                                                                     |                                                         |  |                    |           |                                                                  | PGA VIO 3-0 X 45 (45X3)        | MS300414 | 3-0  |                                                                                       | Non needed                          |  | 45                 | 12        | VIO           |
|                                                                  |            |      |                                                                                     |                                                         |  |                    |           |                                                                  | PGA VIO 3-0 X 50 (50X10)       | MS30500  | 3-0  |                                                                                       | Non needed                          |  | 50                 | 12        | VIO           |
|                                                                  |            |      |                                                                                     |                                                         |  |                    |           |                                                                  | PGA VIO 3-0 X 270 (45X6)       | MS300042 | 3-0  |                                                                                       | Non needed                          |  | 45                 | 12        | VIO           |
|                                                                  |            |      |                                                                                     |                                                         |  |                    |           |                                                                  | PGA VIO 3-0 X 150              | MS300045 | 3-0  |                                                                                       | Non needed                          |  | 150                | 12        | VIO           |
|                                                                  |            |      |                                                                                     |                                                         |  |                    |           |                                                                  | PGA VIO 4-0 X 135 (45X3)       | Ms400284 | 4-0  |                                                                                       | Non needed                          |  | 45                 | 12        | VIO           |

66

Wound Closure | Meril Endo-Surgery Product Catalogue

Meril Endo-Surgery Product Catalogue | Wound Closure

67



# MEGASORB PLUS™

Polyglycolic Acid Suture with Chlorhexidine

Properties

- Retention of 70% knot tensile strength after 14 days
- Complete mass absorption within 60 to 90 days

Indications

- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures


Advantages

- Smooth knot run-down
- Gentle knot repositioning
- High knot holding capacity
- Minimized inflammatory reaction



Also available **MEGASORB T+** Polyglycolic Acid with Triclosan

Product Overview

|                      |                                                                                       |
|----------------------|---------------------------------------------------------------------------------------|
| Structure            | : Multifilament, Braided                                                              |
| Colour               | : Dyed, Undyed                                                                        |
|                      |  |
| Chemical composition | : Polyglycolic acid                                                                   |
| Coating              | : Polycaprolactone + calcium stearate + Antibacterial agent Chlorhexidine             |
| Sizes                | : USP 8/0 (metric 0.4) - USP 2 (metric 5)                                             |
| Type of absorption   | : Hydrolysis                                                                          |
| Sterilization        | : Ethylene oxide                                                                      |

| Description                                                 | FG Code | Size |                                                                                   | Needle Profile<br>Shape, Type &<br>Length | Suture<br>Length<br>(cm) | Pack<br>Size | Suture<br>Colour |
|-------------------------------------------------------------|---------|------|-----------------------------------------------------------------------------------|-------------------------------------------|--------------------------|--------------|------------------|
| MEGASORB PLUS (POLYGLYCOLIC ACID SUTURE WITH CHLORHEXIDINE) |         |      |                                                                                   |                                           |                          |              |                  |
| PGA C PLUS<br>VIO 0 X 90-40MM HC RB                         | MSP2346 | 0    |  | 40 mm 1/2 Circle<br>Round Body            | 90                       | 12           | VIO              |
| PGA C PLUS<br>VIO 1 X 90-40MM HC RB TH                      | MSP2360 | 1    |  | 40 mm 1/2 Circle<br>Round Body            | 90                       | 12           | VIO              |
| PGA C PLUS<br>VIO 1 X 90-40MM HC RB TN                      | MSP2347 | 1    |  | 40 mm 1/2 Circle<br>Round Body            | 90                       | 12           | VIO              |
| PGA C PLUS<br>VIO 1 X 90-40MM HC RC                         | MSP2421 | 1    |  | 40 mm 1/2 Circle<br>Reverse Cutting       | 90                       | 12           | VIO              |
| PGA C PLUS<br>VIO 2-0 X 90-30MM HC RB                       | MSP2317 | 2-0  |  | 30 mm 1/2 Circle<br>Round Body            | 90                       | 12           | VIO              |
| PGA C PLUS<br>VIO 3-0 X 70-20MM HC RB                       | MSP2437 | 3-0  |  | 20 mm 1/2 Circle<br>Round Body            | 70                       | 12           | VIO              |

# MEGASORB RAPID™

Polyglycolic Acid Fast Suture

### Properties

- 50% of the original is retained until initial 6 to 8 Days
- There is a subsequent loss by 17 days post implantation
- Complete absorption usually takes place by 42 Days

### Indications

- Gynaecology/Obstetric procedures (e.g.: episiotomy repair)
- Urology procedures
- Ophthalmic procedures (e.g.: conjunctiva suturing)
- Odontology procedures (e.g.: oral mucosa)
- Paediatric procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

### Advantages

- High initial tensile strength
- Predictable and constant degradation rate
- Good knot security
- Excellent handling properties
- Quick mass absorption



### Product Overview

|                      |                                           |
|----------------------|-------------------------------------------|
| Structure            | : Multifilament, Braided                  |
| Colour               | : Undyed                                  |
| Chemical composition | : Polyglycolic acid                       |
| Sizes                | : USP 6/0 (metric 0.7) - USP 2 (metric 5) |
| Type of absorption   | : Hydrolysis                              |
| Sterilization        | : Gamma radiation                         |





# FILAPRON™

## Poliglecaprone 25 Suture

### Properties

- 80%-60% of the original is retained until initial 7 Days ( DYED )
- 60%-50% of the original is retained until initial 7 Days ( UNDYED )
- 40%-30% of the original is retained until initial 14 Days ( DYED )
- 30%-20% of the original is retained until initial 7 Days ( UNDYED )
- There is a subsequent loss 28 days post implantation
- Complete absorption usually takes place between 90 to 120 Days

### Indications

When short wound support is desirable


- Skinclosure (intracutaneous, subcutaneous)
- Ligatures

### Advantages

- Short-term absorption profile
- High knot pull tensile strength
- Superior handling features



### Product Overview

|                      |                                                                                       |
|----------------------|---------------------------------------------------------------------------------------|
| Structure            | : Monofilament                                                                        |
| Colour               | : Dyed, Undyed                                                                        |
|                      |  |
| Chemical composition | : Poly(glycolide-co-caprolactone)PGCL                                                 |
| Sizes                | : USP 6/0 (metric 0.7) - USP 1 (metric 4)                                             |
| Type of absorption   | : Hydrolysis                                                                          |
| Sterilization        | : Ethylene oxide                                                                      |

Meril Endo-Surgery Product Catalogue | **Wound Closure** 73

Meril Endo-Surgery Product Catalogue | **Wound Closure** 75



Meril Endo-Surgery Product Catalogue | **Wound Closure** 77

# FILAXYN™

## Polydioxanone Suture

### Properties

- 70% of the original is retained until initial 14 Days
- 50%-60% of the original is retained until initial 28 Days
- There is a subsequent minimal absorption until about 90 th post implantation day
- Complete absorption usually takes place between 180 - 210 days

### Indications

- When prolonged wound support is desirable
- Abdominal wall closure
- Orthopedic procedures
- Pediatric sternum closure
- Oncocological procedures

### Advantages

- High initial tensile strength
- Predictable absorption
- Prolonged wound support
- Smooth passage through tissues
- Minimal tissue reaction
- Good handling characteristics
- Standard monofilament knot for security



### Product Overview

|                      |                                           |
|----------------------|-------------------------------------------|
| Structure            | : Monofilament                            |
| Colour               | : Violet                                  |
| Chemical Composition | : Poly (p-dioxanone)                      |
| Sizes                | : USP 6/0 (metric 1.5) - USP 2 (metric 5) |
| Type of absorption   | : Hydrolysis                              |
| Sterilization        | : Ethylene oxide                          |

Meril Endo-Surgery Product Catalogue | **Wound Closure** 81







| Description                          |  |              |     |                                                                                     |                                                    |                    |    | Description |     |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
|--------------------------------------|--|--------------|-----|-------------------------------------------------------------------------------------|----------------------------------------------------|--------------------|----|-------------|-----|---------------|--|---------|--|------|--|-------------------------------------|--|--------------------|--|-----------|--|---------------|--|
| FG Code                              |  | Size         |     | Needle Profile Shape, Type & Length                                                 |                                                    | Suture Length (cm) |    | Pack Size   |     | Suture Colour |  | FG Code |  | Size |  | Needle Profile Shape, Type & Length |  | Suture Length (cm) |  | Pack Size |  | Suture Colour |  |
| FILAXYN (POLYDIOXANONE SUTURE)       |  |              |     |                                                                                     |                                                    |                    |    |             |     |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 5-0 X 70-13MM<br>CU RB       |  | PDX50 1013   | 5-0 |    | 13 mm 3/8 Circle Round Body                        |                    | 70 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 5-0 X 45-13MM<br>CU CT       |  | PDX509863    | 5-0 |    | 13 mm 3/8 Circle Conventional Cutting              |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 5-0 X 70-17MM<br>HC RB DN    |  | PDX509108DN  | 5-0 |    | 17 mm 1/2 Circle Round Body   Round Body           |                    | 70 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 5-0 X 90-17MM<br>HC RB DN    |  | PDX500303DN  | 5-0 |    | 17 mm 1/2 Circle Round Body   Round Body           |                    | 90 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 5-0 X 70-16MM<br>CU CT       |  | PDX509733    | 5-0 |    | 16 mm 3/8 Circle Conventional Cutting              |                    | 70 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 5-0 X 45-19MM<br>CU CT       |  | PDX509869    | 5-0 |    | 19 mm 3/8 Circle Conventional Cutting              |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 5-0 X 70-16MM<br>C RB        |  | PDX5013489   | 5-0 |    | 16 mm 1/2 Circle Round Body                        |                    | 70 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 90-13MM<br>CU RB       |  | PDX609096    | 6-0 |    | 13 mm 3/8 Circle Round Body                        |                    | 90 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 90-8MM<br>CU RB        |  | PDX609091    | 6-0 |    | 8 mm 3/8 Circle Round Body                         |                    | 90 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 45-13MM<br>HC RB       |  | PDX600432    | 6-0 |   | 13 mm 1/2 Circle Round Body                        |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 70-16MM<br>HC RC DN    |  | PDX604323    | 6-0 |  | 16 mm 1/2 Circle Reverse Cutting   Reverse Cutting |                    | 70 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 90-9MM<br>CU TC        |  | PDX601712    | 6-0 |  | 9 mm 3/8 Circle Taper Cut                          |                    | 90 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 45-8MM<br>CU RB DN     |  | PDX60 9091DN | 6-0 |  | 8 mm 3/8 Circle Round Body   Round Body            |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 45-11MM<br>CU RB DN    |  | PDX60 9093DN | 6-0 |  | 11 mm 3/8 Circle Round Body   Round Body           |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 75-13MM<br>CU RB DN    |  | PDX60 127DN  | 6-0 |  | 13 mm 3/8 Circle Round Body   Round Body           |                    | 75 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 45-13MM<br>CU CT       |  | PDX609861    | 6-0 |  | 13 mm 3/8 Circle Conventional Cutting              |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 45-13MM<br>HC SM       |  | PDX607613    | 6-0 |  | 13 mm 1/2 Circle Spatula Micropoint                |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 6-0 X 45-11MM<br>CU RC       |  | PDX600487    | 6-0 |  | 11 mm 3/8 Circle Reverse Cutting                   |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 7-0 X 60-9MM<br>CU RB DN     |  | PDX70 155DN  | 7-0 |  | 9 mm 3/8 Circle Round Body   Round Body            |                    | 60 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 7-0 X 45-13MM<br>CU RB DN    |  | PDX709095DN  | 7-0 |  | 13 mm 3/8 Circle Round Body   Round Body           |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| FILAXYN (POLYDIOXANONE SUTURE)       |  |              |     |                                                                                     |                                                    |                    |    |             |     |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 7-0 X 70-9.3MM<br>CU TC DN   |  | PDX70 1711DN | 7-0 |  | 9.3 mm 3/8 Circle Taper Cut   Taper Cut            |                    | 70 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 3-0 X 70-26MM<br>HC RB DN-3D |  | PDX301316Z   | 3-0 |  | 26 mm 1/2 Circle Round Body   Round Body           |                    | 70 | 36          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 4-0 X 45-19MM<br>CU RC       |  | PDX400422SL  | 4-0 |  | 19 mm 3/8 Circle Reverse Cutting                   |                    | 45 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 3-0 X 70-24MM<br>CU RC       |  | PDX304520    | 3-0 |  | 24mm 3/8 Circle Reverse Cutting                    |                    | 70 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |
| PDO VIO 3-0 X 90-26MM<br>HC RB       |  | PDX300316XL  | 3-0 |  | 26mm 1/2 Circle Round Body                         |                    | 90 | 12          | VIO |               |  |         |  |      |  |                                     |  |                    |  |           |  |               |  |

86

Wound Closure

| Meril Endo-Surgery Product Catalogue

Meril Endo-Surgery Product Catalogue

| Wound Closure

87

# MERISOFT PLAIN<sup>TM</sup>

## Plain Catgut Suture

### Properties

The tensile strength is maintained for not more than 7 to 10 days post implantation, and its absorption is complete within 60 to 90 days.

### Indications

MERISOFT PLAIN<sup>TM</sup> suture is indicated for use in general soft tissue approximation and/or ligation.

### Advantages

- Highest level of collagen guarantees predictable absorption
- Accurate polishing ensures smooth tissue passage
- Packed in liquid to reduce memory & increase pliability
- Rigidly controlled spinning & drying to ensure higher tensile strength



### Product Overview

|                     |                                                                                                         |
|---------------------|---------------------------------------------------------------------------------------------------------|
| Primary packaging   | : Packed in an aluminum foil pouch                                                                      |
| Secondary packaging | : Primary pack is loaded in a secondary pouch made of medical grade paper & laminated polyethylene film |
| Sizes               | : USP 6/0 (metric 1) - USP 2 (metric 6)                                                                 |
| Colour              | : Yellow-tan                                                                                            |
| <div></div>         |                                                                                                         |
| Sterilization       | : Gamma radiation                                                                                       |

Meril Endo-Surgery Product Catalogue | **Wound Closure** 91



# MERISOFT CHROMIC™

## Chromic Catgut Suture

### Properties

MERISOFT CHROMIC™ gut suture is treated with a chromium salt solution to resist the enzymatic action, prolonging the absorption time to 90 to 120 days.

### Indications

MERISOFT CHROMIC™ suture is indicated for use in general soft tissue approximation and/or ligation.

### Advantages

- Highest level of collagen guarantees predictable absorption
- Accurate polishing ensures smooth tissue passage
- Packed in liquid to reduce memory & increase pliability



### Product Overview

|                                                                                       |                                                                                                         |
|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|
| Primary packaging                                                                     | : Packed in an aluminum foil pouch                                                                      |
| Secondary packaging                                                                   | : Primary pack is loaded in a secondary pouch made of medical grade paper & laminated polyethylene film |
| Sizes                                                                                 | : USP 6/0 (metric 1) - USP 2 (metric 6)                                                                 |
| Colour                                                                                | : Brown dyed                                                                                            |
|  |                                                                                                         |
| Sterilization                                                                         | : Gamma radiation                                                                                       |

| Description                              | FG Code  | Size |                                                                                     | Needle Profile Shape, Type & Length | Suture Length (cm) | Pack Size | Suture Colour | Description                                 | FG Code  | Size |                                                                                       | Needle Profile Shape, Type & Length           | Suture Length (cm) | Pack Size | Suture Colour |
|------------------------------------------|----------|------|-------------------------------------------------------------------------------------|-------------------------------------|--------------------|-----------|---------------|---------------------------------------------|----------|------|---------------------------------------------------------------------------------------|-----------------------------------------------|--------------------|-----------|---------------|
| MERISOFT CHROMIC (CHROMIC CATGUT SUTURE) |          |      |                                                                                     |                                     |                    |           |               | MERISOFT CHROMIC (CHROMIC CATGUT SUTURE)    |          |      |                                                                                       |                                               |                    |           |               |
| CHROMIC CATGUT<br>0 X 76-30MM CU RB      | ME4217   | 0    |    | 30 mm 3/8 Circle Round Body         | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>1 X 75-30MM HC RB         | ME014995 | 1    |    | 30 mm 1/2 Circle Round Body                   | 75                 | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-40MM CU RB      | ME4221   | 0    |    | 40 mm 3/8 Circle Round Body         | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>1 X 100-40MM HC TRO       | ME4239   | 1    |    | 40 mm 1/2 Circle Trocar                       | 100                | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 100-45MM HC RB     | ME4229   | 0    |    | 45 mm 1/2 Circle Round Body         | 100                | 12        | CHROMIC       | CHROMIC CATGUT<br>1 X 76-63MM CU RB         | Me4267   | 1    |    | 63 mm 3/8 Circle Round Body                   | 76                 | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-30MM HC RB      | ME4242   | 0    |    | 30 mm 1/2 Circle Round Body         | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>1 X 76-50MM HC RC         | ME4273   | 1    |    | 50 mm 1/2 Circle Reverse Cutting              | 76                 | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-40MM HC RB TN   | ME4246   | 0    |    | 40 mm 1/2 Circle Round Body         | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>1 X 100-40MM HC TC        | ME4297   | 1    |    | 40 mm 1/2 Circle Taper Cut                    | 100                | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-40MM HC RB TH   | ME4258   | 0    |    | 40 mm 1/2 Circle Round Body         | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>2 X 100-45MM HC RB        | ME4228   | 2    |    | 45 mm 1/2 Circle Round Body                   | 100                | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-36MM HC RB      | ME104036 | 0    |    | 36 mm 1/2 Circle Round Body         | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>2 X 100-50MM HC RB        | ME024250 | 2    |    | 50 mm 1/2 Circle Round Body                   | 100                | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 100-40MM HC TC     | ME4279   | 0    |    | 40 mm 1/2 Circle Taper Cut          | 100                | 12        | CHROMIC       | CHROMIC CATGUT<br>2 X 100-40MM HC TRO       | ME4244   | 2    |   | 40 mm 1/2 Circle Trocar                       | 100                | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-45MM ST RB      | ME4205   | 0    |  | 45 mm Straight Round Body           | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>2 X 100-65MM HC RB        | ME4278   | 2    |  | 65 mm 1/2 Circle Round Body                   | 100                | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-45MM CU RB      | ME4225   | 0    |  | 45 mm 3/8 Circle Round Body         | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>2 X 75-40MM HC RC         | ME024244 | 2    |  | 40 mm 1/2 Circle Reverse Cutting              | 75                 | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-45MM HC RB      | ME4250   | 0    |  | 45 mm 1/2 Circle Round Body         | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>2-0 X 76-30MM CU RB       | ME4216   | 2-0  |  | 30 mm 3/8 Circle Round Body                   | 76                 | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-40MM HC TC      | ME4257   | 0    |  | 40 mm 1/2 Circle Taper Cut          | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>2-0 X 76-30MM HC RB       | ME4241   | 2-0  |  | 30 mm 1/2 Circle Round Body                   | 76                 | 12        | CHROMIC       |
| CHROMIC CATGUT<br>0 X 76-40MM HC RC      | ME104040 | 0    |  | 40 mm 1/2 Circle Reverse Cutting    | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>2-0 X 152-36MM HC RCRB DN | ME4277   | 2-0  |  | 36 mm 1/2 Circle Reverse Cutting   Round Body | 152                | 12        | CHROMIC       |
| CHROMIC CATGUT<br>1 X 76-45MM CU RB      | ME4226   | 1    |  | 45 mm 3/8 Circle Round Body         | 76                 | 12        | CHROMIC       | CHROMIC CATGUT<br>2-0 X 75-40MM HC RB       | ME200991 | 2-0  |  | 40 mm 1/2 Circle Round Body                   | 75                 | 12        | CHROMIC       |
| CHROMIC CATGUT<br>1 X 100-45MM HC RB     | ME4227   | 1    |  | 45 mm 1/2 Circle Round Body         | 100                | 12        | CHROMIC       | CHROMIC CATGUT<br>2-0 X 76-26MM HC RB       | ME204225 | 2-0  |  | 26 mm 1/2 Circle Round Body                   | 76                 | 12        | CHROMIC       |
| CHROMIC CATGUT<br>1 X 100-50MM HC RB     | ME4248   | 1    |  | 50 mm 1/2 Circle Round Body         | 100                | 12        | CHROMIC       |                                             |          |      |                                                                                       |                                               |                    |           |               |
| CHROMIC CATGUT<br>1 X 76-40MM HC RB      | ME4259   | 1    |  | 40 mm 1/2 Circle Round Body         | 76                 | 12        | CHROMIC       |                                             |          |      |                                                                                       |                                               |                    |           |               |

Meril Endo-Surgery Product Catalogue | **Wound Closure** 97

# FILASILK™

## Silk Suture

### Indications

- General surgery procedures
- Skin closure
- Oral surgery procedures
- Ophthalmic procedures
- Ligatures

### Advantages

- Excellent handling properties
- Good knot security



### Product Overview

|               |                                           |
|---------------|-------------------------------------------|
| Structure     | : Braided                                 |
| Colour        | : Black                                   |
| Origin        | : Silkworm larval cocoons (Bombyx mori)   |
| Composition   | : Natural silk filaments                  |
| Sizes         | : USP 8/0 (metric 0.4) - USP 2 (metric 5) |
| Sterilization | : Gamma radiation                         |



| Description                 | FG Code    | Size |                                                                                     | Needle Profile Shape, Type & Length   | Suture Length (cm) | Pack Size | Suture Colour | Description                 | FG Code     | Size |                                                                                       | Needle Profile Shape, Type & Length   | Suture Length (cm) | Pack Size | Suture Colour |
|-----------------------------|------------|------|-------------------------------------------------------------------------------------|---------------------------------------|--------------------|-----------|---------------|-----------------------------|-------------|------|---------------------------------------------------------------------------------------|---------------------------------------|--------------------|-----------|---------------|
| FILASILK SILK SUTURE        |            |      |                                                                                     |                                       |                    |           |               | FILASILK SILK SUTURE        |             |      |                                                                                       |                                       |                    |           |               |
| SILK BLK 0 X 90-60MM ST CT  | SLK100789  | 0    |    | 60 mm Straight Conventional Cutting   | 90                 | 12        | BLACK         | SILK BLK 0 X 76-30MM HC RB  | SLK10 5334  | 0    |    | 30 mm 1/2 Circle Round Body           | 76                 | 12        | BLACK         |
| SILK BLK 0 X 90-50MM HC RC  | SLK100552  | 0    |    | 50 mm 1/2 Circle Reverse Cutting      | 90                 | 12        | BLACK         | SILK BLK 0 X 90-40MM HC RB  | SLK105041   | 0    |    | 40 mm 1/2 Circle Round Body           | 90                 | 12        | BLACK         |
| SILK BLK 0 X 90-35MM CU RC  | SLK100371  | 0    |    | 35 mm 3/8 Circle Reverse Cutting      | 90                 | 12        | BLACK         | SILK BLK 0 X 90-40MM HC RC  | SLK102224   | 0    |    | 40 mm 1/2 Circle Reverse Cutting      | 90                 | 12        | BLACK         |
| SILK BLK 0 X 76-45MM CU RC  | SLK10 5037 | 0    |    | 45 mm 3/8 Circle Reverse Cutting      | 76                 | 12        | BLACK         | SILK BLK 0 X 90-75MM ST CT  | SLK100791   | 0    |    | 75 mm Straight Conventional Cutting   | 90                 | 12        | BLACK         |
| SILK BLK 0 X 90-30MM HC RB  | SLK10 5332 | 0    |    | 30 mm 1/2 Circle Round Body           | 90                 | 12        | BLACK         | SILK BLK 0 X 70-39MM CU RC  | SLK101039   | 0    |    | 39 mm 3/8 Circle Reverse Cutting      | 70                 | 12        | BLACK         |
| SILK BLK 0 X 75-30MM HC RC  | SLK102223  | 0    |    | 30 mm 1/2 Circle Reverse Cutting      | 75                 | 12        | BLACK         | SILK BLK 1 X 76-60MM CU CT  | SLK01 5062  | 1    |    | 60 mm 3/8 Circle Conventional Cutting | 76                 | 12        | BLACK         |
| SILK BLK 0 X 45-23MM HC RC  | SLK100373  | 0    |    | 23 mm 1/2 Circle Reverse Cutting      | 45                 | 12        | BLACK         | SILK BLK 1 X 75-30MM CU CT  | SLK012243   | 1    |    | 30 mm 3/8 Circle Conventional Cutting | 75                 | 12        | BLACK         |
| SILK BLK 0 X 75-39MM CU CT  | SLK105039  | 0    |    | 39 mm 3/8 Circle Conventional Cutting | 75                 | 12        | BLACK         | SILK BLK 1 X 75-30MM HC RB  | SLK010723   | 1    |    | 30 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 0 X 75-30MM HC RB  | SLK100334  | 0    |   | 30 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         | SILK BLK 1 X 90-40MM HC RB  | SLK010339   | 1    |   | 40 mm 1/2 Circle Round Body           | 90                 | 12        | BLACK         |
| SILK BLK 0 X 75-24MM CU RC  | SLK105024  | 0    |  | 24 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         | SILK BLK 1 X 75-27MM HC RB  | SLK0115175  | 1    |  | 27 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 0 X 75-37MM HC RB  | SLK100337  | 0    |  | 37 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         | SILK BLK 1 X 100-90MM CU CT | SLK012793   | 1    |  | 90 mm 3/8 Circle Conventional Cutting | 100                | 12        | BLACK         |
| SILK BLK 0 X 75-26MM HC RB  | SLK100383  | 0    |  | 26 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         | SILK BLK 1 X 90-24MM CU RB  | SLK015124   | 1    |  | 24 mm 3/8 Circle Round Body           | 90                 | 12        | BLACK         |
| SILK BLK 0 X 75-45MM CU RC  | SLK105076  | 0    |  | 45 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         | SILK BLK 1 X 90-26MM CU CT  | SLK015126   | 1    |  | 26 mm 3/8 Circle Conventional Cutting | 90                 | 12        | BLACK         |
| SILK BLK 0 X 90-30MM CU RC  | SLK102508  | 0    |  | 30 mm 3/8 Circle Reverse Cutting      | 90                 | 12        | BLACK         | SILK BLK 1 X 75-30MM HC RC  | SLK015130   | 1    |  | 30 mm 1/2 Circle Reverse Cutting      | 75                 | 12        | BLACK         |
| SILK BLK 0 X 90-30MM CU RB  | SLK105030  | 0    |  | 30 mm 3/8 Circle Round Body           | 90                 | 12        | BLACK         | SILK BLK 1 X 75-38MM CU RC  | SLK018138   | 1    |  | 38 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         |
| SILK BLK 0 X 76-30MM CU RB  | SLK10 5017 | 0    |  | 30 mm 3/8 Circle Round Body           | 76                 | 12        | BLACK         | SILK BLK 1 X 75-37MM HC RB  | SLK01 5137N | 1    |  | 37 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 0 X 75-60MM ST CT  | SLK1013113 | 0    |  | 60 mm Straight Conventional Cutting   | 75                 | 12        | BLACK         | SILK BLK 1 X 90-40MM HC RC  | SLK01 5005  | 1    |  | 40 mm 1/2 Circle Reverse Cutting      | 90                 | 12        | BLACK         |
| SILK BLK 0 X 100-90MM CU CT | SLK1013114 | 0    |  | 90 mm 3/8 Circle Conventional Cutting | 100                | 12        | BLACK         | SILK BLK 1 X 75-40MM CU CT  | SLK0113116  | 1    |  | 40 mm 3/8 Circle Conventional Cutting | 75                 | 12        | BLACK         |
| SILK BLK 0 X 75-31MM HC RB  | SLK103341  | 0    |  | 31 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         | SILK BLK 1 X 100-60MM ST CT | SLK010391   | 1    |  | 60 mm Straight Conventional Cutting   | 100                | 12        | BLACK         |
| SILK BLK 0 X 75-40MM CU RC  | SLK105040  | 0    |  | 40 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         |                             |             |      |                                                                                       |                                       |                    |           |               |

| Description                                    | FG Code     | Size |                                                                                     | Needle Profile Shape, Type & Length   | Suture Length (cm) | Pack Size | Suture Colour | Description                     | FG Code     | Size |                                                                                       | Needle Profile Shape, Type & Length   | Suture Length (cm) | Pack Size | Suture Colour |
|------------------------------------------------|-------------|------|-------------------------------------------------------------------------------------|---------------------------------------|--------------------|-----------|---------------|---------------------------------|-------------|------|---------------------------------------------------------------------------------------|---------------------------------------|--------------------|-----------|---------------|
| FILASILK SILK SUTURE                           |             |      |                                                                                     |                                       |                    |           |               | FILASILK SILK SUTURE            |             |      |                                                                                       |                                       |                    |           |               |
| SILK BLK 1 X 75-30MM CU RC                     | SLK010689   | 1    |    | 30 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         | SILK BLK 2-0 X 75-26MM HC RB    | SLK200327   | 2-0  |    | 26 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 2 X 90-90MM CU RC                     | SLK020794   | 2    |    | 90 mm 3/8 Circle Reverse Cutting      | 90                 | 12        | BLACK         | SILK BLK 2-0 X 75-17MM HC RB    | SLK205610   | 2-0  |    | 17 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 2 X 90-40MM HC RB                     | SLK020936   | 2    |    | 40 mm 1/2 Circle Round Body           | 90                 | 12        | BLACK         | SILK BLK 2-0 X 75-25MM HC CT    | SLK205225   | 2-0  |    | 25 mm 1/2 Circle Conventional Cutting | 75                 | 12        | BLACK         |
| SILK BLK 2 X 75-30MM HC RC                     | SLK025230   | 2    |    | 30 mm 1/2 Circle Reverse Cutting      | 75                 | 12        | BLACK         | SILK BLK 2-0 X 76-50MM ST RB    | SLK20 5208  | 2-0  |    | 50 mm Straight Round Body             | 76                 | 12        | BLACK         |
| SILK BLK 2 X 75-40MM HC RC                     | SLK025340   | 2    |    | 40 mm 1/2 Circle Reverse Cutting      | 75                 | 12        | BLACK         | SILK BLK 2-0 X 90-15MM CU RC    | SLK200503   | 2-0  |    | 15 mm 3/8 Circle Reverse Cutting      | 90                 | 12        | BLACK         |
| SILK BLK 2 X 75-40MM HC RB                     | SLK020936SL | 2    |    | 40 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         | SILK BLK 2-0 X 75-22MM HC CT    | SLK2052221  | 2-0  |    | 22 mm 1/2 Circle Conventional Cutting | 75                 | 12        | BLACK         |
| SILK BLK 2 X 100-90MM CU CT                    | SLK020792   | 2    |    | 90 mm 3/8 Circle Conventional Cutting | 100                | 12        | BLACK         | SILK BLK 2-0 X 76-60MM ST CT    | SLK20 5290  | 2-0  |    | 60 mm Straight Conventional Cutting   | 76                 | 12        | BLACK         |
| SILK BLK 2 X 70-36MM CU RC                     | SLK021236   | 2    |    | 36 mm 3/8 Circle Reverse Cutting      | 70                 | 12        | BLACK         | SILK BLK 2-0 X 100-90MM CU CT   | SLK200791   | 2-0  |    | 90 mm 3/8 Circle Conventional Cutting | 100                | 12        | BLACK         |
| SILK BLK 2 X 70-39MM CU RC                     | SLK021239   | 2    |   | 39 mm 3/8 Circle Reverse Cutting      | 70                 | 12        | BLACK         | SILK BLK 2-0 X 45-26MM CU RC    | SLK200685   | 2-0  |    | 26 mm 3/8 Circle Reverse Cutting      | 45                 | 12        | BLACK         |
| SILK BLK 2 X 75-35MM HC RB                     | SLK028235   | 2    |  | 35 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         | SILK BLK 2-0 X 75-30MM HC CT    | SLK205231   | 2-0  |   | 30 mm 1/2 Circle Conventional Cutting | 75                 | 12        | BLACK         |
| SILK BLK 1 X 50-90MM CU CT (2 suture per foil) | SLK012797   | 1    |  | 90 mm 3/8 Circle Conventional Cutting | 50                 | 12        | BLACK         | SILK BLK 2-0 X 90-26MM HC RB    | SLK200327XL | 2-0  |  | 26 mm 1/2 Circle Round Body           | 90                 | 12        | BLACK         |
| SILK BLK 2-0 X 90-20MM CU RC                   | SLK200367   | 2-0  |  | 20 mm 3/8 Circle Reverse Cutting      | 90                 | 12        | BLACK         | SILK BLK 2-0 X 75-26MM CU RC    | SLK200685XL | 2-0  |  | 26 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 90-60MM ST CT                   | SLK200790   | 2-0  |  | 60 mm Straight Conventional Cutting   | 90                 | 12        | BLACK         | SILK BLK 2-0 X 90-24MM CU RB    | SLK205224   | 2-0  |  | 24 mm 3/8 Circle Round Body           | 90                 | 12        | BLACK         |
| SILK BLK 2-0 X 76-45MM CU RC                   | SLK20 5036  | 2-0  |  | 45 mm 3/8 Circle Reverse Cutting      | 76                 | 12        | BLACK         | SILK BLK 2-0 X 90-26MM CU CT    | SLK205226   | 2-0  |  | 26 mm 3/8 Circle Conventional Cutting | 90                 | 12        | BLACK         |
| SILK BLK 2-0 X 76-50MM HC RC                   | SLK20 5052  | 2-0  |  | 50 mm 1/2 Circle Reverse Cutting      | 76                 | 12        | BLACK         | SILK BLK 2-0 X 45-30MM HC RC    | SLK200530   | 2-0  |  | 30 mm 1/2 Circle Reverse Cutting      | 45                 | 12        | BLACK         |
| SILK BLK 2-0 X 90-30MM HC RB                   | SLK20 5331  | 2-0  |  | 30 mm 1/2 Circle Round Body           | 90                 | 12        | BLACK         | SILK BLK 2-0 X 75-30MM CU RC    | SLK205230   | 2-0  |  | 30 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 76-25MM HC TC                   | SLK20 5670  | 2-0  |  | 25 mm 1/2 Circle Taper Cut            | 76                 | 12        | BLACK         | SILK BLK 2-0 X 75-16MM CU CT    | SLK200509   | 2-0  |  | 16 mm 3/8 Circle Conventional Cutting | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 45-23MM HC RC                   | SLK200737   | 2-0  |  | 23 mm 1/2 Circle Reverse Cutting      | 45                 | 12        | BLACK         | SILK BLK 2-0 X 76-25MM HC RB    | SLK205671   | 2-0  |  | 25 mm 1/2 Circle Round Body           | 76                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-22MM HC RB                   | SLK205222   | 2-0  |  | 22 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         | SILK BLK 2-0 X 100-55MM ST CT   | SLK200775   | 2-0  |  | 55 mm Straight Conventional Cutting   | 100                | 12        | BLACK         |
|                                                |             |      |                                                                                     |                                       |                    |           |               | SILK BLK 2-0 X 76-36MM CU RC    | SLK205236   | 2-0  |  | 36 mm 3/8 Circle Reverse Cutting      | 76                 | 12        | BLACK         |
|                                                |             |      |                                                                                     |                                       |                    |           |               | SILK BLK 2-0 X 76-26MM HC RB VB | SLK20 5065  | 2-0  |  | 26 mm 1/2 Circle Round Body           | 76                 | 12        | BLACK         |

| Description                     | FG Code     | Size |                                                                                     | Needle Profile Shape, Type & Length   | Suture Length (cm) | Pack Size | Suture Colour | Description                  | FG Code    | Size |                                                                                       | Needle Profile Shape, Type & Length   | Suture Length (cm) | Pack Size | Suture Colour |
|---------------------------------|-------------|------|-------------------------------------------------------------------------------------|---------------------------------------|--------------------|-----------|---------------|------------------------------|------------|------|---------------------------------------------------------------------------------------|---------------------------------------|--------------------|-----------|---------------|
| FILASILK SILK SUTURE            |             |      |                                                                                     |                                       |                    |           |               | FILASILK SILK SUTURE         |            |      |                                                                                       |                                       |                    |           |               |
| SILK BLK 2-0 X 90-17MM HC TC    | SLK20 5611  | 2-0  |    | 17 mm 1/2 Circle Taper Cut            | 90                 | 12        | BLACK         | SILK BLK 2-0 X 70-36MM HC RB | SLK201336  | 2-0  |    | 36 mm 1/2 Circle Round Body           | 70                 | 12        | BLACK         |
| SILK BLK 2-0 X 76-30MM HC RB VB | SLK20 533VB | 2-0  |    | 30 mm 1/2 Circle Round Body           | 76                 | 12        | BLACK         | SILK BLK 3-0 X 90-60MM ST CT | SLK300622  | 3-0  |    | 60 mm Straight Conventional Cutting   | 90                 | 12        | BLACK         |
| SILK BLK 2-0 X 100-75MM ST CT   | SLK2013102  | 2-0  |    | 75 mm Straight Conventional Cutting   | 100                | 12        | BLACK         | SILK BLK 3-0 X 90-25MM HC RC | SLK300379  | 3-0  |    | 25 mm 1/2 Circle Reverse Cutting      | 90                 | 12        | BLACK         |
| SILK BLK 2-0 X 90-35MM CU RC    | SLK200667   | 2-0  |    | 35 mm 3/8 Circle Reverse Cutting      | 90                 | 12        | BLACK         | SILK BLK 3-0 X 90-16MM CU CT | SLK30 5003 | 3-0  |    | 16 mm 3/8 Circle Conventional Cutting | 90                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-36MM CU RC    | SLK205275   | 2-0  |    | 36 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         | SILK BLK 3-0 X 76-26MM CU RC | SLK30 5028 | 3-0  |    | 26 mm 3/8 Circle Reverse Cutting      | 76                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-35MM CU CT    | SLK202429   | 2-0  |    | 35 mm 3/8 Circle Conventional Cutting | 75                 | 12        | BLACK         | SILK BLK 3-0 X 90-25MM HC RB | SLK30 5085 | 3-0  |    | 25 mm 1/2 Circle Round Body           | 90                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-80MM CU CT    | SLK200261   | 2-0  |    | 80 mm 3/8 Circle Conventional Cutting | 75                 | 12        | BLACK         | SILK BLK 3-0 X 76-20MM HC RB | SLK30 5087 | 3-0  |    | 20 mm 1/2 Circle Round Body           | 76                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-30MM HC RB    | SLK200333   | 2-0  |    | 30 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         | SILK BLK 3-0 X 75-24MM CU RC | SLK300684  | 3-0  |    | 24 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-48MM HC CT    | SLK200552   | 2-0  |   | 48 mm 1/2 Circle Conventional Cutting | 75                 | 12        | BLACK         | SILK BLK 3-0 X 75-17MM HC RB | SLK3054601 | 3-0  |    | 17 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-60MM ST CT    | SLK200790SL | 2-0  |  | 60 mm Straight Conventional Cutting   | 75                 | 12        | BLACK         | SILK BLK 3-0 X 75-22MM HC RB | SLK300587  | 3-0  |   | 22 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 45-24MM CU RC    | SLK200695   | 2-0  |  | 24 mm 3/8 Circle Reverse Cutting      | 45                 | 12        | BLACK         | SILK BLK 3-0 X 75-30MM CU RC | SLK300675  | 3-0  |  | 30 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 100-40MM CU CT   | SLK203959   | 2-0  |  | 40 mm 3/8 Circle Conventional Cutting | 100                | 12        | BLACK         | SILK BLK 3-0 X 75-26MM HC RB | SLK300832  | 3-0  |  | 26 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 45-30MM HC RC C1 | SLK205227   | 2-0  |  | 30 mm 1/2 Circle Reverse Cutting      | 45                 | 12        | BLACK         | SILK BLK 3-0 X 75-30MM HC RB | SLK300571  | 3-0  |  | 30 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-75MM ST CT    | SLK200766   | 2-0  |  | 75 mm Straight Conventional Cutting   | 75                 | 12        | BLACK         | SILK BLK 3-0 X 45-21MM HC RC | SLK300632  | 3-0  |  | 21 mm 1/2 Circle Reverse Cutting      | 45                 | 12        | BLACK         |
| SILK BLK 2-0 X 76-30MM HC RB    | SLK20 5333  | 2-0  |  | 30 mm 1/2 Circle Round Body           | 76                 | 12        | BLACK         | SILK BLK 3-0 X 75-16MM CU RC | SLK300882  | 3-0  |  | 16 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 90-35MM HC CT    | SLK20 5335  | 2-0  |  | 35 mm 1/2 Circle Conventional Cutting | 90                 | 12        | BLACK         | SILK BLK 3-0 X 75-22MM HC CT | SLK302508  | 3-0  |  | 22 mm 1/2 Circle Conventional Cutting | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-39MM CU RC    | SLK200648   | 2-0  |  | 39 mm 3/8 Circle Reverse Cutting      | 75                 | 12        | BLACK         | SILK BLK 3-0 X 45-16MM CU RC | SLK300883  | 3-0  |  | 16 mm 3/8 Circle Reverse Cutting      | 45                 | 12        | BLACK         |
| SILK BLK 2-0 X 100-60MM ST RC   | SLK200723   | 2-0  |  | 60 mm Straight Reverse Cutting        | 100                | 12        | BLACK         | SILK BLK 3-0 X 75-19MM CU RB | SLK305319  | 3-0  |  | 19 mm 3/8 Circle Round Body           | 75                 | 12        | BLACK         |
| SILK BLK 2-0 X 75-40MM HC RB    | SLK200584   | 2-0  |  | 40 mm 1/2 Circle Round Body           | 75                 | 12        | BLACK         | SILK BLK 3-0 X 90-20MM HC RB | SLK300588  | 3-0  |  | 20 mm 1/2 Circle Round Body           | 90                 | 12        | BLACK         |
|                                 |             |      |                                                                                     |                                       |                    |           |               | SILK BLK 3-0 X 90-19MM CU RC | SLK302231  | 3-0  |  | 19 mm 3/8 Circle Reverse Cutting      | 90                 | 12        | BLACK         |
|                                 |             |      |                                                                                     |                                       |                    |           |               | SILK BLK 3-0 X 75-30MM HC RC | SLK305330  | 3-0  |  | 30 mm 1/2 Circle Reverse Cutting      | 75                 | 12        | BLACK         |







| Description                            | FG Code    | Size | Needle Profile<br>Shape, Type &<br>Length | Suture<br>Length<br>(cm) | Pack<br>Size | Suture<br>Colour |
|----------------------------------------|------------|------|-------------------------------------------|--------------------------|--------------|------------------|
| FILASILK SILK SUTURE                   |            |      |                                           |                          |              |                  |
| SILK BLK 1 X 45 (45X12)                | SLS014512  | 1    | Non needed                                | 45                       | 12           | BLACK            |
| SILK BLK 2 X 76 (76X2)                 | SLS02 216  | 2    | Non needed                                | 76                       | 12           | BLACK            |
| SILK BLK 2 X 180                       | SLS02180   | 2    | Non needed                                | 180                      | 12           | BLACK            |
| SILK BLK 2 X 270 (45X6)                | SLS02364   | 2    | Non needed                                | 45                       | 12           | BLACK            |
| SILK BLK 2 X 250                       | SLS02252   | 2    | Non needed                                | 250                      | 12           | BLACK            |
| SILK BLK 2 X 45 (45X12)                | SLS024512  | 2    | Non needed                                | 45                       | 12           | BLACK            |
| SILK BLK 3 X 76 (76X2)                 | SLS03 217  | 3    | Non needed                                | 76                       | 12           | BLACK            |
| SILK BLK 2-0 X 76 (76X2)               | SLS20 213  | 2-0  | Non needed                                | 76                       | 12           | BLACK            |
| SILK BLK 2-0 X 75 (75X12)              | SLS20305   | 2-0  | Non needed                                | 75                       | 12           | BLACK            |
| SILK BLK 2-0 X 180                     | SLS20180   | 2-0  | Non needed                                | 180                      | 12           | BLACK            |
| SILK BLK 2-0 X 75 (75X5)               | SLS20375   | 2-0  | Non needed                                | 75                       | 12           | BLACK            |
| SILK BLK 2-0 X 35 (35X10)              | SLS20 223  | 2-0  | Non needed                                | 35                       | 12           | BLACK            |
| SILK BLK 2-0 X 45 (45X17)              | SLS20203   | 2-0  | Non needed                                | 45                       | 12           | BLACK            |
| SILK BLK 2-0 X 75 (75X10)              | SLS201075  | 2-0  | Non needed                                | 75                       | 12           | BLACK            |
| SILK BLK 2-0 X 250-2D                  | SLS20253Y  | 2-0  | Non needed                                | 250                      | 24           | BLACK            |
| SILK BLK 2-0 X 150 (150X2)             | SLS20605   | 2-0  | Non needed                                | 150                      | 12           | BLACK            |
| SILK BLK 2-0 X 250                     | SLS20253   | 2-0  | Non needed                                | 250                      | 12           | BLACK            |
| SILK BLK 2-0 X 45 (45X12)              | SLS204512  | 2-0  | Non needed                                | 45                       | 12           | BLACK            |
| SILK BLK 3-0 X 76 (76X2)               | SLS30 212  | 3-0  | Non needed                                | 76                       | 12           | BLACK            |
| SILK BLK 3-0 X 35 (35X10)              | SLS30 222  | 3-0  | Non needed                                | 35                       | 12           | BLACK            |
| SILK BLK 3-0 X 75 (75X10)              | SLS30750   | 3-0  | Non needed                                | 75                       | 12           | BLACK            |
| SILK BLK 3-0 X 180                     | SLS30180   | 3-0  | Non needed                                | 180                      | 12           | BLACK            |
| SILK BLK 3-0 X 45 (45X17)              | SLS30202   | 3-0  | Non needed                                | 45                       | 12           | BLACK            |
| SILK BLK 3-0 X 45 (45X12)              | SLS304512  | 3-0  | Non needed                                | 45                       | 12           | BLACK            |
| SILK BLK 4-0 X 76 (76X2)               | SLS40 211  | 4-0  | Non needed                                | 76                       | 12           | BLACK            |
| SILK BLK 0 X 25MTR<br>(6 EA PER BOX)   | SLK10 824R | 0    | Non needed                                | 25 Mtr                   | 6            | BLACK            |
| SILK BLK 1 X 25MTR<br>(6 EA PER BOX)   | SLK01 825R | 1    | Non needed                                | 25 Mtr                   | 6            | BLACK            |
| SILK BLK 2 X 25MTR<br>(6 EA PER BOX)   | SLK02 826R | 2    | Non needed                                | 25 Mtr                   | 6            | BLACK            |
| SILK BLK 2-0 X 25MTR<br>(6 EA PER BOX) | SLK20 823R | 2-0  | Non needed                                | 25 Mtr                   | 6            | BLACK            |
| SILK BLK 3-0 X 25MTR<br>(6 EA PER BOX) | SLK30 822R | 3-0  | Non needed                                | 25 Mtr                   | 6            | BLACK            |

# MERISTEEL™

## Steel Suture

### Indications

- Orthopedic procedures
- Sternum closure

### Advantages

- Exceptional tensile strength
- Excellent tissue compatibility



### Product Overview

|                      |                                       |
|----------------------|---------------------------------------|
| Structure            | : Monofilament                        |
| Colour               | : Silver/Steel                        |
| Chemical composition | : Stainless steel                     |
| Sizes                | : USP 5 (metric 1) - USP 6 (metric 8) |
| Sterilization        | : Gamma radiation                     |

| Description                            | FG Code   | Size |  | Needle Profile Shape, Type & Length   | Suture Length (cm) | Pack Size | Suture Colour | Description                              | FG Code   | Size |  | Needle Profile Shape, Type & Length                    | Suture Length (cm) | Pack Size | Suture Colour |
|----------------------------------------|-----------|------|--|---------------------------------------|--------------------|-----------|---------------|------------------------------------------|-----------|------|--|--------------------------------------------------------|--------------------|-----------|---------------|
| MERISTEEL (STEEL SUTURE)               |           |      |  |                                       |                    |           |               | MERISTEEL (STEEL SUTURE)                 |           |      |  |                                                        |                    |           |               |
| STEEL WIRE 0 X 45-17MM HC BP (45X2)    | STC600    | 0    |  | 17 mm 1/2 Circle Blunt Point          | 45                 | 12        | -             | STEEL WIRE 5 X 45-48MM HC TC (45X4)      | STC05650  | 5    |  | 48 mm 1/2 Circle Taper Cut                             | 45                 | 12        | -             |
| STEEL WIRE 1 X 45-40MM HC TC (45X4)    | STC660    | 1    |  | 40 mm 1/2 Circle Taper Cut            | 45                 | 12        | -             | STEEL WIRE 6 X 45-48MM HC CT (45X4)      | STC654    | 6    |  | 48 mm 1/2 Circle Conventional Cutting                  | 45                 | 12        | -             |
| STEEL WIRE 1 X 45-40MM HC BP (45X4)    | STC661    | 1    |  | 40 mm 1/2 Circle Blunt Point          | 45                 | 12        | -             | STEEL WIRE 6 X 45-44MM HC RB BP (45X4)   | STC9456   | 6    |  | 44 mm 1/2 Circle Blunt Point                           | 45                 | 6         | -             |
| STEEL WIRE 1 X 45-35MM HC RB BLUNT     | STC9435   | 1    |  | 35 mm 1/2 Circle Blunt Point          | 45                 | 12        | -             | STEEL WIRE 6 X 75-44MM HC RB BP (75X2)   | STC9494   | 6    |  | 44 mm 1/2 Circle Blunt Point                           | 75                 | 6         | -             |
| STEEL WIRE 2 X 45-40MM HC TC (45X4)    | STC650    | 2    |  | 40 mm 1/2 Circle Taper Cut            | 45                 | 12        | -             | STEEL WIRE 6 X 45-48MM HC CT (45X1)      | STC9654   | 6    |  | 48 mm 1/2 Circle Conventional Cutting                  | 45                 | 12        | -             |
| STEEL WIRE 2 X 45-40MM HC BP (45X4)    | STC662    | 2    |  | 40 mm 1/2 Circle Blunt Point          | 45                 | 12        | -             | STEEL WIRE 6 X 45-48MM CU CT (45X6)      | STC9654CU | 6    |  | 48 mm 3/8 Circle Conventional Cutting                  | 45                 | 12        | -             |
| STEEL WIRE 2 X 45-40MM HC RB           | STC9540   | 2    |  | 40 mm 1/2 Circle Round Body           | 45                 | 12        | -             | SS WIRE 6X45-48 HCCT (45X1)44 HCBP(45X3) | STC700    | 6    |  | 48/44 mm 1/2 Circle Conventional Cutting   Blunt Point | 45                 | 12        | -             |
| STEEL WIRE 4 X 45-40MM HC BP (45X4)    | STC664    | 4    |  | 40 mm 1/2 Circle Blunt Point          | 45                 | 12        | -             | STEEL WIRE 6 X 75-65MM HC RB BLUNT       | STC9656   | 6    |  | 65 mm 1/2 Circle Blunt Point                           | 75                 | 12        | -             |
| STEEL WIRE 4 X 45-48MM HC CT (45X4)    | STC652    | 4    |  | 48 mm 1/2 Circle Conventional Cutting | 45                 | 12        | -             | STEEL WIRE 7 X 75-50MM HC RC (75X2)      | STC070624 | 7    |  | 50 mm 1/2 Circle Reverse Cutting                       | 75                 | 12        | -             |
| STEEL WIRE 4 X 75-48MM HC TC (75X2)    | STC040651 | 4    |  | 48 mm 1/2 Circle Taper Cut            | 75                 | 12        | -             |                                          |           |      |  |                                                        |                    |           |               |
| STEEL WIRE 5 X 75-55MM HC RC(75X2)     | STC050945 | 5    |  | 55 mm 1/2 Circle Reverse Cutting      | 75                 | 12        | -             |                                          |           |      |  |                                                        |                    |           |               |
| STEEL WIRE 5 X 45-48MM HC CT (45X4)    | STC653    | 5    |  | 48 mm 1/2 Circle Conventional Cutting | 45                 | 12        | -             |                                          |           |      |  |                                                        |                    |           |               |
| STEEL WIRE 5 X 45-44MM HC RB BP (45X4) | STC9455   | 5    |  | 44 mm 1/2 Circle Blunt Point          | 45                 | 6         | -             |                                          |           |      |  |                                                        |                    |           |               |
| STEEL WIRE 5 X 75-44MM HC RB BP (75X2) | STC9495   | 5    |  | 44 mm 1/2 Circle Blunt Point          | 75                 | 6         | -             |                                          |           |      |  |                                                        |                    |           |               |
| STEEL WIRE 5 X 45-48MM HC RB BP (45X4) | STC9497   | 5    |  | 48 mm 1/2 Circle Blunt Point          | 45                 | 12        | -             |                                          |           |      |  |                                                        |                    |           |               |
| STEEL WIRE 5 X 45-48MM HC CT (45X1)    | STC9653   | 5    |  | 48 mm 1/2 Circle Conventional Cutting | 45                 | 12        | -             |                                          |           |      |  |                                                        |                    |           |               |

# FILAMIDE™

## Polyamide Suture

### Indications

- Skin closure (intracutaneous, subcutaneous)
- Plastic and reconstructive procedures
- Microsurgery procedures
- Ophthalmic procedures

### Advantages

- Flexible, easy to handle and tie
- Smooth passage through tissue
- Excellent histocompatibility



### Product Overview

|                      |                                            |
|----------------------|--------------------------------------------|
| Structure            | : Monofilament                             |
| Colour               | : Black                                    |
| Chemical composition | : Polyamide 6 - 6.6                        |
| Sizes                | : USP 10/0 (metric 0.2) - USP 2 (metric 5) |
| Sterilization        | : Gamma radiation                          |



Meril Endo-Surgery Product Catalogue | **Wound Closure** 117

Meril Endo-Surgery Product Catalogue | **Wound Closure** 119

Meril Endo-Surgery Product Catalogue | **Wound Closure** 121

Meril Endo-Surgery Product Catalogue | **Wound Closure** 123





| Description                   | FG Code   | Size | Needle Profile Shape, Type & Length | Suture Length (cm) | Pack Size | Suture Colour |
|-------------------------------|-----------|------|-------------------------------------|--------------------|-----------|---------------|
| FILAMIDE (POLYAMIDE SUTURE)   |           |      |                                     |                    |           |               |
| POLYAMIDE BLK 0 X 76 (76X2)   | NYS10 904 | 0    | Non needled                         | 76                 | 12        | BLACK         |
| POLYAMIDE BLK 1 X 76 (76X2)   | NYS01 905 | 1    | Non needled                         | 76                 | 12        | BLACK         |
| POLYAMIDE BLK 1 X 750 (75X10) | NYS011075 | 1    | Non needled                         | 75                 | 12        | BLACK         |
| POLYAMIDE BLK 1 X 35 (35X5)   | NYS01 705 | 1    | Non needled                         | 35                 | 12        | BLACK         |
| POLYAMIDE BLK 1 X 100 (100X2) | NYS01185  | 1    | Non needled                         | 100                | 12        | BLACK         |
| POLYAMIDE BLK 2 X 76 (76X2)   | NYS02 906 | 2    | Non needled                         | 76                 | 12        | BLACK         |
| POLYAMIDE BLK 2-0 X 76 (76X2) | NYS20 903 | 2-0  | Non needled                         | 76                 | 12        | BLACK         |
| POLYAMIDE BLK 2-0 X 35 (35X5) | NYS20 703 | 2-0  | Non needled                         | 35                 | 12        | BLACK         |

# FILAPROP™

## Polypropylene Suture

### Indications

- Skin closure (intracutaneous, subcutaneous)
- Cardiovascular procedures
- Neurological procedures
- Ophthalmic procedures

### Advantages

- Smooth passage through tissue
- Minimal acute inflammatory reaction in tissues
- Unsurpassed strength
- Optimal elasticity and elongation properties
- Available in Less Memory Tray (LMT)











### Product Overview






|                      |                                            |
|----------------------|--------------------------------------------|
| Structure            | : Monofilament                             |
| Colour               | : Blue                                     |
| Chemical composition | : Polypropylene                            |
| Sizes                | : USP 10/0 (metric 0.2) - USP 2 (metric 5) |
| Sterilization        | : Ethylene oxide                           |

| Description                             | FG Code       | Size |                                                                                     | Needle Profile Shape, Type & Length      | Suture Length (cm) | Pack Size | Suture Colour | Description                              | FG Code      | Size |                                                                                       | Needle Profile Shape, Type & Length      | Suture Length (cm) | Pack Size | Suture Colour |
|-----------------------------------------|---------------|------|-------------------------------------------------------------------------------------|------------------------------------------|--------------------|-----------|---------------|------------------------------------------|--------------|------|---------------------------------------------------------------------------------------|------------------------------------------|--------------------|-----------|---------------|
| FILAPROP (POLYPROPYLENE SUTURE)         |               |      |                                                                                     |                                          |                    |           |               | FILAPROP (POLYPROPYLENE SUTURE)          |              |      |                                                                                       |                                          |                    |           |               |
| POLYPROPYLENE BU 0 X 90-30MM CU RC      | PPL102815     | 0    |    | 30 mm 3/8 Circle Reverse Cutting         | 90                 | 12        | BLUE          | POLYPROPYLENE BU 0 X 70-30MM HC RB HB    | PPL10 842    | 0    |    | 30 mm 1/2 Circle Round Body              | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 70-30MM HC RB      | PPL10 830     | 0    |    | 30 mm 1/2 Circle Round Body              | 70                 | 12        | BLUE          | POLYPROPYLENE BU 0 X 100-26MM HC RB DN   | PPL108026XDN | 0    |    | 26 mm 1/2 Circle Round Body   Round Body | 100                | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 90-30MM HC RB      | PPL10 846     | 0    |    | 30 mm 1/2 Circle Round Body              | 90                 | 12        | BLUE          | POLYPROPYLENE BU 0 X 75-37MM HC TC       | PPL108037    | 0    |    | 37 mm 1/2 Circle Taper Cut               | 75                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 90-30MM HC TC      | PPL10 894     | 0    |    | 30 mm 1/2 Circle Taper Cut               | 90                 | 12        | BLUE          | POLYPROPYLENE BU 0 X 75-38MM CU CT       | PPL108038    | 0    |    | 38 mm 3/8 Circle Conventional Cutting    | 75                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 75-40MM HC RB      | PPL108434     | 0    |    | 40 mm 1/2 Circle Round Body              | 75                 | 12        | BLUE          | POLYPROPYLENE BU 1 X 90-45MM HC TC       | PPL018450    | 1    |    | 45 mm 1/2 Circle Taper Cut               | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 75-37MM HC RB      | PPL108424     | 0    |    | 37 mm 1/2 Circle Round Body              | 75                 | 12        | BLUE          | POLYPROPYLENE BU 1 X 90-50MM HC RC       | PPL018475    | 1    |    | 50 mm 1/2 Circle Reverse Cutting         | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 75-30MM HC RB      | PPL100830     | 0    |    | 30 mm 1/2 Circle Round Body              | 75                 | 12        | BLUE          | POLYPROPYLENE BU 1 X 150-48MM HC RB LOOP | PPL01 834    | 1    |   | 48 mm 1/2 Circle Round Body              | 150                | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 75-24MM CU CT      | PPL108024     | 0    |    | 24 mm 3/8 Circle Conventional Cutting    | 75                 | 12        | BLUE          | POLYPROPYLENE BU 1 X 70-30MM HC RB       | PPL01 840    | 1    |  | 30 mm 1/2 Circle Round Body              | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 75-26MM HC RB      | PPL108026     | 0    |  | 26 mm 1/2 Circle Round Body              | 75                 | 12        | BLUE          | POLYPROPYLENE BU 1 X 70-40MM HC RB       | PPL01 843    | 1    |  | 40 mm 1/2 Circle Round Body              | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 75-45MM HC RC      | PPL108045     | 0    |  | 45 mm 1/2 Circle Reverse Cutting         | 75                 | 12        | BLUE          | POLYPROPYLENE BU 1 X 100-45MM HC RC      | PPL01 883    | 1    |  | 45 mm 1/2 Circle Reverse Cutting         | 100                | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 75-45MM HC RB      | PPL108075     | 0    |  | 45 mm 1/2 Circle Round Body              | 75                 | 12        | BLUE          | POLYPROPYLENE BU 1 X 90-90MM CU CT       | PPL012989    | 1    |  | 90 mm 3/8 Circle Conventional Cutting    | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 90-26MM HC RB DN   | PPL10 824MIPH | 0    |  | 26 mm 1/2 Circle Round Body   Round Body | 90                 | 1         | BLUE          | POLYPROPYLENE BU 1 X 100-40MM HC RC      | PPL01884     | 1    |  | 40 mm 1/2 Circle Reverse Cutting         | 100                | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 100-40MM HC RB     | PPL108434XL   | 0    |  | 40 mm 1/2 Circle Round Body              | 100                | 12        | BLUE          | POLYPROPYLENE BU 1 X 100-36MM HC RB      | PPL01842     | 1    |  | 36 mm 1/2 Circle Round Body              | 100                | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 100-31MM J TYPE TC | PPL10982J     | 0    |  | 31 mm J Shape Taper Cut                  | 100                | 12        | BLUE          | POLYPROPYLENE BU 1 X 90-90MM CU RC       | PPL018190    | 1    |  | 90 mm 3/8 Circle Reverse Cutting         | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 100-31MM HC RB     | PPL108430     | 0    |  | 31 mm 1/2 Circle Round Body              | 100                | 12        | BLUE          | POLYPROPYLENE BU 1 X 75-30MM HCRB TRODN  | PPL018130    | 1    |  | 30 mm 1/2 Circle Round Body   TRO        | 75                 | 12        | BLUE          |
| POLYPROPYLENE BU 0 X 90-40MM HC RC      | PPL108470     | 0    |  | 40 mm 1/2 Circle Reverse Cutting         | 90                 | 12        | BLUE          | POLYPROPYLENE BU 1 X 100-40MM HC RB      | PPL01 843DS  | 1    |  | 40 mm 1/2 Circle Round Body              | 100                | 12        | BLUE          |

Meril Endo-Surgery Product Catalogue | **Wound Closure** 131



| Description                                 | FG Code      | Size |                                                                                     | Needle Profile Shape, Type & Length      | Suture Length (cm) | Pack Size | Suture Colour | Description                                 | FG Code    | Size |                                                                                       | Needle Profile Shape, Type & Length      | Suture Length (cm) | Pack Size | Suture Colour |
|---------------------------------------------|--------------|------|-------------------------------------------------------------------------------------|------------------------------------------|--------------------|-----------|---------------|---------------------------------------------|------------|------|---------------------------------------------------------------------------------------|------------------------------------------|--------------------|-----------|---------------|
| FILAPROP (POLYPROPYLENE SUTURE)             |              |      |                                                                                     |                                          |                    |           |               | FILAPROP (POLYPROPYLENE SUTURE)             |            |      |                                                                                       |                                          |                    |           |               |
| POLYPROPYLENE BU<br>2-0 X 90-31MM HC RB DN  | PPL208526    | 2-0  |    | 31 mm 1/2 Circle Round Body   Round Body | 90                 | 12        | BLUE          | POLYPROPYLENE BU<br>2-0 X 75-70MM ST RB DN  | PPL208400  | 2-0  |    | 70 mm Straight Round Body   Round Body   | 75                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 90-31MM HC TC DN  | PPL208850    | 2-0  |    | 31 mm 1/2 Circle Taper Cut   Taper Cut   | 90                 | 12        | BLUE          | POLYPROPYLENE BU<br>2-0 X 90-30MM HC RB DN  | PPL20844DN | 2-0  |    | 30 mm 1/2 Circle Round Body   Round Body | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 45-26MM CU RC     | PPL208026    | 2-0  |    | 26 mm 3/8 Circle Reverse Cutting         | 45                 | 12        | BLUE          | POLYPROPYLENE BU<br>2-0 X 90-17MM HC TC     | PPL208938  | 2-0  |    | 17 mm 1/2 Circle Taper Cut               | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 90-40MM HC TC DN  | PPL208852DN  | 2-0  |    | 40 mm 1/2 Circle Taper Cut   Taper Cut   | 90                 | 12        | BLUE          | POLYPROPYLENE BU<br>2-0 X 100-40MM CU CT    | PPL208690  | 2-0  |    | 40 mm 3/8 Circle Conventional Cutting    | 100                | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 100-75MM ST CT    | PPL200752    | 2-0  |    | 75 mm Straight Conventional Cutting      | 100                | 12        | BLUE          | POLYPROPYLENE BU<br>3-0 X 70-22MM CU CT     | PPL30 018  | 3-0  |    | 22 mm 3/8 Circle Conventional Cutting    | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 90-55MM ST CT     | PPL208571    | 2-0  |   | 55 mm Straight Conventional Cutting      | 90                 | 12        | BLUE          | POLYPROPYLENE BU<br>3-0 X 90-25MM CU CT     | PPL30 800  | 3-0  |    | 25 mm 3/8 Circle Conventional Cutting    | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 90-40MM HC RB DN  | PPL20 2040DN | 2-0  |  | 40 mm 1/2 Circle Round Body   Round Body | 90                 | 12        | BLUE          | POLYPROPYLENE BU<br>3-0 X 90-25MM HC RB     | PPL30 825  | 3-0  |    | 25 mm 1/2 Circle Round Body              | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 70-30MM HC RB     | PPL20 841    | 2-0  |  | 30 mm 1/2 Circle Round Body              | 70                 | 12        | BLUE          | POLYPROPYLENE BU<br>3-0 X 70-25MM HC RB     | PPL30 838  | 3-0  |    | 25 mm 1/2 Circle Round Body              | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 90-90MM CU CT     | PPL200988    | 2-0  |  | 90 mm 3/8 Circle Conventional Cutting    | 90                 | 12        | BLUE          | POLYPROPYLENE BU<br>3-0 X 70-20MM HC RB     | PPL30 887  | 3-0  |  | 20 mm 1/2 Circle Round Body              | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 90-30MM HC RB-3D  | PPL20844Z    | 2-0  |  | 30 mm 1/2 Circle Round Body              | 90                 | 36        | BLUE          | POLYPROPYLENE BU<br>3-0 X 70-60MM ST CT     | PPL30 888  | 3-0  |  | 60 mm Straight Conventional Cutting      | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 90-26MM HCTCDN-3D | PPL208977Z   | 2-0  |  | 26 mm 1/2 Circle Taper Cut   Taper Cut   | 90                 | 36        | BLUE          | POLYPROPYLENE BU<br>3-0 X 90-26MM HC RB DN  | PPL30 8522 | 3-0  |  | 26 mm 1/2 Circle Round Body   Round Body | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X90-31MM HCRBVB-3D  | PPL20295Z    | 2-0  |  | 31 mm 1/2 Circle Round Body              | 90                 | 36        | BLUE          | POLYPROPYLENE BU<br>3-0 X 90-17.5MM HCTC DN | PPL30 8936 | 3-0  |  | 17.5 mm 1/2 Circle Taper Cut   Taper Cut | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 90-26MM HC TC     | PPL208926    | 2-0  |  | 26 mm 1/2 Circle Taper Cut               | 90                 | 12        | BLUE          | POLYPROPYLENE BU<br>3-0 X 90-17MM HC RB     | PPL308558  | 3-0  |  | 17 mm 1/2 Circle Round Body              | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU<br>2-0 X 90-30MM HC RC     | PPL208931    | 2-0  |  | 30 mm 1/2 Circle Reverse Cutting         | 90                 | 12        | BLUE          | POLYPROPYLENE BU<br>3-0 X 90-16MM HC TC DN  | PPL308667  | 3-0  |  | 16 mm 1/2 Circle Taper Cut   Taper Cut   | 90                 | 12        | BLUE          |
|                                             |              |      |                                                                                     |                                          |                    |           |               | POLYPROPYLENE BU<br>3-0 X 75-22MM HC RB     | PPL308770  | 3-0  |  | 22 mm 1/2 Circle Round Body              | 75                 | 12        | BLUE          |
|                                             |              |      |                                                                                     |                                          |                    |           |               | POLYPROPYLENE BU<br>3-0 X 90-26MM HC TC DN  | PPL308976  | 3-0  |  | 26 mm 1/2 Circle Taper Cut   Taper Cut   | 90                 | 12        | BLUE          |
|                                             |              |      |                                                                                     |                                          |                    |           |               | POLYPROPYLENE BU<br>3-0 X 45-19MM CU RC     | PPL308665  | 3-0  |  | 19 mm 3/8 Circle Reverse Cutting         | 45                 | 12        | BLUE          |

| Description                              |            |                                     |                                                                                     |                                          |               |         |      | Description                               |                    |           |                                                                                       |                                          |    |    |      |
|------------------------------------------|------------|-------------------------------------|-------------------------------------------------------------------------------------|------------------------------------------|---------------|---------|------|-------------------------------------------|--------------------|-----------|---------------------------------------------------------------------------------------|------------------------------------------|----|----|------|
| FG Code                                  | Size       | Needle Profile Shape, Type & Length | Suture Length (cm)                                                                  | Pack Size                                | Suture Colour | FG Code | Size | Needle Profile Shape, Type & Length       | Suture Length (cm) | Pack Size | Suture Colour                                                                         |                                          |    |    |      |
| FILAPROP (POLYPROPYLENE SUTURE)          |            |                                     |                                                                                     |                                          |               |         |      | FILAPROP (POLYPROPYLENE SUTURE)           |                    |           |                                                                                       |                                          |    |    |      |
| POLYPROPYLENE BU 3-0 X 45-24MM CU RC     | PPL308684  | 3-0                                 |    | 24 mm 3/8 Circle Reverse Cutting         | 45            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 90-31MM HC RB DN   | PPL308525DN        | 3-0       |    | 31 mm 1/2 Circle Round Body   Round Body | 90 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 75-26MM HC RB     | PPL308832  | 3-0                                 |    | 26 mm 1/2 Circle Round Body              | 75            | 12      | BLUE |                                           |                    |           |                                                                                       |                                          |    |    |      |
| POLYPROPYLENE BU 3-0 X 75-19MM CU RC     | PPL308687  | 3-0                                 |    | 19 mm 3/8 Circle Reverse Cutting         | 75            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 45-22MM HC CT      | PPL308731          | 3-0       |    | 22 mm 1/2 Circle Conventional Cutting    | 45 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 90-16MM HC RB     | PPL308559  | 3-0                                 |    | 16 mm 1/2 Circle Round Body              | 90            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 45-19MM CU CT      | PPL308878          | 3-0       |    | 19 mm 3/8 Circle Conventional Cutting    | 45 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 90-25MM HC TC     | PPL308977  | 3-0                                 |    | 25 mm 1/2 Circle Taper Cut               | 90            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 90-26MM HC RB DNV  | PPL308355          | 3-0       |    | 26 mm 1/2 Circle Round Body   Round Body | 90 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 90-30MM HC RB     | PPL308525  | 3-0                                 |    | 30 mm 1/2 Circle Round Body              | 90            | 12      | BLUE |                                           |                    |           |                                                                                       |                                          |    |    |      |
| POLYPROPYLENE BU 3-0 X 90-16MM CU RC     | PPL308549  | 3-0                                 |    | 16 mm 3/8 Circle Reverse Cutting         | 90            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 90-31MM HC TC DN   | PPL308849DN        | 3-0       |    | 31 mm 1/2 Circle Taper Cut   Taper Cut   | 90 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 90-26MM CU RC     | PPL308685  | 3-0                                 |    | 26 mm 3/8 Circle Reverse Cutting         | 90            | 12      | BLUE |                                           |                    |           |                                                                                       |                                          |    |    |      |
| POLYPROPYLENE BU 3-0 X 90-17MM HC TC     | PPL308317  | 3-0                                 |    | 17 mm 1/2 Circle Taper Cut               | 90            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 90-40MM HC TC DN   | PPL308851          | 3-0       |   | 40 mm 1/2 Circle Taper Cut   Taper Cut   | 90 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 90-36MM HC RB     | PPL308336  | 3-0                                 |   | 36 mm 1/2 Circle Round Body              | 90            | 12      | BLUE |                                           |                    |           |                                                                                       |                                          |    |    |      |
| POLYPROPYLENE BU 3-0 X 75-22MM HC RB DN  | PPL308771  | 3-0                                 |  | 22 mm 1/2 Circle Round Body   Round Body | 75            | 12      | BLUE | POLYPROPYLENE BU 3-0 X90-26MM HCRB TRODN  | PPL308833          | 3-0       |  | 26 mm 1/2 Circle Round Body   TRO        | 90 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 75-60MM ST CT     | PPL3013413 | 3-0                                 |  | 60 mm Straight Conventional Cutting      | 75            | 12      | BLUE |                                           |                    |           |                                                                                       |                                          |    |    |      |
| POLYPROPYLENE BU 3-0 X 75-16MM CU CT     | PPL3013453 | 3-0                                 |  | 16 mm 3/8 Circle Conventional Cutting    | 75            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 90-20MM HC RBTC DN | PPL309557          | 3-0       |  | 20 mm 1/2 Circle Round Body   Taper Cut  | 90 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 120-26MM HC RB    | PPL308534  | 3-0                                 |  | 26 mm 1/2 Circle Round Body              | 120           | 12      | BLUE |                                           |                    |           |                                                                                       |                                          |    |    |      |
| POLYPROPYLENE BU 3-0 X 75-24MM CU RC     | PPL307694  | 3-0                                 |  | 24 mm 3/8 Circle Reverse Cutting         | 75            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 90-17MM HC RB DN   | PPL308558DN        | 3-0       |  | 17 mm 1/2 Circle Round Body   Round Body | 90 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 90-36MM HC RB DN  | PPL307499  | 3-0                                 |  | 36 mm 1/2 Circle Round Body   Round Body | 90            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 90-16MM CU RC-3D   | PPL308549Z         | 3-0       |  | 16 mm 3/8 Circle Reverse Cutting         | 90 | 36 | BLUE |
| POLYPROPYLENE BU 3-0 X 120-26MM HC TC DN | PPL307778  | 3-0                                 |  | 26 mm 1/2 Circle Taper Cut   Taper Cut   | 120           | 12      | BLUE | POLYPROPYLENE BU 3-0 X 90-60MM ST CT      | PPL300621          | 3-0       |  | 60 mm Straight Conventional Cutting      | 90 | 12 | BLUE |
| POLYPROPYLENE BU 3-0 X 45-26MM CU RC     | PPL308025  | 3-0                                 |  | 26 mm 3/8 Circle Reverse Cutting         | 45            | 12      | BLUE | POLYPROPYLENE BU 3-0 X 45-26MM CU CT EX   | PPL308021          | 3-0       |  | 26 mm 3/8 Circle Conventional Cutting    | 45 | 12 | BLUE |
|                                          |            |                                     |                                                                                     |                                          |               |         |      | POLYPROPYLENE BU 3-0 X 75-36MM CU CT      | PPL300626          | 3-0       |  | 36 mm 3/8 Circle Conventional Cutting    | 75 | 12 | BLUE |
|                                          |            |                                     |                                                                                     |                                          |               |         |      | POLYPROPYLENE BU 4-0 X 90-19MM CU RC      | PPL408015          | 4-0       |  | 19 mm 3/8 Circle Reverse Cutting         | 90 | 12 | BLUE |
|                                          |            |                                     |                                                                                     |                                          |               |         |      | POLYPROPYLENE BU 4-0 X 90-16MM CU CT      | PPL408875          | 4-0       |  | 16 mm 3/8 Circle Conventional Cutting    | 90 | 12 | BLUE |


134

Wound Closure | Meril Endo-Surgery Product Catalogue

Meril Endo-Surgery Product Catalogue |

Wound Closure

135

| Description                              | FG Code    | Size |                                                                                     | Needle Profile Shape, Type & Length                         | Suture Length (cm) | Pack Size | Suture Colour | Description                              | FG Code      | Size |                                                                                       | Needle Profile Shape, Type & Length     | Suture Length (cm) | Pack Size | Suture Colour |
|------------------------------------------|------------|------|-------------------------------------------------------------------------------------|-------------------------------------------------------------|--------------------|-----------|---------------|------------------------------------------|--------------|------|---------------------------------------------------------------------------------------|-----------------------------------------|--------------------|-----------|---------------|
| FILAPROP (POLYPROPYLENE SUTURE)          |            |      |                                                                                     |                                                             |                    |           |               | FILAPROP (POLYPROPYLENE SUTURE)          |              |      |                                                                                       |                                         |                    |           |               |
| POLYPROPYLENE BU 4-0 X 90-25MM CU CT     | PPL408879  | 4-0  |    | 25 mm 3/8 Circle Conventional Cutting                       | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 90-24MM CU RC     | PPL408684    | 4-0  |    | 24 mm 3/8 Circle Reverse Cutting        | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 70-16MM HC RB     | PPL40 849  | 4-0  |    | 16 mm 1/2 Circle Round Body                                 | 70                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 70-19MM CU RC     | PPL408682XL  | 4-0  |    | 19 mm 3/8 Circle Reverse Cutting        | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 70-16MM CU CT     | PPL40 870  | 4-0  |    | 16 mm 3/8 Circle Conventional Cutting                       | 70                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X90-25MM HC TC      | PPL408425    | 4-0  |    | 25 mm 1/2 Circle Taper Cut              | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X90-25MM HC TC DNC1 | PPL40 867  | 4-0  |    | 25 mm 1/2 Circle Taper Cut  Taper Cut                       | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 90-16MM CU RC     | PPL408416    | 4-0  |    | 16 mm 3/8 Circle Reverse Cutting        | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-17.5MM HCTC DN | PPL40 8935 | 4-0  |    | 17.5 mm 1/2 Circle Taper Cut  Taper Cut                     | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 75-22MM HC RB DN  | PPL40 8526DN | 4-0  |    | 22 mm 1/2 Circle Round Body  Round Body | 75                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-26MM HC RB DN  | PPL40 8521 | 4-0  |    | 26 mm 1/2 Circle Round Body  Round Body                     | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 70-17MM HC TC     | PPL40 835    | 4-0  |    | 17 mm 1/2 Circle Taper Cut              | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-16MM HC RB DN  | PPL40 8557 | 4-0  |    | 16 mm 1/2 Circle Round Body  Round Body                     | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 45-16MM CU RC     | PPL408884    | 4-0  |    | 16 mm 3/8 Circle Reverse Cutting        | 45                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-26MM HC RB     | PPL408355  | 4-0  |   | 26 mm 1/2 Circle Round Body                                 | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 45-16MM CU CT     | PPL408007    | 4-0  |    | 16 mm 3/8 Circle Conventional Cutting   | 45                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-17MM HC RB     | PPL408558  | 4-0  |  | 17 mm 1/2 Circle Round Body                                 | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 45-26MM CU CT     | PPL408020SL  | 4-0  |    | 26 mm 3/8 Circle Conventional Cutting   | 45                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 75-13MM HC CT DN  | PPL408666  | 4-0  |  | 13 mm 1/2 Circle Conventional Cutting  Conventional Cutting | 75                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X90-26MM HC TC DNC2 | PPL40 8845   | 4-0  |  | 26 mm 1/2 Circle Taper Cut  Taper Cut   | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-20MM HC RB     | PPL408761  | 4-0  |  | 20 mm 1/2 Circle Round Body                                 | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 45-19MM CU CT     | PPL408011    | 4-0  |  | 19 mm 3/8 Circle Conventional Cutting   | 45                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-20MM HC CT DN  | PPL408840  | 4-0  |  | 20 mm 1/2 Circle Conventional Cutting  Conventional Cutting | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 90-17MM HC TC     | PPL40 857    | 4-0  |  | 17 mm 1/2 Circle Taper Cut              | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-26MM HC CT DN  | PPL408846  | 4-0  |  | 26 mm 1/2 Circle Conventional Cutting  Conventional Cutting | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 45-12MM CU RC     | PPL40 879    | 4-0  |  | 12 mm 3/8 Circle Reverse Cutting        | 45                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 45-19MM CU RC     | PPL408682  | 4-0  |  | 19 mm 3/8 Circle Reverse Cutting                            | 45                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 90-16MM HC TC DN  | PPL408417    | 4-0  |  | 16 mm 1/2 Circle Taper Cut  Taper Cut   | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 45-20MM CU RC     | PPL408683  | 4-0  |  | 20 mm 3/8 Circle Reverse Cutting                            | 45                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 90-17MM HC TC EX  | PPL400835    | 4-0  |  | 17 mm 1/2 Circle Taper Cut              | 90                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-13MM HC RB DN  | PPL408204  | 4-0  |  | 13 mm 1/2 Circle Round Body  Round Body                     | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 70-17MM HC TC C1  | PPL40 899    | 4-0  |  | 17 mm 1/2 Circle Taper Cut              | 70                 | 12        | BLUE          |
| POLYPROPYLENE BU 4-0 X 90-18MM CU RB DN  | PPL408861  | 4-0  |  | 18 mm 3/8 Circle Round Body  Round Body                     | 90                 | 12        | BLUE          | POLYPROPYLENE BU 4-0 X 90-26MM HC RB EX  | PPL408831    | 4-0  |  | 26 mm 1/2 Circle Round Body             | 90                 | 12        | BLUE          |
|                                          |            |      |                                                                                     |                                                             |                    |           |               | POLYPROPYLENE BU 4-0 X 90-20MM HC TC DN  | PPL408841DN  | 4-0  |  | 20 mm 1/2 Circle Taper Cut  Taper Cut   | 90                 | 12        | BLUE          |
|                                          |            |      |                                                                                     |                                                             |                    |           |               | POLYPROPYLENE BU 4-0 X 70-16MM CU CT-3D  | PPL40870Z    | 4-0  |  | 16 mm 3/8 Circle Conventional Cutting   | 70                 | 36        | BLUE          |

Meril Endo-Surgery Product Catalogue | **Wound Closure** 139





Meril Endo-Surgery Product Catalogue | **Wound Closure** 143



Meril Endo-Surgery Product Catalogue | **Wound Closure** 147



# MERICRON XL<sup>TM</sup>

## Polyester Suture

### Indications

- Cardiovascular procedures
- Ophthalmic procedures
- Neurological procedures



### Advantages

- Smoother suture surface
- Less resistance to knot sliding
- Better strength
- Pliable structure









### Product Overview



|                      |                                           |
|----------------------|-------------------------------------------|
| Structure            | : Braided                                 |
| Colour               | : Green, White                            |
| Chemical composition | : Polyethylene terephthalate              |
| Coating              | : Polybutilate                            |
| Sizes                | : USP 6/0 (metric 0.7) - USP 5 (metric 7) |
| Sterilization        | : Gamma Radiation                         |

| Description                                | FG Code   | Size |                                                                                     | Needle Profile Shape, Type & Length    | Suture Length (cm) | Pack Size | Suture Colour | Description                                | FG Code  | Size |                                                                                       | Needle Profile Shape, Type & Length    | Suture Length (cm) | Pack Size | Suture Colour |
|--------------------------------------------|-----------|------|-------------------------------------------------------------------------------------|----------------------------------------|--------------------|-----------|---------------|--------------------------------------------|----------|------|---------------------------------------------------------------------------------------|----------------------------------------|--------------------|-----------|---------------|
| MERICRON XL (POLYESTER SUTURE)             |           |      |                                                                                     |                                        |                    |           |               | MERICRON XL (POLYESTER SUTURE)             |          |      |                                                                                       |                                        |                    |           |               |
| POLYESTER GR<br>0 X 90-50MM HC TC          | Me100979  | 0    |    | 50 mm 1/2 Circle Taper Cut             | 90                 | 12        | GREEN         | POLYESTER GREEN<br>1 X 90-36MM HC TC       | ME014407 | 1    |    | 36 mm 1/2 Circle Taper Cut             | 90                 | 12        | GREEN         |
| POLYESTER GR<br>0 X 90-26MM HC TC DN       | ME6978    | 0    |    | 26 mm 1/2 Circle Taper Cut   Taper Cut | 90                 | 12        | GREEN         | POLYESTER GR<br>2 X 100-45MM HC TC         | ME643    | 2    |    | 45 mm 1/2 Circle Taper Cut             | 100                | 12        | GREEN         |
| POLYESTER GR<br>0 X 75-22MM HC RC          | ME517G    | 0    |    | 22 mm 1/2 Circle Reverse Cutting       | 75                 | 12        | GREEN         | POLYESTER GR<br>2 X 75-45MM HC TC (4 in 1) | ME4843   | 2    |    | 45 mm 1/2 Circle Taper Cut             | 75                 | 6         | GREEN         |
| POLYESTER GR<br>0 X 90-36MM HC RB          | ME844G    | 0    |    | 36 mm 1/2 Circle Round Body            | 90                 | 12        | GREEN         | POLYESTER GR<br>2 X 75-45MM HC TC          | ME4842   | 2    |    | 45 mm 1/2 Circle Taper Cut             | 75                 | 12        | GREEN         |
| POLYESTER GR<br>0 X 75-30MM HC RB          | ME106030  | 0    |    | 30 mm 1/2 Circle Round Body            | 75                 | 12        | GREEN         | POLYESTER GR<br>2 X 90-37MM HC RB          | ME026237 | 2    |    | 37 mm 1/2 Circle Round Body            | 90                 | 12        | GREEN         |
| POLYESTER GR<br>0 X 75-26MM HC RB          | ME10361   | 0    |    | 26 mm 1/2 Circle Round Body            | 75                 | 12        | GREEN         | POLYESTER GR<br>5 X 75-55MM HC TC          | ME646    | 5    |    | 55 mm 1/2 Circle Taper Cut             | 75                 | 12        | GREEN         |
| POLYESTER GR<br>0 X 75-31MM HC RB          | ME100975  | 0    |    | 31 mm 1/2 Circle Round Body            | 75                 | 12        | GREEN         | POLYESTER GR<br>5 X 75-55MM HC TC (4 in 1) | ME4846   | 5    |    | 55 mm 1/2 Circle Taper Cut             | 75                 | 6         | GREEN         |
| POLYESTER GR<br>0 X 90-36MM HC RC          | ME1032049 | 0    |   | 36 mm 1/2 Circle Reverse Cutting       | 90                 | 12        | GREEN         | POLYESTER GR<br>2-0 X 90-17MM HC TC DN C2  | ME617    | 2-0  |   | 17 mm 1/2 Circle Taper Cut   Taper Cut | 90                 | 12        | GREEN         |
| POLYESTER GR<br>0 X 90-30MM CU RC          | ME106680  | 0    |  | 30 mm 3/8 Circle Reverse Cutting       | 90                 | 12        | GREEN         | POLYESTER WH<br>2-0 X 76-17MM HC TC        | ME673W   | 2-0  |  | 17 mm 1/2 Circle Taper Cut             | 76                 | 12        | WHITE         |
| POLYESTER GR<br>0 X 75-36MM CU RC          | ME107562  | 0    |  | 36 mm 3/8 Circle Reverse Cutting       | 75                 | 12        | GREEN         | POLYESTER WH<br>2-0 X 90-25MM HC TC DN     | ME687W   | 2-0  |  | 25 mm 1/2 Circle Taper Cut   Taper Cut | 90                 | 12        | WHITE         |
| POLYESTER GREEN<br>1 X 90-30MM HC RC       | ME0132047 | 1    |  | 30 mm 1/2 Circle Reverse Cutting       | 90                 | 12        | GREEN         | POLYESTER GR<br>2-0 X 90-17MM HC TC DN C1  | ME6937   | 2-0  |  | 17 mm 1/2 Circle Taper Cut   Taper Cut | 90                 | 12        | GREEN         |
| POLYESTER GREEN<br>1 X 75-22MM HC RC       | ME010518  | 1    |  | 22 mm 1/2 Circle Reverse Cutting       | 75                 | 12        | GREEN         | POLYESTER GR<br>2-0 X 90-26MM HC TC DN     | ME6977   | 2-0  |  | 26 mm 1/2 Circle Taper Cut   Taper Cut | 90                 | 12        | GREEN         |
| POLYESTER GREEN<br>X 90-37MM HC RB         | ME016137  | 1    |  | 37 mm 1/2 Circle Round Body            | 90                 | 12        | GREEN         | POLYESTER WH<br>2-0 X 90-26MM HC TC DN     | ME6987   | 2-0  |  | 26 mm 1/2 Circle Taper Cut   Taper Cut | 90                 | 12        | WHITE         |
| POLYESTER GREEN<br>1 X 90-30MM HC RB       | ME016130  | 1    |  | 30 mm 1/2 Circle Round Body            | 90                 | 12        | GREEN         | POLYESTER WH<br>2-0 X 90-17MM HC TC DN     | ME617W   | 2-0  |  | 17 mm 1/2 Circle Taper Cut   Taper Cut | 90                 | 12        | WHITE         |
| POLYESTER GREEN<br>1 X 75-31MM J TYPE RB H | ME010969J | 1    |  | 31 mm J Shape Round Body               | 75                 | 12        | GREEN         |                                            |          |      |                                                                                       |                                        |                    |           |               |
| POLYESTER GREEN<br>1 X 90-30MM HC RB-3D    | ME016130Z | 1    |  | 30 mm 1/2 Circle Round Body            | 90                 | 36        | GREEN         |                                            |          |      |                                                                                       |                                        |                    |           |               |



| Description                                    | FG Code    | Size |                                                                                     | Needle Profile Shape, Type & Length                           | Suture Length (cm) | Pack Size | Suture Colour | Description                                    | FG Code   | Size |                                                                                       | Needle Profile Shape, Type & Length          | Suture Length (cm) | Pack Size | Suture Colour   |
|------------------------------------------------|------------|------|-------------------------------------------------------------------------------------|---------------------------------------------------------------|--------------------|-----------|---------------|------------------------------------------------|-----------|------|---------------------------------------------------------------------------------------|----------------------------------------------|--------------------|-----------|-----------------|
| MERICRON XL (POLYESTER SUTURE)                 |            |      |                                                                                     |                                                               |                    |           |               | MERICRON XL (POLYESTER SUTURE)                 |           |      |                                                                                       |                                              |                    |           |                 |
| POLYESTER GR<br>4-0 X 76-8MM QC SM             | ME682      | 4-0  |    | 8 mm 1/4 Circle<br>Spatula Micropoint                         | 76                 | 12        | GREEN         | POLYESTER GREEN<br>3-0 X 75-17MM HC TC-3D      | ME306935Z | 3-0  |    | 17 mm 3/8 Circle<br>Reverse Cutting          | 75                 | 36        | GREEN           |
| POLYESTER GREEN<br>4-0 X 75-26MM HC RB         | ME406831   | 4-0  |    | 26 mm 1/2 Circle<br>Round Body                                | 75                 | 12        | GREEN         | POLYESTER GREEN<br>4-0 X 45-13MM CU RC-3D      | ME407573Z | 4-0  |    | 13 mm 3/8 Circle<br>Reverse Cutting          | 45                 | 36        | GREEN           |
| POLYESTER GR<br>4-0 X 90-17MM HC TC DN         | ME4013395  | 4-0  |    | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut                  | 90                 | 12        | GREEN         | POLYESTER GREEN<br>4-0 X 45-19MM CU RC-3D      | ME407142Z | 4-0  |    | 19 mm 3/8 Circle<br>Reverse Cutting          | 45                 | 36        | GREEN           |
| POLYESTER GREEN<br>4-0 X 75-20MM HC RB         | ME406761   | 4-0  |    | 20 mm 1/2 Circle<br>Round Body                                | 75                 | 12        | GREEN         | POLYESTER WHGR<br>2-0X75-25MM<br>HC TC DN TP6  | MME6577X  | 2-0  |    | 25 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 75                 | 6         | GREEN/<br>WHITE |
| POLYESTER GREEN<br>4-0 X 75-13MM<br>CU RB DN   | ME6891G    | 4-0  |    | 13 mm 3/8 Circle<br>Round Body  <br>Round Body                | 75                 | 12        | GREEN         | POLYESTER WHGR<br>2-0X75-17MM<br>HC TC DN TP6  | MME6555X  | 2-0  |    | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 75                 | 6         | GREEN/<br>WHITE |
| POLYESTER GREEN<br>4-0 X 90-17.5MM<br>HC RC DN | ME6935G    | 4-0  |    | 17.5 mm 1/2 Circle<br>Reverse Cutting  <br>Reverse Cutting    | 90                 | 12        | GREEN         | POLYESTER WHGR<br>2-0X75-25MM HC TC DN TP5     | MME6577   | 2-0  |    | 25 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 75                 | 6         | GREEN/<br>WHITE |
| POLYESTER GR<br>4-0 X 75-20MM<br>HC RB DN      | ME406761DN | 4-0  |    | 20 mm 1/2 Circle<br>Round Body  <br>Round Body                | 75                 | 12        | GREEN         | POLYESTER WHGR<br>2-0X75-17MM<br>HC TC DN TP5  | MME6555   | 2-0  |    | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 75                 | 6         | GREEN/<br>WHITE |
| POLYESTER GR<br>5-0 X 45-8MM QC SM DN          | ME683      | 5-0  |   | 8 mm 1/4 Circle<br>Spatula Micropoint  <br>Spatula Micropoint | 45                 | 12        | GREEN         | POLYESTER GR<br>2-0 X 90-17MM<br>HC TC DN TP   | PME617    | 2-0  |  | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 90                 | 12        | GREEN           |
| POLYESTER GREEN<br>5-0 X 60-17MM<br>HC RB DN   | ME506757   | 5-0  |  | 17 mm 1/2 Circle<br>Round Body  <br>Round Body                | 60                 | 12        | GREEN         | POLYESTER WH<br>2-0 X 90-17MM<br>HC TC DN TP   | PME617W   | 2-0  |  | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 90                 | 12        | WHITE           |
| POLYESTER GREEN<br>5-0 X 45-8MM<br>HC SM DN    | ME500883DN | 5-0  |  | 8 mm 1/2 Circle<br>Spatula Micropoint  <br>Spatula Micropoint | 45                 | 12        | GREEN         | POLYESTER GR<br>2-0 X 90-25MM<br>HC TC DN TP   | PME687    | 2-0  |  | 25 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 90                 | 12        | GREEN           |
| POLYESTER GREEN<br>5-0 X 75-13MM<br>HC RB DN   | ME6889     | 5-0  |  | 13 mm 1/2 Circle<br>Round Body  <br>Round Body                | 75                 | 12        | GREEN         | POLYESTER WH<br>2-0 X 90-25MM<br>HC TC DN TP   | PME687W   | 2-0  |  | 25 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 90                 | 12        | WHITE           |
| POLYESTER GREEN<br>5-0 X 75-13MM<br>CU RB DN   | ME6890G    | 5-0  |  | 13 mm 3/8 Circle<br>Round Body  <br>Round Body                | 75                 | 12        | GREEN         | POLYSTER XL GR<br>2-0 X 75-26MM<br>HC TC DN TP | ME202613  | 2-0  |  | 26 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 75                 | 12        | GREEN           |
| POLYESTER WH<br>5-0 X 45-8MM<br>QC SM DN       | ME50893DN  | 5-0  |  | 8 mm 1/4 Circle<br>Spatula Micropoint  <br>Spatula Micropoint | 45                 | 12        | WHITE         | POLYSTER XL GR<br>2-0 X 75-17MM<br>HC TC DN TP | ME200553  | 2-0  |  | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 75                 | 12        | GREEN           |
| POLYESTER WH<br>5-0 X 45-8MM<br>HC SM DN       | ME50883DN  | 5-0  |  | 8 mm 1/2 Circle<br>Spatula Micropoint  <br>Spatula Micropoint | 45                 | 12        | WHITE         | POLYESTER XL WHGR<br>2-0X75-17MM<br>HCTCDN TP3 | MME6555MP | 2-0  |  | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut | 75                 | 6         | GREEN/<br>WHITE |
| POLYESTER GR<br>6-0 X 75-11MM<br>CU RB DN      | ME6597     | 6-0  |  | 11 mm 3/8 Circle<br>Round Body  <br>Round Body                | 75                 | 12        | GREEN         |                                                |           |      |                                                                                       |                                              |                    |           |                 |
| POLYESTER GREEN<br>1 X 90-36MM HC RC - 3D      | ME012520Z  | 1    |  | 36 mm 1/2 Circle<br>Reverse Cutting                           | 90                 | 36        | GREEN         |                                                |           |      |                                                                                       |                                              |                    |           |                 |



| Description                                      | FG Code   | Size |                                                                                     | Needle Profile<br>Shape, Type &<br>Length      | Suture<br>Length<br>(cm) | Pack<br>Size | Suture<br>Colour |
|--------------------------------------------------|-----------|------|-------------------------------------------------------------------------------------|------------------------------------------------|--------------------------|--------------|------------------|
| MERICRON XL (POLYESTER SUTURE)                   |           |      |                                                                                     |                                                |                          |              |                  |
| POLYESTER XL WHGR<br>2-0X75-25MM HCTCDN TP3      | MME6577MP | 2-0  |    | 25 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut   | 75                       | 6            | GREEN/<br>WHITE  |
| POLYESTER XL GR<br>2-0 X 90-17MM<br>HC TCDN P3   | PME617MP  | 2-0  |    | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut   | 90                       | 12           | GREEN            |
| POLYESTER XL WH<br>2-0 X 90-17MM<br>HC TCDN P3   | PME617WMP | 2-0  |    | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut   | 90                       | 12           | WHITE            |
| POLYESTER GR<br>2-0 X 90-17MM<br>HC TC DN TP6    | PME6937   | 2-0  |    | 17 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut   | 90                       | 12           | GREEN            |
| POLYESTER GR<br>2-0 X 90-26MM<br>HC TC DN TP6    | PME6977G  | 2-0  |    | 26 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut   | 90                       | 12           | GREEN            |
| POLYESTER WH<br>2-0 X 90-26MM<br>HC TC DN TP6    | PME6987W  | 2-0  |   | 26 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut   | 90                       | 12           | WHITE            |
| POLYESTER XL WH/GR<br>2-0X90-25MM<br>HC TC DNTP5 | SM10B77   | 2-0  |  | 25 mm 1/2 Circle<br>Taper Cut  <br>Taper Cut   | 90                       | 6            | GREEN/<br>WHITE  |
| POLYESTER GR<br>5-0 X 75-13MM<br>HC RB DN TP3    | PME6889G  | 5-0  |  | 13 mm 1/2 Circle<br>Round Body  <br>Round Body | 75                       | 12           | GREEN            |

Note: .....

Note: \_\_\_\_\_

158 **Wound Closure** | Meril Endo-Surgery Product Catalogue

Note: \_\_\_\_\_

Meril Endo-Surgery Product Catalogue | **Wound Closure** 159