

**PAGINA 0**

**CREATA PENTRU A  
PUTEA SEMINA  
ELECTRONIC  
DOCUMENTUL**

## Balloon-expandable Peripheral Stent System

- Broad offering of 6 Fr compatible 0.035" balloon-expandable stent with radiopaque marker technology for optimized visibility
- Low crossing profile
- Minimal shortening for placement confidence


Visi-Pro™ Catheter lengths: 80 cm and 135 cm

Each system includes:  
One stent and delivery catheter system

### Compliance chart

Diameter (mm)	Inflation pressure (atm)				
	8	9	10	11	12
5.0	5.00 <sup>1</sup>	5.09	5.16	5.22	5.28 <sup>2</sup>
6.0	6.00 <sup>1</sup>	6.11	6.22	6.31	6.39 <sup>2</sup>
7.0			7.00 <sup>1</sup>	7.09	7.17 <sup>2</sup>
8.0			8.00 <sup>1</sup>	8.15	8.26 <sup>2</sup>
9.0			9.00 <sup>1</sup>	9.15	9.28 <sup>2</sup>
10.0			10.00 <sup>1</sup>	10.11	10.21 <sup>2</sup>

1. Diameter at nominal pressure
2. Diameter at rated burst pressure



## Order information

Product catalogue number		Stent dimensions		Balloon length (mm)	Recommended sheath size (F)	Recommended guidewire (inch)
Catheter length 80 cm	Catheter length 135 cm	Diameter (mm)	Length (mm)			
PXP35-05-12-080		5	12	15	6*	0.035
PXP35-05-17-080	PXP35-05-17-135	5.0	17	20	6*	0.035
PXP35-05-27-080	PXP35-05-27-135	5.0	27	30	6*	0.035
PXP35-05-37-080	PXP35-05-37-135	5.0	37	40	6*	0.035
PXP35-05-57-080	PXP35-05-57-135	5.0	57	60	6*	0.035
PXP35-06-12-080		6.0	12	15	6*	0.035
PXP35-06-17-080	PXP35-06-17-135	6.0	17	20	6*	0.035
PXP35-06-27-080	PXP35-06-27-135	6.0	27	30	6*	0.035
PXP35-06-37-080	PXP35-06-37-135	6.0	37	40	6*	0.035
PXP35-06-57-080	PXP35-06-57-135	6.0	57	60	6*	0.035
PXP35-07-12-080		7.0	12	15	6*	0.035
PXP35-07-17-080	PXP35-07-17-135	7.0	17	20	6*	0.035
PXP35-07-27-080	PXP35-07-27-135	7.0	27	30	6*	0.035
PXP35-07-37-080	PXP35-07-37-135	7.0	37	40	6*	0.035
PXP35-07-57-080	PXP35-07-57-135	7.0	57	60	6*	0.035
PXP35-08-17-080	PXP35-08-17-135	8.0	17	20	6*	0.035
PXP35-08-27-080	PXP35-08-27-135	8.0	27	30	6*	0.035
PXP35-08-37-080	PXP35-08-37-135	8.0	37	40	6*	0.035
PXP35-08-57-080	PXP35-08-57-135	8.0	57	60	6*	0.035
PXP35-09-17-080	PXP35-09-17-135	9.0	17	20	7	0.035
PXP35-09-27-080	PXP35-09-27-135	9.0	27	30	7	0.035
PXP35-09-37-080	PXP35-09-37-135	9.0	37	40	7	0.035
PXP35-09-57-080	PXP35-09-57-135	9.0	57	60	7	0.035
PXP35-10-17-080	PXP35-10-17-135	10.0	17	20	7	0.035
PXP35-10-27-080	PXP35-10-27-135	10.0	27	30	7	0.035
PXP35-10-37-080	PXP35-10-37-135	10.0	37	40	7	0.035
PXP35-10-57-080	PXP35-10-57-135	10.0	57	60	7	0.035

\* 6F= 0.085" I.D.