


fixLOCK Broad Plate, 4.5 mm

Item Code: 849

Indications

- Temporary internal fixation and stabilization of humerus, femur, and tibial shaft fractures
- Periprosthetic fractures
 - Comminuted fractures
 - Fractures in osteopenic bone
 - Non-unions
 - Mal-unions
 - Osteotomy fixation



Features

- Used over the tensile surface of bones to work on tension band principle
- Uniform stiffness of all segments, prolonging fatigue life of the implant
- Smaller area of post-fixation avascularity
- Dynamic compression unit offers 2 mm of axial compression
- The combination of conventional and locking screws offers optimum fixation regardless of bone density
- Available in both Titanium & Stainless Steel

Sizes Available

No. of Holes	Code	
	S.S.	Titanium
6	849.06	849.06T
7	849.07	849.07T
8	849.08	849.08T
9	849.09	849.09T
10	849.10	849.10T
11	849.11	849.11T
12	849.12	849.12T
14	849.14	849.14T
16	849.16	849.16T