




		F7952	F7952W	F7752	F7752W	F7952P(1)	F7952PW(1)	F7952PL	F7952WPL	F7752P(1)	F7752PW(1)	
<b>Physio-Control</b>	LIFEPAK 9 presa-socket QUICK-COMBO	●	●	●	●	●	●	-	-	●	●	 <p>F7952 F7952W F7752 F7752W F7952P F7952PW F7752P F7752PW</p>
	LIFEPAK 10	●	●	●	●	●	●	-	-	●	●	
	LIFEPAK 12	●	●	●	●	●	●	-	-	●	●	
	LIFEPAK 15	●	●	●	●	●	●	-	-	●	●	
	LIFEPAK 20	●	●	●	●	●	●	-	-	●	●	
	LIFEPAK 500 AED con presa grigia-with grey socket	●	●	●	●	-	-	-	-	-	-	
	LIFEPAK 500 AED con presa rosa-with pink socket	●	●	●	●	-	-	●	●	-	-	
	LIFEPAK 1000 AED	●	●	●	●	-	-	●	●	-	-	
	LIFEPAK CR PLUS AED	-	●	-	●	-	-	●	●	-	-	
LIFEPAK EXPRESS AED	-	●	-	●	-	-	●	●	-	-		
<b>Bexen Cardio</b>	REANIBEX-200	●	●	(*)	(*)	-	-	-	-	-	-	 <p>F7952PL F7952WPL</p>
	REANIBEX-300	●	●	(*)	(*)	-	-	-	-	-	-	
	REANIBEX-500	●	●	(*)	(*)	-	-	-	-	-	-	
	REANIBEX-700	●	●	(*)	(*)	●	●	-	-	(*)	(*)	
	REANIBEX-800	●	●	(*)	(*)	●	●	-	-	(*)	(*)	
<b>Cardioline</b>	ELIFE 700	●	●	●	●	●	●	-	-	●	●	
<b>Mindray</b>	BENEHEART D3	●	●	●	●	●	●	-	-	●	●	
	BENEHEART D6	●	●	●	●	●	●	-	-	●	●	

F7952PL, F7952WPL: riduzione di energia erogata al paziente: 200J sul defibrillatore diventano 50J sul paziente-reduced energy: defibrillator output of 200J is reduced to 50J

(\*) Non per stimolazione – Not suitable for pace

		F7953	F7953W	F7753	F7953P	F7953PW	F7753P	
<b>Agilent Hewlett Packard</b>	CODEMASTER 100	●	●	●	●	●	●	
	CODEMASTER XE	●	●	●	●	●	●	
	CODEMASTER XL	●	●	●	●	●	●	
	CODEMASTER XL+	●	●	●	●	●	●	
	HEARTSTREAM XL, XLT e-and XL+, con cavo-with cable M3507A	●	●	●	●	●	●	
<b>Laerdal Medical</b>	HEARTSTART 4000	●	●	●	●	●	●	