

# SKIN MARKER WITH RULER

## Product

Marking the patient's skin is important in a range of clinical settings and for different reasons but the principal functions can be to:

- ensure surgery is performed on the correct site
- target non-invasive treatments
- plan the surgical approach
- ensure natural realignment of skin flaps during wound closure

## MARKER



- ✓ Double head marker with regular and fine tip
- ✓ Large ink reservoir for longer product life
- ✓ Constant ink flow for clear marking
- ✓ Clip cap for regular tip and standard cap for fine tip

## Gentian violet ink

- ✓ Colour works on all ranges of pigmented skin
- ✓ Quick drying, indelible, biocompatible
- ✓ Non-irritating, non-toxic
- ✓ Non-smearing, does not run



**For quality  
&  
safety**


  
[www.sutures.be](http://www.sutures.be)



## RULER

Flexible 15 cm

In centimetres and inches

Ensures measurement even over contoured surfaces

## Features

- ✓ Single-use, individually packaged in sterile peel pouches
- ✓ Latex-free
- ✓ Tapered tip for ergonomic handling

## SMI tip

For best results, remove oil from the skin using alcohol before marking

## Technical description

Package of 36

14.1 cm long (with caps)

Weight 7 g

High Density Polyethylene (HDPE)

Product	Unit	Order no
Skin marker sterile, with ruler	36	ZSM20

