# **Technical Data Sheet**

Model No.: N62D01 BUDI-DI: 697442484L05QF

Product Name: Disposable Mob Cap

Classification:

Class I, according to Annex VIII of Regulation (EU) 2017/745

Category I, according to annex I of Regulation (EU) 2016/425

Description: Mob Cap is made of PP, with light weight polypropylene fabric for comfort and economy, with soft non-irritating elastic band.

Intended use: The caps are intended to be used in medical institutions, laboratories and other places.

Material: Non-woven (polypropylene)

● Approvals: CE, ISO 13485

Manufacturer: Leboo Healthcare Products Limited

Country of origin: China

Specification: Detail see specification as below.



-		Α			-
	$\gamma \gamma$	$\Upsilon$		$\gamma$	

Material W	Maight/gam)	Size(inch)	Measurement Location(inch)					Chulo	
	Weight(gsm)		Α	В	С	D	Color	Style	
PP	10	21	21	1	1	1			
Tolerance			1	/	/	/	Blue	Single Elastic Pack: 100pcs	

Design: Check: Approve: Date:




⊕ H.O.: A-1506, Lar Valley International, No.168 Guang'anmen Wai St., Xicheng District, Beijing 100055, China.
♠ T: +86 10 63252699
➡ E: leboo@leboo.com.cn
➡ W: www.leboo.com.cn

### 

- (1) Carton Requirements: Five layers corrugated
- (2) Packing Quantity: 100pcs/bag, 1000pcs/ctn5 Layers Carton(Paper)





#### **Storage:**

- (1) Store in dry, clean conditions in original packaging.
- (2) Store away from direct sunlight, source of high temperature and solvent vapors.
- (3) Store with the temperature range -5 $^{\circ}$ C to +45 $^{\circ}$ C and with relative humidity below 65%.

#### **Shelf Life:**

Shelf life is 60 months from the date of manufacture when stored as stated above.

## **Cautions and disposal:**

- (1) Do not contact with heavy oils, sparks or flames;
- (2) Do not use in the environment with conditions of excessive heat;
- (3) Replace if damaged, heavily contaminated or in accordance with local practice;
- (4) Hand and dispose of contaminated products with care and in accordance with national regulations.



Keep away from sunlight



Use only once



Keep dry



Don't use when packing damaged and consult instructions for use



Medical device



Caution



Non-sterile



Consult instructions for use or consult electronic instructions for use



CE Symbol

Design: Check: Approve: Date: