


Lotus 41



Upnmed Equipment

SpO2 Ext.Cable



Brand	Compatible Models	P/N	Description	Male Connector	Female Connector
Mindray (mental 5P)	Mindray 509B, 500A, PM9000 (Mindray tech)	U703-1AL	Spo2 extension cable L=2.4M, 5pin-DB9F require Mindray sensor		
	Mindray MEC-509, 509+, 503-507 (CSI tech)	U703-2AL	Spo2 extension cable L=2.4M, 5pin-DB9F require CSI sensor		
		U703-3AL	Spo2 extension cable L=2.4M, 5pin-DB9F require Ohmeda sensor		
Mindray (Lemo 5P)		U703-4AL	Spo2 extension cable L=2.4M, 5pin-DB9F require Nellcor sensor		
Mindray (Lemo 6P 40°)	Mindray: PM7000, 6201, 8000, 9000, VS800 (mindray module), MEC- 1000, MEC-1200, MEC-100, MEC- 200 Medec: MEDEC-1000, M-1000, M- 1100 Kontron: 7138, 7850, 7855 ChaoLiDe, JinKeWei, QiLan, Biocare: BM9000 Creative: PC9000, UP6000 (Nellcor tech)	U703-5AL	Spo2 extension cable L=2.4M, 6pin 40degree- DB9F require Nellcor sensor		
	PM9000, PM9100, PM9200, PM6201, PM7000, PM8000 (before 2006) (Masimo tech)	U703-8AL	Spo2 extension cable L=2.4M, 6pin 40degree- DB9F require Masimo sensor		
		U703-8C	Spo2 extension cable L=2.4M, 6pin 40 degree- 6J require Masimo LNOPSensor		
		U703-11H	Spo2 extension cable L=2.4M, 6pin 40degree- DB9F require Oximax sensor		
Mindray (6P 60°)	(Nellcor tech)	U703-6A	Spo2 extension cable L=2.4M, 6pin 60degree- DB9F require Nellcor sensor		