




Hefei C&P Nonwoven Products Co.,Ltd  
Add: No.22Park road,Feidong new city development area,Hefei,Anhui,China  
TEL:86-551-67707456 Fax:86-551-67709567

E-mail: [info@cponwoven.com.cn](mailto:info@cponwoven.com.cn) [www.medicaldrape.com](http://www.medicaldrape.com)

## TECHNICAL DATA SHEET

Date:2023-06-01

PRODUCT : REINFORCED SURGICAL GOWN STERILE AND DISPOSABLE



PRODUCT DETAILS:

- REINFORCED SURGICAL GOWN  
THEY ARE ULTRASONIC WEIDING

MATERIAL:


- BLUE 45GSM SMS(3 layers SPUNBOND + MELTBLOWN + SPUNBOND NONWOVEN) REINFORCEMENT 40GSM (POLYETHYLENE +POLYPROPYLENE)
- WITH KNITTED CUFFS AND AJUSTABLE VELCRO NECK FASTNERS

CHARACTERISTICS

- PROTECT DOCTORS AND NURSES DURING SURGERY
- RESISTANT TO LIQUID PENETRATION AND MICROBIAL PENETRATION
- ULTRASONIC WELDING
- EFFECTIVE MICROBIAL BARRIER IN WET AND DRY STATE
- RESISTANT TO BURST, TENSION AND TRACTION
- WITHOUT PARTICAL RELEASE
- ANTI STATIC TREATMENT
- STATIC CAPACITY, LOW NENORY CAPACITY
- RETARTED BURNING
- BACK TIES



- **ODORLESS**
- **COOL AND BREATHEABLE**



	Measuring Part/cm	S	M	L	XL	XXL	Tolerance
<b>COD E</b>		GR000-DE	GR010-DE	GR011-DE	GR012-DE	GR013-DE	
<b>Gown</b>	Length(J)	117	125	132	140	150	
	Cirumference(H)	145	150	157	165	170	
	Chest Width(D)	60	65	71	76	84	
	Arm Length(B)	50	55	59	65	65	
	Sleeve Width(C)	27	28	30	31	31	
	Shoulder(A)	55	57	59	66	73	
	Reinforcement Width(M)	50	50	50	50	50	
	Reinforcement Length(N)	70	70	90	90	90	
Sleeve Reinforcement (O)	40	40	40	40	40		



Hefei C&P Nonwoven Products Co.,Ltd  
Add: No.22Park road,Feidong new city development area,Hefei,Anhui,China  
TEL:86-551-67707456 Fax:86-551-67709567

E-mail: [info@cnonwoven.com.cn](mailto:info@cnonwoven.com.cn) [www.medicaldrape.com](http://www.medicaldrape.com)

## **INTENDED USE**

Destined to be used by surgical team members to prevent the transference of infectious agents during surgical procedures.

## **STERILIZATION**

Products are sterilized by ethylene oxide according to standart EN ISO 11135. They are sterilized with EO sterilizer in house. Sterilization is documented and monitored through necessary forms, biological indicators, chemical indicators, dosimeters and data loggers.

## **PACKAGING**

Individual pouch pack, pouch has necessary peeling opening available.

## **LABELING**

Adhesive labels are identifying the product for traceability and registration purposes. And all labels indicate the information below

- Product series name and product name
- Set type(dimensions and type of drape)
- Set contents with quantities
- Manufacturing date
- Date of expiry
- Lot number
- Reference number
- Sterilization indicator
- Method of sterilization
- Warnings (for ex: do not use open packages)
- Manufacturer contact and company information
- CE mark

## **SHELF LIFE**

Product have a shelf life of three(3) years according to the test of Package Verification Test of Disposable sterilization pouch Accelerated Aged for 3years.

## **STORAGE CONDITIONS**

Keep the product under normal conditions for healthcare materials in a cool and dry place protected from sunlight. The use of medical grade paper and film guarantees its sterility for the duration of its shelf life under normal transport and storage conditions.



Hefei C&P Nonwoven Products Co.,Ltd  
Add: No.22Park road,Feidong new city development area,Hefei,Anhui,China  
TEL:86-551-67707456 Fax:86-551-67709567

E-mail: [info@cponwoven.com.cn](mailto:info@cponwoven.com.cn) [www.medicaldrape.com](http://www.medicaldrape.com)

## QUALITY CERTIFICATION AND COMPLIENCE

- All process runs in compliance with EN 13795 AND EN ISO 13485

### Transport

Gowns are transported in a carton boxes. Basic warnings and information is printed on cartons with manufacturer information and contact details. CE marking and name can be seen on carton

### Comments

The product must be used by specialized personel only. Check the package's integrity before use. Do not use if the package is open or damaged.

