





104 ARN

The product has been designed to provide higher stabilization than the sponge and the Nelson collars following a neck trauma or cervical disc hernia. Made of plastizote and suitable for MR monitoring as it does not contain any metal parts.







104B

Functions like the classical Philadelphia collar. It has a tracheal opening in the front. Made of plastizote and suitable for MR monitoring as it does not contain any metal parts.

## **ARMOR CLASSICS**







Available colors:



This product is actively used for facet joint problems of the lumbago, unbalanced postural muscles of the lumbar spine, ligament problems, osteokondrosis, spondiliosis, spondiliartosis and chronic back pains. Air permeable elastic filamented woven textile. 32 cm back height.







Available colors:



ARC 320K

This product is actively used for facet joint problems of the lumbago, unbalanced postural muscles of the lumbar spine, ligament problems, osteokondrosis, spondiliosis, spondiliartosis and chronic back pains. Air permeable elastic filamented woven textile. 32 cm back height. Comes with built-in additional belt



## **ARMOR CLASSICS**



This product is applicable to correct the dorsal vertebra and reduce kyphosis, in traumas of lumbar and throcic vertebra and postural abnormalities. Contains a pair of steel inserts.



This product has been designed to correct dorsal vertebra, reduce kyphosis and for thoracic vertebral traumas and postural deficiencies. Contains a «Y» shaped light weight durable aluminum insert.