

Save \$650 & get a lab-themed hidden object puzzle [Save and claim >](#)



[Home](#) > [Shop All Products](#) > [Buffers And Standards](#) > [Electrochemistry Buffers](#) > [Anode Buffer Container \(ABC\), Fo...](#)



[Certificates](#)  [SDS](#)

Applied Biosystems™

Anode Buffer Container (ABC), for 3500/SeqStudio™ Flex

The Anode Buffer Container (ABC) contains 1X running buffer to support all electrophoresis applications on the Applied Biosystems 3500 [Read more](#)

Have Questions? [Contact Us](#)

Catalog Number	Quantity
4393927	4 pack

alog number 4393927



Save \$650 & get a lab-themed hidden object puzzle [Save and claim >](#)



Quantity:

4 pack

- Product Overview
- Documents
- FAQ

The Anode Buffer Container (ABC) contains 1X running buffer to support all electrophoresis applications on the Applied Biosystems 3500 and SeqStudio Flex series genetic analyzers. It is made in a ready-to-use, disposable container with a radio frequency identification (RFID) tag incorporated into the label. The top of the ABC is heat sealed with a plastic film, which should be removed prior to direct installation on the instrument. Each package includes four individual containers.

For Research Use Only. Not for use in diagnostic procedures.

Specifications	
Recommended Storage	Contains 4 packs of anode buffer container. Store the Anode Buffer Container at 2-8°C until ready to use. Refer to the expiration date on the label.
For Use With (Equipment)	3500 Series Genetic Analyzers, SeqStudio Flex Series Genetic Analyzers
Quantity	4 pack
Unit Size	Each

Documents & Downloads

Certificates



Search by lot number or partial lot number

Search

Lot #	Certificate Type	Date	Catalog Number(s)
2503021	Certificate of Analysis	May 19, 2025	4393927

Save \$650 & get a lab-themed hidden object puzzle [Save and claim >](#)



 2503017	Certificate of Analysis	May 08, 2025	4393927
 2503019	Certificate of Analysis	May 05, 2025	4393927

5 results displayed, search above for a specific certificate

[Request a Certificate](#)

Safety Data Sheets



Scientific Resources

Brochures

^

-  [CE Running Buffers - Green Fact Sheet](#)
-  [White paper: Shipping Capillary electrophoresis running at ambient temperature](#)

Technical Notes

^

-  [Technical Note: Ambion Shipping Conditions Anode and Cathode Buffer Container](#)

Product Information

Manuals

^

-  [Product Sheet: Anode Buffer Container 3500/Flex Series - SeqStudio Flex and 3500 series instruments](#)

Save \$650 & get a lab-themed hidden object puzzle [Save and claim >](#)



License #481 - Sequencing or Fragment Analysis Intellectual Property

▼

Frequently asked questions (FAQs)

When I started my run on the Applied Biosystems 3500/3500xL Genetic Analyzer, an error occurred: "Unstable Electrophoresis current detected, check for air bubbles" but I don't see any bubbles. What else can cause this?

▼

I have an Applied Biosystems 3500/3500xL Genetic Analyzer. When do my consumables expire?

▼

My Applied Biosystems 3500/3500xL Genetic Analyzer consumables have expired, but how long past expiration can I really use them?

▼

I am only using the Applied Biosystems 3500/3500xL Genetic Analyzer once a week. Should I leave the polymer on the instrument or put it back in the refrigerator?

▼

Do I need to replace the buffer every 14 days if the Applied Biosystems 3500/3500xL Genetic Analyzer is not in use?

▼

[View more](#)