



# **vdx** 4343

#### DR acquisition system with Flat Panel technology

#### Product Data

The VDX 4343TS imaging system is an image acquisition system that can be integrated with general-purpose radiology systems. Designed to improve the general X-ray diagnostic room workflow, it provides high quality images and long-term reliability. The acquisition system uses Canon FDX4343R flat panel detector which provides outstanding image sharpness and detail. The system includes:


- Canon FDX4343R amorphous Silicon Cesium lodide scintillator detector to convert X-ray photons into a digital image (one or two detectors according to the system configuration).
- Integrated Operator Console providing full control of exposure (only with G100C generators), image acquisition, processing and DICOM functions.

## Digital Detector – FDX4343R

Tuno	Flat Danal Datastor
Type	Flat Panel Detector
X-ray conversion layer	Cesium lodide (Csl) with amorphous Silicon (a-Si) photodiode
Nominal Active Image Area	43 x 43 cm (16.9" x 16.9")
Image Resolution (lp/mm)	3.5 lp/mm (Nyquist frequency)
lmage matrix size	3072 x 3072 pixels
lmage depth	14-bit acquired (16384 gray levels)
Output characteristic	Linear response up to the saturation dose
X-ray dose saturation	60 µGy / 6800 µR
Detector pixel size	140 µm
TFT matrix readout time	1.22 s
Acquisition window	Up to 3200 ms in 3 steps (500 ms; 1000 ms; 3200 ms)
MTF (Modulation Transfer Function) RQA5 according to IEC 62220-1-1	1 0.9 0.8 0.7 0.6 E 0.5 0.4 0.3 0.2 0.1 0 0.5 1 1.5 2 2.5 3 3.5 Frequency (Ip/ mm)







### Electrical features

Power supply (synchronizer)	24 Vdc
Absorbed current	
(synchronizer)	0.5 A
Power supply (detector)	100-240 Vac single phase, 50/60 Hz
Absorbed power (detector)	60 VA

## Standards and regulations



CE symbol grants the product compliance to the European Directive for Medical Devices 93/42/FFC and its revised versions

Note: Products are continuously under review in the light of technical advancement. The actual specification may therefore be subject to improvement or modification without notice.

VILLA SISTEMI MEDICALI s.p.a. 20090 BUCCINASCO (MI) – ITALY, Via delle Azalee, 3 Tel. +39-02-488591, Fax +39-02-4881844 Company with Quality System certified by







ISO 9001:2015

ISO13485:2016