

Procedure Name			TOMOS Head, TOMOS Hand, TOMOS Knee, TOMOS Abdomen	TOMOS C-Spine TOMOS L-Spine, TOMOS Hip Joint, TOMOS Chest	TOMO
General	Image Processing	Table Action Mode	Normal	Normal	Normal
		TOUCH Disp. Mode	RF	RF	RF
		Auto Shutter	MED	MED	MED
		VReverse	OFF	OFF	OFF
		HReverse	OFF	OFF	OFF
		Field of View	12 inch	17 inch	HOLD
		Marker	OFF	OFF	OFF
		Position Memory	no select	no select	no select
	Correction Step	RkV(kV)	2	2	2
		Acq.Time (step)	3	3	3
	Destination	Send Dest 1	no select	no select	no select
		Send Dest 2	no select	no select	no select
		Auto Send Timing	Next ACQ	Next ACQ	Next ACQ
		Print Preset	no select	no select	no select