

Sterling[™]

0.018" balloon dilatation catheter

Demand More

The #1 0.018" (0,46mm) balloon dilatation catheter with the most comprehensive size offering available


Sterling™ Balloon Dilatation Catheter

**Boston
Scientific**
Advancing science for life™

Demand More from your 0.018" (0,46mm) balloon catheter

Sterling delivers:

- The most comprehensive 0.018" (0,46mm) size offering
- *Exceptionally low 0.020" (0,51mm) profile on a 0.018" (0,46mm) wire
- Indicated for wider acute use including carotid** and post-stent dilatation




*Tip entry measurements performed at the distal-most portion of the tip. Average measurements taken by Boston Scientific (n=3, 5 mm x 150 mm balloons). Data on file. (TR 90884555). Bench testing may not be representative of clinical performance. All photographs are property of Boston Scientific.

**Monorail™ only

Tip




- Exceptional crossability
- Best in class* 0.020" (0,51mm) tip entry profile on a 0.018" (0,46mm) wire designed to facilitate crossability through tight stenoses
- Seamless, formed tip designed to enhance crossability



*Tip entry measurements performed at the distal-most portion of the tip. Average measurements taken by Boston Scientific (n=3, 5 mm x 150 mm balloons). Data on file. (TR 90884555). Bench testing may not be representative of clinical performance. All photographs are property of Boston Scientific.

Balloon

- Pebax™ semi-compliant balloon material
- Best in class* rapid deflation
- 8x40 mm through 4 Fr** and 7x220 mm through 5 Fr**
- Bioslide™ hydrophilic coating designed to facilitate advancement through complex lesions
- Highly-visible gold RO markers define the edges of the balloon




*Average measurements taken by Boston Scientific (n=3, 5 mm x 150 mm balloons). Data on file. (TR 90884555). Bench testing may not be representative of clinical performance. All photographs are property of Boston Scientific. Illustrations are not necessarily to scale. **4Fr is 1,33mm; 5Fr is 1,67mm

Shaft




- 0.018" (0,46mm) and 0.014" (0,36mm) guidewire compatible
- OTW* and Monorail™ rapid-exchange catheter platforms
- 40, 80, 90, 135, and 150cm shafts
- Shaft technology engineered to promote fast deflation
- Balance of trackability and pushability




*OTW is Over The Wire



Sterling™ OTW*



Sterling Monorail™



*OTW is Over The Wire

Demand More Size Options

		length (mm)											
Balloon		10	15	20	30	40	60	80	100	120	150	200	220
diameter (mm)	2.0 mm												
	2.5 mm												
	3.0 mm												
	3.5 mm												
	4.0 mm												
	4.5 mm												
	5.0 mm												
	5.5 mm												
	6.0 mm												
	6.5 mm												
	7.0 mm												
	8.0 mm												
	9.0 mm												
10.0 mm													

- Available in lengths up to 220 mm to treat long SFA and BTK disease with fewer inflations
- 112 unique balloon sizes, 316 codes
- OTW* and Monorail rapid-exchange platforms on 40, 80, 90, 135, and 150 cm shafts

*OTW is Over The Wire

All cited trademarks are the property of their respective owners. Average measurements taken by Boston Scientific. Bench testing performed by Boston Scientific Corporation.

Sterling™ family of balloon dilatation catheters is STERLING™ MONORAIL™ AND STERLING™ OVER-THE-WIRE PTA BALLOON DILATATION CATHETERS and STERLING™ SL MONORAIL™ AND STERLING™ SL OVER-THE-WIRE PTA BALLOON DILATATION CATHETERS.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations.