BRUCELLOSIS



Rose Bengal (Rapid slide agglutination antigen)



Rapid slide agglutination antigen for the detection of anti-*Brucella abortus, melitensis*, or *suis* antibodies in bovine, ovine, caprine, and porcine species.



www.id-vet.com

With you at every step

BRUCELLOSIS

Rose Bengal (Rapid slide agglutination antigen)

Rapid slide agglutination antigen for the detection of anti-*Brucella abortus, melitensis,* or *suis* antibodies in bovine, ovine, caprine, and porcine species. • Suspension of *Brucella abortus* biovar 1 Weybridge strain No 99, inactivated by heat and phenol and coloured with rose bengal stain in an acidified buffer.

• Calibrated on the OIEISS International *Brucella* standard serum and the national secondary serum to be positive against a 1/45 dilution and negative against a 1/55 dilution, as per the European directive CEE 64/432 and OIE requirements.

Ordering information

Product code	RSA-RB
Kit format	10 ml
Reactions	330 reactions (30 µl per test)

References

(1) OIE Terrestrial Manual 2009; 2.4.3; Bovine brucellosis.

(2) Annex C of the European directive CEE 64/432

Associated Products

Brucellosis positive freeze-dried serum

Product Code	MRI-BRU
Format	1 ml vial (freeze-dried)
Description	Freeze-dried highly positive bovine serum to be used as internal reference material for quality control. Positive when tested with Rose Bengal.

Brucellosis negative freeze-dried serum

Product Code	MRINEG-BRU
Format	1 ml vial (freeze-dried)
Description	Freeze-dried negative bovine serum to be used as internal reference material for quality control. Negative when tested with Rose Bengal.

