

Optional Accessories

for AUTOCON® II 400 SCB and AUTOCON® II 80

			for use with	
			AUTOCON® II 400 SCB	AUTOCON® II 80
	20 0178 31	Three-Pedal Footswitch , for use with AUTOCON® II 400 SCB	●	—
	20 0178 30	Two-Pedal Footswitch , for use with AUTOCON® II 400 SCB and AUTOCON® II 200	●	—
	20 0178 32	One-Pedal Footswitch , for activating coagulation, for use with AUTOCON® II 400 SCB and AUTOCON® II 200	●	—
	20 0178 34	Two-Pedal Footswitch , digital, one-stage, for use with AUTOCON® II 80	—	●
	20 0178 33	One-Pedal Footswitch , digital, one-stage, for use with AUTOCON® II 80	—	●
	27805	Neutral Electrode , of conductive silicone, with 2 rubber ties for fastening, contact surface A = 500 cm², for use with Connecting Cable 27806	●	●
	860021 E	Neutral Electrode , of conductive silicone, with 1 rubber tie for fastening, contact surface A = 187 cm², for use with Connecting Cable 27806	—	●
	27806	Neutral Electrode Connecting Cable , for Neutral Electrodes 27805 and 860021 E, length 400 cm	-111 -115 -122 -125	●
	27806 UR	Neutral Electrode Connecting Cable , for use with Neutral Electrode 27805	-112 -116 -122 -125	●
	27806 US	Neutral Electrode Connecting Cable , for use with Neutral Electrode 27802	-112 -116 -122 -125	●
	27802	Neutral Electrode , for single use, contact surface divided into two, A = 169 cm², package of 50, Connecting Cable 27806 US required	●	●
	27801	Connecting Cable , for connecting Neutral Electrode 27802, length 500 cm	-111 -115 -122 -125	●
	26 5200 43	Electrode Handle , with 2 buttons for activating the unipolar generator, yellow button: unipolar cutting, blue button: unipolar coagulation, Connecting Cable 26 5200 45	●	●
	26 5200 45	High Frequency Cable , for Electrode Handle 26 5200 43, length 400 cm	-111 -115 -122 -125	●
	26 5200 46	Electrode Handle , without buttons, with integrated connecting cable, length 300 cm	-111 -115 -122 -125	●