
Medtronic

Export Advance

▲ Indications, Safety, and Warnings



ADVANCING PERFORMANCE

The Export Advance™ aspiration catheter delivers consistent, high-performing aspiration power when it matters most — restoring flow and protecting patients.

PRODUCT DETAILS†

Superior Deliverability¹

Full-Wall Technology provides variable levels of stiffness without joints, for optimal kink resistance and catheter performance.

Dependable Delivery to the Target Site:

- Preloaded stylet enhances shaft stiffness during catheter delivery, boosting kink resistance, trackability, and pushability²
- Buddy wire compatibility allows for extra support, when needed.†

High-performing Aspiration Power:

- Large extraction lumen (0.044" proximal; 0.043" distal) increases aspiration power.‡
- Optimised hub geometry improves flow.‡
- Soft, short, forward-facing tip design permits excellent particle capture.[†]

ORDERING INFORMATION

Product Code	Guide Compatibility	Wire Compatibility (in)	Length (cm)
ADVANCE	6 F min. guide I.D. 0.070"	0.014	140

- ^{*} ™Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
- † Based on bench test data, 2013. Bench test data may not be indicative of clinical performance.
- [‡] Based on bench test data vs. Export™ AP catheter.
- ¹ Bench testing vs. Terumo Eliminate™* and Eliminate™* 3 aspiration catheters, Vascular Solutions Pronto®* LP, Pronto®* V3 and Pronto®* V4 aspiration catheters, Merit®* ASAP aspiration catheter, Kaneka Thrombuster III GR* aspiration catheter, Spectranetics QuickCat™* aspiration catheter, IHT Polska Hunter®* aspiration catheter, Medrad Possis™* Fetch™* 2 aspiration catheter, and Hexacath Recover* aspiration catheter. 10093906DOC, 2013. Bench test data may not be indicative of clinical performance.
- ² Based on animal studies compared with Export[™] AP aspiration catheter without stylet.

REQUEST MORE INFORMATION

CONTACT US

RELATED PRODUCTS

- Catheters
- Cardiovascular Products

Privacy Statement | Terms of Use | Contact