

Technical Data

Cystine Heart Agar Base

M172

Cystine Heart Agar when enriched with haemoglobin is recommended for the cultivation of *Francisella tularensis*. Without enrichment it supports excellent growth of gram-negative cocci and other pathogenic organisms.

Composition**

Ingredients	Gms / Litre
Beef heart, infusion from	500.000
Proteose peptone	10.000
Dextrose	10.000
Sodium chloride	5.000
L-Cystine	1.000
Agar	15.000
Final pH (at 25°C)	6.8 ± 0.2

^{**}Formula adjusted, standardized to suit performance parameters

Directions

Suspend 51 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. When to be enriched with haemoglobin (2%), suspend 10.2 grams of medium in 100 ml distilled water. Sterilize as above. Cool medium to 50°C and aseptically add 100 ml of 2% sterile haemoglobin solution. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Francisella tularensis is the cause of tularaemia, a plague-like disease of rodents and other small organisms. It was first described in humans in 1907 (1). The organisms are strict aerobes; fresh isolates cannot be cultured on ordinary medium but require a complex medium containing blood, or tissue extracts and cystine. Several media formulations were employed to isolate this microorganism. Blood Dextrose Cystine Agar, described by Francis (2) was found to be satisfactory for cultivating F. tularensis . Addition of 0.05% cystine and 1% dextrose to Heart Infusion Agar can also be employed for cultivation of F. tularensis !(3). Subsequently haemoglobin was added to Cystine Heart Agar Base to develop a satisfactory cultivation medium for F. tularensis (4). This medium is also known as Cystine Glucose Blood Agar and is the most suitable medium for isolating F. tularensis (2). Hemoglobin provides additional nutrients and growth factors. This medium also supports growth of gram-negative cocci and other pathogenic microorganisms without additional enrichment. Cystine Heart Agar Base can be supplemented with Rabbit blood and antimicrobial agents (5).

This medium is a nutritionally rich medium, which may also be used for cultivating many other organisms generally difficult to grow.

Beef heart infusion and proteose peptone are sources of carbon, nitrogen, vitamins and minerals. Dextrose is an energy source. L-Cystine is the source of amino acid. Sodium chloride provides the essential ions. Overgrowth by contaminating organisms can be reduced by incorporating 100-500 units penicillin per ml into the medium (1).

F. tularensis is a Biosafety Level 2 pathogen that can be transmitted by aerosols or by penetration of unbroken skin (5). Wearing of gowns, gloves and masks is recommended for people handling suspected infectious material (6).

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

HiMedia Laboratories Technical Data

Basal medium :Amber coloured clear to slightly opalescent gel After addition of 2% haemoglobin solution: Chocolate brown coloured opaque gel forms in Petri plates

Reaction

Reaction of 5.1% w/v aqueous solution at 25°C. pH: 6.8±0.2

pН

6.60-7.00

Cultural Response

M172: Cultural characteristics observed with added 2% Haemoglobin after an incubation at 35-37°C for 40-48 hours.

Organism Growth

Cultural Response

Francisella tularensis ATCC luxuriant

29684

Neisseria meningitidis ATCC luxuriant

13090

Streptococcus pneumoniae luxuriant

ATCC 6303

Streptococcus pyogenes luxuriant

ATCC 19615

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on label.

Reference

- 1. Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Tenover F. C., Yolken R. H., (Ed.), 1999, Manual of Clinical Microbiology, 7th Ed. American Society for Microbiology, Washington, D.C.
- 2. Francis, 1928, JAMA, 91:1155.
- 3. Shaw, 1930, Zentr. Bakt. I. Abt. Orig., 118:216.
- 4. Rhamy, 1933, Am. J. Clin. Pathol., 3:121.
- 5. Isenberg, (Ed.), 1992, Clinical Microbiology Procedures Handbook, Vol. 1. American Society for Microbiology, Washington, D.C.
- 6. U.S. Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health, 1999, Biosafety in Microbiological and Biomedical Laboratories, 4th Ed., HHS Publication.

Revision: 1 / 2011

CE

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.