

# AXIUM™ DETACHABLE COILS — HIGH PERFORMANCE MEETS EASE OF USE

## Axium™ 3D Detachable Coils

SKU	SECONDARY DIAMETER (MM)	COIL LENGTH (CM)	PRIMARY COIL OD	VOLUME (MM <sup>3</sup> )
QC-2-2-3D	2	2	0.0115 in (0.29 mm)	1.34
QC-2-4-3D		4		2.68
QC-2-6-3D		6		4.02
QC-2.5-2-3D	2.5	2		1.34
QC-2.5-4-3D		4		2.68
QC-2.5-6-3D		6		4.02
QC-2.5-8-3D	2.5	8		5.36
QC-3-4-3D		4		2.68
QC-3-6-3D		6		4.02
QC-3-8-3D	3	8		5.36
QC-3-10-3D		10		6.70
QC-3.5-6-3D		6		4.02
QC-3.5-12-3D	3.5	12	8.04	
QC-3.5-15-3D		15	10.05	
QC-4-6-3D		6	4.75	
QC-4-8-3D	4	8	6.33	
QC-4-10-3D		10	7.92	
QC-4-12-3D		12	9.50	
QC-5-8-3D	5	8	6.33	
QC-5-10-3D		10	7.92	
QC-5-15-3D		15	11.88	
QC-6-10-3D	6	10	7.92	
QC-6-15-3D		15	11.88	
QC-6-20-3D		20	15.83	
QC-7-15-3D	7	15	13.85	
QC-7-20-3D		20	18.47	
QC-7-30-3D		30	27.70	
QC-8-15-3D	8	15	13.85	
QC-8-20-3D		20	18.47	
QC-8-30-3D		30	27.70	
QC-9-20-3D	9	20	18.47	
QC-9-30-3D		30	27.70	
QC-10-20-3D		20	18.47	
QC-10-30-3D	10	30	27.70	
QC-12-30-3D		30	31.96	
QC-12-40-3D		40	42.61	
QC-14-30-3D	14	30	31.96	
QC-14-40-3D		40	42.61	
QC-16-40-3D		40	42.61	
QC-18-40-3D	18	40	42.61	
QC-20-50-3D		50	53.27	
QC-22-50-3D		50	53.27	
QC-25-50-3D	25	50	53.27	

## Axium™ Helix Detachable Coils

SKU	SECONDARY DIAMETER (MM)	COIL LENGTH (CM)	PRIMARY COIL OD	VOLUME (MM <sup>3</sup> )
QC-1.5-1-HELIX	1.5	1	0.0115 in (0.29 mm)	0.67
QC-1.5-2-HELIX		2		1.34
QC-1.5-3-HELIX		3		2.01
QC-1.5-4-HELIX	1.5	4		2.68
QC-2-1-HELIX		1		0.67
QC-2-2-HELIX		2		1.34
QC-2-3-HELIX	2	3		2.01
QC-2-4-HELIX		4		2.68
QC-2-6-HELIX		6		4.02
QC-2-8-HELIX	2	8		5.36
QC-2.5-2-HELIX		2		1.34
QC-2.5-4-HELIX		4		2.68
QC-2.5-6-HELIX	2.5	6	4.02	
QC-2.5-8-HELIX		8	5.36	
QC-3-4-HELIX		4	2.68	
QC-3-6-HELIX	3	6	4.02	
QC-3-8-HELIX		8	5.36	
QC-4-8-HELIX		8	6.33	
QC-4-10-HELIX	4	10	7.92	
QC-4-12-HELIX		12	9.50	
QC-5-15-HELIX		15	11.88	
QC-5-20-HELIX	5	20	15.83	
QC-6-20-HELIX		6	15.83	
QC-7-20-HELIX		7	18.47	
QC-7-30-HELIX	7	30	27.70	
QC-8-20-HELIX		8	18.47	
QC-8-30-HELIX		8	27.70	
QC-9-20-HELIX	9	20	18.47	
QC-9-30-HELIX		30	27.70	
QC-10-20-HELIX		10	18.47	
QC-10-30-HELIX	10	30	27.70	
QC-12-30-HELIX		12	31.96	
QC-12-40-HELIX		12	42.61	
QC-14-30-HELIX	14	30	31.96	
QC-14-40-HELIX		40	42.61	
QC-16-30-HELIX		16	31.96	
QC-16-40-HELIX	16	40	42.61	
QC-18-40-HELIX		18	42.61	
QC-20-40-HELIX		20	42.61	
QC-20-50-HELIX	20	50	53.27	

## NOW AVAILABLE

### Axium™ Prime 3D Detachable Coils (Extra Soft)

SKU	SECONDARY DIAMETER (MM)	COIL LENGTH (CM)	PRIMARY COIL OD	VOLUME (MM <sup>3</sup> )
APB-1-2-3D-ES	1	2	0.0108 in (0.27 mm)	1.18
APB-1-3-3D-ES		3		1.77
APB-1-4-3D-ES		4		2.36
APB-1.5-2-3D-ES	1.5	2		1.18
APB-1.5-3-3D-ES		3		1.77
APB-1.5-4-3D-ES		4		2.36
APB-2-2-3D-ES	2	2		1.18
APB-2-3-3D-ES		3		1.77
APB-2-4-3D-ES		4		2.36
APB-2.5-4-3D-ES	2.5	4		2.36
APB-2.5-6-3D-ES		6		3.55
APB-3-4-3D-ES		4		2.36
APB-3-6-3D-ES	3	6	3.55	
APB-3-8-3D-ES		8	4.73	
APB-3.5-6-3D-ES		6	3.55	
APB-3.5-8-3D-ES	3	8	4.73	
APB-3.5-10-3D-ES		10	5.91	

### Axium™ Prime Helix Detachable Coils (Extra Soft)

SKU	SECONDARY DIAMETER (MM)	COIL LENGTH (CM)	PRIMARY COIL OD	VOLUME (MM <sup>3</sup> )
APB-1-1-HX-ES	1	1	0.0108 in (0.27 mm)	0.59
APB-1-2-HX-ES		2		1.18
APB-1-3-HX-ES		3		1.77
APB-1.5-1-HX-ES	1.5	1		0.59
APB-1.5-2-HX-ES		2		1.18
APB-1.5-3-HX-ES		3		1.77
APB-1.5-4-HX-ES	1.5	4		2.36
APB-2-1-HX-ES		1		0.59
APB-2-2-HX-ES		2		1.18
APB-2-3-HX-ES	2	3		1.77
APB-2-4-HX-ES		4		2.36
APB-2-6-HX-ES		6		3.55
APB-2-8-HX-ES	2	8	4.73	
APB-2.5-3-HX-ES		3	1.77	
APB-2.5-4-HX-ES		4	2.36	
APB-2.5-6-HX-ES	2.5	6	3.55	
APB-2.5-8-HX-ES		8	4.73	
APB-3-4-HX-ES		4	2.36	
APB-3-6-HX-ES	3	6	3.55	
APB-3-8-HX-ES		8	4.73	
APB-3-10-HX-ES		10	5.91	

## Axium™ 3D PGLA MicroFX Detachable Coils

SKU	SECONDARY DIAMETER (MM)	COIL LENGTH (CM)	PRIMARY COIL OD	VOLUME (MM <sup>3</sup> )
PC-2-2-3D	2	2	0.0115 in (0.29 mm)	1.34
PC-2-4-3D		4		2.68
PC-2-6-3D		6		4.02
PC-3-4-3D	3	4		2.68
PC-3-6-3D		6		4.02
PC-3-8-3D		8		5.36
PC-4-6-3D	4	6		4.75
PC-4-8-3D		8		6.33
PC-4-10-3D		10		7.92
PC-4-12-3D	4	12		9.50
PC-5-8-3D		8		6.33
PC-5-10-3D		10		7.92
PC-5-15-3D	5	15	11.88	
PC-6-10-3D		10	7.92	
PC-6-15-3D		15	11.88	
PC-6-20-3D	6	20	15.83	
PC-7-15-3D		15	13.85	
PC-7-20-3D		20	18.47	
PC-7-30-3D	7	30	27.70	
PC-8-15-3D		15	13.85	
PC-8-20-3D		20	18.47	
PC-8-30-3D	8	30	27.70	
PC-9-20-3D		20	18.47	
PC-9-30-3D		30	27.70	
PC-10-20-3D	10	20	18.47	
PC-10-30-3D		30	27.70	
PC-12-30-3D		12	31.96	
PC-12-40-3D	12	40	42.61	
PC-14-30-3D		14	31.96	
PC-14-40-3D		40	42.61	
PC-16-40-3D	16	40	42.61	
PC-18-40-3D		18	42.61	

## Axium™ Helix Nylon MicroFX Detachable Coils

SKU	SECONDARY DIAMETER (MM)	COIL LENGTH (CM)	PRIMARY COIL OD	VOLUME (MM <sup>3</sup> )
NC-2-1-HELIX	2	1	0.0115 in (0.29 mm)	0.67
NC-2-2-HELIX		2		1.34
NC-2-3-HELIX		3		2.01
NC-2-4-HELIX	2	4		2.68
NC-2-6-HELIX		6		4.02
NC-2-8-HELIX		8		5.36
NC-3-4-HELIX	3	4		2.68
NC-3-6-HELIX		6		4.02
NC-3-8-HELIX		8		5.36
NC-4-8-HELIX	4	8		6.33
NC-4-10-HELIX		10		7.92

## Axium™ Helix PGLA MicroFX Detachable Coils

SKU	SECONDARY DIAMETER (MM)	COIL LENGTH (CM)	PRIMARY COIL OD	VOLUME (MM <sup>3</sup> )
PC-2-1-HELIX	2	1	0.0115 in (0.29 mm)	0.59
PC-2-2-HELIX		2		1.18
PC-2-3-HELIX		3		
PC-2-4-HELIX	2	4		
PC-2-6-HELIX		6		
PC-2-8-HELIX		8		1.77
PC-3-4-HELIX	3	4		0.59
PC-3-6-HELIX		6		1.18
PC-3-8-HELIX		8		2.36
PC-4-8-HELIX	4	8		0.59
PC-4-10-HELIX		10		3.55
PC-4-12-HELIX		12		4.73
PC-5-15-HELIX	5	15	1.77	
PC-5-20-HELIX		20	4.73	
PC-6-20-HELIX		6	2.36	
PC-7-20-HELIX	7	20	6.33	
PC-7-30-HELIX		30	7.92	
PC-8-20-HELIX		8	6.33	
PC-8-30-HELIX	8	30	7.92	
PC-9-20-HELIX		20	6.33	
PC-9-30-HELIX		30	7.92	
PC-10-20-HELIX	10	20	6.33	
PC-10-30-HELIX		30	7.92	

## I.D. Instant Detacher

One detacher required per procedure

Product Catalogue Number	Number By Box
ID-1-5	5

## Axium™ Prime 3D Super Soft Detachable Coils

SKU	SECONDARY DIAMETER (MM)	COIL LENGTH (CM)	PRIMARY COIL OD	VOLUME (MM <sup>3</sup> )
APB-4-6-3D-SS	4	6	0.0115 in (0.29 mm)	4.02
APB-4-8-3D-SS		8		5.36
APB-4-10-3D-SS		10		6.70
APB-4-12-3D-SS		12		8.04
APB-5-8-3D-SS	5	8		5.36
APB-5-10-3D-SS		10		6.70
APB-5-15-3D-SS		15		10.05
APB-6-10-3D-SS	6	10		6.70
APB-6-15-3D-SS		15		10.05
APB-6-20-3D-SS		20		13.40

## Axium™ Prime Helix Super Soft Detachable Coils

SKU	SECONDARY DIAMETER (MM)	COIL LENGTH (CM)	PRIMARY COIL OD	VOLUME (MM <sup>3</sup> )
APB-4-6-HX-SS	4	6	0.0115 in (0.29 mm)	4.02
APB-4-8-HX-SS		8		5.36
APB-4-10-HX-SS		10		6.70
APB-4-12-HX-SS		12		8.04
APB-5-10-HX-SS	5	10		6.70
APB-5-15-HX-SS		15		10.05
APB-5-20-HX-SS		20		13.40
APB-6-12-HX-SS	6	12		8.04
APB-6-20-HX-SS		20		13.40

Indications, contraindications, warnings, and instructions for use can be found in the product labeling supplied with each device.

Axium™ and Axium™ Prime detachable coils are intended for the endovascular embolization of intracranial aneurysms. Axium™ and Axium™ Prime detachable coils are also intended for the embolization of other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae.