

To whom it may concern,

By this letter, we **Nipro Medical Europe** on behalf of **Nipro Corporation (Japan)**, would like to provide clarification on the scope of use of NIPRO CF-609N filter:

- The CF-609N filter is compatible and used with SURDIAL X hemodialyses machine, and is used as the component together with the machine;
- The CF-609N filter is recommended to be changed within period of each of 1200 hours of treatment procedure until the machine will warn that filter should be replaced;

Please feel free to contact us with any questions.

Serge Kemps

CEO

