




/ You can find more information about this product in our web-shop.



VACUETTE® SERUM TUBES

VACUETTE® serum tubes are used for determinations in serum for routine clinical chemistry tests and hormones, serology, immunohaematology and TDM.


Serum tubes from Greiner Bio-One are available in three versions: CAT Serum Clot Ac-

tivator, CAT Serum Sep Clot Activator and CAT Serum Fast Separator.

The inner wall of each serum tube is coated with microscopic silica particles (clot activator main component SiO_2) which activate clotting when the tubes are carefully inverted.

VACUETTE® Blood Collection Tubes

VACUETTE® Serum Tubes



IVD

STERILE

CAT Serum Clot Activator

VACUETTE® CAT serum clot activator tubes are coated with micronized silica particles that activate clotting when the tubes are gently turned. Selected tubes are also available with a transparent label.

Item no.	Nominal volume	Cap colour	Ring colour	Thread type	Tube size	Label	Barcode	Qty. inner / outer
454098	1 ml	● red	○ white	PREMIUM	13 x 75	paper	no	50 / 1,200
454096	2 ml	● red	○ white	PREMIUM	13 x 75	paper	no	50 / 1,200
484515	2 ml	● red	○ white	PREMIUM	13 x 75	paper	yes	50 / 1,200
454092	4 ml	● red	● black	PREMIUM	13 x 75	paper	no	50 / 1,200
474092	4 ml	● red	● black	PREMIUM	13 x 75	transparent	no	50 / 1,200
484525	4 ml	● red	● black	PREMIUM	13 x 75	paper	yes	50 / 1,200
456092	6 ml	● red	● black	PREMIUM	13 x 100	paper	no	50 / 1,200
476092	6 ml	● red	● black	PREMIUM	13 x 100	transparent	no	50 / 1,200
486504	6 ml	● red	● black	PREMIUM	13 x 100	paper	yes	50 / 1,200
454236	2 ml	● red	○ white	non-ridged	13 x 75	paper	no	50 / 1,200
454204	4 ml	● red	● black	non-ridged	13 x 75	paper	no	50 / 1,200
456089	6 ml	● red	● black	non-ridged	13 x 100	paper	no	50 / 1,200
456048	6 ml	● red	● black	non-ridged	13 x 100	transparent	no	50 / 1,200
455092	9 ml	● red	● black	non-ridged	16 x 100	paper	no	50 / 1,200