

OXIVIT MED

c/f: 1007600044280; adresa: str. Decebal 82-90, or. Chişinău, Republica Moldova
telefon: + 373 22 808002; fax: + 373 22 808003
web: www.oxivit-med.com; e-mail: info@oxivit-med.com

Date generale despre ofertant

S.C. OXIVIT-MED S.R.L
Administrator: Dmitrii Kojevnikov
Adresa poştală: mun. Chişinău, str. Decebal 82-90
Tel./Fax: 022 808 002, 022 808 003
E-mail: info@oxivit-med.com; oxivit.medical@gmail.com
Cod IBAN: MD09MO2224ASV23488147100
Banca: „Mobiasbanca OTP Group” S.A
Codul băncii: MOBBMD22
Cod fiscal: 1007600044280
Cod TVA: 0306300

Cu respect,

Dmitrii Kojevnikov

Administrator



Benannt durch/Designated by
Zentralstelle der Länder
für Gesundheitsschutz
bei Arzneimitteln und
Medizinprodukten
www.zlg.de
ZLG-BS-244.10.08



Product Service

EC Certificate

Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

(Devices in Class IIa, IIb or III)

No. G1 039709 1263 Rev. 00

Manufacturer:

Medtronic, Inc.

710 Medtronic Parkway
Minneapolis MN 55432
USA

Product Category(ies):

- **Autotransfusion Systems and Associated Disposables**
- **Centrifugal Blood Pumps**
- **Bio-Console Drive Units**
- **Flow Monitoring Systems**
- **Bio-Cal Blood Temperature Controller**
- **Temperature Monitoring Systems and Associated Disposables**
- **Blood Monitoring Systems**
- **Cardioplegia Delivery Systems**
- **Disposable Blood Handling Devices used for Open Heart Surgery**
- **Arterial Filters**
- **Oxygenators including Heat Exchangers, with and without Cardiotomy Reservoirs**
- **Cardiotomy Venous Reservoirs**
- **Venous Reservoir Bags**
- **Perfusion Equipment and Disposable Perfusion Devices**
- **Disposable Medical Devices for Drainage Systems**
- **Disposable Medical Devices for use in Cardiopulmonary Surgery: Cardioplegia, Cannulae, Venting, Suction**
- **Pressure Display System & related accessories of class IIa**
- **Tissue Positioning/Stabilizing Devices**
- **Surgical Site Clearing Devices**
- **Intravascular Shunts**
- **Surgical Retractors**

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final

Page 1 of 3

TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123

TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany



Benannt durch/Designated by
Zentralstelle der Länder
für Gesundheitsschutz
bei Arzneimitteln und
Medizinprodukten
www.zlg.de
ZLG-BS-244.10.08



Product Service

EC Certificate

Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

(Devices in Class IIa, IIb or III)

No. G1 039709 1263 Rev. 00

inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

Report No.: 72150396

Valid from: 2020-02-12

Valid until: 2024-05-26

Date, 2020-02-12

Christoph Dicks
Head of Certification/Notified Body



Benannt durch/Designated by
Zentralstelle der Länder
für Gesundheitsschutz
bei Arzneimitteln und
Medizinprodukten
www.zlg.de
ZLG-BS-244.10.08



Product Service

EC Certificate

Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

(Devices in Class IIa, IIb or III)

No. G1 039709 1263 Rev. 00

Facility(ies):

Medtronic Mexico S.de R.L.de CV
Av. Paseo Cucapah, 10510 El Lago, C.P. 22210 Tijuana, Baja
California, MEXICO

Medtronic Perfusion Systems
7611 Northland Drive, Minneapolis, MN 55428, USA

Medtronic, Inc.
710 Medtronic Parkway, Minneapolis MN 55432, USA

Medtronic

Engineering the extraordinary

Cannula catalog


Medtronic in cardiac surgery

A company
like no other



The largest
cannula portfolio
worldwide, backed
by an even larger
commitment
– **to you.**

The products in this catalog are available for our customers in Europe, Middle East, Africa and Russia



Cortiva™ BioActive Surface

- Durable, non-leaching end point-attached heparin
- Provides thromboresistance and enhanced blood compatibility

All product codes with * are also available with Cortiva™ BioActive Surface. For more information, contact your local Medtronic Cannula Products representative.

Contents

For more information, contact your local Medtronic cannula products representative.

- ▶ **Extracorporeal circuit**
 - ▶ Arterial Cannulae4
 - ▶ Femoral Cannulae 19
 - ▶ Venous Cannulae28
- ▶ **Cardioplegia circuit**
 - ▶ Antegrade Perfusion Cannulae37
 - ▶ Retrograde Perfusion Cannulae47
 - ▶ Cardioplegia Adapters 55
- ▶ **Suction circuit and surgical accessories**
 - ▶ Suction, Sumps, and Vents.....67
 - ▶ Surgical Support Products75
- ▶ **Appendix**
 - ▶ Product Codes..... 86
 - ▶ Index 87
 - ▶ List of products available with Cortiva™ BioActive Surface 93



Alleviate pain. Restore health. Extend life.

Living our mission with you

Our Mission guides our day-to-day work and reminds us that our efforts, large and small, transform the lives of millions of people each year.

We have been working with you to deliver trusted and innovative cannula options for more than 40 years for standard, pediatric, and MICS procedures – more recently ECLS – while offering your favorite cardiac suction, surgical support items, and cardio-kits.

Together, we can deliver the cardiac cannula solutions that matter most to you – and your patients.

Extracorporeal circuit

Arterial cannulae



Arterial cannulae are intended to deliver blood to the patient during cardiopulmonary bypass. Medtronic understands the importance of gentle flow in the often fragile aorta. We offer a number of different product variations to facilitate a variety of surgical techniques and to accommodate varying patient anatomies. Our adult and pediatric arterial cannulae are available in sizes from 6 Fr to 24 Fr in a variety of tip configurations, lengths, and connector sizes to meet your specific needs.

Arterial cannulae



- ▶ EOPA 3D™ 6
- ▶ Select 3D™ II 7
- ▶ EOPA™ Elongated One-Piece 8
- ▶ DLP™ One-Piece Pediatric 9
- ▶ DLP™ Straight Tip 10
- ▶ DLP™ Flexible Arch 11
- ▶ Select Series™ Straight Tip 12
- ▶ Select Series™ Angled Tip 13
- ▶ DLP™ Curved Tip with One-Piece Wirewound Body 14
- ▶ DLP™ Curved Tip Cannulae with Wirewound Body 15
- ▶ DLP™ Curved Tip 16
- ▶ DLP™ Curved Metal Tip 18

Arterial Cannulae

EOPA 3D™ Arterial cannulae



These cannulae offer unique tapered diffuse flow tips and kink-resistant, elongated wirewound bodies. The 3-D tip is designed with three integrated flutes to diffuse and disperse blood flow. Features include adjustable radiopaque suture ring, obturator, hemostasis cap, and depth markings.



Vented 3/8 in (0.95 cm) connector

78220* | 20 Fr (6.7 mm)



78222 | 22 Fr (7.3 mm)

Non-vented 3/8 in (0.95 cm) connector

78320 | 20 Fr (6.7 mm)

78322* | 22 Fr (7.3 mm)

10 per carton



* Also available with Cortiva™ BioActive Surface.

Important Safety Information

EOPA 3D models: Cannula tip must be inserted beyond side ports to minimum depth marks to avoid possibility of tissue dissection of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Arterial Cannulae

Select 3D™ II Arterial cannulae



These cannulae offer beveled tips and tapered, one-piece, kink-resistant wirewound bodies. The tips feature three integrated flutes, which help diffuse and disperse blood flow.



**11.5 in (29.2 cm)
overall length
45° tip**

Vented 3/8 in (0.95 cm) connector


78420	20 Fr (6.7 mm)
78422	22 Fr (7.3 mm)
78424	24 Fr (8.0 mm)



Non-vented 3/8 in (0.95 cm) connector

78520	20 Fr (6.7 mm)
78522	22 Fr (7.3 mm)
78524	24 Fr (8.0 mm)

20 per carton





Arterial Cannulae

EOPA™ Elongated One-Piece Arterial Cannulae



These cannulae offer elongated, one-piece, kink-resistant wirewound bodies with an introducer and hemostasis cap. Features include depth markings and an adjustable radiopaque suture ring.



12 in (30.5 cm) overall length

Blunt tip introducer without guidewire

Vented 3/8 in (0.95 cm) connector

77418	18 Fr (6.0 mm)
77420*	20 Fr (6.7 mm)
77422*	22 Fr (7.3 mm)
77424	24 Fr (8.0 mm)

Non-vented 3/8 in (0.95 cm) connector

77518	18 Fr (6.0 mm)
77520*	20 Fr (6.7 mm)
77522*	22 Fr (7.3 mm)
77524	24 Fr (8.0 mm)

10 per carton



12 in (30.5 cm) overall length

Dilator tip introducer with 0.038 in (0.0965 mm) guidewire

Vented 3/8 in (0.95 cm) connector

77618*	18 Fr (6.0 mm)
77620*	20 Fr (6.7 mm)
77622*	22 Fr (7.3 mm)
77624	24 Fr (8.0 mm)

Non-vented 3/8 in (0.95 cm) connector

77718	18 Fr (6.0 mm)
77720*	20 Fr (6.7 mm)
77722*	22 Fr (7.3 mm)
77724	24 Fr (8.0 mm)

10 per carton

* Also available with Cortiva™ BioActive Surface.

8 Important Safety Information



EOPA models with dilator tip introducer: A guidewire must be used for proper placement of introducer tip while inserting cannula to avoid trauma to the back wall of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgical procedures.

Arterial Cannulae

DLP™ One-Piece Pediatric Arterial Cannulae



These cannulae offer beveled tips with elongated, one-piece, kink-resistant wirewound bodies. Features include an adjustable hemostasis introducer system and depth markings.



9 in (22.9 cm) overall length

Vented 1/4 in (0.64 cm) connector (blue hub)

77006* 6 Fr (2.0 mm)

77008* 8 Fr (2.7 mm)

77010* 10 Fr (3.3 mm)

77012* 12 Fr (4.0 mm)

77014* 14 Fr (4.7 mm)

77016* 16 Fr (5.3 mm)

Non-vented 1/4 in (0.64 cm) connector (blue hub)

77106 6 Fr (2.0 mm)

77108 8 Fr (2.7 mm)

77110 10 Fr (3.3 mm)

77112 12 Fr (4.0 mm)

77114 14 Fr (4.7 mm)

77116 16 Fr (5.3 mm)

5 per carton



9 in (22.9 cm) overall length

Vented 3/16 in (0.48 cm) connector (grey hub)

77206 6 Fr (2.0 mm)

77208 8 Fr (2.7 mm)

77210 10 Fr (3.3 mm)

77212 12 Fr (4.0 mm)

Non-vented 3/16 in (0.48 cm) connector (grey hub)

77306 6 Fr (2.0 mm)

77308 8 Fr (2.7 mm)

77310 10 Fr (3.3 mm)

77312 12 Fr (4.0 mm)

5 per carton






* Also available with Cortiva™ BioActive Surface.

Arterial Cannulae

DLP™ Straight Tip Arterial Cannulae



These pediatric and adult cannulae feature beveled tips and tapered bodies. A variety of sizes are available.



7 in (17.8 cm) overall length

Connection site 1/4 in (0.64 cm)

75008	8 Fr (2.7 mm)
75010	10 Fr (3.3 mm)
75012	12 Fr (4.0 mm)
75014	14 Fr (4.7 mm)

Connection site 3/8 in (0.95 cm)

75318*	18 Fr (6.0 mm)
75320*	20 Fr (6.7 mm)
75322	22 Fr (7.3 mm)
75324	24 Fr (8.0 mm)

20 per carton

8 in (20.3 cm) overall length

Vented 3/8 in (0.95 cm) connector

76020	20 Fr (6.7 mm)
76022	22 Fr (7.3 mm)
76024	24 Fr (8.0 mm)

Non-vented 3/8 in (0.95 cm) connector

76118	18 Fr (6.0 mm)
76120	20 Fr (6.7 mm)
76122	22 Fr (7.3 mm)
76124	24 Fr (8.0 mm)

20 per carton


* Also available with Cortiva™ BioActive Surface.

Arterial Cannulae

DLP™ Flexible Arch Arterial Cannulae



These cannulae offer beveled tips and tapered, one-piece bodies. Features include adjustable ring and depth markings.



11 in (27.9 cm) overall length

Vented 3/8 in (0.95 cm) connector

70420	20 Fr (6.7 mm)
--------------	----------------

70422	22 Fr (7.3 mm)
--------------	----------------

70424	24 Fr (8.0 mm)
--------------	----------------

Non-vented 3/8 in (0.95 cm) connector

71420	20 Fr (6.7 mm)
--------------	----------------

71422	22 Fr (7.3 mm)
--------------	----------------

71424*	24 Fr (8.0 mm)
---------------	----------------

20 per carton


* Also available with Cortiva™ BioActive Surface.

Arterial Cannulae

Select Series™ Straight Tip Arterial Cannulae



These cannulae feature beveled tips and tapered, one-piece, kink-resistant wirewound bodies, and connector peel cap.



12 in (30.5 cm) overall length

Vented 3/8 in (0.95 cm) connector

72120	20 Fr (6.7 mm)
--------------	----------------

72122	22 Fr (7.3 mm)
--------------	----------------

72124*	24 Fr (8.0 mm)
---------------	----------------

Non-vented 3/8 in (0.95 cm) connector

72220	20 Fr (6.7 mm)
--------------	----------------

72222	22 Fr (7.3 mm)
--------------	----------------

72224*	24 Fr (8.0 mm)
---------------	----------------

20 per carton



* Also available with Cortiva™ BioActive Surface.

Arterial Cannulae

Select Series™ Angled Tip Arterial Cannulae



These cannulae feature beveled tips and tapered, one-piece, kink-resistant wirewound bodies, and connector peel cap.



12 in (30.5 cm) overall length

45° tip without side holes

Vented 3/8 in (0.95 cm) connector

72420 20 Fr (6.7 mm)

72422 22 Fr (7.3 mm)

72424 24 Fr (8.0 mm)

Non-vented 3/8 in (0.95 cm) connector

72520 20 Fr (6.7 mm)

72522 22 Fr (7.3 mm)

72524 24 Fr (8.0 mm)

20 per carton



12 in (30.5 cm) overall length

45° tip with side holes

Vented 3/8 in (0.95 cm) connector

73420 20 Fr (6.7 mm)

73422 22 Fr (7.3 mm)


73424 24 Fr (8.0 mm)

Arterial Cannulae

DLP™ Curved Tip Arterial Cannulae with One-piece Wirewound Body



These cannulae feature beveled curved tips and tapered, one-piece, kink-resistant wirewound bodies, tip orientation line, and connector peel cap.



11 in (27.9 cm) overall length

Round flange

Vented 3/8 in (0.95 cm) connector

87220 20 Fr (6.7 mm)

87222* 22 Fr (7.3 mm)

87224 24 Fr (8.0 mm)

Non-vented 3/8 in (0.95 cm) connector

87120 20 Fr (6.7 mm)

87122 22 Fr (7.3 mm)

87124 24 Fr (8.0 mm)

10 per carton

* Also available with Cortiva™ BioActive Surface.

Arterial Cannulae

DLP™ Curved Tip Arterial Cannulae with Wirewound Body



These cannulae feature an overall longer length, beveled curved tips, and wirewound bodies with tip orientation line.




**13 in (33.0 cm)
overall length**

Curved tip with flange

Connection site 3/8 in (0.95 cm)

87024	24 Fr (8.0 mm)
--------------	----------------

10 per carton




**13 in (33.0 cm)
overall length**

Curved tip without flange

Connection site 3/8 in (0.95 cm)

87924	24 Fr (8.0 mm)
--------------	----------------

10 per carton




Arterial Cannulae

DLP™ Curved Tip Arterial Cannulae



These cannulae feature beveled curved tips and tapered bodies with tip orientation line.



7 in (17.8 cm) overall length

Curved tip with flange

Connection site 3/8 in (0.95 cm)

83020	20 Fr (6.7 mm)
83022	22 Fr (7.3 mm)
83024	24 Fr (8.0 mm)

20 per carton

7 in (17.8 cm) overall length

Curved tip without flange

Connection site 3/8 in (0.95 cm)

84020	20 Fr (6.7 mm)
84022	22 Fr (7.3 mm)
84024	24 Fr (8.0 mm)

20 per carton

8 in (20.3 cm) overall length

Curved tip with flange

Vented 3/8 in (0.95 cm) connector

88020	20 Fr (6.7 mm)
88022	22 Fr (7.3 mm)
88024	24 Fr (8.0 mm)

Non-vented 3/8 in (0.95 cm) connector

88120	20 Fr (6.7 mm)
88122	22 Fr (7.3 mm)
88124	24 Fr (8.0 mm)

20 per carton

8 in (20.3 cm) overall length

Curved tip without flange

Vented 3/8 in (0.95 cm) connector

89020	20 Fr (6.7 mm)
89022	22 Fr (7.3 mm)
89024	24 Fr (8.0 mm)

Non-vented 3/8 in (0.95 cm) connector

89120	20 Fr (6.7 mm)
89122	22 Fr (7.3 mm)
89124	24 Fr (8.0 mm)



20 per carton

Arterial Cannulae

DLP™ Curved Tip Arterial Cannulae



These cannulae feature beveled curved tips and tapered bodies with tip orientation line.



10 in (25.4 cm) overall length

Curved tip with flange and soft body

Vented 3/8 in (0.95 cm) connector

81020	20 Fr (6.7 mm)
--------------	----------------

81022	22 Fr (7.3 mm)
--------------	----------------

Non-vented 3/8 in (0.95 cm) connector

81120*	20 Fr (6.7 mm)
---------------	----------------

81122*	22 Fr (7.3 mm)
---------------	----------------

20 per carton

10 in (25.4 cm) overall length

Curved tip without flange and soft body

Vented 3/8 in (0.95 cm) connector

82020	20 Fr (6.7 mm)
--------------	----------------

82022	22 Fr (7.3 mm)
--------------	----------------

Non-vented 3/8 in (0.95 cm) connector

82120	20 Fr (6.7 mm)
--------------	----------------

82122	22 Fr (7.3 mm)
--------------	----------------

20 per carton



* Also available with Cortiva™ BioActive Surface.

Arterial Cannulae

DLP™ Curved Metal Tip Arterial Cannulae



These cannulae feature beveled metal tips and tapered bodies with tip orientation line. A choice of either standard or soft flange is available.



10.5 in (26.7 cm) overall length

Soft flange

Vented 3/8 in (0.95 cm) connector

80020	20 Fr (6.7 mm)
--------------	----------------

80022	22 Fr (7.3 mm)
--------------	----------------

Non-vented 3/8 in (0.95 cm) connector

80120*	20 Fr (6.7 mm)
---------------	----------------

80122	22 Fr (7.3 mm)
--------------	----------------

20 per carton

10.5 in (26.7 cm) overall length

Standard flange

Vented 3/8 in (0.95 cm) connector

80220	20 Fr (6.7 mm)
--------------	----------------

80222	22 Fr (7.3 mm)
--------------	----------------

Non-vented 3/8 in (0.95 cm) connector

80320	20 Fr (6.7 mm)
--------------	----------------

80322	22 Fr (7.3 mm)
--------------	----------------

20 per carton

* Also available with Cortiva™ BioActive Surface.



Extracorporeal circuit

Femoral cannulae



The largest cannula portfolio worldwide.

Our commitment
– to you.

Femoral cannulae provide greater versatility of cannulation techniques to meet the needs of the changing patient population and the desire for smaller thoracic incisions. Designed to be placed in the femoral vein and/or artery, Medtronic femoral cannulae are ideally suited for both percutaneous and direct insertion. The unique flat wirewound, thin-wall construction of the Bio-Medicus™ cannulae maximizes flow rates and cannula flexibility. Our various sizes ensure the availability of a cannula to meet your specific needs.

Femoral cannulae



- ▶ Bio-Medicus™ NextGen Pediatric Arterial and Venous 21
- ▶ Bio-Medicus™ NextGen Arterial or Jugular Venous 22
- ▶ Bio-Medicus™ Percutaneous Kit, Multi-Stage venous 23
- ▶ Bio-Medicus™ NextGen Bi-caval Venous 24
- ▶ Bio-Medicus™ NextGen Insertion Kits 25
- ▶ DLP™ Venous 26
- ▶ DLP™ Arterial 27

Femoral Cannulae

Bio-Medicus™ NextGen Pediatric Arterial and Venous Cannulae




These pediatric cannulae offer a flexible, kink-resistant, thin-wall design that promotes maximum flow rates. The tip design eases insertion and provides a reinforced non-step transition zone along with smooth side holes. Both venous and arterial models are packaged with a tapered solid core and a percutaneous-style, guidewire-capable (0.025 mm dia), non-phthalate PVC introducer, allowing for either a traditional or Seldinger insertion. Each cannula includes a repositionable radiopaque suture ring, an improved friction-fit hemostasis cap, and a retention clip that aids in tubing management. These products are manufactured without DEHP.

Arterial

Non-vented 1/4 in (0.64 cm) connector

		Overall	Tip
96820-108	8 Fr (2.7 mm)	22.9 cm	10.0 cm
96820-110	10 Fr (3.3 mm)	22.9 cm	10.5 cm
96820-112	12 Fr (4.0 mm)	22.9 cm	11.0 cm
96820-114	14 Fr (4.7 mm)	22.9 cm	11.5 cm

1 per carton




Venous


Non-vented 1/4 in (0.64 cm) connector

		Overall	Tip
96830-108	8 Fr (2.7 mm)	22.9 cm	10.0 cm
96830-110	10 Fr (3.3 mm)	22.9 cm	10.5 cm
96830-112	12 Fr (4.0 mm)	22.9 cm	11.0 cm
96830-114	14 Fr (4.7 mm)	22.9 cm	11.5 cm

1 per carton



Bio-Medicus Pediatric Cannulae



Retention Clip



See page 25 for 96553 insertion kit.

Femoral Cannulae

Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannulae



These cannulae offer a flexible, kink-resistant, thin-wall design that promotes maximum flow rates. The tip design eases insertion with a reinforced non-step transition zone, and smooth side holes. A silicone friction-fit cap and non-phthalate PVC introducer provide improved function as well. The percutaneous-style introducer, a repositionable radiopaque suture ring, and a 3/8 inch connection site are standard on all models.

12.5 in (31.8 cm) overall length

7.09 in (18.0 cm) tip

Vented 3/8 in (0.95 cm) connector

Cannula	Cannula Kit	French Size
96570-115	96530-115	15 Fr (5.0 mm)
96570-117	96530-117	17 Fr (5.7 mm)
96570-119	96530-119	19 Fr (6.3 mm)
96570-121	96530-121	21 Fr (7.0 mm)
96570-123	96530-123	23 Fr (7.7 mm)
96570-125	96530-125	25 Fr (8.3 mm)


1 per carton

Order your cannulae as single stand-alone units or with a kit included. For cannula singles, an insertion kit should be ordered.


Cannula kit includes:



White color dilator included with 21 Fr, 23 Fr, and 25 Fr models



Bio-Medicus NextGen Femoral Arterial/Jugular Venous Cannulae



• Bio-Medicus NextGen Femoral Arterial Cannulae
• Bio-Medicus NextGen Jugular Venous Cannulae

➤ See page 25 for 96552 insertion kit.



Femoral Cannulae

Bio-Medicus™ Percutaneous Kit, Multi-Stage Venous Femoral



These adult cannulae feature one-piece, kink-resistant, wirewound bodies with multiple side ports, an introducer with depth markings, and a tip taper location indicator. Insertion kit components are listed below.

Cannula kit includes:




30 in (76.2 cm) overall length
23.6 in (60 cm) tip

Non-vented 3/8 in (0.95 cm) connector

96880-019	19 Fr (6.3 mm)
96880-021	21 Fr (7.0 mm)
96880-025	25 Fr (8.3 mm)

1 per carton




Femoral Cannulae

Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannulae



These cannulae feature a flexible, kink-resistant, thin-wall design that promotes maximum flow rates. The tip design offers an elongated tapered introducer, a reinforced non-step tip transition zone, and smooth side holes. All of these enhancements add up to an easier insertion. Side holes positioned in the SVC and IVC near hepatic veins allow bi-caval drainage and optional isolation of the right atrium. A silicone friction fit cap and non-phthalate PVC introducer provide improved function as well. Each model is standard with a percutaneous-style introducer, a repositionable radiopaque suture ring, and a 3/8 inch connection site.

Cannulae with a kit include:



Bi-caval capability with optimized basket placement

Enables increased venous drainage; ideal for MICS or traditional cases.

Designed to enable right atrium occlusion and access.


Non-vented 3/8 in (0.95 cm) connector

Cannula Singles	Cannula Kit	French Size	Overall Length	Tip Length
96670-115	96600-115	15 Fr (5.0 mm)	64.8 cm	48.9 cm
96670-117	96600-117	17 Fr (5.7 mm)	64.8 cm	48.9 cm
96670-119	96600-119	19 Fr (6.3 mm)	69.9 cm	54.0 cm
96670-121	96600-121	21 Fr (7.0 mm)	69.9 cm	54.0 cm
96670-123†	96600-123†	23 Fr (7.7 mm)	76.2 cm	60.0 cm
96670-125†	96600-125†	25 Fr (8.3 mm)	76.2 cm	60.0 cm
96670-127†	96600-127†	27 Fr (9.0 mm)	76.2 cm	60.0 cm
96670-129†	96600-129†	29 Fr (9.7 mm)	76.2 cm	60.0 cm

1 per carton

Order your cannulae as single stand-alone units or with a kit included. For cannula singles, an insertion kit should be ordered.

† 1/2 in (1.27 cm) to 3/8 in (0.95 cm) tubing adapter included.



See page 25 for 96551 insertion kit.



Femoral Cannulae

Bio-Medicus™ NextGen Insertion Kits



These adult and pediatric insertion kits provide components for peripheral cannulation via Seldinger technique. Designed for percutaneous or direct visualization insertion, these kits provide the following insertion components as listed below.



Venous

96551

- 0.038 in (0.0965 mm) diameter; 180 cm (70.87 in) guidewire
- #11 scalpel blade
- 10 cc syringe
- Stepped vessel dilator, 8 Fr/10 Fr, 12 Fr/14 Fr, and 16 Fr/18 Fr
- Seldinger needle, 18 ga (1.02 mm)

5 kits per carton



Pediatric Arterial/Venous

96553

- 0.025 in (0.0635 mm) diameter; 60 cm (23.62 in) guidewire
- #11 scalpel blade
- 10 cc syringe
- Stepped vessel dilator, 8 Fr/10 Fr, and 12 Fr/14 Fr
- Seldinger needle, 18 ga (1.02 mm)

5 kits per carton



Arterial/Jugular

96552

- 0.038 in (0.0965 mm) diameter; 100 cm (39.37 in) guidewire
- #11 scalpel blade
- 10 cc syringe
- Stepped vessel dilator, 8 Fr/10 Fr, 12 Fr/14 Fr, and 16 Fr/18 Fr
- Seldinger needle, 18 ga (1.02 mm)


5 kits per carton

Femoral Cannulae

DLP™ Femoral Venous Cannulae



These cannulae feature multiple side holes, depth markings, and tapered tip introducer.



25.5 in (64.8 cm) insertion length

Connection site 3/8 in (0.95 cm) and 0.038 in (0.0965 mm) x 180 cm (70.87 in) guidewire

96328	28 Fr (9.3 mm)
--------------	----------------

Non-vented 3/8 in (0.95 cm) connector and 0.038 in (0.0965 cm) x 180 cm (70.87 in) guidewire

96428	28 Fr (9.3 mm)
--------------	----------------

4 per carton





96328

Femoral Cannulae

DLP™ Femoral Arterial Cannulae



These cannulae feature multiple side holes, depth markings, and tapered tip introducer for placement over a 0.038 in (0.0965 cm) guidewire (**not included**).



7 in (17.8 cm) insertion length

Vented 3/8 in (0.95 cm) connector

57414 14 Fr (4.7 mm)

57417 17 Fr (5.7 mm)

57421 21 Fr (7.0 mm)

4 per carton



Extracorporeal circuit

Venous cannulae



Venous cannulae are designed to provide optimal venous drainage while providing flexibility and ease of use. We offer a variety of single stage venous cannulae with kink-resistant, wirewound bodies in sizes from 12 Fr to 40 Fr. These cannulae are offered in both straight and right angle with multiport tips and the option for a thin-wall metal tip design. Our product breadth ensures the availability of a cannula to meet your specific needs.

Venous cannulae





- ▶ DLP™ Single Stage with Right Angle Metal Tip 30
- ▶ DLP™ Single Stage 31
- ▶ DLP™ Single Stage with right angle 32
- ▶ DLP™ Malleable Single Stage 33
- ▶ VC2™ 34
- ▶ MC2™ Two Stage 35
- ▶ MC2X™ Three Stage 36

Venous Cannulae

DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip



These pediatric and adult cannulae feature beveled metal, multiport tips, and kink-resistant wirewound bodies with a tip orientation line. This construction allows for high flow rates with minimal pressure differential.



* Also available with Cortiva™ BioActive Surface.



**14 in (35.6 cm)
overall length**

Connection site 1/4 in (0.64 cm)

67312*	12 Fr (4.0 mm)
67314*	14 Fr (4.7 mm)
67316*	16 Fr (5.3 mm)
67318*	18 Fr (6.0 mm)
67320*	20 Fr (6.7 mm)

Connection site 3/8 in (0.95 cm)

69312	12 Fr (4.0 mm)
69314	14 Fr (4.7 mm)
69316	16 Fr (5.3 mm)
69318	18 Fr (6.0 mm)
69320*	20 Fr (6.7 mm)
69322	22 Fr (7.3 mm)
69324*	24 Fr (8.0 mm)
69328*	28 Fr (9.3 mm)
69331	31 Fr (10.3 mm)

10 per carton



**14 in (35.6 cm)
overall length**

Extended tip

Connection site 3/8 in (0.95 cm)

69428	28 Fr (9.3 mm)
69431	31 Fr (10.3 mm)




10 per carton

Venous Cannulae

DLP™ Single Stage Venous Cannulae



These pediatric and adult cannulae feature tapered, multiport tips and one-piece, kink-resistant wirewound bodies with depth markings.



12-16 in (30.5-40.6 cm) overall length

Connection site 1/4 in (0.64 cm)

66112* 12 Fr (4.0 mm)

66114* 14 Fr (4.7 mm)

66116* 16 Fr (5.3 mm)

Connection site 1/4 in-3/8 in (0.64-0.95 cm)

66118* 18 Fr (6.0 mm)

66120* 20 Fr (6.7 mm)

66122 22 Fr (7.3 mm)

66124* 24 Fr (8.0 mm)

Connection site 3/8 in (0.95 cm)

66126 26 Fr (8.7 mm)

66128* 28 Fr (9.3 mm)

66130 30 Fr (10.0 mm)

66132* 32 Fr (10.7 mm)

66134 34 Fr (11.3 mm)

66236* 36 Fr (12.0 mm)

Connection site 1/2 in (1.27 cm)

66136 36 Fr (12.0 mm)

66140 40 Fr (13.3 mm)

10 per carton



* Also available with Cortiva™ BioActive Surface.

Venous Cannulae

DLP™ Single Stage Venous Cannulae with right angle



These pediatric and adult cannulae feature tapered, multiport tips and one-piece, kink-resistant wirewound bodies with depth markings.



* Also available with Cortiva™ BioActive Surface.

12-15 in (30.5-38.1 cm) overall length

Connection site 1/4 in (0.64 cm)

67512	12 Fr (4.0 mm)
67514	14 Fr (4.7 mm)
67516	16 Fr (5.3 mm)

Connection site 1/4 in-3/8 in (0.64-0.95 cm)

67518	18 Fr (6.0 mm)
67520*	20 Fr (6.7 mm)
67522	22 Fr (7.3 mm)
67524*	24 Fr (8.0 mm)

Connection site 3/8 in (0.95 cm)

67526	26 Fr (8.7 mm)
67528*	28 Fr (9.3 mm)
67530*	30 Fr (10.0 mm)
67532*	32 Fr (10.7 mm)
67534*	34 Fr (11.3 mm)
67636*	36 Fr (12.0 mm)
67638	38 Fr (12.7 mm)





10 per carton

Venous Cannulae

DLP™ Malleable Single Stage Venous Cannulae



These pediatric and adult cannulae feature multiport tips and malleable, one-piece, kink-resistant wirewound bodies with depth markings. This construction allows the surgeon to shape the cannulae. An integrated malleable wire outside the wirewinding optimizes the inner diameter (I.D.) for maximum flow performance.



12-15 in (30.5-38.1 cm) overall length

Connection site 1/4 in (0.64 cm)

68112*	12 Fr (4.0 mm)
68114*	14 Fr (4.7 mm)
68116*	16 Fr (5.3 mm)

Connection site 1/4 in-3/8 in (0.64-0.95 cm)

68118*	18 Fr (6.0 mm)
68120*	20 Fr (6.7 mm)
68122	22 Fr (7.3 mm)
68124*	24 Fr (8.0 mm)

Connection site 3/8 in (0.95 cm)

68126	26 Fr (8.7 mm)
68128*	28 Fr (9.3 mm)
68130	30 Fr (10.0 mm)
68132*	32 Fr (10.7 mm)
68134	34 Fr (11.3 mm)

Connection site 1/2 in (1.27 cm)

68136*	36 Fr (12.0 mm)
68138	38 Fr (12.7 mm)

10 per carton



The malleable single stage venous cannulae allow the surgeon to manipulate the cannula shape into virtually any position.






* Also available with Cortiva™ BioActive Surface.

Venous Cannulae

VC2™ Venous Cannulae



These cannulae offer swirl tips with slotted atrial baskets and kink-resistant, wirewound bodies that feature depth markings. The oval body presents a lower profile in the surgical field.



15 in (38.1 cm) overall length

Wirewound body
Connection site 3/8 in (0.95 cm)

93438	34/38 Fr (11.3/12.7 mm)
--------------	-------------------------

Oval wirewound body
Connection site 1/2 in (1.27 cm)

93464	34/48 Fr (11.3/16.0 mm)
--------------	-------------------------

Wirewound body
Connection site 1/2 in (1.27 cm)

93448	34/48 Fr (11.3/16.0 mm)
--------------	-------------------------

Oval wirewound body, non-vented
1/2 in (1.27 cm) connector

93463C	34/38 Fr (11.3/12.7 mm)
---------------	-------------------------

93464C	34/48 Fr (11.3/16.0 mm)
---------------	-------------------------

10 per carton

Wirewound body, non-vented
1/2 in (1.27 cm) connector

93438C*	34/38 Fr (11.3/12.7 mm)
----------------	-------------------------

93448C	34/48 Fr (11.3/16.0 mm)
---------------	-------------------------

Oval wirewound body
Connection site 3/8 in (0.95 cm)

93463	34/38 Fr (11.3/12.7 mm)
--------------	-------------------------



* Also available with Cortiva™ BioActive Surface.

Venous Cannulae

MC2™ Two Stage Venous Cannulae



These cannulae feature multiport tips with atrial baskets and kink-resistant, wirewound bodies with depth markings. The oval body presents a lower profile in the surgical field.



* Also available with Cortiva™ BioActive Surface.

15 in (38.1 cm) overall length

Connection site 3/8 in (0.95 cm)

91228* 28/36 Fr (9.3/12.0 mm)

Connection site 1/2 in (1.27 cm)

91236 36/46 Fr (12.0/15.3 mm)

91240* 32/40 Fr (10.7/13.3 mm)

91246* 34/46 Fr (11.3/15.3 mm)

91251 36/51 Fr (12.0/17.0 mm)

Non-vented 1/2 in (1.27 cm) connector

91228C* 28/36 Fr (9.3/12.0 mm)

91236C* 36/46 Fr (12.0/15.3 mm)

91240C* 32/40 Fr (10.7/13.3 mm)

91246C* 34/46 Fr (11.3/15.3 mm)

91251C 36/51 Fr (12.0/17.0 mm)

10 per carton

15 in (38.1 cm) overall length

Reduced stage body Connection site 1/2 in (1.27 cm)

91229 29/37 Fr (9.7/12.3 mm)

Non-vented 1/2 in (1.27 cm) connector

91229C 29/37 Fr (9.7/12.3 mm)

Thin-wall body Connection site 3/8 in (0.95 cm)

91329* 29/29 Fr (9.7/9.7 mm)

Non-vented 3/8 in (0.95 cm) connector

91329C* 29/29 Fr (9.7/9.7 mm)

Oval body Connection site 1/2 in (1.27 cm)

91263* 32/40 Fr (10.7/13.3 mm)

91265* 36/46 Fr (12.0/15.3 mm)

Oval body Non-vented 1/2 in (1.27 cm) connector

91263C* 32/40 Fr (10.7/13.3 mm)

91265C* 36/46 Fr (12.0/15.3 mm)


10 per carton

Venous Cannulae

MC2X™ Three Stage Venous Cannulae



These cannulae feature a thin-wall design that provides minimal pressure differential for use with vacuum or kinetically assisted venous drainage.



15 in (38.1 cm) overall length


Connection site 3/8 in (0.95 cm)

91429* | 29/29/29 Fr (9.7/9.7/9.7 mm)

Non-vented
3/8 in (0.95 cm) connector

91429C* | 29/29/29 Fr (9.7/9.7/9.7 mm)

10 per carton



15.5 in (39 cm) overall length

Connection site 3/8 in (0.95 cm)

91437 | 29/46/37 Fr (9.7/15.3/12.3 mm)


Non-vented
3/8 in (0.95 cm) connector

91437C | 29/46/37 Fr (9.7/15.3/12.3 mm)

10 per carton



Thinner Wall Design



MC2X has a 31% thinner wall for a larger inner diameter when compared to competitive cannula²



Larger end hole of the MC2X allows for better flow opportunity³
^{2,3} 2013 Medtronic data on file.

* Also available with Cortiva™ BioActive Surface.



Cardioplegia circuit

Antegrade cannulae



Antegrade cannulae are designed to deliver cardioplegia solution to the heart via the arteries in the normal direction of blood flow (antegrade perfusion). Aortic root cannulae are intended to either aspirate air from the aorta or infuse cardioplegia solution directly into the aortic root near the aortic valve. We offer adult and pediatric aortic root cannulae in varying lengths to accommodate various techniques and patient anatomies. Our product breadth ensures the availability of a cannula to meet your specific needs.

Antegrade cannulae




- ▶ MiAR™ Aortic Root with Flow-Guard 39
- ▶ DLP™ Dual Lumen Cardioplegia 40
- ▶ DLP™ Aortic Root 41
- ▶ DLP™ Aortic Root with Vent Line 42
- ▶ DLP™ Cardioplegia Needles 43
- ▶ DLP™ Coronary Ostial Perfusion 44

Antegrade Cannulae

MiAR™ Aortic Root Cannulae with Flow-Guard



These Minimally Invasive Aortic Root (MiAR) models feature a 12.25 in (31 cm) cannula body to allow access into the chest for small incision techniques. This length also allows positioning of the luer connection to the cardioplegia line out of the surgeon's way, outside the chest incision. Additionally, these MiAR models include the Flow-Guard™ feature to maintain hemostasis during removal of the needle from the cannula.




**12.25 in (31 cm)
overall length**

**Flanged standard tip and
Flow-Guard introducer**

11012L 12 ga (9 Fr)

11014L 14 ga (7 Fr)

10 per carton



Important Safety Information

Care should be taken when inserting the needle to prevent perforation of the back wall of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures. Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques.

Antegrade Cannulae

DLP™ Dual Lumen Cardioplegia Cannulae





The dual lumen cannula features a tip attached to a clear body with a separate vent line. Additional features include left heart venting and simultaneous cardioplegia delivery and aspiration. Each cannula is supplied with a stainless steel introducer needle.

**3.25 in (8.3 cm)
overall length**

Standard tip

30401	12 ga (9 Fr)
--------------	--------------

10 per carton



Dual Lumen Aortic Root
Application Detail

Important Safety Information


Care should be taken when inserting the needle to prevent perforation of the back wall of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures. Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques.

Antegrade Cannulae

DLP™ Aortic Root Cannulae



These pediatric and adult cannulae feature radiopaque tips attached to clear PVC bodies. Additional features available with these cannulae include aortic root pressure monitoring and left heart venting. All cannulae are supplied with a stainless steel introducer needle. All DLP Aortic Root Cannulae can be used to aspirate air emboli as well as administer cardioplegia.



Flow-Guard™ Introducer

2.5 in (6.4 cm) overall length

Flanged standard tip and standard introducer

10218	18 ga (4 Fr) white tip and clear flange
12218	18 ga (4 Fr) blue one-piece tip and flange

20 per carton



5.75 in (14.6 cm) overall length

Slotted long tip and standard introducer

10112	12 ga (9 Fr)
10114	14 ga (7 Fr)

20 per carton



5.5 in (14.0 cm) overall length

Flanged standard tip and standard introducer

10009	9 ga (11 Fr)
10012*	12 ga (9 Fr)
10014	14 ga (7 Fr)
10016	16 ga (5 Fr)
10018	18 ga (4 Fr)

Flanged standard tip and Flow-Guard introducer

11012	12 ga (9 Fr)
11014	14 ga (7 Fr)

20 per carton




5.25 in (13.3 cm) overall length

Flanged pressure monitoring tip and standard introducer

23009	9 ga (11 Fr)
--------------	--------------

20 per carton



Flow-Guard Introducer Application Detail

* Also available with Cortiva™ BioActive Surface.

Important Safety Information




Care should be taken when inserting the needle to prevent perforation of the back wall of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures. Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques.

Antegrade Cannulae

DLP™ Aortic Root Cannulae with Vent Line



These cannulae feature radiopaque tips attached to clear bodies with separate vent lines. Additional features available with these cannulae include aortic root pressure monitoring and left heart venting. All cannulae are supplied with a stainless steel introducer needle.



5.5 in (14.0 cm) overall length

Flanged standard tip and standard introducer

20009	9 ga (11 Fr)
20012*	12 ga (9 Fr)
20012S	12 ga (9 Fr) with two clamps
20014	14 ga (7 Fr)
20014L	14 ga (7 Fr), 8 in (20.3 cm) vent line
20016	16 ga (5 Fr)

Flanged standard tip and Flow-Guard™ introducer

21012	12 ga (9 Fr)
21014	14 ga (7 Fr)

20 per carton



5.75 in (14.6 cm) overall length

Slotted long tip and standard introducer

20112	12 ga (9 Fr)
20114	14 ga (7 Fr)

20 per carton



5.75 in (14.6 cm) overall length

Flanged long tip and standard introducer

20114WF	14 ga (7 Fr)
20114HWF	14 ga (7 Fr) with side holes

20 per carton



5.25 in (13.3 cm) overall length

Flanged pressure monitoring tip and standard introducer

24009	9 ga (11 Fr)
--------------	--------------

20 per carton



* Also available with Cortiva™ BioActive Surface.

Important Safety Information

Care should be taken when inserting the needle to prevent perforation of the back wall of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures. Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques.

Antegrade Cannulae

DLP™ Cardioplegia Needles



These needles feature stainless steel tips with a plastic depth stop. These needles are attached to flexible tubing with drape clamp and female luer.



Pediatric

6.5 in (16.5 cm) overall length

Tip length 1/4 in (0.64 cm) with hub stop

11316	16 ga (5 Fr)
-------	--------------

20 per carton

Adult

10 in (25.4 cm) overall length

Tip length 5/8 in (1.59 cm)
with 4 side holes and flange stop

10313	13 ga (8 Fr)
-------	--------------

20 per carton

Important Safety Information

Care should be taken when inserting the needle to prevent perforation of the back wall of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures. Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques.

Antegrade Cannulae

DLP™ Coronary Ostial Perfusion Cannulae




These cannulae feature silicone bodies with soft bulb-shaped tips. All cannulae have a female luer connection site.




10 in (25.4 cm) overall length

30315	15 Fr (5.0 mm) bulb
30317	17 Fr (5.7 mm) bulb
30320	20 Fr (6.7 mm) bulb

10 per carton



Intracoronary cardioplegia



Important Safety Information



Care should be taken when inserting the needle to prevent perforation of the back wall of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures. Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques.

Antegrade Cannulae

DLP™ Coronary Ostial Perfusion Cannulae



These cannulae feature flanged, radiopaque basket tips or soft tips attached to malleable stainless steel shafts. All cannulae have a female luer connection site.



**6 in (15.2 cm)
overall length**

Basket tip

30010	10 Fr (3.3 mm)
30012	12 Fr (4.0 mm)
30014	14 Fr (4.7 mm)

20 per carton



**6 in (15.2 cm)
overall length**

Soft, concave tip

30050

10 per carton



**6 in (15.2 cm)
overall length**

Soft, convex tip

30055

10 per carton

Important Safety Information


Care should be taken when inserting the needle to prevent perforation of the back wall of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures. Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques.

Antegrade Cannulae


DLP™ Coronary Ostial Perfusion Cannulae



These high flow cannulae feature flanged, radiopaque basket tips or soft silicone tips attached to large-bore, malleable stainless steel shafts. All cannulae have a female luer connection site.



Ostial antegrade cardioplegia



7.5 in (19.1 cm) overall length

90° angle basket tip

30110 | 10 Fr (3.3 mm)

30112 | 12 Fr (4.0 mm)

30114 | 14 Fr (4.7 mm)

10 per carton



7.5 in (19.1 cm) overall length

90° angle soft silicone tip

30155 | 10 Fr (3.3 mm)

10 per carton



7.5 in (19.1 cm) overall length

45° angle basket tip

30212 | 12 Fr (4.0 mm)

10 per carton



7.5 in (19.1 cm) overall length

45° angle soft silicone tip

30255 | 10 Fr (3.3 mm)

10 per carton

Important Safety Information

Care should be taken when inserting the needle to prevent perforation of the back wall of the aorta. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures. Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques.

Cardioplegia circuit

Retrograde cannulae



Retrograde cannulae are designed to deliver cardioplegia solution to the heart via the coronary sinus in the reverse direction of normal blood flow (retrograde perfusion). We offer the retrograde cannulae in pediatric and adult sizes. The pediatric cannulae are offered with silicone bodies, with or without a guidewire stylet, and a manual-inflate cuff. Adult models are offered with a silicone or PVC body, a solid or guidewire stylet with a variety of handle options, and an auto- or manual-inflate cuff. Silicone bodies ensure product malleability under cold cardioplegia temperatures. The purpose of the cuff is to occlude the opening of the coronary sinus, forcing cardioplegia flow to the coronary vessel bed of the heart. Manual-inflate models require the clinician to inflate the cuff prior to delivering cardioplegia. Our product breadth ensures the availability of a cannula to meet your specific needs.

Retrograde cannulae





- ▶ Gundry™ Retrograde Coronary Sinus Perfusion with Manual-Inflate Cuff 49
- ▶ DLP™ Retrograde Coronary Sinus Perfusion with Manual-Inflate Cuff 50
- ▶ DLP™ Retrograde Coronary Sinus Perfusion with Auto-Inflate Cuff 52

Retrograde Cannulae

Gundry™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff



These pediatric and adult cannulae feature silicone bodies and manual-inflate cuffs with pressure monitoring lines. All cannulae have a guidewire stylet and either a male luer handle or a Tru-Touch™ handle. All cannulae come with a syringe for cuff inflation.



9 in (22.9 cm) overall length

Smooth cuff and wirewound body

94110	10 Fr (3.3 mm) guidewire styletw
--------------	----------------------------------

4 per carton

12.5 in (31.8 cm) overall length

Smooth cuff and wirewound body

94115T	15 Fr (5.0 mm) guidewire stylet
---------------	---------------------------------

10 per carton

12.5 in (31.8 cm) overall length

Smooth cuff, wirewound body, and integral stopcock

94113T	13 Fr (4.3 mm) guidewire stylet
---------------	---------------------------------

10 per carton

94615	15 Fr (5.0 mm) guidewire stylet
--------------	---------------------------------

4 per carton

94715	15 Fr (5.0 mm) guidewire stylet with Tru-Touch handle
--------------	---

10 per carton

Important Safety Information



Extreme caution should be exercised while introducing the cannula into the coronary sinus. Do not force the cannula into the coronary sinus as this may cause vessel damage. Do not overinflate the balloon.

Retrograde Cannulae

DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff



These pediatric and adult cannulae feature silicone bodies and manual-inflate cuffs with pressure monitoring lines. All cannulae come with a smooth cuff and a syringe for cuff inflation.



9.5 in (24.1 cm) overall length

Smooth cuff, non-wirewound body, and integral stopcock

94006	6 Fr (2.0 mm) no stylet
94106	6 Fr (2.0 mm) guidewire stylet

4 per carton

12.5 in (31.8 cm) overall length

Smooth cuff and wirewound body

94215T	15 Fr (5.0 mm) solid stylet
---------------	-----------------------------

10 per carton

12.5 in (31.8 cm) overall length

Smooth cuff, wirewound body, and integral stopcock

94725	15 Fr (5.0 mm) solid stylet with Tru-Touch™ handle
--------------	--

10 per carton

Important Safety Information



Extreme caution should be exercised while introducing the cannula into the coronary sinus. Do not force the cannula into the coronary sinus as this may cause vessel damage. Do not overinflate the balloon.

Retrograde Cannulae

DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff



These pediatric and adult cannulae feature silicone bodies and manual-inflate cuffs with pressure monitoring lines. A variety of cuffs and stylets are available. All cannulae come with a syringe for cuff inflation.



12.5 in (31.8 cm) overall length

Elongated cuff, wirewound body, and integral stopcock

94625	15 Fr (5.0 mm) solid stylet
94665	15 Fr (5.0 mm) guidewire stylet

10 per carton



12.5 in (31.8 cm) overall length

Ridged cuff, wirewound body, and integral stopcock

94913	13 Fr (4.3 mm) guidewire stylet
94913L	13 Fr (4.3 mm) guidewire stylet and 6 in (15.2 cm) pressure monitoring and inflation line
94915	15 Fr (5.0 mm) guidewire stylet
94965	15 Fr (5.0 mm) guidewire stylet with Tru-Touch™ handle
94975	15 Fr (5.0 mm) solid stylet with Tru-Touch handle

10 per carton

Important Safety Information


Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques. Due to limitations of direct visualization during minimally invasive techniques, echocardiographic or fluoroscopic imaging is recommended. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Retrograde Cannulae

DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff



These cannulae feature silicone bodies and auto-inflate cuffs with pressure monitoring lines. A variety of cuffs and stylets are available.



12.5 in (31.8 cm) overall length

Smooth preformed cuff and wirewound body

94315T	15 Fr (5.0 mm) solid stylet
94415T	15 Fr (5.0 mm) guidewire stylet
94735	15 Fr (5.0 mm) solid stylet with Tru-Touch™ handle
94745	15 Fr (5.0 mm) guidewire stylet with Tru-Touch handle

10 per carton




12.5 in (31.8 cm) overall length

Ridged preformed cuff, wirewound body, and integral stopcock

94995	15 Fr (5.0 mm) guidewire stylet with Tru-Touch handle
--------------	---

10 per carton



Silicone Auto-inflate Cuff Flow Detail

Silicone auto-inflating cuffs offer the convenience of cuff inflation without the need for a syringe. Their unique flow-through design allows cardioplegia to circulate through the cuff before exiting the cannula tip. These cuffs also feature a dual pressure monitoring system, which normally provides intra-sinus pressure monitoring through the port in the tip. However, should the tip inadvertently become occluded, the second port will monitor intra-cuff pressure. This feature greatly reduces the potential for coronary sinus damage due to overinflation.

Important Safety Information

Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques. Due to limitations of direct visualization during minimally invasive techniques, echocardiographic or fluoroscopic imaging is recommended.

Retrograde Cannulae

DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff



These cannulae feature PVC bodies and polyurethane auto-inflate smooth cuffs with pressure monitoring lines. A variety of stylets are available.



10 in (25.4 cm) overall length

Smooth preformed cuff, short fluted bullet-nosed tip, and integral stopcock

94885K

15 Fr (5.0 mm) solid stylet with Tru-Touch™ handle

10 per carton



11 in (27.9 cm) overall length

Smooth preformed cuff, multiport tip, and integral stopcock


94533

13 Fr (4.3 mm) solid stylet with Tru-Touch handle

94535

15 Fr (5.0 mm) solid stylet with Tru-Touch handle

10 per carton



Important Safety Information

Extreme caution should be exercised while introducing the cannula into the coronary sinus. Do not force the cannula into the coronary sinus as this may cause vessel damage. Do not overinflate the balloon.

Retrograde Cannulae

DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff



These cannulae feature PVC bodies and polyurethane auto-inflate cuffs with pressure monitoring lines. A variety of cuffs and stylets are available.



12.5 in (31.8 cm) overall length

Smooth preformed cuff, fluted bullet-nosed tip, and integral stopcock

94835

15 Fr (5.0 mm) solid stylet with Tru-Touch handle

10 per carton




12.5 in (31.8 cm) overall length

Ridged preformed silicone cuff, short fluted bullet-nosed tip, and integral stopcock

94935

15 Fr (5.0 mm) solid stylet with Tru-Touch handle

10 per carton



PVC Auto-inflate Cuff Flow Detail



94885

Stylet



12 in (30.5 cm) overall length

Smooth preformed cuff, short fluted bullet-nosed tip, and integral stopcock

94885

15 Fr (5.0 mm) solid stylet with Tru-Touch™ handle

94895

15 Fr (5.0 mm) guidewire stylet with Tru-Touch handle

10 per carton

Important Safety Information

Additional care and caution may be necessary due to the unique adaptations required for minimally invasive techniques. Due to limitations of direct visualization during minimally invasive techniques, echocardiographic or fluoroscopic imaging is recommended. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.



Cardioplegia circuit

Cardiograde adapters



Cardioplegia adapters are designed to permit customization of the cardioplegia circuit. Medtronic adapters can be used to connect multiple cardioplegia cannulae and/or vein graft cannulae to a single inlet source, switch between antegrade and retrograde delivery of cardioplegia, and provide the option for attaching a vent line to cardioplegia cannulae. Additional adapters are available to extend the length of the cardioplegia circuit, recirculate the cardioplegia solution, or provide an extension to the pressure monitoring line.

Cardioplegia adapters



- ▶ DLP™ Perfusion 57
- ▶ DLP™ Perfusion/Venting 57
- ▶ DLP™ Extension Line 58
- ▶ DLP™ Pressure Monitoring Extension Line 58
- ▶ DLP™ Straight 59
- ▶ DLP™ "Y" 59
- ▶ DLP™ "Y" Type Recirculating 60
- ▶ DLP™ "Y" Type Coronary Perfusion 60
- ▶ DLP™ "Y" 61
- ▶ DLP™ Vein Graft Cannula 62
- ▶ DLP™ Cardioplegia Adapter with Pressure Port 63
- ▶ ARISS™ Perfusion Switches 63
- ▶ DLP™ Cardioplegia Management Set 64
- ▶ DLP™ Antegrade/Retrograde Stopcock 64
- ▶ DLP™ Multiple Perfusion Sets 65

Cardioplegia Adapters

DLP™ Perfusion Adapters



These adapters feature a barbed connector on one end and either a male luer or female luer on the other. Some models have a length of tubing between the barbed connector and luers.

2 in (5.1 cm) overall length

Male luer

10007	1/4 in (0.64 cm) barbed connector
10008	3/8 in (0.95 cm) barbed connector

40 per carton



3.25 in (8.3 cm) overall length

Tubing and female luer

10022	1/4 in (0.64 cm) barbed connector
--------------	-----------------------------------

40 per carton



DLP™ Perfusion/Venting Adapter

This adapter is a combination of perfusion adapter 10007 with an added high-flow stopcock that allows for venting. This 3-way stopcock has a large-bore inner diameter and has female and male luers on the ends.

Male luer

13002

20 per carton



Important Safety Information

Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Cardioplegia Adapters

DLP™ Extension Line Adapters



These adapters feature soft (silicone-like) PVC tubing with or without a clamp, and male and winged female luers on the ends.

10 in (25.4 cm) overall length

Male and winged female luers

11001	without clamp
--------------	---------------

10 per carton



20 in (50.8 cm) overall length

Male and winged female luers

11001G	with clamp
---------------	------------

20 per carton

DLP™ Pressure Monitoring Extension Line Adapters

These adapters feature combinations of male and female luers on the ends.


6 ft (1.8 m) overall length

Male and female luers

25009	male x female luers
--------------	---------------------

25010	male x male luers
--------------	-------------------

20 per carton



Important Safety Information

Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Cardioplegia Adapters

DLP™ Straight Adapters



These adapters have a male luer on one end and a 1/4 inch (0.64 cm) connector on the other.

5 in (12.7 cm) overall length

Male luers

10001	with vented 1/4 in connector
10001K	with vented 1/4 in connector and clamp

20 per carton



DLP™ "Y" Adapters

These adapters feature a male luer on the single leg. A female luer is on the cardioplegia inlet leg and a 1/4 inch (0.64 cm) connector is on the venting leg.

Male luer on single leg

10005†	7.5 in (19.1 cm) overall length with 3 in (7.6 cm) inlet leg and color-coded clamps
10534	8.5 in (21.6 cm) overall length with one-way valve in vent line and color-coded clamps

20 per carton



†When used with 10005 Venting "Y" Adapter, an Aortic Root Cannula can be used to vent the left heart.

Important Safety Information

Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Cardioplegia Adapters

DLP™ "Y" Type Recirculating Adapters



These adapters feature a soft "Y" that has a two-leg option with 1/4 inch (0.64 cm) connectors and venting capabilities or a three-leg option that terminates with male luers.

8 in (20.3 cm) overall length

Male luer on single leg

10003	color-coded clamps
--------------	--------------------

20 per carton



DLP™ "Y" Type Coronary Perfusion Adapters

These adapters feature a rigid "Y" with two legs in various lengths that terminate with male luers.



Female luer on single leg

10004	3.5 in (8.9 cm) and 10 in (25.4 cm) legs and white clamps
--------------	---

10004D	6 in (15.2 cm) legs and white clamps
---------------	--------------------------------------


10700	3.5 in (8.9 cm) and 10 in (25.4 cm) legs and color-coded clamps
--------------	---

10706	6 in (15.2 cm) soft, silicone-like PVC legs and white clamps
--------------	--

10710	10 in (25.4 cm) soft, silicone-like PVC legs and white clamps
--------------	---

20 per carton

"Y" Type Coronary Perfusion Adapter



Important Safety Information

Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Cardioplegia Adapters

DLP™ "Y" Adapter



This adapter is a combination of Models 10004 and 10005 and features four legs. Two legs terminate in male luers, one terminates with a female luer, and the fourth terminates with a 1/4 inch (0.64 cm) connector.

20.5 in (52.1 cm) overall length

Combination with four legs

10005BH	color-coded clamps
---------	--------------------

20 per carton



This adapter features a female luer on the cardioplegia inlet leg and a 1/4 inch (0.64 cm) connector on the venting leg.

7.5 in (19.1 cm) overall length

Male luer on single leg

10005OS	white clamps with color-coded arrows
---------	--------------------------------------

20 per carton



Important Safety Information

Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Cardioplegia Adapters

DLP™ Vein Graft Cannula



This product features a soft silicone tip with graduated steps to fit any vessel. The tip is attached to a "Y" adapter via tubing with clamps. The "Y" adapter terminates with a male luer.




25.5 in (63.5 cm) overall length

Male luer on single inlet leg

10010

"Y" adapter with female luer and white clamp on male leg

20 per carton



The combination of a Vein Graft Cannula and an Aortic Root Cannula provides a way to perfuse the aortic root and vein graft. In this case the vein graft must be held with the fingers and not tied.

Important Safety Information

Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Cardioplegia Adapters

DLP™ Cardioplegia Adapter with Pressure Port



This adapter features a latex-free rubber septum to allow access for a pressure monitoring needle near the connection to the cardioplegia cannula.



Locking male luer

15004

10 per carton

ARISS™ Perfusion Switches


These table mount models allow the surgeon to switch between antegrade and retrograde cardioplegia delivery. These versions feature pressure monitoring lines that run alongside the inlet and outlet ports. The outlet ports are color coded for easy identification.

13 in (33.0 cm) overall length

Table mount models

13000	unequal legs
13004	equal legs

10 per carton



Important Safety Information


Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Cardioplegia Adapters

DLP™ Cardioplegia Management Set



This product features soft tubing with three legs that terminate in male luer. The legs have varying lengths and are set in a mounting bracket with adhesive pads on the back. The red clamp leg has a length of 10 inches (25.4 cm), the white clamp leg has a length of 20 inches (50.8 cm), and the blue clamp leg has a length of 30 inches (76.2 cm).



Female luer on single leg

13003	color-coded clamps
-------	--------------------

5 per carton

DLP™ Antegrade/Retrograde Stopcock Adapter

This adapter features a large-bore, three-way stopcock and male luers on the outlet ports.

26.5 in (67.3 cm) overall length

Female luer on single leg

13001

10 per carton



Important Safety Information

Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Cardioplegia Adapters

DLP™ Multiple Perfusion Sets



These products allow for simultaneous perfusion of the aortic root and up to three or more vein grafts. They feature inlet ports with male or female luer and clamps attached to an adapter that splits into four or more legs. These legs have color-coded clamps for easy identification and terminate with male luer. Models with a vent line terminate in a 1/4 inch (0.64 cm) connector.

12 in (30.5 cm) overall length

Female luer on single inlet leg with 4 legs

14003	4 - 8 in (20.3 cm) length legs with 3 red clamps and 1 blue clamp
14013	4 - 10 in (25.4 cm) length legs with red, blue, yellow, and white clamps

10 per carton



14000



14001

15 in (38.1 cm) overall length

Female luer on single inlet leg with 4 legs

14000	1 - 3 in (7.6 cm) length leg with blue clamp 3 - 10 in (25.4 cm) length legs with red clamps
--------------	---

10 per carton

Female luer on single inlet leg with 4 legs and vent line

14001	1 - 8 in (20.3 cm) length vented leg with red clamp on vent line 3 - 10 in (25.4 cm) length legs with red clamps
14017	1 - 8 in (20.3 cm) length vented leg with stiff tubing and red clamp on vent line 3 - 10 in (25.4 cm) lengths stiff tubing and red clamps
14031	1 - 12 in (30.5 cm) length vented leg with yellow clamp on vent line 3 - 10 in (25.4 cm) length legs with red clamps

10 per carton



Important Safety Information


Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Cardioplegia Adapters

DLP™ Multiple Perfusion Sets



These products allow for simultaneous perfusion of the aortic root and up to three or more vein grafts. They feature inlet ports with male or female luer and clamps attached to an adapter that splits into four or more legs. These legs have color-coded clamps for easy identification and terminate with male luers.



14007

15 in (38.1 cm) overall length

Female luer on single inlet leg with 4 legs and 3 free-flow vessel cannulae (Model 30003)

14007	1 - 3 in (7.6 cm) length leg with blue clamp
	3 - 10 in (25.4 cm) length legs with red clamps

10 per carton

Female luer on single inlet leg with 4 legs and 3 one-way valve vessel cannulae (Model 30001)


14008	1 - 3 in (7.6 cm) length leg with blue clamp
	3 - 10 in (25.4 cm) length legs with red clamps

10 per carton

Female luer on single inlet leg with 6 legs

14006	1 - 3 in (7.6 cm) length leg with blue clamp
	5 - 10 in (25.4 cm) length legs with red clamps

10 per carton



14006

Important Safety Information

Cardioplegia adapter devices are a combination of products that have their own indications, contraindications, precautions, warnings, and potential adverse events. Refer to those components individually. Cardioplegia adapter devices are not intended for use except as indicated. Use the lowest possible pressure to reduce risk of hemolysis or risk of vessel damage. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.



Products

Suction circuit



We offer suction products that are designed to provide either fine or gross suction to ensure a clear and dry surgical field. The Medtronic handheld cardiac suction tubes are designed for fine suction, while our intracardiac suction tubes and sumps are intended to provide suction within and around the heart. The Medtronic left heart vent catheters are designed for indirect and direct venting of the left ventricle and feature perforated tips. These vents are offered in different lengths, introducer styles, and in both silicone and PVC body construction. Our product breadth ensures the availability of a cannula to meet your specific needs.

Suction circuit




- ▶ DLP™ Cardiac Suction Tubes 69
- ▶ DLP™ Pericardial Sumps 72
- ▶ DLP™ Pericardial/Intracardiac Sumps 72
- ▶ DLP™ Cardiac Sump 73
- ▶ DLP™ Left Heart Vent Catheters 74

Suction circuit

DLP™ Cardiac Suction Tubes



These suction tubes feature comfortable hand grips with a vacuum control port. The tip is made from a malleable stainless steel shaft, which allows the surgeon to gently shape to a desired preference. Models are available with or without a soft silicone tip. All suction tubes are equipped with a length of tubing and clamp terminating with a 1/4 inch (0.64 cm) connector.



Frazier tip



Soft tip



3 in (7.6 cm)

approx. shaft and tip length

Frazier tip

10051	6 Fr (2.0 mm) shaft
--------------	---------------------

Soft tip

10053	6 Fr (2.0 mm) shaft with 10 Fr (3.3 mm) tip
--------------	--

20 per carton

6 in (15.2 cm)

approx. shaft and tip length

Frazier tip

10050	6 Fr (2.0 mm) shaft
--------------	---------------------

Soft tip

10052	6 Fr (2.0 mm) shaft with 10 Fr (3.3 mm) tip
--------------	--


20 per carton

Suction circuit

DLP™ Cardiac Suction Tubes



This suction tube features a comfortable hand grip. The tip is made from a malleable stainless steel shaft, which allows the surgeon to gently shape to a desired preference. The DLP Intracardiac Suction device offers a perforated pool tip to minimize the possibility of tissue occlusion while maintaining an opportunity for drainage. This suction tube is equipped with a length of tubing terminating with a 1/4 inch (0.64 cm) connector.



6 in (15.2 cm) overall length

Perforated pool tip

10060

10 Fr (3.3 mm) shaft with 20 Fr (6.7 mm) tip (28 holes in tip)

20 per carton

Suction circuit

DLP™ Cardiac Suction Tubes



These suction tubes feature a comfortable hand grip. The tip shaft is made from rigid stainless steel. The Macro Sucker offers a fluted pool tip to minimize the possibility of tissue occlusion while maintaining an opportunity for drainage. The Mini Sucker has a soft, vented tip designed to clear blood from the surgical site with gentle suction.



7 in (18.0 cm)

approx. shaft and tip length

Macro Sucker

10061	20 Fr (6.7 mm) tip with 16 Fr (5.3 mm) shaft
--------------	--

20 per carton



4.75 in (12.1 cm)

approx. shaft and tip length

Mini Sucker

10062	9 Fr (3.0 mm) tip with 11 Fr (3.7 mm) shaft
--------------	---

20 per carton

Important Safety Information

Neither Model 10061 nor Model 10062 should be placed through the mitral valve. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Suction circuit DLP™ Pericardial Sump



This sump features a fluted tip that is encased in a stainless steel spring and has a weight at the end. This tip is attached to soft, flexible tubing terminating in a connector.



15 in (38.1 cm) overall length

Weighted tip

12010

20 Fr (6.7 mm), 1/4 in
(0.64 cm) connector

20 per carton

DLP™ Pericardial/Intracardiac Sump

This sump features a weighted perforated pool tip to minimize the possibility of tissue occlusion while maintaining an opportunity for drainage. The perforated tip is attached to tubing that terminates with a 1/4 inch (0.64 cm) connector.



15 in (38.1 cm) overall length

Weighted tip

12112

20 Fr (6.7 mm)

20 per carton


Important Safety Information

The Pericardial Sump (Model 12010) is not intended to be placed through a valve to drain a closed cardiac chamber. Care and caution should be taken to avoid damage to vessels and cardiac tissue during cannulation or other cardiac surgery procedures.

Suction circuit DLP™ Cardiac Sump



These sumps feature a perforated pool tip to minimize the possibility of tissue occlusion while maintaining an opportunity for drainage. These sumps are available in pediatric and adult sizes. The perforated tip is attached to tubing that terminates with a 1/4 inch (0.64 cm) connector.



12 in (30.5 cm) overall length

12013	12 Fr (4.0 mm) (10 holes in tip)
--------------	-------------------------------------

20 per carton

15 in (38.1 cm) overall length

12012	20 Fr (6.7 mm) (28 holes in tip)
--------------	-------------------------------------

20 per carton

Suction circuit

DLP™ Left Heart Vent Catheters



These vent catheters are for direct and indirect venting of the left ventricle and feature perforated tips. A choice of either PVC or silicone is available, along with straight, preformed, or malleable bodies with depth markings. These vents are available in pediatric and adult sizes. Straight body models come with a malleable or stiff guidewire style introducer for easy insertion and placement. All vent catheters terminate with a vented or non-vented 1/4 inch (0.64 cm) connector.



PVC 13 in (33.0 cm) overall length

Preformed 1.5 in (3.8 cm) tip and vented connector

12008*	10 Fr (3.3 mm) (8 holes in tip)
---------------	------------------------------------

20 per carton

Preformed 2 in (5.1 cm) tip and vented connector

12001	13 Fr (4.3 mm) (14 holes in tip)
--------------	-------------------------------------

20 per carton

2 in (5.1 cm) tip, vented connector, and malleable introducer

12101	13 Fr (4.3 mm) (14 holes in tip)
--------------	-------------------------------------

20 per carton



PVC 15 in (38.1 cm) overall length

Malleable body and vented connector

12110	10 Fr (3.3 mm) (9 holes in tip)
--------------	------------------------------------

12113	13 Fr (4.3 mm) (9 holes in tip)
--------------	------------------------------------

12115	15 Fr (5.0 mm) (20 holes in tip)
--------------	-------------------------------------

20 per carton



PVC 15 in (38.1 cm) overall length

Malleable body and non-vented connector

12116	16 Fr (5.3 mm) (11 holes in tip)
--------------	-------------------------------------

12118	18 Fr (6.0 mm) (11 holes in tip)
--------------	-------------------------------------

20 per carton



Silicone 16 in (40.6 cm) overall length

Non-vented connector and malleable introducer


12016	16 Fr (5.3 mm) (20 holes in tip)
--------------	-------------------------------------

12002	20 Fr (6.7 mm) (24 holes in tip)
--------------	-------------------------------------

Non-vented connector and preformed introducer

12220	20 Fr (6.7 mm) (24 holes in tip)
--------------	-------------------------------------

20 per carton



* Also available with Cortiva™ BioActive Surface.



Products

Surgical support



The largest cannula portfolio worldwide.

Our commitment
– **to you.**

Surgical support products are designed to provide timesaving techniques and ease of use – and to facilitate the remainder of the surgical procedure. We offer a wide range of products to ensure that your specific needs are met.

Surgical support



- ▶ DLP™ Disposable Pressure Display Sets 77
- ▶ DLP™ Tourniquet Sets 78
- ▶ TourniKwik™ Tourniquet Sets 79
- ▶ Sure-Snare™ Tourniquet Sets 79
- ▶ DLP™ Vessel Cannulae 80
- ▶ DLP™ Arteriotomy Cannulae 81
- ▶ DLP™ Arteriotomy Cannulae 81
- ▶ DLP™ Electro-Surgical Knife Cleaner 82
- ▶ DLP™ Phrenic Nerve Pads 82
- ▶ DLP™ Coronary Artery Retraction Clips 83
- ▶ DLP™ Pressure Catheter Placement Set 84
- ▶ DLP™ IV-decanter 85
- ▶ DLP™ Rapid Prime Set 86

Surgical support

DLP™ Disposable Pressure Display Sets




These display sets are for use with our Pressure Display Box. They feature a silicone diaphragm inside a plastic dome with a stopcock that attaches to the transducer fitting on the Pressure Display Box. This dome is attached to a coil of tubing that terminates with either a male luer or a stopcock for convenient handling.

45 in (1.14 m) overall length

61000	stopcock
62000	male luer

10 per carton



Surgical support DLP™ Tourniquet Sets




These tourniquet sets consist of color-coded tubes in varying lengths and color combinations. They are available for pediatric and adult uses. Wire snares are included with the tube sets.

Standard sets

Tube length of 4 in (10.2 cm) and 6 in (15.2 cm)

79016	12 Fr [3 - 4 in (10.2 cm) bronze, 3 - 6 in (15.2 cm) bronze, 2 - snares]
--------------	--


40 sets per carton



Tube length of 5.5 in (14.0 cm)

79003	12 Fr (1 - red, 2 - blue, 3 - clear, 1 - snare)
79015	12 Fr (1 - red, 2 - blue, 5 - clear, 2 - snares)

40 sets per carton



Tube length of 7 in (17.8 cm)

79004	12 Fr (1 - red, 1 - blue, 1 - snare)
79006	12 Fr (1 - red, 1 - blue, 3 - clear, 1 - snare)

40 sets per carton


79022	12 Fr (6 - small O.D. white radiopaque, 1 - snare)
79023	19 Fr (3 - large O.D. white radiopaque, 1 - snare)

20 sets per carton

79026	12 Fr (2 - red, 2 - blue, 2 - clear) and 16 Fr (2 - large O.D. clear, 1 - snare)
--------------	--

79027	12 Fr (2 - red, 1 - blue, 1 - snare)
--------------	--------------------------------------

40 sets per carton



Vena cava sets

Tube length of 7 in (17.8 cm)

79005	12 Fr (3 - clear, 1 - snare)
79009	12 Fr (4 - clear, 1 - snare)

40 sets per carton




Surgical support

TourniKwik™ Tourniquet Sets



These tourniquet sets consist of bronze-colored tubes in varying lengths that have a quick tie cap attached to one end. These caps eliminate the use of clamps by wrapping the suture string around the cap. Wire snares are included with each set of tubes.



Tube length of 6 in (15.2 cm)

79010	12 Fr (2 - tubes, 1 - snare)
79011	12 Fr (4 - tubes, 1 - snare)

40 sets per carton


Tube length of 7.5 in (19.1 cm)

79012	12 Fr (2 - tubes, 1 - snare)
79013	12 Fr (4 - tubes, 1 - snare)

40 sets per carton

Sure-Snare™ Tourniquet Sets

These tourniquet sets consist of bronze-colored tubes in varying lengths that have a plug attached to one end. These plugs eliminate the use of clamps by pushing the plug into the tube end where the suture string exits.



Tube length of 5 in (12.7 cm)

79020	22 Fr (2 - tubes, 2 - umbilical tapes, 2 - plastic snares)
--------------	--

20 sets per carton



Tube length of 7.5 in (19.1 cm)

79008	12 Fr (2 - tubes, 1 - wire snare)
--------------	-----------------------------------

20 sets per carton

Surgical support DLP™ Vessel Cannulae



These cannulae feature clear and radiopaque bodies and are available with or without a one-way valve. All the cannulae terminate with a female luer and have tips in various sizes and shapes.

2 in (5.1 cm) overall length

Clear body

30003	3 mm blunt tip
30004	2 mm blunt tip
30005	4 mm acorn tip
30007	3 mm beveled tip

40 per carton



2 in (5.1 cm) overall length

Clear body with radiopaque band

30022	3 mm blunt tip
--------------	----------------

20 per carton



2 in (5.1 cm) overall length

Radiopaque body

30006	3 mm beveled tip
30023	3 mm blunt tip

40 per carton



2.5 in (6.4 cm) overall length

Clear body with one-way valve

30000	3 mm beveled tip
30001	3 mm blunt tip

40 per carton



2.5 in (6.4 cm) overall length

Radiopaque body with one-way valve

30021	3 mm blunt tip
--------------	----------------

40 per carton



Surgical support

DLP™ Arteriotomy Cannulae



This cannula features a stainless steel shaft with a bulb-shaped tip terminating with a female luer.

1.8 in (4.6 cm) overall length

31001	1 mm tip
--------------	----------

20 per carton



These cannulae feature a polyurethane tube with a bulb-shaped radiopaque tip connected to a winged female luer.



1.3 in (3.3 cm) overall length

31002	2 mm tip
--------------	----------

31003	3 mm tip
--------------	----------

20 per carton



2.3 in (5.8 cm) overall length

31102	2 mm tip
--------------	----------

31103	3 mm tip
--------------	----------

31104	4 mm tip
--------------	----------

31105	5 mm tip
--------------	----------

31106	6 mm tip
--------------	----------

20 per carton

Surgical support

DLP™ Electro-Surgical Knife Cleaner



This device keeps electro-surgical knife blades clean with a quick swipe. It features a stainless steel spring attached to a plastic base, which has a self-adhesive pad. This allows the device to be placed at the most convenient site for the surgeon.

2.5 in (6.4 cm) overall length

60001

40 per carton




DLP™ Phrenic Nerve Pads

These pads are designed to help protect the phrenic nerve from the extreme cold during surgery. They are made of a closed-cell foam and contain slots. These slots can hold a pericardial sump in place. A radiopaque marker is tied to the pad.

38002 | medium

38003 | large

10 per carton



Surgical support

DLP™ Coronary Artery Retraction Clips




These clips are designed to help improve exposure to a coronary anastomosis site. Small prongs grab and gently hold tissue away from the vessel to improve vision.

16130	3.0 mm prong length
--------------	---------------------

16150	5.0 mm prong length
--------------	---------------------

10 per carton



Surgical support DLP™ Pressure Catheter Placement Set



These sets are an ideal option for measuring cardiac performance postoperatively. The catheters measure left atrial pressure in real time to accurately display critical preload pressures of the left heart. The catheter tubing has depth markings and is attached to a female luer. There is a self-adhesive anchoring collar for securing the catheter at the desired location.



24 in (61.0 cm) overall length

50010

3 Fr (1.0 mm) catheter with 2 - 17 ga peel-away needles

50011

5 Fr (1.7 mm) catheter with 2 - 14 ga peel-away needles

10 per carton




Catheter Models 50010 and 50011 are inserted through the chest wall and into the right superior pulmonary vein using two peel-away needles. The peel-away needles permit catheter insertion through the chest wall and pulmonary vein without utilizing an integral introducer.

Surgical support DLP™ IV-decanter



This decanter includes several features to facilitate aseptic decanting of fluids from intravenous containers. It includes a large-bore spike that is attached to flexible tubing. This decanter also features a clamp, thereby eliminating spillage, reducing contamination.



8.5 in (21.6 cm) overall length

40040 Flexible tubing and clamp

20 per carton

DLP™ Rapid Prime Set

This set facilitates the transfer of fluids during the priming of the circuit. It features a large-bore spike attached to flexible tubing with a clamp, and terminates with a male luer.



14 in (35.6 cm) overall length

With male luer

10021

10 per carton



Appendix

Product codes



Products codes Index



10001	DL P™ Straight Adapters, with vented 1/4 in connector.....	59	12008	DLP™ Left Heart Vent Catheters, 10 Fr (3.3 mm).....	74
10001K	DL P™ Straight Adapters, with vented 1/4 in connector and clamp.....	59	12010	DLP™ Pericardial Sumps, 1/4 in (0.64 cm) connector	72
10003	DL P™ "Y" Type Recirculating Adapters.....	60	12012	DL P™ Cardiac Sump, 20 Fr (6.7 mm)	73
10004	DL P™ "Y" Type Coronary Perfusion Adapters, 3.5 in (8.9 cm).....	60	12013	DL P™ Cardiac Sump, 12 Fr (4.0 mm).....	73
10004D	DLP™ "Y" Type Coronary Perfusion Adapters, 6 in (15.2 cm).....	60	12016	DLP™ Left Heart Vent Catheters, 16 Fr (5.3 mm).....	74
10005	DL P™ "Y" Adapters, 7.5 in (19.1 cm) overall length.....	59	12101	DLP™ Left Heart Vent Catheters, 13 Fr (4.3 mm).....	74
10005BH	DLP™ "Y" Type Adapters.....	61	12110	DLP™ Left Heart Vent Catheters, 10 Fr (3.3 mm).....	74
10005OS	DLP™ "Y" Type Adapters.....	61	12112	DLP™ Pericardial/Intracardiac Sumps, 20 Fr (6.7 mm).....	72
10007	DLP™ Perfusion Adapters, 1/4 in (0.64 cm) barbed connector.....	57	12113	DLP™ Left Heart Vent Catheters, 13 Fr (4.3 mm).....	74
10008	DLP™ Perfusion Adapters, 3/8 in (0.95 cm) barbed connector.....	57	12115	DLP™ Left Heart Vent Catheters, 15 Fr (5.0 mm).....	74
10009	DLP™ Aortic Root Cannulae, 9 ga (11 Fr).....	41	12116	DLP™ Left Heart Vent Catheters, 16 Fr (5.3 mm).....	74
10010	DLP™ Vein Graft Cannula.....	62	12118	DLP™ Left Heart Vent Catheters, 18 Fr (6.0 mm).....	74
10012	DLP™ Aortic Root Cannulae, 12 ga (9 Fr).....	41	12218	DLP™ Aortic Root Cannulae, 18 ga (4 Fr).....	41
10014	DLP™ Aortic Root Cannulae, 14 ga (7 Fr).....	41	12220	DLP™ Left Heart Vent Catheters, 20 Fr (6.7 mm).....	74
10016	DLP™ Aortic Root Cannulae, 16 ga (5 Fr).....	41	13000	ARISS™ Perfusion Switches.....	63
10018	DLP™ Aortic Root Cannulae, 18 ga (4 Fr).....	41	13001	DL P™ Antegrade/Retrograde Stopcock Adapter	64
10021	DLP™ Rapid Prime Set.....	85	13002	DLP™ Perfusion/Venting Adapter.....	57
10022	DLP™ Perfusion Adapters, 1/4 in (0.64 cm) barbed connector.....	57	13003	DLP™ Cardioplegia Management Set.....	64
10050	DLP™ Cardiac Suction Tubes, 6 Fr (2.0 mm).....	69	13004	ARISS™ Perfusion Switches.....	63
10051	DLP™ Cardiac Suction Tubes, 6 Fr (2.0 mm).....	69	14000	DLP™ Multiple Perfusion Sets	65
10052	DLP™ Cardiac Suction Tubes, 6 Fr (2.0 mm).....	69	14001	DLP™ Multiple Perfusion Sets	65
10053	DLP™ Cardiac Suction Tubes, 6 Fr (2.0 mm).....	69	14003	DLP™ Multiple Perfusion Sets	65
10060	DLP™ Cardiac Suction Tube, 10 Fr (3.3 mm).....	70	14006	DLP™ Multiple Perfusion Sets	66
10061	DLP™ Cardiac Suction Tube, Macro Sucker, 20 Fr (6.7 mm).....	71	14007	DLP™ Multiple Perfusion Sets	66
10062	DLP™ Cardiac Suction Tube, Mini Sucker, 9 Fr (3.0 mm).....	71	14008	DLP™ Multiple Perfusion Sets	66
10112	DLP™ Aortic Root Cannulae, 12 ga (9 Fr).....	41	14013	DLP™ Multiple Perfusion Sets	65
10114	DLP™ Aortic Root Cannulae, 14 ga (7 Fr).....	41	14017	DLP™ Multiple Perfusion Sets	65
10218	DLP™ Aortic Root Cannulae, 18 ga (4 Fr).....	41	14031	DLP™ Multiple Perfusion Sets	65
10313	DLP™ Cardioplegia Needles, 13 ga (8 Fr).....	43	15004	DLP™ Cardioplegia Adapter with Pressure Port	63
10534	DL P™ "Y" Adapters, 8.5 in (21.6 cm) overall length.....	59	16130	DLP™ Coronary Artery Retraction Clips, 3.0 mm prong length	83
10700	DLP™ "Y" Type Coronary Perfusion Adapters, 3.5 in (8.9 cm).....	60	16150	DLP™ Coronary Artery Retraction Clips, 5.0 mm prong length	83
10706	DLP™ "Y" Type Coronary Perfusion Adapters, 6 in (15.2 cm).....	60	20009	DLP™ Aortic Root Cannulae with Vent Line, 9 ga (11 Fr).....	42
10710	DLP™ "Y" Type Coronary Perfusion Adapters	60	20012	DLP™ Aortic Root Cannulae with Vent Line, 12 ga (9 Fr).....	42
11001	DLP™ Extension Line Adapters, 10 in (25.4 cm) length.....	58	20012S	DLP™ Aortic Root Cannulae with Vent Line, 12 ga (9 Fr).....	42
11001G	DLP™ Extension Line Adapters, 20 in (50.8 cm) length.....	58	20014	DLP™ Aortic Root Cannulae with Vent Line, 14 ga (7 Fr).....	42
11012	DLP™ Aortic Root Cannulae, 12 ga (9 Fr).....	41	20014L	DLP™ Aortic Root Cannulae with Vent Line, 14 ga (7 Fr).....	42
11012L	MiAR™ Aortic Root Cannula with Flow-Guard, 12 ga (9 Fr).....	39	20016	DLP™ Aortic Root Cannulae with Vent Line, 16 ga (5 Fr).....	42
11014	DLP™ Aortic Root Cannulae, 14 ga (7 Fr).....	41	20112	DLP™ Aortic Root Cannulae with Vent Line, 12 ga (9 Fr).....	42
11014L	MiAR™ Aortic Root Cannula with Flow-Guard, 14 ga (7 Fr).....	39	20114	DLP™ Aortic Root Cannulae with Vent Line, 14 ga (7 Fr).....	42
11316	DLP™ Cardioplegia Needles, 16 ga (5 Fr).....	43	20114HWF	DLP™ Aortic Root Cannulae with Vent Line, 14 ga (7 Fr).....	42
12001	DLP™ Left Heart Vent Catheters, 13 Fr (4.3 mm)	74	20114WF	DLP™ Aortic Root Cannulae with Vent Line, 14 ga (7 Fr).....	42
12002	DLP™ Left Heart Vent Catheters, 20 Fr (6.7 mm).....	74	21012	DLP™ Aortic Root Cannulae with Vent Line, 12 ga (9 Fr).....	42

Products codes Index



21014	DLP™ Aortic Root Cannulae with Vent Line, 14 ga (7 Fr).....	42	50011	DL P™ Pressure Catheter Placement Set, 5 Fr (1.7 mm).....	84
23009	DLP™ Aortic Root Cannulae, 9 ga (11 Fr).....	41	57414	DLP™ Femoral Arterial Cannulae, 14 Fr (4.7 mm).....	27
24009	DLP™ Aortic Root Cannulae with Vent Line, 9 ga (11 Fr).....	42	57417	DLP™ Femoral Arterial Cannulae, 17 Fr (5.7 mm).....	27
25009	DLP™ Pressure Monitoring Extension Line Adapters, male x female luers.....	58	57421	DLP™ Femoral Arterial Cannulae, 21 Fr (7.0 mm).....	27
25010	DLP™ Pressure Monitoring Extension Line Adapters, male x male luers.....	58	60001	DLP™ Electro-Surgical Knife Cleaner.....	82
30000	DLP™ Vessel Cannulae, 3 mm beveled tip.....	80	61000	DLP™ Disposable Pressure Display Sets.....	77
30001	DLP™ Vessel Cannulae, 3 mm blunt tip.....	80	62000	DLP™ Disposable Pressure Display Sets.....	77
30003	DLP™ Vessel Cannulae, 3 mm blunt tip.....	80	66112	DLP™ Single Stage Venous Cannulae, 12 Fr (4.0 mm).....	31
30004	DLP™ Vessel Cannulae, 2 mm blunt tip.....	80	66114	DLP™ Single Stage Venous Cannulae, 14 Fr (4.7 mm).....	31
30005	DLP™ Vessel Cannulae, 4 mm acorn tip.....	80	66116	DLP™ Single Stage Venous Cannulae, 16 Fr (5.3 mm).....	31
30006	DLP™ Vessel Cannulae, 3 mm beveled tip.....	80	66118	DLP™ Single Stage Venous Cannulae, 18 Fr (6.0 mm).....	31
30007	DLP™ Vessel Cannulae, 3 mm beveled tip.....	80	66120	DLP™ Single Stage Venous Cannulae, 20 Fr (6.7 mm).....	31
30010	DLP™ Coronary Ostial Perfusion Cannulae, 10 Fr (3.3 mm).....	45	66122	DLP™ Single Stage Venous Cannulae, 22 Fr (7.3 mm).....	31
30012	DLP™ Coronary Ostial Perfusion Cannulae, 12 Fr (4.0 mm).....	45	66124	DLP™ Single Stage Venous Cannulae, 24 Fr (8.0 mm).....	31
30014	DLP™ Coronary Ostial Perfusion Cannulae, 14 Fr (4.7 mm).....	45	66126	DLP™ Single Stage Venous Cannulae, 26 Fr (8.7 mm).....	31
30021	DLP™ Vessel Cannulae, 3 mm blunt tip.....	80	66128	DLP™ Single Stage Venous Cannulae, 28 Fr (9.3 mm).....	31
30022	DLP™ Vessel Cannulae, 3 mm blunt tip.....	80	66130	DLP™ Single Stage Venous Cannulae, 30 Fr (10.0 mm).....	31
30023	DLP™ Vessel Cannulae, 3 mm blunt tip.....	80	66132	DLP™ Single Stage Venous Cannulae, 32 Fr (10.7 mm).....	31
30050	DLP™ Coronary Ostial Perfusion Cannulae.....	45	66134	DLP™ Single Stage Venous Cannulae, 34 Fr (11.3 mm).....	31
30055	DLP™ Coronary Ostial Perfusion Cannulae.....	45	66136	DLP™ Single Stage Venous Cannulae, 36 Fr (12.0 mm).....	31
30110	DLP™ Coronary Ostial Perfusion Cannulae, 10 Fr (3.3 mm).....	46	66140	DLP™ Single Stage Venous Cannulae, 40 Fr (13.3 mm).....	31
30112	DLP™ Coronary Ostial Perfusion Cannulae, 12 Fr (4.0 mm).....	46	66236	DLP™ Single Stage Venous Cannulae, 36 Fr (12.0 mm).....	31
30114	DLP™ Coronary Ostial Perfusion Cannulae, 14 Fr (4.7 mm).....	46	67312	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 12 Fr (4.0 mm).....	30
30155	DLP™ Coronary Ostial Perfusion Cannulae, 10 Fr (3.3 mm).....	46	67314	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 14 Fr (4.7 mm).....	30
30212	DLP™ Coronary Ostial Perfusion Cannulae, 12 Fr (4.0 mm).....	46	67316	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 16 Fr (5.3 mm).....	30
30255	DLP™ Coronary Ostial Perfusion Cannulae, 10 Fr (3.3 mm).....	46	67318	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 18 Fr (6.0 mm).....	30
30315	DLP™ Coronary Ostial Perfusion Cannulae, 15 Fr (5.0 mm).....	44	67320	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 20 Fr (6.7 mm).....	30
30317	DLP™ Coronary Ostial Perfusion Cannulae, 17 Fr (5.7 mm).....	44	67512	DLP™ Single Stage Venous Cannulae with right angle, 12 Fr (4.0 mm).....	32
30320	DLP™ Coronary Ostial Perfusion Cannulae, 20 Fr (6.7 mm).....	44	67514	DLP™ Single Stage Venous Cannulae with right angle, 14 Fr (4.7 mm).....	32
30401	DLP™ Dual Lumen Cardioplegia Cannulae, 12 ga (9 Fr).....	40	67516	DLP™ Single Stage Venous Cannulae with right angle, 16 Fr (5.3 mm).....	32
31001	DLP™ Arteriotomy Cannulae, 1 mm tip.....	81	67518	DLP™ Single Stage Venous Cannulae with right angle, 18 Fr (6.0 mm).....	32
31002	DLP™ Arteriotomy Cannulae, 2 mm tip.....	81	67520	DLP™ Single Stage Venous Cannulae with right angle, 20 Fr (6.7 mm).....	32
31003	DLP™ Arteriotomy Cannulae, 3 mm tip.....	81	67522	DLP™ Single Stage Venous Cannulae with right angle, 22 Fr (7.3 mm).....	32
31102	DLP™ Arteriotomy Cannulae, 2 mm tip.....	81	67524	DLP™ Single Stage Venous Cannulae with right angle, 24 Fr (8.0 mm).....	32
31103	DLP™ Arteriotomy Cannulae, 3 mm tip.....	81	67526	DLP™ Single Stage Venous Cannulae with right angle, 26 Fr (8.7 mm).....	32
31104	DLP™ Arteriotomy Cannulae, 4 mm tip.....	81	67528	DLP™ Single Stage Venous Cannulae with right angle, 28 Fr (9.3 mm).....	32
31105	DLP™ Arteriotomy Cannulae, 5 mm tip.....	81	67530	DLP™ Single Stage Venous Cannulae with right angle, 30 Fr (10.0 mm).....	32
31106	DLP™ Arteriotomy Cannulae, 6 mm tip.....	81	67532	DLP™ Single Stage Venous Cannulae with right angle, 32 Fr (10.7 mm).....	32
38002	DLP™ Phrenic Nerve Pads, medium.....	82	67534	DLP™ Single Stage Venous Cannulae with right angle, 34 Fr (11.3 mm).....	32
38003	DLP™ Phrenic Nerve Pads, large.....	82	67636	DLP™ Single Stage Venous Cannulae with right angle, 36 Fr (12.0 mm).....	32
40040	DLP™ IV-decanter.....	85	67638	DLP™ Single Stage Venous Cannulae with right angle, 38 Fr (12.7 mm).....	32
50010	DL P™ Pressure Catheter Placement Set, 3 Fr (1.0 mm).....	84	68112	DLP™ Malleable Single Stage Venous Cannulae, 12 Fr (4.0 mm).....	33

Products codes Index



68114	DLP™ Malleable Single Stage Venous Cannulae, 14 Fr (4.7 mm).....	33	73420	Select Series™ Angled Tip Arterial Cannulae, 20 Fr (6.7 mm).....	13
68116	DLP™ Malleable Single Stage Venous Cannulae, 16 Fr (5.3 mm).....	33	73422	Select Series™ Angled Tip Arterial Cannulae, 22 Fr (7.3 mm).....	13
68118	DLP™ Malleable Single Stage Venous Cannulae, 18 Fr (6.0 mm).....	33	73424	Select Series™ Angled Tip Arterial Cannulae, 24 Fr (8.0 mm).....	13
68120	DLP™ Malleable Single Stage Venous Cannulae, 20 Fr (6.7 mm).....	33	75008	DLP™ Straight Tip Arterial Cannulae, 8 Fr (2.7 mm).....	10
68122	DLP™ Malleable Single Stage Venous Cannulae, 22 Fr (7.3 mm).....	33	75010	DLP™ Straight Tip Arterial Cannulae, 10 Fr (3.3 mm).....	10
68124	DLP™ Malleable Single Stage Venous Cannulae, 24 Fr (8.0 mm).....	33	75012	DLP™ Straight Tip Arterial Cannulae, 12 Fr (4.0 mm).....	10
68126	DLP™ Malleable Single Stage Venous Cannulae, 26 Fr (8.7 mm).....	33	75014	DLP™ Straight Tip Arterial Cannulae, 14 Fr (4.7 mm).....	10
68128	DLP™ Malleable Single Stage Venous Cannulae, 28 Fr (9.3 mm).....	33	75318	DLP™ Straight Tip Arterial Cannulae, 18 Fr (6.0 mm).....	10
68130	DLP™ Malleable Single Stage Venous Cannulae, 30 Fr (10.0 mm).....	33	75320	DLP™ Straight Tip Arterial Cannulae, 20 Fr (6.7 mm).....	10
68132	DLP™ Malleable Single Stage Venous Cannulae, 32 Fr (10.7 mm).....	33	75322	DLP™ Straight Tip Arterial Cannulae, 22 Fr (7.3 mm).....	10
68134	DLP™ Malleable Single Stage Venous Cannulae, 34 Fr (11.3 mm).....	33	75324	DLP™ Straight Tip Arterial Cannulae, 24 Fr (8.0 mm).....	10
68136	DLP™ Malleable Single Stage Venous Cannulae, 36 Fr (12.0 mm).....	33	76020	DLP™ Straight Tip Arterial Cannulae, 20 Fr (6.7 mm).....	10
68138	DLP™ Malleable Single Stage Venous Cannulae, 38 Fr (12.7 mm).....	33	76022	DLP™ Straight Tip Arterial Cannulae, 22 Fr (7.3 mm).....	10
69312	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 12 Fr (4.0 mm).....	30	76024	DLP™ Straight Tip Arterial Cannulae, 24 Fr (8.0 mm).....	10
69314	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 14 Fr (4.7 mm).....	30	76118	DLP™ Straight Tip Arterial Cannulae, 18 Fr (6.0 mm).....	10
69316	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 16 Fr (5.3 mm).....	30	76120	DLP™ Straight Tip Arterial Cannulae, 20 Fr (6.7 mm).....	10
69318	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 18 Fr (6.0 mm).....	30	76122	DLP™ Straight Tip Arterial Cannulae, 22 Fr (7.3 mm).....	10
69320	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 20 Fr (6.7 mm).....	30	76124	DLP™ Straight Tip Arterial Cannulae, 24 Fr (8.0 mm).....	10
69322	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 22 Fr (7.3 mm).....	30	77006	DLP™ One-Piece Pediatric Arterial Cannulae, 6 Fr (2.0 mm).....	9
69324	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 24 Fr (8.0 mm).....	30	77008	DLP™ One-Piece Pediatric Arterial Cannulae, 8 Fr (2.7 mm).....	9
69328	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 28 Fr (9.3 mm).....	30	77010	DLP™ One-Piece Pediatric Arterial Cannulae, 10 Fr (3.3 mm).....	9
69331	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 31 Fr (10.3 mm).....	30	77012	DLP™ One-Piece Pediatric Arterial Cannulae, 12 Fr (4.0 mm).....	9
69428	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 28 Fr (9.3 mm).....	30	77014	DLP™ One-Piece Pediatric Arterial Cannulae, 14 Fr (4.7 mm).....	9
69431	DLP™ Single Stage Venous Cannulae with Right Angle Metal Tip, 31 Fr (10.3 mm).....	30	77016	DLP™ One-Piece Pediatric Arterial Cannulae, 16 Fr (5.3 mm).....	9
70420	DLP™ Flexible Arch Arterial Cannulae, 20 Fr (6.7 mm).....	11	77106	DLP™ One-Piece Pediatric Arterial Cannulae, 6 Fr (2.0 mm).....	9
70422	DLP™ Flexible Arch Arterial Cannulae, 22 Fr (7.3 mm).....	11	77108	DLP™ One-Piece Pediatric Arterial Cannulae, 8 Fr (2.7 mm).....	9
70424	DLP™ Flexible Arch Arterial Cannulae, 24 Fr (8.0 mm).....	11	77110	DLP™ One-Piece Pediatric Arterial Cannulae, 10 Fr (3.3 mm).....	9
71420	DLP™ Flexible Arch Arterial Cannulae, 20 Fr (6.7 mm).....	11	77112	DLP™ One-Piece Pediatric Arterial Cannulae, 12 Fr (4.0 mm).....	9
71422	DLP™ Flexible Arch Arterial Cannulae, 22 Fr (7.3 mm).....	11	77114	DLP™ One-Piece Pediatric Arterial Cannulae, 14 Fr (4.7 mm).....	9
71424	DLP™ Flexible Arch Arterial Cannulae, 24 Fr (8.0 mm).....	11	77116	DLP™ One-Piece Pediatric Arterial Cannulae, 16 Fr (5.3 mm).....	9
72120	Select Series™ Straight Tip Arterial Cannulae, 20 Fr (6.7 mm).....	12	77206	DLP™ One-Piece Pediatric Arterial Cannulae, 6 Fr (2.0 mm).....	9
72122	Select Series™ Straight Tip Arterial Cannulae, 22 Fr (7.3 mm).....	12	77208	DLP™ One-Piece Pediatric Arterial Cannulae, 8 Fr (2.7 mm).....	9
72124	Select Series™ Straight Tip Arterial Cannulae, 24 Fr (8.0 mm).....	12	77210	DLP™ One-Piece Pediatric Arterial Cannulae, 10 Fr (3.3 mm).....	9
72220	Select Series™ Straight Tip Arterial Cannulae, 20 Fr (6.7 mm).....	12	77212	DLP™ One-Piece Pediatric Arterial Cannulae, 12 Fr (4.0 mm).....	9
72222	Select Series™ Straight Tip Arterial Cannulae, 22 Fr (7.3 mm).....	12	77306	DLP™ One-Piece Pediatric Arterial Cannulae, 6 Fr (2.0 mm).....	9
72224	Select Series™ Straight Tip Arterial Cannulae, 24 Fr (8.0 mm).....	12	77308	DLP™ One-Piece Pediatric Arterial Cannulae, 8 Fr (2.7 mm).....	9
72420	Select Series™ Angled Tip Arterial Cannulae, 20 Fr (6.7 mm).....	13	77310	DLP™ One-Piece Pediatric Arterial Cannulae, 10 Fr (3.3 mm).....	9
72422	Select Series™ Angled Tip Arterial Cannulae, 22 Fr (7.3 mm).....	13	77312	DLP™ One-Piece Pediatric Arterial Cannulae, 12 Fr (4.0 mm).....	9
72424	Select Series™ Angled Tip Arterial Cannulae, 24 Fr (8.0 mm).....	13	77418	EOPA™ Elongated One-Piece Arterial Cannulae, 18 Fr (6.0 mm).....	8
72520	Select Series™ Angled Tip Arterial Cannulae, 20 Fr (6.7 mm).....	13	77420	EOPA™ Elongated One-Piece Arterial Cannulae, 20 Fr (6.7 mm).....	8
72522	Select Series™ Angled Tip Arterial Cannulae, 22 Fr (7.3 mm).....	13	77422	EOPA™ Elongated One-Piece Arterial Cannulae, 22 Fr (7.3 mm).....	8
72524	Select Series™ Angled Tip Arterial Cannulae, 24 Fr (8.0 mm).....	13	77424	EOPA™ Elongated One-Piece Arterial Cannulae, 24 Fr (8.0 mm).....	8

Products codes Index



77518	EOPA™ Elongated One-Piece Arterial Cannulae, 18 Fr (6.0 mm).....	8	80122	DLP™ Curved Metal Tip Arterial Cannulae, 22 Fr (7.3 mm).....	18
77520	EOPA™ Elongated One-Piece Arterial Cannulae, 20 Fr (6.7 mm).....	8	80220	DLP™ Curved Metal Tip Arterial Cannulae, 20 Fr (6.7 mm).....	18
77522	EOPA™ Elongated One-Piece Arterial Cannulae, 22 Fr (7.3 mm).....	8	80222	DLP™ Curved Metal Tip Arterial Cannulae, 22 Fr (7.3 mm).....	18
77524	EOPA™ Elongated One-Piece Arterial Cannulae, 24 Fr (8.0 mm).....	8	80320	DLP™ Curved Metal Tip Arterial Cannulae, 20 Fr (6.7 mm).....	18
77618	EOPA™ Elongated One-Piece Arterial Cannulae, 18 Fr (6.0 mm).....	8	80322	DLP™ Curved Metal Tip Arterial Cannulae, 22 Fr (7.3 mm).....	18
77620	EOPA™ Elongated One-Piece Arterial Cannulae, 20 Fr (6.7 mm).....	8	81020	DLP™ Curved Tip Arterial Cannulae, 20 Fr (6.7 mm).....	17
77622	EOPA™ Elongated One-Piece Arterial Cannulae, 22 Fr (7.3 mm).....	8	81022	DLP™ Curved Tip Arterial Cannulae, 22 Fr (7.3 mm).....	17
77624	EOPA™ Elongated One-Piece Arterial Cannulae, 24 Fr (8.0 mm).....	8	81120	DLP™ Curved Tip Arterial Cannulae, 20 Fr (6.7 mm).....	17
77718	EOPA™ Elongated One-Piece Arterial Cannulae, 18 Fr (6.0 mm).....	8	81122	DLP™ Curved Tip Arterial Cannulae, 22 Fr (7.3 mm).....	17
77720	EOPA™ Elongated One-Piece Arterial Cannulae, 20 Fr (6.7 mm).....	8	82020	DLP™ Curved Tip Arterial Cannulae, 20 Fr (6.7 mm).....	17
77722	EOPA™ Elongated One-Piece Arterial Cannulae, 22 Fr (7.3 mm).....	8	82022	DLP™ Curved Tip Arterial Cannulae, 22 Fr (7.3 mm).....	17
77724	EOPA™ Elongated One-Piece Arterial Cannulae, 24 Fr (8.0 mm).....	8	82120	DLP™ Curved Tip Arterial Cannulae, 20 Fr (6.7 mm).....	17
78220	EOPA 3D™ Arterial Cannulae, 20 Fr (6.7 mm).....	6	82122	DLP™ Curved Tip Arterial Cannulae, 22 Fr (7.3 mm).....	17
78222	EOPA 3D™ Arterial Cannulae, 22 Fr (7.3 mm).....	6	83020	DLP™ Curved Tip Arterial Cannulae, 20 Fr (6.7 mm).....	16
78320	EOPA 3D™ Arterial Cannulae, 20 Fr (6.7 mm).....	6	83022	DLP™ Curved Tip Arterial Cannulae, 22 Fr (7.3 mm).....	16
78322	EOPA 3D™ Arterial Cannulae, 22 Fr (7.3 mm).....	6	83024	DLP™ Curved Tip Arterial Cannulae, 24 Fr (8.0 mm).....	16
78420	Select 3D™ II Arterial Cannulae, 20 Fr (6.7 mm).....	7	84020	DLP™ Curved Tip Arterial Cannulae, 20 Fr (6.7 mm).....	16
78422	Select 3D™ II Arterial Cannulae, 22 Fr (7.3 mm).....	7	84022	DLP™ Curved Tip Arterial Cannulae, 22 Fr (7.3 mm).....	16
78424	Select 3D™ II Arterial Cannulae, 24 Fr (8.0 mm).....	7	84024	DLP™ Curved Tip Arterial Cannulae, 24 Fr (8.0 mm).....	16
78520	Select 3D™ II Arterial Cannulae, 20 Fr (6.7 mm).....	7	87024	DLP™ Curved Tip Arterial Cannulae with Wirewound Body, 24 Fr (8.0 mm).....	15
78522	Select 3D™ II Arterial Cannulae, 22 Fr (7.3 mm).....	7	87120	DLP™ Curved Tip Arterial Cannulae with One-piece Wirewound Body, 20 Fr (6.7 mm).....	14
78524	Select 3D™ II Arterial Cannulae, 24 Fr (8.0 mm).....	7	87122	DLP™ Curved Tip Arterial Cannulae with One-piece Wirewound Body, 22 Fr (7.3 mm).....	14
79003	DLP™ Tourniquet Sets.....	78	87124	DLP™ Curved Tip Arterial Cannulae with One-piece Wirewound Body, 24 Fr (8.0 mm).....	14
79004	DLP™ Tourniquet Sets.....	78	87224	DLP™ Curved Tip Arterial Cannulae with One-piece Wirewound Body, 24 Fr (8.0 mm).....	14
79005	DLP™ Tourniquet Sets.....	78	87924	DLP™ Curved Tip Arterial Cannulae with Wirewound Body 24 Fr (8.0 mm).....	15
79006	DLP™ Tourniquet Sets.....	78	88020	DLP™ Curved Tip Arterial Cannulae 20 Fr (6.7 mm).....	16
79008	Sure-Snare™ Tourniquet Sets.....	79	88022	DLP™ Curved Tip Arterial Cannulae 22 Fr (7.3 mm).....	16
79009	DLP™ Tourniquet Sets.....	78	88024	DLP™ Curved Tip Arterial Cannulae 24 Fr (8.0 mm).....	16
79010	TourniKwik™ Tourniquet Sets.....	79	88120	DLP™ Curved Tip Arterial Cannulae 20 Fr (6.7 mm).....	16
79011	TourniKwik™ Tourniquet Sets.....	79	88122	DLP™ Curved Tip Arterial Cannulae 22 Fr (7.3 mm).....	16
79012	TourniKwik™ Tourniquet Sets.....	79	88124	DLP™ Curved Tip Arterial Cannulae 24 Fr (8.0 mm).....	16
79013	TourniKwik™ Tourniquet Sets.....	79	89020	DLP™ Curved Tip Arterial Cannulae 20 Fr (6.7 mm).....	16
79015	DLP™ Tourniquet Sets.....	78	89022	DLP™ Curved Tip Arterial Cannulae 22 Fr (7.3 mm).....	16
79016	DLP™ Tourniquet Sets.....	78	89024	DLP™ Curved Tip Arterial Cannulae 24 Fr (8.0 mm).....	16
79020	Sure-Snare™ Tourniquet Sets.....	79	89120	DLP™ Curved Tip Arterial Cannulae 20 Fr (6.7 mm).....	16
79022	DLP™ Tourniquet Sets.....	78	89122	DLP™ Curved Tip Arterial Cannulae 22 Fr (7.3 mm).....	16
79023	DLP™ Tourniquet Sets.....	78	89124	DLP™ Curved Tip Arterial Cannulae 24 Fr (8.0 mm).....	16
79026	DLP™ Tourniquet Sets.....	78	91228	MC2™ Two Stage Venous Cannulae, 28/36 Fr (9.3/12.0 mm).....	35
79027	DLP™ Tourniquet Sets.....	78	91228C	MC2™ Two Stage Venous Cannulae, 28/36 Fr (9.3/12.0 mm).....	35
80020	DLP™ Curved Metal Tip Arterial Cannulae, 20 Fr (6.7 mm).....	18	91229	MC2™ Two Stage Venous Cannulae, 29/37 Fr (9.7/12.3 mm).....	35
80022	DLP™ Curved Metal Tip Arterial Cannulae, 22 Fr (7.3 mm).....	18	91229C	MC2™ Two Stage Venous Cannulae, 29/37 Fr (9.7/12.3 mm).....	35
80120	DLP™ Curved Metal Tip Arterial Cannulae, 20 Fr (6.7 mm).....	18	91236	MC2™ Two Stage Venous Cannulae, 36/46 Fr (12.0/15.3 mm).....	35

Products codes Index



91236C	MC2™ Two Stage Venous Cannulae, 36/46 Fr (12.0/15.3 mm).....	35	94535	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm)	53
91240	MC2™ Two Stage Venous Cannulae, 32/40 Fr (10.7/13.3 mm).....	35	94615	Gundry™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm)	49
91240C	MC2™ Two Stage Venous Cannulae, 32/40 Fr (10.7/13.3 mm).....	35	94625	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm)	51
91246	MC2™ Two Stage Venous Cannulae, 34/46 Fr (11.3/15.3 mm).....	35	94665	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm)	51
91246C	MC2™ Two Stage Venous Cannulae, 34/46 Fr (11.3/15.3 mm).....	35	94715	Gundry™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm)	49
91251	MC2™ Two Stage Venous Cannulae, 36/51 Fr (12.0/17.0 mm).....	35	94725	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm)	50
91251C	MC2™ Two Stage Venous Cannulae, 36/51 Fr (12.0/17.0 mm).....	35	94735	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm)	52
91263	MC2™ Two Stage Venous Cannulae, 32/40 Fr (10.7/13.3 mm).....	35	94745	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm).....	52
91263C	MC2™ Two Stage Venous Cannulae, 32/40 Fr (10.7/13.3 mm).....	35	94835	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm).....	54
91265	MC2™ Two Stage Venous Cannulae, 36/46 Fr (12.0/15.3 mm).....	35	94885	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm).....	54
91265C	MC2™ Two Stage Venous Cannulae, 36/46 Fr (12.0/15.3 mm).....	35	94885K	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm).....	53
91329	MC2™ Two Stage Venous Cannulae, 29/29 Fr (9.7/9.7 mm).....	35	94895	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm).....	54
91329C	MC2™ Two Stage Venous Cannulae, 29/29 Fr (9.7/9.7 mm).....	35	94913	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 13 Fr (4.3 mm).....	51
91429	MC2X™ Three Stage Venous Cannulae, 29/29/29 Fr (9.7/9.7/9.7 mm).....	36	94913L	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 13 Fr (4.3 mm).....	51
91429C	MC2X™ Three Stage Venous Cannulae, 29/29/29 Fr (9.7/9.7/9.7 mm).....	36	94915	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm).....	51
91437	MC2X™ Three Stage Venous Cannulae, 29/46/37 Fr (9.7/15.3/12.3 mm).....	36	94935	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm).....	54
91437C	MC2X™ Three Stage Venous Cannulae, 29/46/37 Fr (9.7/15.3/12.3 mm).....	36	94965	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm).....	51
93438	VC2™ Venous Cannulae, 34/38 Fr (11.3/12.7 mm).....	34	94975	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm).....	51
93438C	VC2™ Venous Cannulae, 34/38 Fr (11.3/12.7 mm).....	34	94995	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm).....	52
93448	VC2™ Venous Cannulae, 34/48 Fr (11.3/16.0 mm).....	34	96328	DLP™ Femoral Venous Cannulae, 28 Fr (9.3 mm).....	26
93448C	VC2™ Venous Cannulae, 34/48 Fr (11.3/16.0 mm).....	34	96428	DLP™ Femoral Venous Cannulae, 28 Fr (9.3 mm).....	26
93463	VC2™ Venous Cannulae, 34/38 Fr (11.3/12.7 mm).....	34	96530-115	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannula Kits, 15 Fr (5.0 mm).....	22
93463C	VC2™ Venous Cannulae, 34/38 Fr (11.3/12.7 mm).....	34	96530-117	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannula Kits, 17 Fr (5.7 mm).....	22
93464	VC2™ Venous Cannulae, 34/48 Fr (11.3/16.0 mm).....	34			
93464C	VC2™ Venous Cannulae, 34/48 Fr (11.3/16.0 mm).....	34			
94006	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 6 Fr (2.0 mm).....	50			
94106	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 6 Fr (2.0 mm).....	50			
94110	Gundry™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 10 Fr (3.3 mm).....	49			
94113T	Gundry™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 13 Fr (4.3 mm).....	49			
94115T	Gundry™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm).....	49			
94215T	DLP™ Retrograde Coronary Sinus Perfusion Cannulae with Manual-Inflate Cuff, 15 Fr (5.0 mm).....	50			
94315T	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm).....	52			
94415T	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 15 Fr (5.0 mm).....	52			
94533	DL P™ Retrograde Coronary Sinus Perfusion Cannulae with Auto-Inflate Cuff, 13 Fr (4.3 mm)	53			

Products codes Index



96530-119	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannula Kits, 19 Fr (6.3 mm).....	22	96820-108	Bio-Medicus™ NextGen Pediatric Arterial Cannulae, 8 Fr (2.7 mm).....	21
96530-121	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannula Kits, 21 Fr (7.0 mm).....	22	96820-110	Bio-Medicus™ NextGen Pediatric Arterial Cannulae, 10 Fr (3.3 mm).....	21
96530-123	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannula Kits, 23 Fr (7.7 mm).....	22	96820-112	Bio-Medicus™ NextGen Pediatric Arterial Cannulae, 12 Fr (4.0 mm).....	21
96530-125	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannula Kits, 25 Fr (8.3 mm).....	22	96820-114	Bio-Medicus™ NextGen Pediatric Arterial Cannulae, 14 Fr (4.7 mm).....	21
96551	Bio-Medicus™ NextGen Venous Insertion Kit.....	25	96830-108	Bio-Medicus™ NextGen Pediatric Venous Cannulae, 8 Fr (2.7 mm).....	21
96552	Bio-Medicus™ NextGen Arterial/Jugular Venous Insertion Kit.....	25	96830-110	Bio-Medicus™ NextGen Pediatric Venous Cannulae, 10 Fr (3.3 mm).....	21
96553	Bio-Medicus™ NextGen Pediatric Arterial/Venous Insertion Kit.....	25	96830-112	Bio-Medicus™ NextGen Pediatric Venous Cannulae, 12 Fr (4.0 mm).....	21
96570-115	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannulae, 15 Fr (5.0 mm).....	22	96830-114	Bio-Medicus™ NextGen Pediatric Venous Cannulae, 14 Fr (4.7 mm).....	21
96570-117	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannulae, 17 Fr (5.7 mm).....	22	96880-019	Bio-Medicus™ Percutaneous Kit, Multi-Stage Venous Femoral, 19 Fr (6.3 mm).....	23
96570-119	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannulae, 19 Fr (6.3 mm).....	22	96880-021	Bio-Medicus™ Percutaneous Kit, Multi-Stage Venous Femoral, 21 Fr (7.0 mm).....	23
96570-121	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannulae, 21 Fr (7.0 mm).....	22	96880-025	Bio-Medicus™ Percutaneous Kit, Multi-Stage Venous Femoral, 25 Fr (8.3 mm).....	23
96570-123	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannulae, 23 Fr (7.7 mm).....	22			
96570-125	Bio-Medicus™ NextGen Femoral Arterial or Jugular Venous Cannulae, 25 Fr (8.3 mm).....	22			
96600-115	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannula Kits, 15 Fr (5.0 mm).....	24			
96600-117	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannula Kits, 17 Fr (5.7 mm).....	24			
96600-119	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannula Kits, 19 Fr (6.3 mm).....	24			
96600-121	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannula Kits, 21 Fr (7.0 mm).....	24			
96600-123	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannula Kits, 23 Fr (7.7 mm).....	24			
96600-125	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannula Kits, 25 Fr (8.3 mm).....	24			
96600-127	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannula Kits, 27 Fr (9.0 mm).....	24			
96600-129	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannula Kits, 29 Fr (9.7 mm).....	24			
96670-115	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannulae, 15 Fr (5.0 mm).....	24			
96670-117	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannulae, 17 Fr (5.7 mm).....	24			
96670-119	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannulae, 19 Fr (6.3 mm).....	24			
96670-121	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannulae, 21 Fr (7.0 mm).....	24			
96670-123†	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannulae, 23 Fr (7.7 mm).....	24			
96670-125†	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannulae, 25 Fr (8.3 mm).....	24			
96670-127†	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannulae, 27 Fr (9.0 mm).....	24			
96670-129†	Bio-Medicus™ NextGen Femoral Bi-caval Venous Cannulae, 29 Fr (9.7 mm).....	24			

Products codes

List of products available with Cortiva™ BioActive Surface



Extracorporeal Circuit - Arterial Cannulae

CB71424	Flexible Arch Arterial, 11" (28cm), 24 Fr. (8.0mm), 3/8" Non-Vented Connector with Cortiva™ BioActive Surface
CB72124	Arterial Cannula, Wirewound, Straight, 24 Fr., Vented connector with Cortiva™ BioActive Surface
CB72224	Arterial Cannula, Wirewound, Straight, Non-Vented Connector, 24 Fr. with Cortiva™ BioActive Surface
CB75318	Arterial Cannula, Adult Straight Tip with 3/8" conn. Site, 18 Fr. (6.0mm) with Cortiva™ BioActive Surface
CB75320	Arterial Cannula, Adult Straight Tip with 3/8" conn. Site, 20 Fr. (6.7mm) with Cortiva™ BioActive Surface
CB77006	Arterial, One Piece Pediatric 9" (23cm), 1/4" Vented Connector, 6 Fr. (2.0mm) with Cortiva™ BioActive Surface
CB77008	Arterial, One Piece Pediatric 9" (23cm), 1/4" Vented Connector, 8 Fr. (2.7mm) with Cortiva™ BioActive Surface
CB77010	Arterial, One Piece Pediatric 9" (23cm), 1/4" Vented Connector, 10 Fr. (3.3mm) with Cortiva™ BioActive Surface
CB77012	Arterial, One Piece Pediatric 9" (23cm), 1/4" Vented Connector, 12 Fr. (4.0mm) with Cortiva™ BioActive Surface
CB77014	Arterial, One Piece Pediatric 9" (23cm), 1/4" Vented Connector, 14 Fr. (4.7mm) with Cortiva™ BioActive Surface
CB77016	Arterial, One Piece Pediatric 9" (23cm), 1/4" Vented Connector, 16 Fr. (5.3mm) with Cortiva™ BioActive Surface
CB77420	Arterial Elongated One Piece 12" (30.5 cm), Blunt Tip introducer w/o guidewire, 20 Fr. (6.7 mm), vented with Cortiva™ BioActive Surface
CB77422	Arterial Elongated One Piece 12" (30.5 cm), Blunt Tip introducer w/o guidewire, 22 Fr. (7.3 mm), vented with Cortiva™ BioActive Surface
CB77520	Arterial Elongated One Piece 12" (30.5 cm), Blunt Tip introducer w/o guidewire, 20 Fr. (6.7 mm), non-vented with Cortiva™ BioActive Surface
CB77522	Arterial Elongated One Piece 12" (30.5 cm), Blunt Tip introducer w/o guidewire, 22 Fr. (7.3 mm), non-vented with Cortiva™ BioActive Surface
CB77618	Arterial Elongated One Piece 12" (30.5 cm), Dilator Tip introducer w/guidewire, 18 Fr. (6.0 mm), vented with Cortiva™ BioActive Surface
CB77620	Arterial Elongated One Piece 12" (30.5 cm), Dilator Tip introducer w/guidewire, 20 Fr. (6.7 mm), vented with Cortiva™ BioActive Surface
CB77622	Arterial Elongated One Piece 12" (30.5 cm), Dilator Tip introducer w/guidewire, 22 Fr. (7.3 mm), vented with Cortiva™ BioActive Surface
CB77720	Arterial Elongated One Piece 12" (30.5 cm), Dilator Tip introducer w/guidewire, 20 Fr. (6.7 mm), non-vented with Cortiva™ BioActive Surface
CB77722	Arterial Elongated One Piece 12" (30.5 cm), Dilator Tip introducer w/guidewire, 22 Fr. (7.3 mm), non-vented with Cortiva™ BioActive Surface
CB78222	EOPA 3D 22FR VENTED with Cortiva™ BioActive Surface
CB78322	EOPA 3D 22FR NON-VENTED with Cortiva™ BioActive Surface

CB80120	Arterial, Metal Tip, 8.5" (21.6cm), Soft Flange, 20 Fr. (6.7mm), 3/8" Non-Vent Conn with Cortiva™ BioActive Surface
CB81120	Arterial Cannula, Ultraflex Curved Tip, 10" (25cm) Soft Body, 3/8" Non-Vented Connector, 20 Fr. with Cortiva™ BioActive Surface
CB81122	Arterial Cannula, Ultraflex Curved Tip, 10" (25cm) Soft Body, 3/8" Non-Vented Connector, 22 Fr. with Cortiva™ BioActive Surface
CB87222	WireWound curved Tip Arterial Cannula, 22 Fr., 3/8" Vented Connector with Cortiva™ BioActive Surface

Extracorporeal Circuit - Femoral Cannulae

CB96345-023	Adult Venous Femoral Cannula with 1/2" Non-Vented Connector, 30" (76.2 cm), 23 Fr. with Cortiva™ BioActive Surface
CB96345-025	Adult Venous Femoral Cannula with 1/2" Non-Vented Connector, 30" (76.2 cm), 25 Fr. with Cortiva™ BioActive Surface
CB96345-027	Adult Venous Femoral Cannula with 1/2" Non-Vented Connector, 30" (76.2 cm), 27 Fr. with Cortiva™ BioActive Surface
CB96345-029	Adult Venous Femoral Cannula with 1/2" Non-Vented Connector, 30" (76.2 cm), 29 Fr. with Cortiva™ BioActive Surface
CB96570-015	Femoral Arterial Cannula w/introducer 3/8" Connector with Lock-Vent, 17" (43.2 cm), 15 Fr. with Cortiva™ BioActive Surface
CB96570-017	Femoral Arterial Cannula w/introducer 3/8" Connector with Lock-Vent, 17" (43.2 cm), 17 Fr. with Cortiva™ BioActive Surface
CB96570-019	Femoral Arterial Cannula w/introducer 3/8" Connector with Lock-Vent, 17" (43.2 cm), 19 Fr. with Cortiva™ BioActive Surface
CB96570-021	Femoral Arterial Cannula w/introducer 3/8" Connector with Lock-Vent, 17" (43.2 cm), 21 Fr. with Cortiva™ BioActive Surface
CB96670-015	Femoral Venous Cannula w/introducer 3/8" Non-Vented Connector, 30" (76.2 cm), 15 Fr. with Cortiva™ BioActive Surface
CB96670-017	Femoral Venous Cannula w/introducer 3/8" Non-Vented Connector, 30" (76.2 cm), 17 Fr. with Cortiva™ BioActive Surface
CB96670-019	Femoral Venous Cannula w/introducer 3/8" Non-Vented Connector, 30" (76.2 cm), 19 Fr. with Cortiva™ BioActive Surface
CB96670-021	Femoral Venous Cannula w/introducer 3/8" Non-Vented Connector, 30" (76.2 cm), 21 Fr. with Cortiva™ BioActive Surface
CB96825-008	Pediatric Arterial Cannula with 1/4" Non-Vented Connector, 0.38mm, 7.5" (19 cm), 8 Fr. with Cortiva™ BioActive Surface
CB96825-010	Pediatric Arterial Cannula with 1/4" Non-Vented Connector, 0.38mm, 7.5" (19 cm), 10 Fr. with Cortiva™ BioActive Surface
CB96825-012	Pediatric Arterial Cannula with 1/4" Non-Vented Connector, 0.38mm, 7.5" (19 cm), 12 Fr. with Cortiva™ BioActive Surface
CB96825-014	Pediatric Arterial Cannula with 1/4" Non-Vented Connector, 0.38mm, 7.5" (19 cm), 14 Fr. with Cortiva™ BioActive Surface
CB96835-008	Pediatric Venous Cannula with 1/4" Non-Vented Connector, 0.38mm, 8 Fr. 7.5" (19 cm) with Cortiva™ BioActive Surface

Products codes

List of products available with Cortiva™ BioActive Surface



- CB96835-010 Pediatric Venous Cannula with 1/4" Non-Vented Connector, 0.38mm, 10 Fr. 7.5" (19 cm) with Cortiva™ BioActive Surface
- CB96835-012 Pediatric Venous Cannula with 1/4" Non-Vented Connector, 0.38mm, 12 Fr. 7.5" (19 cm) with Cortiva™ BioActive Surface
- CB96835-014 Pediatric Venous Cannula with 1/4" Non-Vented Connector, 0.38mm, 14 Fr. 7.5" (19 cm) with Cortiva™ BioActive Surface

Extracorporeal Circuit - Venous Cannulae

- CB66112 Single Stage Venous Cannula Straight Tip, 12 Fr. (4.0mm), 1/4" Conn Site with Cortiva™ BioActive Surface
- CB66114 Single Stage Venous Cannula Straight Tip, 14 Fr. (4.7mm), 1/4" Conn Site with Cortiva™ BioActive Surface
- CB66116 Single Stage Venous Cannula Straight Tip, 16 Fr. (5.3mm), 1/4" Conn Site with Cortiva™ BioActive Surface
- CB66118 Single Stage Venous Cannula Straight Tip, 18 Fr. (6.0mm), 1/4" - 3/8" Conn Site with Cortiva™ BioActive Surface
- CB66120 Single Stage Venous Cannula Straight Tip, 20 Fr. (6.7mm), 1/4" - 3/8" Conn Site with Cortiva™ BioActive Surface
- CB66124 Single Stage Venous Cannula Straight Tip, 24 Fr. (8.0mm), 1/4" - 3/8" Conn Site with Cortiva™ BioActive Surface
- CB66128 Single Stage Venous Cannula Straight Tip, 28 Fr. (9.3mm), 3/8" Conn Site with Cortiva™ BioActive Surface
- CB66132 Single Stage Venous Cannula Straight Tip, 32 Fr. (10.7mm), 3/8" Conn Site with Cortiva™ BioActive Surface
- CB66236 Single Stage Venous Cannula Straight Tip, 36 Fr. (12mm), 3/8" Conn Site with Cortiva™ BioActive Surface
- CB67312 Single Stage Venous Cannula, Right Angle, Metal Tip, Pediatric Size, 12 Fr., 1/4" conn. Site with Cortiva™ BioActive Surface
- CB67314 Single Stage Venous Cannula, Right Angle, Metal Tip, Pediatric Size, 14 Fr., 1/4" conn. Site with Cortiva™ BioActive Surface
- CB67316 Single Stage Venous Cannula, Right Angle, Metal Tip, Pediatric Size, 16 Fr., 1/4" conn. Site with Cortiva™ BioActive Surface
- CB67318 Single Stage Venous Cannula, Right Angle, Metal Tip, Pediatric Size, 18 Fr., 1/4" conn. Site with Cortiva™ BioActive Surface
- CB67320 Single Stage Venous Cannula, Right Angle, Metal Tip, Pediatric Size, 20 Fr., 1/4" conn. Site with Cortiva™ BioActive Surface
- CB67520 Single Stage Venous Cann., Angled Tip, 20 Fr. (6.7mm), 1/4" - 3/8" Conn Site with Cortiva™ BioActive Surface
- CB67524 Single Stage Venous Cann., Angled Tip, 24 Fr. (8.0mm), 1/4" - 3/8" Conn Site with Cortiva™ BioActive Surface
- CB67528 Single Stage Venous Cannula, Angled Tip, 28 Fr. (9.3mm), 3/8" Conn Site with Cortiva™ BioActive Surface
- CB67530 Single Stage Venous Cann., Angled Tip, 30 Fr. (10.0mm), 3/8" Conn Site with Cortiva™ BioActive Surface
- CB67532 Single Stage Venous Cannula, Angled Tip, 32 Fr. (10.7mm), 3/8" Conn Site with Cortiva™ BioActive Surface
- CB67534 Single Stage Venous Cann., Angled Tip, 34 Fr. (11.3mm), 3/8" Conn Site with Cortiva™ BioActive Surface
- CB67636 Single Stage Venous Cannula, Angled Tip, 36 Fr. (12mm), 3/8" conn. Site with Cortiva™ BioActive Surface
- CB68112 Malleable Single Stage Venous Cannula, 12 Fr. with Cortiva™ BioActive Surface
- CB68114 Malleable Single Stage Venous Cannula, 14 Fr. with Cortiva™ BioActive Surface
- CB68116 Malleable Single Stage Venous Cannula, 16 Fr. with Cortiva™ BioActive Surface
- CB68118 Malleable Single Stage Venous Cannula, 18 Fr. with Cortiva™ BioActive Surface
- CB68120 Malleable Single Stage Venous Cannula, 20 Fr. with Cortiva™ BioActive Surface
- CB68124 Malleable Single Stage Venous Cannula, 24 Fr. with Cortiva™ BioActive Surface
- CB68128 Malleable Single Stage Venous Cannula, 28 Fr. with Cortiva™ BioActive Surface
- CB68132 Malleable Single Stage Venous Cannula, 32 Fr. with Cortiva™ BioActive Surface
- CB68136 Malleable Single Stage Venous Cannula, 36Fr. with Cortiva™ BioActive Surface
- CB69320 Single Stage Venous Cannula, Right Angle, Metal Tip, High Flow, 3/8" conn site, 20 Fr. with Cortiva™ BioActive Surface
- CB69324 Single Stage Venous Cannula, Right Angle, Metal Tip, High Flow, 3/8" conn site, 24 Fr. with Cortiva™ BioActive Surface
- CB69328 Single Stage Venous Cannula, Right Angle, Metal Tip, High Flow, 3/8" conn site, 28 Fr. with Cortiva™ BioActive Surface
- CB91228 MC2 Two Stage Venous Cannula 28/36 Fr. with 3/8" conn. Site with Cortiva™ BioActive Surface
- CB91228C MC2 Two Stage Venous Cannula 28/36 Fr. with 1/2" Connector with Cortiva™ BioActive Surface
- CB91236C MC2 Two Stage Venous Cannula 36/46 Fr. with 1/2" Connector with Cortiva™ BioActive Surface
- CB91240 MC2 Two Stage Venous Cannula 32/40 Fr. with 1/2" conn. Site with Cortiva™ BioActive Surface
- CB91240C MC2 Two Stage Venous Cannula 32/40 Fr. with 1/2" Connector with Cortiva™ BioActive Surface
- CB91246 MC2 Two Stage Venous Cannula 34/46 Fr. with 1/2" conn. Site with Cortiva™ BioActive Surface
- CB91246C MC2 Two Stage Venous Cannula 34/46 Fr. with 1/2" Connector with Cortiva™ BioActive Surface
- CB91263 Oval MC2 Venous Cannula 32/40 Fr. with 1/2" conn. Site with Cortiva™ BioActive Surface
- CB91263C Oval MC2 Venous Cannula 32/40 Fr. with 1/2" Connector with Cortiva™ BioActive Surface
- CB91265 Oval MC2 Venous Cannula 36/46 Fr. with 1/2" conn. Site with Cortiva™ BioActive Surface
- CB91265C Oval MC2 Venous Cannula 36/46 Fr. with 1/2" Connector with Cortiva™ BioActive Surface

Products codes

List of products available with Cortiva™ BioActive Surface



- CB91329 Two Stage Venous Cannula 29/29 Fr. with 3/8" conn. Site with Cortiva™ BioActive Surface
- CB91329C Two Stage Venous Cannula 29/29 Fr. with 3/8" Connector with Cortiva™ BioActive Surface
- CB91429 Three Stage Venous Cannula 29/ 29/ 29 Fr. with 3/8" conn. Site with Cortiva™ BioActive Surface
- CB91429C Three Stage Venous Cannula 29/ 29/ 29 Fr. with 3/8" barbed Connector with Cortiva™ BioActive Surface
- CB93438C VC2 JR Atrial Caval, 34 Fr./38 Fr. (11.3mm/12mm) Wirewound 1/2" Connector with Cortiva™ BioActive Surface

Cardioplegia Circuit - Antegrade Perfusion Cannulae

- CB10012 Standard Aortic Root Cannula, 12ga (9 Fr.) with Cortiva™ BioActive Surface
- CB20012 Standard Aortic Root Cannula with Vent Line, 12ga (9 Fr.) with Cortiva™ BioActive Surface

Suction Circuit and Surgical Accessories - Suction Circuit Products

- CB12008 Left Heart Vent Catheter, PVC - Pediatric 10 Fr. (3.3mm), preformed 1/5" tip, vented slip conn., 13" (33.0cm) with Cortiva™ BioActive Surface



Medtronic

Europe
Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

medtronic.eu

Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. ®, ™*Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

UC202211218EE © Medtronic 2022.
All rights reserved.

Medtronic.eu

For more information, contact your local Medtronic Cannula Products representative.
The products in this catalog are available for our customers in Europe, Middle East, Africa and Russia.

See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at www.medtronic.eu.

For applicable products, consult instructions for use on www.medtronic.com/manuals. Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

This information is intended only for users in markets where Medtronic products and therapies are approved or available for use as indicated within the respective product manuals. Content on specific Medtronic products and therapies is not intended for users in markets that do not have authorization for use.

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

No. CE 598775
Issued To: Osypka AG
Earl-H.-Wood-Straße 1
Rheinfelden
79618
Germany

In respect of:

The design, development and manufacture of RF ablation catheters, multipolar steerable diagnostic catheters, diagnostic catheters for intracardiac recording and temporary pacing, valvuloplasty balloon catheters, endovascular grasping and snare/loop-catheters, TMA temporary pacing leads, external cardiac pacemaker and cardioversion devices and those aspects related to obtaining and maintaining sterility for temporary pacing lead adapters, cables DX, PK and TX, cables for EP catheters and cables for pacing leads.

on the basis of our examination of the quality assurance system under the requirements of Council Directive 93/42/EEC, Annex II excluding section 4. The quality assurance system meets the requirements of the directive. For the placing on the market of class III products an Annex II section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):



Albert Roossien, Regulatory Lead

First Issued: 2016-07-22

Date: 2019-02-27

Expiry Date: 2023-12-18

...making excellence a habit.™


Page 1 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

EndoStim™

TEMPORARY PACING LEADS



Temporary Pacing Leads

EndoStim™

Pacing leads for transvenous applications to temporarily stimulate the ventricle or the atrium. Based on the application of external pacemakers suitable for emergency, pre- and postoperative stimulation as well as for rapid-pacing during TAVI procedures.

Characteristics


The EndoStim™ portfolio includes bipolar and quadripolar pacing leads. A floating balloon catheter completes the product range.

The following French sizes, shapes, lengths, inter-electrode distances and connectors are available:

	Bipolar EndoStim™ Leads	Quadripolar EndoStim™ Leads	Pacing Catheter with Balloon
French Size	4F, 5F, 6F	4F, 5F, 6F	5F
Length	90 cm, 125 cm	90 cm, 125 cm	110 cm
Shape	Straight, Cournand, Josephson	Straight, Josephson	Cournand
Inter-Electrode Distance	2 mm, 5 mm, 10 mm	5 mm, 10 mm, 2-5-2 mm	10 mm
Connector	Exposed 2 mm pin with protective cap	Shrouded 2 mm pin	Exposed 2 mm pin with protective cap

Shape

EndoStim™ pacing leads are available in three shapes: Straight, Cournand and Josephson.



Connector

EndoStim™ pacing leads are either available with shrouded pins or with exposed pins with protective cap.



shrouded pins
(Bipolar pacing leads)



exposed pins with protective cap
(Quadripolar pacing leads, pacing lead with balloon)

Order Information Temporary Pacing Leads

Bipolar EndoStim™ Temporary Pacing Leads				(Packaging Unit: 5 pieces)		
Article No.	Number of Poles	French Size	Length	Shape	Inter-Electrode Distance	Connector
27039-F	2	4F	90 cm	Straight	2 mm	Exposed Pins
27003-F	2	4F	90 cm	Straight	5 mm	Exposed Pins
27010-F	2	4F	90 cm	Straight	10 mm	Exposed Pins
27040-F	2	4F	90 cm	Cournand	2 mm	Exposed Pins
27018-F	2	4F	90 cm	Cournand	5 mm	Exposed Pins
27015-F	2	4F	90 cm	Cournand	10 mm	Exposed Pins
27038-F	2	4F	90 cm	Josephson	2 mm	Exposed Pins
27006-F	2	4F	90 cm	Josephson	5 mm	Exposed Pins
27009-F	2	4F	90 cm	Josephson	10 mm	Exposed Pins
27051-F	2	5F	90 cm	Straight	2 mm	Exposed Pins
27045-F	2	5F	90 cm	Straight	5 mm	Exposed Pins
27048-F	2	5F	90 cm	Straight	10 mm	Exposed Pins
27052-F	2	5F	90 cm	Cournand	2 mm	Exposed Pins
27046-F	2	5F	90 cm	Cournand	5 mm	Exposed Pins
27049-F	2	5F	90 cm	Cournand	10 mm	Exposed Pins
27050-F	2	5F	90 cm	Josephson	2 mm	Exposed Pins
27044-F	2	5F	90 cm	Josephson	5 mm	Exposed Pins
27047-F	2	5F	90 cm	Josephson	10 mm	Exposed Pins
27042-F	2	5F	125 cm	Straight	2 mm	Exposed Pins
27004-F	2	5F	125 cm	Straight	5 mm	Exposed Pins
27012-F	2	5F	125 cm	Straight	10 mm	Exposed Pins
27043-F	2	5F	125 cm	Cournand	2 mm	Exposed Pins
27019-F	2	5F	125 cm	Cournand	5 mm	Exposed Pins
27016-F	2	5F	125 cm	Cournand	10 mm	Exposed Pins
27041-F	2	5F	125 cm	Josephson	2 mm	Exposed Pins
27007-F	2	5F	125 cm	Josephson	5 mm	Exposed Pins
27011-F	2	5F	125 cm	Josephson	10 mm	Exposed Pins
27063-F	2	6F	90 cm	Straight	2 mm	Exposed Pins
27057-F	2	6F	90 cm	Straight	5 mm	Exposed Pins
27060-F	2	6F	90 cm	Straight	10 mm	Exposed Pins
27055-F	2	6F	90 cm	Cournand	2 mm	Exposed Pins
27058-F	2	6F	90 cm	Cournand	5 mm	Exposed Pins
27061-F	2	6F	90 cm	Cournand	10 mm	Exposed Pins
27062-F	2	6F	90 cm	Josephson	2 mm	Exposed Pins
27056-F	2	6F	90 cm	Josephson	5 mm	Exposed Pins
27059-F	2	6F	90 cm	Josephson	10 mm	Exposed Pins
27054-F	2	6F	125 cm	Straight	2 mm	Exposed Pins
27005-F	2	6F	125 cm	Straight	5 mm	Exposed Pins
27014-F	2	6F	125 cm	Straight	10 mm	Exposed Pins
27055-F	2	6F	125 cm	Cournand	2 mm	Exposed Pins
27020-F	2	6F	125 cm	Cournand	5 mm	Exposed Pins
27017-F	2	6F	125 cm	Cournand	10 mm	Exposed Pins
27053-F	2	6F	125 cm	Josephson	2 mm	Exposed Pins
27008-F	2	6F	125 cm	Josephson	5 mm	Exposed Pins
27013-F	2	6F	125 cm	Josephson	10 mm	Exposed Pins

Temporary Pacing Leads

Order Information Temporary Pacing Leads

Quadripolar EndoStim™ Temporary Pacing Leads				(Packaging Unit: 5 pieces)		
Article No.	Number of Poles	French Size	Length	Shape	Inter-Electrode Distance	Connector
27022-F	4	4F	90 cm	Straight	5 mm	Shrouded Pins
27028-F	4	4F	90 cm	Straight	10 mm	Shrouded Pins
27033-F	4	4F	90 cm	Straight	2-5-2 mm	Shrouded Pins
27021-F	4	4F	90 cm	Josephson	5 mm	Shrouded Pins
27027-F	4	4F	90 cm	Josephson	10 mm	Shrouded Pins
27032-F	4	4F	90 cm	Josephson	2-5-2 mm	Shrouded Pins
27024-F	4	5F	125 cm	Straight	5 mm	Shrouded Pins
27030-F	4	5F	125 cm	Straight	10 mm	Shrouded Pins
27035-F	4	5F	125 cm	Straight	2-5-2 mm	Shrouded Pins
27023-F	4	5F	125 cm	Josephson	5 mm	Shrouded Pins
27029-F	4	5F	125 cm	Josephson	10 mm	Shrouded Pins
27034-F	4	5F	125 cm	Josephson	2-5-2 mm	Shrouded Pins
27026-F	4	6F	125 cm	Straight	5 mm	Shrouded Pins
27031-F	4	6F	125 cm	Straight	10 mm	Shrouded Pins
27037-F	4	6F	125 cm	Straight	2-5-2 mm	Shrouded Pins
27025-F	4	6F	125 cm	Josephson	5 mm	Shrouded Pins
27031-F	4	6F	125 cm	Josephson	10 mm	Shrouded Pins
27036-F	4	6F	125 cm	Josephson	2-5-2 mm	Shrouded Pins

Bipolar EndoStim™ Temporary Pacing Lead with Balloon								(Packaging Unit: 1 piece)	
Article No.	Number of Poles	French Size	Length	Shape	Inter-Electrode Distance	Ballon Diameter	Ballon Volume	Connector	
29207-BS	2	5F	110 cm	Cournand	10 mm	9 mm	1.0 cc	Exposed Pins	



Latex balloon

EasyStim™ Adapter – Accessory for Temporary Pacing Leads

The EasyStim™ adapter allows to connect shrouded 2mm pins to collets for exposed (naked) 2mm pins.

EasyStim™ Adapter		(Packaging Unit: 20 double packs)
Article No.	Product Name	Sterility
94005	EasyStim	no-sterile
94006	EasyStim.S	sterile



Osypka Medical GmbH

Albert-Einstein-Strasse 3
 12489 Berlin, Germany
 Phone: +49 (30) 6392 8300
 Fax: +49 (30) 6392 8301
 E-Mail: mail@osypkamed.com
 www.osypkamed.com



EasyStim and EndoStim are trademarks of Osypka Medical,
 SPIKE FLOW is a trademark of FIAB
 © Osypka Medical 2022



0344



Benannt durch/Designated by
Zentralstelle der Länder
für Gesundheitsschutz
bei Arzneimitteln und
Medizinprodukten
www.zlg.de
ZLG-BS-244.10.08



Product Service

EC Certificate

Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

(Devices in Class IIa, IIb or III)

No. G1 077790 0060 Rev. 00

Manufacturer:

Covidien LLC

15 Hampshire Street
Mansfield MA 02048
USA

**Product Category(ies): Oximetry and Capnography Monitor Systems
Temperature Monitor Systems, Patient Warming
Device Systems, Disposable Airway Management
Devices, Tracheal Tubes, Tracheostomy Tubes,
Speaking Valves, and Intubating Stylets, Ventilator
Systems and Patient Interface Circuit Systems,
EEG Monitoring Systems, Breathing Therapy and
Humidification, Heated Inspiratory Line
Humidifiers, Multi-patient Physiologic Monitoring
System and Data Analytics Software,
Gastrointestinal Measurement and Dilation System,
Electrosurgical Diathermy System Electrode.**

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

Report No.: 72145607

Valid from: 2020-06-29

Valid until: 2024-05-26

Date, 2020-06-29

Christoph Dicks
Head of Certification/Notified Body

SOLUȚII PENTRU TOT CEEA CE ESTE ESENȚIAL PENTRU INTUBAȚIE.

Catalog de produse cu tuburi endotraheale Shiley™



Medtronic
Further. Together™



CUPRINS

Ghid de informare rapidă	04
Soluții generale pentru sala de operații	06
Soluții de specialitate sala de operații	10
Dificultăți de intubație și soluții în caz de urgență	24
Accesorii pentru tuburile de intubație	26

GHID DE INFORMARE RAPIDĂ
TUBURI ENDOTRAHEALE

DESCRIERE	CANTIT. PER CUTIE	MĂRIME (mm)				
		2.0	2.5	3.0	3.5	4.0
Tub traheal oral/nazal TaperGuard™, cu orificiu de tip „ochiul lui Murphy”	10					
Tub traheal TaperGuard™ cu stilet, cu orificiu de tip „ochiul lui Murphy”	10					
Tub traheal oral evacuare TaperGuard™, cu orificiu de tip „ochiul lui Murphy”	10					
Tub traheal oral/nazal SealGuard™, cu orificiu de tip „ochiul lui Murphy”	10					
Tub traheal oral evacuare SealGuard™, cu orificiu de tip „ochiul lui Murphy”	10					
Tub traheal oral/nazal, fără manșetă, cu orificiu de tip „ochiul lui Murphy”	10	9320E	9325E	9336E	9335E	9342E
Tub traheal oral/nazal, cu manșetă, cu orificiu de tip „ochiul lui Murphy”	10			9430E		9440E
Tub traheal oral/nazal Hi-Contour cu manșetă ce asigură volum mare și presiune scăzută, cu orificiu de tip „ochiul lui Murphy”	10			107-30G	107-35G	107-40G
Tub traheal oral/nazal Hi-Contour cu manșetă ce asigură cu volum mare și presiune scăzută, cu orificiu de tip „ochiul lui Murphy”, sistem BRANDT™	10					
Tub traheal oral/nazal, fără manșetă, cu orificiu de tip „ochiul lui Murphy”	10	86232	86233	86234	86235	86236
Tub traheal oral/nazal cu formă adaptabilă anatomic, fără manșetă, cu orificiu de tip „ochiul lui Murphy”	10		411-25	411-30	411-35	411-40
Tub traheal RAE oral, cu manșetă, cu orificiu de tip „ochiul lui Murphy”	10					115-40
Tub traheal RAE oral, cu manșetă, cu orificiu de tip „ochiul lui Murphy”, sistem BRANDT™	10					315-40
Tub traheal RAE oral, fără manșetă, cu orificiu de tip „ochiul lui Murphy”	10			113-30	113-35	113-40
Tub traheal RAE nazal, cu manșetă, cu orificiu de tip „ochiul lui Murphy”	10					
Tub traheal RAE nazal, cu manșetă, sistem BRANDT™, cu orificiu de tip „ochiul lui Murphy”	10					
Tub traheal RAE nazal, cu manșetă, din material moale	10					
Tub traheal RAE nazal, fără manșetă, cu orificiu de tip „ochiul lui Murphy”	10			114-30	114-35	114-40
Tub traheal oral/nazal Lo-Contour cu manșetă ce asigură volum mare și presiune scăzută, ranforsat	10			118-30	118-35	118-40
Tub traheal oral/nazal, cu manșetă, ranforsat, cu orificiu de tip „ochiul lui Murphy”	10					
Tub traheal oral/nazal Lo-Contour cu manșetă ce asigură volum mare și presiune scăzută, ranforsat, sistem BRANDT™	10					
Tub traheal oral/nazal, cu manșetă, ranforsat, cu stilet	10					
Tub traheal oral/nazal, cu manșetă, ranforsat, cu stilet	10					
Tub traheal oral/nazal, fără manșetă, ranforsat	10			127-30-2	127-35-2	127-40-2
Tub traheal oral/nazal, fără manșetă, ranforsat, foarte moale	10		127-25	127-30	127-35	127-40
Tub traheal oral/nazal cu laser, cu manșetă dublă, cu orificiu de tip „ochiul lui Murphy”	1					
Tub traheal oral/nazal cu laser, fără manșetă, cu pensă Magill	1			161-30	161-35	161-40
Tub traheal microlaringeal oral/nazal, cu manșetă, cu orificiu de tip „ochiul lui Murphy”	10					121-40
Tub traheal oral Oxford, cu manșetă	10					
Tub traheal pentru laringectomie, cu manșetă, cu orificiu de tip „ochiul lui Murphy”	10					

TUBURI PENTRU INTUBAȚIA ENDOBRONȘICĂ

DESCRIERE	CANTIT. PER CUTIE	MĂRIME (Fr)				
		24	26	28	30	32
Tub endobronșic, pentru partea stângă	1			125028		125032
Tub endobronșic, pentru partea stângă	5			125528		125532
Tub endobronșic, pentru partea dreaptă	1					
Tub endobronșic, pentru partea dreaptă	5					
Tub endobronșic, pentru partea stângă, cu pinten de oprire la nivelul carinei	1					
Tub endobronșic, pentru partea stângă, cu sistem CPAP	1			125128		125132
Tub endobronșic, pentru partea dreaptă, cu sistem CPAP	1					
Sondă esofago-traheală cu lumen dublu Folosiți mărimea 37 pentru pacienții cu înălțimea de la 4' (122 cm) până la 5' 6" (183 cm) Folosiți mărimea 41 pentru pacienții cu înălțimea de la 5' (152 cm) și mai înalți						
Sondă esofago-traheală cu lumen dublu, cu linii de raze X	4					
Sondă esofago-traheală cu lumen dublu — cu o singură linie	4					
Sondă esofago-traheală cu lumen dublu — cu linie Rollup Kit	4					
Model didactic Trainer de sondă esofago-traheală cu lumen dublu, nu pentru uz uman, ci doar cu rol demonstrativ, de instruire pe manechin	1					
Tub de intubare nazofaringian	10		8888-247023	8888-247031	8888-247049	8888-247056

MÄRIME (mm)											
4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0	9.5	10
	18750	18755	18760	18765	18770	18775	18780	18785	18790	18795	18710
	18750S	18755S	18760S	18765S	18770S	18775S	18780S	18785S	18790S		18710S
			18860	18865	18870	18875	18880	18885	18890		
	109850	109855	109860	109865	109870	109875	109880	109885	109890		109810
			110860	110865	110870	110875	110880	110885	110890		
9345E	9350E	9360E	9366E	9365E	9370E						
	9450E	9555E	9460E	9465E	9570E	9475E	9480E	9485E	9590E	9495E	9500E
107-45G	107-50G	107-55G	107-60	107-65	107-70	107-75	107-80	107-85	107-90	107-95	107-10
	123-50	123-55	123-60	123-65	123-70	123-75	123-80	123-85	123-90	123-95	
86237	86238	86239	86240	86241	86242						
411-45	411-50	411-55	411-60	411-65	411-70						
115-45	115-50	115-55	115-60	115-65	115-70	115-75	115-80	115-85	115-90		
315-45	315-50	315-55	315-60	315-65	315-70	315-75	315-80	315-85	315-90		
113-45	113-50	113-55	113-60	113-65	113-70						
			119-60	119-65	119-70	119-75	119-80				
			319-60	319-65	319-70	319-75	319-80				
			119760	119765	119770	119775	119780				
114-45	114-50	114-55	114-60	114-65	114-70						
118-45	118-50	118-55	118-60	118-65	118-70	118-75	118-80	118-75	118-90	118-95	
			118-60M	118-65M	118-70M	118-75M	118-80M	118-75M	118-90M	118-95M	
	318-50	318-55	318-60	318-65	318-70	318-75	318-80	318-85	318-90	318-95	
	117050	117055	117060	117065	117070	117075	117080	117085	117090		
	117150	117155	117160	117165	117170	117175	117180	117185	117190		
127-45-2	127-50-2	127-55-2	127-60-2	127-65-2	127-70-2						
127-45	127-50	127-55									
160-45	160-50	160-55	160-60								
	121-50		121-60								
			128-60		128-70		128-80		128-90		128-10
					129-70		129-80				

MÄRIME (Fr)				
34	35	37	39	41
	125035	125037	125039	125041
	125535	125537	125539	125541
	126035	126037	126039	126041
	126535	126537	126539	126541
	125235	125237	125239	125241
	125135	125137	125139	125141
	126135	126137	126139	126141
		5-18537		5-18541
		5-18237		5-18241
		5-18437		5-18441
				5-18141
8888-247064				

SOLUȚII GENERALE PENTRU SALA DE OPERAȚII



Mască laringiană Shiley™

Caracteristici și beneficii

Pentru pacienți, mască laringiană Shiley™

- Tub cu formă adaptabilă anatomic
- Manșetă moale cu linie de umflare integrată
- De unică folosință, sterilă
- Fără latex
- Fără DEHP

Pentru furnizorii de servicii medicale, mască laringiană Shiley™ oferă

- Tub transparent cu formă anatomică
- Marcaje informaționale ușor de citit
- Vârf moale rotunjit pentru introducere facilă
- Gamă completă de mărimi de la copii la adulți
- Model durabil din PVC



Informații referitoare la comandă

REF	DESCRIERE	GREUTATEA PACIENTULUI (kg)	I.D. (mm)	O.D. (mm)	MANȘETĂ Ø (mm)
33710	Mască laringiană 1.0 Shiley™	<5	6.0	10.0	4
33715	Mască laringiană 1.5 Shiley™	5–10	7.0	11.6	7
33720	Mască laringiană 2.0 Shiley™	10–20	9.0	13.5	10
33725	Mască laringiană 2.5 Shiley™	20–30	9.0	13.5	14
33730	Mască laringiană 3.0 Shiley™	30–50	10.0	15.8	20
33740	Mască laringiană 4.0 Shiley™	50–70	10.0	15.8	30
33750	Mască laringiană 5.0 Shiley™	70–100	11.5	17.5	40
33760	Mască laringiană 6.0 Shiley™	>100	11.5	17.5	50

FLEXIBIL. EFICIENT.

Tub traheal oral/nazal Hi-Contour cu manșetă de etanșare ce asigură volum mare și presiune mică

Caracteristici și beneficii

- Manșetă ce asigură volum mare și presiune mică
- Manșetă flexibilă, cu pereți subțiri
- Marcaje de profunzime deasupra manșetei
- Linie de raze x pe pe toată lungimea
- Orificiu de tip „ochiul lui Murphy” încorporat
- Potrivit atât pentru intubația orală cât și pentru cea nazală

Informații referitoare la comandă

	REF	Sistem REF BRANDT™	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	MANȘETĂ Ø (mm)
PEDIATRIC	107-30G	-	3.0	4.3	190	6
	107-35G	-	3.5	4.9	213	7
	107-40G	-	4.0	5.6	232	8
	107-45G	-	4.5	6.2	253	11
	107-50G	123-50	5.0	6.9	274	18
	107-55G	123-55	5.5	7.5	304	21
	107-60	123-60	6.0	8.2	317	22
	107-65	123-65	6.5	8.8	327	23
	107-70	123-70	7.0	9.6	339	25
	107-75	123-75	7.5	10.2	350	26
107-80	123-80	8.0	10.9	360	30	
107-85	123-85	8.5	11.5	370	32	
107-90	123-90	9.0	12.1	375	33	
107-95	123-95	9.5	12.8	375	33	
107-10	-	10.0	13.5	375	34	

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.



BRANDT™



Tub traheal oral/nazal Curity™, cu manșetă/fără manșetă

Tuburi de bază standard

Pentru utilizare în timpul intubării pe termen scurt în sala de operații și în urgențe.

Caracteristici și beneficii

- Tub traheal cu manșetă cu valvă Luer și balon pilot
- Manșetă ce asigură volum mare și presiune mică
- Marcaje ce indică profunzimea intubației și conector de 15 mm premontat
- Linie de raze x pe toată lungimea, pentru verificarea siguranței poziționării
- Orificiu de tip „ochiul lui Murphy” încorporat

Informații referitoare la comandă

REF	DESCRIERE	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	MANȘETĂ Ø (mm)
-----	-----------	--------------	--------------	-----------------	-------------------

Fără manșetă


PEDIATRIC	9320E	Tub endotraheal 2.0 mm UNCUF	2.0	2.9	130	-
	9325E	Tub endotraheal 2.5 mm UNCUF	2.5	3.6	140	-
	9336E	Tub endotraheal 3.0 mm UNCUF	3.0	4.2	160	-
	9335E	Tub endotraheal 3.5 mm UNCUF	3.5	4.9	180	-
	9342E	Tub endotraheal 4.0 mm UNCUF	4.0	5.5	200	-
	9345E	Tub endotraheal 4.5 mm UNCUF	4.5	6.2	220	-
	9350E	Tub endotraheal 5.0 mm UNCUF	5.0	6.8	240	-
	9360E	Tub endotraheal 5.5 mm UNCUF	5.5	7.5	270	-
	9366E	Tub endotraheal 6.0 mm UNCUF	6.0	8.2	280	-
	9365E	Tub endotraheal 6.5 mm UNCUF	6.5	8.8	290	-
9370E	Tub endotraheal 7.0 mm UNCUF	7.0	9.6	310	-	

Cu manșetă

PEDIATRIC	9430E	Tub endotraheal 3.0 mm CUFDF	3.0	4.2	160	8
	9440E	Tub endotraheal 4.0 mm CUFDF	4.0	5.5	200	11
	9450E	Tub endotraheal 5.0 mm CUFDF	5.0	6.8	240	16
	9555E	Tub endotraheal 5.5 mm CUFDF	5.5	7.5	270	17
	9460E	Tub endotraheal 6.0 mm CUFDF	6.0	8.2	280	22
	9465E	Tub endotraheal 6.5 mm CUFDF	6.5	8.8	290	22
	9570E	Tub endotraheal 7.0 mm CUFDF	7.0	9.6	310	25
	9475E	Tub endotraheal 7.5 mm CUFDF	7.5	10.2	320	25
	9480E	Tub endotraheal 8.0 mm CUFDF	8.0	10.9	320	27
	9485E	Tub endotraheal 8.5 mm CUFDF	8.5	11.5	320	27
	9590E	Tub endotraheal 9.0 mm CUFDF	9.0	12.1	320	29
	9495E	Tub endotraheal 9.5 mm CUFDF	9.5	12.8	320	29
	9500E	Tub endotraheal 10.0 mm CUFDF	10.0	13.5	320	32

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.



Tub traheal oral/nazal, fără manșetă

Caracteristici și beneficii


- Capăt cu contur delicat și înveliș protector
- Gradații pe ambele părți ale tubului traheal și marcaje speciale la capăt
- Linie de raze x pe toată lungimea
- Orificiu de tip „ochiul lui Murphy” încorporat
- Potrivit atât pentru intubarea orală cât și pentru cea nazală



Tub traheal oral/nazal Contour, fără manșetă

Caracteristici și beneficii

- Material moale — Poate fi potrivit pentru intubarea nazală
- Tuburi semitransparente pentru vizualizarea directă a condensării aerului expirat
- Capăt cu înveliș protector special
- Gradații pe ambele părți ale tubului traheal și marcaje speciale la capăt
- Linie de raze x pe toată lungimea
- Orificiu de tip „ochiul lui Murphy” încorporat



Informații referitoare la comandă

	REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)
PEDIATRIC	86232	2.0	2.9	171
	86233	2.5	3.6	171
	86234	3.0	4.3	190
	86235	3.5	4.9	213
	86236	4.0	5.6	232
	86237	4.5	6.2	253
	86238	5.0	6.9	274
	86239	5.5	7.5	304
	86240	6.0	8.2	317
	86241	6.5	8.8	327
86242	7.0	9.6	339	

Cantitatea cutiei: 10
Toate tuburile sunt livrate în ambalaj steril, de unică folosință.

Informații referitoare la comandă

	REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)
PEDIATRIC	411-25	2.5	3.6	171
	411-30	3.0	4.3	190
	411-35	3.5	5.1	213
	411-40	4.0	6.0	232
	411-45	4.5	6.5	253
	411-50	5.0	6.9	274
	411-55	5.5	7.5	304
	411-60	6.0	8.2	317
	411-65	6.5	8.8	327
	411-70	7.0	9.6	339

Cantitatea cutiei: 10
Toate tuburile sunt livrate în ambalaj steril, de unică folosință.



SOLUȚII DE SPECIALITATE PENTRU SALA DE OPERAȚII

AJUTĂ LA REDUCEREA IMPACTULUI INTUBAȚIEI ASUPRA TRAHEEI.¹

Tuburi traheale orale/nazale TaperGuard™

Caracteristici și beneficii

- Manșetă conică
- Marcaje de profunzime deasupra manșetei
- Vârf cu înveliș protector special, cu orificiu de tip „ochiul lui Murphy”
- Potrivit pentru intubația orală și cea nazală
- Linie de raze x pe toată lungimea
- Disponibil și cu stilet de intubație preîncărcat

Informații referitoare la comandă

REF	REF with stylet	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	MANȘETĂ Ø (mm)
18750	18750S	5.0	6.9	274	17.4
18755	18755S	5.5	7.5	304	17.4
18760	18760S	6.0	8.2	317	20.6
18765	18765S	6.5	8.9	327	20.6
18770	18770S	7.0	9.5	339	25.4
18775	18775S	7.5	10.2	350	25.4
18780	18780S	8.0	10.8	360	25.4
18785	18785S	8.5	11.4	370	28.6
18790	18790S	9.0	12.1	375	28.6
18795		9.5	12.8	375	28.6
18710	18710S	10.0	13.5	375	28.6

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.

1. Lichtenthal PR, Borg UB. Endotracheal cuff pressure: role of tracheal size and cuff volume. *Crit Care*. 2011;15:147.



Tub traheal oral de evacuare TaperGuard™

Caracteristici și beneficii:


- Tubul endotraheal oral de evacuare Shiley™ TaperGuard încorporează tehnologia de evacuare care, în numeroase studii efectuate de-a lungul ultimului deceniu, s-a demonstrat că reduce VAP cu o medie de 50%.¹
- Orificiu de tip „ochiul lui Murphy” cu înveliș protector încorporat
- Linie de raze x pe toată lungimea

Informații referitoare la comandă

REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	MANȘETĂ Ø
Oral				
18860	6.0	9	308	20.6
18865	6.5	9.8	317	20.6
18870	7.0	10.4	330	25.4
18875	7.5	11.2	350	25.4
18880	8.0	11.8	360	25.4
18885	8.5	12.6	370	28.6
18890	9.0	13.1	375	28.6

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.



1. Muscedere J, Rewa O, McKechnie K, Jiang X, Laporta D, Heyland DK. Subglottic secretion drainage for the prevention of ventila-tor-associated pneumonia: a systematic review and meta-analysis. *Crit Care Med.* 2011;39(8):1985-1991.

Tub traheal oral de evacuare SealGuard™

Caracteristici și beneficii

- Tub endotraheal cu lumen pentru drenajul secrețiilor din regiunea subglotică și orificiu de aspirație
- Manșetă din material poliuretanic dur
- Poate reduce riscul infecției respiratorii în timpul ventilației pe termen lung prin drenarea secrețiilor din regiunea subglotică
- Linie de raze x pe toată lungimea
- Orificiu de tip „ochiul lui Murphy” încorporat

Informații referitoare la comandă

REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	MANȘETĂ Ø
Oral				
110860	6.0	9	308	22
110865	6.5	9.8	317	22
110870	7.0	10.4	330	24
110875	7.5	11.2	350	26
110880	8.0	11.8	360	27
110885	8.5	12.6	370	29
110890	9.0	13.1	375	30

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință



1. Lorente L, Lecuona M, Jiménez A, Mora ML, Sierra A. Influence of an endotracheal tube with polyurethane cuff and subglottic drainage on pneumonia. *Am J Respir Crit Care Med.* 2007 Dec 1;176(11):1079-83. Epub 2007 Sep 13.

Tub traheal oral/nazal SealGuard™

Caracteristici și beneficii

- Manșetă conică
- Manșetă din material poliuretanic dur
- Forma special retractabilă a vârfului sprijină intubația nazală
- Orificiu de tip „ochiul lui Murphy” încorporat
- Adecvat pentru intubația orală și cea nazală
- Linie de raze x pe toată lungimea

Informații referitoare la comandă

REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	MANȘETĂ Ø
109850	5.0	6.9	274	20
109855	5.5	7.5	304	20
109860	6.0	8.2	317	22
109865	6.5	8.8	327	22
109870	7.0	9.6	339	26
109875	7.5	10.2	350	27
109880	8.0	10.9	360	28
109885	8.5	11.5	370	29
109890	9.0	12.1	375	31
109810	10.0	13.5	375	33

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință



Tub traheal oral/nazal Lo-Contour, cu manșetă, ranforsat

Caracteristici și beneficii

- Spirală cu ranforsare metalică și încapsulare uniformă, securizată în peretele tubului, pentru a preveni încovoierea, comparativ cu tuburile endotraheale tradiționale
- Manșetă ce asigură volum mare și presiune mică
- Două marcaje de profunzime deasupra manșetei
- Ranforsare pentru a susține controlul razelor x și verificarea poziției tubului
- Conector de 15mm atașat de tubul traheal
- Potrivit pentru intubația orală și cea nazală
- Disponibil cu sau fără orificiul de tip „ochiul lui Murphy”
- Disponibil cu conector neatașat
- Disponibil cu stilet de intubație preîncărcat
- La folosirea tuburilor traheale ranforsate se recomandă un dispozitiv antimușcătură sau unul de tip Guedel opțional. Tubul traheal ranforsat este disponibil și cu sistemul BRANDT™



Informații referitoare la comandă

	REF	REF cu gaură Murphy Eye	REF cu sistem BRANDT™	REF cu stilet	REF cu racord demontabil de 15 mm	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	MANȘETĂ Ø (mm)	
PEDIATRIC	118-30	-	-	-	-	3.0	5.0	195	8	
	118-35	-	-	-	-	3.5	5.2	212	9	
	118-40	-	-	-	-	4.0	6.2	236	9	
	118-45	-	-	-	-	4.5	6.7	250	11	
	118-50	-	-	318-50	117050	117150	5.0	6.9	308	17
	118-55	-	-	318-55	117055	117155	5.5	7.5	320	19
	118-60	118-60M	318-60	117060	117160	6.0	8.2	320	21	
	118-65	118-65M	318-65	117065	117165	6.5	8.8	330	22	
	118-70	118-70M	318-70	117070	117170	7.0	9.6	340	24	
	118-75	118-75M	318-75	117075	117175	7.5	10.2	350	26	
	118-80	118-80M	318-80	117080	117180	8.0	10.8	360	28	
	118-85	118-85M	318-85	117085	117185	8.5	11.5	365	29	
	118-90	118-90M	318-90	117090	117190	9.0	12.1	370	31	
	118-95	118-95M	318-95	-	-	9.5	12.8	370	32	

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.

BRANDT™



Gaură Murphy Eye



Tub traheal oral/nazal, fără manșetă, ranforsat, foarte moale

Caracteristici și beneficii

- Material siliconic foarte flexibil și ranforsare a peretelui interior al tubului, pentru a preveni încovoierea, comparativ cu tuburile endotraheale tradiționale
- Marcare pe ambele laturi ale vârfului
- Ranforsare pentru a susține controlul razelor x și verificarea poziției tubului
- Conector de 15 mm inserat complet
- Pentru o amplasare mai comodă poate fi recomandat stiletul de intubație
- Potrivit pentru intubația orală și cea nazală
- La folosirea tuburilor traheale ranforsate se poate recomanda un dispozitiv antimușcătură sau unul de tip Guedel



Tub traheal oral/nazal SafetyFlex fără manșetă, ranforsat

Caracteristici și beneficii

- Spirală cu ranforsare metalică și încapsulare uniformă, securizată în peretele tubului, pentru a preveni încovoierea, comparativ cu tuburile endotraheale tradiționale
- Marcarea vârfului poate ajuta la poziționarea corectă a tubului traheal dacă este posibilă vizualizarea corzilor vocale
- Ranforsare pentru a susține controlul razelor x și verificarea poziției tubului
- Conector de 15 mm
- Potrivit pentru intubația orală și cea nazală
- La folosirea tuburilor traheale ranforsate se poate recomanda un dispozitiv antimușcătură sau unul de tip Guedel

Informații referitoare la comandă

	REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)
PEDIATRIC	127-30-2	3.0	5.0	192
	127-35-2	3.5	5.2	212
	127-40-2	4.0	6.2	232
	127-45-2	4.5	6.7	252
	127-50-2	5.0	6.9	290
	127-55-2	5.5	7.5	298
	127-60-2	6.0	8.2	308
	127-65-2	6.5	8.8	318
	127-70-2	7.0	9.6	328

Informații referitoare la comandă

	REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)
PEDIATRIC	127-25	2.5	4.0	159
	127-30	3.0	4.7	179
	127-35	3.5	5.3	202
	127-40	4.0	6.0	222
	127-45	4.5	6.7	242
	127-50	5.0	7.3	262
	127-55	5.5	8.0	294

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.


Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.

Tub traheal RAE oral și nazal, cu manșetă

Caracteristici și beneficii

- Tub preformat pentru a diminua riscul de încovoiere¹
- Acces îmbunătățit la câmpul operator, conexiunile fiind situate în afara acestei zone¹
- Manșetă ce asigură volum mare și presiune mică
- Vârf cu înveliș protector special
- Orificiu de tip „ochiul lui Murphy” încorporat
- Linie de raze x pe toată lungimea



Informații referitoare la comandă

	REF	REF cu sistem BRANT™	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	MANȘETĂ Ø (mm)
PEDIATRIC	Oral					
	115-40	315-40	4.0	5.6	216	8
	115-45	315-45	4.5	6.2	237	11
	115-50	315-50	5.0	6.9	251	19
	115-55	315-55	5.5	7.6	263	21
	115-60	315-60	6.0	8.2	287	22
	115-65	315-65	6.5	8.7	297	24
	115-70	315-70	7.0	9.5	317	28
	115-75	315-75	7.5	10.1	327	30
	115-80	315-80	8.0	10.8	345	31
	115-85	315-85	8.5	11.1	355	33
	115-90	315-90	9.0	11.9	374	34
	Nazal					
	119-60	319-60	6.0	8.2	386	22
	119-65	319-65	6.5	8.7	396	24
	119-70	319-70	7.0	9.5	406	28
	119-75	319-75	7.5	10.1	416	30
	119-80	319-80	8.0	10.8	425	31

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.



BRANDT™



1. Nagelhout, John J., and Karen L. Plaus. „38/Anesthesia for Ear, Nose, Throat, Maxillofacial Surgery. *Nurse Anesthesia*. St. Louis: Elsevier, 2005.

Tub traheal RAE oral și nazal, fără manșetă

Caracteristici și beneficii

- Tub traheal preformat pentru diminuarea riscului de încovoiere¹
- Acces îmbunătățit la câmpul operator, conexiunile fiind situate în afara acestei zone¹
- Două orificii distale de tip „ochiul lui Murphy”, câte unul pe fiecare parte
- Linie de raze x pe toată lungimea

Informații referitoare la comandă

REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	
Oral				
PEDIATRIC	113-30	3.0	4.2	174
	113-35	3.5	4.8	195
	113-40	4.0	5.6	210
	113-45	4.5	6.2	232
	113-50	5.0	6.9	246
	113-55	5.5	7.6	263
	113-60	6.0	8.2	273
	113-65	6.5	8.7	285
113-70	7.0	9.5	300	
Nazal				
PEDIATRIC	114-30	3.0	4.2	190
	114-35	3.5	4.8	220
	114-40	4.0	5.6	240
	114-45	4.5	6.2	272
	114-50	5.0	6.9	286
	114-55	5.5	7.6	293
	114-60	6.0	8.2	303
	114-65	6.5	8.7	310
114-70	7.0	9.5	315	

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.



1. Nagelhout, John J., and Karen L. Plaus. „38/Anesthesia for Ear, Nose, Throat, Maxillofacial Surgery. *Nurse Anesthesia*. St. Louis: Elsevier, 2005.

Tub traheal RAE nazal cu manșetă și finisaj moale

Caracteristici și beneficii

- Tub traheal moale, cu manșetă ai cărei pereți subțiri asigură un volum mare și o presiune scăzută
- Material mai moale, cu o formulă specială pentru intubația nazală
- Linie de raze x pe toată lungimea
- Două marcaje de profunzime circulare pot ajuta la poziționarea corectă
- Curbură dublă, doar pentru intubația nazală

Informații referitoare la comandă

REF	I.D. (mm)	O.D. (mm)	LUNGIME De la vârf la semnul nărilor (mm)	MANȘETĂ Ø (mm)
119760	6.0	8.8	28	24
119765	6.5	9.4	29	28
119770	7.0	10.2	30	30
119775	7.5	10.9	31	31
119780	8.0	11.6	32	32

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.



Tub traheal oral/nazal microlaringeal cu manșetă

Caracteristici și beneficii


- Tubul cu diametru mic poate oferi acces maxim la câmpul operator
- Lungimea și diametrul manșetei comparabile cu cele ale unui tub standard cu diametrul interior de 8 mm
- Manșetă specială ce asigură presiune scăzută și volum mare
- Orificiu de tip „ochiul lui Murphy” încorporat
- Potrivit pentru intubația nazală
- Linie de raze x pe toată lungimea

Informații referitoare la comandă

REF	I.D. (mm)	O.D. (mm)	LUNGIME Cu racord (mm)	MANȘETĂ Ø (mm)
121-40	4.0	5.6	368	25
121-50	5.0	6.9	368	27
121-60	6.0	8.2	368	27

Cantitatea cutiei: 10

Toate tuburile sunt livrate în ambalaj steril, de unică folosință.



Tub traheal pentru laringectomie cu manșetă

Caracteristici și beneficii


- Două curbe preformate separate pot facilita amplasarea tubului
- Intubare realizată prin traheostomie
- Manșetă ce asigură volum mare și presiune scăzută
- Vârf cu înveliș protector special
- Orificiu de tip „ochiul lui Murphy” încorporat
- Linie de raze x pe toată lungimea

Informații referitoare la comandă

REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)	MANȘETĂ Ø (mm)
129-70	7.0	9.5	358	28
129-80	8.0	10.8	358	31

Cantitatea cutiei: 10


Toate tuburile sunt livrate în ambalaj steril, de unică folosință.



Tub traheal oral cu laser, cu manșetă dublă și fără manșetă

Caracteristici și beneficii

- Tub din oțel inoxidabil, flexibil, neinflamabil, rezistent la laser, pentru ventilație în timpul intervenției chirurgicale cu laser CO₂ și KTP a supra regiunii laringiene/traheale
- La contactul accidental dintre axul tubului și fasciculul laser are loc reflexia unui fascicul defocalizat¹
- Profil spiralat din oțel pentru etanșeizare ermetică
- Conector de 15 mm atașat de tubul traheal
- Manșete duble cu diametru foarte mare și linii de umflare individuale
- Baloane pilot cu valve autoetanșante prevăzute cu marcajul “proximal” și “depărtat”
- Orificu de tip „ochiul lui Murphy” încorporat
- Ambele manșete pot fi umplute cu soluție izotonică pentru asigurarea unei protecții adecvate în timpul intervenției chirurgicale cu laser



Tub traheal cu laser fără manșetă

Caracteristici și beneficii

- Tub din oțel inoxidabil, flexibil, neinflamabil, rezistent la laser, pentru ventilație în timpul intervenției chirurgicale cu laser CO₂ și KTP a supra regiunii laringiene/traheale
- La contactul accidental dintre axul tubului și fasciculul laser are loc reflexia unui fascicul defocalizat¹
- Profil spiralat din oțel pentru etanșeizare ermetică
- Conector de 15 mm atașat de tubul traheal
- Vârf metalic șlefuit fin



Informații referitoare la comandă

REF	I.D. (mm)	O.D. (mm)	LUNGIME (mm)
160-45	4.5	7.0	349
160-50	5.0	7.5	349
160-55	5.5	7.9	349
160-60	6.0	8.5	349
PEDIATRIC 161-30	3.0	5.2	238
161-35	3.5	5.7	238
161-40	4.0	6.1	238

Cantitatea cutiei: 1

Toate tuburile sunt livrate în ambalaj steril, de unică folosință

1. Lai HC, Juang SE, Lui TJ, Ho WM. Fires of endotracheal tubes of three different materials during carbon dioxide laser surgery. *Acta Anesthesiology Sin.* 2002; 40(1): 47-51.

Tub endobronșic, pentru partea stângă, cu manșetă din poliuretan

Caracteristici și beneficii

- Pentru intubarea endobronșică pe partea stângă, atunci când se indică ventilația diferențială sau pe un singur plămân în chirurgia toracică sau în terapia intensivă
 - Materialul poliuretanic al manșetei asigură o forță de punșionare superioară a celei obținute prin folosirea de materiale convenționale¹
 - Manșeta bronșică cu profil specific poate contribui la localizarea vârfului distal atunci când verificarea este confirmată de bronhoscopul din fibră optică
 - Manșeta bronșică albastră poate contribui la identificarea vârfului distal când se verifică amplasarea corectă printr-un bronhoscop din fibră optică
- Vârful tubului bronșic este ușor curbat pentru a ajuta la îmbunătățirea poziționării în bronhia principală stângă
- Markeri de radioopacitate la vârful distal, deasupra manșetei bronșice și la deschiderea traheală, pentru a ajuta la localizarea și verificarea poziției tubului
- Pachetul cu accesorii conține 2 conectoare Opti-Port™, 2 catetere de aspirație a secrețiilor endobronșice
- Gentle-Flo™ și 1 conector "Y"
- Disponibil și cu pinten de oprire la nivelul carinei a radioopacității, pentru localizarea exactă a tubului endobronșic în interiorul carinei



Informații referitoare la comandă

Standard REF	REF Cu sistem CPAP	REF Cu cârlig bronșic	REF cantitate: 5 (fără accesorii)	I.D. (mm)*	O.D. (mm)	O.D. (Fr)	LUNGIME (mm)
Stânga							
125028	125128	-	125528	3.1	9.3	28	420
125032	125132	-	125532	3.4	10.7	32	420
125035	125135	125235	125535	4.8	11.7	35	420
125037	125137	125237	125537	5.1	12.3	37	420
125039	125139	125239	125539	5.3	13.0	39	420
125041	125141	125241	125541	5.4	13.7	41	420

Accesorii: Conține REF 125-10 și 2 catetere de aspirație Gentle-Flo™ de mari dimensiuni

* Limita aplicabilă atunci când se folosește bronhoscopul din fibră optică


1. Pe baza testării pe banc

Tub endobronșic, pentru partea dreaptă, cu manșetă din poliuretan

Pentru intubarea endobronșică pe partea dreaptă, atunci când atunci când se indică ventilația diferențială sau pe un singur plămân în chirurgia toracică sau în terapia intensivă.

Caracteristici și beneficii

- Materialul poliuretanic al manșetei asigură o forță de puncționare superioară celei obținute prin folosirea de materiale convenționale¹
- Manșetă bronșică marcată cu cod albastru, lumen proximal bronșic și balon pilot bronșic, pentru identificare facilă
- Manșeta bronșică specială în formă de S poate oferi o marjă crescută de siguranță la poziționarea tubului endobronșic în apropierea lobului superior drept
- Manșeta bronșică albastră poate contribui la identificarea vârfului distal când se verifică amplasarea corectă cu ajutorul unui bronhoscop din fibră optică
- Markeri de radioopacitate la nivelul capătului distal, în jurul orificiului lateral, deasupra manșetei bronșice și la nivelul orificiului traheal ajută la localizarea și verificarea poziției tubului
- Pachetul cu accesorii cuprinde 2 conectoare Opti-Port™, 2 catetere de aspirație a secrețiilor endobronșice Gentle-Flo™ și 1 racord "Y"



Informații referitoare la comandă

Standard de REF	REF Cu sistem CPAP	REF Cu cârlig bronșic	Cantitatea REF : 5 (fără accesorii)	I.D. (mm)*	O.D. (mm)	O.D. (Fr)	LUNGIME (mm)
Dreapta							
126035	126135	-	126535	4.8	11.7	35	420
126037	126137	-	126537	5.1	12.3	37	420
126039	126139	-	126539	5.3	13.0	39	420
126041	126141	-	126541	5.4	3.7	41	420

Accesorii: Conține REF 125-10 și 2 catetere de aspirație Gentle-Flo™ de mari dimensiuni

* Limita aplicabilă atunci când se folosește bronhoscopul din fibră optică

Pachet cu accesorii pentru tuburile endobronșice

Caracteristici și beneficii

- Conectorul pivotant în unghi drept, rotativ, permite racordarea și poate reduce riscul de dislocare a tuburilor endobronșice
- Poate facilita utilizarea instrumentelor din fibră optică și introducerea cateterelor de aspirație (e.g., Gentle-Flo™)
- Conținut: 2 conectoare Opti-Port™, 1 conector Y
- Inclus în pachetele individuale ale tuburilor endobronșice

Informații referitoare la comandă

REF	Descriere	Cutie /Cantitate
125-10	Pachet cu accesorii	10



Sistemul CPAP¹

Pentru evitarea sau tratarea hipoxemiei din timpul anesteziei pe un singur plămân.

Caracteristici și beneficii

- Valva CPAP permite ajustarea nivelurilor de presiune de la 1 la 10 cm H₂O
- Conectorul pivotant în unghi drept, rotativ, permite atașarea la tubul endobronșic
- Compact și ușor de folosit, chiar și atunci când spațiul este limitat
- Nu conține latex

Informații referitoare la comandă

REF	Cutie / Cantitate
125-20	5

¹Supapa CPAP este autorizată de Dr. David Alfery.

DIFICULTĂȚI DE INTUBAȚIE ȘI SOLUȚII ÎN CAZ DE URGENȚĂ



Dispozitiv de intubație esofago-traheală cu lumen dublu

La stabilirea unei intubații dificile, timpul și precizia sunt cruciale. Spre deosebire de tuburile endotraheale, acest tub este proiectat pentru a stabili o permeabilitate liberă a căilor respiratorii, fără intervenție chirurgicală, când este plasat fie în trahee, fie în esofag.

Caracteristici și beneficii

- Lumenul dublu special poate stabili o intubație eficientă fie prin esofag, fie prin trahee
- Balonul faringian gonflabil, proiectat pentru a menține ferm dispozitivul pe loc și a contribui la prevenirea evacuării gazului pe nas sau pe gură
- Lumenul de mărime naturală poate facilita aspirația conținutului gastric fără a întrerupe ventilația pacientului, dacă dispozitivul de intubație esofago-traheală este amplasat în esofag
- Manșetă esofagiană gonflabilă pentru a etanșa esofagul ca să împiedice pătrunderea gazului în stomac și aspirarea conținutului gastric
- Potrivit pentru utilizarea înainte de spitalizare, în intervenții chirurgicale și în cazuri de urgență
- Linie de raze x pe toată lungimea
- Dispozitiv de intubație esofago-traheală cu lumen dublu:
 - Utilizați mărimea 37 pentru pacienții cu înălțimea de la 4' (122 cm) până la 5' 6" (183 cm)
 - Utilizați mărimea 41 pentru pacienții cu înălțimea de 5' (152 cm) și mai înalți



Informații referitoare la comandă

REF	DESCRIERE	MĂRIMEA (Fr)	CUTIE /CANTITATE
5-18537	Tăvi standard cu dispozitive pentru intubație esofago-traheală	37	4
5-18541	Tăvi standard cu dispozitive pentru intubație esofago-traheală	41	4
5-18237	Dispozitive pentru intubație esofago-traheală ambalate individual	37	4
5-18241	Dispozitive pentru intubație esofago-traheală ambalate individual	41	4
5-18437	Role cu dispozitive pentru intubație esofago-traheală	37	4
5-18441	Role cu dispozitive pentru intubație esofago-traheală	41	4
5-18141	Model didactic de dispozitiv pentru intubație esofago-traheală — dispozitiv demonstrativ	41	1

ACCESORII PENTRU TUBURI DE INTUBAȚIE



Tub de intubație Guedel

Produs de unică folosință

Caracteristici și beneficii

- Dispozitiv antimușcatură pentru intubația orală
- Formă adaptabilă anatomic
- Codat cromatic



Informații referitoare la comandă

REF	MĂRIME	LUNGIME (mm)*	COD CROMATIC
steril			
287/7839	000	30	roz
287/7838	00	40	portocaliu
287/7636	0	50	albastru
287/7635	1	60	violet
287/7634	2	70	alb
287/7633	3	80	verde
287/7632	4	90	galben
287/7631	5	100	roșu
287/7843	6	110	Vernil
Ambalat curat			
287P7839SP	000	30	roz
287P7838SP	00	40	portocaliu
287P7636SP	0	50	albastru
287P7635SP	1	60	violet
287P7634SP	2	70	alb
287P7633SP	3	80	verde
287P7632SP	4	90	galben
287P7631SP	5	100	roșu
287P7843SP	6	110	light green

* Măsurarea axială de sub capătul proximal la centrul orificiului distal
Cantitatea cutiei: 10
Ambalate individual, gata pentru utilizare

Tub de intubație Guedel

Reutilizabil

Caracteristici și beneficii

- Opac, moale, material care nu conține latex
- Dispozitiv antimușcatură pentru intubația orală
- Formă adaptabilă anatomic
- Gata pentru utilizare
- Codat cromatic

Informații referitoare la comandă

REF	MĂRIME	LUNGIME (mm)*	COD CROMATIC
287/25000	000	30	roz
287/25001	00	40	portocaliu
287/25002	0	50	albastru
287/25003	1	60	violet
287/25004	2	70	alb
287/25005	3	80	verde
287/25006	4	90	galben
287/25007	5	100	roșu
287/25008	6	110	vernil

* Măsurarea axială de sub capătul proximal la centrul orificiului distal
Cantitatea cutiei: 10
Ambalate individual, gata pentru utilizare

Tub de intubație Guedel

Steril, gata pentru utilizare, fără codare cromatică

Caracteristici și beneficii

- Opac, moale, materialul nu conține latex
- Dispozitiv antimușcătură pentru intubația orală
- Formă adaptabilă anatomic

Informații referitoare la comandă

REF	MĂRIME	LUNGIME (mm)*
287/25018	6	110
287/25017	5	100
287/25016	4	90
287/25015	3	80
287/25014	2	70
287/25013	1	60
287/25012	0	50
287/25011	00	40
287/25010	000	30

Măsurarea axială de sub capătul proximal la centrul orificiului distal

Cantitatea cutiei: 10

Ambalate individual, gata pentru utilizare

Strângerea tubului traheal

Caracteristici și beneficii

- Stabilizează tuburile traheale introduse pe cale orală și nazală fără să necesite bandă adezivă
- Fabricat din Velfoam™** pentru a reduce iritarea pielii
- Ambalate curat, doar pentru unică folosință



Informații referitoare la comandă

Strângere tub traheal REF	LUNGIME (mm)
330	700

Cantitatea cutiei: 25

Tub de intubație nazofaringian

Caracteristici și beneficii

Permite deschiderea căii respiratorii în timpul aspirării și bronhoscopiei. Materialul moale ajută la reducerea senzației de tuse și a reflexului de sufocare a pacientului. Protejează mucoasa în timpul manevrelor repetate de trecere a cateterului.

Informații referitoare la comandă

REF	I.D. (mm)	O.D. (Fr)	(mm)	LUNGIME (mm)
8888247023	6.5	26	8.7	114
8888247031	7.0	28	9.3	127
8888247049	7.5	30	10.0	140
8888247056	8.0	32	10.7	152
	8.5	34	11.3	159

Cantitatea cutiei: 10

Ambalate individual, gata pentru folosire

Stilet de intubație

Caracteristici și beneficii

- Aluminiiu maleabil cu înveliș din plastic fin pentru a contribui la intubația traheală în condiții dificile
- Vârf distal moale pentru a reduce posibilitatea producerii unor leziuni în timpul intubației
- Suprafață specială pentru a facilita introducerea și scoaterea stiletului
- Ambalat curat, doar pentru unică folosință
- Disponibil și steril

Informații referitoare la comandă

REF	TUB RECOMANDAT I.D. (mm)	LUNGIME (mm)
BANDĂ DIN SATIN		
Sondă de intubație ambalată curat		
116-06	2.5 până la 4.5	280
116-10	până la	350
116-14	> 5.0	350
Sondă de intubație ambalată individual		
116-06-S	2.5 până la 4.5	280
116-10-S	4.0 până la 6.0	350
116-14-S	> 5.0	350

Cantitatea cutiei: 10



Detector CO₂ colorimetric pentru adulți și copii

Detectoarele CO₂ colorimetrice pot verifica dacă pacientul este intubat corect chiar de la început —și dacă rămâne așa în timpul transportului.

Detectoarele noastre CO₂ vă oferă un mod facil, economic de a evalua CO₂ expirat în timpul procedurilor medicale de menținere a funcției respiratorii.

Caracteristici și beneficii

- Se atașează direct la tubul endotraheal pentru a indica nivelurile de CO₂ expirat în cursul intubației timp de până la 2 ore
- Răspunde rapid la CO₂ expirat printr-o simplă modificare a culorii de la violet galben
- Facilitează vizualizarea printr-un afișaj constant al răspunsului la fiecare respirație
- Ajută la evaluarea eficacității CPR

Informații referitoare la comandă

REF	DESCRIERE	MĂRIME RECOMANDATĂ ÎN FUNCȚIE DE GREUTATEA PACIENTULUI (kg)	VOLUM INTERN (cc)	REZISTENȚA LA FLUX	GREUTATEA DETECTORULUI (g)	ORIFICIILE RACORDULUI Capătul pacientului (mm)	Capătul circuitului (mm)	TIMPUL DE FOLOSIRE (ore)
EASYCAP II PEDICAP	Casetă cu 24	peste 15	25	3.0cm H ₂ O±1.0cm Flux de 60 L/min	Sub 20 g	22 OD/15 ID	15 OD/13 ID	Până la 2 ore
EASYCAP II 6 PEDICAP 6	Casetă cu 6	1 până la 15	3	2.5cm H ₂ O±0.5cm Flux de 10 L/min	Sub 5 g	18 OD/15 ID	15 OD/5 ID	Până la 2 ore



Materialele promoționale sunt doar pentru distribuirea în afara Statelor Unite. Nu toate produsele sunt disponibile în toate țările. Pentru informații suplimentare, contactați reprezentanții de vânzări.

IMPORTANT: Citiți instrucțiunile, contraindicațiile, avertizările și măsurile de siguranță de pe ambalaj.

Proveniența ilustrațiilor: Thinkstock

© 2016 Medtronic. Toate drepturile rezervate. Medtronic, Medtronic logo și Further, Together sunt mărci înregistrate ale Medtronic.™* Mărcile terței părți sunt mărcile de fabrică ale proprietarilor lor. Toate celelalte mărci sunt mărcile de fabrică ale companiei Medtronic.

08/2016-16-AW-0026-[Wf#1038973]

6135 Gunbarrel Avenue 800.635.5267 medtronic.com/covidien
Boulder, CO 80301

Medtronic