

TYA01~TYA04 >>>

Code 编号	Description 描述	Packing 包装
TYA01	Centrifuge Tube 50ml, Free Standing, Bulk Pack 离心管50ml, 可立式, 散装	50pcs/bag, 500pcs/ctn
TYA01-D	Centrifuge Tube 50ml, Free Standing, Individual Packing, E.O. Sterile 离心管50ml, 可立式, 独立包装, EO消毒	1pc/bag, 50pcs/bag, 500pcs/ctn
TYA02	Centrifuge Tube 50ml, Conical Bottom, Bulk Pack 离心管50ml, 锥形底, 散装	50pcs/bag, 500pcs/ctn
TYA02-D	Centrifuge Tube 50ml, Conical Bottom, Individual packing, E.O. Sterile 离心管50ml, 锥形底, 独立包装, EO消毒	1pc/bag, 50pcs/bag, 500pcs/ctn
TYA03	Centrifuge Tube 15ml, Conical Bottom, Bulk Pack 离心管15ml, 锥形底, 散装	100pcs/bag, 2000pcs/ctn
TYA03-D	Centrifuge Tube 15ml, Conical Bottom, Individual packing, E.O. Sterile 离心管15ml, 锥形底, 独立包装, EO消毒	1pc/bag, 100pcs/bag, 1000pcs/ctn
TYA04	Centrifuge Tube 10ml, Conical Bottom, Bulk Pack 离心管10ml, 锥形底, 散装	100pcs/bag, 2000pcs/ctn



TYA05~TYA08 >>>

Code 编号	Description 描述	Packing 包装
TYA05	Centrifuge Tube 10ml, Bulk Pack 离心管10ml, 散装	200pcs/bag, 3000pcs/ctn
TYA06	Centrifuge Tube 7ml, Bulk Pack 离心管7ml, 散装	200pcs/bag, 4000pcs/ctn
TYA07	Centrifuge Tube 5ml, Bulk Pack 离心管5ml, 散装	200pcs/bag, 6000pcs/ctn
TYA08	Centrifuge Tube 2ml, Bulk Pack 离心管2ml, 散装	500pcs/bag, 10000pcs/ctn



TYA09~TYA10 >>>

Code 编号	Description 描述	Packing 包装
TYA09	Centrifuge Tube 1.5ml, Shield Lid, Bulk Pack 离心管1.5ml, 有齿, 散装	500pcs/bag, 10000pcs/ctn
TYA10	Centrifuge Tube 1.5ml, Bulk Pack 离心管1.5ml, 散装	500pcs/bag, 10000pcs/ctn




TYA11~TYA13 >>>

Code 编号	Description 描述	Packing 包装
TYA11	Centrifuge Tube 0.5ml, Shield Lid, Bulk Pack 离心管0.5ml, 有齿, 散装	1000pcs/bag, 25000pcs/ctn
TYA12	Centrifuge Tube 0.2ml, Without Graduation, Flat Lid, Bulk Pack 离心管0.2ml, 无刻度, 平盖, 散装	1000pcs/bag, 50000pcs/ctn
TYA13	Centrifuge Tube 0.2ml, Without Graduation, Domed Lid, Bulk Pack 离心管0.2ml, 无刻度, 鼓盖, 散装	1000pcs/bag, 50000pcs/ctn

TYA14~TYA15 >>>

Code 编号	Description 描述	Packing 包装
TYA14	8-pc PCR Tube Strip for 0.2ml, Tube and Lid Separate Packing PCR 8连管, 单孔容量0.2ml, 管盖分开包装	200pcs/bag, 5000pcs/ctn
TYA15	12-pc PCR Tube Strip for 0.2ml, Tube and Lid Separate Packing PCR 12连管, 单孔容量0.2ml, 管盖分开包装	200pcs/bag, 5000pcs/ctn





TYB10~TYB11 >>>

Code 编号	Description 描述	Packing 包装
TYB10	81-well 1.8ml Cryotube Box 81孔1.8ml冷冻管盒	100pcs/ctn
TYB11	100-well 1.8ml Cryotube Box 100孔1.8ml冷冻管盒	50pcs/ctn



TYB12 >>>

Code 编号	Description 描述	Packing 包装
TYB12	100-well 1.8ml Cryotube Box, with Hinged Lid 100孔1.8ml冷冻管盒, 搭扣式	50pcs/ctn



TYB13 >>>

Code 编号	Description 描述	Packing 包装
TYB13	36-well 5ml Cryotube Box 36孔5ml冷冻管盒	50pcs/ctn



TYB14 >>>

Code 编号	Description 描述	Packing 包装
TYB14	96-well Freezing Box 96孔冰盒	50pcs/ctn



TYB15 >>>

Code 编号	Description 描述	Packing 包装
TYB15	81-well 1.8ml Cryotube Paper Box 81孔1.8ml冷冻管纸盒	100pcs/ctn

TYB16 >>>

Code 编号	Description 描述	Packing 包装
TYB16	100-well 1.8ml Cryotube Paper Box 100孔1.8ml冷冻管纸盒	100pcs/ctn



TYC01~TYC05 >>>

Code 编号	Description 描述	Packing 包装
TYC01	Gilson Pipette Tip 10ul 吉尔森吸头10ul	1000pcs/bag, 50000pcs/ctn
TYC02	Gilson Pipette Tip 200ul 吉尔森吸头200ul	1000pcs/bag, 50000pcs/ctn
TYC03	Gilson Graduated Tip 200ul 吉尔森刻度吸头200ul	1000pcs/bag, 50000pcs/ctn
TYC04	Gilson Pipette Tip 1000ul 吉尔森吸头1000ul	500pcs/bag, 15000pcs/ctn
TYC05	Gilson Graduated Tip 1000ul 吉尔森刻度吸头1000ul	500pcs/bag, 15000pcs/ctn



TYC06~TYC08 >>>

Code 编号	Description 描述	Packing 包装
TYC06	Eppendorf Pipette Tip 10ul 爱频道夫吸头10ul	1000pcs/bag, 50000pcs/ctn
TYC07	Eppendorf Pipette Tip 200ul 爱频道夫吸头200ul	1000pcs/bag, 40000pcs/ctn
TYC08	Eppendorf Pipette Tip 1000ul 爱频道夫吸头1000ul	500pcs/bag, 15000pcs/ctn



TYC09 >>>

Code 编号	Description 描述	Packing 包装
TYC09	Gilson Pipette Tip 1250ul 吉尔森吸头1250ul	500pcs/bag, 10000pcs/ctn



TYC10~TYC11 >>>

Code 编号	Description 描述	Packing 包装
TYC10	Fixed Pipette 定量移液器	1pcs/box, 50pcs/ctn
TYC11	Adjustable Pipette 可调节移液器	1pc/box, 50pcs/ctn



TYC12~TYC13 >>>

Code 编号	Description 描述	Packing 包装
TYC12	Filter Tip 10ul & Box 吉尔森滤芯吸头10ul和吸头盒	96/rack, 10racks/unit, 5units/ctn
TYC13	Filter Tip 200ul & Box 吉尔森滤芯吸头200ul和吸头盒	96/rack, 10racks/unit, 5units/ctn

TYC14 >>>

Code 编号	Description 描述	Packing 包装
TYC14	Filter Tip 1000ul & Box 吉尔森滤芯吸头1000ul和吸头盒	96/rack, 10racks/unit, 5units/ctn

TYC15~TYC17 >>>

Code 编号	Description 描述	Packing 包装
TYC15	Gilson Pipette tip 5ml 吉尔森吸头5ml	250pcs/bag, 7500pcs/ctn
TYC16	Dalong Pipette tip 5ml 大龙吸头5ml	250pcs/bag, 4000pcs/ctn
TYC17	Finland Pipette tip 5ml 芬兰吸头5ml	250pcs/bag, 5000pcs/ctn

TYC18 >>>

Code 编号	Description 描述	Packing 包装
TYC18	96-well 10ul Tip Box 96孔10ul吸头盒	150pcs/ctn



TYC19 >>>

Code 编号	Description 描述	Packing 包装
TYC19	96-well 200ul Tip Box 96孔200ul吸头盒	100pcs/ctn



TYC20 >>>

Code 编号	Description 描述	Packing 包装
TYC20	60-well 1000ul Tip Box 60孔1000ul吸头盒	100pcs/ctn




TYC21 >>>

Code 编号	Description 描述	Packing 包装
TYC21	28-well 5ml Tip Box (suit for Gilson or Dalong) 28孔5ml吸头盒(适合吉尔森或大龙款)	50pcs/ctn

TYC22~TYC24-1 >>>

Code 编号	Description 描述	Packing 包装
TYC22	1ml Transfer Pipette 1ml巴式吸管	500pcs/bag, 5000pcs/ctn
TYC22-1	1ml Transfer Pipette, Individual Packing, EO Sterile 1ml巴式吸管, 独立包装, EO消毒	1pc/blister bag, 200pcs/box, 2000pcs/ctn
TYC23	2ml Transfer Pipette 2ml巴式吸管	500pcs/bag, 5000pcs/ctn
TYC23-1	2ml Transfer Pipette, Individual Packing, EO Sterile 2ml巴式吸管, 独立包装, EO消毒	1pc/blister bag, 200pcs/box, 2000pcs/ctn
TYC24	3ml Transfer Pipette 3ml巴式吸管	500pcs/bag, 5000pcs/ctn
TYC24-1	3ml Transfer Pipette, Individual Packing, EO Sterile 3ml巴式吸管, 独立包装, EO消毒	1pc/blister bag, 200pcs/box, 2000pcs/ctn





TYD01~TYD14 >>>

Code 编号	Description 描述	Packing 包装
TYD01	Test Tube 12*75mm, PS 试管12*75mm, PS	500pcs/bag, 5000pcs/ctn
TYD02	Test Tube 12*75mm, PP 试管12*75mm, PP	500pcs/bag, 5000pcs/ctn
TYD03	Test Tube 13*75mm, PS 试管13*75mm, PS	500pcs/bag, 5000pcs/ctn
TYD04	Test Tube 13*75mm, PP 试管13*75mm, PP	500pcs/bag, 5000pcs/ctn
TYD05	Test Tube 12*100mm, PS 试管12*100mm, PS	300pcs/bag, 3000pcs/ctn
TYD06	Test Tube 12*100mm, PP 试管12*100mm, PP	300pcs/bag, 3000pcs/ctn
TYD07	Test Tube 13*100mm, PS 试管13*100mm, PS	400pcs/bag, 3000pcs/ctn
TYD08	Test Tube 13*100mm, PP 试管13*100mm, PP	400pcs/bag, 3000pcs/ctn
TYD09	Test Tube 15*100mm, PS 试管15*100mm, PS	250pcs/bag, 2500pcs/ctn
TYD10	Test Tube 15*100mm, PP 试管15*100mm, PP	250pcs/bag, 2500pcs/ctn
TYD11	Test Tube 16*100mm, PS 试管16*100mm, PS	150pcs/bag, 2400pcs/ctn
TYD12	Test Tube 16*100mm, PP 试管16*100mm, PP	150pcs/bag, 2400pcs/ctn
TYD13	Test Tube 16*102mm, PS 试管16*102mm, PS	250pcs/bag, 2500pcs/ctn
TYD14	Test Tube 16*102mm, PP 试管16*102mm, PP	250pcs/bag, 2500pcs/ctn




TYD15 >>>

Code 编号	Description 描述	Packing 包装
TYD15	10ml Test Tube with Screw Cap, PS 10ml试管配螺旋盖, PS	100pcs/bag, 2500pcs/ctn



TYD16 >>>

Code 编号	Description 描述	Packing 包装
TYD16	10ml Test Tube with Screw Cap, PS 10ml试管配螺旋盖, PS	100pcs/bag, 2000pcs/ctn



TYD17~TYD19 >>>

Code 编号	Description 描述	Packing 包装
TYD17	Tube Stopper 12mm, PE 试管塞12mm, PE	1000pcs/bag, 20000pcs/ctn
TYD18	Tube Stopper 13mm, PE 试管塞13mm, PE	1000pcs/bag, 25000pcs/ctn
TYD19	Tube Stopper 16mm, PE 试管塞16mm, PE	1000pcs/bag, 10000pcs/ctn



TYD20~TYD22 >>>

Code 编号	Description 描述	Packing 包装
TYD20	Tube Rack 90 Well, 12-13mm 试管架90孔, 12-13mm	1pcs/bag, 50pcs/ctn
TYD21	Tube Rack 60 Well, 15-17mm 试管架60孔, 15-17mm	1pcs/bag, 50pcs/ctn
TYD22	Tube Rack 40 Well, 20-21mm 试管架40孔, 20-21mm	1pcs/bag, 50pcs/ctn

TYG01~TYG10 >>>



Code 编号	Description 描述	Packing 包装
TYG01	Glass Test Tube 12*75mm 玻璃试管12x75mm	250pcs/box, 2000pcs/ctn
TYG02	Glass Test Tube 12*100mm 玻璃试管12x100mm	250pcs/box, 2000pcs/ctn
TYG03	Glass Test Tube 13*100mm 玻璃试管13x100mm	250pcs/box, 2000pcs/ctn
TYG04	Glass Test Tube 15*100mm 玻璃试管15x100mm	250pcs/box, 1000pcs/ctn
TYG05	Glass Test Tube 16*100mm 玻璃试管16x100mm	250pcs/box, 1000pcs/ctn
TYG06	Glass Test Tube 16*150mm 玻璃试管16x150mm	250pcs/box, 1000pcs/ctn
TYG07	Glass Test Tube 18*180mm 玻璃试管18x180mm	100pcs/box, 600pcs/ctn
TYG08	Glass Test Tube 20*150mm 玻璃试管20x150mm	100pcs/box, 400pcs/ctn
TYG09	Glass Test Tube 20*200mm 玻璃试管20x200mm	100pcs/box, 400pcs/ctn
TYG10	Glass Test Tube 25*200mm 玻璃试管25x200mm	100pcs/box, 400pcs/ctn

TYG11~TYG18 >>>



Code 编号	Description 描述	Packing 包装
TYG11	Glass Test Tube 13*100mm with Screw Cap 罗口玻璃管13x100mm	100pcs/box, 800pcs/ctn
TYG12	Glass Test Tube 15*100mm with Screw Cap 罗口玻璃管15x100mm	100pcs/box, 1200pcs/ctn
TYG13	Glass Test Tube 15*150mm with Screw Cap 罗口玻璃管15x150mm	100pcs/box, 800pcs/ctn
TYG14	Glass Test Tube 16*150mm with Screw Cap 罗口玻璃管16x150mm	100pcs/box, 1000pcs/ctn
TYG15	Glass Test Tube 18*180mm with Screw Cap 罗口玻璃管18x180mm	100pcs/box, 338pcs/ctn
TYG16	Glass Test Tube 20*150mm with Screw Cap 罗口玻璃管20x150mm	100pcs/box, 800pcs/ctn
TYG17	Glass Test Tube 20*200mm with Screw Cap 罗口玻璃管20x200mm	100pcs/box, 400pcs/ctn
TYG18	Glass Test Tube 25*200mm with Screw Cap 罗口玻璃管25x200mm	100pcs/box, 400pcs/ctn



TYG19~TYG19-1 >>>

Code 编号	Description 描述	Packing 包装
TYG19	Blue Micro Hematocrit Capillary Tube, Glass 蓝标玻璃毛细管, 100支/桶	100pcs/barrel, 50000pcs/ctn
TYG19-1	Blue Micro Hematocrit Capillary Tube, Glass 蓝标玻璃毛细管, 500支/桶	500pcs/barrel, 50000pcs/ctn



TYG20~TYG20-1 >>>

Code 编号	Description 描述	Packing 包装
TYG20	Red Micro Hematocrit Capillary Tube, Glass 红标玻璃毛细管, 100支/桶	100pcs/barrel, 50000pcs/ctn
TYG20-1	Red Micro Hematocrit Capillary Tube, Glass 红标玻璃毛细管, 500支/桶	500pcs/barrel, 50000pcs/ctn




TYH01 >>>

Code 编号	Description 描述	Packing 包装
TYH01	Embedding Cassette Stripe Holes, POM 包埋盒, 长条型, POM材质	500pcs/bag, 5000pcs/ctn



TYH02 >>>

Code 编号	Description 描述	Packing 包装
TYH02	Embedding Cassette Square Holes, POM 包埋盒, 方格型, POM材质	500pcs/bag, 5000pcs/ctn



TYH03 >>>

Code 编号	Description 描述	Packing 包装
TYH03	Embedding Cassette Round Holes without Cover, POM 包埋盒, 圆孔无盖, POM材质	500pcs/bag, 5000pcs/ctn

TYH04~TYH09 >>>


Code 编号	Description 描述	Packing 包装
TYH04	Histology Tissue Molds 7*7*6, Stainless Steel 包埋盒底模, 7x7x6, 不锈钢材质	50pcs/bag
TYH05	Histology Tissue Molds 15*15*6, Stainless Steel 包埋盒底模, 15x15x6, 不锈钢材质	50pcs/bag
TYH06	Histology Tissue Molds 24*24*6, Stainless Steel 包埋盒底模, 24x24x6, 不锈钢材质	50pcs/bag
TYH07	Histology Tissue Molds 30*24*6, Stainless Steel 包埋盒底模, 30x24x6, 不锈钢材质	50pcs/bag
TYH08	Histology Tissue Molds 37*24*6, Stainless Steel 包埋盒底模, 37x24x6, 不锈钢材质	50pcs/bag
TYH09	Cover for Embedding Cassette, Stainless Steel 包埋盒盖子, 不锈钢材质	100pcs/bag



TYI-7101~TYI-7105 >>>

Code 编号	Description 描述	Packing 包装
TYI-7101	Polished Edges, Plain 抛光边	50pcs/box, 50box/ctn
TYI-7102	Cut Edges, Plain 未磨边	50pcs/ box, 50box/ctn
TYI-7105	Polished Edges, Frosted One End On One Side 抛光边, 单头单面磨砂	50pcs/ box, 50box/ctn

Note: 72pcs/box for slide is available base on client' s request



TYI01 >>>

Code 编号	Description 描述	Packing 包装
TYI01	Adhesion Microscope Slide 粘附载玻片 (液基细胞子, 冷冻病理切片, 免疫组化专用)	50pcs/box, 50box/ctn



TYI-1818C~TYI-2460C >>>

Code 编号	Description 描述	Packing 包装
TYI-1818C	18*18mm	100pcs/box, 500box/ctn
TYI-2020C	20*20mm	100pcs/box, 500box/ctn
TYI-2222C	22*22mm	100pcs/box, 500box/ctn
TYI-2424C	24*24mm	100pcs/box, 500box/ctn
TYI-2432C	24*32mm	100pcs/box, 250box/ctn
TYI-2450C	24*50mm	100pcs/box, 250box/ctn
TYI-2460C	24*60mm	100pcs/box, 250box/ctn




TYI02~TYI07 >>>

Code 编号	Description 描述	Packing 包装
TYI02	1 Place, Paper 1片装, 纸质	1000pcs/ctn
TYI03	2 Place, Paper 2片装, 纸质	500pcs/ctn
TYI04	3 Place, Paper 3片装, 纸质	500pcs/ctn
TYI05	4 Place, Paper 4片装, 纸质	200pcs/ctn
TYI06	10 Place, Paper 10片装, 纸质	100pcs/ctn
TYI07	20 Place, Paper 20片装, 纸质	50pcs/ctn




TYI08 >>>

Code 编号	Description 描述	Packing 包装
TYI08	Slide Tray for Microscope Slide, without Cover 晾片板, 20片装, 无盖	100pcs/ctn



TYI09 >>>

Code 编号	Description 描述	Packing 包装
TYI09	5 Place Mailer, PP 5片装邮寄夹, PP	100pcs/bag, 1200pcs/ctn



TYI10 >>>

Code 编号	Description 描述	Packing 包装
TYI10	25 Place Slide Box 25片装载玻片盒	1pc/bag, 200pcs/ctn



TYJ06 Blood Collection Pen
采血笔

Length: 12.8cm
长度: 12.8厘米
Diameter: 1.6cm
宽度: 1.6厘米


Packing: 1pc/box, 100pcs/ctn



TYJ07 Safety Pressure Activated Blood Lancet
安全式方针

Material: stainless steel
材料: 不锈钢
Sterilization: gamma radiation
消毒: 伽玛消毒
Gauge: 21g, 23g, 26g, 28g, 30g
尺寸: 21g, 23g, 26g, 28g, 30g

Packing: 100pcs/box, 20box/ctn



TYJ08 Safety Bean Sprout Blood Lancet
安全豆芽式采血针

Material: stainless steel
材料: 不锈钢
Sterilization: gamma radiation
消毒: 伽玛消毒
Gauge: 21g, 23g, 26g, 28g, 30g
尺寸: 21g, 23g, 26g, 28g, 30g


Packing: 100pcs/box, 20box/ctn



TYJ09 Blood Collection Needle
真空采血针

Sterilization: gamma radiation
消毒: 伽玛消毒
Gauge: 18g-24g
尺寸: 18g-24g
Color: green, black, yellow, blue, pink
颜色: 绿色, 黑色, 黄色, 蓝色, 粉色

Packing: 100pcs/box, 5000pcs/ctn



TYJ10 Cotton Swab
棉签


Sterilization: EO or gamma radiation
消毒: EO或者伽玛消毒
Stick Material: wooden or plastic
材质: 木棒或者塑料棒
Size: 2.2*15mm
尺寸: 2.2x15mm

Packing: 100pcs/box, 10000pcs/ctn

TYJ11 Female Swab
女性拭子

Sterilization: EO or gamma radiation
消毒: EO或者伽玛消毒
Stick Material: wooden or plastic
材质: 木棒或者塑料棒

Packing: 100pcs/box, 2000pcs/ctn



TYJ12 Transport Swab with Gel
凝胶拭子

Medium: amies (cir), stuart (cir), cary blair
Sterilization: EO or gamma radiation
消毒: EO或者伽玛消毒
Stick Material: wooden or plastic
材质: 木棒或者塑料棒

Packing: 100pcs/box, 1000pcs/ctn





Plain Tube
普通管

Usage: Applying to blood specimen collection and storage for biochemical, immunological and serum clinical examination. Whether it is glass or plastic tube, the wall is extremely smooth, and can effectively prevent blood cells attaching to the wall, which can minimize the specimens in vitro hemolysis.

用途：用于生化、免疫、血清等临床医学检验中血液样本的采集和贮藏。无论是玻璃管还是塑料管，其内壁都极其光滑，能有效避免血细胞附壁，最大限度减少标本体外溶血。

Code 编号	Material 材质	Size(mm) 规格	Vol 采血量	Packing 包装	Ctn size(cm) 箱子尺寸	Gw (kg) 毛重
TYK01	PET	13*75	3-5ml	100pcs/rack 1800pcs/ctn	52*38*31	13.5
TYK01-1	Glass					18.2
TYK02	PET	13*100	5-7ml	100pcs/rack 1800pcs/ctn	52*38*38	15.8
TYK02-1	Glass					22
TYK03	PET	16*100	8-10ml	100pcs/rack 1200pcs/ctn	51*41*25	14
TYK03-1	Glass					17.5



Gel and Clot Activator Tube
分离胶管

Usage: Applying to the acquisition of high-quality serum, specimen for biochemical, immunological and serum clinical examination. Serum separation gel is a kind of semi-solid inert gel, composed of several compounds and synthesized by medical polymeric material. After blood centrifugal, the specimen can be separated into serum top, separation gel middle and blood corpuscle lower. In the middle of the tube, a thick layer of separation gel can completely isolate serum from blood corpuscle, so as to ensure the purity of serum.

用途：用于生化、免疫、血清等临床医学检验中获得高品质血清样本。血清分离胶是一种由数种化合物所组成之半固体惰性胶，该分离胶全部采用医用高分子材料合成，采用离心后，将分离成血清（上层）、血清分离胶（中层）及血球（下层）。试管中间厚厚一层血清分离胶能完全隔离血清及血球，让血球无法接触到血清，从而保证了血清的纯净度。

Code 编号	Material 材质	Additive 添加剂	Size(mm) 规格	Vol 采血量	Packing 包装	Ctn size(cm) 箱子尺寸	Gw (kg) 毛重
TYK04	PET	Gel & Clot Activator	13*75	3-5ml	100pcs/rack 1800pcs/ctn	52*38*31	16
TYK04-1	Glass						24.2
TYK05	PET	Gel & Clot Activator	13*100	5-7ml	100pcs/rack 1800pcs/ctn	52*38*38	17.2
TYK05-1	Glass						26
TYK06	PET	Gel & Clot Activator	16*100	8-10ml	100pcs/rack 1200pcs/ctn	51*41*25	18
TYK06-1	Glass						22.5



Pro-Coagulation Tube
促凝管

Usage: Applying to fast collection of blood specimen for biochemical and immunological clinical examination. A kind of coagulant, which is different from other chemical coagulants having effects on the test results, was sprayed evenly on the inner wall of the tube to speed the blood solidification without any heating devices (such as water container). After blood collection, let the tube stand for 15 minutes (make sure the ambient temperature over 16°C), and then you can get the ideal and high purity serum.

用途：用于生化、免疫等临床医学检验中要求快速血液样本的采集。该储血管内添加的促凝剂，有别于其他化学促凝法对检测结果的影响。促凝剂均匀喷涂在储血管内壁上，无需借助于任何加热设施（如水温箱等），完全采用纯物理化，致使血液快速凝固。采血完毕后，静置立放15分钟（离心环境温度保证在16°C），即可得到理想、高纯净度的血清。

Code 编号	Material 材质	Additive 添加剂	Size(mm) 规格	Vol 采血量	Packing 包装	Ctn size(cm) 箱子尺寸	Gw (kg) 毛重
TYK07	PET	pro-coagulation	13*75	3-5ml	100pcs/rack 1800pcs/ctn	52*38*31	13.5
TYK07-1	Glass						18.2
TYK08	PET	pro-coagulation	13*100	5-7ml	100pcs/rack 1800pcs/ctn	52*38*38	15.8
TYK08-1	Glass						22
TYK09	PET	pro-coagulation	16*100	8-10ml	100pcs/rack 1200pcs/ctn	51*41*25	14
TYK09-1	Glass						17.5



PT Tube
血凝管

Usage: Applying to blood clotting (PT, APTT, coagulation factor) test for clinical examination. The concentration of sodium citrate is 109mol/L and the accurate ratio of blood and anticoagulant is 9:1.

用途：用于临床医学检验中凝血机制（PT、APTT、血凝因子）检验。柠檬酸钠浓度109mol/L，采集血液与抗凝剂的精确比为9：1。

Code 编号	Material 材质	Additive 添加剂	Size(mm) 规格	Vol 采血量	Packing 包装	Ctn size(cm) 箱子尺寸	Gw (kg) 毛重
TYK10	PET	Sodium Citrate 1:9	13*75	1.8-4.5ml	100pcs/rack 1800pcs/ctn	52*38*38	14.2
TYK10-1	Glass						18.2



Edta Tube
血常规管

Usage: applying to clinical haematology test, such as blood routine examination, all kinds of blood cells analysis (such as blood cell morphology examination, blood cell count, hematocrit, hemoglobin, sicking test, etc.)

用途：用于临床血液学检验，如血常规等。适用各种血细胞分析（如血细胞形态镜检、血细胞计数、红细胞比积、血红蛋白、镰变实验等）。

Code 编号	Material 材质	Additive 添加剂	Size(mm) 规格	Vol 采血量	Packing 包装	Ctn size(cm) 箱子尺寸	Gw (kg) 毛重
TYK11	PET	EDTA K2/K3	13*75	1-5ml	100pcs/rack 1800pcs/ctn	52*38*31	13.5
TYK11-1	Glass						18.2
TYK12	PET	EDTA K2/K3	13*100	5-7ml	100pcs/rack 1800pcs/ctn	52*38*38	15.8
TYK12-1	Glass						22
TYK13	PET	EDTA K2/K3	16*100	8-10ml	100pcs/rack 1200pcs/ctn	51*41*25	14
TYK13-1	Glass						17.5



Heparin Tube
肝素管

Usage: applying to clinical biochemistry, blood specimen collection for emergency biochemical tests and anticoagulation(blood rheology measurement, biochemical test, genetic diagnosis, enzymology, enzyme immunoassay, etc.)

用途：用于临床生化、急诊生化检验中血液样本的采集和抗凝（血变学测定、生化、基因诊断、酶学、酶免法等检验）。

Code 编号	Material 材质	Additive 添加剂	Size(mm) 规格	Vol 采血量	Packing 包装	Ctn size(cm) 箱子尺寸	Gw (kg) 毛重
TYK14	PET	Heparin lithium	13*75	3-5ml	100pcs/rack 1800pcs/ctn	52*38*31	13.5
TYK14-1	Glass						18.2
TYK15	PET	Heparin lithium	13*100	5-7ml	100pcs/rack 1800pcs/ctn	52*38*38	15.8
TYK15-1	Glass						22
TYK16	PET	Heparin lithium	16*100	8-10ml	100pcs/rack 1200pcs/ctn	51*41*25	14
TYK16-1	Glass						17.5




Glucose Tube
血糖管

Usage: Applying to clinical determination of blood glucose. The original characters of blood specimen can be maintained for a long time.

用途：用于临床血糖值的测定，能较长时间保持血液样本的原始性状。

Code 编号	Material 材质	Additive 添加剂	Size(mm) 规格	Vol 采血量	Packing 包装	Ctn size(cm) 箱子尺寸	Gw (kg) 毛重
TYK17	PET	Sodium fluoride /potassium oxalate	13*75	2-5ml	100pcs/rack 1800pcs/ctn	52*38*31	13.5
TYK17-1	Glass						18.2
TYK18	PET	Sodium fluoride /potassium oxalate	13*100	5-7ml	100pcs/rack 1800pcs/ctn	52*38*38	15.8
TYK18-1	Glass						22
TYK19	PET	Sodium fluoride /potassium oxalate	16*100	8-10ml	100pcs/rack 1200pcs/ctn	51*41*25	14
TYK19-1	Glass						17.5



ESR Vacuum Tube
血沉管

Usage: 1.28ml tube with a total capacity of 1.6ml, is equipped with 0.129mol/L sodium citrate 0.32ml. Collect blood 1.28ml to the central scale and after blood collection, we immediately invert the tube gently for 4-5 times so as to mix blood completely with anticoagulant, then put the tube into a D2010 automatic ESR analyzer and after 30 minutes read the data. This product can save time and effort and at the same time prevent cross infection.

用途：1.28ml管装有0.129mol/L柠檬酸钠0.32ml，采血1.28ml，至中央刻度处，总容量为1.6ml。采血后，立即轻轻反复颠倒4-5次，使抗凝剂充分与血液混匀。然后放入D2010型自动血沉分析仪内，30分钟后即可读数。该产品省时、省力，从源头杜绝了交叉感染。

Code 编号	Material 材质	Additive 添加剂	Size(mm) 规格	Vol 采血量	Packing 包装	Ctn size(cm) 箱子尺寸	Gw (kg) 毛重
TYK20	PET	Sodium citrate1:4	13*75	1.6-2.4ml	100pcs/rack 1800pcs/ctn	52*38*31	13.5
TYK20-1	Glass						18.2

Code 编号	Material 材质	Additive 添加剂	Size(mm) 规格	Vol 采血量	Packing 包装	Ctn size(cm) 箱子尺寸	Gw (kg) 毛重
TYK21	Glass	Sodium citrate1:4	9*120	1.28ml	100pcs/rack 1800pcs/ctn	48*37*45	18.2



Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109

7109 Polished Edges Single Frosted End Glass Microscope Slide For Laboratory

Microscope slide is a glass or quartz sheet that be used to hold specimen when observing things with a microscope. When making a sample, cell or tissue section is placed on the slide and a cover glass is placed on it for observation. In optics, which is a thin sheet of a glass material used to produce a phase difference.

Product Description:

Microscope slide is a glass or quartz sheet that be used to hold specimen when observing things with a microscope. When making a sample, cell or tissue section is placed on the slide and a cover glass is placed on it for observation. In optics, which is a thin sheet of a glass material used to produce a phase difference.

Product Features:

【Material】 : High quality clear optical glass with good light transmission

【Ready to use】 : Pre-cleaned and ready to use. A thorough automatic cleaning process remove all particles, debris, surface; contaminants and ensures that slides are ready to use

【Safety Corners】 : Polished ground edges, Single side color frosted end

【Universal Use】 : Suitable for various Microscope

【Cleaning method】 : Wash with clean water or alcohol, and then wipe with cotton or gauze (it is better to use mirror wiping paper, fingers to avoid touching the slide, so as not to leave fingerprints on it, affecting the use of next observation)

Product Specification:

Model	Size	Thickness	Description	Packing
TYI-7101-A	25.4x76.2 (1"x3")	1.0-1.2mm	Polished Edges,Plain	50pcs/box, 50box/ctn
TYI-7101-B			Polished Edges,Plain	
TYI-7102			Cut Edges, Plain	
TYI-7105			Polished Edges,Single Frosted End	
TYI-7107			Polished Edges,Double Frosted End	
TYI-7109			Clear glass, Polished Edges, Single Frosted End	
TYI01			Adhension Slide, Positive Charged	



Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109

Product Details:



Packaging Details

50pcs in one box,50box in one carton

Carton size: 37x 18x 16cm

Gross Weight: 13KG

Net Weight: 12KG

Also can customize package as per your request.



Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109



Autoclavable 8 Channel 12 Channel Adjustable Mechanical Pipette

This pipette is a general purpose pipette for sampling and dispensing accurate liquid volumes. The pipettes operate on the air displacement principle and use disposable tips. The black tip cone can be sterilized 121°C.

This pipette is a general purpose pipette for sampling and dispensing accurate liquid volumes. The pipettes operate on the air displacement principle and use disposable tips. The black tip cone can be sterilized 121°C.

The pipettes cover a volume range of 0.5ul to 5ml.

All pipettes have been quality tested according to our company's demands. The quality control according to our company's involves gravimetric testing of each pipette with distilled water at 22°C using the manufacture's original tips.

Method of Application

1. Choose the right pipettor
2. Set the fluid volume
3. Install the suction head
4. Moving fluid
5. Placement of pipette

Method of Getting Liquor

1. Forward pipetting

Press the pipetting lever to the first stop position, then slowly release the button to return to the original point; then press the pipetting lever to the first stop position to discharge the liquid, stop for a while and continue to press the pipetting lever to the second stop Drain the remaining liquid at the point, and finally slowly loosen the pipetting lever.



Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109

2. Reverse pipetting

First press the button to the second stop position, slowly release the pipetting lever to return to the original point, press the pipetting lever to the first stop position when discharging liquid to discharge the set volume of liquid, continue to hold the pipetting lever. The operating lever is located at the first stop position. Remove the tip with residual liquid and discard it.

Product Specifications

Pipette 8 channel			
Item	Description	Increment	Packing Unit
NW3033	0.5-10ul	0.1ul	Independent packing
	5-50ul	0.5ul	Independent packing
	50-300ul	5ul	Independent packing

Pipette 12 channel			
Item	Description	Increment	Packing Unit
NW3034	0.5-10ul	0.1ul	Independent packing
NW3043	5-50ul	0.5ul	Independent packing
NW3044	50-300ul	5ul	Independent packing

Product Details:





Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109

Autoclavable Single Channel High Temperature-Resistant Accurate Fixed Adjustable Pipette

Single channel Pipette, also called pipette gun, is a measuring tool that can transfer liquid from the original container to another container within a certain range.It is widely used in biology, chemistry, Medical diagnosis and Laboratory testing etc.

Product Description:

Single channel Pipette, also called pipette gun, is a measuring tool that can transfer liquid from the original container to another container within a certain range.It is widely used in biology, chemistry, Medical diagnosis and Laboratory testing etc.

Product Features:

【Volume】 : This Pipette has two types, Fixed or Adjustable volume, a variety of ranges can be selected, you can choose according to your application.

Fixed volume: 1ul,2ul,5ul,10ul,20ul,30ul,50ul,100ul,200ul,500ul,1000ul,5000ul,10000ul etc

Adjustable volume: 0.1-2.5ul,0.5-10ul,5-50ul,50-200ul,200-1000ul, 1000-5000ul, 2-10ml etc

【Tolerable temperature】 : Resistant to high temperature, the whole product can be autoclaved and high temperature disinfection.

【Lightweight and Ergonomic】 : Ergonomic design provides excellent operating experience

Lightweight handle design, Large display window allows for easy volume identification

【Highly accuracy】 : Every pipette is tested to ensure precision and accuracy. If the volume during use is not accurate, you can re-calibrate.

【Simple maintenance】 : Easy calibration and maintenance

Product Specification:

Item	Volume	Type	Packing unit
6015001	1 ul	Fixed	Independent packing
6015002	2 ul		
6015003	2.5 ul		
6015004	5 ul		
6015005	7 ul		
6015006	10 ul		
6015007	15 ul		
6015008	20 ul		
6015009	25 ul		
6015010	30ul		
6015011	40ul		
6015012	50ul		
6015013	75 ul		
6015014	80 ul		
6015015	100 ul		
6015016	200 ul		
6015017	250 ul		
6015018	500 ul		
6015019	1000 ul		



Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109

6015020	2000 ul		
6015021	2500 ul		
6015022	5000 ul		
6015023	10000 ul		

Item	Volume range	Increment	Type	Packing Unit
6015024	0.1-2.5 ul	0.05 μ l	Adjustable	Independent packing
6015025	0.5-10 ul	0.1 μ l		
6015026	2-10 ul	0.5 μ l		
6015027	5-50 ul	0.5 μ l		
6015028	10-100 ul	1 μ l		
6015029	20-200ul	1 μ l		
6015030	50-200 ul	1 μ l		
6015031	100-1000 ul	5 μ l		
6015032	200-1000 ul	5 μ l		
6015033	1000-5000 ul	50 μ l		
6015034	2-10ml	0.1ml		

Product Details:



Packaging Details

One piece in independent packing, then suitable quantity in cartons.

Single package size: 33X13X5 cm

Single gross weight: 0.300 kg

Also can customize package as per your request.



Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109

Universal Filter Pipette Tips With Box

Introduction :

Filter Pipette Tips are racked, sterile and universal-fit. They are the recommended filtered pipette tips for different brand pipettors. They are also compatible with Thermo, Brand, Eppendorf, Gilson and Rainin micropipettes. Our filter tips offer extra protection from contamination by use of a barrier between the tipcone of the pipette and the sample.

Specifications :

Made of high quality polypropylene material


Unique long tube mouth design makes fluid transfer more accurate

With super hydrophobicity, the hydrophobic filter element forms a barrier to aerosols, eliminating the risk of cross-contamination;

Suitable for commonly used pipettes on the market, with optimized aperture to ensure smooth sample extraction; (at 121°C for 20 min);

Code	Volume	Type	Material	Color	Packing
TYM-10	10ul	DNASE RANSE FREE	PP	Transparent	96pcs/box, 10boxes/middle box, 50boxes/ctn
TYM-50	50ul				
TYM-100	100ul				
TYM-200	200ul				
TYM-1000	1000ul				
TYM-1250	1250ul				

Details





Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109

Feature:

1. Manufactured from high-density, inert Ultra-high-molecular-weight polyethylene
2. Provides a strong barrier against contamination while not impeding airflow
3. Hydrophobic properties prevent absorption of samples
4. Certified free of RNase, DNase, DNA, PCR inhibitors and pyrogen-free



Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109

Sterile Plastic Urine Container Urine Cup 30ml 40ml 60ml 100ml 120m

Product Description:

Material: made of polypropylene (PP) and PS, With screw cap

Color: transparent

Cap color: Red, Blue, Yellow, etc.

Usage: The specimen containers are mainly used for urine collection

Sterile: E.O. Sterile or Non-Sterile

Product Features:

Capacity: 30ml, 40ml, 60ml, 100ml, 120ml

Urine Container material: PP or PS

Caps are available in different colors

Availabel with or without label, individually packed or bulk packed

Product Specifications:

Code	Description	Diameter (cm)	Height (cm)	Package	Carton size (cm)	Weight (kg)
TYF01	30ml specimen cup Ps, no graduation, bulk pack, Non-sterile Leak proof	2.3	9.4	1000pcs/ctn	52x45x34	12
TYF01-S	30ml specimen cup Ps, no graduation, single bag, Leak proof Sterile			1pc/bag, 800pcs/ctn	56x36x42	8
TYF02	60ml specimen cup PP, molded graduation, bulk pack, Non-sterile Leak proof	3.7	6.5	500pcs/carton	60x29x38	6
TYF02-S	60ml specimen cup PP, molded graduation,			1pc/bag 500pcs/carton	56x36x42	5.8



Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China

Tel: 0086-023-67728212

Fax:0086-023-67728109

	individual pack, sterile, Leakproof					
TYF03	120ml specimen cup PP, molded graduation, bulk pack, Non-sterile Leak proof	4.7	7, 2	500pcs/carton	60x30x35	8
TYF03-S	120ml specimen cup PP, molded- graduation, individual pack, sterile, Leakproof			1pc/polybag 500pcs/carton	60x50x55	8.5

Product Advantages:

Water-proof design can prevent the specimen from leaking out.

Product Details:





Chongqing New World Trading Co.,Ltd
Address: Suite 2-4, Block 30, No. 11 Tianlan Road,
Jiangbei District, Chongqing, China
Tel: 0086-023-67728212 Fax:0086-023-67728109

Lab Disposable 1ml 2ml 5ml 10ml Plastic Sterile Serological Pipette

The serological pipette is a measuring device used to accurately pipette a certain volume of solution. The pipette is a measuring instrument, which is only used to measure the volume of the solution it releases. It is a slender glass tube with an enlarged part in the middle. The lower end is sharp-mouthed, and a marking line is engraved on the upper end of the neck, which is a mark of the exact volume to be removed.

Send Inquiry

Product Description:

The serological pipette is a measuring device used to accurately pipette a certain volume of solution. The pipette is a measuring instrument, which is only used to measure the volume of the solution it releases. It is a slender glass tube with an enlarged part in the middle. The lower end is sharp-mouthed, and a marking line is engraved on the upper end of the neck, which is a mark of the exact volume to be removed.

Product Features:

1. Single Wrapped in Paper&Plastic Bag
2. Graduations Inaccuracy < 2%
3. E-Beam Sterile
4. Certified with DNase/RNase Free, Non-Pyrogenic

Product Specifications:

Code	Specification	Volume	Packing	Carton Size (cm)	Weight (kg)
TYM01	Serological Pipette PS Material Clear Graduated	1ml	50pcs/bag, 10bags/box, 6box/ctn	67.5*36*47	15
TYM02		2ml	50pcs/bag, 10bags/box, 6box/ctn	67.5*36*47	22.5
TYM03		5ml	50pcs/bag, 10bags/box, 3box/ctn	67*43*34	18
TYM04- 1		10ml	50pcs/bag, 8bags/box, 3box/ctn	67*43*34	16.5

Product Advantages:

1. Available in six capacity: 1ml, 2ml, 5ml, 10ml, 25ml and 50ml.
2. Smooth and clean without drawing lines.
3. Excellent double-side graduation design for liquid adding or reducing.
4. All serological pipettes are fixed with high quality filters

Product Details:



Smooth suction mouth
Not easy to adhered



The surface is neat and smooth, which greatly reduces the adhesion to the pipe wall so the sampling accuracy is higher

Different Color Signs





Chongqing New World Trading Co.,Ltd

Suite 2-4, Block 30, No. 11 Tianlan Road, Jiangbei District, Chongqing, China
Tel: 0086-023 67728212 Fax: 0086-023 67728109 whatsapp:008613677626206 Email:coco@cnwtc.com

Specification Sheet

Product Name	Glass Test Tube
Model No	TYG04-1
Type	15X120mm
Size	Diameter: 15mm +/- 1mm Inner diameter:12.8+/-0.1mm Length: 120mm +/- 1mm
Color :	transparent
Material	Glass
Packing	Packing method: 250pcs/box ,1000pcs/ctn Carton size:41.5x29x29.5cm Gross weight: 13kgs

Chongqing New World Trading Co.,Ltd

Coco Jiang

